

# Spirit Select<sup>®</sup> Bed

## control panel guide

Caregivers' guide to quickly reference control panel notifications on Spirit Select beds.

**Scale**  
Press and hold. Follow instructions on screen to zero the scale.

**Watchdog<sup>®</sup> Bed Monitoring System**  
Press and hold to reset Watchdog. Press once to view configuration changes on the screen.

**Lights menu**  
Press and hold to access Lights Menu.

**Nurse reminder**  
Press and select desired time or Nurse Call option to set Nurse Reminder.

**View weight**  
Press to wake up the screen and view patient weight

**CPR release**  
Press and hold to engage CPR release. Manual release is located at the head of bed.

**Bed exit alarm**  
When patient is in bed, Press to set Bed Exit Alarm. To turn off, press again. When patient sets off the alarm, press once to pause. Press twice to turn off.

**Sound**  
Press until desired alarm volume is lit. To change alarm tone, press and hold and select desired tone.

The control panel features a central LCD display showing '0.0 ZERO' and a 'SCALE' icon. To the right of the display are 'MENU' and 'ON ENTER' buttons. Below the display is a grid of function buttons: 'WATCHDOG', 'LIGHT', 'REMEMBER', 'HEAD', 'BED', 'KNEE', 'TREND', '100% TREND', 'CHAIR', 'LEVEL ALL', 'CPR', 'AC DC ATTN BATTERY', 'RESET SET', 'BED EXIT ALARM', and 'CHIME'. A 'VOLUME' control is located at the bottom right of the panel.

# Spirit Select<sup>®</sup> Bed

## Lights guide

stryker

Caregivers' guide to quickly reference Watchdog<sup>®</sup> bed and status indicator notification lights on Spirit Select beds.


**Watchdog**  
light is  
turned off



**Watchdog**  
is set and has  
not changed  
**Brake** is set



**Parameters**  
have changed  
(bed height, siderails,  
optional parameters:  
head, knee, bed exit  
alarm, lock-outs)



**Brake**  
is not set



**Nurse  
reminder**  
time is up



**Bumper  
nightlight**  
color options

