Performance-PRO™ XT

REF Model 6086

# stryker®

## Operations/Maintenance Manual

Manuel d'utilisation et d'entretien  
Bedienungs- und wartungshandbuch  
Gebruiks-/onderhoudshandleiding  
Manuale d'uso e manutenzione  
Manual de uso y mantenimiento  
Manual de funcionamento/manutenção  
Användar-/underhållshandbok  
Betjenings-/vedligeholdelsesmanual  
Käyttö- ja huolto-ohjekirja  
Instrukcja obsługi/konserwacji  
Návod k obsluze a údržbě

CE

# International Addresses

## EUROPE HEADQUARTERS

**Stryker SA**  
Cite-Centre, Grand-Rue 90  
P.O. Box 1568  
1820 Montreux, Switzerland  
Phone: 011-41-21-966-12-01  
Fax: 011-41-21-966-12-00

## EASTERN EUROPE

**Stryker SA - Export Business**  
Grand-Rue 90  
P.O. Box 1567  
1820 Montreux, Switzerland  
Phone: 011-41-21-966-14-00  
Fax: 011-41-21-966-14-01

## AUSTRALIA

**Stryker Australia**  
8 Herbert Street  
St Leonards NSW 2065  
Australia  
Phone: 61-2-9467-1000  
Fax: 61-2-9467-1010

## AUSTRIA

**Stryker-Howmedica Osterreich GmbH**  
Millenium Tower  
Handelskai 94-96  
120 Wien  
Austria  
Phone: 011-43-1-240-27-6400  
Fax: 011-43-1-240-27-6410

## BELGIUM

**NV Stryker SA**  
Ikaros Business Park- Fase III  
Ikaroslaan 12  
1930 Zaventem, Brussels  
Belgium  
Phone: 011-32-2-717-9210  
Fax: 011-32-2-717-9249

## DENMARK

**Stryker Denmark**  
Filial of Stryker AB  
Postbox 772  
1532 Copenhagen  
Denmark  
Phone: 011-45-33-93-6099  
Fax: 011-45-33-93-2069

## UK/IRELAND

**Stryker UK Limited**  
Stryker House  
Hambridge Road  
Newbury, Berkshire  
RG14 5EG, England  
Phone: 011-44-1635-556-500  
Phone: 011-44-1635-262-400  
Fax: 011-44-1635-580-300

## FINLAND

**Stryker AB, Finland**  
PO 80 Makelankatu 2  
00501 Helsinki  
Finland  
Phone: 011-35-89-774-4680  
Fax: 011-35-89-774-46820

## FRANCE

**Stryker France S.A.S.**  
ZAC - Avenue de Satolas Green  
69330 Pusignan  
France  
Phone: 011-33-472-45-36-00  
Fax: 011-33-472-45-36-99

## GERMANY

**Stryker Howmedica GmbH**  
Dr. Homer Stryker Platz 1  
47228 Duisburg  
Germany  
Phone: 011-49-2065-837-0  
Fax: 011-49-2065-837-837

## GREECE

**Stryker Hellas EPE**  
455 Messogion ave  
153 43 Agia Paraskevi  
Athens, Greece  
Phone: 011-30-2-10-600-32-22  
Fax: 011-30-2-10-600-48-12

## ITALY

**Stryker Italia Srl**  
Via Ghisalba 15B  
00188 Roma  
Italy  
Phone: 011-39-06-33-05-41  
Fax: 011-39-06-33-614-067

## MIDDLE EAST / NORTH AFRICA

**Stryker Osteonics SA**  
Twin Towers  
11th Floor, Suite 1101 & 1102  
P.O. Box 41446  
Baniyas Road  
Dubai, Deira, UAE  
Phone: 011-97-14-222-2842  
Fax: 011-97-14-224-7381

## NETHERLANDS

**Stryker Nederlands BV**  
(P.O. Box 13, 4180 BA Waardenburg)  
4181 CD Waardenburg  
The Netherlands  
Phone: 011-31-418-569-700  
Fax: 011-31-418-569-777

## NORWAY

**Stryker Norway**  
Norsk Fillial  
Nedre Vollgate 3  
0158 Oslo  
Norway  
Phone: 011-47-22-42-22-44  
Fax: 011-47-22-42-22-54

## POLAND

**Stryker Polska Sp. ZO.O**  
Kolejowa 5/7  
01-217 Warszawa  
Poland  
Phone: 011-48-22-434-88-50  
Fax: 011-48-22-434-88-60

## PORTUGAL

**Stryker Portugal Produtos Medicos, LTDA.**  
Avenida Marechal Gomes da Costa, 35  
1800-255 Lisboa  
Portugal  
Phone: 011-35-1-21-839-49-10  
Fax: 011-35-1-21-839-49-19

## ROMANIA

**Stryker Osteonics Romania S.R.L.**  
19, Leonida Str.  
District 2  
7000 Bucharest  
Romania  
Phone: 011-40-2-12-12-11-22  
Fax: 011-40-2-12-12-11-33

## SOUTH AFRICA

**Stryker Osteonics PTY. LTD.**  
3 Susan Street  
Strydom Park  
Johannesburg, 2194 - South Africa  
Mailing Address:  
P.O. Box 48039  
2129 Roosevelt Park  
Johannesburg, South Africa  
Phone: 011-27-11-791-4644  
Fax: 011-27-11-791-4696

## SPAIN

**Stryker Iberia SL**  
c/Sepulveda n17  
28108 Alcobendas  
Madrid  
Spain  
Phone: 011-34-91-728-35-00  
Phone: 011-34-91-358-20-44  
Fax: 011-34-91-358-07-48

## SWEDEN

**Stryker AB/Scandinavia**  
Box 50425  
SE-204 14 Malmo  
Sweden  
Phone: 011-46-40-691-81-00  
Fax: 011-46-40-691-81-91

# Table of Contents

Symbols and Definitions . . . . .	<a href="#">1-3</a>
Symbols . . . . .	<a href="#">1-3</a>
Warning/Caution/Note Definition . . . . .	<a href="#">1-3</a>
Introduction . . . . .	<a href="#">1-4</a>
Product Description . . . . .	<a href="#">1-4</a>
Intended Use of Product. . . . .	<a href="#">1-4</a>
Specifications . . . . .	<a href="#">1-5</a>
Contact Information . . . . .	<a href="#">1-6</a>
Serial Number Location . . . . .	<a href="#">1-6</a>
Product Illustration . . . . .	<a href="#">1-7</a>
Summary of Safety Precautions . . . . .	<a href="#">1-8</a>
Setup Procedures. . . . .	<a href="#">1-12</a>
Cot Fastener Installation . . . . .	<a href="#">1-13</a>
Vehicle Safety Hook Selection . . . . .	<a href="#">1-15</a>
Vehicle Safety Hook Installation . . . . .	<a href="#">1-16</a>
Vehicle Configuration . . . . .	<a href="#">1-16</a>
Required Hardware for Installation of the Safety Hook (Not Supplied) . . . . .	<a href="#">1-16</a>
Front to Back Positioning of the Safety Hook . . . . .	<a href="#">1-17</a>
Side to Side Positioning of the Safety Hook . . . . .	<a href="#">1-18</a>
Installing the Safety Hook. . . . .	<a href="#">1-18</a>
Adjusting Cot Load Height. . . . .	<a href="#">1-19</a>
Cot Positions . . . . .	<a href="#">1-20</a>
Operation Guide . . . . .	<a href="#">1-21</a>
Operating Guidelines . . . . .	<a href="#">1-21</a>
Proper Lifting Techniques. . . . .	<a href="#">1-21</a>
Transferring the Patient to the Cot . . . . .	<a href="#">1-22</a>
Rolling the Cot. . . . .	<a href="#">1-22</a>
Adjusting the Height of the Cot with Two Operators. . . . .	<a href="#">1-23</a>
Adjusting the Height of an Empty Cot with One Operator. . . . .	<a href="#">1-24</a>
Loading or Unloading the Cot. . . . .	<a href="#">1-25</a>
Loading or Unloading the Cot with the Power-LOAD Option . . . . .	<a href="#">1-25</a>
Loading the Cot into a Vehicle with Two Operators . . . . .	<a href="#">1-26</a>
Loading an Empty Cot into a Vehicle with One Operator . . . . .	<a href="#">1-27</a>
Unloading the Cot from a Vehicle with Two Operators. . . . .	<a href="#">1-28</a>
Unloading an Empty Cot from a Vehicle with One Operator. . . . .	<a href="#">1-29</a>
Using Additional Assistance . . . . .	<a href="#">1-30</a>
Operating the Siderails. . . . .	<a href="#">1-31</a>
Operating the Backrest . . . . .	<a href="#">1-31</a>
Operating the Retractable Head Section. . . . .	<a href="#">1-32</a>
Adjusting the Footrest . . . . .	<a href="#">1-33</a>
Adjusting the Optional Knee Gatch . . . . .	<a href="#">1-34</a>
Operating the Optional Wheel Lock(s) . . . . .	<a href="#">1-35</a>
Using Restraint Straps . . . . .	<a href="#">1-36</a>
Using the Restraint Belt Extension . . . . .	<a href="#">1-39</a>

# Table of Contents

Optional Accessories . . . . .	<a href="#">1-40</a>
Installing the Base Storage Net . . . . .	<a href="#">1-41</a>
Using the Defibrillator Platform . . . . .	<a href="#">1-41</a>
Using the Equipment Hook . . . . .	<a href="#">1-42</a>
Using the Head Extension with Pillow . . . . .	<a href="#">1-42</a>
Operating the Optional Two-Stage I.V. Pole . . . . .	<a href="#">1-43</a>
Operating the Optional Three-Stage I.V. Pole . . . . .	<a href="#">1-44</a>
Using the Kickstand for Dialysis Scale . . . . .	<a href="#">1-45</a>
Attaching an Oxygen Bottle to an Oxygen Bottle Holder . . . . .	<a href="#">1-46</a>
Using the Retractable Head Section Oxygen Bottle Holder . . . . .	<a href="#">1-47</a>
Attaching the Pedi-Mate® Infant Restraint System . . . . .	<a href="#">1-48</a>
Installing the Backrest Storage Pouch . . . . .	<a href="#">1-50</a>
Installing the Head End Storage Flat . . . . .	<a href="#">1-51</a>
Using the Transfer Flat . . . . .	<a href="#">1-51</a>
Cleaning . . . . .	<a href="#">1-52</a>
Washing Procedure . . . . .	<a href="#">1-52</a>
Washing Limitations . . . . .	<a href="#">1-52</a>
Removal of Iodine Compounds . . . . .	<a href="#">1-53</a>
Preventative Maintenance . . . . .	<a href="#">1-54</a>
Checklist . . . . .	<a href="#">1-54</a>
Regular Inspection and Adjustments . . . . .	<a href="#">1-55</a>
Maintenance Record . . . . .	<a href="#">1-57</a>
Training Record . . . . .	<a href="#">1-58</a>
Quick Reference Replacement Parts List . . . . .	<a href="#">1-59</a>
Service Information . . . . .	<a href="#">1-60</a>
Backrest Adjustment . . . . .	<a href="#">1-60</a>
Wheel Locking Force Adjustment . . . . .	<a href="#">1-61</a>
Cot Retaining Post Adjustment . . . . .	<a href="#">1-62</a>
Cot Retaining Post Replacement . . . . .	<a href="#">1-63</a>
Cot Retaining Post Screw Replacement . . . . .	<a href="#">1-63</a>
Headsection Replacement . . . . .	<a href="#">1-64</a>
Backrest Gas Cylinder Replacement . . . . .	<a href="#">1-64</a>
Inner, Inner Tube Replacement . . . . .	<a href="#">1-65</a>
Outer, Inner Tube Replacement . . . . .	<a href="#">1-66</a>
Outer, Outer Tube Replacement . . . . .	<a href="#">1-67</a>
Inner, Outer Tube Replacement . . . . .	<a href="#">1-68</a>
Siderail Assembly Replacement . . . . .	<a href="#">1-69</a>
Warranty . . . . .	<a href="#">1-70</a>
Stryker EMS Return Policy . . . . .	<a href="#">1-71</a>
Return Authorization . . . . .	<a href="#">1-71</a>
Damaged Merchandise . . . . .	<a href="#">1-71</a>
International Warranty Clause . . . . .	<a href="#">1-71</a>
Patent Information . . . . .	<a href="#">1-71</a>

## SYMBOLS

	Attention, consult accompanying documents
	Safe Working Load Symbol
	Pinch Point

## WARNING/CAUTION/NOTE DEFINITION

The words WARNING, CAUTION and NOTE carry special meanings and should be carefully reviewed.

### **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

### **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

### **NOTE**

Provides special information to make maintenance easier or important instructions clearer.

# Introduction

---

This manual is designed to assist you with the operation and maintenance of the Stryker **Performance-PRO™ XT** cot. Read this manual thoroughly before using the equipment or beginning maintenance on it. To ensure safe operation of this equipment, it is recommended that methods and procedures be established for educating and training staff on the safe operation of this cot.

## PRODUCT DESCRIPTION

The Stryker Model 6086 **Performance-PRO™ XT** is a manual ambulance cot that consists of a platform mounted on a wheeled X-frame designed to support and transport a maximum weight of 700 pounds in pre-hospital and hospital environments. The device is collapsible for use in emergency vehicles and has an adjustable load height feature to allow the device to be set to different ambulance deck heights for proper body mechanics during loading and unloading. Duplicate foot-end controls on the upper and lower lift bars accommodate different operator positions or sizes and the side release handle allows a single operator to raise and lower an unoccupied cot. The device is equipped with the following: a retractable head section for 360-degree mobility in any height position, side rails, patient securement straps, an adjustable pneumatic backrest and various optional accessories that assist with transport of the patient. Maximum patient comfort is attainable with the three different litter positions of shock, flat leg and optional knee gatch positioning.

## INTENDED USE OF PRODUCT

The Stryker Model 6086 **Performance-PRO™ XT** is a non-powered wheeled stretcher, which is intended to support and transport the entire body of a traumatized, ambulatory or non-ambulatory human patient (includes infants and adults). The device is designed to support patients in a supine (horizontal) or sitting position and facilitate the transportation of associated medical equipment (i.e. oxygen bottles, monitors, and/or pumps) in emergency/transport vehicles. This ambulance cot is intended to be used in pre-hospital and hospital environments, in emergency and non-emergency applications. It is rated to a maximum capacity of 700 pounds (sum of the patient, mattress and accessory weight) and the intended operators of the device are trained professionals including: emergency medical service and medical care center personnel, as well as medical first responders. The expected service life of the product is 7 years. Ambulance cots are intended for transportation purposes. They are not intended for extended stay or to be used as hospital beds. They are also not intended to be used in devices which modify air pressure, such as hyperbaric chambers.

## SPECIFICATIONS

 <b>Safe Working Load</b> <b>Note:</b> Safe Working Load indicates the sum of the patient mattress and accessory weight.		700 pounds	317,5 kg
Backrest Articulation/Shock Position		0° to 73° / +15°	
Overall Length/Minimum Length/Width		80" / 64" / 23"	
Height <sup>1</sup>	Position 1	13.8"	35,1 cm
	Position 2	22"	55,9 cm
	Position 3	25.8"	65,5 cm
	Position 4	28.1"	71,4 cm
	Position 5	31.9"	81 cm
	Position 6	34.6"	87,9 cm
	Position 7 (LOW)	37.3"	94,7 cm
	Position 8 (MID)	40"	101,6 cm
	Position 9 (HIGH)	42.2"	107,2 cm
Weight <sup>2</sup>		89 lb	40,37 kg
Caster Diameter/Width		6" / 2"	
Minimum Operators Required for an Occupied Cot		2	
Minimum Operators Required for an Unoccupied Cot		1	
Recommended Fastener Systems		Model 6370 or 6377 Floor Mount Type Model 6371 Wall Mount Type	
Maximum Loading Height <sup>3</sup>		Up to 34"	Up to 86,4 cm
Single Wheel Lock / Double Wheel Lock		Optional	

<sup>1</sup> Height is measured from the bottom of the mattress at seat section to ground level.

<sup>2</sup> Cot is weighed without mattress and restraints.

<sup>3</sup> Load wheel height can be set between 27.25" (69,2 cm) and 34" (86,4).

Stryker reserves the right to change specifications without notice.

The **Performance-PRO™ XT** is designed to conform to the Federal Specification for the Star-of-Life Ambulance (KKK-A-1822).

The **Performance-PRO™ XT** is designed to be compatible with competitive cot fastener systems.

Patents pending.

**The yellow and black color scheme is a proprietary trademark of Stryker Corporation.**

# Introduction

## CONTACT INFORMATION

Contact Stryker Customer Service or Technical Support at: (800) 327-0770 or (269) 324-6500.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

## SERIAL NUMBER LOCATION

Please have the serial number (Figure 1) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.



Figure 1: Cot Serial Number & Location

## PRODUCT ILLUSTRATION

**Figure 2: Cot Components**

[Return To Table of Contents](#)

# Summary of Safety Precautions

---

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

## WARNING

- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
  - Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
  - Do not modify the cot or any components of the cot. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty (see [page 1-70](#)).
  - It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Cot Fastener system meets the installation specifications listed on [page 1-14](#). Injury may result if a non-compatible cot is used in the Stryker Cot Fastener system.
  - Have the vehicle safety hook installed by a certified mechanic. Improper safety hook installation can cause injury to the patient or operator and/or damage to the cot.
  - Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-15](#).
  - The face of the safety hook that engages the safety bar should be located at least 3-3/4" from the leading edge of the door sill. After installation, verify that the cot legs lock into the load position without contacting the vehicle bumper.
  - To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.
  - Verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle or injury to the patient or operator and/or damage to the cot may occur.
  - The cot must have at least 5/8" of clearance between the vehicle bumper and the cot to disengage the safety bar when unloading the cot from the vehicle. Verify that the cot legs lock into the load position before disengaging the safety bar from the safety hook. Failure to properly lock the cot height into position can cause injury to the patient or operator and/or damage to the cot.
  - Before placing the cot into service, confirm that the cot load height is set correctly for your vehicle.
  - Practice changing height positions and loading the cot until operation of the product is fully understood. Improper use can cause injury.
  - Do not allow untrained assistants to assist in the operation of the cot. Untrained technicians/assistants can cause injury to the patient or themselves.
  - Do not ride on the base of the cot. Damage to the product could occur, resulting in injury to the patient or operator.
  - Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. Transporting the cot in a lowered position, head or foot end first, minimizes the potential of a cot tip.
  - Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.
  - Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.
  - Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the product.
  - Never apply the optional wheel lock(s) while a patient is on the cot. Tipping could occur if the cot is moved while the wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
  - Siderails are not intended to serve as a patient restraint device. See [page 1-36](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.
  - High obstacles such as curbing, steps or rough terrain can cause the cot to tip, possibly causing injury to the patient or operator.
  - If the cot is equipped with the optional kickstand, make sure that the kickstand remains in the retracted position and does not engage during transport.
-

## WARNING

- Transporting the cot in lower positions reduces the potential of a cot tip. If possible, obtain additional assistance or take an alternate route.
- When operating the side release handle, keep hands away from the foot end release handle to avoid injury.
- If lowering the cot to the lowest position (position 1), remove your foot from the base tube or injury could result.
- Power-LOAD is designed to be compatible with the 6085/6086 Performance-PRO XT, 6500/6506 Power-PRO XT, and 6510/6516 Power-PRO IT cots with the Power-LOAD option only. In certain situations, you can use Power-LOAD as a standard antler for most X-frame cots, but a rail clamp assembly is required for all cots without the Power-LOAD option.
- It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Model 6390 Power-LOAD system is a Power-LOAD compatible cot. Injury may result if a non-compatible cot is used in the Stryker Model 6390 Power-LOAD system.
- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading/unloading a patient. Injury to the patient or operator could result.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- To avoid injury, always verify that the head section is locked into place prior to operating the cot.
- Be sure that the undercarriage has engaged and is locked before removing the loading wheels from the patient compartment floor of the vehicle. An unlocked undercarriage will not support the cot and injury to the patient or operator could result.
- Siderails are not intended to serve as a patient restraint device. See [page 1-36](#) for proper restraint strap usage. Failure to utilize the siderails properly could result in patient injury.
- Do not attempt to load the cot into the patient compartment with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the product.
- Never install or use wheel locks on a cot with excessively worn wheels. Installing or using wheel locks on wheels with less than a 6" diameter could compromise the holding ability of the wheel lock, resulting in injury to the patient or operator and/or damage to the cot or other equipment.
- Do not attach restraints to the base tubes, cross tubes, or fowler skin. Improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.
- Stryker recommends a two person operation when using the kickstand.
- Make sure that the patient weight is centered on the cot before using the kickstand.
- Engage the kickstand with your foot only.
- Lower cot height prior to engaging kickstand for increased stability.
- Make sure that the kickstand remains in the retracted position and does not engage during transport.
- Do not use the kickstand as a brake.
- Do not engage kickstand on a sloped surface.
- If the cot is equipped with the optional retractable head section oxygen bottle holder, use caution while the oxygen bottle holder is installed to avoid pinching your fingers between the fowler bracket and the oxygen bottle.
- To avoid accidental release of the Pedi-Mate® and possible injury to the infant, ensure that the buckle is located away from obstructions on the cot or accessories.
- When the optional head end storage flat is being used, ensure that it does not interfere with the operation of the retractable head section, safety bar and safety hook. Injury to the patient or operator could result.

# Summary of Safety Precautions

---

## WARNING

- When cleaning, use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the cot.
  - **SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY.** If the products described above are used to clean Stryker patient care equipment, measures must be taken to ensure that the cots are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the cots leaves a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.
  - Failure to properly clean or dispose of contaminated mattress or cot components increases the risk of exposure to bloodborne pathogens and may cause injury to the patient or the operator.
- 

## CAUTION

- Set the cot load height to the proper stop height prior to operation.
  - Installation of the safety hook should be done by a certified mechanic familiar with ambulance vehicle construction. Consult the vehicle manufacturer before installing the safety hook and be sure that the installation of the safety hook does not damage or interfere with the brake lines, oxygen lines, fuel lines, fuel tank or electrical wiring of the vehicle.
  - Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.
  - Do not allow the cot undercarriage to drop unassisted (commonly known as a “hot drop”) when removing the cot from the vehicle. Repeated hot dropping causes premature wear or damage to the cot.
  - Wheel lock(s) are only intended to help prevent the cot from rolling while unattended. Wheel locks may not provide sufficient resistance on all surfaces or under loads.
  - Ensure that the restraints are not entangled in the base frame when raising and lowering the cot.
  - The weight of the equipment in the base storage net (if equipped) must not exceed 20 lb (9 kg).
  - Be careful when retracting the base to avoid damaging items stored in the base storage net.
  - To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 35 lb (15.9 kg).
  - To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 40 lb (18 kg).
  - To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 40 lb (18 kg).
  - Do not use two head end oxygen bottle holders at the same time.
  - Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.
  - The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 20 lb (9 kg).
  - The weight of the equipment in the head end storage flat (if equipped) must not exceed 40 lb (18 kg).
  - **DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT.**
  - Maximum water temperature should not exceed 180°F/82°C.
  - Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 24 inches (61 cm) from the unit.
  - Allow cot to air dry.
  - Towel dry all casters and interface points.
  - Failure to comply with these instructions may invalidate any/all warranties.
-


Figure 3: Potential Pinch Points

 **WARNING**

Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

# Setup Procedures

---

Ensure that all shipping and packaging materials have been removed from the product(s) prior to use.

Unpack the cartons and check all items for proper operation. It is important that the cot is working properly before it is put into service. See “Figure 2: Cot Components” on [page 1-7](#) to identify all of the cot components.

Before placing the cot into service, check these components:

- Check for loose fasteners. Replace if loose. Reference all assembly drawings.
- All welds are intact (not cracked or broken)
- No bent or broken tubing or sheet metal
- No debris in wheels
- All wheels are secure, and rolling and swivelling properly
- Both siderails move and latch properly
- Backrest operates properly
- Optional accessories are intact and operate properly
- Height positioning latch function operates properly
- Cot secure in each height position (see [page 1-20](#))
- Undercarriage folds properly
- Retractable head section operates properly
- Safety bar operates properly
- Foot rest operates properly
- No rips or cracks in the mattress cover
- Body restraints are intact and operate properly
- Wheel lock(s) operate properly (optional equipment)
- Vehicle safety hook engages the safety bar so that the cot loads and unloads properly from the vehicle (see [page 1-15](#))
- Approved cot fastener (Stryker Model 6370/6377/6378/6379 or 6371 Cot Fastener – Not included) installed in the vehicle (see [page 1-13](#))
- Adjust cot load height (see [page 1-19](#))

The patient compartment of the vehicle in which the cot will be used must have a:

- Smooth rear edge for cot loading.
- Level floor large enough for the folded cot.
- Stryker Model 6370/6377/6378/6379 or 6371 Cot Fastener System or Stryker Model 6390 Power-LOAD (not included)
- 34” (86,4 cm) maximum loading height.
- Space to properly install the safety hook.

**Note:** Loose items or debris on the patient compartment floor can interfere with the operation of the safety hook and cot fastener. Keep the patient compartment floor clear.

When necessary, modify the vehicle to fit the cot. Do not modify the cot.

---

 **WARNING**

- Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
  - Do not modify the cot or any components of the cot. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty (see [page 1-70](#)).
-

# Cot Fastener Installation

---

**Note:** The Cot Fastener Installation instructions on [page 1-13](#) through [page 1-18](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

The Stryker Cot Fastener Systems are designed to be compatible only with cots which conform to the installation specifications listed on [page 1-14](#).

---

 **WARNING**

It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Cot Fastener Systems meets the installation specifications listed on [page 1-14](#). Injury may result if a non-compatible cot is used in the Stryker Fastener System.


---

**Note:** Adjustment of the rail clamp assembly may be required in order to compensate for any variation in the cot retaining post position depending on the cot manufacturer and model number.


For more information about the Stryker Cot Fastener Systems, see the Cot Fastener Operations/Maintenance Manual.

# Cot Fastener Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual (6390-009-001) for installation instructions.

**Figure 4: Installation Specifications - Floor Mount Fastener**

**Figure 5: Installation Specifications - Wall Mount Fastener**

# Vehicle Safety Hook Selection

**Note:** The Vehicle Safety Hook Selection and Installation instructions on [page 1-15](#) through [page 1-18](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions. Power-LOAD ships and is installed with its own safety hook, thus no additional hook is needed.

The vehicle safety hook is a device that ships with the cot. The cot safety bar and vehicle safety hook are designed to keep the cot from being accidentally removed from the vehicle and to provide increased operator assurance and confidence when loading and unloading. The safety hook was designed for compatibility and proper operation when loading and unloading the cot from a vehicle that is compliant with Federal Regulation KKK-A-1822.

Stryker offers three different types of safety hooks that are ordered and shipped with your cot. These safety hook types are designed to meet the needs of various emergency vehicle configurations, specifically the length and location of the floor structure support that is located in the rear of the vehicle.

Consider the following information when selecting which safety hook is appropriate for your vehicle configuration:

- Determine the location of the floor structure support where there is adequate room to mount the safety hook.
- Ensure that the safety hook can be securely mounted into the back of the vehicle while providing adequate bumper clearance to allow the cot to be loaded and unloaded from the vehicle.
- Note the differences in vehicle design. Each safety hook provides a different mounting location option to maintain the appropriate distance between the face of the safety hook and the edge of the door sill.


Figure 6: Safety Hook Types

Due to the differences in vehicle dimensions and the floor structure support locations, each safety hook requires a different mounting location. See “Vehicle Safety Hook Installation” to determine the correct positioning for safety hook installation.

**Note:** When replacing an existing safety hook with a new style, adjust the mounting location to maintain the proper position of the safety hook face.

# Vehicle Safety Hook Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## VEHICLE CONFIGURATION

According to federal regulations (reference KKK-A-1822), the bumper height of the vehicle shall be installed equidistant  $\pm 5$  cm (2 inches) from the vehicle floor to the ground level, which is defined as the vehicle deck height. The bumper step shall have a minimum depth of 13 cm (5 inches) and a maximum depth of 25 cm (10 inches). If the bumper depth is greater than 18 cm (7 inches), then the bumper must be able to fold. Installation of the safety hook into any vehicle compliant with this federal specification provides adequate clearance for the cot base to lower to its fully extended position. The cot is compatible with all vehicle deck heights (see specifications for maximum load height) as long as the vehicle meets the federal specifications that are outlined in KKK-A-1822.


Figure 7: Vehicle Deck Height


Figure 8: Vehicle Deck Height

### CAUTION

- Set the cot load height to the proper stop height prior to operation.
- Installation of the safety hook should be done by a certified mechanic familiar with ambulance vehicle construction. Consult the vehicle manufacturer before installing the safety hook and be sure that the installation of the safety hook does not damage or interfere with the brake lines, oxygen lines, fuel lines, fuel tank or electrical wiring of the vehicle.

## REQUIRED HARDWARE FOR INSTALLATION OF THE SAFETY HOOK (NOT SUPPLIED)

- (2) Grade 5, Minimum 1/4"-20 Socket Head Cap Screws\* for the short or long safety hook
- (2) Grade 5, Minimum 1/4"-20 Flat Socket Head Cap Screws\* for the J hook
- (2) Flat Washers
- (2) Lock Washers
- (2) 1/4"-20 Nuts

\* The length of the socket head cap screws depends on the thickness of the vehicle floor. Use screws that are long enough to go completely through the patient compartment floor, washer and nut by at least two full threads.

# Vehicle Safety Hook Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## WARNING

- Have the vehicle safety hook installed by a certified mechanic. Improper safety hook installation can cause injury to the patient or operator and/or damage to the cot.
- Failure to install the safety hook can cause injury to the patient or operator.
- The face of the safety hook that engages the safety bar should be located at least 3-3/4" from the leading edge of the door sill. After installation, verify that the cot legs lock into the load position without contacting the vehicle bumper.
- To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.

**Note:** Stryker recommends that, prior to installation, the certified mechanic plan the placement of the safety hook in the rear of the vehicle.

Before installing the safety hook into your vehicle, check the front to back and side to side positioning when unloading and loading the cot to ensure that the safety hook will be installed properly. The cot safety bar must engage the safety hook every time, regardless of cot position.

## FRONT TO BACK POSITIONING OF THE SAFETY HOOK

1. Select the appropriate safety hook for your vehicle configuration.
2. Position the safety hook at least 3-3/4" from the leading edge of the door sill.
3. Ensure that the safety hook can be securely mounted into the back of the vehicle while providing adequate bumper clearance to allow the cot to be loaded and unloaded from the vehicle.
4. See "Side to Side Positioning of the Safety Hook" to confirm the side to side placement.


Figure 9: Safety Hook Placement

[Return To Table of Contents](#)

# Vehicle Safety Hook Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## SIDE TO SIDE POSITIONING OF THE SAFETY HOOK


1. Remove the cot from the fastener and unload it from the vehicle.
2. While the cot is being removed, note the position of the load wheels and the safety bar.
3. Mark the center of the cot safety bar on the vehicle floor.
4. Verify that the position marked in Step 3 is where the safety bar engages the safety hook every time when unloading the cot in a variety of positions (all the way to the left and all the way to the right), regardless of cot position.
  - If the cot safety bar does not engage the safety hook in any of these positions (left, center, or right), modify the vehicle, not the cot or safety hook.
  - If the cot safety bar engages the safety hook every time, install the safety hook.

## INSTALLING THE SAFETY HOOK

1. Determine the correct safety hook front to back and side to side positioning, so the cot safety bar engages the safety hook every time.
2. Drill the holes for the screws.
3. Fasten the safety hook to the patient compartment floor and verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle.

### WARNING

- Verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle or injury to the patient or operator and/or damage to the cot may occur.
- The cot must have at least 5/8" of clearance between the vehicle bumper and the cot to disengage the safety bar when unloading the cot from the vehicle. Verify that the cot legs lock into the load position before disengaging the safety bar from the safety hook. Failure to properly lock the cot height into position can cause injury to the patient or operator and/or damage to the cot.

**Figure 10: Safety Hook Placement  
(For Reference Only)**

**Figure 11: Safety Bar Engaging Safety Hook**

# Adjusting Cot Load Height

Before placing the cot into service, confirm that the cot load height is set correctly for your vehicle. The cot load height can be adjusted to match the height of the vehicle deck. If the cot does not line up correctly, adjustments may need to be made to the cot load height.

## WARNING

Before placing the cot into service, confirm that the cot load height is set correctly for your vehicle.

To check the load height of the **Performance-PRO™ XT** cot:

1. Roll the cot up to the loading area of your vehicle.
2. Compare the difference between the deck height of the vehicle and the load height of the cot.
3. Select HIGH, MID or LOW depending on the cot load height requirements (see Figure 12). For example:
  - The HIGH marking on the rack is recommended for vehicle deck heights above 32 inches.
  - The MID marking on the rack is recommended for vehicle deck heights between 30 and 32 inches.
  - The LOW marking on the rack is recommended for vehicle deck heights below 30 inches.
4. Verify that the safety hook always engages the cot safety bar, regardless of how the cot is unloaded from the vehicle. If the safety bar misses the safety hook, select the next lower height setting.




Figure 12: Cot Load Height


# Cot Positions


Position 1 - Use for patient transfer


Position 2 - Use for patient transfer/cot rolling


Position 3 - Use for patient transfer/cot rolling


Position 4 - Use for patient transfer/cot rolling


Position 5 - Use for patient transfer/cot rolling


Position 6 - Use for patient transfer/cot rolling

Position 7 - Use for patient transfer/cot rolling (LOW)

Position 8 - Use for patient transfer/cot rolling (MID)

Position 9 - Use for patient transfer/cot rolling (HIGH)

## OPERATING GUIDELINES

- Use the cot only as described in this manual.
- Read all labels and instructions on the cot before using the cot.
- Loading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator. If additional assistance is needed, see “Using Additional Assistance” on [page 1-30](#).
- Do not adjust, roll or load the cot into a vehicle without advising the patient. Stay with the patient and control the cot at all times.
- The cot can be transported in any position. Stryker recommends that the operators transport the patient in the lowest comfortable position to maneuver the cot.
- Only use the wheel lock(s) during patient transfer or without a patient on the cot.
- Do not leave wheel lock(s) engaged while transporting the cot. Failure to do so may cause wheel damage.
- Always use the restraint straps.
- Use properly trained helpers, when necessary, to control the cot.

---

### **WARNING**

- Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
- Practice changing height positions and loading the cot until operation of the product is fully understood. Improper use can cause injury.
- Do not allow untrained assistants to assist in the operation of the cot. Untrained technicians/assistants can cause injury to the patient or themselves.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- Do not ride on the base of the cot. Damage to the product could occur, resulting in injury to the patient or operator.
- Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. Transporting the cot in a lowered position, head or foot end first, minimizes the potential of a cot tip.
- Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.

---

### **CAUTION**

Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.

---

## PROPER LIFTING TECHNIQUES

When lifting the cot and patient, there are five basic guidelines to help you avoid injury:

- Keep your hands close to your body.
- Keep your back straight.
- Coordinate your movements with your partner and lift with your legs.
- Avoid twisting.
- Always operate the cot as described in this manual.

## TRANSFERRING THE PATIENT TO THE COT

### To transfer the patient to the cot:

1. Roll the cot to the patient.
2. Place the cot beside the patient and raise or lower the cot to the level of the patient.
3. Lower the siderails and open the restraint straps.
4. Transfer the patient to the cot using accepted EMS procedures.
5. Use all the restraint straps to secure the patient to the cot (see [page 1-36](#)).
6. Adjust the backrest and foot rest as necessary.

**Note:** When transferring larger patients, use of the Transfer Flat (6005-001-001) is recommended.

---

### **WARNING**

- Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.
  - Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the product.
  - Never apply the optional wheel lock(s) while a patient is on the cot. Tipping could occur if the cot is moved while the wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
  - Siderails are not intended to serve as a patient restraint device. See [page 1-36](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.
- 

## ROLLING THE COT

### When rolling the cot:

- Make sure that all of the restraint straps are securely buckled around the patient (see [page 1-36](#)).
- Position an operator at the foot end and one at the head end of the cot **at all times** when rolling the cot with a patient on it.
- Approach door sills and/or other low obstacles squarely and lift each set of wheels over the obstacle separately.

---

### **WARNING**

- High obstacles such as curbing, steps or rough terrain can cause the cot to tip, possibly causing injury to the patient or operator.
  - If the cot is equipped with the optional kickstand, make sure that the kickstand remains in the retracted position and does not engage during transport.
  - Transporting the cot in lower positions reduces the potential of a cot tip. If possible, obtain additional assistance or take an alternate route.
-

## ADJUSTING THE HEIGHT OF THE COT WITH TWO OPERATORS

### WARNING

- Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- When operating the side release handle, keep hands away from the foot end release handle to avoid injury.

You can raise or lower an unoccupied cot with one operator. If a patient is on the cot, a minimum of **two (2) trained operators** (one located at each end of the cot) are required to raise or lower the cot.

### To raise or lower the cot from the ends:

1. The operator at the foot end of the cot squeezes the release handle (A or B) while a secure grip is maintained on the lifting bars (see Figure 13).
2. Both operators must lift the cot until the weight is off the latching mechanism (approximately 1/4").
3. The operator at the foot end squeezes and holds the release handle and both operators then raise or lower the cot together. The handle is released when the desired position is reached. Both operators should maintain a secure grip on the litter frame until the latching mechanism is securely locked into position.

### To raise or lower the cot from the sides:

1. Check the cot to determine if the side release handle is on the patient left or right side.
2. The operator on the patient's right or left (depending on the location of the release handle) reaches the release handle at the midpoint of the litter (C). Both operators must lift the cot until the weight is off the latching mechanism (approximately 1/4") (see Figure 13).
3. The operator at the patient's right or left (depending on the location of the release handle) squeezes and holds the release handle. Both operators then raise or lower the cot together. The handle is released when the desired position is reached. Both operators should maintain a secure grip on the litter frame until the latching mechanism is securely locked into position.


Figure 13: Adjusting the Cot Height

# Cot Operation

## ADJUSTING THE HEIGHT OF AN EMPTY COT WITH ONE OPERATOR

### To raise or lower the cot from the foot end:

1. Standing at the foot end of the cot, grasp the lower foot end lift tube.
2. Tip the cot up onto the load wheels (see Figure 14).
3. Squeeze and hold the release handle and raise or lower the foot end to the desired position. The handle is released when the desired position is reached.
4. Lower the cot back onto the four base wheels (see Figure 15).

Figure 14: Cot Tipped on Load Wheels

Figure 15: Cot Lowered to Ground

### To raise or lower the cot from the side:

1. Place one foot on the outer base tube.
2. Grasp the side release handle with one hand. Place your other hand on the outer support rail to help stabilize the cot (see Figure 16).
3. Squeeze the side release handle and raise or lower the cot to the desired position. The handle is released when the desired position is reached (see Figure 17).

### WARNING

If lowering the cot to the lowest position (position 1), remove your foot from the base tube or injury could result.

Figure 16: Holding Outer Support Rail

Figure 17: Lowering Cot from Side

## LOADING OR UNLOADING THE COT

The cot loading and unloading instructions on [page 1-26](#) through [page 1-29](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6086 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for loading and unloading instructions.

## LOADING OR UNLOADING THE COT WITH THE POWER-LOAD OPTION

The Model 6086 **Performance-PRO™ XT** cot is fully compatible with the Model 6390 Power-LOAD system if it is ordered with the Power-LOAD option or compatibility kit.

For more information about using your Power-LOAD compatible cot, see the Power-LOAD Operations/Maintenance Manual.

---

### **WARNING**

- Power-LOAD is designed to be compatible with the 6085/6086 Performance-PRO XT, 6500/6506 Power-PRO XT, and 6510/6516 Power-PRO IT cots with the Power-LOAD option only. In certain situations, you can use Power-LOAD as a standard antler for most X-frame cots, but a rail clamp assembly is required for all cots without the Power-LOAD option.
  - It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Model 6390 Power-LOAD system is a Power-LOAD compatible cot. Injury may result if a non-compatible cot is used in the Stryker Model 6390 Power-LOAD system.
-

## LOADING THE COT INTO A VEHICLE WITH TWO OPERATORS

### WARNING

- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame. (See [page 1-16](#) for safety hook installation instructions.)
- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-16](#).

Figure 18: Safety Bar Engaging the Safety Hook

### To load the cot into a vehicle with two operators:

1. Place the cot in a loading position (any position where the loading wheels meet the vehicle floor height). Roll the cot to the open door of the patient compartment. Lift the vehicle bumper to the raised position (if equipped).
2. Push the cot forward until the loading wheels are on the patient compartment floor and the safety bar passes the safety hook as shown in Figure 18.
3. For maximum clearance to lift the base, pull the cot back until the safety bar engages the safety hook. Operator 2 should verify that the bar engages the safety hook.
4. **Operator 1** – Grasp the cot frame at the foot end. Lift the foot end of the cot until the weight is off of the latching mechanism. Squeeze and hold the release handle (A or B, as shown in Figure 13 on [page 1-23](#)).
5. **Operator 2** – Stabilize the cot by placing your hand on the outer rail (C). Grasp the base frame where indicated (D). After the foot end operator lifts the cot and squeezes the release handle, raise the undercarriage until it stops in the uppermost position and hold it there (see Figure 19). The foot end operator should release the handle to lock the base in the retracted position.
6. **Both Operators** – Push the cot into the patient compartment (see Figure 20), engaging the cot fastener (not included).

Figure 19: 2 Operators with One Lifting the Base

Figure 20: 2 Operators with Base Full Up

## LOADING AN EMPTY COT INTO A VEHICLE WITH ONE OPERATOR

### WARNING

- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading or unloading a patient. Injury to the patient or operator could result.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

### To load an empty cot into a vehicle with one operator:

1. Place the cot in a loading position (any position in which the load wheels meet the vehicle floor height).
2. Lift the vehicle bumper to the raised position (if equipped).
3. Roll the cot to the open door of the patient compartment.
4. Push the cot forward until the loading wheels are on the compartment floor and the safety bar passes the safety hook.
5. Pull the cot back until the safety bar engages the safety hook.
6. Grasp the cot frame at the foot end and squeeze and hold the release handle (see Figure 21).
7. Lower the foot end of the cot to the ground, making sure that the cot locks in position 1 (see Figure 22).
8. Lift the foot end of the cot until it is level with the compartment floor (see Figure 23).
9. Grasp the base of the cot with one hand and pull up the base of the cot towards the litter, reducing the space between the base and the litter.
10. Push the cot into the patient compartment by guiding it into the cot fastener.

### WARNING

Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.

Figure 21: Squeeze the Release Handle

Figure 22: Lower the Foot End of the Cot

Figure 23: Pull Up the Base of the Cot

## UNLOADING THE COT FROM A VEHICLE WITH TWO OPERATORS

### WARNING

- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-16](#).
- To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

### To unload the cot from a vehicle with two operators:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-13](#)).
3. Operator 1 – Grasp the cot frame. Pull the cot out of the patient compartment until the safety bar engages the safety hook (see [Figure 24](#)).
4. Operator 2 – Grasp the base frame where indicated, lift slightly, and lower the base frame to its fully extended position while operator 1 squeezes and holds the release handle (see [Figure 25](#)).
5. Operator 1 – Let go of the release handle and make sure that the undercarriage locks into place. Set the cot onto the ground.
6. Operator 2 – Disengage the safety bar from the safety hook by pushing the safety bar release lever forward.
7. Remove the cot loading wheels from the vehicle. Place the cot in any position, except full down for rolling.

### CAUTION

Do not allow the cot undercarriage to drop unassisted (commonly known as a “hot drop”) when removing the cot from the vehicle. Repeated hot dropping causes premature wear or damage to the cot.

### WARNING

Be sure that the undercarriage has engaged and is locked before removing the loading wheels from the patient compartment floor of the vehicle. An unlocked undercarriage will not support the cot and injury to the patient or operator could result.

Figure 24: 2 Operators with Base Full Up

Figure 25: 2 Operators with One Lowering the Base

## UNLOADING AN EMPTY COT FROM A VEHICLE WITH ONE OPERATOR

### WARNING

- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading or unloading a patient. Injury to the patient or operator could result.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

### To unload an empty cot from a vehicle with one operator:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-13](#)). Grasp the cot frame at the foot end; pull the cot from the vehicle until the safety bar engages the safety hook (see Figure 26).
3. Lower the foot end of the cot to the ground (see Figure 27).
4. Squeeze and hold the release handle (see Figure 28) and raise the foot end of the cot back to a level position with the compartment floor.
5. Disengage the safety bar from the safety hook by pushing the safety bar release lever forward and roll the cot out of the vehicle.

Figure 26: Pull the Base of the Cot

Figure 27: Lower the Foot End of the Cot

Figure 28: Squeeze the Release Handle

# Operation Guide

English

## USING ADDITIONAL ASSISTANCE

### IF EQUIPPED WITH THE RIGHT HAND RELEASE OPTION

	Changing Levels	Rolling	Loading/Unloading
Two Operators Two Helpers			
Two Operators Four Helpers			

### IF EQUIPPED WITH THE LEFT HAND RELEASE OPTION

	Changing Levels	Rolling	Loading/Unloading
Two Operators Two Helpers			
Two Operators Four Helpers			

[Return To Table of Contents](#)

## OPERATING THE SIDERAILS

To raise the siderails, as shown in Figure 29, lift up on the siderail until the latch clicks and the siderail locks into place.

To lower the siderails, squeeze handle (B) to release the siderail latch. Guide the siderail down toward the foot end until flat. Ensure that the siderails are lowered when a patient is being transferred to or from the cot.

### WARNING

Siderails are not intended to serve as a patient restraint device. See [page 1-36](#) proper restraint strap usage. Failure to use the siderails properly could result in patient injury.


Figure 29: Backrest Elevated and Siderails Raised

## OPERATING THE BACKREST

To raise the backrest, as shown in Figure 29, squeeze handle (A) for pneumatic assist in lifting the backrest to the desired height.

To lower the backrest, squeeze handle (A) and push down on the backrest frame until the backrest has reached the desired height.

## OPERATING THE RETRACTABLE HEAD SECTION

The head section telescopes from a first position suitable for loading the cot into an emergency vehicle to a second position retracted within the litter frame. When retracted, the cot can roll in any direction on the caster wheels even in the lowest position, allowing improved mobility and maneuverability.

### To extend the head section:

1. Grasp the outer rail with one hand for support and pull the handle (A), rotating the handle towards the head end of the cot to release the head section from the locked position
2. While holding the handle (A) in the released position, pull the head section away from the litter frame, lengthening the head section until it engages in the fully extended position.
3. Release handle (A) to lock the head section in the extended position.

### To retract the head section:

1. Grasp the outer rail with one hand for support and release the handle (A), rotate the handle towards the head end of the cot to release the head section from the locked position.
2. While holding the handle (A) in the released position, push the head section toward the litter frame, retracting the head section until it engages in the retracted position.
3. Release handle (A) to lock the head section in the retracted position.

### WARNING

- To avoid injury, always verify that the head section is locked into place prior to operating the cot.
- Do not attempt to load the cot into the patient compartment with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the product.


Figure 30: Head Section Extended


Figure 31: Head Section Retracted


Figure 32: Head Section Release Handles

## ADJUSTING THE FOOTREST

The footrest is adjustable to allow for elevation of the patient's legs (see Figure 33).

**To raise the footrest,** lift the foot rest frame (A) as high as possible until it locks into place. The support bracket engages automatically when released.

**To lower the footrest,** lift the foot rest frame (A) and, while holding the frame, lift up on the release handle (B) until the bracket disengages. Carefully lower the footrest until it rests flat.


Figure 33: Footrest Elevated

## ADJUSTING THE OPTIONAL KNEE GATCH

### To raise the knee gatch (see Figure 34):

1. Lift either of the red lifting loops (A) until the knee gatch is in its fully raised position.
2. Slowly lower the knee gatch to allow the support bracket to engage in the locking mechanism.
3. Check to be sure the lock is fully engaged before releasing the lifting loop.

### To lower the knee gatch:

1. Lift either of the red lifting loops to relieve pressure on the locking mechanism and while holding the loop, push on the red release handle (B) until the bracket disengages.
2. Carefully lower the knee gatch to the flat position.


Figure 34: Raised Knee Gatch


Figure 35: Raised Knee Gatch in Trend

### To raise the knee gatch in trend (see Figure 35):

1. Lift the foot rest frame (C) as high as possible until it locks into place.
2. The support bracket engages automatically when released.

### To lower the knee gatch in trend:

1. Lift the foot rest frame (C) and, while holding the frame, lift up on the release handle (B) until the bracket disengages.
2. Carefully lower the footrest until it rests flat.

## OPERATING THE OPTIONAL WHEEL LOCK(S)

To activate the optional wheel lock(s), press fully down on the pedal (A) as shown in Figure 36 until it stops and is resting firmly against the surface of the wheel.

To release the optional wheel lock(s), depress the upper face of the pedal with your foot or lift up with your toe under the pedal. The upper portion of the pedal will rest against the caster frame when the wheel lock is released.


Figure 36: Wheel Lock

### WARNING

- Never apply the optional wheel lock(s) while a patient is on the cot. Tipping could occur if the cot is moved while a wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
- Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the cot.
- Never install or use wheel locks on a cot with excessively worn wheels. Installing or using wheel locks on wheels with less than a 6" diameter could compromise the holding ability of the wheel lock, resulting in injury to the patient or operator and/or damage to the cot or other equipment.

### CAUTION

Wheel lock(s) are only intended to help prevent the cot from rolling while unattended and to aid in patient transfer. A wheel lock may not provide sufficient resistance on all surfaces or under loads.

## USING RESTRAINT STRAPS

### WARNING

Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.

Always secure the patient on the cot with all of the restraint straps.

To attach the restraint strap to the cot:

1. Wrap the restraint strap around the cot frame as shown in Figure 37.
2. Push the restraint strap buckle through the loop as shown in Figure 38.
3. Pull the buckle through the loop to secure the restraint strap to the cot as shown in Figure 39.

Figure 37: Wrap strap around cot

Figure 38: Push buckle through loop

Figure 39: Pull strap to tighten

4. Repeat steps 1-3 until all restraint straps are securely attached to the cot in the required attachment locations as shown in Figure 40.

Buckle the restraints across the patient's chest/shoulders, waist and legs. Note the attachment locations in Figure 41.

Keep the restraint straps buckled (as shown in Figure 40) when the cot is not being used with a patient to avoid damage to the buckles and straps.

Figure 40: Restraint Strap Locations

## USING RESTRAINT STRAPS (CONTINUED)


When attaching the restraint straps to the cot, the attachment points should provide both strong anchorage and proper restraint position while not interfering with equipment and accessories.

**⚠ WARNING**

Do not attach restraints to the base tubes, cross tubes, or fowler skin. Improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.

(Back View)

Shoulder/Chest Restraint Straps

Knee Restraint Straps

Foot Restraint Straps


Figure 41: All Straps (Front View)

[Return To Table of Contents](#)

## USING RESTRAINT STRAPS (CONTINUED)

### CAUTION

Ensure that the restraints are not entangled in the base frame when raising and lowering the cot.

When the cot is put into service, open the restraints and place them at either side of the cot until the patient is positioned on the cot mattress. Lengthen the restraints, buckle them around the patient and shorten them until the required tightness is achieved.

- **To open the restraint**, press the red button (A) on the front of the buckle “receiver”. This releases the buckle latch plate (B) which can then be pulled out of the receiver (Figure 42).
- **To close the restraint**, push the latch plate into the receiver until a “click” is heard. When fastening the chest restraint ensure that the latch plate passes through both links (C) on the shoulder strap (Figure 42).
- **To lengthen the restraint**, grasp the buckle latch plate, turn it at an angle to the webbing, then pull it out (Figure 43). A hemmed tab at the end of the webbing prevents the latch plate from coming off of the strap.
- **To shorten the restraint**, grasp the hemmed tab and pull the webbing back through the latch plate until the required tightness is achieved (Figure 44).


Figure 42: Buckling the Safety Restraints


Figure 43: Lengthening the Safety Restraint



Figure 44: Shortening the Safety Restraint

Whenever a restraint is buckled on a patient, verify that the latch plate is fully engaged and any extra webbing is not tangled in the cot or hanging loose.

Inspect the restraints **at least** once a month (more frequently if used heavily). Inspection should include checking for a bent or broken receiver or latch plate, torn or frayed webbing, etc. Any restraint showing wear or not operating properly **must** be replaced immediately.

## USING THE RESTRAINT BELT EXTENSION

Use the restraint belt extension, as shown in Figure 45, for extra length when buckling the lap belt around large patients.

**Figure 45: Attaching the Restraint Belt Extension**

# Optional Accessories

The accessories listed below can be purchased and installed on the **Performance-PRO™ XT** cot.

Accessory	Part Number	Operation Guide Page Number
Base Storage Net	6500-160-000	<a href="#">page 1-41</a>
Defibrillator Platform	6500-170-000	<a href="#">page 1-41</a>
Equipment Hook	6500-147-000	<a href="#">page 1-42</a>
Head Extension with Pillow	6100-044-000	<a href="#">page 1-42</a>
I.V. Pole Assembly, Two-Stage, Right	6500-210-000	<a href="#">page 1-43</a>
I.V. Pole Assembly, Three-Stage, Right	6500-215-000	<a href="#">page 1-44</a>
I.V. Pole Assembly, Two-Stage, Left	6500-211-000	<a href="#">page 1-43</a>
I.V. Pole Assembly, Three-Stage, Left	6500-216-000	<a href="#">page 1-44</a>
I.V. Pole Assembly, Two-Stage, Dual	6500-212-000	<a href="#">page 1-43</a>
I.V. Pole Assembly, Three-Stage, Dual	6500-217-000	<a href="#">page 1-44</a>
Kickstand Assembly	6085-002-000	<a href="#">page 1-45</a>
Oxygen Bottle Holder, Foot End	6500-140-000	<a href="#">page 1-46</a>
Oxygen Bottle Holder, Head End	6500-141-000	
Oxygen Bottle Holder, Removable	6080-140-000	
Oxygen Bottle Holder, Retractable Head Section	6085-046-000	<a href="#">page 1-47</a>
Pedi-Mate Restraint Package	6091-300-010	<a href="#">page 1-48</a>
Pocketed Back Rest Pouch	6500-130-000	<a href="#">page 1-50</a>
Storage Flat, Head End	6085-035-000	<a href="#">page 1-51</a>
Transfer Flat	6005-001-001	<a href="#">page 1-51</a>

## INSTALLING THE BASE STORAGE NET

To install the base storage net, wrap the Velcro® straps around the base tubes.

---

### CAUTION

- The weight of the equipment in the base storage net (if equipped) must not exceed 20 lb (9 kg).
- Be careful when retracting the base to avoid damaging items stored in the base storage net.

---

**Note:** The kickstand (p/n 6085-002-000) is not compatible with the optional base storage net (p/n 6500-160-000).

## USING THE DEFIBRILLATOR PLATFORM

See the Defibrillator Tray Operations/Maintenance manual for operation guide, safety precautions, cleaning, preventative maintenance, assembly drawings, and warranty information.

## USING THE EQUIPMENT HOOK

Use the equipment hook (A) (see Figure 46) to hang additional accessories or equipment, such as defibrillators or monitors.

### CAUTION

To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 35 lb (15.9 kg).


Figure 46: Equipment Hook

## USING THE HEAD EXTENSION WITH PILLOW

You can install the head extension assembly over the fowler to provide head end support.

Attach the pillow to the head extension by placing the support into the flap on the bottom of the pillow. Secure the pillow to the Velcro® on the bottom of the support.

**Note:** The head extension with pillow (p/n 6100-044-000) is not compatible with the optional equipment hook (p/n 6500-147-000) or optional fowler oxygen bottle holder (p/n 6500-141-000).

## OPERATING THE OPTIONAL TWO-STAGE I.V. POLE

### To use the two-stage I.V. pole (see Figure 48):

1. Lift and pivot the pole from the storage position and push down until it is locked into the receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the telescoping portion in place.
4. Hang the I.V. bags on the I.V. hook (D).
5. Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
6. Turn the lock actuator (B) clockwise to tighten.
7. Lift up and pivot the pole down into the storage position (see Figure 47).

### CAUTION

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 40 lb (18 kg).

**Note:** The dual two-stage I.V. poles (p/n 6500-212-000) are not compatible with either the patient right (6500-210-000) or patient left (6500-211-000) two-stage I.V. pole options.


Figure 47: Two-Stage I.V. Pole Storage Position


Figure 48: Two-Stage I.V. Pole

## OPERATING THE OPTIONAL THREE-STAGE I.V. POLE

**To use the three-stage I.V. pole (see Figure 50):**

1. Lift and pivot the pole from the storage position and push down until it is locked into the receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the bottom telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the bottom telescoping portion in place.
4. For a higher I.V. pole, pull up on section (D) until the spring clip (E) engages.
5. Hang I.V. bags on the I.V. hook (F).
6. To lower the I.V. pole, push in on the spring clip (E) and slide section (D) down into section (C). Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
7. Turn the lock actuator (B) clockwise to tighten.
8. Lift up and pivot the pole down into the storage position (see Figure 49).

### CAUTION

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 40 lb (18 kg).

**Note:** The three-stage I.V. pole options (p/n 6500-217-000) are not compatible with either the patient right (6500-215-000) or patient left (6500-216-000) two-stage I.V. pole options.


Figure 49: Three-Stage I.V. Pole Storage Position


Figure 50: Three-Stage I.V. Pole

## USING THE KICKSTAND FOR DIALYSIS SCALE

The kickstand is intended for weighing patients on a scale.

**Note:**

- The kickstand assembly is configured for an X-frame cot retention system only.
- The kickstand (p/n 6085-002-000) is not compatible with the optional base storage net (p/n 6500-160-000).

**⚠ WARNING**

- Stryker recommends a two person operation when using the kickstand.
- Make sure that the patient weight is centered on the cot before using the kickstand.
- Engage the kickstand with your foot only.
- Lower cot height prior to engaging kickstand for increased stability.
- Make sure that the kickstand remains in the retracted position and does not engage during transport.
- Do not use the kickstand as a brake.
- Do not engage kickstand on a sloped surface.

**To use the kickstand:**

1. Operator 1 engages the kickstand with their foot as shown in Figure 51.
2. Operator 2 lifts the foot end of the cot at a height sufficient to actuate the kickstand.
3. Both operators must make sure that the kickstand is in the forward locked position as shown in Figure 52.


Figure 51

Figure 52

**To release the kickstand:**

1. Operator 1 lifts the foot end of the cot until both wheels are off of the floor.
2. Operator 2 rolls the cot slightly forward to make sure that the kickstand retracts on its own as shown in Figure 53.


Figure 53

## ATTACHING AN OXYGEN BOTTLE TO AN OXYGEN BOTTLE HOLDER

### To attach an oxygen bottle:

1. Place an oxygen bottle in the holder.
2. Insert the lower strap through the buckle and affix the strap onto itself to secure the oxygen bottle to the holder.

**Note:** Inspect the straps and clips for wear between use and replace the strap if it is no longer holding the oxygen bottle.

---

### CAUTION

- To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 40 lb (18 kg).
- Do not use two head end oxygen bottle holders at the same time.

---

**Note:** The optional fowler oxygen bottle holder (p/n 6500-141-000) is not compatible with the optional retractable head section oxygen bottle holder (p/n 6085-046-000).

## USING THE RETRACTABLE HEAD SECTION OXYGEN BOTTLE HOLDER

To attach an oxygen bottle to the retractable head section oxygen bottle holder:

1. Center the oxygen bottle on the cradled surface of item (A) as shown in Figure 54.
2. Tighten both straps (B) around the oxygen bottle.
3. Secure the slack on the straps to the Velcro® on the straps.


Figure 54: Retractable Head Section Oxygen Bottle Holder

**Note:** Inspect the straps and clips for wear between use and replace the strap if it is no longer holding the oxygen bottle.

---

**⚠ WARNING**

If the cot is equipped with the optional retractable head section oxygen bottle holder, use caution while the oxygen bottle holder is installed to avoid pinching your fingers between the fowler bracket and the oxygen bottle.

---

**⚠ CAUTION**

- To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 40 lb (18 kg).
  - Do not use two head end oxygen bottle holders at the same time.
-

## ATTACHING THE PEDI-MATE® INFANT RESTRAINT SYSTEM

See the Pedi-Mate® users manual for the manufacturer's recommendations for the use, operation and care of the Pedi-Mate® Infant Restraint System.

### To secure the Pedi-Mate® to the cot:

1. Remove any restraints that are already attached to the cot.
2. Raise the cot backrest to the full upright position.
3. Position the Pedi-Mate® pad flat on the backrest with the black backrest straps out (see Figure 55).

Figure 55: Positioning the Pedi-Mate®

4. Wrap the straps around the backrest and insert the ends of the straps through the brackets. Securely fasten the buckle (see Figure 56).


Figure 56: Fastening the Pedi-Mate® Buckle

### WARNING

To avoid accidental release of the Pedi-Mate®, and possible injury to the infant, ensure that the buckle is located away from obstructions on the cot or accessories.

## ATTACHING THE PEDI-MATE® INFANT RESTRAINT SYSTEM (CONTINUED)

5. Pull firmly on the end of the adjustable backrest strap and tighten it securely.
6. Insert the mainframe straps between the cot frame and the mattress. To ensure that the release button is located toward the foot end of the cot, insert the buckle behind the litter cross brace and bring it up in front of the cross brace. Secure the buckle around the cross brace, leaving a little slack in the strap for final adjustment (see Figure 57).

Figure 57: Securing the Safety Restraints on a Cot

### **WARNING**

To avoid accidental release of the Pedi-Mate®, and possible injury to the infant, ensure that the buckle is located away from obstructions on the cot or accessories.

7. Verify that all of the straps are snug and fastened securely (see Figure 58).

Figure 58: Pedi-Mate® Strapped to a Cot

**Note:** These are general instructions for installation of the Pedi-Mate®. Safe and proper use of the Pedi-Mate® is solely at the discretion of the user. Stryker recommends that all users be trained on the proper use of the Pedi-Mate® before using it in an actual situation. Retain these instructions for future reference. Include them with the product in the event of transfer to new users.

Pedi-Mate® is a registered trademark of Ferno-Washington, Inc.

## INSTALLING THE BACKREST STORAGE POUCH

Install the optional backrest storage pouch using the Velcro® straps as shown in Figure 59. Insert each strap through a hole in the backrest skin and mount the pouch flat against the backrest.

### CAUTION

- Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.
- The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 20 lb (9 kg).


Figure 59: Backrest Storage Pouch

## INSTALLING THE HEAD END STORAGE FLAT

### WARNING

When the optional head end storage flat is being used, ensure that it does not interfere with the operation of the retractable head section, safety bar and safety hook. Injury to the patient or operator could result.


Figure 60: Head End Storage Flat

### To install the optional head end storage flat (see Figure 60):

1. Install the Velcro® straps (A) near the pneumatic cylinder and around the cross bar of the retractable head section.
2. Buckle the restraint straps (B) around the outer rails of the retractable head section.

### CAUTION

- Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.
- The weight of the equipment in the head end storage flat (if equipped) must not exceed 40 lb (18 kg).

## USING THE TRANSFER FLAT

When transferring larger patients, use of the Transfer Flat (6005-001-001) is recommended.

[Return To Table of Contents](#)

# Cleaning

---

The **Performance-PRO™** XT cot is designed to be power washable. The unit may show some signs of oxidation or discoloration from continuous washing, however, no degradation of the cot's performance characteristics or functionality will occur due to power washing as long as the proper procedures are followed.

Thoroughly clean the cot once a month. Clean Velcro® **AFTER EACH USE**. Saturate Velcro® with disinfectant and allow disinfectant to evaporate. Appropriate disinfectant for nylon Velcro® should be determined by the service.

## WASHING PROCEDURE

- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- The preferred method Stryker Medical recommends for power washing the cot is with the standard hospital surgical cart washer or hand held wand unit.

## WASHING LIMITATIONS

---

### **WARNING**

When cleaning, use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the cot.

---

### **CAUTION**

- DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT.
  - Maximum water temperature should not exceed 180°F/82°C.
  - Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 24 inches (61 cm) from the unit.
  - Allow cot to air dry.
  - Towel dry all casters and interface points.
  - Failure to comply with these instructions may invalidate any/all warranties.
-

In general, when used in those concentrations recommended by the manufacturer, either phenolic type or quaternary (**excluding Virex® TB**) type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may result.

Suggested cleaners for the cot surfaces and restraint straps:

- Quaternary Cleaners (active ingredient - ammonium chloride)
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)

**Note:** Do not immerse the restraint strap metal buckle components in water. Rinse with clean water. Allow to air dry.

Avoid over saturation and ensure that the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

---

 **WARNING**

SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products described above are used to clean Stryker EMS equipment, measures must be taken to ensure the cots are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the cots will leave a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.

**Note:** Failure to follow the above directions when using these types of cleaners may void this product's warranty (see [page 1-70](#)).

## REMOVAL OF IODINE COMPOUNDS

Use a solution of 1/2 Tablespoon Sodium Thiosulfate in a pint of warm water to clean the stained area. Clean as soon as possible after staining occurs. If stains are not immediately removed, allow solution to soak or stand on the surface. Rinse surfaces which have been exposed to the solution in clear water before returning unit to service.

---

 **WARNING**

Failure to properly clean or dispose of contaminated mattress or cot components will increase the risk of exposure to bloodborne pathogens and may cause injury to the patient or the operator.

---

# Preventative Maintenance

Preventative maintenance should be performed at a minimum of annually. A preventative maintenance program should be established for all Stryker Medical equipment. Preventative maintenance may need to be performed more frequently based on the usage level of the product.

Operation	Schedule	Procedure
Cleaning and Disinfecting	Each Use.	See <a href="#">page 1-22</a> .
Inspection	For 1-25 calls per month, inspect cot every 6 months  For 26-200 calls per month, inspect cot every 3 months  For 201+ calls per month, inspect cot monthly.	See checklist below.

**Note:** Use the Maintenance Record form on [page 1-27](#) to keep up to date maintenance records.

## CHECKLIST

- All fasteners secure (reference assembly drawings).
- All welds intact, not cracked or broken.
- No bent or broken tubing or sheet metal.
- No wear on the bumper and bumper housing.
- No debris in wheels.
- All wheels secure, rolling and swiveling properly.
- Optional wheel lock holds wheel securely when on and clears wheel when off.
- Siderails move and latch properly.
- Backrest operating properly.
- Optional accessories intact and operating properly.
- Height positioning latch functioning properly.
- Cot secure in each height position.
- Undercarriage folds properly.
- Retractable head section operating properly.
- Safety bar operating properly.
- Footrest operating properly.
- No rips or cracks in mattress cover.
- Body restraints intact and working properly.
- Lubricate base tubes (optional).

Product Serial Number:		

Completed by: \_\_\_\_\_ Date: \_\_\_\_\_

# Preventative Maintenance

## REGULAR INSPECTION AND ADJUSTMENTS

### Maintenance Intervals

Use this schedule as a general guide for maintenance. Factors such as weather, terrain, geographical location, and individual usage will alter the required maintenance schedule. If you are unsure as to how to perform these checks, contact your Stryker service technician. If you are in doubt as to what intervals to follow when maintaining your product, consult your Stryker service technician. Check each routine and replace damaged or worn parts if necessary.

Item	Routine	Every (whichever comes first)			
		One month	Three months	Six months	12 Months
Settings	Verify the cot and fastener fit and function properly	X			
	Verify the safety bar engages the vehicle safety hook properly	X			
Litter	Inspect the cot frame/litter	X			
	Verify all welds intact, not cracked or broken				X
	Verify no bent, broken or damaged components			X	
	Verify all fasteners secure (reference assembly drawings)		X		
	Verify warning labels present, legible (reference assembly drawings)				X
	Verify no damage or tears on cot grips			X	
	Verify the siderails operate and latch properly			X	
	Verify the backrest cylinder operates properly		X		
	Adjust pneumatic cylinder for full range of motion, if required		X		
	Verify the footrest operates properly			X	
Mattress	Inspect the fowler bumper and screws				X
	Inspect the dead stop (6085-001-094) for wear. Replace if necessary		X		
Restraints	Verify no cracks or tears on cot mattress			X	
	Inspect patient restraints for proper function and no excessive wear (bent or broken receiver or latch plate, torn or frayed webbing, etc.)	X			
Base	Inspect the cot frame/base	X			
	Verify all welds intact, not cracked or broken				X
	Verify no bent, broken, or damaged components			X	
	Verify all fasteners secure (reference assembly drawings)		X		
	Verify that the cot retaining post is secure. If not secure, replace the screw. See "Cot Retaining Post Screw Replacement" on page 1-63.			X	
Verify no excessive damage to X-frame guards			X		

[Return To Table of Contents](#)

# Preventative Maintenance

Item	Routine	Every (whichever comes first)			
		One month	Three months	Six months	12 Months
Wheels	Verify wheels are free of debris			X	
	Verify rubber is in good condition				X
	Verify all wheels secure, rolling and swiveling properly	X			
	Check and adjust optional wheel locks as necessary				X
X-Frame	Verify smooth operation of X-frame		X		
Head Section	Verify all fasteners secure (reference assembly drawings)		X		
	Verify no bent, broken, or damaged components			X	
	Verify the head section extends and locks properly		X		
	Verify the grip bar has no excessive damage or tears			X	
	Verify load wheels are secure and roll properly			X	
	Verify the safety bar operates properly. Pull toward the head section to ensure that it swings and rotates freely and pulls back to home position.	X			
Accessories	Inspect the straps and clips on the retractable head section oxygen bottle holder (optional) for wear				X
	Verify the I.V. pole (optional) operates properly		X		
	Verify the head extension & pillow (optional) operates properly		X		
	Verify the Pedi-Mate® restraint package (optional) operates properly		X		
	Verify the restraint extender (optional) operates properly		X		
	Verify the oxygen bottle holder (optional) operates properly		X		
	Verify the transfer flat (optional) operates properly		X		
	Verify the equipment hook (optional) operates properly		X		
	Verify the pocketed back rest pouch (optional) operates properly		X		
	Verify the base storage net (optional) operates properly		X		
	Verify the kickstand (optional) retracts fully to the transport position		X		
	Verify that the kickstand (optional) bolts are tightened properly		X		
	Lubricate the kickstand spring and internal spring housing (optional) using Tri-Flow® lubrication			X	

# Quick Reference Replacement Parts List

The parts and accessories listed on these pages are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call Stryker Customer Service USA: 1-800-327-0770 (Option 2) for availability and pricing.

Part Name	Part Number
Dead Stop	6085-001-094
Gas Spring	1010-031-077
Head Extension - Pillow Only	6100-045-000
Head Section	6500-002-020
I.V. Pole, Two-Stage, Right	6500-210-000
I.V. Pole, Two-Stage, Left	6500-211-000
I.V. Pole, Two-Stage, Dual	6500-212-000
I.V. Pole, Three-Stage, Right	6500-215-000
I.V. Pole, Three-Stage, Left	6500-216-000
I.V. Pole, Three-Stage, Dual	6500-217-000
Kit, Retractable Head Section Oxygen Bottle Holder	6085-700-003
Label, "Lift Here"	6080-090-108
Label, Side Release	6085-001-159
Mattress, Bolster	6090-041-010
Mattress, Flat	6090-042-010
Mattress, Gatch Compatible Option	6550-001-084
Restraint Belt Extension	6082-160-050
Restraint Package, Domestic	6082-260-010
Safety Hook, J	6092-936-018
Safety Hook, Long	6060-936-018
Safety Hook, Short	6060-936-017
Sensor Housing Cover	6500-001-199
Label, Sensor Housing	6085-001-156
Siderail Assembly	6082-026-010
Touch-Up Paint (Yellow)	6060-199-010
Touch-Up Paint (Black)	6060-199-011
Velcro Adhesive Loop Pile, Litter	6060-032-046
Wheel Lock, Single Adjustable	6082-501-010
Wheel Lock, Dual Adjustable	6082-502-010

## BACKREST ADJUSTMENT

### Tools Required:

- 1/2" Combination Wrench
- 5/32" Hex Wrench
- 3/32" Hex Wrench
- Small Slotted Screwdriver

### Procedure:

1. For easier access, move the backrest to 73 degrees.

**Note:** Before continuing with the backrest adjustment procedure, be sure that the cylinder (A) is completely threaded into the yoke (B) so no threads are showing on the shaft of the cylinder. If threads are showing, use a 3/32" hex wrench to remove the set screw (C) in the center of the yoke. Using a small slotted screwdriver, remove the E-clip and pin (D & E) that holds the bottom of the pneumatic cylinder. Thread the cylinder shaft (A) completely into the yoke (B). Replace the E-clip and pin (D & E) and replace the set screw (C) using Loctite®. (Figure 61 and Figure 62)


Figure 61

2. Using a 1/2" combination wrench, loosen the hex nut (F) on the backrest pivot (J) while holding the set screw (H) fixed in the pivot (Figure 62).
3. Using a 5/32" hex wrench, turn the set screw (H) until there is no play between the backrest release handle (K) and the pneumatic cylinder release button (Figure 62).

**Note:** Make sure that the backrest travels from flat to at least 73 degrees. If it does not, turn the set screw clockwise 1/2 turn. Repeat until at least 73 degrees of travel is achieved.

4. Lower the backrest to a 5–10 degree angle and release the handle. Apply approximately 50 lb of downward force to the end of the backrest. If the backrest drifts down, turn the set screw counterclockwise. Repeat until the backrest does not drift downward.
5. Using the 1/2" combination wrench, tighten the hex nut (F) while holding the set screw fixed in the pivot (Figure 62).
6. Verify proper operation of the unit before returning it to service.


Figure 62

## WHEEL LOCKING FORCE ADJUSTMENT

### Tools Required:

- 5/32" Hex Wrench
- 7/16" Combination Wrench or Socket

### Procedure:

1. Using the 5/32" hex wrench and 7/16" combination wrench or socket, remove the socket screw from the center of the lock pedal. The wheel lock is initially assembled with the pedal set at the minimum locking force. The marker on the pedal (A) is aligned with the marker on the octagonal sleeve (B) (Figure 63).
2. Remove the sleeve (B). Rotate the sleeve counterclockwise to increase the pedal locking force and clockwise to decrease the locking force. Insert the sleeve into the pedal (Figure 63).
3. Using the 5/32" hex wrench and 7/16" combination wrench or socket, reinstall the socket screw.
4. Test the pedal locking force and verify that the pedal holds properly before returning it to service.


Figure 63: Wheel Locking Force Adjustment

## COT RETAINING POST ADJUSTMENT

### Tools Required:

- 3/16" Hex Wrench

### CAUTION

The cot retaining post comes preconfigured for an X-frame cot, if the fastener has been configured for an H-frame style cot, the cot retaining post must be adjusted to accommodate the fastener.

### Procedure:

1. Using a 3/16" hex wrench, remove the two socket head cap screws (A) that hold the pin brackets (B) to the base frame (C) (Figure 64).
2. Turn the bottom pin bracket 180°.
3. Using a 3/16" hex wrench, reinstall the two socket head cap screws that were removed in step 1.
4. Verify proper operation of the unit before returning it to service.

**Note:** If the arrow (D) on the bottom bracket of the retaining post (E) points toward the head end of the cot, the retaining post is set for an X-frame style cot. If the arrow points toward the foot end of the cot the post is set for an H-frame style cot (Figure 64).


Figure 64: Cot Retaining Post

## COT RETAINING POST REPLACEMENT

### Tools Required:

- T30 Torx Driver
- 5/32" Hex Wrench
- Torque Wrench (in-lb)

### Procedure:

1. Raise the cot to the full upright position.
2. Turn the cot onto the patient left side.

**Note:** Locate the arrow that is located on the bottom bracket. The replacement retaining post bracket will need to be assembled in the same orientation.


Figure 65

3. Using a T30 Torx driver, remove the two socket head cap screws (A) that secure the current cot retaining post to the base tube (B) (Figure 65). Discard the screws and cot retaining post.
4. Insert the button head cap screw (C) through the retaining post cap (D) and post tube (E), and then into the top pin bracket (F) (Figure 65).
5. Using a 5/32" hex wrench, tighten the button head cap screw (C) completely to secure the retaining post cap (D) and post tube (E) to the top pin bracket (F) (Figure 65). Using a torque wrench, torque the screw to 100-140 in-lb.
6. Assemble the cot retaining post across the base tube. Align the holes of the retaining post halves and insert the two socket head cap screws into the threaded holes of the retaining post top.
7. Using a T30 Torx driver, tighten the two socket head cap screws completely.
8. Verify proper operation of the unit before returning it to service.

**Note:** Adjustment of the rail clamp assembly may be required in order to compensate for any variation in cot retaining post position depending on the ambulance cot manufacturer and model number.

## COT RETAINING POST SCREW REPLACEMENT

### Tools Required:

- T25 Torx Driver
- 5/32" Hex Wrench
- Torque Wrench (in-lb)

### Procedure:

1. Using a T25 Torx driver or 5/32" hex wrench, remove the button head cap screw that secures the retaining post cap and post tube to the top pin bracket. Discard the screw.
2. Using a 5/32" hex wrench, install and tighten the button head cap screw (p/n 0004-503-000) completely to secure the retaining post cap and tube to the top portion of the retaining post assembly. Using a torque wrench, torque the screw to 100-140 in-lb.

**Note:** If you cannot torque the screw to 100-140 in-lb, then you must replace the entire cot retaining post. See "Cot Retaining Post Replacement".

3. Verify proper operation of the unit before returning it to service.

## HEADSECTION REPLACEMENT

### Tools Required:

- 7/16" Combination Wrench
- 3/16" Hex Wrench

### Procedure:

1. Raise the cot and the backrest to the full upright position.
2. Using a 7/16" combination wrench and a 3/16" hex wrench, remove the two screws (A) that secure the cap bearings to the base litter interface bracket (one on each side) (Figure 66).
3. Squeeze the head release handles and slowly remove the head section assembly.
4. Reverse steps to reinstall.
5. Verify proper operation of the unit before returning it to service.


Figure 66

## BACKREST GAS CYLINDER REPLACEMENT

### Tools Required:

- 3/32" Hex Wrench
- Slotted Screwdriver

### Procedure:

1. Raise the cot and the backrest to the full upright position.
2. Using a 3/32" hex wrench, loosen the set screw (A) that holds the gas shaft to the yoke (Figure 67).
3. Using a slotted screwdriver, remove the e-clip from the clevis pin that secures the bottom of the gas cylinder.
4. Unscrew the gas cylinder shaft from the yoke.
5. Reverse the above procedures to install the new gas cylinder. See "Backrest Adjustment" on [page 1-60](#).
6. Verify proper operation of the unit before returning it to service.


Figure 67

## INNER, INNER TUBE REPLACEMENT

### Tools Required:

- (2) Saw Horse
- 3/8" Combination Wrench
- 7/16" Combination Wrench
- T25 Torx Driver
- (2) 9/16" Combination Wrench
- Dead Blow Hammer

### Procedure:

1. Using a T25 Torx driver, remove the four button head cap screws (A) from the base stiffener (B) (Figure 68).
2. Lower the cot to the full down position.
3. Using two saw horses, flip the cot upside down on the saw horses to support the cot.
4. Using a 3/8" and 7/16" combination wrench, remove all four caster mount bolts (C) that secure the outer lift tube assemblies to the foot base tubes (Figure 69).
5. Using two 9/16" combination wrenches and a dead blow hammer, remove the base connecting rod (D) and centerlock hex jam nut (E).

**Note:** Save the bearings (F) that fall out. (Figure 71).

6. Remove the outer base tube weldment (G) on the side that is damaged. Leave the opposite side on to support the X-frame (Figure 70).
7. Slide the foot base tubes (H) through the X-frame legs to loosen the X-frame (Figure 70).

**Note:** The opposite side will still have the outer base tube and foot base tubes attached to the X-frame.

8. Remove the inner tube and save the bushings (J) to reuse on the new tube.
9. Reverse the above procedures to install the new inner tube.
10. Verify proper operation of the unit before returning it to service.


Figure 68


Figure 69


Figure 70


Figure 71

[Return To Table of Contents](#)

## OUTER, INNER TUBE REPLACEMENT

### Tools Required:

- (2) Saw Horse
- 3/8" Combination Wrench
- 7/16" Combination Wrench
- T25 Torx Driver
- (2) 9/16" Combination Wrench
- Dead Blow Hammer
- T27 Torx Driver
- 

### Procedure:

1. Complete steps 1-7 from "Inner, Inner Tube Replacement" on [page 1-65](#).
2. Using a T27 Torx driver, remove the truss head screw (A) that secures the X-frame guard (B) to the inner tube (see Figure 72).
3. Using a T25 Torx driver, remove the two button head cap screws (C) that secure the tube bearing.
4. Remove the inner tube and save the tube bushings to reuse on the new tube.
5. Reverse the above procedures to install the new inner tube.

**Note:** Make sure that the X-frame guard screw hole on the new inner tube is facing the top of the cot when reassembling.

6. Verify proper operation of the unit before returning it to service.


Figure 72


Figure 73

## OUTER, OUTER TUBE REPLACEMENT

### Tools Required:

- (2) Saw Horse
- 3/8" Combination Wrench
- 7/16" Combination Wrench
- T25 Torx Driver
- (2) 9/16" Combination Wrench
- Dead Blow Hammer
- 3/16" Allen Wrench
- 1/2" Combination Wrench
- Needle Nose Pliers

### Procedure:

1. Complete steps 2-7 from "Inner, Inner Tube Replacement" on [page 1-65](#).
2. Using a T25 Torx driver, remove the two button head cap screws that secure the tube bearing (see "Figure 73" on [page 1-66](#)).
3. Remove the inner tube.
4. Remove the outer tube X-frame guard and set aside to reuse on the new outer tube (see "Figure 72" on page 1-66).
5. Using a 3/16" Allen wrench and 1/2" combination wrench, remove the bolt that holds the timing link to the outer tube and remove the timing link.
6. Depending on where your side release handle is located, you may need to remove the height adjustment rack springs. Using needle nose pliers, remove the two return springs for the height adjustment rack.
7. Using a 3/16" Allen wrench, remove the four bolts that secure the slider housing and set aside.
8. Remove the outer tube and save the tube bushings to reuse on the new tube.
9. Reverse the above procedures to install the new outer tube.

**Note:** Make sure that the X-frame guard screw hole on the new inner tube is facing the top of the cot when reassembling.

10. Verify proper operation of the unit before returning it to service.

## INNER, OUTER TUBE REPLACEMENT

### Tools Required:

- (2) Saw Horse
- 3/8" Combination Wrench
- 7/16" Combination Wrench
- T25 Torx Driver
- (2) 9/16" Combination Wrench
- Dead Blow Hammer
- 3/16" Allen Wrench

### Procedure:

1. Complete steps 1-7 from "Inner, Inner Tube Replacement" on [page 1-65](#).
2. Remove the inner tube.
3. Using a T25 Torx driver, remove the screw that secures the base dead stop from the outer tube.
4. Using a 7/16" combination wrench and a 3/16" Allen wrench, remove the two screws that secure the cap bearings to the base litter interface bracket (one on each side).
5. Squeeze the head release handles and slowly remove the head section assembly.
6. Using a 7/16" combination wrench and 3/16" Allen wrench, remove the bolt that secures the litter interface bracket.
7. Using a 3/16" Allen wrench, remove the two screws that hold the litter interface bracket together.
8. For the trend option, using a T25 driver, remove the two screws that secure the trend support bracket.
9. For the gatch option, using a T25 driver, remove the four screws that secure the gatch support brackets.
10. Pull outward on the outer rail until the litter interface bracket is off of the litter crosstube.
11. Remove the outer tube from the litter crosstube and save the tube bushings to reuse on the new tube.
12. Reverse the above procedures to install the new outer tube.

**Note:** Make sure that the X-frame guard screw hole on the new inner tube is facing the top of the cot when reassembling.

13. Verify proper operation of the unit before returning it to service.

## SIDERAIL ASSEMBLY REPLACEMENT

### Tools Required:

- T25 Torx Driver

### Procedure:

1. Raise the cot to the full upright position.
2. Raise the siderail to the up and locked position.
3. Using a T25 driver, remove the three spindle screws that secure the siderail assembly.
4. Remove the siderail.
5. Reverse the above procedures to install the new siderail assembly.
6. Verify proper operation of the unit before returning it to service.

# Warranty

---

**Stryker EMS**, a division of the Stryker Corporation, offers two distinct warranty options in the United States:

**One (1) year parts and labor.** Under this option, Stryker EMS warrants to the original purchaser that its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective.

**Two (2) year parts.** Under this option, Stryker EMS warrants to the original purchaser that non-expendable components of its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of two (2) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts for, or replacing, at its option, any product which is, in the sole discretion of Stryker, found to be defective. Expendable components, i.e. mattresses, restraints, I.V. poles, storage nets, storage pouches, oxygen straps, and other soft goods, have a one (1) year limited warranty with this option.

Under either warranty option, the Stryker **Performance-PRO™ XT** is designed for a 7 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual. Stryker warrants to the original purchaser that the welds on the **Performance-PRO™ XT** will be free from structural defects for the expected 7 year life of the product as long as the original purchaser owns the product. Original purchasers will also obtain a three (3) year limited parts warranty for the X-frame components.

Upon Stryker's request, purchaser shall return to Stryker's factory any product or part (freight prepaid by Stryker) for which an original purchaser makes a warranty claim.

Any improper use or alteration or repair by unauthorized service providers in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

This statement constitutes Stryker EMS's entire warranty with respect to the aforesaid equipment. STRYKER MAKES NO OTHER WARRANTY OR REPRESENTATION EITHER EXPRESSED OR IMPLIED, EXCEPT AS SET FORTH HEREIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL STRYKER BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

## STRYKER EMS RETURN POLICY

Cots, Stair Chairs, Evacuation Chairs, Cot Fasteners and Aftermarket Accessories may be returned up to 180 days of receipt if they meet the following guidelines:

### Prior to 30 Days

- 30 day money back guarantee in effect
- Stryker EMS is responsible for all charges
- Returns will not be approved on modified items

### Prior to 90 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 10% restocking fee

### Prior to 180 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 25% restocking fee

## RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. **SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.**

## DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Claims for any short shipment must be made within thirty (30) days of invoice.

## INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

## PATENT INFORMATION

The Stryker **Performance-PRO™ XT** cot is covered by one or more of the following patents:

United States	5,575,026	6,276,010	6,648,343	6,908,133	6,796,757
	5,537,700	6,125,485	6,735,794	7,100,224	7,398,571
	D527,103				

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