M-1® Roll-In System

Operations Manual

REF 6100

EN	FR	PL
CS	HR	PT
DA	HU	RO
DE	IT	RU
EL	LT	SK
ES	LV	SL
ET	NL	SV
FI	NO	TR

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	CE mark
	UK Conformity Assessment mark
	Importer
	Unique device identifier
	Authorized representative in the European Community
	Authorized representative in Switzerland
	European medical device
	Catalogue number
	Lot (batch) code
	Serial number
	For US Patents see www.stryker.com/patents
	Manufacturer
	Date of manufacture

	Safe working load
	Type BF applied part

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Indications for use	4
Clinical benefits	5
Contraindications	5
Expected service life	5
Disposal/recycle	5
Specifications	5
Schnitzler Pediatric Restraint System compatibility	6
Product illustration	7
Contact information	8
Serial number location	9
Setup	10
Installation	11
Installing the cot fastener	11
Installing the safety hook option	12
Installing the catch plate	14
Operation	16
Operating guidelines	16
Proper lifting techniques	16
Transferring the patient to the cot	16
Rolling the cot with a patient	17
Changing the cot height	17
Loading the cot into a vehicle	18
Unloading the cot from a vehicle	18
Unloading the cot with the safety hook option	19
Positioning operators and helpers	19
Raising or lowering the backrest	20
Raising or lowering the siderails	20
Raising or lowering the footrest	20
Raising or lowering the knee Gatch option	20
Applying or releasing a wheel lock	21
Applying or releasing the caster Steer-Lock	21
Connecting the litter to the base	22
Removing the litter from the base	23
Changing cot height with the optional base head end retractable lifting handles	23
Extending or retracting the retractable lifting handles	24
Securing the patient with the restraint straps	25
Adjusting restraint straps	26
Adding a restraint strap extension	27
Securing the patient with the Pedi-Mate® infant restraint system	27
Attaching the defibrillator platform	28
Attaching the head extension with pillow	30
Positioning the three-stage IV pole option	30
Attaching an oxygen bottle to the oxygen bottle holder	31
Accessories and parts	32
Cleaning	33
Disinfecting	34
Preventive maintenance	35
Regular inspection and adjustments	35
Every month or two hours	35
Every three months or six hours	36
Every six months or 12 hours	36
Every 12 months or 24 hours	36

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always replace the cot fastener if it has been involved in an accident to avoid the risk of injury due to product damage.
- Always make sure that all cots meet the installation specifications for the Stryker cot fastener system.
- Do not allow untrained assistants to assist in the operation of the product.
- Always follow proper hand placement on hand grips. Keep all hands clear of the red safety bar pivots when you load or unload the cot or change cot height position.
- Do not ride on the base of the cot.
- Do not transport the cot sideways to avoid the risk of tipping. Always transport the cot in a lowered position, head end or foot end first to minimize the risk of tipping.
- Always keep hands, fingers, and feet away from moving parts. Use caution when you place your hands and feet near the base tubes while you raise or lower the cot.
- Always use all restraint straps to secure the patient on the product. An unrestrained patient may fall from the product and be injured.
- Do not leave a patient unattended. Hold the product while a patient is on the product.
- Do not apply a wheel lock when a patient is on the product or when you move the product to avoid the risk of tipping.
- Do not use siderails as a patient restraint device.
- Always transport the cot at a lower height to reduce the risk of a cot tip. If possible, obtain additional assistance or take an alternate route.
- Always avoid high obstacles, such as curbing, steps, or rough terrain to avoid the risk of the product tipping over.
- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
- Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.
- Always lock the base frame before you remove the loading wheels from the vehicle patient compartment floor or loading tray system. An unlocked base frame will not support the cot and injury to the patient or operator could result.
- Always practice changing height positions and loading the cot until you fully understand the operation of the product.
- Always support the load of the patient, cot, and accessories after the weight is off of the ground.
- Always follow proper hand placement on hand grips. Keep all hands clear of the red safety bar pivots when you load, unload, or change the cot height position.
- Do not install or apply a wheel lock on a product with worn wheels that are less than 6 in. diameter.
- Do not leave a patient or occupant unattended. Hold the product while a patient or occupant is on the product.
- Do not attach restraint straps to the base tubes or cross tubes.

- Always locate the buckle away from obstructions or accessories on the cot to avoid the risk of accidental release of **Pedi-Mate®** infant restraint system and injury to the infant.
 - Do not use the oxygen bottle holder to hold an oxygen bottle holder when the transport vehicle is in motion. Always place the oxygen bottle holder in an appropriate storage location when the transport vehicle is in motion.
 - Always inspect the straps and clips for wear between use. Replace the strap if it no longer holds the oxygen bottle.
 - Always use any appropriate personal protective equipment while power washing to avoid inhaling contagion. Power washing equipment may aerate contamination.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not use the product with incompatible restraints.
 - Always install the cot fastener by qualified personnel only to be in conformance to this standard. 1/4" steel backing plates (or equivalent) are required beneath the mounting surface (ambulance floor or tray system) at all mounting points.
 - Always clear any obstacles that may interfere and cause injury to the operator or patient before operating the product.
 - Do not apply the **Steer-Lock** and attempt to push the product sideways.
 - Do not entangle the restraint straps in the base frame when you raise or lower the cot.
 - Always secure the defibrillator platform to the product when you use the defibrillator platform.
 - Always use and adjust the straps that are provided with the defibrillator platform to secure the defibrillator.
 - Always change the attachment location or adjust the straps for your specific defibrillator size or shape.
 - Do not load the defibrillator platform above the safe working load of 30 lb (13.6 kg).
 - Do not load the IV pole above the safe working load of 25 lb (11.3 kg).
 - Do not load the oxygen bottle holder above the safe working load of 15 lb (6.8 kg).
 - Do not use two oxygen bottle holders at the same time.
 - Do not steam clean or ultrasonically clean the product.
 - Do not exceed 180° F (82° C) as the maximum water temperature.
 - Do not exceed 1500 psi (130.5 bar) as the maximum water pressure. If you are using a hand held wand to wash the product, keep the pressure nozzle at a minimum of 24 in. (61 cm) from the product.
 - Always allow to air dry.
 - Do not clean, disinfect, service, or perform maintenance while the product is in use.
 - Always use authorized parts to avoid the risk of product damage.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 6100 **M-1**® Roll-In System is a manual ambulance cot that consists of a platform on a wheeled frame designed to support and transport a maximum weight of 500 lb (227 kg) in pre-hospital and hospital environments. The device is collapsible for use in emergency vehicles. The removable litter and patented **Steer-Lock System™** provide enhanced maneuverability. To assist with patient transport, the product features:

- Lift handles at four corners
- Patient restraint straps
- Adjustable pneumatic backrest
- Optional accessories

For maximum patient comfort, you can select from one of three different litter positions, including:

- Shock
- Flat leg
- Optional knee Gatch

Indications for use

The Stryker Model 6100 **M-1** Roll-In System is a non-powered wheeled stretcher, which is intended to support and transport the entire body of a traumatized, ambulatory, or non-ambulatory human patient (includes infants and adults). In emergency or transport vehicles, the product:

- Supports patients in a supine (horizontal) or sitting position
- Facilitates the transportation of associated medical equipment (i.e. oxygen bottles or IV bags)

This ambulance cot is intended to be used in pre-hospital and hospital environments, in emergency and non-emergency applications. It is rated to a maximum capacity of 500 lb (227 kg), which is the sum of the patient, mattress, and accessory weight. Intended operators of the product include:

- Trained professionals
- Emergency medical service
- Medical care center personnel
- Medical first responders

Ambulance cots are for transportation purposes. They are not intended for extended stay or for use as hospital beds. They are also not intended for use in devices that change air pressure, such as hyperbaric chambers.

Clinical benefits

Cot: patient transport

Fastener: support cot for transport

Cot and fastener system: support and transport patients

Contraindications

None known.

Expected service life

The Stryker Model 6100 **M-1** Roll-In System has a five year expected service life under normal use conditions and with appropriate periodic maintenance.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

For the Stryker Model 6100 **M-1** Roll-In System base unit (6100-003-000), knee Gatch litter (6100-031-000), and flat foot section litter (6100-033-000).

	Note - Safe working load indicates the sum of the patient, mattress, and accessory weight.	500 lb	227 kg
Backrest articulation/shock position (6100-031-000 and 6100-033-000)		0° - 75°	
Length/width (6100-003-000)		77 in. / 21 in.	195.6 cm / 53.3 cm
Length/width (6100-031-000)		74.5 in. / 22.1 in.	189.2 cm / 56.1 cm
Length/width (6100-033-000)		74.8 in. / 22.1 in.	189.9 cm / 56.1 cm
Height ¹ (6100-003-000)	Position 1	14.4 in.	36.6 cm
	Position 2	24.4 in.	62.0 cm
	Position 3	30.0 in.	76.2 cm
	Position 4	34.6 in.	87.9 cm
	Position 5	37.6 in.	95.5 cm
	Position 6	39.3 in.	99.8 cm
	Position 7	13.0 in.	33.0 cm

	Position 8	Not applicable	Not applicable
Height ¹ (6100-031-000 and 6100-033-000)		7.3 in.	18.5 cm
Weight ² (6100-003-000)		77.7 lb	35.2 kg
Weight ² (6100-031-000)		44.6 lb	20.2 kg
Weight ² (6100-033-000)		40.2 lb	19.0 kg
Maximum weight capacity (6100-031-000 and 6100-033-000)		500 lb	227 kg
Wheel diameter/ width (6100-031-000 and 6100-033-000)		4.0 in. / 0.8 in.	10.2 cm / 2.0 cm
Handle extension (6100-031-000 and 6100-033-000)		7.75 in.	19.7 cm
Caster diameter/width (6100-003-000)		6 in. / 2 in.	15.2 cm / 5.1 cm
Shock position (6100-031-000)		+17°	
Knee Gatch position (6100-031-000)		30°	
Minimum operators required for loading/ unloading		1	
Recommended fastener systems		Model 6381 floor mounted - center mount, straight guide Model 6373 tray mounted - center mount, din guide Model 6376 floor mounted - center mount, straight guide - no side restraints	
Recommended floor/tray height		27 in. - 31 in.	68.6 cm - 78.7 cm

¹Height is measured from bottom of mattress at seat section to ground level.

²Standard base unit is weighed without optional accessories.

Stryker reserves the right to change specifications without notice.

The yellow and black color scheme is a proprietary trademark of Stryker Corporation.

Environmental conditions	Operation	Storage and transportation
Temperature		
Relative humidity		
Atmospheric pressure		


Schnitzler Pediatric Restraint System compatibility


CAUTION - Do not use the product with incompatible restraints.

XPS mattress (6500-003-130 or 6506-003-130)	Schnitzler Pediatric Restraint System (0058-384-000) with serial number less than 720 XPSNPR-2123	Schnitzler Pediatric Restraint System (0058-384-000) with serial number equal to or greater than 720 XPSNPR-2123
Mattress with lot code less than 18001001	Compatible	Compatible
Mattress with lot code equal to or greater than 18001001	Not compatible	Compatible

The Schnitzler Pediatric Restraint System (0058-385-000) is compatible for use with the standard mattress on Model 6100 M-1 Roll-In System, Model 6506 Power-PRO XT Cot, and Model 6550 Power-PRO TL Cot.

Product illustration

A	Base frame
B	Caster Steer-Lock knob
C	Foot end
D	Foot end release handle (green)
E	Foot end release handle (red)
F	Hand grip
G	Head end
H	Head end release handle (red)
I	Litter release button (blue)
J	Litter retention bracket
K	Loading wheel
L	Release button (green)
M	Retractable base lifting handle option

N	Transport wheel
P	Wheel lock
R	Retractable litter lifting handle
T	Knee Gatch release lever
U	Knee Gatch lifting handle
V	Roller wheel
W	DIN retention pin option
X	Litter frame
Y	Siderail release handle
Z	Perimeter bumper
AA	Backrest
AB	Backrest adjustment release handle

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
 3800 E. Centre Avenue
 Portage, MI 49002
 USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location


Figure 1 – Serial number location

Setup

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Remove all shipping and packaging materials from the product prior to use.

Unpack the cartons and check all items for proper operation. Make sure that the product operates properly before it is put into service.

Make sure that the patient compartment of the vehicle where you will be using the product has a:

- Smooth rear edge for product loading
- Level floor that is large enough for the folded product and cot fastener system
- Stryker Model 6373, 6376, or 6381 cot fastener system (not included)
- Floor or tray system loading height between 27.0 in. (68.6 cm) and 31.0 in. (78.7 cm)

Note - Always remove loose items or debris from the patient compartment floor. Loose items or debris may interfere with the cot rolling smoothly or the operation of the cot fastener.

When necessary, modify the vehicle to fit the cot. Do not modify the cot.

This manual should be considered a permanent part of the product and should remain with the product even if the product is subsequently sold.

Stryker continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Installation

Installing the cot fastener

The Stryker Model 6373, 6376, and 6381 cot fastener systems are compatible only with cots that conform to the installation specifications. Cots that meet these specifications include:

- Stryker Model 6100 **M-1** Roll-In System
- Ferno-Washington Model X-2¹

WARNING

- Always replace the cot fastener if it has been involved in an accident to avoid the risk of injury due to product damage.
- Always make sure that all cots meet the installation specifications for the Stryker cot fastener system.

CAUTION - Always install the cot fastener by qualified personnel only to be in conformance to this standard. 1/4" steel backing plates (or equivalent) are required beneath the mounting surface (ambulance floor or tray system) at all mounting points.

For more information about the Stryker cot fastener systems, see the **M-1** Cot Fastener Operations and Maintenance Manuals.

¹ 2000 model year or earlier. Stryker is not responsible for changes in specifications to other manufacturer's cots.


Figure 2 – Model 6376: View 2A


Figure 3 – Model 6376: View 2B


Figure 4 – Model 6373: View 2C


Figure 5 – Model 6373: View 2D

Installing the safety hook option

Safety hook kit (610000010001):

- (1) Hook arm assembly (610000010016)
- (1) Guide bracket (610000010008)
- (1) Roller assembly (610000010011)

Tools required (metric or imperial):

- Blue **Loctite**® 242
- T27H star key
- M4 Allen key or 5/32" hex wrench
- M6 Allen key or 1/4" hex wrench
- M11 spanner or 7/16" wrench
- M13 spanner or 1/2" wrench
- M16 spanner or 5/8" wrench
- Torque wrench

Procedure:

1. Apply the brakes. See *Applying or releasing a wheel lock* (page 21).
2. Lower the base frame to any mid height position. See *Changing the cot height* (page 17).
3. Remove the litter from the cot.
4. Lay the cot base frame on its side for easy access to components.

Note - Place a towel or sheet under the cot to protect from scratches or damage.

5. To install the supplied hook arm assembly:

- a. Using a T27H star key and M16 spanner, remove the button head cap screw (A) and flange nut (B) (Figure 6). Save the screw and nut.
- b. Apply Blue **Loctite**® 242 to the button head cap screw (A) or flange nut (B).
- c. Install the supplied hook arm assembly (C) (610000010016) with the screw and nut removed in step 5a (Figure 6).


Figure 6 – Hook arm

6. To install the supplied guide bracket:

- a. Using an M4 Allen key and M11 spanner, remove the two Fiberlock hex nuts (D) from the two hex socket head cap screws (E) at the foot end of the cot (Figure 7). Save the nuts and screws.
- b. Place the supplied guide bracket (F) (610000010008) over the two hex socket head cap screws (saved in step 3a) (Figure 7). Using a torque wrench, secure with the Fiberlock hex nuts (D) (removed in step 5a).


Figure 7 – Guide bracket

7. Unscrew the nuts from the supplied roller assembly (610000010011).
8. Using an M6 Allen key and an M13 spanner, install the supplied roller assembly (H) and nuts (G) on the slide tube support (Figure 8).


Figure 8 – Roller assembly

9. Return the cot to its upright position.
10. Replace the litter (removed in step 3).
11. Verify proper operation.
 - a. Use one finger to raise and lower the hook arm assembly. Confirm that the hook arm assembly moves freely on the cot.
 - b. Swing the head end legs back and forth to make sure that the safety hook and guide bracket do not cause interference.
12. Complete *Installing the catch plate* (page 14).

Installing the catch plate

The catch plate is intended for use with the **M-1** safety hook option (610000010001).

Catch plate kit (610000010022)

- (1) Catch plate (610000010015)

Additional part used (not included):

- (2) M6 Allen bolt

Tools required (metric or imperial):

- M5 Allen key or 3/16" hex wrench
- M6 Allen key or 1/4" hex wrench

Stryker recommends that you follow the specifications listed below for the catch plate installation:

- Angle (A) of the loading tray does not exceed 16° (Figure 9)
- Height (B) from the top of the catch plate to the ground is between 28.0" and 29.4" (Figure 9)
- Horizontal clearance (C) between the safety hook and catch plate is sufficient for unloading (Figure 10)


Figure 9 – Catch plate and loading tray specifications


Figure 10 – Safety hook and catch plate clearance

Procedure:

1. Using an M6 Allen key, install the supplied catch plate (610000010015) (D) and two M6 Allen bolts (E) (not included) directly behind the cot fastener latch (Figure 11).


Figure 11 – Catch plate

Note - Before you drill the two holes for the M6 Allen bolts, make sure that there is free space under the table to avoid product or vehicle damage.

2. Test unloading the cot from a vehicle to confirm that the safety hook connects with the catch plate. Using an M5 Allen key, adjust the hook tip (F), as needed, to make sure that the safety hook and catch plate connect while unloading (Figure 12).


Figure 12 – Hook tip adjustment

Note - Push the cot legs against the loading tray when you test unloading.

Operation

Operating guidelines

WARNING

- Do not allow untrained assistants to assist in the operation of the product.
 - Always follow proper hand placement on hand grips. Keep all hands clear of the red safety bar pivots when you load or unload the cot or change cot height position.
 - Do not ride on the base of the cot.
 - Do not transport the cot sideways to avoid the risk of tipping. Always transport the cot in a lowered position, head end or foot end first to minimize the risk of tipping.
 - Always keep hands, fingers, and feet away from moving parts. Use caution when you place your hands and feet near the base tubes while you raise or lower the cot.
-

CAUTION - Always clear any obstacles that may interfere and cause injury to the operator or patient before operating the product.

- Operate the product only as described in this manual.
- Read all labels and instructions on the product before use.
- Practice changing height positions and loading the cot until you fully understand the operation of the product.
- Always load or unload an occupied cot with a minimum of two trained operators. Two operators must be present when a cot is occupied. Stryker recommends that both operators are at the foot end to reduce the load on each operator. One or two operators can lift from the foot end of the cot.
- Do not adjust, roll, or load the cot into a vehicle without advising the patient. Stay with the patient and control the product at all times.
- You can transport the cot in any position. Stryker recommends that the operators transport the patient in the lowest comfortable position to maneuver the cot.
- Only use the wheel locks during patient transfer or without a patient on the product.
- Do not apply a wheel lock when an occupant is on the product or when you move the product to avoid the risk of tipping.
- Always use all restraint straps to secure the patient on the product. An unrestrained patient may fall from the product.
- Use trained helpers to control the cot, when necessary.

Proper lifting techniques

When you lift the product and patient, follow these proper lifting techniques to avoid the risk of injury:

- Keep your hands close to your body
- Keep your back straight
- Coordinate all movement with your partner
- Lift with your legs
- Avoid twisting

Transferring the patient to the cot

WARNING

- Always use all restraint straps to secure the patient on the product. An unrestrained patient may fall from the product and be injured.
- Do not leave a patient unattended. Hold the product while a patient is on the product.

- Do not apply a wheel lock when a patient is on the product or when you move the product to avoid the risk of tipping.
 - Do not use siderails as a patient restraint device.
-

To transfer the patient to the product:

1. Roll the product to the patient. See *Rolling the cot with a patient* (page 17).
2. Place the product beside the patient and raise or lower the product to the level of the patient.
3. Lower the siderails and open the restraint straps.
4. Transfer the patient to the product. Follow accepted EMS procedures.
5. Secure the patient to the product with all of the restraint straps.
6. Raise the siderails and adjust the backrest and footrest as necessary.

Rolling the cot with a patient

WARNING

- Always transport the cot at a lower height to reduce the risk of a cot tip. If possible, obtain additional assistance or take an alternate route.
 - Always avoid high obstacles, such as curbing, steps, or rough terrain to avoid the risk of the product tipping over.
-

To roll the cot with a patient:

1. Position one operator at the foot end and one operator at the head end of the cot.
2. Lift each set of wheels over the door sill or obstacle separately.

Changing the cot height

WARNING

- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
 - Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.
 - Always lock the base frame before you remove the loading wheels from the vehicle patient compartment floor or loading tray system. An unlocked base frame will not support the cot and injury to the patient or operator could result.
 - Always practice changing height positions and loading the cot until you fully understand the operation of the product.
 - Always keep hands, fingers, and feet away from moving parts. Use caution when you place your hands and feet near the base tubes while you raise or lower the cot.
-

Note

- Always make sure that both operators communicate to achieve the intended operation.
- One operator can lower the cot alone by releasing alternate ends of the cot.

Always use additional assistance if needed. A minimum of one trained operator is required for this task. You must be able to lift the total weight of the patient, cot, and any other items on the cot. The higher you must lift the cot, the more difficult it becomes to hold the weight. You may need help to load the cot into a vehicle if you are too short or if the patient is too heavy to safely lift.

See *Positioning operators and helpers* (page 19).

To change the cot height:

1. Position one operator at the foot end and one at the head end of the cot.
2. Grasp the hand grips on the base unit.
3. Lift the cot up approximately 1/4 in. until the weight is off of the latching mechanism.

4. Squeeze and hold the red release handle to raise or lower the cot to the desired height.
5. Release grip on the release handle to stop the cot in the next available cot position.

Loading the cot into a vehicle

WARNING

- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
 - Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.
 - Always keep hands, fingers, and feet away from moving parts. Use caution when you place your hands and feet near the base tubes while you raise or lower the cot.
 - Always support the load of the patient, cot, and accessories after the weight is off of the ground.
-

Note - Always remove loose items or debris from the patient compartment floor. Loose items or debris may interfere with the cot rolling smoothly or the operation of the cot faster.

Always use additional assistance if needed. A minimum of one trained operator is required for this task. You must be able to lift the total weight of the patient, cot, and any other items on the cot. The higher you must lift the cot, the more difficult it becomes to hold the weight. You may need help to load the cot into a vehicle if you are too short or if the patient is too heavy to safely lift.

See *Positioning operators and helpers* (page 19).

To load the cot into a vehicle:

1. Place the cot in a loading position. Make sure that load wheels meet the vehicle floor height.
2. Lift the vehicle bumper to the raised position, if equipped.
3. Roll the cot to the open door of the patient compartment or loading tray system.
4. Push the cot forward until the load wheels are on the patient compartment floor or the loading tray system.
5. Guide the cot forward until the front legs touch the rear of the vehicle or loading tray system.
6. At the foot end of the cot, push and hold the green release button while you squeeze and hold the green release handle.
7. Guide the cot forward until the rear legs touch the rear of the vehicle or loading tray system. The front legs will fold back.
8. Relax your grip on the green release handle and let go of the green release button.
9. Lift the foot end of the cot until the weight is off of the cot base.
10. Squeeze and hold the red release handle to guide the cot into the vehicle or loading tray system. The rear legs will fold back.
11. Make sure that the cot is secure in the cot fastener system.

Unloading the cot from a vehicle

WARNING

- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
 - Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.
 - Always lock the base frame before you remove the loading wheels from the vehicle patient compartment floor or loading tray system. An unlocked base frame will not support the cot and injury to the patient or operator could result.
 - Always keep hands, fingers, and feet away from moving parts. Use caution when you place your hands and feet near the base tubes while you raise or lower the cot.
 - Always support the load of the patient, cot, and accessories after the weight is off of the ground.
-

Always use additional assistance if needed. A minimum of one trained operator is required for this task. You must be able to lift the total weight of the patient, cot, and any other items on the cot. The higher you must lift the cot, the more difficult it becomes to hold the weight. You may need help to load the cot into a vehicle if you are too short or if the patient is too heavy to safely lift.

See *Positioning operators and helpers* (page 19).

To unload the cot from a vehicle:

1. Lift the vehicle bumper to the raised position, if equipped.
2. Release the cot from the cot fastener.
3. At the foot end of the cot, hold the hand grips while you squeeze and hold the red release handle.
4. Guide the cot out of the vehicle or loading tray system until the rear legs of the cot lower and fully extend.
5. Release the red release handle after the legs are fully extended.
6. Make sure that the rear legs are locked.
7. Guide the cot backward out of the vehicle or loading tray system until the front legs of the cot lower and fully extend.
8. Make sure that the front legs are locked. If required, slightly lift the cot to allow the legs to fully extend and lock.
9. Remove the load wheels from the patient compartment of the vehicle.
10. Roll the cot away from the vehicle.







Unloading the cot with the safety hook option

To unload the cot with the safety hook option, make sure that the cot legs fully extend and lock in place. The safety hook will not release the cot if the cot legs are not fully extended and locked.

Note - Using an M5 Allen key, adjust the hook tip (located on the hook arm) as needed to increase engagement.

Positioning operators and helpers

WARNING - Always follow proper hand placement on hand grips. Keep all hands clear of the red safety bar pivots when you load, unload, or change the cot height position.

	Changing levels	Rolling	Loading and unloading
Two operators (O) Two helpers (H)			
Two operators (O) Four helpers (H)			

Raising or lowering the backrest

To raise the backrest, squeeze the backrest adjustment handle to move the backrest to the desired height.

To lower the backrest, squeeze the backrest adjustment handle while you push down on the backrest to the desired height.

Raising or lowering the siderails

Always lower the siderails, if equipped, when you transfer a patient to the cot or from the cot.

WARNING - Do not use siderails as a patient restraint device.

To raise the siderails, lift up on the siderail until the latch clicks and the siderail locks into place. Always keep the siderails in the raised position unless you are transferring the patient.

To lower the siderails, squeeze the siderail release handle to release the siderail latch. Guide the siderail down toward the foot end of the cot until the siderail lays flat. Always lower the siderails when you transfer a patient to or from the cot.

Raising or lowering the footrest

You can adjust the footrest to elevate the patient's legs.

To raise the footrest, lift the frame as high as possible. The support bracket automatically secures the frame when you release the footrest.

To lower the footrest, lift the frame and lift up on the footrest release handle until the frame releases the support bracket. Lower the footrest until it lays flat.

Raising or lowering the knee Gatch option

To raise the knee Gatch:

1. Lift either of the red lifting loops (A) until the knee Gatch is in the highest position (Figure 13).
2. Lower the knee Gatch to secure the support bracket to the locking mechanism.
3. Make sure that the locking mechanism is secure before you release the lifting loop.

To lower the knee Gatch, lift either of the red lifting loops (A) to relieve pressure on the locking mechanism. While you hold the loop, push on the red footrest release handle (B) until the support bracket is released (Figure 13). Lower the knee Gatch until it lays flat.

To raise the knee Gatch in Trendelenburg, lift the footrest frame (C) as high as possible until the frame locks into place (Figure 13). The support bracket automatically connects when you release the frame.

To lower the knee Gatch in Trendelenburg, lift the footrest frame (C) and, while you hold the frame, lift up on the red footrest release handle (B) until the frame releases the support bracket (Figure 13). Lower the footrest until it lays flat.

Note - The foot end oxygen bottle holder option (6500-240-000) is not compatible with the knee Gatch option (6500-082-000).


Figure 13 – Gatch

Applying or releasing a wheel lock

WARNING

- Do not apply a wheel lock when a patient is on the product or when you move the product to avoid the risk of tipping.
- Do not install or apply a wheel lock on a product with worn wheels that are less than 6 in. diameter.
- Do not leave a patient or occupant unattended. Hold the product while a patient or occupant is on the product.

To apply a wheel lock, press down on the pedal until it stops and is resting against the surface of the wheel.

To release a wheel lock, press down on the top of the pedal with your foot or lift up on the pedal with your toe. The top of the pedal will rest against the caster frame when you release the wheel lock.

Note - Wheel locks help prevent the product from rolling while unattended. Wheel locks may not provide sufficient resistance on all surfaces or under loads.

Applying or releasing the caster Steer-Lock

The **Steer-Lock** feature applies the head end casters to improve steering control, when needed. You can release the **Steer-Lock** to allow for free swiveling for greater mobility.

WARNING

- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
- Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.

CAUTION - Do not apply the **Steer-Lock** and attempt to push the product sideways.

Note

- You can apply or release the caster **Steer-Lock** from either side of the base frame.
- Always release the caster **Steer-Lock** when you load the product into the cot fastener. This allows the head end casters to swivel freely so that the head end of the cot aligns with the fastener.

To apply the **Steer-Lock**, rotate the red locking knob into the locked position. Push the cot forward so the head end casters swivel and lock into position.

To release the **Steer-Lock**, rotate the red locking knob into the unlocked position. Push the cot in any direction. Make sure that the casters swivel freely.


Figure 14 – Steer-Lock knob location

Connecting the litter to the base

You can attach the Model 6100-031-000 knee Gatch litter and the Model 6100-033-000 flat foot section litter in either direction.

To connect the removable litter to the base (Figure 15):

1. Lower the litter onto the base.
2. Align the roller guides on the litter with the center base tube.
3. Push the litter forward toward the head end to attach to the litter retention brackets.
4. Make sure that you lock the litter onto the base before you release your grip on the litter.


Figure 15 – Connecting and removing the litter

Removing the litter from the base

To remove the removable litter from the base (Figure 15):

1. Push and hold the blue litter release button that is located at the foot end of the base.
2. Pull the litter toward the food end of the cot to remove the litter from the litter retention brackets.
3. Lift to remove the litter from the base.

Changing cot height with the optional base head end retractable lifting handles

WARNING

- Always lock the caster **Steer-Lock** knob and casters when you load or unload the cot from a vehicle or loading tray system to avoid the risk of the product tipping.
 - Always lock the caster **Steer-Lock** knob and casters when you change the cot height to avoid the risk of the product tipping.
 - Always lock the base frame before you remove the loading wheels from the vehicle patient compartment floor or loading tray system. An unlocked base frame will not support the cot and injury to the patient or operator could result.
 - Always practice changing height positions and loading the cot until you fully understand the operation of the product.
-

Note

- Always make sure that both operators communicate to achieve the intended operation.
- One operator can lower the cot alone by releasing alternate ends of the cot.

Always use additional assistance if needed. A minimum of one trained operator is required for this task. You must be able to lift the total weight of the patient, cot, and any other items on the cot. The higher you must lift the cot, the more difficult it becomes to hold the weight. You may need help to load the cot into a vehicle if you are too short or if the patient is too heavy to safely lift.

See *Positioning operators and helpers* (page 19).

To change the cot height:

1. Pivot the retractable lifting handles (A) into position (Figure 16).
2. Position an operator at the foot end and one at the head end of the cot.
3. Grasp the hand grips on the base unit.
4. Lift the cot up approximately 1/4 in. until the weight is off of the latching mechanism.
5. Squeeze and hold the red release handle to raise or lower the cot to the desired height.
6. Release grip on the release handle to stop the cot in the next available cot position.


Figure 16 – Positioning the lifting handles

Extending or retracting the retractable lifting handles

Note - Always make sure that you lock the retractable lifting handles before you lift the product.

To extend or retract the retractable lifting handles (Figure 17):

1. Press the release buttons to pull out or push in the retractable lifting handles.
2. Release the release buttons when the handles lock into the intermediate position.
3. Press the release buttons to pull out or push in the retractable lifting handles again.
4. Release the release buttons when the handles lock into the fully extended or retracted position.


Figure 17 – Extending and retracting the lifting handles

Securing the patient with the restraint straps

WARNING

- Always use all restraint straps to secure the patient on the product. An unrestrained patient may fall from the product and be injured.
 - Do not attach restraint straps to the base tubes or cross tubes.
-

CAUTION - Do not entangle the restraint straps in the base frame when you raise or lower the cot.

Secure restraint straps to the cot in the required attachment locations (Figure 18 and Figure 19). Restraint strap attachment locations should provide strong anchorage and proper restraint position. Do not allow restraint straps to interfere with equipment or accessories. Buckle the restraint straps across the patient's shoulders, waist, and legs. Buckle restraint straps when the cot is not in use.

To install the restraints:

1. Wrap the restraint strap around the cot frame.
2. Push the restraint strap buckle through the loop.
3. Pull the buckle around the loop to secure the restraint strap to the cot.


Figure 18 – Restraint strap attachment points


Figure 19 – Restraint strap attachment points - BS EN 1789:2007+A1:2010

Adjusting restraint straps

Open the restraint straps and place them at either side of the cot while you position the patient on the cot mattress. Lengthen the restraint straps, buckle them around the patient, and shorten them to tighten.

- To open the restraint strap, press the red button on the front of the buckle receiver. This allows you to release the buckle latch plate and pull it out of the receiver.
- To close the restraint strap, push the latch plate into the receiver until you hear a click.
- To lengthen the restraint strap, grasp the buckle latch plate, turn it at an angle to the webbing, then pull it out. A hemmed tab at the end of the webbing prevents the latch plate from coming off of the strap.
- To shorten the restraint strap, grasp the hemmed tab and pull the webbing back through the latch plate to tighten.

When you buckle a restraint strap around a patient, secure the latch plate and remove any loose webbing from the cot.

Inspect the restraint straps at least once a month (more if used often). Check for a bent or broken receiver or latch plate, or torn or frayed webbing. Replace any worn or inoperable restraint strap.


Figure 20 – Buckle the straps around the patient


Figure 21 – Lengthen the restraint strap


Figure 22 – Shorten the restraint strap

Adding a restraint strap extension

Add a restraint strap extension (6082-160-050) for extra length when you buckle the lap belt around larger patients.


Figure 23 – Restraint strap extension

Securing the patient with the Pedi-Mate® infant restraint system

WARNING - Always locate the buckle away from obstructions or accessories on the cot to avoid the risk of accidental release of **Pedi-Mate®** infant restraint system and injury to the infant.

To secure the patient with the **Pedi-Mate®** infant restraint system:

1. Remove any restraints from the cot.
2. Raise the cot backrest to the full upright position.
3. Position the **Pedi-Mate®** pad flat on the cot backrest with the black backrest straps facing out.
4. Wrap the straps around the cot backrest and insert the ends of the straps through the brackets.
5. Fasten the buckle.
6. Pull on the end of the adjustable backrest strap to tighten.
7. Insert the mainframe straps between the product frame and the mattress.
8. Insert the buckle behind the litter cross brace and bring it up in front of the litter cross brace.
9. Fasten the buckle around the litter cross brace. Leave slack in the strap for final adjustment.
10. Tighten all straps.

Note - For more information about the **Pedi-Mate®** infant restraint system, see the manufacturer's instructions for use, operation, and care. Safe and proper use of the **Pedi-Mate®** infant restraint system is at the discretion of the user. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product. These are general instructions only. Keep these instructions for future reference. These instructions are a permanent part of the product and should remain with the product even if the product is sold.

Attaching the defibrillator platform

CAUTION

- Always secure the defibrillator platform to the product when you use the defibrillator platform.
 - Always use and adjust the straps that are provided with the defibrillator platform to secure the defibrillator.
 - Always change the attachment location or adjust the straps for your specific defibrillator size or shape.
 - Do not load the defibrillator platform above the safe working load of 30 lb (13.6 kg).
-

To attach the defibrillator platform:

1. Place the defibrillator platform in the stored position (Figure 24).
2. Open and expand the defibrillator platform legs (Figure 25).


Figure 24 – Stored position


Figure 25 – Open defibrillator platform legs

3. If the product is equipped with an IV pole, raise the IV pole (A) to the up position (Figure 26).
4. Place the defibrillator platform on the product frame.
5. Position the inside legs (B) of the defibrillator platform toward the head end of the product (Figure 26).
6. Place the latch hook (C) under the litter frame or foot end fastener and push the tab (D) up until it locks into place on one side with an audible click (Figure 27). Repeat on the other side. For **Power-LOAD** compatible cots, if equipped, you must lengthen and attach the straps to the foot end fastener (Figure 28).


Figure 26 – Raise IV pole and position defibrillator tray (Model 6506 shown)


Figure 27 – Latch hook placement (Model 6506 shown)

7. Make sure that the defibrillator platform is secured to the product.
8. Place the defibrillator on the defibrillator platform.
9. Secure the straps (E) to hold the defibrillator on the defibrillator platform (Figure 29).


Figure 28 – Tab (Model 6506 shown)


Figure 29 – Attach defibrillator platform to the cot

Note - If the defibrillator platform is not attached to the product after you connect both latch hooks, unlatch the tab. To adjust, loosen or tighten the strap and then push the tab up until it locks into place with an audible click.

Attaching the head extension with pillow

You can attach the head extension assembly on the backrest to provide head end support.

To attach the pillow to the head extension, place the support into the flap on the bottom of the pillow. Secure the pillow to the **Velcro®** on the bottom of the support.

Note

- The head extension with pillow (6100-044-000) is not compatible with the equipment hook option (6500-147-000), Fowler oxygen bottle holder option (6500-241-000), or retractable head section oxygen bottle holder option (6085-046-000).
- The head extension pillow is a Type BF applied part.

Positioning the three-stage IV pole option

CAUTION - Do not load the IV pole above the safe working load of 25 lb (11.3 kg).

To position the IV pole (Figure 30):

1. Lift and pivot the IV pole from the storage position and push down until the IV pole locks into the receptacle (A).
2. To raise the height of the pole, turn the locking collar (B) counterclockwise and pull up on the telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the locking collar (B) clockwise to lock the telescoping portion in place.
4. For a higher IV pole, pull up on section (D) until the spring clip (E) connects.
5. Hang the IV bags on the IV hook (F).
6. To lower the IV pole, push in on the spring clip (E) and slide section (D) down into section (C). Turn the locking collar (B) counterclockwise and slide section (C) into the bottom tube.
7. Turn the locking collar (B) clockwise to tighten.
8. Lift up and pivot the pole down into the storage position.

Note - The dual three-stage IV poles (6500-317-000 or 6550-317-000) are not compatible with either the patient right (6500-315-000 or 6550-315-000) or the patient left (6500-316-000 or 6550-316-000) three-stage IV pole options.


Figure 30 – IV pole raised position

Attaching an oxygen bottle to the oxygen bottle holder

WARNING

- Do not use the oxygen bottle holder to hold an oxygen bottle holder when the transport vehicle is in motion. Always place the oxygen bottle holder in an appropriate storage location when the transport vehicle is in motion.
 - Always inspect the straps and clips for wear between use. Replace the strap if it no longer holds the oxygen bottle.
-

CAUTION

- Do not load the oxygen bottle holder above the safe working load of 15 lb (6.8 kg).
 - Do not use two oxygen bottle holders at the same time.
-

To attach an oxygen bottle to the oxygen bottle holder:

1. Place an oxygen bottle in the holder.
2. Insert the lower strap through the buckle and affix the strap onto itself to secure the oxygen bottle to the holder.

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Name	Number
Base option, non AS/NZS 4535	6100-003-950
Base option, AS/NZS 4535	6100-018-000
Catch plate kit	610000010022
Defibrillator platform	6100-170-010
DIN lock down pin, option	6100-042-000
HAVASU™ Euro IV pole, three-stage, right	6100-115-020
HAVASU™ Euro IV pole, three-stage, left	6100-116-020
Head extension with pillow	6100-044-000
Int. head end retractable lift handles	6100-010-000
IV pole option, right	6100-115-000
IV pole option, left	6100-116-000
Litter, flat foot section	6100-033-000
Litter, knee Gatch	6100-031-000
Mattress, bolster	6090-041-010
Mattress, flat	6090-042-010
Mattress, premium bolster	6100-041-010
No siderail, option	6100-021-000
Safety hook kit, option	610000010001
Pillow (premium mattress only)	6100-041-030
Removable oxygen bottle holder	6080-140-000
Siderail, fold down	6100-020-000

Cleaning

WARNING - Always use any appropriate personal protective equipment while power washing to avoid inhaling contagion. Power washing equipment may aerate contamination.

CAUTION

- Do not steam clean or ultrasonically clean the product.
 - Do not exceed 180° F (82° C) as the maximum water temperature.
 - Do not exceed 1500 psi (130.5 bar) as the maximum water pressure. If you are using a hand held wand to wash the product, keep the pressure nozzle at a minimum of 24 in. (61 cm) from the product.
 - Always allow to air dry.
 - Do not clean, disinfect, service, or perform maintenance while the product is in use.
-

The product is power washable. The product may show some signs of oxidation or discoloration from continuous washing. No degradation of the product's performance will occur from power washing as long as you follow the proper procedures.

- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- The preferred method Stryker Medical recommends for power washing the product is with the standard hospital surgical cart washer or hand held wand unit.
- Clean the cot once a month.
- Clean **Velcro®** after each use. Saturate **Velcro®** with disinfectant and allow disinfectant to evaporate. Appropriate disinfectant for nylon **Velcro®** should be determined by the service.
- Immersing restraint strap metal buckles can cause buckle corrosion and is not recommended. Rinse with clean water and allow to air dry to reduce chance of corrosion. Replace restraints if metal buckles are corroded.
- Direct skin contact with visibly soiled, permeable material may increase the risk of infection.

Note - Water that gets into the **M-1** Roll-In System will drain through the drain tube to the underside of the vehicle.

Disinfecting

CAUTION - Do not clean, disinfect, service, or perform maintenance while the product is in use.

In general, when used in concentrations recommended by the manufacturer, you can use either phenolic or quaternary (excluding **Virex® TB**) type disinfectants. Iodophor type disinfectants are not recommended for use because staining may occur.

The recommended disinfectants for this product's surfaces include the following:

- Quaternary (active ingredient - ammonium chloride)
- Phenolic (active ingredient - o-phenylphenol)
- Chlorinated bleach solution (use up to UK disinfecting 10,000 ppm available chlorine (941 mL of a 5.25% sodium hypochlorite solution per 4000 mL of water)
- Alcohol (active ingredient - 70% isopropyl alcohol)

To wipe down the product with disinfectant between uses:

1. Follow the manufacturer's dilution recommendations exactly.
2. Apply the recommended disinfectant solution by spray or pre-soaked wipes.
3. Hand wash all surfaces of the product with the recommended disinfectant.
4. Disinfect all exposed surfaces. Pay attention to high contact areas.
5. Follow the disinfecting solution manufacturer's instructions for appropriate contact time and rinsing requirements.
6. Dry the product thoroughly before returning the product to service.

Avoid over saturation. Do not allow the product to remain wet.

Follow the manufacturer's dilution recommendations for appropriate contact time and rinsing requirements. Follow the chemical manufacturer's guidelines for proper disinfecting.

Note

- Failure to follow the above directions when using these types of cleaners may void this product's warranty.
- Always wipe the product with clean water and dry after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature corrosion of critical components.

Preventive maintenance

CAUTION

- Always use authorized parts to avoid the risk of product damage.
- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Establish and follow a maintenance schedule and keep records of the maintenance activity. Remove product from service before you perform the preventive maintenance inspection. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

When using maintenance products, follow the directions of the manufacturer and reference all material safety data sheets (MSDS).

Operation	Schedule	Procedure
Cleaning and disinfecting	<i>Cleaning</i> (page 33) and <i>Disinfecting</i> (page 34)	
Inspection	For 1-25 calls per month, inspect the cot every 6 months For 26-200 calls per month, inspect the cot every 3 months For 201+ calls per month, inspect the cot monthly	See checklist below

Regular inspection and adjustments

The following schedule is a general guide to maintenance. Factors such as weather, terrain, geographical location, and individual usage will alter the required maintenance schedule. If you are unsure how to perform these checks, contact your Stryker service technician. If you are in doubt as to what intervals to follow to maintain your product, consult your Stryker service technician. Check each routine and replace worn parts if necessary.

Every month or two hours

Inspect these items every month or two hours, whichever comes first.

Item	Inspect
Settings	Cot and fastener fit and function
Litter	Frame and litter
Restraints	Restraints function with no excessive wear (such as bent or broken receiver or latch plate or torn or frayed webbing)
Base	Frame and base
Wheels	All wheels are secure, roll, and swivel
	Load wheels are secure and roll

Every three months or six hours

Inspect these items every three months or six hours, whichever comes first.

Item	Inspect
Litter	All fasteners are secure
	Backrest cylinder operates
	Adjust pneumatic cylinder for full range of motion, if required
Mattress	No cracks or tears
Base	All fasteners are secure
	Plastic wear strips are intact on the front and rear leg assemblies
	Cot front legs drop to a locked position under its own weight
Accessories and parts	All accessories and parts operate (such as IV pole, head extension pillow, defibrillator platform, and Pedi-Mate® restraint package)
	Velcro® is not peeling off of the guide bracket, replace if required
	Catch plate is not deformed, replace if required

Every six months or 12 hours

Inspect these items every six months or 12 hours, whichever comes first.

Item	Inspect
Litter	No bent, broken, or damaged components
	No damage or tears on cot grips
	Siderails operate and latch
	Footrest operates
Base	No bent, broken, or damaged components
Wheels	Free of debris

Every 12 months or 24 hours

Inspect these items every 12 months or 24 hours, whichever comes first.

Item	Inspect
Litter	All welds are intact, not cracked, or broken
	Warning labels present, legible
	Fowler bumper and rivets intact and secure
	Backrest operates

	Extension lift handles operate
	Siderails operate and latch
	Knee Gatch/Trendelenburg operates
	Litter retention pins are secure and not bent or broken
Base	All welds are intact, not cracked, or broken
Wheels	Rubber is in good condition
	Check and adjust wheel locks
	Check and adjust caster Steer-Lock System
Accessories and parts	Straps and clips on the oxygen bottle holder for wear
	Hook tip connects with catch plate while unloading, adjust hook tip location as needed
	Catch plate is secure and not deformed, replace if required
	Hook arm hinge is secure and operates
	Velcro® is not peeling off of the guide bracket, replace if required

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