Power-PRO™ XT

REF 6506

# stryker®

## Operations/Maintenance Manual

Manuel d'utilisation et d'entretien  
Bedienungs- und Wartungshandbuch  
Gebruiks-/onderhoudshandleiding  
Manuale d'uso e manutenzione  
Manual de uso y mantenimiento  
Manual de funcionamento/manutenção  
Användar-/underhållshandbok  
Betjenings-/vedligeholdelsesmanual  
Käyttö- ja huolto-ohjekirja  
Instrukcja obsługi/konserwacji  
Návod k obsluze a údržbě



CE

# Symbols

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## ROMANIA

**Stryker Osteonics Romania S.R.L.**  
19, Leonida Str.  
District 2  
7000 Bucharest  
Romania  
Phone: 40-2-12-12-11-22  
Fax: 40-2-12-12-11-33

## SOUTH AFRICA

**Stryker Osteonics PTY. LTD.**  
3 Susan Street  
Strydom Park  
Johannesburg, 2194 - South Africa  
Mailing Address:  
P.O. Box 48039  
2129 Roosevelt Park  
Johannesburg, South Africa  
Phone: 27-11-791-4644  
Fax: 27-11-791-4696

## SPAIN

**Stryker Iberia SL**  
c/Sepulveda n17  
28108 Alcobendas  
Madrid  
Spain  
Phone: 34-91-728-35-00  
Phone: 34-91-358-20-44  
Fax: 34-91-358-07-48

## SWEDEN

**Stryker AB/Scandinavia**  
Box 50425  
SE-204 14 Malmo  
Sweden  
Phone: 46-40-691-81-00  
Fax: 46-40-691-81-91



# Symbols

	Refer to instruction manual/booklet
	Operating instructions
	CE Mark
	Manufacturer
	Safe working load
	Dangerous voltage
	General warning
	Caution
	Warning; crushing of hands
	No pushing
	Do not lubricate
	SMRT™ Power System
	Extend
	Retract
	Warning; non-ionizing radiation
<b>IPX6</b>	Protection from powerful water jets



# Symbols

English

	Type B Applied Part
	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1: 2005 and CAN/CSA-C22.2 No. 60601-1:08.
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.
	Keep dry
	450 lb/205 kg weight capacity
	Do not stack more than three high
	Dangerous voltage
	This device complies with Part 15 of the FCC rules
	Date of manufacture
	Serial number
	Lot (batch) code
	Catalog/model number
	Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)
	Lift here
	This way up
	Fragile

# Symbols

English

	Two person lift
	Direct current
	European authorized representative
	English text below this symbol is intended for USA audiences only
	Approved by independent communications authority of South Africa
	Registered in United Arab Emirates by the Telecommunications Regulatory Authority
	Product complies with applicable EMC standards in Australia/New Zealand
	Box manufacturer's certificate - this packaging box has minimum bursting test value of 500 lbs per sq inch

[Return To Table of Contents](#)

# Table of Contents

English

Symbols . . . . .	<a href="#">1-1</a>
Warning/Caution/Note Definition . . . . .	<a href="#">1-8</a>
Introduction . . . . .	<a href="#">1-9</a>
Product Description . . . . .	<a href="#">1-9</a>
Intended Use of Product. . . . .	<a href="#">1-9</a>
Expected Service Life . . . . .	<a href="#">1-9</a>
Contraindications . . . . .	<a href="#">1-9</a>
Specifications . . . . .	<a href="#">1-10</a>
Contact Information . . . . .	<a href="#">1-12</a>
Serial Number Location . . . . .	<a href="#">1-12</a>
Serial Number Key . . . . .	<a href="#">1-12</a>
Product Illustration . . . . .	<a href="#">1-13</a>
Summary of Safety Precautions . . . . .	<a href="#">1-14</a>
Pinch Points . . . . .	<a href="#">1-19</a>
Mechanical Stability. . . . .	<a href="#">1-20</a>
Setup Procedures. . . . .	<a href="#">1-21</a>
Setting Cot Load Height and “Jog” Function . . . . .	<a href="#">1-22</a>
Cot Fastener Installation . . . . .	<a href="#">1-23</a>
Installing the In-Fastener Shut-Off . . . . .	<a href="#">1-25</a>
Vehicle Safety Hook Selection . . . . .	<a href="#">1-26</a>
Vehicle Safety Hook Installation . . . . .	<a href="#">1-27</a>
Vehicle Configuration . . . . .	<a href="#">1-27</a>
Required Hardware for Installation of the Safety Hook (Not Supplied) . . . . .	<a href="#">1-27</a>
Front to Back Positioning of the Safety Hook . . . . .	<a href="#">1-28</a>
Side to Side Positioning of the Safety Hook . . . . .	<a href="#">1-29</a>
Installing the Safety Hook. . . . .	<a href="#">1-29</a>
Power-PRO™ Cot User Controls . . . . .	<a href="#">1-30</a>
Using the Cot Control Switches . . . . .	<a href="#">1-30</a>
Checking the Cot Battery Power Level . . . . .	<a href="#">1-31</a>
Checking the Hour Meter/LCD Error Display . . . . .	<a href="#">1-32</a>
Operation Guide . . . . .	<a href="#">1-33</a>
Operating Guidelines . . . . .	<a href="#">1-33</a>
Proper Lifting Techniques. . . . .	<a href="#">1-33</a>
Transferring the Patient to the Cot . . . . .	<a href="#">1-34</a>
Rolling the Cot. . . . .	<a href="#">1-34</a>
Adjusting The Height of the Cot . . . . .	<a href="#">1-35</a>
Loading or Unloading the Cot. . . . .	<a href="#">1-36</a>
Loading or Unloading the Cot with the Power-LOAD Option . . . . .	<a href="#">1-36</a>
High Speed Retract/Extend . . . . .	<a href="#">1-36</a>
Loading the Cot into a Vehicle with Two Operators - Powered Method. . . . .	<a href="#">1-37</a>
Loading an Empty Cot into a Vehicle with One Operator - Powered Method. . . . .	<a href="#">1-39</a>
Unloading the Cot from a Vehicle with Two Operators - Powered Method . . . . .	<a href="#">1-40</a>
Unloading an Empty Cot from a Vehicle with One Operator - Powered Method . . . . .	<a href="#">1-42</a>
Using the Manual Override. . . . .	<a href="#">1-43</a>

# Table of Contents

English

Loading the Cot into a Vehicle with Two Operators - Manual Method . . . . .	<a href="#">1-44</a>
Unloading the Cot from a Vehicle with Two Operators - Manual Method . . . . .	<a href="#">1-46</a>
Unloading an Empty Cot from a Vehicle with One Operator - Manual Method . . . . .	<a href="#">1-48</a>
Using Additional Assistance . . . . .	<a href="#">1-49</a>
Operating the Backrest . . . . .	<a href="#">1-50</a>
Raising and Lowering the Siderails (Standard) . . . . .	<a href="#">1-50</a>
Raising and Lowering the Siderails (XPS Option) . . . . .	<a href="#">1-51</a>
Operating the Retractable Head Section . . . . .	<a href="#">1-52</a>
Adjusting the Footrest . . . . .	<a href="#">1-53</a>
Raising and Lowering the Optional Knee Gatch . . . . .	<a href="#">1-54</a>
Operating the Optional Wheel Locks . . . . .	<a href="#">1-55</a>
Operating the Optional Steer-Lock . . . . .	<a href="#">1-56</a>
Using the Optional Kickstand for Dialysis Scale . . . . .	<a href="#">1-57</a>
Using Restraint Straps . . . . .	<a href="#">1-58</a>
Using X-Restraint Straps . . . . .	<a href="#">1-60</a>
Adjusting Restraint Straps . . . . .	<a href="#">1-62</a>
Using the Restraint Extender . . . . .	<a href="#">1-63</a>
Attaching the Pedi-Mate® Infant Restraint System . . . . .	<a href="#">1-64</a>
Removing and Replacing a SMRT™ Pak . . . . .	<a href="#">1-66</a>
Using the Defibrillator Platform . . . . .	<a href="#">1-67</a>
Using the Equipment Hook . . . . .	<a href="#">1-69</a>
Using the Head Extension with Pillow . . . . .	<a href="#">1-69</a>
Operating the Optional Two-Stage I.V. Pole . . . . .	<a href="#">1-70</a>
Operating the Optional Three-Stage I.V. Pole . . . . .	<a href="#">1-71</a>
Attaching an Oxygen Bottle to an Oxygen Bottle Holder . . . . .	<a href="#">1-72</a>
Using the Retractable Head Section Oxygen Bottle Holder . . . . .	<a href="#">1-73</a>
Installing the Base Storage Net . . . . .	<a href="#">1-74</a>
Installing the Backrest Storage Pouch . . . . .	<a href="#">1-75</a>
Installing the Head End Storage Flat . . . . .	<a href="#">1-76</a>
Using the Transfer Flat . . . . .	<a href="#">1-76</a>
Attaching the Mattress . . . . .	<a href="#">1-77</a>
Cleaning . . . . .	<a href="#">1-78</a>
Washing Procedure . . . . .	<a href="#">1-78</a>
Washing Limitations . . . . .	<a href="#">1-78</a>
Removal of Iodine Compounds . . . . .	<a href="#">1-79</a>
Preventive Maintenance . . . . .	<a href="#">1-80</a>
Lubrication . . . . .	<a href="#">1-80</a>
Regular Inspection and Adjustments . . . . .	<a href="#">1-81</a>
Maintenance Record . . . . .	<a href="#">1-84</a>
Training Record . . . . .	<a href="#">1-85</a>
Troubleshooting Guide . . . . .	<a href="#">1-86</a>
Electronics and Hydraulics Locator . . . . .	<a href="#">1-86</a>
Hydraulic Assembly . . . . .	<a href="#">1-87</a>
Hydraulic Assembly Wiring Schematics . . . . .	<a href="#">1-87</a>

# Table of Contents

English

Electrical System Block Diagram . . . . .	<a href="#">1-88</a>
Troubleshooting Guide . . . . .	<a href="#">1-90</a>
LCD Error Codes . . . . .	<a href="#">1-93</a>
Main Cable Assembly . . . . .	<a href="#">1-94</a>
Main Cable Assembly Wiring Schematics . . . . .	<a href="#">1-94</a>
Control Board Assembly . . . . .	<a href="#">1-95</a>
Control Board Wiring Schematics . . . . .	<a href="#">1-95</a>
Quick Reference Replacement Parts List . . . . .	<a href="#">1-96</a>
Accessories . . . . .	<a href="#">1-98</a>
Service Information . . . . .	<a href="#">1-99</a>
Backrest Adjustment . . . . .	<a href="#">1-99</a>
Headsection Replacement . . . . .	<a href="#">1-100</a>
Backrest Gas Cylinder Replacement . . . . .	<a href="#">1-100</a>
Manual Release Cable Adjustment . . . . .	<a href="#">1-101</a>
Filling the Hydraulics Assembly Reservoir . . . . .	<a href="#">1-102</a>
Wheel Locking Force Adjustment . . . . .	<a href="#">1-103</a>
Steer-Lock Mechanism Adjustment . . . . .	<a href="#">1-104</a>
Cot Retaining Post Adjustment . . . . .	<a href="#">1-105</a>
Cot Retaining Post Replacement . . . . .	<a href="#">1-106</a>
Cot Retaining Post Screw Replacement . . . . .	<a href="#">1-106</a>
Hydraulic A Valve or B Valve Replacement . . . . .	<a href="#">1-107</a>
Hydraulic Manual Release Valve Replacement . . . . .	<a href="#">1-108</a>
Hydraulic Cylinder Replacement . . . . .	<a href="#">1-109</a>
Hydraulic Hose Replacement . . . . .	<a href="#">1-110</a>
Terminal Block Replacement . . . . .	<a href="#">1-111</a>
Siderail Assembly Replacement (Standard) . . . . .	<a href="#">1-112</a>
Siderail Assembly Replacement (XPS Option) . . . . .	<a href="#">1-113</a>
Ratchet Assembly Replacement (XPS Option) . . . . .	<a href="#">1-114</a>
Release Handle Assembly Replacement (XPS Option) . . . . .	<a href="#">1-115</a>
Spring Handle Assembly Replacement (XPS Option) . . . . .	<a href="#">1-115</a>
Child Restraint with X-Restraint Package Installation . . . . .	<a href="#">1-116</a>
Warranty . . . . .	<a href="#">1-117</a>
Stryker EMS Return Policy . . . . .	<a href="#">1-118</a>
Return Authorization . . . . .	<a href="#">1-118</a>
Damaged Merchandise . . . . .	<a href="#">1-118</a>
International Warranty Clause . . . . .	<a href="#">1-118</a>
Patent Information . . . . .	<a href="#">1-118</a>
EMC Information . . . . .	<a href="#">1-119</a>

# Warning/Caution/Note Definition

---

English

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

---

## **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

---

## **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

---

## **NOTE**

Provides special information to make maintenance easier or important instructions clearer.

# Introduction

---

This manual is designed to assist you with the operation and maintenance of the Stryker **Power-PRO™** XT cot. Read this manual thoroughly before using the equipment or beginning maintenance on it. To ensure safe operation of this equipment, it is recommended that methods and procedures be established for educating and training staff on the safe operation of this cot.

## PRODUCT DESCRIPTION

The Stryker Model 6506 **Power-PRO™** XT is a powered ambulance cot that consists of a platform mounted on a wheeled X-frame designed to support and transport a maximum weight of 700 lb (318 kg) in pre-hospital and hospital environments. The device is collapsible for use in emergency vehicles and has an adjustable load height feature to allow the device to be set to different ambulance deck heights for proper body mechanics during loading and unloading. The NiCd battery-powered hydraulic lift system allows operators to raise and lower the cot using the powered controls, while duplicate foot-end controls on the upper and lower lift bars accommodate different operator positions or sizes. The cot is equipped with a manual back-up release handle to allow the operation of cot functions in the event of power loss. The device is equipped with the following: a retractable head section for 360-degree mobility in any height position, side rails, patient securement straps, an adjustable pneumatic backrest and various optional accessories that assist with transport of the patient. Maximum patient comfort is attainable with the three different litter positions of shock, flat leg and optional knee gatch positioning.

## INTENDED USE OF PRODUCT

The Stryker **Power-PRO™** is a powered wheeled stretcher, which is intended to support and transport the entire body of a traumatized, ambulatory or non-ambulatory human patient (includes infants and adults). The battery-powered hydraulic lift system is intended to help reduce the effort required by the operator to raise and lower the cot. The device is designed to support patients in a supine (horizontal) or sitting position and facilitate the transportation of associated medical equipment (i.e. oxygen bottles, monitors, and/or pumps) in emergency/transport vehicles. This ambulance cot is intended to be used in pre-hospital and hospital environments, in emergency and non-emergency applications. It is rated to a maximum capacity of 700 lb (318 kg) (sum of the patient, mattress and accessory weight) and the intended operators of the device are trained professionals including emergency medical service and medical care center personnel, as well as medical first responders.

## EXPECTED SERVICE LIFE

- 7 years for **Power-PRO™** XT cot
- 7 years for **SMRT™** charger
- 2 years for **SMRT™** Pak battery

## CONTRAINDICATIONS

- **Power-PRO™** XT is not intended for extended stay or to be used as a hospital bed.
- **Power-PRO™** XT is not intended to be used in devices which modify air pressure, such as hyperbaric chambers.

# Introduction

English

## SPECIFICATIONS

	<p>Safe Working Load  <b>Note:</b> Safe Working Load indicates the sum of the patient, mattress and accessory weight.</p>	700 lb	318 kg
Maximum Unassisted Lift Capacity <sup>1</sup>		500 lb	227 kg
Backrest Articulation/Shock Position (Standard Fowler - 6506-012-003)		0° to 73° / +15°	
Backrest Articulation/Shock Position (1865 Fowler Option - 6506-012-004)		0° to 75° / +15°	
Overall Length/Minimum Length/Width		81 in / 63 in / 23 in	206 cm / 160 cm / 58 cm
Height <sup>2</sup>		Adjustable from 14 in to 41.5 in	Adjustable from 36 cm to 105 cm
Weight <sup>3</sup>		125 lb	57 kg
Caster Diameter/Width		6 in / 2 in	15 cm / 5 cm
Minimum Operators Required for Loading/ Unloading an Occupied Cot		2	
Minimum Operators Required for Loading/ Unloading an Unoccupied Cot		1	
Recommended Fastener Systems		Model 6370 or 6377 Floor Mount Type Model 6371 Wall Mount Type Model 6390 Power-LOAD	
Recommended Loading Height <sup>4</sup>		Up to 36 in	Up to 91 cm
Recommended Working Height (excluding mattress)		15.75 in.	40 cm
Single Adjustable Wheel Lock/ Double Adjustable Wheel Lock		Optional	
Hydraulic Oil		Stryker Part Number 6500-001-293	
Power System			
Battery		24V $\approx$ NiCd - <b>SMRT™</b> Power System	
Charger		100-240V $\sim$ 1.20 A, 50/60Hz or 12V $\approx$ 4.16 A - <b>SMRT™</b> Power System	
Cot Duty Cycle		16.7% (1 Min. On / 5 Min. Off)	
Standards (Cots and Chargers) <sup>5</sup>		ANSI/AAMI ES60601-1: 2005, CAN/CSA-C22.2 No. 60601-1:08, BS EN 1789:2007+A2:2014, AS/NZS-4535, KKK-A-1822, SAE J3027, BS EN 1865-2:2010 (for 1865 fowler option), BS EN 1865-3:2012 (for XPS option)	

<sup>1</sup> Cot loads over 300 lb (136 kg) may require additional assistance to meet the set cot load height.

<sup>2</sup> Height measured from bottom of mattress at seat section to ground level.







<sup>3</sup> Cot is weighed with one battery and without mattress and restraints.

<sup>4</sup> Cot may be set to any ambulance deck height ranging from 26" to 36" (66 cm to 91 cm).

<sup>5</sup> To meet SAE J3027 and AS/NZS-4535 crash-test standards with the use of a crash-rated fastener, such as Power-LOAD (Model 6390) (when equipped with the XPS option siderails), you must install the X-restraint package (6500-001-430) and knee gatch bolster mattress (6500-002-150/6506-002-150). To meet BS EN 1789:2007+A2:2014 crash-test standards with the use of a crash-rated fastener, such as Power-LOAD (Model 6390), you must install the G-rated restraint package (6500-002-030) and knee gatch bolster mattress (6500-002-150/6506-002-150). Cot is compliant to the BS EN 1865-3:2012+A1:2015 standard with the XPS option (6506-040-000). Cot is compliant to the BS EN 1865-2:2010+A1:2015 standard with the 1865 fowler option (6506-012-004). The Britax Meridian SICT Series No. 7200/A/2010 Convertible Child Restraint with the X-restraint package (6500-001-430) has been dynamically crash tested with a 10 kg crash dummy to 18.2 G forward and 10 G sideward per AS/NZS-4535: 1999 crash-test standards.

# Introduction

## SPECIFICATIONS (CONTINUED)

Environmental Conditions	Operation	Storage and Transportation
Temperature		
Relative Humidity		
Atmospheric Pressure		

Stryker reserves the right to change specifications without notice.

The **Power-PRO™** XT is designed to conform to the Federal Specification for the Star-of-Life Ambulance (KKK-A-1822).

The **Power-PRO™** XT is designed to be compatible with competitive cot fastener systems.

Patents pending.

**The yellow and black color scheme is a proprietary trademark of Stryker Corporation.**

Stryker hereby declares that this **Power-PRO™** XT ambulance cot (model 6506) is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. A copy of the original declaration of conformity can be obtained by contacting Stryker Medical at 3800 E. Centre Ave. Portage, MI 49002 Attn. Regulatory Affairs.

### CAUTION

- Changes or modifications to the unit not expressly approved by Stryker could void the user's authority to operate the system.
- This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their expense.

# Introduction

English

## CONTACT INFORMATION

Contact Stryker Customer Service or Technical Support at: (800) 327-0770.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

Please have the serial number (A) of your Stryker product available (as shown in Figure 1) when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

## SERIAL NUMBER LOCATION


Figure 1

## SERIAL NUMBER KEY

See Figure 1 and the following key for additional serial number information:

(1)	2 digit month
(2)	2 digit year
(3)	5 digit sequence that starts with 39000 each month

# Introduction

English

## PRODUCT ILLUSTRATION


Figure 2

[Return To Table of Contents](#)

# Summary of Safety Precautions

---

English

Carefully read and strictly follow the warnings and cautions listed on these pages. Service only by qualified personnel.

---

## WARNING

- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
  - To avoid the risk of patient or operator injury, use both hands when transporting the cot.
  - Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
  - Do not modify the cot or any components of the cot. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty (see [page 1-117](#)).
  - Any emergency vehicle to be used with this cot **must** have the in-fastener shut-off system installed (if not using Power-LOAD) (see [page 1-25](#)).
  - It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Cot Fastener Systems meets the installation specifications listed on [page 1-24](#). Injury may result if a non-compatible cot is used in the Stryker Fastener System.
  - The in-fastener shut-off must be positioned properly before placing the cot into service. Failure to install the in-fastener shut-off may cause injury to the patient or operator and/or damage to the vehicle.
  - Do not attempt to operate the cot when it is loaded into a cot fastener.
  - The in-fastener shut-off is **only** a means for disabling the electronic functionality. Damage to the product and/or injury to the patient and/or operator may occur if used for any other purpose.
  - Have the vehicle safety hook installed by a certified mechanic. Improper safety hook installation can cause injury to the patient or operator and/or damage to the cot.
  - Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-27](#).
  - The face of the safety hook that engages the safety bar should be located at least 3-3/4" from the leading edge of the door sill. After installation, verify that the cot legs lock into the load position without contacting the vehicle bumper.
  - To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.
  - Verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle or injury to the patient or operator and/or damage to the cot may occur.
  - The cot must have at least 5/8" of clearance between the vehicle bumper and the cot to disengage the safety bar when unloading the cot from the vehicle. Verify that the cot legs lock into the load position before disengaging the safety bar from the safety hook. Failure to properly lock the cot height into position can cause injury to the patient or operator and/or damage to the cot.
  - To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
  - Do not remove the battery when the cot is activated.
  - Avoid direct contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.
  - Entanglement in powered cot mechanisms can cause serious injury. Operate the cot only when all persons are clear of the mechanisms.
  - Inspect **SMRT™** Paks for damage before every use.
  - Practice changing height positions and loading the cot until operation of the product is fully understood. Improper use can cause injury.
  - Do not allow untrained assistants to assist in the operation of the cot. Untrained technicians/assistants can cause injury to the patient or themselves.
  - Do not ride on the base of the cot. Damage to the product could occur, resulting in injury to the patient or operator.
  - Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. Transporting the cot in a lowered position, head or foot end first, minimizes the potential of a cot tip.
  - Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.
-

# Summary of Safety Precautions

## WARNING (CONTINUED)

- Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.
- Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the product.
- Never apply the optional wheel locks while a patient is on the cot. Tipping could occur if the cot is moved while the wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
- Siderails are not intended to serve as a patient restraint device. See [page 1-58](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.
- Hydraulically raising or lowering the cot may temporarily affect electronic patient monitoring equipment. For best results, patient monitoring should be conducted when the cot is idle.
- High obstacles such as curbing, steps or rough terrain can cause the cot to tip, possibly causing injury to the patient or operator.
- If the cot is equipped with the optional kickstand, make sure that the kickstand remains in the retracted position and does not engage during transport.
- Transporting the cot in lower positions reduces the potential of a cot tip. If possible, obtain additional assistance or take an alternate route.
- Power-LOAD is designed to be compatible with the 6085/6086 Performance-PRO XT, 6500/6506 **Power-PRO™** XT, and 6510/6516 **Power-PRO™** IT cots with the Power-LOAD option only. In certain situations, you can use Power-LOAD as a standard antler for most X-frame cots, but a rail clamp assembly is required for all cots without the Power-LOAD option.
- It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Model 6390 Power-LOAD system is a Power-LOAD compatible cot. Injury may result if a non-compatible cot is used in the Stryker Model 6390 Power-LOAD system.
- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- When using a cot fastener, do not load the cot into the vehicle with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the cot.
- Whenever the weight of the cot and patient is no longer supported by the wheels, the cot will **automatically** enter the high speed retract mode if the retract (-) button is pressed.
- After the weight is off of the ground, the operators must support the load of the patient, cot and any accessories. Failure to support the load properly may cause injury to the patient or operator.
- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading/unloading a patient. Injury to the patient or operator could result.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Do not press the extend (+) button until the safety bar has engaged the safety hook.
- To avoid injury, always verify that the head section is locked into place prior to operating the cot.
- Do not attempt to load the cot into the patient compartment with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the product.
- Never install or use a wheel lock on a cot with excessively worn wheels. Installing or using a wheel lock on a wheel with less than a 6" diameter could compromise the holding ability of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the cot or other equipment.
- Do not attach restraints to the base tubes, cross tubes, or fowler skin. Improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.

[Return To Table of Contents](#)

# Summary of Safety Precautions

---

English

## WARNING (CONTINUED)

- To avoid the risk of patient injury or equipment damage, ensure that you properly mount and secure the defibrillator platform to the cot.
- To avoid the risk of patient injury or equipment damage, you must use the provided straps to secure the defibrillator to the defibrillator platform.
- Due to the difference in sizes and shapes of available defibrillators, you may have to change the location and adjustment of the straps that secure the defibrillator to the defibrillator platform. To avoid the risk of patient injury or equipment damage, use and adjust all straps properly to ensure the security of the defibrillator.
- To avoid the risk of patient injury or equipment damage, the weight placed on the defibrillator platform must not exceed 30 lb (13.6 kg).
- Stryker recommends a two person operation when using the kickstand.
- Make sure that the patient weight is centered on the cot before using the kickstand.
- Engage the kickstand with your foot only.
- Lower cot height prior to engaging kickstand for increased stability.
- Make sure that the kickstand remains in the retracted position and does not engage during transport.
- Do not use the kickstand as a brake.
- Do not engage kickstand on a sloped surface.
- If the cot is equipped with the optional retractable head section oxygen bottle holder, use caution while the oxygen bottle holder is installed to avoid pinching your fingers between the fowler bracket and the oxygen bottle.
- To avoid accidental release of the Pedi-Mate®, and possible injury to the infant, ensure that the restraint buckle is located away from obstructions on the cot or accessories.
- When the optional head end storage flat is being used, ensure that it does not interfere with the operation of the retractable head section, safety bar and safety hook. Injury to the patient or operator could result.
- When cleaning, use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the cot.
- **SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY.** If the products described above are used to clean Stryker EMS equipment, measures must be taken to ensure the cots are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the cots will leave a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.
- Failure to properly clean or dispose of contaminated mattress or other cot components will increase the risk of bloodborne pathogens and may cause injury to the patient or operator.
- Escaping fluid under pressure can penetrate the skin causing serious injury. Avoid the hazard by relieving pressure before disconnecting hydraulic or other lines. Tighten all connections before applying pressure. If an accident occurs, see a doctor immediately. Any fluid injected into the skin must be surgically removed within a few hours or gangrene may result. Doctors unfamiliar with this type of injury should reference a knowledgeable medical source.
- To avoid the risk of injury, do not use bare hands to check for hydraulic leaks.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like **Power-PRO™**. Install and place **Power-PRO™** into service according to the EMC information in this manual. Portable and mobile RF communications equipment can affect the function of **Power-PRO™**.
- The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Stryker as replacement parts for internal components may result in increased emissions or decreased immunity of the **Power-PRO™** cot.
- The Power-LOAD system and the **Power-PRO™** cot should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, observe the **Power-PRO™** cot to verify normal operation in the configuration in which it will be used.
- **Power-PRO™** operates at the following frequencies: 70 - 125 kHz for inductive charging and 13.56 MHz±7 kHz, Amplitude Modulated (OOK), ERP: -79.57 dBm. The **Power-PRO™** cot may be interfered with by other equipment, even if that other equipment complies with CISPR emission requirements.

# Summary of Safety Precautions

---

## CAUTION

- Changes or modifications to the unit not expressly approved by Stryker could void the user's authority to operate the system.
- This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their expense.
- The cot can be set at any cot load height position. Establish the required cot load height before placing the cot into service.
- Set the cot load height to the proper stop height prior to operation.
- Installation of the safety hook should be done by a certified mechanic familiar with ambulance vehicle construction. Consult the vehicle manufacturer before installing the safety hook and be sure that the installation of the safety hook does not damage or interfere with the brake lines, oxygen lines, fuel lines, fuel tank or electrical wiring of the vehicle.
- Only use the battery and charger as specified in the **SMRT™** Power System Operations/Maintenance Manual.
- The cot is not for use with an AC adapter.
- When charging a battery in an ambulance, locate the charger in an enclosed cabinet and out of patient reach during transport.
- Ensure that the battery is fully charged prior to placing into service. An uncharged or depleted battery may cause poor cot performance.
- Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.
- Do not "jog" the cot past the established cot load height of the product when the safety bar engages the vehicle safety hook or damage may occur to the product.
- When unloading the cot from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.
- Do not use the XPS option with a standard mattress. Use the wider gatch bolster mattress (6500-003-130) with the XPS option.
- Do not sit or stand on the siderails (XPS option).
- Do not use the siderails (XPS option) as a patient transfer device or surface (for example, to slide a patient from the cot to another surface).
- Do not position patients with full weight on the siderails (XPS option).
- Do not use the siderails (XPS option) as a push/pull device or to steer the unit.
- Remove the battery if the cot is not going to be used for an extended period of time (more than 24 hours).
- Wheel locks are only intended to help prevent the cot from rolling while unattended and to aid in patient transfer. A wheel lock may not provide sufficient resistance on all surfaces or under loads.
- Ensure that the restraints are not entangled in the base frame when raising and lowering the cot.
- The weight of the equipment in the base storage net (if equipped) must not exceed 20 lb (9 kg).
- Be careful when retracting the base to avoid damaging items stored in the base storage net.
- To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 35 lb (15.8 kg).
- To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 25 lb (11.3 kg).
- To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 15 lb (6.8 kg).
- Do not use two head end oxygen bottle holders at the same time.
- Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.
- The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 20 lb (9 kg).
- The weight of the equipment in the head end storage flat (if equipped) must not exceed 40 lb (18 kg).

# Summary of Safety Precautions

---

English

## CAUTION (CONTINUED)

- DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT.
  - Maximum water temperature should not exceed 180°F/82°C.
  - Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 24 inches (61 cm) from the unit.
  - Allow cot to air dry.
  - Towel dry all casters and interface points.
  - Failure to comply with these instructions may invalidate any/all warranties.
  - Always remove the battery before washing the cot.
  - A preventive maintenance program should be established for all Stryker EMS equipment. Preventive maintenance may need to be performed more frequently based on the usage level of the product. Close attention should be given to safety features including, but not limited to:
    - Hydraulic power mechanism
    - All electrical controls return to off or neutral position when released.
- For additional maintenance information, see the preventive maintenance information.
- Improper maintenance can cause injury or damage to the product. Maintain the cot as described in this manual. Use only Stryker approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty (see [page 1-117](#)).
  - Failure to use authorized parts, lubricants, etc. could cause damage to the cot and will void the warranty of the product.
  - Hydraulic lines, hoses, and connections can fail or loosen due to physical damage, kinks, age, and environment exposure. Check hoses and lines regularly to avoid damage to the cot. Check and tighten loose connections.
  - Do not tip the cot onto its load wheels and actuate the product as this will allow air to enter the hydraulic system.
  - Do not lubricate the bearings in the X-frame as it will degrade the performance of the cot and may void its warranty (see [page 1-117](#)).
  - The cot retaining post is shipped preconfigured for an X-frame cot. If the cot fastener has been configured for an H-frame cot, you must adjust the cot retaining post to accommodate the cot fastener.
-

# Summary of Safety Precautions

---

## PINCH POINTS


Figure 3

---

### **WARNING**

Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

---

# Summary of Safety Precautions

---

English

## MECHANICAL STABILITY


Figure 4

---

### WARNING

To avoid the risk of patient or operator injury, use both hands when transporting the cot.

---

#### Notes:

- If the cot is on a plane steeper than five degrees, place the cot in the fully lowered position.
- The defibrillator option and the foot end oxygen bottle holder option can not be used at the same time.

# Setup Procedures

---

Ensure that all shipping and packaging materials have been removed from the products prior to use.

Unpack the cartons and check all items for proper operation. It is important that the cot is working properly before it is put into service. See Figure 2 on [page 1-13](#) to identify all of the cot components.

The patient compartment of the vehicle in which the cot will be used must have a:

- Smooth rear edge for cot loading
- Level floor large enough for the folded cot
- Stryker Model 6370/6377/6378/6379 or 6371 Cot Fastener System or Stryker Model 6390 Power-LOAD (not included)
- In-fastener shut-off module installed and positioned properly (if not using Power-LOAD) (see [page 1-25](#))
- Space to properly install the safety hook

**Note:** Loose items or debris on the patient compartment floor can interfere with the operation of the safety hook and cot fastener. Keep the patient compartment floor clear.

When necessary, modify the vehicle to fit the cot. Do not modify the cot.

---

 **WARNING**

- Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
- Do not modify the cot or any components of the cot. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty (see [page 1-117](#)).
- Any emergency vehicle to be used with this cot **must** have the in-fastener shut-off system installed (if not using Power-LOAD) (see [page 1-25](#)).

---

**Note:**

- This manual should be considered a permanent part of the cot and should remain with the product even if the cot is subsequently sold.
- Stryker continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your cot and this manual. If you have any questions, please contact Stryker Customer Service or Technical Support at (800) 327-0770.

# Setup Procedures

English

## SETTING COT LOAD HEIGHT AND “JOG” FUNCTION

The cot control mechanism uses height sensors to set the load height stop for the cot. These height sensors match the load wheel height for a specific ambulance deck height.

The cot load height can be set from 26” to 36” (66 cm to 91 cm) as measured from the ground to the bottom of the load wheel. Determine the cot load height before placing the cot into service. You can modify the cot load height at any time, but you must determine and set the cot load height before the cot is placed into service.

### To set the cot load height:

1. Locate the sensor housing on the patient right side of the cot as shown in Figure 5.
2. Using a T27 Torx driver, remove the sensor housing cover by loosening the two (2) screws (one on each end) as shown in Figure 6.
3. Adjust the left height sensor only as shown in Figure 7.
  - a. Move the sensor to the left to increase the set load height or move the sensor to the right to decrease the set load height.
  - b. Press the retract (-) button to lower the cot to its lowest position, then press the extend (+) button to raise the cot to its highest set load height.
  - c. Measure the cot height from the bottom of the load wheels to the floor.

**Note:** Add an additional 1/2" (1,3 cm) to your deck height measurement to allow for variations with patient height and other equipment added to the cot.

- d. Repeat steps 3a and 3b until the desired cot load height is reached.
4. After the proper load wheel height is set, ensure that all of the height sensor cables are secure and lying flat inside of the housing between the sensors as shown in Figure 8.
  5. Using a T27 Torx driver, replace the sensor housing cover by reinstalling the two screws that were removed in step 2.
  6. Following completion of the sensor height adjustment, verify that the cot properly engages the safety hook.

Figure 5

Figure 6

Figure 7

Figure 8

### CAUTION

The cot can be set at any cot load height position. Establish the required cot load height before placing the cot into service.

# Cot Fastener Installation

---

**Note:** The Cot Fastener Installation instructions on [page 1-23](#) through [page 1-25](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

The Stryker Cot Fastener Systems are designed to be compatible only with cots which conform to the installation specifications listed on [page 1-24](#).

---

 **WARNING**

It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Cot Fastener Systems meets the installation specifications listed on [page 1-24](#). Injury may result if a non-compatible cot is used in the Stryker Fastener System.

---

**Note:** Adjustment of the rail clamp assembly may be required in order to compensate for any variation in the cot retaining post position depending on the cot manufacturer and model number.

For more information about the Stryker Cot Fastener Systems, see the Cot Fastener Operations/Maintenance Manual.

# Cot Fastener Installation

English

**Note:** These installation instructions are intended for cots with cot fastener systems (NOT Power-LOAD). For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.


Figure 9


Figure 10

# Cot Fastener Installation

## INSTALLING THE IN-FASTENER SHUT-OFF

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

### **WARNING**

The in-fastener shut-off must be positioned properly before placing the cot into service. Failure to install the in-fastener shut-off may cause injury to the patient or operator and/or damage to the vehicle.

The cot and fastener system have an integrated in-fastener shut-off function that disables the cot motor when the cot is secured into the cot fastener. Securely tighten the bolts on the fastener before installing the shut-off bracket. Install the shut-off bracket onto the rail clamp assembly before putting the cot into service.

1. Place the cot into a loading position (any position where the load wheels of the head section meet the vehicle floor height).
2. Lift the vehicle bumper to the raised position (if equipped).
3. Roll the cot to the open door of the patient compartment.
4. Push the cot forward until the load wheels are on the compartment floor and the safety bar passes the safety hook.
5. For maximum clearance to lift the base, pull the cot back until the safety bar engages the safety hook.
6. Raise the base and push the cot into the patient compartment following the appropriate loading instructions.
7. Engage the extended head section of the cot into the cot fastener antler and secure the cot post into the fastener rail clamp.
8. Adjust the shut-off bracket along the rail clamp until the “diamond” (A) on the sensor housing is lined up with the pop rivet head (B) as shown in Figure 11.
9. Using a T27 Torx driver, securely fasten the bolts to attach the shut-off bracket to the rail clamp assembly.
10. Press the retract (–) button to ensure that the motor does not turn on while the cot is in the fastener. The battery indicator will still illuminate. If the motor turns on, readjust the shut-off bracket.


Figure 11


Figure 12

### **WARNING**

- Do not attempt to operate the cot when it is loaded into a cot fastener.
- The in-fastener shut-off is **only** a means for disabling the electronic functionality. Damage to the product and/or injury to the patient and/or operator may occur if used for any other purpose.
- Any emergency vehicle to be used with this cot **must** have the in-fastener shut-off system installed (if not using Power-LOAD).

[Return To Table of Contents](#)

# Vehicle Safety Hook Selection

**Note:** The Vehicle Safety Hook Selection and Installation instructions on [page 1-26](#) through [page 1-29](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions. Power-LOAD ships and is installed with its own safety hook, thus no additional hook is needed.

The vehicle safety hook is a device that ships with the cot. The cot safety bar and vehicle safety hook are designed to keep the cot from being accidentally removed from the vehicle and to provide increased operator assurance and confidence when loading and unloading. The safety hook was designed for compatibility and proper operation when loading and unloading the cot from a vehicle that is compliant with Federal Regulation KKK-A-1822.

Stryker offers three different types of safety hooks that are ordered and shipped with your cot. These safety hook types are designed to meet the needs of various emergency vehicle configurations, specifically the length and location of the floor structure support that is located in the rear of the vehicle.

Consider the following information when selecting which safety hook is appropriate for your vehicle configuration:

- Determine the location of the floor structure support where there is adequate room to mount the safety hook.
- Ensure that the safety hook can be securely mounted into the back of the vehicle while providing adequate bumper clearance to allow the cot to be loaded and unloaded from the vehicle.
- Note the differences in vehicle design. Each safety hook provides a different mounting location option to maintain the appropriate distance between the face of the safety hook and the edge of the door sill.


Figure 13

Due to the differences in vehicle dimensions and the floor structure support locations, each safety hook requires a different mounting location. See "Vehicle Safety Hook Installation" to determine the correct positioning for safety hook installation.

**Note:** When replacing an existing safety hook with a new style, adjust the mounting location to maintain the proper position of the safety hook face.

# Vehicle Safety Hook Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## VEHICLE CONFIGURATION

According to federal regulations (reference KKK-A-1822), the bumper height of the vehicle shall be installed equidistant  $\pm 5$  cm (2 inches) from the vehicle floor to the ground level, which is defined as the vehicle deck height. The bumper step shall have a minimum depth of 13 cm (5 inches) and a maximum depth of 25 cm (10 inches). If the bumper depth is greater than 18 cm (7 inches), then the bumper must be able to fold. Installation of the safety hook into any vehicle compliant with this federal specification provides adequate clearance for the cot base to lower to its fully extended position. The cot is compatible with all vehicle deck heights (see specifications for maximum load height) as long as the vehicle meets the federal specifications that are outlined in KKK-A-1822.


Figure 14


Figure 15

### CAUTION

- Set the cot load height to the proper stop height prior to operation.
- Installation of the safety hook should be done by a certified mechanic familiar with ambulance vehicle construction. Consult the vehicle manufacturer before installing the safety hook and be sure that the installation of the safety hook does not damage or interfere with the brake lines, oxygen lines, fuel lines, fuel tank or electrical wiring of the vehicle.

## REQUIRED HARDWARE FOR INSTALLATION OF THE SAFETY HOOK (NOT SUPPLIED)

- (2) Grade 5, Minimum 1/4"-20 Socket Head Cap Screws\* for the short or long safety hook
- (2) Grade 5, Minimum 1/4"-20 Flat Socket Head Cap Screws\* for the J hook
- (2) Flat Washers
- (2) Lock Washers
- (2) 1/4"-20 Nuts

\* The length of the socket head cap screws depends on the thickness of the vehicle floor. Use screws that are long enough to go completely through the patient compartment floor, washer and nut by at least two full threads.

# Vehicle Safety Hook Installation

English

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## WARNING


- Have the vehicle safety hook installed by a certified mechanic. Improper safety hook installation can cause injury to the patient or operator and/or damage to the cot.
- Failure to install the safety hook can cause injury to the patient or operator.
- The face of the safety hook that engages the safety bar should be located at least 3-3/4" from the leading edge of the door sill. After installation, verify that the cot legs lock into the load position without contacting the vehicle bumper.
- To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.

**Note:** Stryker recommends that, prior to installation, the certified mechanic plan the placement of the safety hook in the rear of the vehicle.

Before installing the safety hook into your vehicle, check the front to back and side to side positioning when unloading and loading the cot to ensure that the safety hook will be installed properly. The cot safety bar must engage the safety hook every time, regardless of cot position.

## FRONT TO BACK POSITIONING OF THE SAFETY HOOK

1. Select the appropriate safety hook for your vehicle configuration.
2. Position the safety hook at least 3-3/4" from the leading edge of the door sill.
3. Ensure that the safety hook can be securely mounted into the back of the vehicle while providing adequate bumper clearance to allow the cot to be loaded and unloaded from the vehicle.
4. See "Side to Side Positioning of the Safety Hook" to confirm the side to side placement.

[Return To Table of Contents](#)

# Vehicle Safety Hook Installation

**Note:** These installation instructions are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for installation instructions.

## SIDE TO SIDE POSITIONING OF THE SAFETY HOOK

1. Remove the cot from the fastener and unload it from the vehicle.
2. While the cot is being removed, note the position of the load wheels and the safety bar.
3. Mark the center of the cot safety bar on the vehicle floor.
4. Verify that the position marked in Step 3 is where the safety bar engages the safety hook every time when unloading the cot in a variety of positions (all the way to the left and all the way to the right), regardless of cot position.
  - If the cot safety bar does not engage the safety hook in any of these positions (left, center, or right), modify the vehicle, not the cot or safety hook.
  - If the cot safety bar engages the safety hook every time, install the safety hook.

## INSTALLING THE SAFETY HOOK

1. Determine the correct safety hook front to back and side to side positioning, so the cot safety bar engages the safety hook every time.
2. Drill the holes for the screws.
3. Fasten the safety hook to the patient compartment floor and verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle.

### WARNING

- Verify that the safety hook always engages the cot safety bar regardless of how the cot is unloaded from the vehicle or injury to the patient or operator and/or damage to the cot may occur.
- The cot must have at least 5/8" of clearance between the vehicle bumper and the cot to disengage the safety bar when unloading the cot from the vehicle. Verify that the cot legs lock into the load position before disengaging the safety bar from the safety hook. Failure to properly lock the cot height into position can cause injury to the patient or operator and/or damage to the cot.


Figure 17


Figure 18

# Power-PRO™ Cot User Controls

English

## USING THE COT CONTROL SWITCHES

There are two identical cot control switches located on the **Power-PRO™** cots. Press the buttons on either of these switches to extend the cot, retract the cot, or release the cot from Power-LOAD (if applicable).

This Figure 19 and table highlight the three buttons located on the cot control switch.


Figure 19

Ref	Name	Description	Description (with use of Power-LOAD)
1	Release	Not applicable	Press to unlock the cot from Power-LOAD
2	Retract (-)	Press and hold to lower the litter or retract the cot undercarriage when loading	Press and hold to fully retract the cot undercarriage
3	Extend (+)	Press and hold to raise the litter or extend the cot undercarriage when loading	Press and hold to fully extend the cot undercarriage

# Power-PRO™ Cot User Controls

English

## CHECKING THE COT BATTERY POWER LEVEL

To check the battery power level, press the retract (–) button (A) as shown in Figure 20 on the cot control switch to activate the cot battery LED indicator (B) as shown in Figure 21 on [page 1-32](#).

The cot battery LED indicator is located at the **Power-PRO™** foot end control enclosure (shown as a battery symbol).

- The LED is solid green when the battery is fully charged or has adequately charged battery power.
- The LED flashes amber when the battery needs to be recharged or replaced.
- The LED is solid amber to indicate a battery error.


Figure 20

See the **SMRT™** Power System Operations/Maintenance Manual for additional **SMRT™** Pak and **SMRT™** Charger information.

### Notes:

- Automatic charging will only occur with **SMRT™** Pak batteries.
- Only use Stryker-approved batteries with **Power-PRO™**.
- If applicable, Power-LOAD automatically charges the **Power-PRO™ SMRT™** Pak battery when the cot is locked into Power-LOAD in the transport position (no cable or connectors required). The cot battery LED indicator momentarily flashes green to signify that it is charging.

### WARNING

- To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
- Do not remove the battery when the cot is activated.
- Avoid direct contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.

### CAUTION

- Only use the battery and charger as specified in the **SMRT™** Power System Operations/Maintenance Manual.
- The cot is not for use with an AC adapter.
- When charging a battery in an ambulance, locate the charger in an enclosed cabinet and out of patient reach during transport.
- Ensure that the battery is fully charged prior to placing into service. An uncharged or depleted battery may cause poor cot performance.

[Return To Table of Contents](#)

# Power-PRO™ Cot User Controls

English

## CHECKING THE HOUR METER/LCD ERROR DISPLAY

The hour meter (C), located on the foot end control enclosure, indicates the amount of time (HHH.H hours) that the hydraulics have been activated as shown in Figure 21. You can use the hour meter to determine the frequency for preventive maintenance procedures as listed on [page 1-80](#).

The error display (C), located on the foot end control enclosure, provides error code information for troubleshooting. See “LCD Error Codes” on [page 1-93](#).


Figure 21

# Operation Guide

---

English

## OPERATING GUIDELINES

- Use the cot only as described in this manual.
- Read all labels and instructions on the cot before using the cot.
- Before first and every use, inspect the **SMRT™** Pak housing and terminal area for cracks and/or damage.
- Loading or unloading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator. If additional assistance is needed, see “Using Additional Assistance” on [page 1-49](#).
- Do not adjust, roll or load the cot into a vehicle without advising the patient. Stay with the patient and control the cot at all times.
- The cot can be transported in any position. Stryker recommends that the operators transport the patient in the lowest comfortable position to maneuver the cot.
- Only use the wheel locks during patient transfer or without a patient on the cot.
- Do not leave wheel locks engaged while transporting the cot. Failure to do so may cause wheel damage.
- Always use the restraint straps.
- Use properly trained helpers, when necessary, to control the cot.

---

### **WARNING**

- Improper usage of the cot can cause injury to the patient or operator. Operate the cot only as described in this manual.
- Entanglement in powered cot mechanisms can cause serious injury. Operate the cot only when all persons are clear of the mechanisms.
- Inspect **SMRT™** Paks for damage before every use.
- Practice changing height positions and loading the cot until operation of the product is fully understood. Improper use can cause injury.
- Do not allow untrained assistants to assist in the operation of the cot. Untrained technicians/assistants can cause injury to the patient or themselves.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- Do not ride on the base of the cot. Damage to the product could occur, resulting in injury to the patient or operator.
- Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. Transporting the cot in a lowered position, head or foot end first, minimizes the potential of a cot tip.
- Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.
- Any emergency vehicle to be used with this cot **must** have the in-fastener shut-off system installed (if not using Power-LOAD) (see [page 1-25](#)).

---

### **CAUTION**

Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.

---

## PROPER LIFTING TECHNIQUES

When lifting the cot and patient, there are five basic guidelines to help you avoid injury:

- Keep your hands close to your body.
- Keep your back straight.
- Coordinate your movements with your partner and lift with your legs.
- Avoid twisting.
- Always operate the cot as described in this manual.

[Return To Table of Contents](#)

# Operation Guide

---

English

## TRANSFERRING THE PATIENT TO THE COT

### To transfer the patient to the cot:

1. Roll the cot to the patient.
2. Place the cot beside the patient and raise or lower the cot to the level of the patient.
3. Lower the siderails and open the restraint straps.
4. Transfer the patient to the cot using accepted EMS procedures.
5. Use all the restraint straps to secure the patient to the cot (see [page 1-58](#)).
6. Adjust the backrest and foot rest as necessary.

**Note:** When transferring larger patients, use of the Transfer Flat (6005-001-001) is recommended.

---

### WARNING

- Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.
  - Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the product.
  - Never apply the optional wheel locks while a patient is on the cot. Tipping could occur if the cot is moved while the wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
  - Siderails are not intended to serve as a patient restraint device. See [page 1-58](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.
  - Hydraulically raising or lowering the cot may temporarily affect electronic patient monitoring equipment. For best results, patient monitoring should be conducted when the cot is idle.
- 

## ROLLING THE COT

### When rolling the cot:

- Make sure that all of the restraint straps are securely buckled around the patient (see [page 1-58](#)).
- Position an operator at the foot end and one at the head end of the cot **at all times** when rolling the cot with a patient on it.
- Approach door sills and/or other low obstacles squarely and lift each set of wheels over the obstacle separately.

---

### WARNING

- High obstacles such as curbing, steps or rough terrain can cause the cot to tip, possibly causing injury to the patient or operator.
  - If the cot is equipped with the optional kickstand, make sure that the kickstand remains in the retracted position and does not engage during transport.
  - Transporting the cot in lower positions reduces the potential of a cot tip. If possible, obtain additional assistance or take an alternate route.
-

# Operation Guide

---

English

## ADJUSTING THE HEIGHT OF THE COT

---

### **WARNING**

- Grasping the cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the cot.
  - Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- 

You can raise or lower an unoccupied cot with one operator. If a patient is on the cot, a minimum of **two (2) trained operators** (one located at each end of the cot) are required to raise or lower the cot.

### **To raise or lower an unoccupied cot:**

1. **Operator 1 (Foot End)** – Grasp the cot frame at the foot end and press either the extend (+) button on the control switch to raise the litter or the retract (-) button on the control switch to lower the litter to the desired position.

### **To raise or lower the cot with a patient:**

1. **Operator 1 (Foot End)** – Grasp the cot frame at the foot end and press either the extend (+) button on the control switch to raise the litter or the retract (-) button on the control switch to lower the litter to the desired position.
2. **Operator 2 (Head End)** – Maintain a firm grip on the outer rail until the cot is securely in the desired position.

**Note:** If the extend (+) button on the control switch remains activated after reaching the set load height, the motor will remain halted until the operator releases the button. After the button is released, press the extend (+) button again to “jog” the cot height up further.

---

### **CAUTION**

Do not “jog” the cot past the established cot load height of the product when the safety bar engages the vehicle safety hook or damage may occur to the product.

---

[Return To Table of Contents](#)

# Operation Guide

---

English

## LOADING OR UNLOADING THE COT

The cot loading and unloading instructions on [page 1-36](#) through [page 1-49](#) are intended for cots that you will NOT use with Power-LOAD. For Model 6506 cots with the Power-LOAD option, see the Power-LOAD Operations/Maintenance Manual for loading and unloading instructions.

## LOADING OR UNLOADING THE COT WITH THE POWER-LOAD OPTION

The Model 6506 **Power-PRO™** XT cot is fully compatible with the Model 6390 Power-LOAD system if it is ordered with the Power-LOAD option or compatibility kit (6506-700-001).

For more information about using your Power-LOAD compatible cot, see the Power-LOAD Operations/Maintenance Manual.

---

### **WARNING**

- Power-LOAD is designed to be compatible with the 6085/6086 Performance-PRO XT, 6500/6506 **Power-PRO™** XT, and 6510/6516 **Power-PRO™** IT cots with the Power-LOAD option only. In certain situations, you can use Power-LOAD as a standard antler for most X-frame cots, but a rail clamp assembly is required for all cots without the Power-LOAD option.
  - It is the responsibility of the cot operator to ensure that the cot being used in the Stryker Model 6390 Power-LOAD system is a Power-LOAD compatible cot. Injury may result if a non-compatible cot is used in the Stryker Model 6390 Power-LOAD system.
- 

## HIGH SPEED RETRACT/EXTEND

The cot is equipped with a high-speed retract mode to expedite loading/unloading the cot into and out of a vehicle.

- The undercarriage **rapidly** retracts toward the highest position once the weight of the cot and patient is no longer supported by the wheels. Press the retract (–) button to actuate the control switch.
  - The undercarriage **rapidly** extends toward the lowest position once the weight of the cot and patient is no longer supported by the wheels. Press the extend (+) button to actuate the control switch.
- 

### **WARNING**

- Whenever the weight of the cot and patient is no longer supported by the wheels, the cot will **automatically** enter the high speed retract mode if the retract (–) button is pressed.
  - After the weight is off of the ground, the operators must support the load of the patient, cot and any accessories. Failure to support the load properly may cause injury to the patient or operator.
-

# Operation Guide

English

## LOADING THE COT INTO A VEHICLE WITH TWO OPERATORS - POWERED METHOD

Loading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator.

### WARNING

- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-27](#).

### To load the cot into a vehicle with two operators:

1. Ensure that the retractable head section is fully extended and locked.
2. Place the cot in a loading position (any position where the load wheels meet the vehicle floor height).
3. Lift the vehicle bumper to the raised position (if equipped).
4. Roll the cot to the open door of the patient compartment.
5. Push the cot forward until the load wheels are on the compartment floor and the safety bar passes the safety hook as shown in Figure 22.
6. For maximum clearance to lift the base, pull the cot back until the safety bar engages the safety hook.
7. **Operator 2** – Verify that the safety bar engages the safety hook.

Figure 22

[Return To Table of Contents](#)

# Operation Guide

English

## LOADING THE COT INTO A VEHICLE WITH TWO OPERATORS - POWERED METHOD (CONTINUED)

8. Load the cot either from the foot end or with one operator at the foot end and one on the side:

### With both operators at the foot end (preferred method):

- **Both Operators** – Grasp the cot frame at the foot end (Figure 23).
- **Operator 1** – Press the retract (-) button until the undercarriage of the cot retracts fully (Figure 24).

Figure 23

Figure 24

### With one operator at the foot end and one on the side:

- **Operator 1** – Grasp the cot frame at the foot end and press the retract (-) button (Figure 25) until the undercarriage of the cot retracts fully (Figure 25).

Figure 25

Figure 26

- **Operator 2** – Securely grasp the cot outer rail to stabilize the cot during retraction.

9. **Both Operators** – Push the cot into the patient compartment as shown in Figure 24 or Figure 26 until the cot engages the cot fastener (not included).

### WARNING

When using a cot fastener, do not load the cot into the vehicle with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the cot.

[Return To Table of Contents](#)

# Operation Guide

## LOADING AN EMPTY COT INTO A VEHICLE WITH ONE OPERATOR - POWERED METHOD

Loading an **unoccupied** cot into the emergency vehicle can be accomplished by a single operator.

### **WARNING**

- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading/unloading a patient. Injury to the patient or operator could result.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

Figure 27

### To load an empty cot into a vehicle with one operator:

1. Place the cot into a loading position (any position where the load wheels of the head section meet the vehicle floor height).
2. Lift the vehicle bumper to the raised position (if equipped).
3. Roll the cot to the open door of the patient compartment.
4. Push the cot forward until the load wheels are on the patient compartment floor (Figure 27) and the safety bar passes the safety hook.
5. For maximum clearance to lift the base, pull the cot back until the safety bar engages the safety hook.
6. Grasp the cot frame at the foot end and press the retract (-) button, until the undercarriage of the cot retracts into its highest position as shown in Figure 28.
7. Push the cot into the patient compartment until the cot engages the cot fastener (not included) as shown in Figure 29.

Figure 28

### **WARNING**

When using a cot fastener, do not load the cot into the vehicle with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the cot.

Figure 29

# Operation Guide

English

## UNLOADING THE COT FROM A VEHICLE WITH TWO OPERATORS - POWERED METHOD

Unloading an occupied cot from a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator.

### WARNING

- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-27](#).
- To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Do not press the extend (+) button until the safety bar has engaged the safety hook.

### To unload the cot from a vehicle with two operators:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-24](#)).
3. Unload the cot either from the foot end or with one operator at the foot end and one on the side:

#### With both operators at the foot end (preferred method):

- **Both Operators** – Grasp the cot frame at the foot end. Pull the cot out of the patient compartment until the safety bar engages the safety hook as shown in Figure 30.
- **Both Operators** – Verify that the safety bar engages the safety hook.
- **Operator 1** – Depress the extend (+) button to lower the undercarriage to its fully extended position (Figure 31).

**Note:** You can use the manual release or a combination of the manual release followed by the extend (+) button. If the extend (+) button is used, you must ensure that the manual release is fully engaged before pressing the extend (+) button.

Figure 30

Figure 31

# Operation Guide

## UNLOADING THE COT FROM A VEHICLE WITH TWO OPERATORS - POWERED METHOD (CONTINUED)

With one operator at the foot end and one on the side (Figure 32):

- **Operator 1** – Grasp the cot frame at the foot end. Pull the cot out of the patient compartment until the safety bar engages the safety hook.
- **Operator 2** – Verify that the safety bar engages the safety hook.
- **Operator 2** – Stabilize the cot during the unloading operation by securely grasping the outer rail.
- **Operator 1** – Depress the extend (+) button to lower the undercarriage to its fully extended position as shown in Figure 33.

**Note:** You can use the manual release or a combination of the manual release followed by the extend (+) button. If the extend (+) button is used, you must ensure that the manual release is fully engaged before pressing the extend (+) button.

Figure 32

Figure 33

4. **Operator 2** – Pull the safety bar release lever forward to disengage the safety bar from the safety hook in the patient compartment as shown in Figure 34.
5. Remove the load wheels from the patient compartment of the vehicle.

### CAUTION

- When unloading the cot from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.
- Do not “jog” the cot past the established cot load height of the product when the safety bar engages the vehicle safety hook or damage may occur to the product.


Figure 34

# Operation Guide

English

## UNLOADING AN EMPTY COT FROM A VEHICLE WITH ONE OPERATOR - POWERED METHOD

Unloading an **unoccupied** cot from a vehicle can be accomplished by a single operator.

### WARNING

- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading or unloading a patient. Injury to the patient or operator could result.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

### To unload an empty cot from a vehicle with one operator:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-24](#)).
3. Grasp the cot frame at the foot end.
4. Pull the cot from the vehicle until the safety bar engages the safety hook as shown in Figure 35.
5. Depress the extend (+) button to lower the undercarriage to its fully extended position as shown in Figure 36.
6. Disengage the safety bar from the safety hook by pulling the safety bar release lever forward and roll the cot out of the vehicle as shown in Figure 37.
7. Remove the load wheels from the patient compartment of the vehicle.

### CAUTION

- When unloading the cot from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.
- Do not “jog” the cot past the established cot load height of the product when the safety bar engages the vehicle safety hook or damage may occur to the product.

Figure 35

Figure 36

Figure 37

# Operation Guide

English

## USING THE MANUAL OVERRIDE

In the event of loss of electrical function, the cot is equipped with a manual override to allow manual operation of the product until electrical functionality is restored. You can use the red manual back-up release handle to raise or lower the cot.

The **red** manual back-up release handle (A) is located along the patient left side of the lower lift bar at the foot end of the cot as shown in Figure 38.


Figure 38

### To raise or lower the cot with the manual back-up release handle:

1. **Both Operators** – Lift the cot during the raise/lower operation to support the weight of the cot at each end.
2. **Operator 1 (Foot End)** – Pull the manual back-up release handle toward the lift bar. While the manual back-up release handle is pulled, raise or lower the cot to the desired position and then release the handle to lock the cot into position.

### Notes:

- The operators must lift the cot weight slightly off of the wheels to use the manual extend or retract while a patient is on the cot.
- Activation of the manual back-up release handle may cause the cot to lower slowly if less than 50 lb (23 kg) are on the cot.
- Hydraulic fluid will become more viscous when the cot is used for extended periods in cold temperatures. When using the manual back-up release function to extend the base during unloading in cold weather conditions, hold the release handle for approximately one second after the cot wheels touch the ground to minimize sagging of the litter as the cot is removed from the ambulance.

[Return To Table of Contents](#)

# Operation Guide

---

English

## LOADING THE COT INTO A VEHICLE WITH TWO OPERATORS - MANUAL METHOD

Loading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator.

---

### **WARNING**

- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-27](#).

---

### **To load the cot into a vehicle with two operators using the manual back-up release handle:**

1. Place the cot in a loading position (any position where the load wheels meet the vehicle floor height).
2. Lift the vehicle bumper to the raised position (if equipped).
3. Roll the cot to the open door of the patient compartment.
4. Push the cot forward until the load wheels are on the patient compartment floor and the safety bar passes the safety hook.
5. For maximum clearance to lift the base, pull the cot back until the safety bar engages the safety hook.

# Operation Guide

## LOADING THE COT INTO A VEHICLE WITH TWO OPERATORS - MANUAL METHOD (CONTINUED)

6. **Operator 2** – Verify that the safety bar engages the safety hook.
7. **Operator 1** – Grasp the cot frame at the foot end. Lift the foot end of the cot until the weight is off of the cot base. Squeeze and hold the release handle as shown in Figure 39.
8. **Operator 2** – Stabilize the cot by placing your hand on the outer rail. Grasp the base frame as shown in Figure 40. After the foot end operator has lifted the cot and squeezed the release handle, raise the undercarriage until it stops in the highest position and hold it there.
9. **Both Operators** – Push the cot into the patient compartment, engaging the cot fastener (not included) as shown in Figure 41.

**Note:** When operating the manual back-up release handle, avoid rapid lifting or lowering of the base or movement may appear sluggish; lift with a slow constant motion.

Figure 39

Figure 40

Figure 41

# Operation Guide

English

## UNLOADING THE COT FROM A VEHICLE WITH TWO OPERATORS - MANUAL METHOD

Unloading an occupied cot from a vehicle requires a minimum of **two (2) trained operators**. One or two operators can lift from the foot end of the cot. Stryker recommends that both operators are at the foot end to reduce the load on each operator.

### WARNING

- Two operators must be present when the cot is occupied.
- Operators must be able to lift the total weight of the patient, cot and any items on the cot.
- The higher an operator must lift the cot, the more difficult it becomes to hold the weight. An operator may need help loading the cot if he/she is too short or if the patient is too heavy to lift safely. The operator must be able to lift the cot high enough for the cot legs to unfold completely and lock when the cot is unloaded. A shorter operator needs to raise their arms higher to enable the undercarriage to unfold.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.
- There must be a safety hook properly installed in the vehicle so that the bumper does not interfere with the front legs of the base frame.
- Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described on [page 1-27](#).
- To avoid injury, verify that the safety bar has engaged the safety hook before removing the cot from the patient compartment.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Do not press the extend (+) button until the safety bar has engaged the safety hook.

### To unload the cot from a vehicle with two operators:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-24](#)).
3. Unload the cot either from the foot end or with one operator at the foot end and one on the side:

#### With both operators at the foot end (preferred method):

- **Both Operators** – Grasp the cot frame at the foot end as shown in Figure 42.
- **Operator 1** – Pull the manual back-up release handle to lower the undercarriage to its fully extended position (Figure 43). Pull the cot out of the patient compartment until the safety bar engages the safety hook.
- **Operator 2** – Verify that the safety bar engages the safety hook.

Figure 42

Figure 43

# Operation Guide

## UNLOADING THE COT FROM A VEHICLE WITH TWO OPERATORS - MANUAL METHOD (CONTINUED)

With one operator at the foot end and one on the side (Figure 44):

- **Operator 1** – Grasp the cot frame at the foot end. Pull the manual back-up release handle to lower the undercarriage to its fully extended position (Figure 45). Pull the cot out of the patient compartment until the safety bar engages the safety hook.
- **Operator 2** – Verify that the safety bar engages the safety hook.
- **Operator 2** – Stabilize the cot during the unloading operation by securely grasping the outer rail.

Figure 44

Figure 45

4. **Operator 2** – Pull the safety bar release lever forward to disengage the safety bar from the safety hook in the patient compartment (Figure 46).
5. Remove the load wheels from the patient compartment of the vehicle.

### CAUTION

When unloading the cot from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.


Figure 46

# Operation Guide

English

## UNLOADING AN EMPTY COT FROM A VEHICLE WITH ONE OPERATOR - MANUAL METHOD

Unloading an **unoccupied** cot from a vehicle can be accomplished by a single operator.

### WARNING

- The one person loading and unloading procedures are for use only with an empty cot. Do not use the procedures when loading or unloading a patient. Injury to the patient or operator could result.
- Do not pull or lift on the safety bar when unloading the cot. Damage to the safety bar could result and injury to the patient or operator could occur.
- Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

### To unload an empty cot from a vehicle with one operator:

1. Lift the vehicle bumper to the raised position (if equipped).
2. Disengage the cot from the cot fastener. (For more information about the cot fastener, see [page 1-24](#)).
3. Grasp the cot frame at the foot end.
4. Pull the cot from the vehicle until the safety bar engages the safety hook as shown in Figure 47.
5. Pull the manual back-up release handle to lower the undercarriage to its fully extended position as shown in Figure 48.
6. Disengage the safety bar from the safety hook by pulling the safety bar release lever forward and roll the cot out of the vehicle.
7. Remove the load wheels from the patient compartment of the vehicle.

Figure 47

Figure 48

### CAUTION

When unloading the cot from the patient compartment, ensure that the caster wheels are safely set on the ground or damage to the product may occur.

# Operation Guide

## USING ADDITIONAL ASSISTANCE

	Changing Levels	Rolling	Loading/Unloading
Two Operators Two Helpers			
Two Operators Four Helpers			

Figure 49

**⚠ WARNING**

Ensure proper hand placement on hand grips. Hands should be clear of red safety bar pivots while loading and unloading the cot or whenever changing height position of the cot with two or more operators.

# Operation Guide

English

## OPERATING THE BACKREST

**To raise the backrest**, as shown in Figure 50, squeeze handle (A) for pneumatic assist in lifting the backrest to the desired height.

**To lower the backrest**, squeeze handle (A) and push down on the backrest frame until the backrest has reached the desired height.


Figure 50

## RAISING AND LOWERING THE SIDERAILS (STANDARD)

**To raise the siderails**, as shown in Figure 50, lift up on the siderail until the latch clicks and the siderail locks into place.

**To lower the siderails**, squeeze handle (B) to release the siderail latch. Guide the siderail down toward the foot end until flat. Ensure that the siderails are lowered when a patient is being transferred to or from the cot.

### **WARNING**

Siderails are not intended to serve as a patient restraint device. See [page 1-58](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.

[Return To Table of Contents](#)

# Operation Guide

## RAISING AND LOWERING THE SIDERAILS (XPS OPTION)

### WARNING

Siderails are not intended to serve as a patient restraint device. See [page 1-58](#) for proper restraint strap usage. Failure to use the restraint straps properly could result in patient injury.

You can order your cot with the eXpandable patient surface (XPS) option or upgrade your cot to add the XPS option. Siderails (XPS option) are attached to the cot and are always available for your use. The siderails (XPS option) adjust according to patient size and lock into seven positions. The siderails also adjust to fit through standard doorways or elevators.

**To raise the siderails**, lift up on the siderail until it locks into the desired position.

**To lower the siderails**, lift up to relieve the weight, then pull the red lever (A) (Figure 51).

The XPS option is not a primary patient support surface. It includes a wider mattress and is intended to enhance patient comfort.

### CAUTION

- Do not use the XPS option with a standard mattress. Use the wider gatch bolster mattress (6500-003-130) with the XPS option.
- Do not sit or stand on the siderails (XPS option).
- Do not use the siderails (XPS option) as a patient transfer device or surface (for example, to slide a patient from the cot to another surface).
- Do not position patients with full weight on the siderails (XPS option).
- Do not use the siderails (XPS option) as a push/pull device or to steer the unit.


Figure 51

## OPERATING THE RETRACTABLE HEAD SECTION

The head section telescopes from a first position suitable for loading the cot into an emergency vehicle to a second position retracted within the litter frame. When retracted, the cot can roll in any direction on the caster wheels even in the lowest position, allowing for improved mobility and maneuverability.

### To extend the head section:

1. Grasp the outer rail with one hand for support and pull the handle (A), rotating the handle toward the head end of the cot to release the head section from the locked position.
2. While holding the handle (A) in the released position, pull the head section away from the litter frame, lengthening the head section until it engages in the fully extended position.
3. Release handle (A) to lock the head section in the extended position.

### To retract the head section:

1. Grasp the outer rail with one hand for support and release the handle (A), rotate the handle toward the head end of the cot to release the head section from the locked position.
2. While holding the handle (A) in the released position, push the head section toward the litter frame, retracting the head section until it engages in the retracted position.
3. Release handle (A) to lock the head section in the retracted position.

### WARNING

- To avoid injury, always verify that the head section is locked into place prior to operating the cot.
- Do not attempt to load the cot into the patient compartment with the head section retracted. Loading the cot with the head section retracted may cause the product to tip or not engage properly in the cot fastener, possibly causing injury to the patient or operator and/or damage to the product.


Figure 52


Figure 53


Figure 54

# Operation Guide

## ADJUSTING THE FOOTREST

The footrest is adjustable to elevate the patient's legs (see Figure 55).

**To raise the footrest**, lift the footrest frame (A) as high as possible. The support bracket engages automatically when released.

**To lower the footrest**, lift the footrest frame (A) and while holding the frame, lift up on the release handle (B) until the bracket disengages. Carefully lower the footrest until flat.


Figure 55

# Operation Guide

English

## RAISING AND LOWERING THE OPTIONAL KNEE GATCH

### To raise the knee gatch (see Figure 56):

1. Lift either of the red lifting loops (A) until the knee gatch is in its fully raised position.
2. Slowly lower the knee gatch to allow the support bracket to engage in the locking mechanism.
3. Check to be sure the lock is fully engaged before releasing the lifting loop.

### To lower the knee gatch:

1. Lift either of the red lifting loops to relieve pressure on the locking mechanism and while holding the loop, push on the red release handle (B) until the bracket disengages.
2. Carefully lower the knee gatch to the flat position.


Figure 56

### To raise the knee gatch in trend (see Figure 56):

1. Lift the foot rest frame (C) as high as possible until it locks into place.
2. The support bracket engages automatically when released.

### To lower the knee gatch in trend:

1. Lift the foot rest frame (C) and, while holding the frame, lift up on the release handle (B) until the bracket disengages.
2. Carefully lower the footrest until it rests flat.

# Operation Guide

English

## OPERATING THE OPTIONAL WHEEL LOCKS

To activate the optional wheel locks, press fully down on the pedal (A) as shown in Figure 57 until it stops and is resting firmly against the surface of the wheel.

To release the optional wheel locks, depress the upper face of the pedal with your foot or lift up with your toe under the pedal. The upper portion of the pedal will rest against the caster frame when the wheel lock is released.


Figure 57

### WARNING

- Never apply the optional wheel locks while a patient is on the cot. Tipping could occur if the cot is moved while a wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
- Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the cot.
- Never install or use a wheel lock on a cot with excessively worn wheels. Installing or using a wheel lock on a wheel with less than a 6" diameter could compromise the holding ability of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the cot or other equipment.

### CAUTION

Wheel locks are only intended to help prevent the cot from rolling while unattended and to aid in patient transfer. A wheel lock may not provide sufficient resistance on all surfaces or under loads.

[Return To Table of Contents](#)

# Operation Guide

English

## OPERATING THE OPTIONAL STEER-LOCK

### To activate steer-lock from the cot foot or head end:

- From the cot foot end, press the red (lock) side of the foot pedal as shown in Figure 58 or from the cot head end, press down on either red pedal as shown in Figure 59.
- Rotate the cot until at least one head end caster is locked.

Figure 58

Figure 59

### To deactivate steer-lock from the cot foot or head end:

- From the cot foot end, press the green (unlock) side of the foot pedal as shown in Figure 60 or from the cot head end, lift up on either red pedal at the head end as shown in Figure 61.

Figure 60

Figure 61

# Operation Guide

## USING THE OPTIONAL KICKSTAND FOR DIALYSIS SCALE

The kickstand is intended for weighing patients on a scale.

### Note:

- The optional kickstand assembly is configured for an X-frame cot retention system only.
- The optional kickstand (6085-102-000) is not compatible with the optional base storage net (6500-160-000).

### WARNING

- Stryker recommends a two person operation when using the kickstand.
- Make sure that the patient weight is centered on the cot before using the kickstand.
- Engage the kickstand with your foot only.
- Lower cot height prior to engaging kickstand for increased stability.
- Make sure that the kickstand remains in the retracted position and does not engage during transport.
- Do not use the kickstand as a brake.
- Do not engage kickstand on a sloped surface.

### To use the kickstand:

1. Operator 1 engages the kickstand with their foot as shown in Figure 62.
2. Operator 2 lifts the foot end of the cot at a height sufficient to actuate the kickstand.
3. Both operators must make sure that the kickstand is in the forward locked position as shown in Figure 63.

Figure 62

Figure 63

### To release the kickstand:

1. Operator 1 lifts the foot end of the cot until both wheels are off of the floor.
2. Operator 2 rolls the cot slightly forward to make sure that the kickstand retracts on its own as shown in Figure 64.


Figure 64

# Operation Guide

English

## USING RESTRAINT STRAPS

### WARNING

Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.

Always secure the patient on the cot with all of the restraint straps.

To attach the restraint strap to the cot:

1. Wrap the restraint strap around the cot frame as shown in Figure 65.
2. Push the restraint strap buckle through the loop as shown in Figure 66.
3. Pull the buckle through the loop to secure the restraint strap to the cot as shown in Figure 67.
4. Repeat steps 1-3 until all restraint straps are securely attached to the cot in the required attachment locations as shown in Figure 68.

Figure 65

Figure 66

Figure 67

Buckle the restraints across the patient's chest/shoulders, waist and legs. Note the attachment locations in Figure 68.

Keep the restraint straps buckled (as shown in Figure 68) when the cot is not being used with a patient to avoid damage to the buckles and straps.

**Note:** Restraint straps are a Type B applied part.

Figure 68

# Operation Guide

## USING RESTRAINT STRAPS (CONTINUED)

When attaching the restraint straps to the cot, the attachment points should provide both strong anchorage and proper restraint position while not interfering with equipment and accessories.

**⚠ WARNING**

Do not attach restraints to the base tubes, cross tubes, or fowler skin. Improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.

(Back View)

Shoulder/Chest Restraint Straps

Knee Restraint Straps

Foot Restraint Straps


Figure 69

# Operation Guide

English

## USING X-RESTRAINT STRAPS

### WARNING

Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot.

To attach the restraint strap at the shoulder, thigh, or ankle:

1. Wrap the restraint strap around the cot frame (Figure A).
2. Push the restraint strap buckle through the loop (Figure B).
3. Pull the buckle through the loop to secure the restraint strap to the cot (Figure C).

Figure A: Wrap strap around cot

Figure B: Push buckle through loop

Figure C: Pull strap to tighten

To attach the restraint strap at the waist:

1. Wrap the restraint strap around the cot frame.
2. Push both buckles through the loop for the double buckle strap (Figure D).
3. Push the buckle and the tongue through the loop for the tongue and buckle strap (Figure E).
4. Pull the buckle through the loop to secure the restraint strap to the cot.

Figure D

Figure E

**Note:** The restraint strap at the waist must form an X with the restraint strap at the shoulder (Figure F).

Figure F: Front view

Figure G: Back view

[Return To Table of Contents](#)

# Operation Guide

## USING X-RESTRAINT STRAPS (CONTINUED)

Secure restraint straps to the cot in the required attachment locations (Figure H). Restraint strap attachment locations should provide strong anchorage and proper restraint position. Do not allow restraint straps to interfere with equipment or accessories. Buckle the restraints across the patient's shoulders, waist, thighs, and ankles. Buckle restraint straps when the cot is not in use.

### WARNING

Do not attach restraint straps to the base tubes or cross tubes.


Figure H: Attachment locations

# Operation Guide

English

## ADJUSTING RESTRAINT STRAPS

### ⚠ CAUTION

Ensure that the restraints are not entangled in the base frame when raising and lowering the cot.

When the cot is put into service, open the restraints and place them at either side of the cot until the patient is positioned on the cot mattress. Lengthen the restraints, buckle them around the patient and shorten them until the required tightness is achieved.

- **To open the restraint**, press the red button (A) on the front of the buckle “receiver”. This releases the buckle latch plate (B) which can then be pulled out of the receiver (Figure 70).
- **To close the restraint**, push the latch plate into the receiver until a “click” is heard. For G-restraint straps only, ensure that the latch plate passes through both links (C) on the shoulder strap when fastening the chest restraint (Figure 70).
- **To lengthen the restraint**, grasp the buckle latch plate, turn it at an angle to the webbing, then pull it out (Figure 71). A hemmed tab at the end of the webbing prevents the latch plate from coming off of the strap.
- **To shorten the restraint**, grasp the hemmed tab and pull the webbing back through the latch plate until the required tightness is achieved (Figure 72).


Figure 70

Figure 71


Figure 72

Whenever a restraint is buckled on a patient, verify that the latch plate is fully engaged and any extra webbing is not tangled in the cot or hanging loose.

Inspect the restraints **at least** once a month (more frequently if used heavily). Inspection should include checking for a bent or broken receiver or latch plate, torn or frayed webbing, etc. Any restraint showing wear or not operating properly **must** be replaced immediately.

# Operation Guide

---

## USING THE RESTRAINT EXTENDER

Use the restraint extender, as shown in Figure 73, for extra length when buckling the lap belt around large patients.


Figure 73

# Operation Guide

English

## ATTACHING THE PEDI-MATE® INFANT RESTRAINT SYSTEM

See the Pedi-Mate® users manual for the manufacturer's recommendations for the use, operation and care of the Pedi-Mate® Infant Restraint System.

**Note:** The Pedi-Mate® Infant Restraint System is a Type B applied part.

### To secure the Pedi-Mate® to the cot:

1. Remove any restraints that are already attached to the cot.
2. Raise the cot backrest to the full upright position.
3. Position the Pedi-Mate® pad flat on the backrest with the black backrest straps out (see Figure 74).

Figure 74

4. Wrap the straps around the backrest and insert the ends of the straps through the brackets. Securely fasten the buckle (see Figure 75).


Figure 75

### **WARNING**

To avoid accidental release of the Pedi-Mate®, and possible injury to the infant, ensure that the buckle is located away from obstructions on the cot or accessories.

[Return To Table of Contents](#)

# Operation Guide

## ATTACHING THE PEDI-MATE® INFANT RESTRAINT SYSTEM (CONTINUED)

5. Pull firmly on the end of the adjustable backrest strap and tighten it securely.
6. Insert the mainframe straps between the cot frame and the mattress. To ensure that the release button is located toward the foot end of the cot, insert the buckle behind the litter cross brace and bring it up in front of the cross brace. Secure the buckle around the cross brace, leaving a little slack in the strap for final adjustment (see Figure 76).

Figure 76

### **WARNING**

To avoid accidental release of the Pedi-Mate®, and possible injury to the infant, ensure that the buckle is located away from obstructions on the cot or accessories.

7. Verify that all of the straps are snug and fastened securely (see Figure 77).

Figure 77

**Note:** These are general instructions for installation of the Pedi-Mate®. Safe and proper use of the Pedi-Mate® is solely at the discretion of the user. Stryker recommends that all users be trained on the proper use of the Pedi-Mate® before using it in an actual situation. Retain these instructions for future reference. Include them with the product in the event of transfer to new users.

Pedi-Mate® is a registered trademark of Ferno-Washington, Inc.

[Return To Table of Contents](#)

# Operation Guide

English

## REMOVING AND REPLACING A SMRT™ PAK

The cot is supplied with two removable 24V **SMRT™** Paks as the power source.

See the **SMRT™** Power System Operations/Maintenance Manual for additional **SMRT™** Pak and **SMRT™** Charger information.

### **WARNING**

- To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
- Do not remove the battery when the cot is activated.
- Avoid direct contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.

### To remove the **SMRT™** Pak:

1. Press the red one hand release button (C) or press the battery release button (A) to release the **SMRT™** Pak (B) from the cot as shown in Figure 78.
2. Slide the released **SMRT™** Pak out of the enclosure.


Figure 78

### To reinstall or replace the **SMRT™** Pak:

1. Align the tabs in the battery enclosure.
2. Push the **SMRT™** Pak into the enclosure until the latch clicks into place.
  - The cot power indicator LED is solid green if the **SMRT™** Pak is charged and ready.
  - The cot power indicator LED flashes amber if the **SMRT™** Pak needs to be recharged or replaced.

**Note:** Batteries slowly lose power when not on the charger.

### **CAUTION**

Remove the battery if the cot is not going to be used for an extended period of time (more than 24 hours).

## USING THE DEFIBRILLATOR PLATFORM

### To install the defibrillator platform:

1. Place the defibrillator platform in the stored position as shown in Figure 79.
2. Open and expand the defibrillator platform legs as shown in Figure 80.


Figure 79


Figure 80

3. If the cot is equipped with an I.V. pole, raise the I.V. pole (A) to the up position as shown in Figure 81.
4. Place the defibrillator platform on the cot frame as shown in Figure 81. Position the inside legs (B) of the defibrillator platform toward the head end of the cot to properly fit the defibrillator platform on the cot frame.


Figure 81

# Operation Guide

English

## USING THE DEFIBRILLATOR PLATFORM (CONTINUED)

- Place the latch hook (C) under the cot litter frame or foot end fastener and push the tab (D) up until it locks into place with an audible click as shown in Figure 82. Repeat on the other side of the defibrillator platform.

### Notes:

- For Power-LOAD compatible cots, you must lengthen and attach the straps to the foot end fastener as indicated in Figure 83.
- If the defibrillator platform is not securely attached to the cot when both latch hooks are engaged or if you cannot secure the latch hooks around the cot litter frame, unlatch the tab (E), loosen or tighten the strap (F) until the proper adjustment is achieved, and then push the tab (E) up until it locks into place with an audible click (Figure 84).


Figure 82

### WARNING

To avoid the risk of patient injury or equipment damage, ensure that you properly mount and secure the defibrillator platform to the cot.

- Ensure that the defibrillator platform is properly secured to the cot.
- Place the defibrillator on the defibrillator platform and secure the straps.


Figure 83

### WARNING

- To avoid the risk of patient injury or equipment damage, you must use the provided straps to secure the defibrillator to the defibrillator platform.
- Due to the difference in sizes and shapes of available defibrillators, you may have to change the location and adjustment of the straps that secure the defibrillator to the defibrillator platform. To avoid the risk of patient injury or equipment damage, use and adjust all straps properly to ensure the security of the defibrillator.
- To avoid the risk of patient injury or equipment damage, the weight placed on the defibrillator platform must not exceed 30 lb (13.6 kg).


Figure 84

# Operation Guide

English

## USING THE EQUIPMENT HOOK

Use the equipment hook (A) (see Figure 85) to hang additional accessories or equipment, such as defibrillators or monitors.

### CAUTION

To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 35 lb (15.8 kg).


Figure 85

## USING THE HEAD EXTENSION WITH PILLOW

You can install the head extension assembly over the fowler to provide head end support.

Attach the pillow to the head extension by placing the support into the flap on the bottom of the pillow. Secure the pillow to the Velcro® on the bottom of the support.

**Note:** The head extension with pillow (6100-044-000) is not compatible with the optional equipment hook (6500-147-000) or optional fowler oxygen bottle holder (6500-241-000).

**Note:** The head extension pillow is a Type B applied part.

[Return To Table of Contents](#)

# Operation Guide

English

## OPERATING THE OPTIONAL TWO-STAGE I.V. POLE

To use the two-stage I.V. pole (see Figure 87):

1. Lift and pivot the pole from the storage position and push down until it is locked into the receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the telescoping portion in place.
4. Hang the I.V. bags on the I.V. hook (D).
5. Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
6. Turn the lock actuator (B) clockwise to tighten.
7. Lift up and pivot the pole down into the storage position (see Figure 86).

### CAUTION

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 25 lb (11.3 kg).

**Note:** The dual two-stage I.V. poles (6500-312-000) are not compatible with either the patient right (6500-310-000) or patient left (6500-311-000) two-stage I.V. pole options.


Figure 86


Figure 87

# Operation Guide

## OPERATING THE OPTIONAL THREE-STAGE I.V. POLE

To use the three-stage I.V. pole (see Figure 89):

1. Lift and pivot the pole from the storage position and push down until it is locked into the receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the bottom telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the bottom telescoping portion in place.
4. For a higher I.V. pole, pull up on section (D) until the spring clip (E) engages.
5. Hang I.V. bags on the I.V. hook (F).
6. To lower the I.V. pole, push in on the spring clip (E) and slide section (D) down into section (C). Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
7. Turn the lock actuator (B) clockwise to tighten.
8. Lift up and pivot the pole down into the storage position (Figure 88).

### CAUTION

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 25 lb (11.3 kg).

**Note:** The dual three-stage I.V. poles (6500-317-000) are not compatible with either the patient right (6500-315-000) or patient left (6500-316-000) two-stage I.V. pole options.


Figure 88


Figure 89

# Operation Guide

---

English

## ATTACHING AN OXYGEN BOTTLE TO AN OXYGEN BOTTLE HOLDER

### To attach an oxygen bottle:

1. Place an oxygen bottle in the holder.
2. Insert the lower strap through the buckle and affix the strap onto itself to secure the oxygen bottle to the holder.

**Note:** Inspect the straps and clips for wear between use and replace the strap if it is no longer holding the oxygen bottle.

---

 **CAUTION**

- To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 15 lb (6.8 kg).
- Do not use two head end oxygen bottle holders at the same time.

---

**Note:** The optional fowler oxygen bottle holder (6500-241-000) is not compatible with the optional retractable head section oxygen bottle holder (6085-046-000).

# Operation Guide

## USING THE RETRACTABLE HEAD SECTION OXYGEN BOTTLE HOLDER

To attach an oxygen bottle to the retractable head section oxygen bottle holder:

1. Center the oxygen bottle on the cradled surface of item (A) as shown in Figure 90.
2. Tighten both straps (B) around the oxygen bottle.
3. Secure the slack on the straps to the Velcro® on the straps.


Figure 90

**Note:** Inspect the straps and clips for wear between use and replace the strap if it is no longer holding the oxygen bottle.

### **WARNING**

If the cot is equipped with the optional retractable head section oxygen bottle holder, use caution while the oxygen bottle holder is installed to avoid pinching your fingers between the fowler bracket and the oxygen bottle.

### **CAUTION**

- To avoid damage to the oxygen bottle holder (if equipped), the weight of the equipment must not exceed 15 lb (6.8 kg).
- Do not use two head end oxygen bottle holders at the same time.

# Operation Guide

English

## INSTALLING THE BASE STORAGE NET

To install the base storage net, wrap the Velcro® straps around the base tubes.

### ⚠ CAUTION

- The weight of the equipment in the base storage net (if equipped) must not exceed 20 lb (9 kg).
- Be careful when retracting the base to avoid damaging items stored in the base storage net.

**Note:** The optional kickstand (6085-102-000) is not compatible with the optional base storage net (6500-160-000).

# Operation Guide

## INSTALLING THE BACKREST STORAGE POUCH

Install the optional backrest storage pouch using the Velcro® straps as shown in Figure 91. Insert each strap through a hole in the backrest skin and mount the pouch flat against the backrest.

### CAUTION

- Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.
- The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 20 lb (9 kg).


Figure 91

# Operation Guide

English

## INSTALLING THE HEAD END STORAGE FLAT

### WARNING

When the optional head end storage flat is being used, ensure that it does not interfere with the operation of the retractable head section, safety bar and safety hook. Injury to the patient or operator could result.

**To install the optional head end storage flat (see Figure 92):**

1. Install the Velcro® straps (A) near the pneumatic cylinder and around the cross bar of the retractable head section.
2. Buckle the restraint straps (B) around the outer rails of the retractable head section.


Figure 92

### CAUTION

The weight of the equipment in the head end storage flat (if equipped) must not exceed 40 lb (18 kg).

## USING THE TRANSFER FLAT

When transferring larger patients, use of the Transfer Flat (6005-001-001) is recommended.

# Operation Guide

---

English

## ATTACHING THE MATTRESS

Use the standard gatch bolster mattress (6500-002-150/6506-002-150) with the standard siderail. Use the wider gatch bolster mattress (6500-003-130) with the optional eXpandable patient system (XPS).

### To attach the mattress to the cot:

1. Align the Velcro® on the back of the mattress with the Velcro® on the cot litter.
2. Attach the strap at the foot end of the mattress through the two holes in the foot end skin on the cot litter.
3. Pull the strap through the buckle and attach the Velcro® to secure the strap.

**Note:** The optional gatch bolster mattress for XPS (6500-003-130) is not compatible with the standard siderail (6506-031-000).

**Note:** The mattress is a Type B applied part.

---

### CAUTION

Do not store items under the cot mattress. Storing items under the mattress can interfere with the operation of the cot.

---

# Cleaning

---

English

The **Power-PRO™** XT cot is designed to be power washable. The unit may show some signs of oxidation or discoloration from continuous washing, however, no degradation of the cot's performance characteristics or functionality will occur due to power washing as long as the proper procedures are followed.

Thoroughly clean the cot once a month. Clean Velcro® **AFTER EACH USE**. Saturate Velcro® with disinfectant and allow disinfectant to evaporate. Appropriate disinfectant for nylon Velcro® should be determined by the service.

## WASHING PROCEDURE

- **Always remove the battery!** Never wash the cot with the battery installed.
- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- The preferred method Stryker Medical recommends for power washing the cot is with the standard hospital surgical cart washer or hand held wand unit.

## WASHING LIMITATIONS

---

### **WARNING**

When cleaning, use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the cot.

---

### **CAUTION**

- **DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT.**
  - Maximum water temperature should not exceed 180°F/82°C.
  - Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 24 inches (61 cm) from the unit.
  - Allow cot to air dry.
  - Towel dry all casters and interface points.
  - Failure to comply with these instructions may invalidate any/all warranties.
  - Always remove the battery before washing the cot.
-

# Cleaning

---

In general, when used in those concentrations recommended by the manufacturer, either phenolic type or quaternary (**excluding Virex® TB**) type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may result.

Suggested cleaners for the cot surfaces and restraint straps:

- Quaternary Cleaners (active ingredient - ammonium chloride)
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)

**Note:** Do not immerse the restraint strap metal buckle components in water. Rinse with clean water. Allow to air dry.

Avoid over saturation and ensure that the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

---

## WARNING

SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products described above are used to clean Stryker EMS equipment, measures must be taken to ensure the cots are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the cots will leave a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.

**Note:** Failure to follow the above directions when using these types of cleaners may void this product's warranty (see [page 1-117](#)).

## REMOVAL OF IODINE COMPOUNDS

Use a solution of 1/2 Tablespoon Sodium Thiosulfate in a pint of warm water to clean the stained area. Clean as soon as possible after staining occurs. If stains are not immediately removed, allow solution to soak or stand on the surface. Rinse surfaces which have been exposed to the solution in clear water before returning unit to service.

---

## WARNING

Failure to properly clean or dispose of contaminated mattress or cot components will increase the risk of exposure to bloodborne pathogens and may cause injury to the patient or the operator.

---

# Preventive Maintenance

---

English

A preventive maintenance program should be established for all Stryker Medical equipment. Preventive maintenance may need to be performed more frequently based on the usage level of the product. The cot requires regular maintenance. Establish and follow a maintenance schedule and keep records of maintenance activity (see [page 1-84](#) for a form).

---

## CAUTION

A preventive maintenance program should be established for all Stryker EMS equipment. Preventive maintenance may need to be performed more frequently based on the usage level of the product. Close attention should be given to safety features including, but not limited to:

- Hydraulic power mechanism
- All electrical controls return to off or neutral position when released

For additional maintenance information, see the preventive maintenance information on [page 1-84](#).

---

## WARNING

- Do not modify the cot or any components of the cot. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product will also void its warranty (see [page 1-117](#)).
  - Escaping fluid under pressure can penetrate the skin causing serious injury. Avoid the hazard by relieving pressure before disconnecting hydraulic or other lines. Tighten all connections before applying pressure. If an accident occurs, see a doctor immediately. Any fluid injected into the skin must be surgically removed within a few hours or gangrene may result. Doctors unfamiliar with this type of injury should reference a knowledgeable medical source.
  - To avoid the risk of injury, do not use bare hands to check for hydraulic leaks.
- 

When using maintenance products, follow the directions of the manufacturer and reference all material safety data sheets.

---

## CAUTION

- Improper maintenance can cause injury or damage to the product. Maintain the cot as described in this manual. Use only Stryker approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty (see [page 1-117](#)).
  - Failure to use authorized parts, lubricants, etc. could cause damage to the cot and will void the warranty of the product.
  - Hydraulic lines, hoses, and connections can fail or loosen due to physical damage, kinks, age, and environment exposure. Check hoses and lines regularly to avoid damage to the cot. Check and tighten loose connections.
  - Do not tip the cot onto its load wheels and actuate the product as this will allow air to enter the hydraulic system.
- 

## LUBRICATION

The cot has been designed to operate without the need for lubrication.

---

## CAUTION

Do not lubricate the bearings in the X-frame as it will degrade the performance of the cot and may void its warranty (see [page 1-117](#)).

---

You can apply Tri-Flow® lubrication to the optional kickstand during periodic maintenance (see [page 1-81](#)).

---

# Preventive Maintenance

## REGULAR INSPECTION AND ADJUSTMENTS

### Maintenance Intervals

The following schedule is intended as a general guide to maintenance. Bear in mind that such factors as weather, terrain, geographical location, and individual usage will alter the required maintenance schedule. If you are unsure as to how to perform these checks please contact your Stryker service technician. If you are in doubt as to what intervals to follow in maintaining your product, consult your Stryker service technician. Use the hour meter ([page 1-32](#)) to determine the frequency for preventive maintenance procedures. Check each routine and replace damaged or worn parts if necessary.

Item	Routine	Every (whichever comes first)			
		1 Month or 2 hours	3 Months or 6 hours	6 Months or 12 hours	12 Months or 24 hours
Settings	Verify the in-fastener shut-off is configured properly	X			
	Verify the cot and fastener fit and function properly				X
	Verify the safety bar engages the vehicle safety hook properly				X
Cylinder	All fasteners are secure (reference all assembly drawings)		X		
	Verify the cylinder is adjusted so the lock nut is tight and the cot stops moving when it hits the dead stops				X
	Inspect for and verify that there are no hydraulic fluid (red) leaks; inspect the fittings and tighten as necessary		X		
	Extend cylinder rod completely and wipe down rod with soft cloth and household cleaner	X			
Hydraulics	Inspect motor mount and verify that all fasteners are secure		X		
	Verify that there are no hydraulic fluid leaks		X		
	Inspect the reservoir and verify that there are no leaks		X		
	Inspect hoses and fittings for damage or wear			X	
	Verify the hydraulic velocity fuse - Place a weight of approximately 50 lb on the cot, raise the cot, lift the cot with two operators, pull the manual back-up release handle, rapidly set the cot down, verify that the cot does not drop				X
Electronic Controls	Extend cot to raised position, measure and check load height			X	
	Verify "jog" function is operating			X	
	Verify high speed retract is working			X	
Switches	Verify there is no damage or wear to either switch			X	
	Verify both switches operate correctly			X	
Cables/Wires	Verify there is no damage or pinching of wiring harness, cables or lines		X		
	Check routings and connections, verify there are no hanging wires	X			
	Verify there are no damaged connectors		X		

[Return To Table of Contents](#)

# Preventive Maintenance

English

Item	Routine	Every (whichever comes first)					
		1 Month or 2 hours	3 Months or 6 hours	6 Months or 12 hours	12 Months or 24 hours		
Manual Back-up Release Handle	Verify that the manual back-up release handle functions properly	X					
	Verify the manual back-up release handle returns to the stowed position				X		
	Verify the base extends/retracts smoothly when the manual back-up release handle is engaged		X				
	With 100 lb or more on the cot, verify the cot does not lower when the manual back-up release handle is pulled		X				
	Inspect the cot frame/litter	X					
	Verify all welds intact, not cracked or broken				X		
	Verify no bent, broken or damaged components			X			
	Verify all fasteners secure (reference all assembly drawings)		X				
	Verify warning labels present, legible (reference assembly drawings)				X		
	Verify no damage or tears on cot grips			X			
Litter	Verify the siderails operate and latch properly			X			
	Verify the backrest cylinder operates properly		X				
	Adjust pneumatic cylinder for full range of motion, if required		X				
	Verify the footrest operates properly			X			
	Verify no cracks or tears on cot mattress			X			
	Mattress	Inspect patient restraints for proper function and no excessive wear (bent or broken receiver or latch plate, torn or frayed webbing, etc.)	X				
		Base	Inspect the cot frame/base	X			
			Verify all welds intact, not cracked or broken				X
			Verify no bent, broken, or damaged components			X	
	Wheels	Verify all fasteners secure		X			
Verify that the cot retaining post is secure. If not secure, then the screw must be replaced. See "Cot Retaining Post Screw Replacement" on page 1-106.				X			
Verify no excessive damage to X-frame guards				X			
Verify wheels are free of debris				X			
X-Frame	Verify all wheels secure, rolling and swiveling properly	X					
	Check and adjust optional wheel locks as necessary				X		
	Verify smooth operation of X-frame		X				

# Preventive Maintenance

Item	Routine	Every (whichever comes first)			
		1 Month or 2 hours	3 Months or 6 hours	6 Months or 12 hours	12 Months or 24 hours
Head Section	Verify all fasteners secure		X		
	Verify no bent, broken, or damaged components			X	
	Verify the head section extends and locks properly		X		
	Verify the grip bar has no excessive damage or tears			X	
	Verify load wheels are secure and roll properly			X	
	Verify the safety bar operates properly. Pull toward the head section to ensure that it swings and rotates freely and pulls back to home position.	X			
Battery	Inspect the <b>SMART™</b> Pak housing and terminal area for cracks or damage	X			
Accessories	Inspect the straps and clips on the retractable head section oxygen bottle holder (optional) for wear				X
	Verify the I.V. pole (optional) operates properly		X		
	Verify the head extension & pillow (optional) operates properly		X		
	Verify the restraint extender (optional) operates properly		X		
	Verify that the defibrillator platform straps (optional) are intact - not frayed or torn				X
	Verify that both defibrillator platform (optional) latch hooks are intact and secure				X
	Verify the kickstand (optional) retracts fully to the transport position		X		
	Verify that the kickstand (optional) bolts are tightened properly		X		
	Lubricate the kickstand spring and internal spring housing (optional) using Tri-Flow® lubrication			X	


# Troubleshooting Guide

English

## ELECTRONICS AND HYDRAULICS LOCATOR

**Note:** Some components have been removed for clarity.

**Figure 93**

# Troubleshooting Guide

## HYDRAULIC ASSEMBLY



Figure 94

## HYDRAULIC ASSEMBLY WIRING SCHEMATICS


[Return To Table of Contents](#)

# Troubleshooting Guide

English

## ELECTRICAL SYSTEM BLOCK DIAGRAM

Lift and Extend (Unload) Functions


# Troubleshooting Guide

## Electrical System Block Diagram Lower and Retract (Load) Functions

English

[Return To Table of Contents](#)

# Troubleshooting Guide

English

## TROUBLESHOOTING GUIDE

Check for proper operation after each step. When the problem is fixed, return the cot to service. If assistance is needed at any time during troubleshooting, please contact a service technician at (800) 327-0770 or (269) 324-6500.

PROBLEM	SOLUTION	PAGES
Litter drifts (without patient weight)	<ol style="list-style-type: none"> <li>1. Flush the hydraulic system by squeezing the manual release handle while simultaneously pressing the (+) power button for approximately 15 seconds. Repeat if necessary.</li> <li>2. Check the manual release cable adjustment.</li> <li>3. Change the 'locking' manual valve.</li> <li>4. Change the 'B' valve.</li> </ol>	
Base drifts (without patient weight)	<ol style="list-style-type: none"> <li>1. Flush the hydraulic system by squeezing the manual release handle while simultaneously pressing the (+) power button for approximately 15 seconds. Repeat if necessary.</li> <li>2. Check the manual release cable adjustment.</li> <li>3. Change the 'non-locking' manual valve.</li> <li>4. Change the 'A' valve.</li> </ol>	
Litter does not lower in the powered mode	<ol style="list-style-type: none"> <li>1. Check the power indicator LED.               <ol style="list-style-type: none"> <li>a. If blinking constant amber, change the battery.</li> </ol> </li> <li>2. Check for error on LCD.</li> <li>3. Check for broken or disconnected wires.</li> <li>4. Check for 24V DC at connector (C) on the main cable by the motor while pressing the retract (-) button. If voltage is present, replace (in order) the hall effects sensor, solenoid, and or 'B' valve. If voltage is not present, go to step 5.</li> <li>5. Check for 24V DC on electronics assembly pins 1 blue and 5 orange on (F) while pressing the retract (-) button. If voltage is not present, replace the electronics assembly. If voltage is present, replace the wire harness.               <ol style="list-style-type: none"> <li>a. If the green light turns on, but does not lower, try the other switch. If the other switch works, replace the bad switch.</li> </ol> </li> </ol>	<p><a href="#">page 1-94</a>  <a href="#">page 1-95</a></p>

# Troubleshooting Guide

## TROUBLESHOOTING GUIDE (CONTINUED)

PROBLEM	SOLUTION	PAGES
<p>Litter does not extend in the powered mode</p>	<ol style="list-style-type: none"> <li>1. Check the power indicator LED.               <ol style="list-style-type: none"> <li>a. If blinking constant amber, change the battery.</li> </ol> </li> <li>2. Check for error on LCD.</li> <li>3. Check for broken or disconnected wires.</li> <li>4. Check for 24V DC at connector (C) on the main cable by the motor while pressing the retract (-) button. If voltage is present, replace (in order) the hall effects sensor, solenoid, and or 'B' valve. If voltage is not present, go to step 5.</li> <li>5. Check for 24V DC on electronics assembly pins 1 blue and 5 orange on (F) while pressing the retract (-) button. If voltage is not present, replace the electronics assembly. If voltage is present, replace the wire harness.               <ol style="list-style-type: none"> <li>a. If the green light turns on, but does not lower, try the other switch. If the other switch works, replace the bad switch.</li> </ol> </li> <li>6. Check the motor.               <ol style="list-style-type: none"> <li>a. If the motor runs, but does not raise the cot:                   <ol style="list-style-type: none"> <li>i. Check the manual release cable for too much tension.</li> <li>ii. Lightly tap the manual locking valve.</li> <li>iii. Replace the manual locking valve.</li> </ol> </li> <li>b. If the motor is stalled, replace the 'A' valve.</li> <li>c. If the light is green, but the motor does not run:                   <ol style="list-style-type: none"> <li>i. Check for 24V DC at connector (E) on the main cable. If voltage is present, replace the hall effects sensor. If the hall effects sensor is replaced, and the motor still does not run, replace the hydraulic sub assembly. If voltage is not present, go to step ii.</li> <li>ii. Check for 24V DC on electronics assembly connection (H) (-) lead on black (+) lead on green while pressing the extend (+) button. If voltage is not present, replace the electronics assembly. If voltage is present, replace the main cable.</li> </ol> </li> </ol> </li> </ol>	<p><a href="#">page 1-94</a>  <a href="#">page 1-95</a></p>

# Troubleshooting Guide

English

## TROUBLESHOOTING GUIDE (CONTINUED)

PROBLEM	SOLUTION	PAGES
Base does not retract in the powered mode	<ol style="list-style-type: none"> <li>1. Check the power indicator LED.               <ol style="list-style-type: none"> <li>a. If blinking constant amber, change the battery.</li> </ol> </li> <li>2. Check for error on LCD.</li> <li>3. Check for broken or disconnected wires.</li> <li>4. Check for 24V DC at connector (C) on the main cable by the motor while pressing the retract (-) button. If voltage is present, replace (in order) the hall effects sensor, solenoid, and or 'B' valve. If voltage is not present, go to step 5.</li> <li>5. Check for 24V DC on electronics assembly pins 1 blue and 5 orange on (F) while pressing the retract (-) button. If voltage is not present, replace the electronics assembly. If voltage is present, replace the wire harness.</li> </ol>	<a href="#">page 1-94</a> <a href="#">page 1-95</a>
Base does not extend in the manual mode	<ol style="list-style-type: none"> <li>1. Check the manual cable adjustment.</li> <li>2. Change the 'non-locking' manual valve.</li> </ol>	
Base does not retract in the manual mode	<ol style="list-style-type: none"> <li>1. Check the manual release cable adjustment.</li> <li>2. Change 'locking' manual valve.</li> </ol>	
Litter does not retract in the manual mode (with patient weight)	<ol style="list-style-type: none"> <li>1. Make sure that the weight is off of the casters before lowering the cot.</li> <li>2. Check the manual cable adjustment.</li> <li>3. Replace the 'locking' manual valve.</li> </ol>	
Litter does not extend in the manual mode	<ol style="list-style-type: none"> <li>1. Check the manual cable adjustment.</li> <li>2. Change the 'non-locking' manual valve.</li> </ol>	
High speed retract does not engage	<ol style="list-style-type: none"> <li>1. Check that the weight is off of the casters.</li> <li>2. Change the pressure switch.</li> <li>3. Change the hall effect cable.</li> </ol>	

# Troubleshooting Guide

## LCD ERROR CODES

### MAIN MICRO

LCD DISPLAY	ERROR DESCRIPTION	DETECTION PERIOD
ERR 01	RAM diagnostic failure	Initialization
ERR 02	Program memory failure	Initialization
ERR 03	EE diagnostic failure	Initialization
ERR 04	EEPROM type and hardware type incompatible	Initialization
ERR 10	Valves diagnostic failure	Initialization
ERR 61	EEPROM rev and firmware rev incompatible	Initialization
ERR 21	Motor shorted	Initialization
ERR 22	Motor open	Initialization
ERR 23	High power gating relay shorted	Initialization
ERR 51	Motor drive FET shorted - Q15	Initialization
ERR 52	Motor drive FET shorted - Q11	Initialization
ERR 55	Motor drive FET shorted - Q16	Initialization
ERR 56	Motor drive FET shorted - Q12	Initialization
ERR 62	Main Micro and ASIC current limit mismatch	Initialization
ERR 80	Extend (+) or retract (-) button detected without key	Run Time
ERR 31	Electronics board over temp (280.22 °F +/- 5%)	Run Time
ERR 81	Bad hall effect sensor combination	Run Time
ERR 93	Safety Micro non-responsive	Run Time

### SAFETY MICRO

LCD DISPLAY	ERROR DESCRIPTION	DETECTION PERIOD
ERR 05	RAM diagnostic failure	Initialization
ERR 06	Program memory diagnostic failure	Initialization
ERR 08	EEPROM type and hardware type incompatible	Initialization
ERR 40	Data error	Run Time
ERR 41	Charging failed battery voltage	Run Time
ERR 42	Charging failed read battery	Run Time
ERR 43	Charging failed battery charging time or over voltage limit	Run Time
ERR 44	Charging failed charging current	Run Time
ERR 45	Charging failed delta temp	Run Time
ERR 63	EEPROM rev and firmware rev incompatible	Initialization
ERR 83	Extend (+) or retract (-) button detected without key	Run Time
ERR 90	ASIC driving without microprocessor instruction	Run Time

# Troubleshooting Guide

English

## MAIN CABLE ASSEMBLY



Figure 95

## MAIN CABLE ASSEMBLY WIRING SCHEMATICS

# Troubleshooting Guide

## CONTROL BOARD ASSEMBLY



Figure 96

## CONTROL BOARD WIRING SCHEMATICS

ITEM	CONNECTION	CONNECTION	ITEM
------	------------	------------	------

MOTOR CABLE

BATTERY POSITIVE

BATTERY NEGATIVE

[Return To Table of Contents](#)

# Quick Reference Replacement Parts List

English

The parts listed on these pages are all currently available for purchase. Some of the parts identified on the assembly drawing in this manual may not be individually available for purchase. Please call Stryker Customer Service USA: 1-800-327-0770 for availability and pricing.

Part Name	Part Number
Cable, Hall Effect Sensor	6500-001-160
DC Battery Charger, 110V, Domestic	6500-070-000
DC Battery Charger 12V/24V, In-Ambulance	6500-072-000
Electronics Assembly	6500-002-014
Gas Cylinder, Backrest	1010-031-077
Headsection (not Power-LOAD compatible)	6506-001-020
Headsection (Power-LOAD compatible)	6506-001-021
Hydraulic Oil	6500-001-293
Kit, Base Tube Protector	6500-700-015
Kit, Battery Pack, <b>SMRT™</b> Pak	6500-700-046
Kit, <b>SMRT™</b> Power System 12V DC (Car Charger), includes charger, 2 paks, and power cord	6500-700-040
Kit, <b>SMRT™</b> Power System 120V AC (Wall Charger), includes charger, 2 paks, and power cord	6500-700-041
Kit, Brake Adjuster	6080-700-018
Kit, Brake - Single Wheel, Foot End Right	6080-201-000
Kit, Brake - Dual Wheel, Foot End Left and Right	6080-202-000
Kit, Equipment Hook	6500-700-003
Kit, Head Extension	6100-700-012
Kit, Lift Capable Safety Bar	6082-700-031
Kit, Oxygen Bottle Holder, Foot End	6500-700-017
Kit, Oxygen Bottle Holder, Fowler	6500-700-011
Kit, Oxygen Bottle Holder, Retractable Head Section	6085-700-003
Mounting Bracket, <b>SMRT™</b> Charger	6500-201-100
Restraint Belt Extension	6082-160-050
Restraint, Leg	6500-001-395
Restraint, Waist	6500-001-393
Restraint Package, Domestic	6500-002-030
Restraint, Shoulder Harness	6500-001-391
Restraint Strap Plastic Cap (Short)	6082-160-051
Restraint Strap Plastic Cap (Tall)	6082-160-055

# Quick Reference Replacement Parts List

English

Part Name	Part Number
Safety Hook, J	6092-036-018
Safety Hook, Long	6060-036-017
Safety Hook, Short	6060-036-018
Siderail, Standard	6082-026-010
Siderail, XPS Option (XPS Siderail and XPS Mattress Kit)	6506-700-004
Touch-Up Paint (Yellow)	6060-199-010
Touch-Up Paint (Black)	7000-001-322
Valve, "A"	6500-001-286
Valve, "B"	6500-001-287
Valve, Locking	6500-001-288
Valve, Non-Locking	6500-001-289
Velcro® 4.6" Adhesive Loop Pile, Litter	6060-032-046
Velcro® 10"	6082-001-082
Wheel Lock	6086-200-010

[Return To Table of Contents](#)

# Accessories

English

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Part Name	Part Number
Defibrillator platform	6506-170-000
Equipment hook	6500-147-000
Head extension with pillow	6100-044-000
I.V. pole, two-stage, right	6500-310-000
I.V. pole, two-stage, left	6500-311-000
I.V. pole, two-stage, dual	6500-312-000
I.V. pole, three-stage, right	6500-315-000
I.V. pole, three-stage, left	6500-316-000
I.V. pole, three-stage, dual	6500-317-000
Mattress, knee gatch bolster	6500-002-150
Mattress, knee gatch bolster, grey	6506-002-150
Mattress, knee gatch bolster, XPS	6500-003-130
Oxygen bottle holder, foot end	6500-240-000
Oxygen bottle holder, head end	6500-241-000
Oxygen bottle holder, removable	6080-140-000
Oxygen bottle holder, retractable head section	6085-046-000
Pedi-Mate restraint package	6091-300-010
Storage flat, head end	6500-128-000
Storage net, base	6500-160-000
Storage pouch, back rest	6500-130-000

## BACKREST ADJUSTMENT

### Tools Required:

- 1/2" Combination Wrench
- 5/32" Hex Wrench
- 3/32" Hex Wrench
- Small Slotted Screwdriver

### Procedure:

1. For easier access, move the backrest to 73 degrees.

**Note:** Before continuing with the backrest adjustment procedure, be sure that the cylinder (A) is completely threaded into the yoke (B) so no threads are showing on the shaft of the cylinder. If threads are showing, use a 3/32" hex wrench to remove the set screw (C) in the center of the yoke. Using a small slotted screwdriver, remove the E-clip and pin (D & E) that holds the bottom of the pneumatic cylinder. Thread the cylinder shaft (A) completely into the yoke (B). Replace the E-clip and pin (D & E) and replace the set screw (C) using Loctite®. (Figure 97 and Figure 98)

2. Using a 1/2" combination wrench, loosen the hex nut (F) on the backrest pivot (J) while holding the set screw (H) fixed in the pivot Figure 98).
3. Using a 5/32" hex wrench, turn the set screw (H) until there is no play between the backrest release handle (K) and the pneumatic cylinder release button (Figure 98).

**Note:** Make sure that the backrest travels from flat to at least 73 degrees. If it does not, turn the set screw clockwise 1/2 turn. Repeat until at least 73 degrees of travel is achieved.

4. Lower the backrest to a 5–10 degree angle and release the handle. Apply approximately 50 lb of downward force to the end of the backrest. If the backrest drifts down, turn the set screw counterclockwise. Repeat until the backrest does not drift downward.
5. Using the 1/2" combination wrench, tighten the hex nut (F) while holding the set screw fixed in the pivot (Figure 98).
6. Verify proper operation of the unit before returning it to service.


Figure 97


Figure 98

# Service Information

English

## HEADSECTION REPLACEMENT

### Tools Required:

- 7/16" Combination Wrench
- 3/16" Hex Wrench

### Procedure:

1. Raise the cot and the backrest to the full upright position.
2. Using a 7/16" combination wrench and a 3/16" hex wrench, remove the two screws (A) that secure the cap bearings to the base litter interface bracket (one on each side) (Figure 99).
3. Squeeze the head release handles and slowly remove the head section assembly.
4. Reverse steps to reinstall.
5. Verify proper operation of the unit before returning it to service.


Figure 99

## BACKREST GAS CYLINDER REPLACEMENT

### Tools Required:

- 3/32" Hex Wrench
- Slotted Screwdriver

### Procedure:

1. Raise the cot and the backrest to the full upright position.
2. Using a 3/32" hex wrench, loosen the set screw (A) that holds the gas shaft to the yoke (Figure 100).
3. Using a slotted screwdriver, remove the e-clip from the clevis pin that secures the bottom of the gas cylinder.
4. Unscrew the gas cylinder shaft from the yoke.
5. Reverse the above procedures to install the new gas cylinder. See "Backrest Adjustment" on [page 1-99](#).
6. Verify proper operation of the unit before returning it to service.


Figure 100

# Service Information

English


## MANUAL RELEASE CABLE ADJUSTMENT

### Tools Required:

- 8 mm Combination Wrench
- 10 mm Combination Wrench
- (2) Weight (50 lb each)

### Procedure:

1. Support the litter so no weight is on the base.
2. Ensure that the manual release cable is intact (A) (Figure 101).
3. Using a 10 mm combination wrench, loosen the cable lock nut (B) (Figure 101).
4. Using a 8 mm combination wrench, adjust the tension on the manual release cable so it just starts to touch the manual release dual pull bracket (C) (Figure 101).

**Note:** The manual release dual pull bracket should not be tight against the manual valve nuts.

5. Tighten the cable lock nut.
6. Test for proper adjustment by following steps A-D:
  - A. Place 50 lb of weight on the hydraulic skin.
  - B. Load height must read 34-1/2" to 35-1/2".
  - C. Place 100 lb of weight on the hydraulic skin, raise cot to full height, pull the manual release handle and ensure that the cot does not drop.
  - D. Remove 100 lb of weight, raise cot to full height, pull the manual release handle, and ensure that the cot drops.

**Note:** If steps A-D do not work properly, repeat steps 3-6.

7. Verify proper operation of the unit before returning it to service.

[Return To Table of Contents](#)

# Service Information

English

## FILLING THE HYDRAULICS ASSEMBLY RESERVOIR

Use only Mobil Mercon® V Synthetic Blend Oil (6500-001-293)

**Note:** Any time you work with the hydraulics you may lose some oil.

### Tools Required:

- 3/16" Hex Wrench

### Procedure:

1. Raise the cot to the full up position.
2. Ensure that the fill port is horizontal and lined up with the hole in the motor mount.
3. Remove the port plug (A) using a 3/16" hex wrench (Figure 102).
4. Fill the reservoir up to the bottom of the fill port.
5. Replace the plug and run the cot up and down a few times.
6. Verify proper operation of the unit before returning it to service.


Figure 102

## WHEEL LOCKING FORCE ADJUSTMENT

### Tools Required:

- 5/32" Hex Wrench
- 7/16" Combination Wrench or Socket

### Procedure:

1. Using the 5/32" hex wrench and 7/16" combination wrench or socket, remove the socket screw from the center of the lock pedal. The wheel lock is initially assembled with the pedal set at the minimum locking force. The marker on the pedal (A) is aligned with the marker on the octagonal sleeve (B) (Figure 103).
2. Remove the sleeve (B). Rotate the sleeve counterclockwise to increase the pedal locking force and clockwise to decrease the locking force. Insert the sleeve into the pedal (Figure 103).
3. Using the 5/32" hex wrench and 7/16" combination wrench or socket, reinstall the socket screw.
4. Test the pedal locking force and verify that the pedal holds properly before returning it to service.


Figure 103

# Service Information

English

## STEER-LOCK MECHANISM ADJUSTMENT

### Tools Required:

- 9/16" Combination Wrench

### Procedure:

#### If your steer-lock mechanism will not engage:

- Using a 9/16" combination wrench, adjust the barrel nuts toward the foot end of the cot (Figure 104).

**Note:** After adjustment, make sure that a minimum of one full thread is exposed on each side of the barrel nut.

#### If your steer-lock mechanism will not disengage:

- Using a 9/16" combination wrench, adjust the barrel nuts toward the head end of the cot (Figure 105).

**Note:** After adjustment, make sure that a minimum of one full thread is exposed on each side of the barrel nut.

Figure 104

Figure 105

# Service Information

English

## COT RETAINING POST ADJUSTMENT

### Tools Required:

- T30 Torx Driver

### CAUTION

The cot retaining post is shipped preconfigured for an X-frame cot. If the cot fastener has been configured for an H-frame cot, you must adjust the cot retaining post to accommodate the cot fastener.

### Procedure:

1. Using a T30 Torx driver, remove the two socket head cap screws (A) that hold the brackets (B) to the base frame (C) (Figure 106). Save both screws for reinstallation.
2. Turn the bottom bracket 180°.
3. Using a T30 Torx driver, reinstall the two socket head cap screws that were removed in step 1.
4. Verify proper operation of the unit before returning it to service.

To determine if your cot is an X-frame or H-frame cot, look for an arrow or groove on the bottom bracket of the cot retaining post.

- The cot retaining post is set for an X-frame cot if the arrow on the bottom bracket of the retaining post points toward the head end of the cot or if the groove in the bottom bracket is located on the inside of the patient left side of the base tube.
- The cot retaining post is set for an H-frame cot if the arrow on the bottom bracket of the retaining post points toward the foot end of the cot or if the groove in the bottom bracket is located on the outside of the patient left side of the base tube.


Figure 106

[Return To Table of Contents](#)

# Service Information

English

## COT RETAINING POST REPLACEMENT

### Tools Required:

- T30 Torx Driver
- 5/32" Hex Wrench
- Torque Wrench (in-lb)

### Procedure:

1. Raise the cot to the full upright position.
2. Turn the cot onto the patient left side.

**Note:** Locate the arrow or groove on the bottom bracket. The replacement retaining post bracket will need to be assembled in the same orientation.


Figure 107

3. Using a T30 Torx driver, remove the two socket head cap screws (A) that secure the current cot retaining post to the base tube (B) (Figure 107). Discard the screws and cot retaining post.
4. Insert the button head cap screw (C) through the retaining post cap (D) and post tube (E), and then into the top pin bracket (F) (Figure 107).
5. Using a 5/32" hex wrench, tighten the button head cap screw (C) completely to secure the retaining post cap (D) and post tube (E) to the top pin bracket (F) (Figure 107). Using a torque wrench, torque the screw to 100-140 in-lb.
6. Assemble the cot retaining post across the base tube. Align the holes of the brackets and insert two socket head cap screws (A) into the threaded holes of the bottom pin bracket (G) (Figure 107).
7. Using a T30 Torx driver, tighten the two socket head cap screws completely.
8. Verify proper operation of the unit before returning it to service.

**Note:** Adjustment of the rail clamp assembly may be required in order to compensate for any variation in cot retaining post position depending on the ambulance cot manufacturer and model number.

## COT RETAINING POST SCREW REPLACEMENT

### Tools Required:

- T25 Torx Driver
- 5/32" Hex Wrench
- Torque Wrench (in-lb)

### Procedure:

1. Using a T25 Torx driver or 5/32" hex wrench, remove the button head cap screw that secures the retaining post cap and post tube to the top pin bracket. Discard the screw.
2. Using a 5/32" hex wrench, install and tighten the button head cap screw (0004-503-000) completely to secure the retaining post cap and tube to the top portion of the retaining post assembly. Using a torque wrench, torque the screw to 100-140 in-lb.

**Note:** If you cannot torque the screw to 100-140 in-lb, then you must replace the entire cot retaining post. See "Cot Retaining Post Replacement".

3. Verify proper operation of the unit before returning it to service.

# Service Information

English

## HYDRAULIC A VALVE OR B VALVE REPLACEMENT

### Tools Required:

- T25 Torx Driver
- 3/4" Combination Wrench
- 7/8" Combination Wrench
- (2) Saw Horse

### Procedure:

1. Raise the cot to the full upright position.
2. Using two saw horses, support the litter and activate the manual release handle to relieve any hydraulic oil pressure.
3. Using a T25 Torx driver, remove the seat pan from the litter to access the hydraulic assembly.
4. Disconnect all connections to the main cable assembly.
5. Using a 3/4" combination wrench, remove the nut that secures the solenoid to the A valve (A) or B valve (B) (Figure 108). Save the nut for reinstallation.
6. Remove the solenoid from the valve. Save the solenoid for reinstallation.
7. Using a 7/8" combination wrench, remove the A valve or B valve from the hydraulic subassembly.

**Note:** Hydraulic oil will leak from the valve and manifold. Lay down towels to catch the oil.

8. Reverse steps to reinstall.
9. Check functionality by running the cot up and down several times. If necessary, add hydraulic oil. See "Filling the Hydraulics Assembly Reservoir" on [page 1-102](#).
10. Verify proper operation of the unit before returning it to service.


Figure 108

[Return To Table of Contents](#)

# Service Information

English

## HYDRAULIC MANUAL RELEASE VALVE REPLACEMENT

### Tools Required:

- T27 Torx Driver
- 7/16" Combination Wrench
- 1/8" Hex Wrench
- 7/8" Hex Wrench
- (2) Saw Horse

### Procedure:

1. Raise the cot to the full upright position.
2. Using two saw horses, support the litter and activate the manual release handle to relieve any hydraulic oil pressure.
3. Using a T27 Torx driver, remove the two button head cap screws (A) that secure the manual release cable bracket to the bottom of the hydraulic subassembly (Figure 109).
4. Using a 1/8" hex wrench, place the hex wrench through the stem at the groove in the valve body to hold the valve stem in position.
5. Using a 7/16" combination wrench, remove the Nylock hex nut (B) from each of the valve stems (Figure 109).
6. Using a 7/8" combination wrench, remove the valve (C or D) to be replaced (Figure 110).
7. Reverse steps to reinstall.
8. Check functionality by running the cot up and down several times. If necessary, add hydraulic oil. See "Filling the Hydraulics Assembly Reservoir" on [page 1-102](#).
9. Verify proper operation of the unit before returning it to service.


Figure 109


Figure 110

## HYDRAULIC CYLINDER REPLACEMENT

### Tools Required:

- 9/16" Combination Wrench
- 3/4" Combination Wrench
- 11/16" Combination Wrench
- 13/16" Combination Wrench
- 3/8" Combination Wrench
- 1/8" Hex Wrench
- (2) Saw Horse

### Procedure:

1. Raise the cot to the full upright position.
2. Using two saw horses, support the litter and activate the manual release handle and manually compress ram to remove the tension on the base cross tube connecting bolt.
3. Using a 3/4" and 9/16" combination wrench, remove the rod attachment pin (A), washer (B), and Nylock hex nut (C) that secure the hydraulic cylinder to the base (Figure 111).
4. Activate the manual release handle and fully compress the hydraulic cylinder.
5. Using a 11/16" and 13/16" combination wrench, remove both hoses (D) from the hydraulic cylinder (E) (Figure 112).

**Note:** Hydraulic oil will leak from the hoses and cylinder. Lay down towels to catch the oil.

6. Keep the hose ends high and upright to minimize the amount of fluid lost.
7. Using a 1/8" hex wrench and 3/8" combination wrench, remove the two socket head set screws (F) and Fiberlock hex nuts (G) that secure the hydraulic cylinder to the base (Figure 113).
8. Reverse steps to reinstall.
9. Check functionality by running the cot up and down several times. If necessary, add hydraulic oil. See "Filling the Hydraulics Assembly Reservoir" on [page 1-102](#).
10. Verify proper operation of the unit before returning it to service.


Figure 111


Figure 112


Figure 113

[Return To Table of Contents](#)

# Service Information

English

## HYDRAULIC HOSE REPLACEMENT

### Tools Required:

- 13/16" Combination Wrench
- 11/16" Combination Wrench
- (2) Saw Horse

### Procedure:

1. Raise the cot to the full upright position.
2. Using two saw horses, support the litter and activate the manual release handle to relieve any hydraulic oil pressure.
3. Using 11/16" and 13/16" combination wrenches, remove the damaged hose (A or B) (Figure 114).

### Notes:

- Pay attention to the routing of the hydraulic hose for reinstallation.
  - Hydraulic oil will leak from the hoses and cylinder. Lay down towels to catch the oil.
4. Reverse steps to reinstall.
  5. Check functionality by running the cot up and down several times. If necessary, add hydraulic oil. See "Filling the Hydraulics Assembly Reservoir" on [page 1-102](#).
  6. Verify proper operation of the unit before returning it to service.


Figure 114

## TERMINAL BLOCK REPLACEMENT

### Tools Required:

- T20 Torx Driver
- T25 Torx Driver

### Procedure:

1. Raise the cot to full up position.
2. Remove the battery and save for reinstallation.
3. Using a T25 Torx driver, remove the six outer button head cap screws (A) from the face plate (Figure 115). Save all screws for reinstallation.
4. Using a T20 Torx driver, remove the four inner delta screws (B) from the face plate to remove the face plate (Figure 115). Save all screws and the face plate for reinstallation.
5. Using a T20 Torx driver, remove the four delta screws (C) that secure the electronics assembly to the foot end enclosure and pull the electronics assembly out (Figure 116). Save all parts for reinstallation.
6. Unplug the black and red wires that connect the cot connector cable assembly (D) to the control board (E) (Figure 117).
7. Using a T20 Torx driver, remove the two delta screws (F) from the bottom plate of the foot end enclosure (G) to remove the enclosure (Figure 117). Save all parts for reinstallation.
8. Remove the terminal block and discard.
9. Reverse steps to reinstall.
10. Check functionality by running the cot up and down several times.
11. Verify proper operation of the unit before returning it to service.


Figure 115


Figure 116


Figure 117

# Service Information

---

English

## SIDERAIL ASSEMBLY REPLACEMENT (STANDARD)

### Tools Required:

- T25 Torx Driver

### Procedure:

1. Raise the cot to the full upright position.
2. Raise the siderail to the up and locked position.
3. Using a T25 driver, remove the three spindle screws that secure the siderail assembly.
4. Remove the siderail.
5. Reverse the above procedures to install the new siderail assembly.
6. Verify proper operation of the unit before returning it to service.

# Service Information

## SIDERAIL ASSEMBLY REPLACEMENT (XPS OPTION)

### Tools Required:

- T25 Torx Driver
- 1/4" Hex Wrench
- 3/16" Hex Wrench
- Slotted Screwdriver
- Deadblow Hammer
- Torque Wrench (ft-lb)

### Procedure:

1. Raise the cot and the fowler to the full upright position.
2. Remove the mattress.
3. Using a T25 Torx driver, remove the button head cap screw (A) and black bumper (B) on the side where you are replacing the siderail (Figure 118). Save the screw and bumper for reinstallation.
4. Using a slotted screwdriver, remove the outer rail bumper.

**Note:** When removing the outer rail bumper, hold on to the siderail main assembly, so it does not fall off. Also, note that the head end and middle siderail pivots may be loose and could fall off of the main assembly.


Figure 118

5. Using a 1/4" hex wrench, remove the socket head cap screws (C) that secure the siderail clamp (D) to the ratchet assembly at the foot end of the main assembly (Figure 119). During installation, torque both new screws to  $22 \pm 3.3$  ft-lb.

**Note:** The siderail will be loose, so do not operate or pull on the siderail.


Figure 119

6. Using a 3/16" hex wrench, remove the two socket head cap screws (E) that secure the middle siderail clamp (F) to the outer rail assembly (Figure 119).
7. Using a 3/16" hex wrench, remove the two socket head cap screws (G) that secure the top and bottom of the base/litter interface bracket to the outer rail assembly (Figure 119).
8. Reverse steps to reinstall. Use a deadblow hammer to reinstall the outer rail bumper.
9. Verify proper operation of the unit before returning it to service.

# Service Information

English

## RATCHET ASSEMBLY REPLACEMENT (XPS OPTION)

### Tools Required:

- 3/32" Hex Wrench
- 1/4" Hex Wrench
- 3/16" Hex Wrench
- Torque Wrench (ft-lb)

### Procedure:

1. Raise the cot and the fowler to the full upright position.
2. Raise the siderail to the up and locked position.
3. Using a 3/32" hex wrench, remove the two screws (A) that secure the ratchet cover (B) to the ratchet assembly. Remove the cover (Figure 120).
4. Using a 1/4" hex wrench, remove the socket head cap screws (C) that secure the siderail clamp (D) to the ratchet assembly at the foot end of the main assembly (Figure 120). During installation, torque both new screws to  $22 \pm 3.3$  ft-lb.
5. Using a 3/16" hex wrench, remove the four screws (E) that secure the ratchet assembly (F) to the overmold assembly and discard, then remove the ratchet assembly (Figure 120). During installation, torque the four new screws to  $9.5 \pm 1.5$  ft-lb.
6. Grasp the ratchet assembly and pull it toward the head end of the cot to remove.
7. Reverse steps to reinstall.
8. Verify proper operation of the unit before returning it to service.


Figure 120

# Service Information

English

## RELEASE HANDLE ASSEMBLY REPLACEMENT (XPS OPTION)

### Tools Required:

- 3/32" Hex Wrench
- Small Slotted Screwdriver

### Procedure:

1. Raise the cot and the fowler to the full upright position.
2. Raise the siderail to the up and locked position.
3. Using a 3/32" hex wrench, remove the four screws (A) that secure the release cover (B) to the overmold assembly to remove the release cover (Figure 121).
4. Using a small slotted screwdriver, pry the release handle return spring (C) up to remove the spring (Figure 121).
5. Grasp the release handle assembly (D), and hinge it upward on the spring side to remove it from the cover (Figure 121).
6. Reverse steps to reinstall.
7. Verify proper operation of the unit before returning it to service.


Figure 121

## SPRING HANDLE ASSEMBLY REPLACEMENT (XPS OPTION)

### Tools Required:

- 3/32" Hex Wrench
- Small Slotted Screwdriver

### Procedure:

1. Raise the cot and the fowler to the full upright position.
2. Raise the siderail to the up and locked position.
3. Using a 3/32" hex wrench, remove the four screws (A) that secure the release cover (B) to the overmold assembly to remove the release cover (Figure 121).
4. Using a small slotted screwdriver, pry the release handle return spring (C) up to remove the spring (Figure 121).
5. Reverse steps to reinstall.
6. Verify proper operation of the unit before returning it to service.

[Return To Table of Contents](#)

# Child Restraint with X-Restraint Package Installation

English

To attach the Britax Meridian SICT Series No. 7200/A/2010 Convertible Child Restraint to the cot with the X-Restraint package (6500-001-430):

1. Lay the cot flat.
2. Open the waist restraint strap. Lay the straps to the side and out of the way.
3. Position the child restraint to face the foot end of the cot.
4. Place the child restraint into the reclined position.
5. Raise the cot fowler up to align with the back of the child restraint.
6. Loop the upper anchorage strap with the anchor fitting and attachment clip from the child restraint around the cot fowler (Figure 122). Pull the adjustment strap to tighten and remove slack.
7. Pull the waist restraint under the siderail and through the bottom of the child restraint (Figure 123).
8. Firmly push down on the seat with your hand, while pulling the waist restraint with the other hand to tighten the waist restraint.
9. If your ankle restraint is long enough to reach the foot end of the child restraint, pull the ankle restraint through the foot end of the child restraint (Figure 124).
10. If your ankle restraint is too short to reach the foot end of the child restraint:
  - Remove the thigh restraint from your cot.
  - Relocate the thigh restraint to the foot end lift bars (Figure 125).
  - Pull the thigh restraint through the foot end of the child restraint (Figure 124).
11. Firmly push down on the seat with your hand, while pulling the leg restraint with the other hand to tighten the leg restraint.
12. While holding the tension, close the rear facing lock-off on the buckle side to hold the tension (Figure 126).
13. Close the rear facing lock-off on the opposite side.
14. Place the baby into the child restraint and secure per manufacturer's instructions.

Figure 122

Figure 123


Figure 124

Figure 125

Figure 126

[Return To Table of Contents](#)

# Warranty

---

Stryker EMS, a division of the Stryker Corporation, offers one warranty option in the United States:

**Two (2) year parts and labor.** Stryker EMS warrants to the original purchaser that its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of two (2) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective. Expendable components, i.e. mattresses, restraints, I.V. poles, storage nets, storage pouches, oxygen straps, and other soft goods, have a one (1) year limited warranty.

The Stryker **Power-PRO™** XT is designed for a 7 year expected service life under normal use conditions, and with appropriate periodic maintenance as described in the maintenance manual. Stryker warrants to the original purchaser that the welds on the **Power-PRO™** XT will be free from structural defects for the expected 7 year life of the product as long as the original purchaser owns the product. Original purchasers will also obtain a three (3) year limited parts warranty for the X-frame components of the **Power-PRO™** cot and a three (3) year limited power train warranty covering the motor pump assembly and hydraulic cylinder assembly. Stryker's obligation under this three (3) year limited warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any part that is, in the sole discretion of Stryker, found to be defective.

**SMRT™ Power Warranties.** Stryker EMS warrants the **SMRT™** Charger for the same duration as the Stryker product for which it is furnished. All **SMRT™** Paks are warranted to be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of one (1) year.

Upon Stryker's request, purchaser shall return to Stryker's factory any product or part (freight prepaid by Stryker) for which an original purchaser makes a warranty claim.

## WARRANTY EXCLUSIONS AND OTHER CONDITIONS

- Any improper use or alteration or repair by unauthorized service providers in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty.
- Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty.
- There are extenuating circumstances and events that may alter the performance of the products such as an ambulance accident. In Stryker's discretion, certain circumstances may allow for evaluation of the product post ambulance accident which could allow for continued use of the product. If products recommended to be removed from service are put back into service, Stryker will consider the product as being subject to unusual stress and improperly maintained. Products which are subject to unusual stress and improper maintenance are not subject to Stryker's warranty as noted above. In addition, Stryker will not indemnify any customer for any third-party claims related to injuries caused by products that have been involved in accidents.
- This warranty is void if the label bearing the serial number of the product has been removed or defaced.
- This warranty is void if the product is not purchased from an authorized Stryker dealer.

This statement constitutes Stryker EMS's entire warranty with respect to the aforesaid equipment. STRYKER MAKES NO OTHER WARRANTY OR REPRESENTATION EITHER EXPRESSED OR IMPLIED, EXCEPT AS SET FORTH HEREIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL STRYKER BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

No employee or representative of Stryker is authorized to change this warranty in any way.

# Warranty

---

English

## STRYKER EMS RETURN POLICY

Cots, Stair Chairs, Evacuation Chairs, Cot Fasteners and Aftermarket Accessories may be returned up to 180 days of receipt if they meet the following guidelines:

### Prior to 30 Days

- 30 day money back guarantee in effect
- Stryker EMS is responsible for all charges
- Returns will not be approved on modified items

### Prior to 90 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 10% restocking fee

### Prior to 180 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 25% restocking fee

## RETURN AUTHORIZATION

Stryker customer service department must approve any merchandise return and will provide an authorization number to be printed on any returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.

## DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Claims for any short shipment must be made within thirty (30) days of invoice.

## INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

## PATENT INFORMATION

The Stryker **Power-PRO™** XT cot is covered by one or more of the following patents:

United States	5,537,700	5,575,026	6,908,133	7,398,571	7,540,047
---------------	-----------	-----------	-----------	-----------	-----------

Other patents pending

The Stryker **SMRT™** Power System is covered by one or more of the following patents:

United States	5,977,746	6,018,227
---------------	-----------	-----------

Other patents pending

# EMC Information

## POWER-PRO XT

Guidance and manufacturer's declaration - electromagnetic emissions		
The Model 6506 <b>Power-PRO™</b> XT cot is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 6506 <b>Power-PRO™</b> XT cot should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Model 6506 <b>Power-PRO™</b> XT cot uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group 2	The Model 6506 <b>Power-PRO™</b> XT cot must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Cot: Class A	The Model 6506 <b>Power-PRO™</b> XT cot is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
	<b>SMRT™</b> Charger (6500-201-010): Class A	The <b>SMRT™</b> Charger is suitable for use in all establishments other than domestic establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Cot: N/A	Not applicable
	<b>SMRT™</b> Charger (6500-201-010): Class A	
Voltage fluctuations flicker emissions IEC 61000-3-3	Cot: N/A	Not applicable
	<b>SMRT™</b> Charger (6500-201-010): Complies	

# EMC Information

English

## POWER-PRO™ XT (CONTINUED)

Guidance and manufacturer's declaration - electromagnetic immunity			
The Model 6506 Power-PRO™ XT cot is intended for use in the electromagnetic environment specified below. The customer or user of the Model 6506 Power-PRO™ XT cot should assure that it is used in such an environment.			
IMMUNITY test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN/IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Applies to: • Cot • <b>SMRT™</b> Charger (6500-201-010)
Electrostatic fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines  $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines  $\pm 1$ kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. Applies to: • <b>SMRT™</b> Charger (6500-201-010)
Surge IEC 61000-4-5	$\pm 8$ kV differential mode  $\pm 2$ kV common mode	$\pm 8$ kV differential mode  $\pm 2$ kV common mode	Mains power quality should be that of a typical commercial or hospital environment. Applies to: • <b>SMRT™</b> Charger (6500-201-010)
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 0,5 cycle $40\% U_T$ ( $60\%$ dip in $U_T$ ) for 5 cycles $70\% U_T$ ( $30\%$ dip in $U_T$ ) for 25 cycles $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 5 seconds	$<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 0,5 cycle $40\% U_T$ ( $60\%$ dip in $U_T$ ) for 5 cycles $70\% U_T$ ( $30\%$ dip in $U_T$ ) for 25 cycles $<5\% U_T$ ( $>95\%$ dip in $U_T$ ) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the charger requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60Hz) magnetic field EN/IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. Applies to: • Cot • <b>SMRT™</b> Charger (6500-201-010)
<b>Note:</b> $U_T$ is the alternating current mains voltage prior to application of the test level.			

# EMC Information

## POWER-PRO XT (CONTINUED)

Guidance and manufacturer's declaration - electromagnetic immunity			
The Model 6506 <b>Power-PRO™</b> XT cot is intended for use in the electromagnetic environment specified below. The customer or user of the Model 6506 <b>Power-PRO™</b> XT cot should assure that it is used in such an environment.			
IMMUNITY test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model 6506 <b>Power-PRO™</b> XT cot, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p><b>Recommended separation distance</b>  <math>d=1.2\sqrt{P}</math></p> <p>Applies to:</p> <ul style="list-style-type: none"> <li>• <b>SMRT™</b> Charger (6500-201-010)</li> </ul>

# EMC Information

## POWER-PRO™ XT (CONTINUED)

English

### Guidance and manufacturer's declaration - electromagnetic immunity

The Model 6506 **Power-PRO™** XT cot is intended for use in the electromagnetic environment specified below. The customer or user of the Model 6506 **Power-PRO™** XT cot should assure that it is used in such an environment.

IMMUNITY test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF EN/IEC 61000-4-3	20 V/m 80 MHz to 2,5 GHz	20 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Model 6506 <b>Power-PRO™</b> XT cot, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $D=(1.2)(\sqrt{P})$ $D=(0.18)(\sqrt{P})$ 80 MHz to 800 MHz $D=(0.35)(\sqrt{P})$ 800 MHz to 2,5 GHz

where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and  $d$  is the recommended separation distance in metres (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup>

Interference may occur in the vicinity of equipment with the following symbol:

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 6506 **Power-PRO™** XT cot is used exceeds the applicable RF compliance level above, the Model 6506 **Power-PRO™** XT cot should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model 6506 **Power-PRO™** XT cot.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths are less than 20 V/m.

[Return To Table of Contents](#)

# EMC Information

## POWER-PRO™ XT (CONTINUED)

Recommended separations distances between portable and mobile RF communication equipment and the Model 6506 Power-PRO™ XT cot			
The Model 6506 Power-PRO™ XT cot is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model 6506 Power-PRO™ XT cot can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 6506 Power-PRO™ XT cot as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $D=(1.2)(\sqrt{P})$	80 MHz to 800 MHz $D=(0.18)(\sqrt{P})$	800 MHz to 2,5 GHz $D=(0.35)(\sqrt{P})$
0.01	0.12	0.018	0.035
0.1	0.38	0.57	0.11
1	1.2	0.18	0.35
10	3.8	0.57	1.1
100	12	1.8	3.5
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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