M-1® Cot Fastener

Operations Manual

REF 6373

REF 6376

REF 6381

EN
CS
DA
DE
EL
ES
ET
FI

FR
HR
HU
IT
LT
LV
NL
NO

PL
PT
RO
RU
SK
SL
SV
TR

Symbols

	Consult instructions for use
	General warning
	Caution
	CE mark
	UK Conformity Assessment mark
	Importer
	Unique device identifier
	Authorized representative in the European Community
	Authorized representative in Switzerland
	European medical device
	Catalogue number
	Lot (batch) code
	Serial number
	For US patents see www.stryker.com/patents
	Manufacturer
	Date of manufacture
	Safe working load

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	3
Product description	3
Indications for use	3
Clinical benefits	3
Contraindications	3
Expected service life	3
Disposal/recycle	3
Specifications	4
Product illustration	5
Contact information	5
Serial number location	6
Date of manufacture	6
Cot compatibility	6
Installation	7
Installing the cot fastener	7
Installing the catch plate	10
Installation checklist	13
Operation	14
Operating the cot fastener	14
Testing the cot fastener	15
Cleaning	16
Suggested cleaners	16

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Install the cot fastener by qualified personnel only. Improper installation could result in injury to the patient or operator.
 - Always make sure that all cots meet the installation specifications for the Stryker cot fastener system.
 - Always replace the cot fastener if it has been involved in an accident to avoid the risk of injury due to product damage.
 - Do not allow untrained assistants to assist in the operation of the product.
 - Always make sure that the cot is locked in the fastener before you transport the patient. Injury may occur if the cot is not secure.
 - Always use any appropriate personal protective equipment while power washing to avoid inhaling contagion. Power washing equipment may aerate contamination.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not steam clean or ultrasonically clean the product.
 - Do not exceed 180 °F (82 °C) as the maximum water temperature.
 - Do not exceed 240 °F (115 °C) as the maximum air dry temperature (cart washers).
 - Do not exceed 1500 psi (103.4 bar) as the maximum water pressure. If you are using a hand held wand to wash the product, keep the pressure nozzle at a minimum of 24 in. (61 cm) from the product.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 6373/6376/6381 **M-1**® Cot Fastener holds a compatible ambulance cot in an emergency vehicle for patient transport purposes. The cot fastener performs this function after the user manually loads the compatible cot into the emergency vehicle and guides it into the cot fastener.

Indications for use

The Stryker Model 6373/6376/6381 is a manual, ambulance cot fastener that secures a compatible ambulance cot during transport. The cot fastener is intended to be used in pre-hospital environments, in emergency and non-emergency applications. Operate the product by trained professionals only. This includes emergency medical service and medical first responders.

Clinical benefits

Cot: patient transport

Fastener: support cot for transport

Cot and fastener system: support and transport patients

Contraindications

None known.

Expected service life

Cot fasteners have a seven year expected service life under normal use conditions and with appropriate periodic maintenance.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.






Specifications

Model	6373 - tray mounted cot fastener (6373-000-000)
	6376 - floor mounted cot fastener (6376-000-000)
	6381 - floor mounted cot fastener (6381-000-000)
Compatibility	See <i>Cot compatibility</i> (page 6)

The Model 6100 **M-1**® Roll-In System cot, in combination with a Model 6373/6376/6381 **M-1** Cot Fastener System, is designed in accordance to BS EN 1789, Medical Vehicles and their Equipment - Road Ambulances. The cot fastener system must be installed by qualified personnel to be in conformance to this standard. 1/4 in. steel backing plates (or equivalent) are required beneath the mounting surface (ambulance floor or tray system) at all mounting points. Before putting the system into service, test the installation using BS EN 1789 as a minimum.

The Model 6100 **M-1** Roll-In System cot, equipped with the AS/NZS 4535 upgrade, in combination with a Model 6373/6376/6381 **M-1** Cot Fastener System, is designed in accordance to AS/NZS 4535 Australian/New Zealand Standard - Ambulance Restraint Systems. The cot fastener system must be installed by qualified personnel to be in conformance to this standard. 1/4 in. steel backing plates (or equivalent) are required beneath the mounting surface (ambulance floor or tray system) at all mounting points. Before putting the system into service, test the installation using AS/NZS 4535 as a minimum.

Side restraint belts are required to conform to BS EN 1789 Medical Vehicles and their Equipment - Road Ambulances. Side restraint belts are not required for applications where this standard does not apply.

Environmental conditions	Operation	Storage and transportation
Temperature		
Relative humidity		
Atmospheric pressure		

Product illustration


Figure 1 – M-1® Cot Fastener

A	Guide rail
B	Front hook
C	Bumper
D	Head end
E	Foot end
F	Red release button
G	Latch housing

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Figure 2 – Serial number location

Date of manufacture

The year of manufacture is the first 2 digits of the serial number.

Cot compatibility

The Stryker Model 6373, 6376, and 6381 **M-1** Cot Fastener systems are compatible only with cots that conform to the installation specifications.

WARNING

- Install the cot fastener by qualified personnel only. Improper installation could result in injury to the patient or operator.
 - Always make sure that all cots meet the installation specifications for the Stryker cot fastener system.
 - Always replace the cot fastener if it has been involved in an accident to avoid the risk of injury due to product damage.
-

Stryker

- Model 6100 **M-1** Roll-In System

Ferno-Washington

- Model X-2*

Note - * 2000 model year or earlier. Stryker is not responsible for changes in specifications to other manufacturer's cots.

Installation

Installing the cot fastener

The Stryker **M-1** Cot Fastener systems are compatible only with cots that conform to the installation specifications.

WARNING

- Install the cot fastener by qualified personnel only. Improper installation could result in injury to the patient or operator.
 - Always make sure that all cots meet the installation specifications for the Stryker cot fastener system.
-

Note - Federal ambulance specifications regarding aisle width and attendant seating may affect positioning of the mounting brackets. If specifications do not apply, the installer should make sure that there is adequate clearance for proper patient access, traction splint projection, rear door closure, etc.

To install the cot fastener in a vehicle:

1. Determine the desired position of the cot inside the vehicle allowing for proper clearance around the cot.

Note - If DIN hardware is already in place, use the dimensions in Figure 3, Figure 4, Figure 5, and Figure 6 as a reference to locate the fastener assembly.


Figure 3 – Reference dimension - View 2A


Figure 4 – Reference dimension - View 2B


Figure 5 – Reference dimension - View 2C


Figure 6 – Reference dimension - View 2D

2. Drill the 18 mounting holes using the assembled fastener as a template.

Note - All dimensions are for reference only. All holes should provide adequate clearance for 8mm screws. A 9mm drill size is recommended.

3. Anchor the foot end latch housing (A) to the floor or tray using four M8 flat head cap screws (B) (not supplied).

Note - A steel backing plate (C) (or equivalent) is required for added support (not supplied) (Figure 7 or Figure 8).


Figure 7 – Floor mount foot end - View 3A


Figure 8 – Tray mount foot end - View 3B

4. Anchor the head end hooks (A) to the floor or tray using eight M8 flat head cap screws (B) (not supplied).

Note

- A steel backing plate (C) (or equivalent) is required for added support (not supplied) (Figure 9 or Figure 10).
- The red release knobs must be oriented toward the center of the fastener (Figure 11).


Figure 9 – Floor mount head end - View 4A


Figure 10 – Tray mount head end - View 4B


Figure 11 – Proper orientation of head end hooks - View 4C

5. Anchor the guide rails to the floor or tray using two M8 button head cap screws (A), two M8 flat washers (B), and two M8 nuts (C) (not provided) (Figure 12 or Figure 13).


Figure 12 – Floor mount guide rails - View 5A


Figure 13 – Tray mount guide rails - View 5B

Installing the catch plate

The catch plate is intended for use with the M-1 safety hook option (610000010001).

Catch plate kit (610000010022)

- (1) Catch plate (610000010015)

Additional part used (not included):

- (2) M6 Allen bolt

Tools required (metric or imperial):

- M5 Allen key or 3/16" hex wrench
- M6 Allen key or 1/4" hex wrench

Stryker recommends that you follow the specifications listed below for the catch plate installation:

- Angle (A) of the loading tray does not exceed 16° (Figure 14)
- Height (B) from the top of the catch plate to the ground is between 28.0" and 29.4" (Figure 14)
- Horizontal clearance (C) between the safety hook and catch plate is sufficient for unloading (Figure 15)



Figure 14 – Catch plate and loading tray specifications

Figure 15 – Safety hook and catch plate clearance

Procedure:

1. Using an M6 Allen key, install the supplied catch plate (610000010015) (D) and two M6 Allen bolts (E) (not included) directly behind the cot fastener latch (Figure 16).


Figure 16 – Catch plate

Note - Before you drill the two holes for the M6 Allen bolts, make sure that there is free space under the table to avoid product or vehicle damage.

2. Test unloading the cot from a vehicle to confirm that the safety hook connects with the catch plate. Using an M5 Allen key, adjust the hook tip (F), as needed, to make sure that the safety hook and catch plate connect while unloading (Figure 17).


Figure 17 – Hook tip adjustment

Note - Push the cot legs against the loading tray when you test unloading.

Installation checklist

Follow this checklist with a Stryker Model 6373/6376/6381 **M-1** Cot Fastener compatible cot.

- _____ Make sure that you do not have any unused components after installation. Your **M-1** Cot Fastener does not ship with any extra components. If you have any unused components after installation, call Stryker service.
- _____ Visually check that all bolts and screws are tight with no signs of protruding or missing fasteners.
- _____ Lift the vehicle bumper to the raised position, if equipped.
- _____ Load the compatible cot into the cot fastener. See *Cot compatibility* (page 6).
- _____ Make sure that the cot is locked into the cot fastener by pulling in and out and side to side on the foot end of the cot.
- _____ Test the cot for proper operation. See *Testing the cot fastener* (page 15).
- _____ Check that the angle of the loading tray does not exceed 16° (**M-1** safety hook option).
- _____ Check that the height from the top of the catch plate to the ground is between 28.0” and 29.4” (**M-1** safety hook option).
- _____ Check that the horizontal clearance between the safety hook and catch plate is sufficient for unloading (**M-1** safety hook option).

Product serial number:			
Installed by:		Date:	
Inspected by:		Date:	

Note - Maintain a copy of this record for at least seven years.

Operation

Operating the cot fastener

WARNING

- Do not allow untrained assistants to assist in the operation of the product.
 - Always make sure that the cot is locked in the fastener before you transport the patient. Injury may occur if the cot is not secure.
-

CAUTION - Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.

1. Before you operate the fastener, press the red release button (F) on the latch housing (G) to make sure that the latch is in the released position (Figure 1, *Product illustration* (page 5)).
2. Make sure that the foot end latch housing is between the loading wheels when you load the cot into the fastener.

Note - The fastener guide rails (A) will center the head end of the cot (Figure 1, *Product illustration* (page 5)).

3. Guide the cot into the vehicle or loading tray system until the rear latch post locks into the latch housing.

Note - The foot end latch housing is equipped with a two-stage locking device (Figure 18). When the cot is locked in the fastener, the red release button will return to its original position.


Figure 18 – Two-stage locking device

A	Flipper
B	Red release knob
C	Front retaining hook

The Stryker Model 6373, 6376, and 6381 cot fastener systems are equipped with a two-stage locking device. In the unlikely event the vehicle is involved in an accident, the two-stage locking device at the head end of the fastener will deploy to lock the head end of the cot in the cot fastener system until you manually release the cot.

To release the cot from the two-stage locking device, pull both of the red release knobs (B) on the inside of the front retaining hooks (C) (Figure 18).

Testing the cot fastener

To test the cot fastener:

1. Guide the cot into the cot fastener.
2. Make sure that the head end retaining hooks secure the loading wheel cross tube on the cot.
3. Make sure that the foot end latch housing secures the retaining pin on the cot.
4. Make sure that the two-stage locking device functions and is free of debris.
5. Test the cot fastener using the appropriate standard as a minimum.
 - BS EN 1789 Medical Vehicles and Their Equipment - Road Ambulances
 - AS NZS 4535 Australian/New Zealand Standard - Ambulance Restraint Systems

Cleaning

WARNING - Always use any appropriate personal protective equipment while power washing to avoid inhaling contagion. Power washing equipment may aerate contamination.

CAUTION

- Do not steam clean or ultrasonically clean the product.
 - Do not exceed 180 °F (82 °C) as the maximum water temperature.
 - Do not exceed 240 °F (115 °C) as the maximum air dry temperature (cart washers).
 - Do not exceed 1500 psi (103.4 bar) as the maximum water pressure. If you are using a hand held wand to wash the product, keep the pressure nozzle at a minimum of 24 in. (61 cm) from the product.
-

The product is power washable. The product may show some signs of oxidation or discoloration from continuous washing. No degradation of the product's performance will occur from power washing as long as you follow the proper procedures.

- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- Power wash with recommended cleaners. Hose down the product and towel dry the guide.
- Power wash the rails and interface plates with a hand held wand unit or wipe the product with a clean cloth and recommended cleaners.

Suggested cleaners

In general, when used in concentrations recommended by the manufacturer, either phenolic type or quaternary (excluding **Virex® TB**) type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may occur.

Suggested cleaners include:

- Quaternary cleaners (active ingredient - ammonium chloride)
- Phenolic cleaners (active ingredient - o-phenylphenol)
- Chlorinated bleach solution (5.25% - less than 1 part bleach to 100 parts water)

Avoid oversaturation. Do not allow the product to stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

Note

- Failure to follow the above directions when using these types of cleaners may void this product's warranty.
- Always wipe the product with clean water and dry after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature corrosion of critical components.

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