**Transfer Flat
Drap de transfert
Transfermatte
Tilzeil
Telo per trasferimento
Överföringsdyna
Transportdug
Siirtoalusta**

REF

Model/Modèle/Modell/Modello/Malli 6005

stryker®

Operations Manual

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Table of Contents

Symbols and Definitions	1-4
Symbols	1-4
Warning/Caution/Note Definition	1-4
Introduction	1-5
Product Description	1-5
Intended Use of Product	1-5
Specifications	1-6
Contact Information	1-7
Lot Number Location	1-7
Summary of Safety Precautions	1-8
Setup Procedures	1-10
Product Inspection	1-10
Operation Guide	1-11
Transferring and Transporting a Patient	1-11
Transferring a Patient from a Cot to a Patient Support Platform	1-13
Storing the Transfer Flat	1-14
Cleaning	1-15
Recommended Routine Care/Cleaning Method	1-15
Washing Procedure	1-15
Removal of Iodine Compounds	1-16
Warranty	1-17
Stryker EMS Return Policy	1-18
Return Authorization	1-18
Damaged Merchandise	1-18
International Warranty Clause	1-18

Symbols and Definitions

English

SYMBOLS

	Warning, consult accompanying documentation
	Safe Working Load
	No lifting with equipment
	Maximum weight
	Manufacturer

WARNING/CAUTION/NOTE DEFINITION

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

NOTE

Provides special information to make maintenance easier or important instructions clearer.

[Return To Table of Contents](#)

Introduction

This manual is designed to assist you with the operation of the Transfer Flat. Read this manual thoroughly before using the unit. To ensure safe operation of the Transfer Flat, it is recommended that methods and procedures be established for educating and training staff on the safe operation of this unit.

PRODUCT DESCRIPTION

The model 6005 Transfer Flat is a support surface with perimeter handles to facilitate the transfer and transport of patients during EMS service calls. It is a comfortable, flexible, durable and cleanable support. A retaining feature to facilitate folding and storage has also been integrated into the device.

INTENDED USE OF PRODUCT

The model 6005 Transfer Flat is intended to be used as a manual patient transfer device. The device shall be used by emergency medical service personnel to transfer patients, with minimal disturbance, to another surface.


WARNING

- Do not modify the Transfer Flat. Modifying the Transfer Flat can cause unpredictable operation resulting in injury to the patient or operator. Modifying the Transfer Flat will also void its warranty.
 - The Transfer Flat is not designed or recommended for use with patients with suspected cervical, spinal, or fracture injuries.
-

Introduction

English

SPECIFICATIONS

	Safe Working Load Note: Safe Working Load indicates the sum of the patient and accessory weight.	1600 lb	725 kg
Length		108"	274 cm
Width		60"	152 cm
Weight		10.0 lb	4.5 kg
Operation and Storage Temperature			

Stryker reserves the right to change specifications without notice.

The yellow and black color scheme is a proprietary trademark of Stryker Corporation.

Introduction

CONTACT INFORMATION

Contact Stryker Customer Service or Technical Support at (800) 327-0770 or (269) 324-6500.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Please have the lot number, located on the specification label (A) of your Stryker product available (as shown in Figure 1), when calling Stryker Customer Service or Technical Support. Include the lot number in all written communication.

LOT NUMBER LOCATION


Figure 1: Lot Number & Location

Summary of Safety Precautions

Carefully read and strictly follow the warnings and cautions listed on these pages.

WARNING

- Do not modify the Transfer Flat. Modifying the Transfer Flat can cause unpredictable operation resulting in injury to the patient or operator. Modifying the Transfer Flat will also void its warranty.
 - The Transfer Flat is not designed or recommended for use with patients with suspected cervical, spinal, or fracture injuries.
 - Improper usage of the Transfer Flat can cause injury to the patient or operator. Operate the unit only as described in this manual.
 - Do not operate the Transfer Flat without proper training.
 - Do not exceed the 1600 pound weight capacity for the Transfer Flat.
 - Do not use the Transfer Flat if any bodily fluids are present. Inspect the Transfer Flat for fluids and clean prior to use.
 - Inspect the Transfer Flat before and after each use for wear, fraying, or other damage. If device is damaged, discontinue use.
 - To assist with proper operation and to avoid injury, use industry approved techniques and guidelines for lifting when using the Transfer Flat to transfer a patient.
 - To avoid injury to the patient and/or operator, operators should never attempt to transfer or transport patient loads greater than what they can safely lift. Evaluate the load prior to attempting the transfer or transport and use additional assistance, if necessary.
 - Always use all restraint straps to secure the patient to the cot. An unrestrained patient may fall from the cot and be injured.
 - The patient must be secured on the cot before transport.
 - When transferring the patient with the Transfer Flat, ensure that there is enough clearance underneath the Transfer Flat to avoid impact with any obstacles.
 - When transferring or transporting the patient, verify that the Transfer Flat handles do not interfere with the cot functionality.
 - When loading the cot into an ambulance, ensure that the Transfer Flat does not interfere with the operation of the safety hook, cot fastener, or the in-fastener shut-off system.
 - The Transfer Flat is not designed for use with winch/tow straps, lifts, or hoist mechanisms. To avoid injury or damage to the equipment, do not attach any type of hoist mechanism to the Transfer Flat.
 - The Transfer Flat is not radiolucent and it cannot be used during x-ray.
 - Water, ice, and debris on the stairs can affect operator footing and proper operation of the Transfer Flat. To reduce the risk of injury, clear the path or consider an alternate route.
-

Summary of Safety Precautions

WARNING - CONTINUED

- Only use the Transfer Flat for patient transfers between fixed patient support surfaces that are level with one another.
 - Do not power wash. The Transfer Flat is not designed to be power washable. Risk of airborne contamination may cause injury to users.
 - **SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY.** If the products described on [page 1-15](#) are used to clean Stryker EMS equipment, measures must be taken to ensure the products are wiped with clean water and air dried thoroughly following cleaning.
 - Failure to properly clean or dispose of contaminated Transfer Flats increases the risk of exposure to bloodborne pathogens and may cause injury to the patient or the operator.
-

CAUTION

- The Transfer Flat should not be used as a permanent surface on any patient support platform.
 - Do not store the Transfer Flat in direct sunlight or in a moist environment. Keep device and equipment out of direct contact with flames.
 - To avoid the risk of damage, the Transfer Flat must be stored where the temperature range is between -22 to 158 degrees Fahrenheit (-30 to 70 degrees Celsius).
 - Using the transfer flat to drag a patient could compromise patient comfort and risk further injury.
 - Always abide to industry approved lifting and patient handling techniques with every patient.
 - Avoid surfaces that could damage the transfer flat.
 - Always inspect the product for damage before and after each use.
 - Discard any transfer flat that has been damaged due to excessive wear.
 - Read this manual carefully before using the unit. Always follow correct operation procedures as outlined in this manual.
-

Setup Procedures

English

Ensure that all shipping and packaging materials have been removed from the product(s) prior to use.

Unpack the cartons and check all items for proper functionality. It is important that the Transfer Flat is working properly before it is put into service. Use the Product Inspection checklist below and the operation instructions to check the unit before it is put into service. Thoroughly inspect for shipping damage.

WARNING

- Improper usage of the Transfer Flat can cause injury to the patient or operator. Operate the unit only as described in this manual.
 - Do not modify the Transfer Flat. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product also voids its warranty.
-

Note:

- This manual should be considered a permanent part of the unit and should accompany the product even if the unit is subsequently sold.
- Stryker continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor modifications between your Transfer Flat and this manual. If you have any questions, please contact Stryker Customer Service or Technical Support at (800) 327-0770 or (269) 324-6500.

Product Inspection

The condition of the Transfer Flat is the responsibility of the owner. It is important that the Transfer Flat is not damaged before use. Use the following checklist and the operation instructions to check the unit before each use.

- _____ Inspect unit for wear, fraying, tears, missing stitches, loose handles or other damage
- _____ Inspect unit for bodily fluids to prevent biocontamination
- _____ Inspect each handle for any visible damage and tug on each handle to ensure that all threads are secure
- _____ Verify that snap buttons work correctly before storing the unit

Operation Guide

TRANSFERRING AND TRANSPORTING A PATIENT

These instructions outline the proper techniques for safely securing and transporting a patient using the Transfer Flat.

WARNING

- Do not operate the Transfer Flat without proper training.
- Do not exceed the 1600 pound weight capacity for the Transfer Flat.
- Do not use the Transfer Flat if any bodily fluids are present. Inspect the Transfer Flat for fluids and clean prior to use.
- Inspect the Transfer Flat before and after each use for wear, fraying, or other damage. If device is damaged, discontinue use.
- The Transfer Flat is not designed or recommended for use with patients with suspected cervical, spinal, or fracture injuries
- To assist with proper operation and to avoid injury, use industry approved techniques and guidelines for lifting when using the Transfer Flat to transfer a patient.
- To avoid injury to the patient and/or operator, operators should never attempt to transfer or transport patient loads greater than what they can safely lift. Evaluate the load prior to attempting the transfer or transport and use additional assistance, if necessary.

When the patient is lying down on the floor or on a bed:

1. Unfold the Transfer Flat.
2. Use appropriate industry techniques to center the patient on the Transfer Flat as shown in Figure 2.


Figure 2

Operation Guide

English

- Using the handles, lift and carry the patient to a cot as shown in Figure 3. The Transfer Flat can accommodate up to twelve operators.


Figure 3

- Wrap the flat around the patient to ensure that no handles are hanging off the side, front, or back of the unit to prevent any interference with the functionality of other devices.
- Ensure that the patient is properly secured on the cot.

WARNING

- Always use all restraint straps to secure the patient to the cot. An unrestrained patient may fall from the cot and be injured.
- The patient must be secured on the cot before transport.
- When transferring the patient with the Transfer Flat, ensure that there is enough clearance underneath the Transfer Flat to avoid impact with any obstacles.
- When transferring or transporting the patient, verify that the Transfer Flat handles do not interfere with the cot functionality.
- When loading the cot into an ambulance, ensure that the Transfer Flat does not interfere with the operation of the safety hook, cot fastener, or the in-fastener shut-off system.
- The Transfer Flat is not designed for use with winch/tow straps, lifts, or hoist mechanisms. To avoid injury or damage to the equipment, do not attach any type of hoist mechanism to the Transfer Flat.
- The Transfer Flat is not radiolucent and it cannot be used during x-ray.
- Water, ice, and debris on the stairs can affect operator footing and proper operation of the Transfer Flat. To reduce the risk of injury, clear the path or consider an alternate route.

Operation Guide

TRANSFERRING A PATIENT FROM A COT TO A PATIENT SUPPORT PLATFORM

To transfer the patient from one patient support platform (for example, bed, stretcher, gurney, operating table) to another:

1. Position one patient support platform alongside the other patient support platform as closely as possible.
2. Secure both patient support platforms. When transferring to another platform, refer to the transfer and safety requirements of the specified platform.
3. Raise the receiving patient support platform outer siderails, if applicable to prevent overtravel (not shown in Figure).

WARNING

Only use the Transfer Flat for patient transfers between fixed patient support surfaces that are level with one another.

4. Adjust the patient support platform heights, so that they are level with one another.
5. Use the Transfer Flat to move the patient across from one surface to the other as shown in Figure 4.


Figure 4

CAUTION

The Transfer Flat should not be used as a permanent surface on any patient support platform.

Operation Guide

English

STORING THE TRANSFER FLAT

The Transfer Flat rolls up and forms its own carrying case. Make sure that the transfer flat is clean and dry before storage.

To store the Transfer Flat:

1. Fold each long side over toward the center of the Transfer Flat as shown in Figures 5 and 6.


Figure 5


Figure 6

2. Starting at either end, roll up the Transfer Flat as shown in Figure 7.
3. Slip the strap through the D-ring and snap it securely as shown in Figure 8. The handle can be used to carry the Transfer Flat.


Figure 7


Figure 8

CAUTION

- Do not store the Transfer Flat in direct sunlight or in a moist environment. Keep device and equipment out of direct contact with flames.
- To avoid the risk of damage, the Transfer Flat must be stored where the temperature range is between -22 to 158 degrees Fahrenheit (-30 to 70 degrees Celsius).

Cleaning

Thoroughly clean the Transfer Flat after each use.

RECOMMENDED ROUTINE CARE/CLEANING METHOD

- Apply disinfectant solution either by spray, solution or pre-impregnated wipes (do not soak Transfer Flats).
- **DO NOT STEAM CLEAN, PRESSURE WASH, OR ULTRASONICALLY CLEAN** Transfer Flats.
- Allow to air dry.

WARNING

Do not power wash. The Transfer Flat is not designed to be power washable. Risk of airborne contamination may cause injury to users.

WASHING PROCEDURE

IMPORTANT: DILUTE ALL DISINFECTANTS IN ACCORDANCE WITH MANUFACTURER'S DIRECTIONS.

Follow the cleaning solution manufacturer's dilution recommendations exactly.

In general, when used in those concentrations recommended by the manufacturer, either phenolic type or quaternary (**excluding Virex® TB**) type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may result.

Suggested cleaners for the Transfer Flat:

- Quaternary Cleaners (active ingredient - ammonium chloride)
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)

WASHING PROCEDURE - CONTINUED

These products are NOT considered mild detergents. They are corrosive in nature and may cause damage to your Transfer Flat if used improperly. Transfer Flats must be rinsed with clean water and dried thoroughly after using corrosives such as quaternary or chlorine bleach. Failure to properly rinse and dry the Transfer Flat leaves a corrosive residue on the surface, likely causing premature corrosion.

Note: Chlorine bleach and quaternary cleaners may cause discoloration or fading of the Transfer Flat.

Avoid over saturation and ensure that the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting. Extended exposure to these chemicals may damage the device.

WARNING

SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products described above are used to clean Stryker EMS equipment, measures must be taken to ensure the products are wiped with clean water and air dried thoroughly following cleaning.

REMOVAL OF IODINE COMPOUNDS

Use a solution of 1/2 Tablespoon Sodium Thiosulfate in a pint of warm water to clean the stained area. Clean as soon as possible after staining occurs. If stains are not immediately removed, allow solution to soak or stand on the surface. Rinse surfaces which have been exposed to the solution in clear water before returning unit to service.

WARNING

Failure to properly clean or dispose of contaminated Transfer Flats increases the risk of exposure to bloodborne pathogens and may cause injury to the patient or the operator.

Note: Staining or discoloration may occur from constant cleaning. Always inspect the Transfer Flat prior to use to ensure proper functionality.

Warranty

Stryker EMS, a division of the Stryker Corporation, offers the following warranty option in the United States:

One (1) year parts. Stryker EMS warrants to the original purchaser that its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective.

The Transfer Flat is designed for a 5 year expected service life under normal use conditions.

If Stryker requests, products or parts for which an original purchaser makes a warranty claim, the purchaser shall return the product or part prepaid freight to Stryker's factory.

Any improper use or alteration or repair by unauthorized service providers in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

This statement constitutes Stryker EMS's entire warranty with respect to the aforesaid equipment. STRYKER MAKES NO OTHER WARRANTY OR REPRESENTATION EITHER EXPRESSED OR IMPLIED, EXCEPT AS SET FORTH HEREIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL STRYKER BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

Warranty

English

STRYKER EMS RETURN POLICY

Cots, Stair Chairs, Evacuation Chairs, Evacuation Sleds, Cot Fasteners, Transfer Flats, and Aftermarket Accessories may be returned up to 180 days of receipt if they meet the following guidelines:

Prior to 30 Days

- 30 day money back guarantee in effect
- Stryker EMS is responsible for all charges
- Returns will not be approved on modified items

Prior to 90 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 10% restocking fee

Prior to 180 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 25% restocking fee

RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. **SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.**

DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full. Claims for any short shipment must be made within thirty (30) days of invoice.

INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

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