Altrix® Temperature Management Hose

Operations Manual  
REF 8003-002-007  
REF 8003-002-107

R<sub>x</sub> ONLY  
CE  
EN NL  
CS NO  
DA PL  
DE PT  
ES RU  
FR SV  
IT TR  
KO ZH

8003-009-200 Rev AA.0 2023-11

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Altrix**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

Stryker EMEA Supply Chain  
Services B. V.  
Frans Maasweg 2  
Venlo 5928 SB, The  
Netherlands

2023-11 8003-009-200 Rev AA.0 WCR:

EN

Symbols

|                     |                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                     | Refer to instruction manual/booklet                                                                                                                                                                                                                                                                                                                                      |
|                     | Consult instructions for use                                                                                                                                                                                                                                                                                                                                             |
|                     | General warning                                                                                                                                                                                                                                                                                                                                                          |
|                     | Caution                                                                                                                                                                                                                                                                                                                                                                  |
| REF                 | Catalogue number                                                                                                                                                                                                                                                                                                                                                         |
| LOT                 | Batch code                                                                                                                                                                                                                                                                                                                                                               |
| CE                  | CE mark                                                                                                                                                                                                                                                                                                                                                                  |
| UK CA               | UK Conformity Assessment mark                                                                                                                                                                                                                                                                                                                                            |
|                     | Importer                                                                                                                                                                                                                                                                                                                                                                 |
| UDI                 | Unique device identifier                                                                                                                                                                                                                                                                                                                                                 |
|                     | Manufacturer                                                                                                                                                                                                                                                                                                                                                             |
|                     | Date of manufacture                                                                                                                                                                                                                                                                                                                                                      |
| CH REP              | Authorized representative in Switzerland                                                                                                                                                                                                                                                                                                                                 |
| EC REP              | Authorized representative in the European Community                                                                                                                                                                                                                                                                                                                      |
| MD                  | European medical device                                                                                                                                                                                                                                                                                                                                                  |
| US Patents          | For US Patents see <a href="http://www.stryker.com/patents">www.stryker.com/patents</a>                                                                                                                                                                                                                                                                                  |
| R <sub>x</sub> ONLY | <b>CAUTION</b> - US Federal Law restricts this device to sale by or on the order of a licensed physician.                                                                                                                                                                                                                                                                |
|                     | In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling. |
|                     | Temperature controller                                                                                                                                                                                                                                                                                                                                                   |

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The **Altrix** Temperature Management Hose provides an interface from the thermal transfer devices that connect to the **Altrix** Precision Temperature Management System (referred to as the controller).

Expected life

This **Altrix** Temperature Management Hose has a one year use expected life under normal use conditions.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

|                         |                                                                    |                                                                                    |
|-------------------------|--------------------------------------------------------------------|------------------------------------------------------------------------------------|
| Patient facing material | Silcryn specialty thermoplastic elastomer                          |                                                                                    |
| Length                  | 10 ft                                                              | 3.05 m                                                                             |
| Compatible controller   | Model 8001 <b>Altrix</b> ® Precision Temperature Management System | See Model 8001 Operations Manual for cleaning, disinfection, and more information. |

Stryker reserves the right to change specifications without notice.

European REACH

In accordance with the European REACH regulation and other environmental regulatory requirements the components that contain declarable substances are listed.

| Description                               | Number                       | Substance of very high concern (SVHC) chemical name |
|-------------------------------------------|------------------------------|-----------------------------------------------------|
| <b>Altrix</b> Temperature Management Hose | 8003-002-007<br>8003-002-107 | bis (2-ethylhexyl) phthalate                        |

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

**Note** - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the lot batch code (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Setup

- Connect the hose (A) to the controller hose connection port (Figure 1).


Figure 1 –

- Connect the hose end (B) to the thermal transfer device (C).

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