Model 6507 Power-PRO 2 post-crash inspection instructions

Follow these instructions for cots with the Stryker X-restraint package (6500-001-430) or the G-rated restraint package (6500-002-030).

Tools required:

- | | | |
|-------------------|---------------------------------|--------------------|
| • T30 Torx driver | • Camera | • Flat screwdriver |
| • 3/4" hex wrench | • Weight (approximately 250 lb) | • 5/8" hex wrench |

Procedure:

1. Remove the product from the vehicle to perform the inspection.

Note - Inspect the product for general damage (such as misalignment, broken or missing parts, cracks, and loose fasteners) as you perform the inspection.

2. Make sure that the cot wheels roll smoothly and the casters are tight yet rotate (Figure 1 and Figure 2). A crushed caster horn indicates that the product was involved in a major crash.

Figure 1 –

Figure 2 –

3. Using a flat screwdriver, remove the brake/steer pedals at both ends of the product. Inspect for cracks or deformation.
4. Using a T30 Torx driver, remove the caster mount covers. Inspect the caster mount castings for cracks or deformation (Figure 3).
5. Make sure that the rail clamp post is secure and the base tube is not damaged (Figure 4).

Figure 3 –

Figure 4 –

6. Inspect the base tube welds for cracks (Figure 5).
7. Make sure that the head end base tube is not bent (Figure 6). Superficial damage to the outer plastic cover is acceptable. A bent base tube indicates that the product was involved in a major crash.

Figure 5 –

Figure 6 –

8. Inspect for local damage (crushing or a notch) at the middle of the outer plastic cross tube cover.
9. Apply an approximate 250 lb weight to the cot to inspect the X-frame. Check the castings that interface with the base and litter for cracks or damage (Figure 7 and Figure 8). Cracked castings indicate that the product was involved in a major crash.

Figure 7 –

Figure 8 –

10. Inspect the welds on the X-frame legs for cracks (Figure 9 and Figure 10).

Figure 9 –


Figure 10 –

11. Inspect the inner length of the base tube legs for dents.
12. Support both ends of the litter and remove the weight from the cot.
13. Using a 3/4" hex wrench and a 5/8" hex wrench, detach the rod end of the hydraulic cylinder from the base (Figure 11). Save the rod end.

14. Manually raise and lower the base to make sure that the X-frame moves through its range of motion.
15. Using a 3/4" hex wrench and a 5/8" hex wrench, reattach the hydraulic cylinder (detached in step 13). Check for damage to the foot end lift bars.
16. Inspect the cross tubes at both ends of the hydraulic cylinder for damage (Figure 11 and Figure 12).

Figure 11 –

Figure 12 –

17. Make sure that you can raise and lower the cot smoothly.
 - a. Support both ends of the litter and make sure that the **high speed retract** feature functions.
 - b. Check the adjustment of the manual release cable.
 - c. Operate the cot with both sets of switches.
 18. Inspect the foot end hitch hooks for damage (Figure 13). Squeeze the hooks together and make sure that the spring returns the hooks to their original position.
 19. Check the electronics enclosure for damage.
 - a. Make sure that the battery slides in and out and the latch mechanism locks the battery into the enclosure.
 - b. Check for contact between the enclosure and the foot end hitch.
- Note** - The battery indicator light should turn green when you insert a fully charged battery.
20. Inspect the underside of the patient left and patient right slider block housings for cracking (Figure 14).


Figure 13 –

Figure 14 –

21. Check for deformation of the outer rails at the restraint attachment locations (Figure 15). Deformation of the outer rail(s) indicates that the product was involved in a major crash.

Note - A telescoping head and foot section is a good indication of deformation.
22. Inspect the motor mounts for deformation or damage. The motor mount attachment to the V-tube should be centered on the V. Pay attention to deformation at the cylinder pivot (Figure 16). Deformation of the motor mounts, which prevent the cot from locking into the **Performance-LOAD** fastener, indicates that the product was involved in a major crash.

Figure 15 –

Figure 16 –

23. Make sure that the head section telescopes in and out. The head section should lock in the fully extended and fully retracted positions. Check for damage to the head section release handles and guards (Figure 17).

24. Inspect the load wheel horns. Check for cracks or damage, the horns are tight, and the wheels roll (Figure 18).

Figure 17 –

Figure 18 –

25. Make sure that the safety bar rotates and returns to the neutral position.

26. Raise and lower the Fowler to make sure that it moves smoothly and the Fowler cylinder release is adjusted (Figure 19). Check the Fowler cylinder for leaks or damage.

27. Inspect the litter pivots for cracks or damage (Figure 20).

Figure 19 –

Figure 20 –

28. Check the Gatch/Trendelenburg attachment for cracks or damage (Figure 21).

29. Inspect the litter pivots for cracks or damage (Figure 22).

Figure 21 –

Figure 22 –

30. Make sure that the Gatch/Trendelenburg move smoothly and that the release mechanism functions.

31. Replace the cot restraints (Figure 23).

Figure 23 –

32. If the cot is equipped with XPS, make sure that the wings pivot and the ratcheting mechanisms lock and release in all seven positions.

33. If the cot has standard siderails, check for cracks and that the siderail moves, locks, and unlocks.

34. Make sure that **Power-PRO 2** locks into and releases from its intended fastener.

35. Verify proper operation before you return the product to service. Follow all preventive maintenance procedures in the Operations/Maintenance Manual.

36. Complete the attached *Model 6507 Power-PRO 2 post-crash inspection report* (page 6) and send to ambulanceaccidents@stryker.com. Make sure that you include all requested pictures.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**, **Performance-LOAD**, **Power-LOAD**, **Power-PRO**. All other trademarks are trademarks of their respective owners or holders.

stryker

3800 E. Centre Avenue, Portage, MI 49002 USA
1-800-327-0770 | 2025-01

Model 6507 Power-PRO 2 post-crash inspection report

Perform this inspection report for the following cot configuration: **Power-LOAD** and **Performance-LOAD** compatible, XPS, XPS mattress, X-restraints, Gatch, and single (patient right) wheel lock.

Inspector name:	Inspection date:
Fastener serial number:	Service report number:

Make sure that you provide all required pictures and document details thoroughly.

Cot serial number Note - Send a picture of your product's serial number with this report (Figure 24).	 Figure 24 – Serial number example picture
---	---

Crash event details

Crash damage indicators - inspect the following parts for damage	
Crushed cot caster horns Crushed caster horns are more likely to occur after a side crash. Note - Send a picture of all four caster horns with this report (Figure 25).	 Figure 25 – Caster horn example picture
Details:	
Dented inner base tube legs Inspect the inner length of the base tube legs for dents. Note - Send a picture of the inner base tube legs with this report (Figure 26).	 Figure 26 – Inner base tube legs example picture

Crash damage indicators - inspect the following parts for damage

Details:

Damaged outer rail

- Check for a bent or crushed outer rail at the restraint attachment locations. Loosen or remove the restraints to inspect. The telescoping head section release should operate properly, and the head section should be able to extend and retract smoothly.

Note - Send a picture of the outer rail with this report (Figure 27).

Figure 27 – Outer rail example picture

Details:

Deformed motor mount

- Damage to the motor mount may appear on cots crashed with a **Performance-LOAD** fastener. Front crashes create damage under the seat section, near the motor. Side crashes can twist the motor mount near the Fowler cylinder pivot. Make sure that the attachment to the v-tube is centered.

Note - Send a picture of the motor mount with this report (Figure 28).

Figure 28 – Motor mount example picture

Details:

XPS locking mechanism

- Make sure that the XPS siderails lock and unlock properly. The XPS locking mechanism is designed to break away during a side crash.

Note - Send a picture of the XPS siderails with this report (Figure 29).

Figure 29 – XPS locking mechanism example picture

Details:

Crash damage indicators - inspect the following parts for damage

Cracked litter pivot casting

- Check for cracks on the castings that secure the inner legs to the litter at the v-tube. This damage occurs after front and side crashes with a **Performance-LOAD** fastener. If a casting cover is present, remove and inspect the casting.

Note - Send a picture of the litter pivot casting with this report (Figure 30).

Figure 30 – Litter pivot casting example picture

Details:

Bent head end base tube

- Check the cross tube at the head end side of the cot base for deformation. This damage may occur in a front crash with the **Power-LOAD** fastener. Damage to the plastic cover is not a failure; rotating the plastic cover is a good way to check if the tube is bent.

Note - Send a picture of the head end base tube with this report (Figure 31).

Figure 31 – Head end base tube example picture

Details:

Damaged hooks

- Inspect both patient left and patient right hooks for damage and dents along the front face and edges.

Note - Send a picture of the hook(s) with this report (Figure 32).

Figure 32 – Damaged hook example picture

Crash damage indicators - inspect the following parts for damage

Details:

Provide details in the space below to describe other damage (exclude normal wear):

Note

- Record any damage found during the inspection procedures. Please complete an additional service report for the repair of these items.
 - If any of the damage indicators are present, Stryker will recommend that you remove the product from service.
- ☐ Check the box if an additional service report is required. Service report number: _____
- ☐ Check the box if the other damage described indicates that the product should be recommended for removal from service.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Stryker**, **Performance-LOAD**, **Power-LOAD**, **Power-PRO**. All other trademarks are trademarks of their respective owners or holders.

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