# Model 3009 ProCuity NAWI scale verification of conformity instructions

For individuals authorized by their region's legal metrology only.

**Note** - Recommended procedure for field verification of conformity does not supersede local requirements.

## Tools required:

- (1) 5 kg calibrated weight
- (1) 10 kg calibrated weight
- (12) 20 kg calibrated weight
- Hole punch
- T25 Torx driver
- T30 Torx driver
- Ambient temperature sensor
- (1) Scale calibration counter label (300900370004)
- Verification sticker

**Note** - If you cannot complete scale verification, the product must remain out of service until you resolve the scale issue and complete scale verification.

3009 ProCuity NAWI scale specifications	
Accuracy class	(III)
Min	20 kg
Max	250 kg
e	2 kg
Verification period	1 year
Tare type	Subtractive
Maximum tare	60 kg

## Access the service menu

1. Lift the footboard until the screen turns off.
2. Restart the footboard.
3. When **ProCuity** appears on the screen, press and hold **ProCuity** for five seconds.
4. Release **ProCuity**, tap the bottom left corner of the screen, and tap the bottom right corner of the screen.


Figure 1 – Service menu screen

## NAWI test mode

1. Access the service menu (page 1).
2. Select Diagnostic Info.
3. Select Scale.
4. Select the information icon in the upper right corner.
5. Record the weight values rounded to the nearest tenth.


Figure 2 – NAWI test mode screen

## Verification instructions

1. Record the following in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17) or equivalent:
  - Operator name or ID
  - Product serial number
  - Time and date
  - Temperature and humidity (humidity optional)
  - Scale software version

**Note** - Access the scale information on the ProCuity service menu screen.

2. Confirm that the product complies with the Type Approval Certificate.
3. Confirm that the zone of operation is set to the correct country.

- a. Access the service menu (page 1).
- b. Select Configuration→Bed Configuration→Country.
4. Confirm that the product specification label (A) at the head end of the product is readable and not damaged (Figure 3).
5. Confirm that the NAWI certification label (B) is readable and not damaged.

**Figure 3 – Labels on head of the product**

6. Push down on the brake pedal to apply the brake.
7. Lower all siderails to the lowest height position.
8. Remove the support surface from the litter deck.
9. Raise the product to the highest height position (Figure 4).


**Figure 4 – Product at highest height position**

10. Raise the Fowler backrest to the highest height position.
11. Raise the Gatch to the highest height position and fold the foot section toward the head end of the product (Figure 5).


**Figure 5 – Raised Gatch and folded foot section**

12. Unplug the power cord from the wall outlet.
13. Turn off the battery switch at the head end of the product to turn the product off.



**Figure 6 – Battery switch**

14. Using a T25 Torx driver, remove the two screws (BJ) that secure the litter electronics box cover (E) to the electronics box assembly (Figure 7). Save the screws.

**Figure 7 – Removing the cover**

15. Using a T30 Torx driver, remove the two screws (AN) that secure the litter electronics box cover (E) to the electronics box assembly (Figure 7). Save the screws and cover.
16. Confirm that the anti-tamper label (D) is readable and not damaged on the main control board and the head end load cell locations (Figure 8).

**Figure 8 – Anti-tamper label**

17. Using a T25 Torx driver, reinstall the two screws (BJ) to secure the litter electronics box cover (E) to the electronics box assembly (Figure 7).
18. Using a T30 Torx driver, reinstall the two screws (AN) to secure the litter electronics box cover (E) to the electronics box assembly (Figure 7).
19. Extend the bed extender to the out and locked position. Fold the bed extender pan toward the foot end of the product (Figure 9).


**Figure 9 – Bed extender out and locked**

**Figure 10 – Foot end anti-tamper labels**

20. Confirm that the anti-tamper labels (F) are readable and not damaged on the foot end load cell locations (Figure 10).
21. Lower the bed extender litter skin and retract the bed extender litter to the in and locked position.
22. Plug the power cord into the wall outlet.
23. Turn on the battery switch at the head end of the product to turn the product on.
24. Lower the foot section into the normal use position.
25. Lower the Gatch to the flat position.
26. Lower the Fowler backrest to the flat position.
27. Using calibrated weights, pre-load the product to 250 kg (Figure 11).

**Figure 11 – Pre-load - 250 kg**

28. Remove all calibrated weights from the product (Figure 12).


Figure 12 – Pre-load - no weight

## NAWI scale calibration

1. Access scale maintenance in the service menu to start the NAWI scale calibration process.
2. Select Configuration→Calibration→Scale.
3. Follow the on-screen prompts to calibrate the scale. When prompted to select load cells replaced, if you did not replace any load cells press continue. Perform the factory zero reset as the last step.  
**Note** - Contact Stryker Customer Service or Technical Support at 1-800-327-0770 to receive the controlled NAWI code.
4. Remove all calibrated and uncalibrated weight sources from the product.

## Eccentricity test

Record the values in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17). Take all values from *NAWI test mode* (page 2).

1. Using calibrated weights, place 85 kg at the head, left of the product (Figure 13).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 13 – Eccentricity test - 85 kg head left

2. Using calibrated weights, place 85 kg at the foot, left of the product (Figure 14).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 14 – Eccentricity test - 85 kg foot left

3. Using calibrated weights, place 85 kg at the foot, right of the product (Figure 15).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 15 – Eccentricity test - 85 kg foot right

4. Using calibrated weights, place 85 kg at the head, right of the product (Figure 16).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 16 – Eccentricity test - 85 kg head right

5. Remove all calibrated weights from the product (Figure 17).


Figure 17 – Eccentricity test - no weight

## Weighing performance test

Record the values in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17). Take all values from *NAWI test mode* (page 2).

1. Using calibrated weights, place 20 kg on the product (Figure 18).

**Note** - The weight indication must be  $\pm 1.0$  kg.



Figure 18 – Weigh performance test - 20 kg

2. Using calibrated weights, place 60 kg on the product (Figure 19).


**Note** - The weight indication must be  $\pm 1.0$  kg.

**Figure 19 – Weigh performance test - 60 kg**

3. Using calibrated weights, place 100 kg on the product (Figure 20).


**Note** - The weight indication must be  $\pm 1.0$  kg.

**Figure 20 – Weigh performance test - 100 kg**

4. Using calibrated weights, place 160 kg on the product (Figure 21).


**Note** - The weight indication must be  $\pm 2.0$  kg.

**Figure 21 – Weigh performance test - 160 kg**

5. Using calibrated weights, place 200 kg on the product (Figure 22).


**Note** - The weight indication must be  $\pm 2.0$  kg.

**Figure 22 – Weigh performance test - 200 kg**

6. Using calibrated weights, place 250 kg on the product (Figure 23).

**Note** - The weight indication must be  $\pm 2.0$  kg.

**Figure 23 – Weigh performance test - 250 kg**

7. Using calibrated weights, place 200 kg on the product (Figure 22).

**Note** - The weight indication must be  $\pm 2.0$  kg.

8. Using calibrated weights, place 160 kg on the product (Figure 21).

**Note** - The weight indication must be  $\pm 2.0$  kg.

9. Using calibrated weights, place 100 kg on the product (Figure 20).

**Note** - The weight indication must be  $\pm 1.0$  kg.

10. Using calibrated weights, place 60 kg on the product (Figure 19).

**Note** - The weight indication must be  $\pm 1.0$  kg.

11. Using calibrated weights, place 20 kg on the product (Figure 18).

**Note** - The weight indication must be  $\pm 1.0$  kg.

12. Remove all calibrated weights from the product (Figure 24).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 24 – Weight performance test - 0 kg

## Repeatability test

The absolute value for the range of repeatability must not exceed 2.0 kg for tests 1, 2, and 3.

Record the values in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17). Take all values from *NAWI test mode* (page 2).

1. Using calibrated weights, place 200 kg on the product (Figure 25).

**Note** - The weight indication must be  $\pm 2.0$  kg.



Figure 25 – Repeatability test - 200 kg

2. Remove all calibrated weights from the product (Figure 26).

**Figure 26 – Repeatability test - no weight**


3. Repeat steps 1-2 two more times for a total of three readings.
4. Subtract the minimum value recorded from the maximum value recorded. Record the results.

**Note** - The repeatability range must be  $\leq 2.0$  kg.

## Tare accuracy

Record the value in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17). Take all values from *NAWI test mode* (page 2).

1. Using calibrated weights, place 40 kg on the product (Figure 27).

**Figure 27 – Tare test - 0 kg net**

2. Tare the scale with 40 kg on the product.
- Note** - The accuracy of tare must be  $0.0$  kg  $\pm 0.5$  kg.
3. Using calibrated weights, place 20 kg on the product (Figure 28).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 28 – Tare test - 20 kg net

4. Using calibrated weights, place 60 kg on the product (Figure 29).

**Note** - The weight indication must be  $\pm 1.0$  kg.


Figure 29 – Tare test - 60 kg net

5. Using calibrated weights, place 100 kg on the product (Figure 30).

**Note** - The weight indication must be  $\pm 1.0$  kg.



Figure 30 – Tare test - 100 kg net

6. Using calibrated weights, place 160 kg on the product (Figure 31).


**Note** - The weight indication must be  $\pm 2.0$  kg.

**Figure 31 – Tare test - 160 kg net**

7. Using calibrated weights, place 210 kg on the product (Figure 32).

**Note** - The weight indication must be  $\pm 2.0$  kg.

**Figure 32 – Tare test - 210 kg net**

8. Using calibrated weights, place 160 kg on the product (Figure 31).

**Note** - The weight indication must be  $\pm 2.0$  kg.

9. Using calibrated weights, place 100 kg on the product (Figure 30).

**Note** - The weight indication must be  $\pm 1.0$  kg.

10. Using calibrated weights, place 60 kg on the product (Figure 29).

**Note** - The weight indication must be  $\pm 1.0$  kg.

11. Using calibrated weights, place 20 kg on the product (Figure 28).

**Note** - The weight indication must be  $\pm 1.0$  kg.

12. Using calibrated weights, place 0 kg on the product (Figure 27).

**Note** - The weight indication must be  $\pm 1.0$  kg.

13. Remove all calibrated weights from the product (Figure 33).


Figure 33 – Tare test - no weight

## Post test

1. Record the following values in the *Model 3009 ProCuity NAWI scale verification of conformity form* (page 17) in the post-calibration section.
  - a. Time and date
  - b. Temperature

**Note** - The temperature must not change  $\pm 5^{\circ}\text{C}$  from the start of the procedure.
2. Record the calibration counter.
  - a. Access the service menu (page 1).
  - b. Select **Diagnostic Information→Scales**.
3. Record the country.
  - a. Access the service menu (page 1).
  - b. Select **Configuration→Bed Configuration→Country**.
4. Remove the current scale calibration counter label (Figure 34).


Figure 34 – Calibration counter label

5. Using a hole punch, punch the calibration counter and zone on the new scale calibration counter label (300900370004) to match the scale calibration menu.
6. Affix the punched scale calibration counter label to the product (Figure 34).
7. Affix the verification sticker next to the NAWI certification label as prescribed by the local legal metrology.

# Model 3009 ProCuity NAWI scale verification of conformity form

Use this form to record data from the Model 3009 ProCuity NAWI scale verification of conformity instructions. Record the weight values rounded to the nearest tenth.

**Note** - If you cannot complete scale verification, the product must remain out of service until the scale issue is resolved and scale verification is complete.

Operator name or ID	
Product serial number	
Time and date (pre-calibration)	
Temperature (pre-calibration)	
Humidity (pre-calibration)	
<b>Note</b> - Humidity optional	
Scale software version	
Calibrated weight IDs	

Eccentricity test		
Operation	Recorded value	Acceptance criteria
Head, left		85.0 kg ± 1.0 kg
Foot, left		85.0 kg ± 1.0 kg
Foot, right		85.0 kg ± 1.0 kg
Head, right		85.0 kg ± 1.0 kg

Weighing performance test		
Operation/scale performance	Recorded value	Acceptance criteria
20 kg		20.0 kg ± 1.0 kg
60 kg		60.0 kg ± 1.0 kg
100 kg		100.0 kg ± 1.0 kg
160 kg		160.0 kg ± 2.0 kg
200 kg		200.0 kg ± 2.0 kg
250 kg		250.0 kg ± 2.0 kg
200 kg		200.0 kg ± 2.0 kg
160 kg		160.0 kg ± 2.0 kg
100 kg		100.0 kg ± 1.0 kg
60 kg		60.0 kg ± 1.0 kg
20 kg		20.0 kg ± 1.0 kg
0 kg		0.0 kg ± 1.0 kg

Repeatability test		
Operation	Recorded value	Acceptance criteria
Test 1: 200 kg		200.0 kg ± 2.0 kg
Test 2: 200 kg		200.0 kg ± 2.0 kg
Test 3: 200 kg		200.0 kg ± 2.0 kg
Repeatability range (maximum value - minimum value)		≤ 2.0 kg

Tare accuracy (40 kg subtractive tare)		
Operation/tare accuracy	Recorded value	Acceptance criteria
0 kg load		0.0 kg ± 0.5 kg
20 kg load		20.0 kg ± 1.0 kg
60 kg load		60.0 kg ± 1.0 kg
100 kg load		100.0 kg ± 1.0 kg
160 kg load		160.0 kg ± 2.0 kg
210 kg load		210.0 kg ± 2.0 kg
160 kg load		160.0 kg ± 2.0 kg
100 kg load		100.0 kg ± 1.0 kg
60 kg load		60.0 kg ± 1.0 kg
20 kg load		20.0 kg ± 1.0 kg
0 kg load		0.0 kg ± 1.0 kg

Post test		
Post-calibration	Recorded value	
Time and date		
Temperature		
Calibration counter		
Calibration country		

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Eksantriklik testi		
Kullanma	Kaydedilen değer	Kabul kriterleri
Baş, sol		85,0 kg ± 1,0 kg
Ayak, sol		85,0 kg ± 1,0 kg
Ayak, sağ		85,0 kg ± 1,0 kg
Baş, sağ		85,0 kg ± 1,0 kg

Tartma performans testi		
Kullanma/tartı performansı	Kaydedilen değer	Kabul kriterleri
20 kg		20,0 kg ± 1,0 kg
60 kg		60,0 kg ± 1,0 kg
100 kg		100,0 kg ± 1,0 kg
160 kg		160,0 kg ± 2,0 kg
200 kg		200,0 kg ± 2,0 kg
250 kg		250,0 kg ± 2,0 kg
200 kg		200,0 kg ± 2,0 kg
160 kg		160,0 kg ± 2,0 kg
100 kg		100,0 kg ± 1,0 kg
60 kg		60,0 kg ± 1,0 kg
20 kg		20,0 kg ± 1,0 kg
0 kg		0,0 kg ± 1,0 kg

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# stryker

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **ProCuity, Stryker**. All other trademarks are trademarks of their respective owners or holders.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA