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English

CF302 Series Inflation Unit Instructions for Use

Figure 1



Figure 2

CONTRAINDICATIONS

Air support therapy is not recommended when spinal stability is a concern.

INDICATIONS

The Sof-Care® Care-Free® Companion Pump (CF302) Series is for use with only Gaymar Sof-Care® and duo-gard® overlays. These products are designed to assist in the prevention and treatment of pressure ulcers.

The Sof-Care® Care-Free® Companion Pump (CF302) unit will inflate and maintain Sof-Care® overlays at therapeutic pressure levels, even when minor punctures occur. The Sof-Care® overlay redistributes the patient's weight over the surface and effectively reduces tissue interface pressure.

SAFETY PRECAUTIONS

DANGER

If a bulge in the top surface of the overlay is noticed, discontinue use and call your dealer.

A bulge may cause the patient to roll over when not desired, potentially resulting in serious injury or death.


CAUTION

- Check patient's skin regularly. Consult physician if any redness or skin break occurs. Serious injury could result if the patient's skin condition is left untreated.
- Do not place the pump in the patient's bed, in contact with the patient, or under sheets or other coverings. This could cause serious injury or could affect pump performance.

INSTRUCTIONS

1. Set up overlay and pump as shown in Figure 1. Place corner straps around bed mattress.
2. Place the pump near the bed. If placing it on footboard, place the unit on the outside (not the patient side), with the rubber feet against the footboard as shown in Figure 2.
3. Turn unit on using illuminated on/off switch. Unit will take approximately 40 minutes to inflate the overlay.
4. Place the patient on the overlay and set the comfort control knob to the appropriate setting based on patient comfort.
5. Perform Hand Check.
A Hand Check must be performed every 8 hours to verify proper operation.
To perform the Hand Check:
With the patient on his or her back, slide hand, flat and palm up, between the overlay and the mattress. Hand should be directly under the air cell that is under the patient's buttocks (or other bony area).


HAND CHECK

If the patient's body is not in direct contact with hand, the system is operating correctly. The patient's body is in direct contact with flat hand, the system is not operating properly. Adjust the pressure control to a higher setting. Wait 10

minutes and repeat the Hand Check. If the Hand Check fails, check that the hoses are not kinked or pinched. If repeated Hand Check fails and hoses are not kinked, contact Dealer or Gaymar for further instruction.

CPR DEFLATE.

WARNING

Deflate before CPR or CPR will be ineffective.

CLEANING

The pump housing should be cleaned between patients.

Before cleaning, unplug the pump. Clean housing with a mild detergent and damp cloth. Wipe pump dry with a clean cloth before operating. Do not submerge. Do not autoclave.

TECHNICAL DATA

Voltage: 230 VAC

Current: 0.05A

Wattage: 10 Watts

Overcurrent Protection: Fused, Rated 1A 250 VAC

Frequency: 50 Hz

Operating Temperature: 15°C to 32°C, (60°F to 90°F)

Classification:

- EN 60601-1, EN 60601-1-2
- IPX0, enclosed equipment without protection against ingress of water
- Not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen (O₂) or nitrous oxide (N₂O).
- Continuous operation

The control unit is warranted for one (1) year.

SYMBOLS

Caution, consult accompanying documents

Dangerous voltage

Type BF equipment

Double Insulation

Consult Operating Instructions

Disposal: Contact local distributor who will take the necessary steps according to your national market.

CUSTOMER SERVICE

If you have purchased this pump through a local dealer and are experiencing problems, please contact that local dealer for assistance. For additional assistance, contact our Customer Service Department at:

Gaymar Industries, Inc.
10 Centre Drive
Orchard Park, NY 14127-2295 USA
Telephone: (716) 662-8636
FAX: (716) 662-0730

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©2009. Gaymar Industries, Inc.
Sof-Care and Duo-Gard are registered trademarks of Gaymar
Stryker France S.A.S.