SofCare®
Stretcher Overlay
Recouvrement de brancard
Revestimento de maca
担架覆蓋垫
ストレッチャーパッド

REF 2780-400-000
 2780-500-000

Operations Manual
Manuel d'utilisation
Manual de utilização
使用手册
操作マニュアル

ENGLISH (EN)
 FRANÇAIS (FR)
 PORTUGUÊS (PT)
 简体中文 (ZH-CN)
 日本語 (JA)

2014/12 B.0 2780-209-002 Rev B www.stryker.com

Introduction

This manual assists you with the operation of the Stryker Model 2780-400-000 and 2780-500-000 SofCare® Stretcher Overlay. Read this manual before operating this product and keep a copy on file. Set methods and procedures to educate and train your staff on the safe operation of this product.

Intended use

The SofCare Overlay assists in the prevention and treatment of existing pressure ulcers and is for single patient use. The product is for use by trained healthcare professionals or family members. Use this product in an acute or long term care or home healthcare settings. Use settings also include critical care, step-down, med/surg, sub-acute, and post anesthesia care unit (PACU).

Product description

The SofCare Overlay combines three layers of air cells to create a contoured support surface. An inflated SofCare Overlay redistributes the pressure on the patient's skin when placed on a mattress. The SofCare Overlay is inflatable with an optional pump.

Expected life

SofCare Overlay has a 30 day single patient use expected life under normal use conditions.

Contraindications

- Air support therapy is not recommended when spinal stability is a concern.
- This product is not intended to support a patient in the prone position.

Specifications

Model	2780-400-000		2780-500-000	
Length	75 in.	190,5 cm	75 in.	190,5 cm
Width	27 in.	68,6 cm	30 in.	76,2 cm
Therapeutic working load	100 lb to 350 lb (45,3 kg to 158,8 kg)			
Safe working load	350 lb (158,8 kg)			
Material	Vynyl			
Operating ambient temperature range and humidity limits				
Storage temperature and humidity limits				
Optional pumps	Companion pump 2.0 (2781-100-000)		Hand Inflation Pump (2781-200-000)	

Symbols

	General warning		Caution		Operating instructions
	Foot end		Batch code		Temperature limitation
	Safe working load		Humidity limitation		No sharp objects
	Model number		Manufacturer		Plastic #3- Polyvinyl Chloride (PVC or V)

Stryker reserves the right to change specifications without notice.

Operation

Setting up the product

- Remove the product from the plastic bag. Place the product across the stretcher with the hose connection on the right side of a supine patient.
- Standing at the foot end of the stretcher, unfold the cushion until you see the air valve on the right side.
- Place the top strap around the top lower section of the stretcher (Figure 2). Secure both side straps to the stretcher of the tie down points as (Figure 3).


Figure 2: Strap location for top straps

Figure 3: Strap location for side straps

Connecting an optional non-electric pump to the product

- Attach the screw cap nozzle of the hand pump to the air valve on the cushion.
- Use approximately 50 complete strokes to inflate the Overlay (1 complete stroke = 1 in and out motion).
- Remove the hand pump. In a rapid movement, replace the blue cap until sealed.
- Check the pressure setting by hand.

Note: Overinflation may damage the overlay.

Adding or changing liners

- Apply the liners loose to avoid a hammock effect.
- Use minimal layers of sheeting and incontinence pads.
- Change the liners per your hospital procedures.

Maintaining patient comfort and protection

The Overlay has an under inflated appearance and may not need inflation. The best way to determine proper inflation is to check by hand.

Checking pressure setting by hand

Check the pressure setting by hand to make sure that the patient receives the proper therapy. Perform a hand check with the patient in the worst case position, which is the side-lying position. This position maximizes the pressure at the sacrum.

Figure 4: Under inflation

Figure 5: Proper inflation

- Slide your hand, palm up with fingers flat, between the product and the mattress at the patient's sacrum.
 - Do not lean on the product or lift at the side as these actions can lead to false readings.
 - Adjust the pressure setpoint until you can no longer feel the patient's sacrum.
- If you can feel the bony prominence (Figure 4), inflate the product.
- Wait two minutes and repeat the hand check until the you have adequate inflation (Figure 5).

Note: Check the patient at regular intervals as directed by hospital protocol.

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WARNING/CAUTION/NOTE DEFINITION

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note: Provides special information to make maintenance easier or important instructions clearer.

WARNING

- Do not use this product with the bed siderails down. Siderails are to remind the patient of the edge of the unit. Do not use the siderails as a patient restraining device.
- Do not reuse this product on another patient to avoid the risk of cross-contamination and infection.
- Always deflate the product before performing CPR.
- Do not use this product as a transfer device.
- Always keep the bedding dry. Failure to keep the patient's skin dry may lead to maceration.
- Do not continue to use the product if a bulge in the top surface is present as this may reduce the reliability of the product function.
- Always secure the product to the mattress using the corner straps.
- Always inspect the product for tears, cuts, holes, or any other damage before use. Always check the patient's skin at regular intervals for redness or breaks. If redness or breaks in the skin occur, consult the attending physician.
- Always consult the pump operations manual for use before connecting the system to this product.

CAUTION

- Always use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the product lowers the systems ability to aid in pressure reduction.
- Do not pull liners too tight over the product to avoid the hammock effect and reduce the effectiveness of the product.
- Always latch the siderails when the bed is in a raised position.
- Always be aware that the use of a mattress overlay may reduce the effectiveness and coverage of a siderail.

Contact information

Contact Stryker Customer Service or Technical Support at 1-800-327-0770.

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Have the Lot batch code (E) (Figure 1) available when you call Customer Service. Include the Lot batch code in all written communication.

Deflating for CPR

- Unscrew the hose connector cap (Figure 6) from the product.
- For slow or partial deflation, insert a pen point or other device into one of the four holes in the valve disc (Figure 7). Air will escape at a slow pace.
- For fast or complete deflation, pull the tab to remove the disc from the air valve (Figure 8).
- Proceed per your hospital CPR procedures.


Figure 6: Unscrew cap

Figure 7: Puncture

Figure 8: Remove disc

Cleaning

- Clean the product with a clean damp cloth, soap, and water.
- Wipe dry.
- Dispose of product according to your local waste management policy.

Notes:

- Disinfectant cleaners will degrade the product. Do not use disinfectants, such as solvents, alcohol based products, or Virex TB germicidal products.
- Iodine based disinfectants will stain the product. Staining does not affect the patient nor the use of the product.

Warranty

Stryker Medical, a division of Stryker Corporation ("Stryker") warrants that its models 2780-400-000 and 2780-500-000 SofCare Stretcher Overlay products will be free from defects in manufacturing and workmanship for a maximum of 30 day