This manual assists you with the operation of the Stryker Model 8001-061-701 and 8001-061-702 **MuT-Blanket®**. Read this manual before operating this product and keep a copy on file. Set methods and procedures to educate and train your staff on the safe operation of this product.

Intended use and Contraindications

Always consult the controller operations manual before using this product.

Product description

The **MuT-Blanket** is a lightweight polymer. The **Click-Tite** connectors provide the means to connect to the temperature controller.

Expected life
These **MuT-Blanket** products have a 90 day reusable expected life upon first use and adhere to normal use conditions.

Specifications

Material	vinyl, polypropylene	
Model	8001-061-701	8001-061-702
Patient	adult	pediatric
Length	69 in. 175 cm	36 in. 91 cm
Width	25 in. 63 cm	25 in. 63 cm
Compatible controllers	Medi-Therm series temperature controller	
Environmental conditions	Operation	Storage and transportation
Temperature		
Relative humidity		

Stryker reserves the right to change specifications without notice.

In accordance with the European REACH regulation, this article contains bis(2-ethylhexyl) phthalate.

Symbol	Meaning
	General warning
	Refer to instruction manual/ booklet
	No sharp objects
	Batch code
	Quantity
	Date of manufacture
	Authorized representative in the European Community
	CAUTION US Federal Law restricts this device to sale by or on the order of a licensed physician.

English (EN), Español (ES), Deutsch (DE), Français (FR), Italiano (IT), Nederlands (NL), Svenska (SV), Dansk (DA), Português (PT), Norsk (NN, Pøisi (PL), Русский (RU), Türkçe (TR), Český (CS), 中文 (ZH-CN), 한국어 (KO)

Contact Information

Contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

To view your product manual online, see <https://techweb.stryker.com>.

How the Lot batch code (A) (Figure 1) available when you call Customer Service. Include the Lot batch code in all written communication.

Safety Precautions

- WARNING**
 - Always check the product for leaks, tears, cuts, holes, stains, or any other damage prior to and during use.
 - Always check the patient's skin condition in areas of contact with the blanket according to hospital protocol or as directed by a physician.
 - Always turn or reposition the patient over the duration of therapy according to hospital protocol.
 - Always keep the area between the patient and the thermal transfer device dry.
 - Application of heat or cold are most effective when you apply the body wrap directly to the skin.
 - Do not allow hoses to fold, kink, or to wrap around a patient.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

CAUTION

- Do not use if the package is damaged.
- Always use **MuT-Blanket** products with an approved temperature controller. Using devices with other temperature controllers has not been tested.
- Do not use **quaternaries** that contain glycol ethers as they may damage the product.
- Do not use phenolic cleaners as they may damage the product.
- Do not use this product to position or transfer a patient.
- Do not use this product to position or transfer a patient.
- Always pre-fill the product before applying to the patient to avoid the risk of pressure ulcers.
- Do not use three or more adult **MuT-Blanket's** at the same time to avoid the risk of water overflow when you power off the controller.
- Always consult the controller operations manual before connecting this product.
- Do not use sharp objects or pins with this product.
- Always use a dry absorbent sheet between the patient and the plastic side of any thermal transfer device.
- Do not block fluid pathways of the thermal transfer device with restraints.
- Always use minimal layers of sheeting. Too many layers between the patient's skin and the blanket will reduce the cooling or warming capabilities of the system.

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Recheque el paciente's skin condition

- Check the patient at regular intervals as directed by hospital protocol. Note any change in the skin integrity that relates to:
 - Excessive moisture - dry the skin surface by wiping away the moisture
 - Color of the epidermis
 - Skin texture
 - Patient's skin condition is acceptable to continue therapy

Check the blanket

- Recheck the blanket for the following:
 - position
 - water status
 - dry surface
 - leaks
 - cracks


Desinfection considerations

- Remove the blanket to expose patient's chest.
- Remove excess moisture.

Desconectar the blanket from the hoses

- Power off the controller. Allow most of the water to drain back into the machine. (Approximately 10 minutes).
- Close blanket clamps (Figure 2).
- Close the connector hose clamps.
- Unsnap the fitting. Pull the male coupling out to disconnect (Figure 3).
- Disconnect the power cord into a wall outlet.
- Discard blanket after use per facility protocol.

Note: Do not discard the controller connector hose after use, save for next use.

Cleaning

- Clean the product with a clean damp cloth, soap, and water.
- Wipe dry.

Disinfecting

- Suggested disinfectants:
- Chlorinated bleach solution (1 part bleach solution (5.25% sodium hypochlorite) to 10 parts water (10% bleach cleaner))
 - 70% Isopropyl alcohol

To disinfect the product after each patient use, follow these steps in order:

- Clean and dry the product (see Cleaning) before disinfectants are applied.
- Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes (do not soak the product).

Note: Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.

- Rinse the product with a clean, damp cloth to remove excess disinfectant.
- Wipe the product with a clean, dry cloth to remove any excess liquid or disinfectant.
- Allow the product to air dry before returning to service.

Disposing of the product

Dispose of the product in accordance with your local waste management policy.

Warranty

Stryker Medical, a division of Stryker Corporation ("Stryker") warrants that its models 8001-061-701 and 8001-061-702 **MuT-Blanket** products will be free from defects in manufacturing and workmanship for a maximum of 90 day use period. "Use" period begins when such devices individual package is opened. Stryker's obligation under this warranty is expressly limited to supplying a product replacement, at its option, any product which is, in the sole discretion of Stryker, found to be defective. If requested by Stryker, products for which a claim is made shall be returned prepaid to the factory. Any improper use or any alterations or repair by others in such manner as in Stryker's judgment affects the product's reliability and adversely shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

Warranty exclusion and damage limitations
This express warranty set forth herein is the only warranty applicable to the product. Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by Stryker. In no event shall Stryker be liable for incidental or consequential damages.

Return authorization

Product cannot be returned without prior approval from the Stryker Customer Service Department and without a return authorization number. Will be provided on the date of the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged merchandise

ICC Regulations require that claims for damaged product must be made with fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

International warranty clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for extra information.

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