Rapr•Round®

Hyper/Hypothermia Vest and Leg Wrap

REF DHV530 Vest (S/M)

REF DHV535 Vest (L)

REF DHL540 Leg Wrap

Instructions for Use

Read and understand the operating instructions and all safety precautions prior to use.

Symbols

Attention, consult accompanying documents

Instructions for Use

Single patient use

Do not puncture

Do not use if package is damaged

LATEX FREE

Latex free

SAFETY PRECAUTIONS

WARNINGS

The Rapr•Round body wraps are used with a device that generates heat. Excessive heat could cause thermal stress or skin lesions. Failure to follow these precautions could result in death or serious injury:

- US Federal Law restricts this device to sale by or on the order of a licensed health care professional.
- The body wraps have not been tested with other manufacturer's machines. The vest and leg wraps are only for use with Stryker Medi-Therm® series Hyper/Hypothermia machines.
- If the patient's temperature is not responding as expected, notify the attending physician promptly.
- Failure to monitor core temperature could result in abnormal elevation or reduction of body temperature resulting in serious injury or death.
- Do not apply heat to the lower extremities during aortic cross-clamping. Ischemic injury may occur.
- Do not use the body wraps to position or transfer patients.
- Keep hoses free and clear of patient. Do not allow hoses to become folded, kinked or wrapped around the patient.
- To insure proper water flow through the body wrap, follow the Flow Check Procedure.
- To prevent patient injury, follow hospital/departmental protocol to secure patient on table/bed.
- Patients who have peripheral vascular disease, inadequate heart function, reduction in blood volume or who are temperature-sensitive, may deviate from the normal response to the external application of heat and cold. This patient population needs to be monitored closely for signs of intolerance.
- Electrosurgical pads may be used when conductive warming body wraps are applied. In the unlikely event that the body wrap may leak, do not place electrosurgical pads between the patient and the conductive body wrap.
- Prep solutions pooling beneath the patient may cause chemical injury.
- Check the integrity of the skin according to unit protocol when regulating temperature with external devices. Frequency of assessment and documentation will vary depending upon the individual response of the patient.
- Risk of cross-contamination (e.g. infection, allergic reaction) if product were to be reused.

NOTE

- Application of heat and cold is most effective when body wraps are applied directly to patient's skin.
- Following application of the vest and/or leg wrap(s), you may see the pattern of these body wraps on the patient's skin. This pattern on the skin may be reddened (hyperemic). This occurrence is normal and will disappear within the same length of time that the wraps had been applied. This phenomenon is often mistaken for injury but is in fact reactive hyperemia.^{1,2}

References:


- 1 Stewart, T.P., Magnano, S.J. Burns or Pressure Ulcers in the Surgical Patient? Decubitus 1988 1: 3640
- 2 Stewart, T.P. When a burn is not a burn. International Journal of Dermatology 1998 37: 729732

Sizing

The leg wrap, DHL540, is available in one size. The vest is available in two sizes. The DHV535 (large—L) vest is designed to fit a chest size from 117cm (46") to 137cm (54"). The DHV530 (small/medium—S/M) vest is for chest sizes from 81cm (32") to 117cm (46").

General Instructions


Figure 1

1. Although not necessary, body wraps may be pre-filled prior to patient application. Adjust temperature as necessary.
2. Apply wraps directly to skin.
3. Place the side marked "GAYMAR" away from the patient's skin. (When vest is installed properly, the hoses will be on the patient's left side. See Figure 1.
4. Fasten wraps closed with Velcro-like closure.
5. Do not use pins or sharp objects with the body wraps.
6. Patients should not be lying on hoses.
7. Hoses should not be wrapped around patient.
8. Any folds or creases may restrict or prevent the flow of water.
9. Leg Wrap: Make sure the cutout portion of the leg wrap is located over the knee.
10. After patient positioning and/or repositioning, perform the Flow Check Procedure.

Flow Check Procedure

Figure 2

When using two gray hose sets (see Figure 2), one connected to a vest and one connected to the leg wrap(s), check for blocked water flow as follows:

1. Disconnect one gray hose set from the machine connection. See Figure 2. If CHECK WATER FLOW light illuminates, check for kinks, folds or closed clamps in the other body wraps connected to the machine.
2. Reconnect hose set to machine.
3. Disconnect second hose from the machine. If CHECK WATER FLOW light illuminates, check for kinks, folds or closed clamps in the other body wraps connected to the machine.
4. Reconnect hose set to the machine.

NOTE When using only one gray hose set, connected to only the vest or leg wrap(s), it is not necessary to follow this procedure. The CHECK WATER FLOW light will illuminate if the flow is interrupted.

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Symboles

Attention, consultez les documents joints

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Stryker France S.A.S.