stryker[®]

Operations manual

stryker®

Operations Manual

Symbols

English
EN

	General warning
	Caution
	Warning; crushing of hands
	This product is equipped with a hospital-grade attachment (receptacle) plug. You can only achieve grounding reliability when you use a hospital-grade receptacle.
	Warning: Non-protectively earthed, potential for risk of electric shock
	In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return or collection systems available in your country.
	Safe working load
	Maximum patient weight
	Alternating current
	Dangerous voltage
	Type B applied part
IPX4	Protection from liquid splash
	Protective earth terminal
	Product provides terminal for connection of a potential equalization conductor. The potential equalization conductor provides direct connection between the product and potential equalization busbar of the electrical installation.

Symbols

English
EN

	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1: 2005 and CANCSA-C22.2 No. 60601-1:08.
	Caution, electrostatic sensitive
	Catalogue number
	Serial number
	Manufacturer
	Date of manufacture

Table of Contents

English
EN

Warning/Caution/Note Definition	1-3
Summary of safety precautions	1-4
Introduction	1-6
Product description	1-6
Intended use	1-6
Expected service life	1-6
Contraindications	1-7
Specifications	1-7
Product illustration	1-9
Contact information	1-9
Serial number location	1-10
Setup	1-11
Equipping the optional nurse call communication	1-11
Operation	1-13
Applying or releasing the brakes	1-13
Applying or releasing the steer pedal	1-14
Shutting down the Birthing Bed/emergency stop	1-14
Storing the power cord	1-15
Transporting the Birthing Bed	1-15
Positioning the Lock-Rite foot section	1-16
Positioning the patient on the Birthing Bed	1-16
Operating the optional nurse call communication	1-17
Activating the CPR release	1-19
Raising or lowering the siderails	1-19
Positioning the labor grips	1-21
Adjusting the optional Skoocher system	1-21
Adjusting for Trendelenburg Positioning	1-22
Adjusting for Pelvic Tilt	1-22
Adjusting the optional lumbar pillow	1-22
Illuminating the room with the night light	1-22
Positioning the Glideaway footrests	1-22
Adjusting the Glideaway Footrests	1-23
Attaching or removing the optional calf supports	1-25
Adjusting the Glideaway footrests with optional calf supports	1-26
Removing the Lock-Rite foot section	1-28
Attaching the Lock-Rite foot section	1-28
Attaching the Mattress	1-29
Removing the fluid basin	1-31
Positioning the Birthing Bed manually with the manual hand crank	1-31
Removing or replacing the headboard	1-32
Installing the optional labor bar	1-32
Operating the three-stage permanently attached IV pole	1-33
Nurse control panel (Outside siderail)	1-34

Table of Contents

English
EN

Patient control panel (Inside siderail).....	1-36
Patient control panel with optional Smart TV (Inside siderail).....	1-37
Patient lumbar control functions with optional Smart TV and optional Air Mattress (Inside siderail).....	1-38
Motion pendant functions (optional) - 5000-081-021.....	1-39
Motion pendant with nurse call functions (optional) - 5000-081-023.....	1-40
Motion and communication pendant functions (optional) - 5000-081-020.....	1-41
Preventive maintenance	1-43
Cleaning and disinfecting with wipes	1-44
Cleaning.....	1-45
Cleaning the mattress	1-45
EMC Information	1-47

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note: Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

English
EN

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Only operate the product when all operators are clear of the mechanisms.
 - Do not remove, sit, or place any load on the cover panels to avoid the risk of possible electrical shock hazard and product damage. Service should only be performed by qualified personnel.
 - An operator should not touch connectors and the patient simultaneously to avoid the risk of patient or operator injury.
 - Do not use the product in the presence of flammable anesthetics.
 - Always keep feet clear from the area above the base cover or below the base cover when lowering the product or when applying the brakes or releasing the brakes.
 - Always apply the brakes when a patient is getting in the product or out of the product to avoid instability.
 - Always apply the brakes when the patient is unattended.
 - Always allow enough clearance between the head end of the product and the adjacent wall, so you can easily unplug the power cord from the wall outlet during an emergency.
 - Always properly handle the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, immediately remove the product from service and contact the appropriate maintenance personnel.
 - Always unplug the product power cord from the wall when using oxygen administering equipment. Possible fire hazard exists when this product is used with oxygen administering equipment other than nasal, mask type, or half bed-length tent type.
 - Always store the power cord before transporting the Birthing Bed.
 - Do not use the foot section locking bar as a grasping bar or other patient assist device.
 - Always leave the product in the lowest position when the patient is unattended to help reduce the number and severity of patient falls.
 - Always make sure that all persons and equipment are away from the area below and around the fowler before you activate the CPR release. The CPR release is for emergency use only.
 - Always set the siderail position necessary to make sure that a patient is safely in the product.
 - Always lock the nurse control panel and patient control panel when the patient is unattended.
 - Always keep the siderails outside of the oxygen tent.
 - Always keep hands clear of potential pinch points while positioning the labor grips.
 - Always lower and lock the locking bar before placing weight on the foot section.
 - Do not extend the oxygen tent below the mattress platform.
 - Do not clean, service, or perform maintenance while the product is in use.
 - When large spills occur in the area of the circuit boards, 120VAC cables, and motors, always immediately unplug the power cord from the wall outlet. Remove the patient from the product, clean up the fluid, and have service personnel completely inspect the product. Fluids can cause unpredictable operation and decreased functionality of any electrical product. Do not return the product to service until it is completely dry and has been thoroughly tested for safe operation.
 - Always completely dry the mattress after cleaning. Failure to thoroughly rinse and dry mattresses after cleaning may cause damage to the mattress and may void this product's warranty.
-

Summary of safety precautions

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always plug the product directly into a properly grounded, three-prong receptacle. You can only achieve grounding reliability when you use a hospital-grade receptacle. This product is equipped with a hospital-grade plug for protection against electric shock hazard.
 - Perform preventive maintenance annually to make sure that all product features are functioning properly.
 - Maintain the cleanliness and integrity of both ground chains to minimize static build-up and discharge.
 - Always apply the brake to prevent unintended movement.
 - Do not apply the brake pedal to stop a moving product.
 - The maximum safe working load for each calf support is 86 lb (39 kg).
 - The maximum safe working load of the foot section is 300 lb.
 - Only use compatible mattresses that are authorized by Stryker.
 - Inspect the mattress cover after each use. Openings in the mattress cover may allow fluids to enter the mattress. Discontinue using the mattress and safely dispose the mattress if any rips or cracks are found in the mattress cover.
 - Always disconnect the power cord when you use the manual hand crank.
 - Do not hang IV bags that exceed 40 lb onto the IV pole.
 - Always make sure that the IV pole is at a low height to allow it to pass safely through door openings.
 - Do not use the IV pole as a push/pull device.
 - Always place the pendant safely on the mattress while the pendant is in use.
 - Do not squeeze or pinch cables in the bed frame when routing cables from other equipment.
 - Always unplug the product before cleaning or servicing.
 - Always unplug the product, set the brakes, and place blocks under the litter frame for support when working under the Birthing Bed.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature corrosion of critical components. Failure to follow these cleaning instructions may void your warranty.
 - Do not steam clean, pressure wash, ultrasonically clean, or immerse any part of the product in water. Exposure to water may damage the internal electric parts. These methods of cleaning are not recommended and may void this product's warranty.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Notes

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Birthing Bed is an AC-powered, adjustable obstetrical hospital bed designed to aid in the birthing process. The Birthing Bed provides comfort for both the patient and operator with an ergonomic design. The Birthing Bed offers a trend-like position, a nurse call feature, labor bar, calf supports, and foot supports to assist in the labor and delivery process. The Birthing Bed offers the **Skoocher**® option to allow for customized positioning and assists in moving the patient closer to the perineal edge of the product.

Intended use

The Birthing Bed is intended to be used by patients during the labor and delivery process, and those requiring intra-hospital transfers to and from delivery rooms, in a health care environment. The safe working load (the sum of the patient, the mattress, and accessory weight) for the Birthing Bed is 500 lb (227 kg).

The Birthing Bed is intended to support a human patient during the labor and delivery process. The Birthing Bed (frame and mattress) can come in contact with the patient's skin.

The Birthing Bed is intended to be used in an acute health care environment, including hospitals and other institutions where obstetric procedures are provided for labor and delivery. Intended operators of the Birthing Bed are health care professionals and maintenance technicians (for maintenance related issues). Patients and bystanders can operate limited functions, which shall be limited to bed motion features that are related to comfort only and the nurse call feature.

The Birthing Bed is compatible with a mattress with the following dimensions:

- Head end mattress – 49.5 in. (126 cm) x 33 in. (84 cm) x 5 in. (13 cm)
- Foot section mattress – 30 in. (76 cm) x 30 in. (76 cm) x 3 in. (8 cm)

Expected service life

The Birthing Bed has a 10 year expected service life under normal use, conditions, and with appropriate periodic maintenance.

The Birthing Bed mattress has a two year expected service life under normal use, conditions, and with appropriate periodic maintenance.

Contraindications

The Birthing Bed is not intended to:

- use its siderails as a patient restraint device
- be used without a mattress
- support more than one individual at a time
- be used on patients less than two years old
- be a sterile product
- be used for pediatric use

Specifications

	Safe working load Note: Safe working load indicates the sum of the patient, mattress, and accessory weight	500 lb	227 kg
	Maximum patient weight	396 lb	180 kg
Weight of product	Standard configuration	475 lb	215 kg
	All options/accessories	525 lb	238 kg
Product length/width	Siderails up	93 in. x 41 in.	236 cm x 104 cm
	Siderails down	93 in. x 37 in.	236 cm x 94 cm
Product height (to top of seat) (with 6 in. casters)	Low	17.5 in.	44 cm
	High	35.5 in.	90 cm
Mattress size	Head end	49.5 in. x 33 in. x 5 in.	126 cm x 84 cm x 13 cm
	Foot section	30 in. x 30 in. x 3 in.	76 cm x 76 cm x 8 cm
Foot section Travel height	Up	0 in.	0 cm
	Down	7 in.	17.8 cm
Break-away point from wall		60 in.	152 cm
Litter position		Maximum Elevation: Head 70°, Trend-Like 8°	
Caster size		6 in. Standard; 8 in. Optional	
Fluid basin		7 quarts	7 liters
Electrical	Standard 4 motor function: Head-Unit-Foot-Trend-Like 120VAC, 60 Hz, 10 Amp Current leakage less than 300 microamperes (per UL 60601-1) Hospital grade plug and 3-wire heavy duty cord Compatible with non-flammable anesthetic agents and oxygen by nasal catheter or mask type		







Introduction

English
EN

Specifications (Continued)

	Over voltage category — II (per IEC 60601-1) 120VAC - 1500V
Circuit breaker, Thermal	250VAC, 10 Amp (0059-196-000)
Illumination	Range: 100 lx to 1500 lx
Battery voltage	9VDC — Nurse Call Option (Alkaline Battery)
Rated duty cycle	3 minutes / 57 minutes — 3 minutes ON and 57 minutes OFF

Class I Equipment: Equipment that protects against electrical shock and does not solely rely on basic insulation, but which includes an additional safety precaution that is provided for the connection of the equipment to the protective earth conductor in the fixed wiring of the installation that accessible metal parts cannot become live in the event of a failure of basic insulation.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (Non-Condensing)		
Atmospheric pressure		

Specifications listed are approximate and may vary slightly from product to product or by power supply fluctuations.

Stryker reserves the right to change specifications without notice.

Product illustration


Figure 1-1: Birthing Bed

A	Brake/Steer pedal	K	Head end mattress
B	Calf support (optional)	L	Labor grip
C	Caster	M	Patient lumbar control functions with Smart TV and Air Mattress (optional)
D	CPR release	N	Perineal Edge
E	Fluid basin	O	Lock-Rite™ foot section
F	Footrest	P	Siderail
G	Foot section locking bar	Q	Nurse control panel (Outside siderail)
H	Foot section mattress	R	Patient control panel (Inside siderail)
I	Foot section release handle	S	Seat section
J	Headboard	T	Skoocher system (optional)

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Introduction

English
EN

Contact information (Continued)

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location


Figure 1-2: Serial Number Location

Setup

To unpack your product, see the unpacking instructions that are attached to the product inside of the shipping crate.

WARNING

Only operate the product when all operators are clear of the mechanisms.

CAUTION

Always plug the product directly into a properly grounded, three-prong receptacle. You can only achieve grounding reliability when you use a hospital-grade receptacle. This product is equipped with a hospital-grade plug for protection against electric shock hazard.

Before placing the product into service, make sure that these components are working properly:

1. Plug the power cord into a properly grounded hospital-grade wall receptacle.
 2. Depress the brake pedal to make sure that the brakes hold.
 3. Raise and lower the siderails to make sure that they move smoothly and lock into position.
 4. Press each button on the inside and outside siderail control panels to make sure that each function is working properly.
 5. Press each button on the optional hand pendant to make sure that each function is working properly.
-

WARNING

- Do not remove, sit, or place any load on the cover panels to avoid the risk of possible electrical shock hazard and product damage. Service should only be performed by qualified personnel.
 - An operator should not touch connectors and the patient simultaneously to avoid the risk of patient or operator injury.
 - Do not use the product in the presence of flammable anesthetics.
-

CAUTION

- Perform preventive maintenance annually to make sure that all product features are functioning properly.
 - Maintain the cleanliness and integrity of both ground chains to minimize static build-up and discharge.
-

Equipping the optional nurse call communication

If your product is equipped with the optional nurse call communication, follow the setup instructions below before putting the product into service.

1. Plug the interface cable into the 37-pin connector in the litter frame at the head end of the product (B) ([Figure 1-3 on page 1-12](#)) and into the applicable connection (patient station, head wall, or docker station).
Note: Only connect the 37-pin connector to the Head Wall Output Configuration (B) ([Figure 1-3 on page 1-12](#)) or product Communications Tester (sold separately).
2. Test the interface cable to verify connectivity.
3. Plug the power cord into the wall outlet.
4. Push in the pendant port switch and turn the switch to the **ON** position (90° clockwise) (A) ([Figure 1-3 on page 1-12](#)).
5. Press the **Nurse call** button (see [Nurse control panel \(Outside siderail\) on page 1-34](#)) to verify the connection between the product's nurse call signal and the hospital's nurse call system.

Note: A 9V battery, located in the head end of the product, powers the nurse call signal. If the product is unplugged from a wall outlet while the nurse call pendant port switch is set to the **ON** position, the charge of the 9V battery will begin to drain.

Setup

English
EN

Equipping the optional nurse call communication (Continued)

To operate the optional nurse call communication, see [Operating the optional nurse call communication on page 1-17](#).


Figure 1-3: Nurse call pendant port and 37-pin connector

Applying or releasing the brakes

WARNING

- Always keep feet clear from the area above the base cover or below the base cover when lowering the product or when applying the brakes or releasing the brakes.
- Always apply the brakes when a patient is getting in the product or out of the product to avoid instability.
- Always apply the brakes when the patient is unattended.

CAUTION

- Always apply the brake to prevent unintended movement.
- Do not apply the brake pedal to stop a moving product.

The brake pedals are located on both the left and right sides of the product.


Figure 1-4: Brake pedal (Left side shown)

To apply the brakes, fully depress the pedal toward the head end. The brake pedal locks all four caster wheel brakes and holds the product in place.


Figure 1-5: Apply the brakes

To release the brakes, depress the pedal toward the foot end until the pedal is in the neutral position. This releases all four caster wheels and allows you to freely move the product forward, backward, and from side to side.


Figure 1-6: Release the brakes

Applying or releasing the steer pedal

The steer pedals are located on both the left and right sides of the product.


Figure 1-7: Steer pedal (Left side shown)

To apply the steer pedal, fully depress the pedal toward the foot end. The steer pedal locks the right side caster on the foot end and allows you to easily maneuver the product forward and backward in a straight path.


Figure 1-8: Apply the steer pedal

To release the steer pedal, depress the pedal toward the head end until the pedal is in the neutral position. This releases the right side caster on the foot end and allows you to freely move the product forward, backward, and from side to side.


Figure 1-9: Release the steer pedal

Shutting down the Birthing Bed/emergency stop

WARNING

Always allow enough clearance between the head end of the product and the adjacent wall, so you can easily unplug the power cord from the wall outlet during an emergency.

To shut down the Birthing Bed:

1. Unplug the power cord from the wall outlet.
2. Store the power cord (see [Storing the power cord on page 1-15](#)).

Storing the power cord

WARNING

- Always properly handle the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, immediately remove the product from service and contact the appropriate maintenance personnel.
- Always unplug the product power cord from the wall when using oxygen administering equipment. Possible fire hazard exists when this product is used with oxygen administering equipment other than nasal, mask type, or half bed-length tent type.
- Always store the power cord before transporting the Birthing Bed.

To store the power cord, wrap the power cord around the roller bumpers that are located between the headboard and the IV pole at the head end of the product.


Figure 1-10: Storing the power cord

Transporting the Birthing Bed

WARNING

Always store the power cord before transporting the Birthing Bed.

To transport the Birthing Bed:

1. Press the **Siderail Bed Motion Control Lockout** button to lock all siderail bed motion controls (see [Nurse control panel \(Outside siderail\) on page 1-34](#)).
2. Unplug the power cord from the wall outlet.
3. Store the power cord (see [Storing the power cord on page 1-15](#)).
4. Lower the IV pole to allow it to safely pass through door openings.
5. Raise and lock the siderails in the highest height (see [Raising or lowering the siderails on page 1-19](#)).
6. Depress the pedal toward the foot end until the pedal is in the neutral position to release the brakes (see [Applying or releasing the brakes on page 1-13](#)).
7. Push the product from the head end.

Operation

English
EN

Positioning the Lock-Rite foot section

WARNING

Do not use the foot section locking bar as a grasping bar or other patient assist device.

To position the Lock-Rite foot section:

1. Position the patient's perineum out and over the edge of the perineal edge.
2. Lower the foot section to its lowest position.
3. Place and position the patient's feet on the footrests or legs in the optional calf supports.

Positioning the patient on the Birthing Bed

WARNING

Always leave the product in the lowest position when the patient is unattended to help reduce the number and severity of patient falls.

CAUTION

The maximum safe working load of the foot section is 300 lb.

To position the patient on the Birthing Bed:

1. Place and position the patient's feet on the footrests.
2. Raise or lower the footrests to a position that is comfortable for the patient.
3. Remove the foot section.
4. Press the **Litter Up** button (see [Nurse control panel \(Outside siderail\) on page 1-34](#)) to raise the bed to a comfortable height. Release the button when you achieve the desired height.
5. Position the patient's perineum out and over the edge of the perineal edge.

Other possible patient positions:

Upright Sitting

You can position the footrests or optional calf supports for the patient with the foot section still intact. This allows the operator to provide hands-free care from the foot end of the bed.

Squatting

By supporting themselves on the labor bar, the patient can fully stretch the pelvis and vaginal opening, allowing the baby to descend.

Resting/Squatting

The use of a high semi-fowler and of the Skoocher system allows for a comfortable resting position.

Positioning the patient on the Birthing Bed (Continued)

McRobert's Maneuver

For expected shoulder dystocia, the adjustable fowler and length-adjustable footrests facilitate pelvic opening. The operator's hands are free to apply suprapubic pressure.

Lithotomy

The optional calf supports and optional Skoocher system help position the patient for better visibility and access, whether for an episiotomy, forcep application, or if the patient is anesthetized.

Side Lying

The optional calf supports are accessible and easy to position for appropriate flexion. The optional calf supports can hold the patient's leg to increase pelvic stretch.

Knee/Chest

You can lower the footrest to increase back stretch, help relieve back labor, and aid in fetal rotation.

Kneeling, Leaning Forward

Adjusting the Skoocher system and the upright fowler position puts the patient in an increased pelvic tilt, allowing for better fetal descent and rotation in the posterior position.

Operating the optional nurse call communication

Nurse call allows the patient to send a signal to the nurse station.

To activate nurse call, press the **Nurse call** button (see [Nurse control panel \(Outside siderail\) on page 1-34](#)).

Communication between the patient and the nurse station is established at the moment when the nursing staff responds to the nurse call signal.

Note: A 9V battery, located in the head end of the product, powers the nurse call signal. If the product is unplugged from a wall outlet while the nurse call pendant port switch is set to the **ON** position, the charge of the 9V battery will begin to drain.

If the Power LED is blinking, you need to replace the 9V nurse call battery. See the Birthing Bed Maintenance Manual for replacement procedures.

Operation

English
EN

Operating the optional nurse call communication (Continued)


Figure 1-11: Nurse call switch

Activating the CPR release

WARNING

Always make sure that all persons and equipment are away from the area below and around the fowler before you activate the CPR release. The CPR release is for emergency use only.

When the fowler is raised and quick access to the patient is needed, you can position the Birthing Bed quickly to 0° by activating the CPR release.

The instant CPR release lever is located at the head end section on both the left and right sides of the fowler (Figure 1-12 on page 1-19).

To activate the CPR release, grasp the lever on either side of the fowler and squeeze tightly. The fowler will instantly descend to the lowest position.


Figure 1-12: CPR release

Raising or lowering the siderails

WARNING

- Always set the siderail position necessary to make sure that a patient is safely in the product.
- Always lock the nurse control panel and patient control panel when the patient is unattended.
- Always keep the siderails outside of the oxygen tent.

Note: Siderails are not intended to be used as a patient restraint device.

Raising or lowering the siderails (Continued)

When raising the siderails, listen for the "click" that indicates that the siderail has locked in the raised position. Pull firmly on the siderail to make sure it is locked into position.

To raise the siderail to its highest position, rotate the siderail forward to the full up position (Figure 1-13 on page 1-20).

Note: The siderail does not lock into the intermediate position when you raise the siderail.

Figure 1-13: Siderail full up

To lower the siderail to the intermediate position, push in the release handle and rotate the siderail backward until it locks into the intermediate position (Figure 1-14 on page 1-20).

Figure 1-14: Siderail intermediate height

To lower the siderail to its lowest position, push in the release handle again and rotate the siderail backward until it is completely lowered (Figure 1-15 on page 1-20).

Note: Make sure that the siderail is in the lowest position before raising the siderail directly to the full up position. If you do not completely lower the siderail, the siderail will lock into the intermediate position when you raise the siderail.

Figure 1-15: Siderail lowered

Positioning the labor grips

WARNING

Always keep hands clear of potential pinch points while positioning the labor grips.

Labor grips offer the patient a point of leverage while in the Birthing Bed.

To raise the labor grips, pull the labeled lever forward and rotate the grip out and above the seat section until the labor grip locks into place ([Figure 1-16 on page 1-21](#)).

To lower the labor grips, pull the labeled lever forward to unlock the labor grip and then rotate the grip in and below the seat section until the labor grip locks into place ([Figure 1-17 on page 1-21](#)).

Figure 1-16: Labor grip up

Figure 1-17: Labor grip down

Adjusting the optional Skoocher system

Skoocher is a powered system that shifts the perineal edge toward the operator.

You can adjust the seat depth from 12.5 in. to 8.5 in. to accommodate patients of different sizes or to bring the patient closer to or farther from the operator. When the fowler is lowered, the seat depth will automatically adjust to 12.5 in.

To adjust the Skoocher:

1. Raise the fowler to a minimum of 35°.
Note: When the fowler is elevated to 35° or greater, the Skoocher LED will illuminate on the siderail to indicate that Skoocher is active.
2. To move the patient toward the perineal edge and the operator ([Figure 1-18 on page 1-22](#)), press the **Skooch In** button to move the fowler toward the foot end of the product (see [Nurse control panel \(Outside siderail\) on page 1-34](#)).
3. To move the patient away from the perineal edge and the operator ([Figure 1-19 on page 1-22](#)), press the **Skooch Out** button to move the fowler toward the head end of the product (see [Nurse control panel \(Outside siderail\) on page 1-34](#)).

Adjusting the optional Skoocher system (Continued)

Figure 1-18: Skooch in

Figure 1-19: Skooch out

Adjusting for Trendelenburg Positioning

To adjust for Trendelenburg positioning, press the **Trend/Pelvic Tilt** button on the siderail (see [Nurse control panel \(Outside siderail\) on page 1-34](#)) or on the optional pendant (see [Motion pendant functions \(optional\) - 5000-081-021 on page 1-39](#)). Release the button when you achieve the desired angle.

Adjusting for Pelvic Tilt

Pelvic tilt provides additional tilt and comfort for the patient when you raise the fowler.

To adjust for pelvic tilt, press the **Trend/Pelvic Tilt** button on the siderail (see [Nurse control panel \(Outside siderail\) on page 1-34](#)) or on the optional pendant (see [Motion pendant functions \(optional\) - 5000-081-021 on page 1-39](#)). Release the button when you achieve the desired angle.

Adjusting the optional lumbar pillow

The inflatable, full-width lumbar pillow supports the patient for back labor, and the seat pillow enhances perineal visibility.

To adjust the lumbar pillow, see [Patient lumbar control functions with optional Smart TV and optional Air Mattress \(Inside siderail\) on page 1-38](#).

Note: The lumbar controls are only compatible with the Air mattress.

Illuminating the room with the night light

The Birthing Bed is equipped with a night light that illuminates the floor area. The night light automatically activates when the light in the room becomes dim.

Positioning the Glideaway footrests

To position the Glideaway footrest:

1. Grasp the footrest and pull it up and around the foot section ([Figure 1-20 on page 1-23](#)).
2. Grasp and attempt to move the footrest to make sure that the footrest is locked in the upright position.

Positioning the Glideaway footrests (Continued)

3. To stow the footrest, pull the purple lever that is located at the bottom of the footrest mechanism (Figure 1-22 on page 1-23).
4. Swing the entire footrest around and under the foot section (Figure 1-23 on page 1-23).

Note: You cannot stow the footrest unless it is in the original upright position (Figure 1-21 on page 1-23).

Figure 1-20: Accessing the footrest

Figure 1-21: Original upright position

Figure 1-22: Purple lever


Figure 1-23: Stowing the footrest

Adjusting the Glideaway Footrests

Adjusting the Footrest Width

You can pivot the footrests in and out to adjust the patient's feet and legs to the desired width.

Adjusting the Glideaway Footrests (Continued)

To pivot the footrest:

1. Squeeze the gray release handles on the abductors that are located at the foot end of the Birthing Bed (Figure 1-25 on page 1-24).
2. Pivot the footrest to the desired position (Figure 1-26 on page 1-24).
3. Release the handle to lock the footrest in the desired position.

Notes

- The footrest will not lock into any position beyond 36°.
- For emergency situations, the footrests will pivot to 90°.

Figure 1-24: Abductor Location

Figure 1-25: Abductor Release Handle

Figure 1-26: Pivoted Footrests

Adjusting the Footrest Proximity

You can adjust the proximity of the footrests to accommodate the different heights of patients.

To adjust the proximity of the footrest:

1. Squeeze either of the two purple handles at the opening on the footrest mechanism (Figure 1-27 on page 1-25).
2. Slide the footrest to the desired position (Figure 1-28 on page 1-25).
3. Release the handle to lock the footrest in the desired position.

Adjusting the Glideaway Footrests (Continued)

Figure 1-27: Purple Handles


Figure 1-28: Footrest Proximity

Adjusting the Lock-Rite Foot Section Height

You can adjust the foot section up to 7 in. below the top of the perineal edge to accommodate different sized patients.

To adjust the height of the foot section, press the **Foot Section Down** button on the siderail (see [Nurse control panel \(Outside siderail\) on page 1-34](#)) or on the optional pendant (see [Motion pendant functions \(optional\) - 5000-081-021 on page 1-39](#)). Release the button when you achieve the desired height.


Figure 1-29: Adjusting the Lock-Rite Foot Section Height

Attaching or removing the optional calf supports

CAUTION

The maximum safe working load for each calf support is 86 lb (39 kg).

1. Grasp the footrest and pull it up and around the foot section ([Figure 1-20 on page 1-23](#)).
2. Grasp and attempt to move the footrest to make sure that the footrest is locked in the upright position.

Attaching or removing the optional calf supports (Continued)

3. Push the purple button and pivot the footrest into the full upright position (A) ([Figure 1-30 on page 1-26](#)).

To attach the calf support to the footrest:

1. Pull the pull knob (C) ([Figure 1-30 on page 1-26](#)).
2. While pulling the pull knob, center the screws on the calf support into the large holes in the footrest plate.
3. Lower the calf support into the footrest plate (B) ([Figure 1-30 on page 1-26](#)).
4. Release the pull knob into the footrest plate.

To remove the calf support from the footrest:

1. Pull the pull knob (C) ([Figure 1-30 on page 1-26](#)).
2. While pulling the pull knob, slide the screws on the calf support into the large holes in the footrest plate.
3. Lift the calf support off of the footrest plate (B) ([Figure 1-30 on page 1-26](#)).
4. Release the pull knob.


Figure 1-30: Attaching the calf supports

Adjusting the Glideaway footrests with optional calf supports

Note: You can make all footrest adjustments (see [Adjusting the Glideaway Footrests on page 1-23](#)) with the optional calf supports attached.

1. Grasp the footrest with the optional calf support attached and pull it up and around the foot section ([Figure 1-38 on page 1-27](#)).
2. Grasp and attempt to move the footrest to make sure that the footrest is locked in the upright position.
3. Push the purple button and pivot the footrest into the full upright position (A) ([Figure 1-30 on page 1-26](#)).

Adjusting the Glideaway footrests with optional calf supports (Continued)

To adjust the calf support:

1. Release the calf support by grasping the gray calf support handle and rotating it to the left (Figure 1-33 on page 1-27).
2. Adjust the calf support to the desired position.
3. Lock the calf support by rotating the gray calf support handle to the right (Figure 1-33 on page 1-27).

Figure 1-31: Stowed footrest with calf support

Figure 1-32: Accessing the calf support


Figure 1-33: Calf support handle

Figure 1-34: Positioned footrest with calf support

To stow the calf support:

1. Adjust the calf support into the divot housing (Figure 1-35 on page 1-27).
2. Lock the calf support in the divot housing by rotating the gray calf support handle to the right (Figure 1-33 on page 1-27).
3. Push the purple button to restore the footrest to its original position (Figure 1-21 on page 1-23).
Note: You cannot stow the footrest unless it is in the original upright position (Figure 1-21 on page 1-23).
4. To stow the footrest, pull the purple lever that is located at the bottom of the footrest mechanism (Figure 1-22 on page 1-23).
5. Swing the entire footrest around and under the foot section (Figure 1-31 on page 1-27).


Figure 1-35: Divot housing

Figure 1-36: Restoring the calf support


Figure 1-37: Restoring the footrest

Figure 1-38: Stowing the footrest with calf support

Operation

English
EN

Removing the Lock-Rite foot section

WARNING

Do not use the foot section locking bar as a grasping bar or other patient assist device.

CAUTION

The maximum safe working load of the foot section is 300 lb.

Note: Before removing the foot section, put the footrests into position above the foot section and place the patient's feet in the footrests (see [Positioning the patient on the Birthing Bed on page 1-16](#)).

To remove the foot section:

1. Walk into the curved end of the foot section.
2. Squeeze the foot section release handle (B) into the foot section locking bar (A) ([Figure 1-47 on page 1-29](#)) and raise both bars up until they latch into the full up position ([Figure 1-40 on page 1-28](#)).
3. Release the foot section release handle ([Figure 1-40 on page 1-28](#)).
4. While holding the foot section locking bar, leverage the foot section against your torso and slowly walk the foot section away from the Birthing Bed ([Figure 1-41 on page 1-28](#)).

Note: If the foot section is placed with the perineal edge on the floor, it will stand independently ([Figure 1-42 on page 1-28](#)).


Figure 1-39: Squeezing the foot section release handle and locking bar

Figure 1-40: Releasing the foot section release handle

Figure 1-41: Removing the foot section

Figure 1-42: Perineal Edge

Attaching the Lock-Rite foot section

To attach the foot section:

1. If the footrests are abducted, return them to their original positions ([Figure 1-43 on page 1-29](#)).
2. Grasp the foot section locking bar and hold the foot section against your torso.
3. Place the foot section on top of the abduction arms ([Figure 1-44 on page 1-29](#)).
4. While squeezing both the foot section release handle (B) and the foot section locking bar (A) ([Figure 1-47 on page 1-29](#)), leverage the foot section against your torso and slowly walk the foot section into the Birthing Bed ([Figure 1-45 on page 1-29](#)).
5. To lock the foot section onto the Birthing Bed, squeeze and lower both the foot section release handle and foot section locking bar into the locked position ([Figure 1-46 on page 1-29](#)). As both bars are lowered, the gap between the foot section and the Birthing Bed will narrow.

Note: You can only lower the foot section locking bar when the foot section is in the "locking zone", which is located approximately 2 in. from the head end mattress.

Operation

Attaching the Lock-Rite foot section (Continued)

6. Before placing weight on the foot section, make sure that the foot section is securely locked onto the Birthing Bed by grabbing the foot section locking bar and attempting to move the foot section.

WARNING

Always lower and lock the locking bar before placing weight on the foot section.

Figure 1-43: Returning abductors to original positions


Figure 1-44: Attaching the foot section

Figure 1-45: Locking the foot section


Figure 1-46: Locking the foot section onto the Birthing Bed


Figure 1-47: Foot section

Attaching the Mattress

CAUTION

Only use compatible mattresses that are authorized by Stryker.

To attach the foot section mattress:

1. Snap the foot section mattress to the edges of the Lock-Rite foot section.

Operation

English
EN

Attaching the Mattress (Continued)

Compatible Foot Section Mattress (Detachable Part):

4701-045-005	Foot Mattress
4701-045-015	Foot Mattress with Flame Barrier

Enhanced Comfort Mattress: This Non-Air mattress contains 5 in. of memory foam integrated with an egg crate topper.

To attach the Enhanced mattress to the Birthing Bed:

1. Insert the locating pins, located on the underside of the mattress, into the holes in the metal seat support.
2. Snap the mattress to the fowler of the Birthing Bed.

Compatible Enhanced Comfort Mattress (Detachable Part):

4701-045-000	Enhanced Comfort Mattress
4701-045-010	Enhanced Comfort Mattress with Flame Barrier

Air Mattress: The Air mattress contains two internal powered air pillows that the patient controls with the Patient Lumbar Control Functions. These air pillows can firm or soften the back and seat sections of the mattress for patient comfort (see [Patient lumbar control functions with optional Smart TV and optional Air Mattress \(Inside siderail\) on page 1-38](#)).

Note: The lumbar controls are only compatible with the Air mattress.

To attach the Air mattress to the Birthing Bed:

1. Insert the locating pins, located on the underside of the mattress, into the holes in the metal seat support.
2. Connect all the hoses from the Air mattress to the Birthing Bed.
3. Snap the mattress to the fowler of the Birthing Bed.

Compatible Air Mattress (Detachable Part):

4701-048-000	Air Mattress
4701-048-020	Air Mattress with Flame Barrier

WARNING

Do not extend the oxygen tent below the mattress platform.

CAUTION

Inspect the mattress cover after each use. Openings in the mattress cover may allow fluids to enter the mattress. Discontinue using the mattress and safely dispose the mattress if any rips or cracks are found in the mattress cover.

Removing the fluid basin

You can use standard Bio-Med disposal bags (red bags) or any other type of fluid basin liner with the fluid basin.

To remove the fluid basin:

1. Remove the Lock-Rite foot section (see [Removing the Lock-Rite foot section on page 1-28](#)).
2. Gently push the basin upward from the bottom ([Figure 1-48 on page 1-31](#)).
3. When the pegs on the product are even with the large circular holes of the fluid basin, slowly remove the basin from the pegs ([Figure 1-49 on page 1-31](#)).

Figure 1-48: Lifting the fluid basin

Figure 1-49: Removing the fluid basin

Positioning the Birthing Bed manually with the manual hand crank

CAUTION

Always disconnect the power cord when you use the manual hand crank.

A manual override is available to position the product when the power cord is not plugged into the wall socket.

To manually position the Birthing Bed, insert the manual hand crank into either of the sockets at the head end of the product and rotate the crank until you achieve the desired position ([Figure 1-50 on page 1-31](#)).


Figure 1-50: Manual hand crank locations

Removing or replacing the headboard

To remove the headboard, grasp the handles and lift the headboard straight up and off the product (Figure 1-51 on page 1-32).

To replace the headboard, align the plastic inserts on the bottom of the headboard with the slots at the head end of the product, and then lower the headboard until it completely seats in the slots (Figure 1-52 on page 1-32).


Figure 1-51: Removing the headboard


Figure 1-52: Replacing the headboard

Installing the optional labor bar

The patient can use the labor bar to aid with various birthing positions, such as squatting or kneeling. The optional labor bar is designed to support a 250 lb load.

To install the labor bar, insert the ends of the labor bar into the sockets that are located at the pivot point on top of the footrest abductors (Figure 1-53 on page 1-32).


Figure 1-53: Installing the labor bar


Figure 1-54: Labor bar

Operation

Operating the three-stage permanently attached IV pole

⚠ CAUTION

- Do not hang IV bags that exceed 40 lb onto the IV pole.
- Always make sure that the IV pole is at a low height to allow it to pass safely through door openings.
- Do not use the IV pole as a push/pull device.

Note: You can only install the three-stage permanently attached IV pole at the head end of the product.

To use the three-stage permanently attached IV pole:

1. Lift and pivot the pole from the storage position.
2. To raise the height of the pole, pull up on the telescoping portion (A) until it locks into place at its fully raised position (Figure 1-55 on page 1-33).
3. For a higher IV pole, pull up on section (B) (Figure 1-55 on page 1-33). Release section (B) at any desired height and it will lock into place.
4. Rotate the IV hangers (C) to the desired position and hang the IV bags (Figure 1-55 on page 1-33).
5. To lower the IV pole, push up on the black portion of the grip (D) while holding on to section (B) until it lowers (Figure 1-55 on page 1-33).
6. Turn the latch (E) until the telescoping portion (A) lowers (Figure 1-55 on page 1-33).



Figure 1-55: Three-Stage Permanently Attached IV Pole

Operation

English
EN

Nurse control panel (Outside siderail)

	Button name	Function
A	Fowler up	Raises the fowler
B	Fowler down	Lowers the fowler
C	Litter up	Raises the litter
D	Litter down	Lowers the litter
E	Foot section up	Raises the foot section
F	Foot section down	Lowers the foot section
<p>Note: When the product is at a low height and you press Foot Section Down, the litter may raise automatically to provide adequate clearance. The fowler up LED on the siderail flashes to indicate activation.</p>		
G	Skooch In (optional)	Moves the fowler toward the foot end of the product
H	Skooch Out (optional)	Moves the fowler toward the head end of the product
<p>Note: The fowler will rise to 35° before the seat depth decreases when you press Decrease Seat Depth. After the fowler is elevated to 35° or greater, the Skoocher LED illuminates to indicate that Skoocher is active. You can then change the seat depth.</p>		
I	Skoocher LED (optional)	Illuminates to indicate that Skoocher is active, and that you can change seat depth
J	Trend/Pelvic Tilt	Lowers the head end of the product into a Trend-Like position. Can also be used for Pelvic Tilt.
<p>Note: When the product is at a low height and you press Trend/Pelvic Tilt, the litter may raise automatically to provide adequate clearance. The fowler up LED on the siderail will flash to indicate activation.</p>		

Operation


Nurse control panel (Outside siderail) (Continued)

	Button name	Function
K	Nurse call (optional)	Activates nurse call
<p>Note: The nurse call light on the siderails show the patient which button to push to contact the nursing station in a darkened room.</p>		
L	Level	Raises the head end of the product or returns the product to level
M	Inside siderail light	<p>Illuminates the inside siderail buttons</p> <p>Five settings are available for the inside siderail lights:</p> <ul style="list-style-type: none"> • Press the button once to turn on the lights at low intensity • Press the button twice for medium intensity • Press the button three times for high intensity • Press the button four times to only leave the nurse call light on • Press the button five times to turn all siderail panel lights off
N	Siderail bed motion control lockout	<p>Locks all siderail bed motion controls.</p> <hr/> <p>⚠ WARNING</p> <p>Always lock the nurse control panel and patient control panel when the patient is unattended.</p> <hr/> <p>To lock the patient control panel, press the lockout button once. The lockout button illuminates when activated.</p> <p>To unlock the patient control panel, press the lockout button a second time.</p>
O	Power LED	<p>Illuminates when you plug the product into a hospital-grade receptacle.</p> <p>Blinks when you need to replace the 9V nurse call battery.</p>

Operation

English
EN


Patient control panel (Inside siderail)

	Button name	Function
A	Foot section up	Raises the foot section
B	Foot section down	Lowers the foot section
C	Fowler up	Raises the fowler
D	Fowler down	Lowers the fowler
E	Nurse call LED (optional)	Illuminates amber when the patient presses the Nurse Call button
F	Nurse call (optional)	Activates nurse call
G	Nurse call answer LED (optional)	Illuminates green when a nurse answers a call
H	TV/Radio power (optional)	Turns on the TV or the radio
I	TV/Radio volume up (optional)	Increases the volume
J	TV/Radio volume down (optional)	Decreases the volume
K	Room light (optional)	Turns the room light on or off
L	Product overhead light (optional)	Turns the product overhead light on or off

Operation

Patient control panel with optional Smart TV (Inside siderail)


	Button name	Function
A	Foot section up	Raises the foot section
B	Foot section down	Lowers the foot section
C	Fowler up	Raises the fowler
D	Fowler down	Lowers the fowler
E	Nurse call LED (optional)	Illuminates amber when the patient presses the Nurse Call button
F	Nurse call (optional)	Activate nurse call
G	Nurse call answer LED (optional)	Illuminates green when a nurse answers a call
H	Smart TV power (optional)	Turns on Smart TV
I	Smart TV volume up (optional)	Increases Smart TV volume
J	Smart TV volume down (optional)	Decreases Smart TV volume
K	Channel up (optional)	Changes the TV channel up
L	Channel down (optional)	Changes the TV channel down

Operation

English
EN

Patient lumbar control functions with optional Smart TV and optional Air Mattress (Inside siderail)


	Button name	Function
A	Product overhead light	Turns the product overhead light on or off
B	Room light	Turns the room light on or off
C	Firm seat	Increases support to the seat section
D	Soften seat	Decreases support to the seat section
E	Firm back	Increases support to the patient's lower back
F	Soften back	Decreases support to the patient's lower back

Operation

Motion pendant functions (optional) - 5000-081-021

CAUTION

- Always place the pendant safely on the mattress while the pendant is in use.
- Always unplug the pendant from the product and store it when the pendant is not in use.
- Always unplug the pendant from the product and store it when using oxygen administering equipment to avoid a possible fire hazard.
- Do not squeeze or pinch cables in the bed frame when routing cables from other equipment.

	Button name	Function
A	Trend/Pelvic Tilt	Lowers the head end of the product into a Trend-Like position. Can also be used for Pelvic Tilt.
B	Reverse Trend	Raises the head end of the product or returns the product to level
C	Fowler up	Raises the fowler
D	Fowler down	Lowers the fowler
E	Foot section up	Raises the foot section
F	Foot section down	Lowers the foot section
G	Litter up	Raises the litter
H	Litter down	Lowers the litter


Operation

English
EN

Motion pendant with nurse call functions (optional) - 5000-081-023

CAUTION

- Always place the pendant safely on the mattress while the pendant is in use.
- Always unplug the pendant from the product and store it when the pendant is not in use.
- Always unplug the pendant from the product and store it when using oxygen administering equipment to avoid a possible fire hazard.
- Do not squeeze or pinch cables in the bed frame when routing cables from other equipment.

	Button name	Function
A	Nurse call	Activates nurse call
B	Trend/Pelvic Tilt	Lowers the head end of the product into a Trend-Like position. Can also be used for Pelvic Tilt.
C	Reverse Trend	Raises the head end of the product or returns the product to level
D	Fowler up	Raises the fowler
E	Fowler down	Lowers the fowler
F	Foot section up	Raises the foot section
G	Foot section down	Lowers the foot section

Operation


Motion pendant with nurse call functions (optional) - 5000-081-023 (Continued)

	Button name	Function
H	Litter up	Raises the litter
I	Litter down	Lowers the litter

Motion and communication pendant functions (optional) - 5000-081-020

⚠ CAUTION

- Always place the pendant safely on the mattress while the pendant is in use.
- Always unplug the pendant from the product and store it when the pendant is not in use.
- Always unplug the pendant from the product and store it when using oxygen administering equipment to avoid a possible fire hazard.
- Do not squeeze or pinch cables in the bed frame when routing cables from other equipment.

	Button name	Function
A	TV/Radio	Turns on the TV or the radio and changes the channels
B	TV/Radio volume up	Increases the volume
C	TV/Radio volume down	Decreases the volume
D	Nurse call	Activates nurse call
E	Room light	Turns the room light on or off

Operation

English
EN

Motion and communication pendant functions (optional) - 5000-081-020 (Continued)

	Button name	Function
F	Overhead light	Turns the product overhead light on or off
G	Fowler up	Raises the fowler
H	Fowler down	Lowers the fowler
I	Foot section up	Raises the foot section
J	Foot section down	Lowers the foot section

Preventive maintenance

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage.

Remove product from service before performing preventive maintenance. Preventive maintenance should only be performed by trained or certified personnel.

Note: Clean and disinfect the exterior of the mattress before inspection, if applicable.

Inspect the following items:

- _____ All welds and all fasteners are secure
- _____ Tubing or sheet metal for bends or breaks
- _____ Casters are free of debris
- _____ Casters are secure and swivel properly
- _____ Casters lock securely by depressing the brake pedal
- _____ Steer caster latches properly
- _____ Siderails move and latch properly
- _____ Fowler operates properly
- _____ Fowler (back rest) slide operates properly
- _____ Litter Up/Down operates properly
- _____ Foot section operates properly
- _____ Footrests operate properly
- _____ Trend operates properly
- _____ IV pole is intact and operating properly
- _____ Mattress cover after each use
- _____ Mattress cover for rips or cracks
- _____ Headboard and siderail panels for cracks or splits
- _____ All covers are not damaged and are not producing sharp edges
- _____ Night light operating properly
- _____ CPR release operating properly
- _____ All siderail motion functionality
- _____ Nurse call functionality (Optional equipment)
- _____ 9V Nurse call battery (replace annually) (Optional equipment)
- _____ Lubricate where required
- _____ Pendant for any physical damage
- _____ Power cord not frayed
- _____ Cables not worn or pinched
- _____ All electrical connections tight
- _____ All grounds secure to the frame
- _____ Ground impedance not more than 100 mΩ (milliohms)
- _____ Current leakage not more than 300 μA (microamps)
- _____ Ground chains are clean, intact, and have at least two links touching the floor
- _____ Enclosure is free from wear, tear, stresses and mechanical damage
- _____ No rust or corrosion of parts
- _____ Labels for legibility, proper adherence, and integrity

Product Serial Number:		

Completed by: _____ Date: _____

Cleaning and disinfecting with wipes

English
EN

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred wipes (2060-000-001 6" x 10" or 2060-000-002 9" x 12") include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether - less than 3%

Note: For safety information, read the product label.

To clean or disinfect the external product surface:

1. To clean, wipe external surfaces with a fresh, clean wipe to remove all visible soils. Repeat as necessary until the product is clean.

Notes

- Use as many wipes as necessary.
 - Complete step 1 before you disinfect.
2. To disinfect, wipe external surfaces with a fresh, clean wipe until wet. Allow the external surface to remain wet for two minutes at room temperature.
 3. Allow the product to dry before you return it to service.

Cleaning

WARNING

- Do not clean, service, or perform maintenance while the product is in use.
 - When large spills occur in the area of the circuit boards, 120VAC cables, and motors, always immediately unplug the power cord from the wall outlet. Remove the patient from the product, clean up the fluid, and have service personnel completely inspect the product. Fluids can cause unpredictable operation and decreased functionality of any electrical product. Do not return the product to service until it is completely dry and has been thoroughly tested for safe operation.
-

CAUTION

- Always unplug the product before cleaning or servicing.
 - Always unplug the product, set the brakes, and place blocks under the litter frame for support when working under the Birthing Bed.
-

The recommended cleaners for this product's surfaces include the following:

- Quaternary Cleaners (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)
- ≤ 21% propanol alcohol

Hand wash all surfaces of the product with warm water and mild detergent. Dry thoroughly.

Avoid oversaturation and make sure the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

CAUTION

- Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature corrosion of critical components. Failure to follow these cleaning instructions may void your warranty.
 - Do not steam clean, pressure wash, ultrasonically clean, or immerse any part of the product in water. Exposure to water may damage the internal electric parts. These methods of cleaning are not recommended and may void this product's warranty.
-

Cleaning the mattress

For mattress cleaning and disinfecting, use warm water and a neutral detergent. You can also use a sodium hypochlorite solution for cleaning. If you use a chlorinated bleach solution (5.25%), dilute the solution to less than 2 parts bleach to 100 parts water.

WARNING

Always completely dry the mattress after cleaning. Failure to thoroughly rinse and dry mattresses after cleaning may cause damage to the mattress and may void this product's warranty.

CAUTION

Inspect the mattress cover after each use. Openings in the mattress cover may allow fluids to enter the mattress. Discontinue using the mattress and safely dispose the mattress if any rips or cracks are found in the mattress cover.

Cleaning

English
EN

Cleaning the mattress (Continued)

Virex® TB is not an approved cleaning agent for this product. Do not use for product cleaning.

Removing and cleaning the mattress:

To remove the head end mattress, pull the mattress back and wipe clean.

To remove the foot section mattress, unsnap along the edges, pull up on the mattress, and wipe clean. Additional snaps are located along the perineal edge of the foot section for additional cleaning access.

To clean the perineal edge, unsnap the mattress, pull up and wipe clean.

Figure 1-56: Cleaning behind the head end mattress

Figure 1-57: Cleaning the Perineal Edge

EMC Information

Guidance and manufacturer's declaration - electromagnetic emissions		
The Model 4701 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 4701 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The Model 4701 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Model 4701 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class A	
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Model 4701 is suitable for use in the electromagnetic environment specified below. The customer or the user of the Model 4701 should assure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) EN/IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast Transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 8 kV differential mode ± 2 kV common mode	± 8 kV differential mode ± 2 kV common mode	Main power quality is that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	Main power quality should be that of a typical commercial or hospital environment. If the user of the Model 4701 requires continued operation during power main interruptions, it is recommended that the device be powered from an

EMC Information

English
EN

(Continued)

Guidance and manufacturer's declaration - electromagnetic immunity			
			uninterrupted power supply or a battery.
Power frequency (50/60Hz) magnetic field EN/IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage before applications of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The Model 4701 is suited for use in the electromagnetic environment specified below. The customer or the user of the Model 4701 should assure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Model 4701, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended separation distance</p> $D=(1.2) (\sqrt{P})$ $D=(1.2) (\sqrt{P})$ <p>80 MHz to 800 MHz</p> $D=(2.3) (\sqrt{P})$ <p>800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

EMC Information

(Continued)

Guidance and manufacturer's declaration - electromagnetic immunity			
			
<p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 4701 is used exceeds the applicable RF compliance level above, the Model 4701 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Model 4701.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.</p>			

Recommended separations distances between portable and mobile RF communications equipment and the Model 4701			
<p>The Model 4701 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model 4701 can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 4701 as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $D=(1.2) (\sqrt{P})$	80 MHz to 800 MHz $D=(1.2) (\sqrt{P})$	800 MHz to 2.5 GHz $D=(2.3) (\sqrt{P})$
0.01	1.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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