

ELECTRIC ACUTE CARE BED FL20E Model



GOBED+™

TECHNICAL ASSISTANCE AND PARTS

1 888 233-6888 (Canada)

1 800 327-0770 (United States)

Outside Canada and the United States: Contact your local representative

72-0439 REV A
F15-43-E-A

January 2006
Printed in Canada

TABLE OF CONTENTS

1. INTRODUCTION.....	4
1.1 Bed Specifications	4
1.2 Technical Support.....	5
1.3 Warning / Caution / Note Definitions	5
1.4 Safety Tips and Guidelines	5
1.5 Warranty.....	7
Limited Warranty	7
To Obtain Service and/or Parts.....	7
Return Authorization	7
Damaged Merchandise.....	7
1.6 Symbols.....	7
1.7 Set-Up Procedure	8
1.8 Cleaning and Preventative Maintenance.....	9
Bed Cleaning and Mattress Care	9
Preventative Maintenance	10
1.9 Bed Position Pictograms.....	11
1.10 Bed Illustration	11
2. OPERATION GUIDE	12
2.1 Switching On Power	12
2.2 Applying the Brakes.....	12
2.3 Moving the bed	12
2.4 Foley Bag Hook Usage.....	13
2.5 Foot Support Arm Operation.....	13
2.6 Patient Restraint Strap Locations.....	13
2.7 CPR Emergency Release (Optional).....	13
Lowering the Fowler During a Power Failure	14
2.8 Night Light Usage (Optional).....	14
2.9 Nurse Call Usage (Optional)	14
2.10 Auxiliary Power Outlet Usage (Option Available w/120V~ Bed Model).....	14
2.11 Positioning Siderails	15
2.12 Head and Foot Board Operation	15
2.13 Head Siderail Function Guide	16
2.14 Foot Board Control Panel Guide	17
2.15 Function Lockout Usage	17
2.16 Weigh System Control Panel Guide (Optional)	18
Operating (Zeroing) the Scale Before Putting a New Patient on the Bed	18
Adding or Removing Equipment when a Patient is on the Bed	19
Changing Manually the Numerical Value of Displayed Weight.....	19
Changing the Weight Measure Unit	19
2.17 Chaperone Bed Exit (Optional)	20
Zeroing the Bed Exit System	20
Adding or Removing Equipment	20
2.18 Chaperone Bed Exit with Zone Control (Optional).....	21
Zone Settings	21
Zeroing the Bed Exit System	21
Adding or Removing Equipment	22
3. ACCESSORIES	22

1. INTRODUCTION

This manual is designed to assist in the operation of Stryker's GOBED⁺ Acute Care bed. Read it thoroughly before operating the bed. Hospital staff should be able to refer to this manual at all time when using the bed.

This Operations Manual is an integral part of the bed and should be included if the bed is sold or transferred.

1.1 BED SPECIFICATIONS *

Safe Working Load	579 lb (263 kg)
Weigh System Capacity	Patients weighing up to 500 lb (227 kg)
Weigh System Accuracy	± 2 % of patient weight from 100 lb (45.3 kg) to 500 lb (227 kg) ± 2 lb of patient weight under 100 lb (45.3 kg)
Overall Bed Length/Width	92 7/8 x 42 9/16" (235 x 108.1 cm) - siderails up 92 7/8 x 38 3/4" (235 x 98.4 cm) - siderails down
Weight w/Boards	417 lb (189 kg)
Patient Sleep Surface	35 x 80" (89 x 203 cm) extendable to 82" (208 cm) & 84" (213 cm)
Minimum/Maximum Bed Height	14.5 to 28.9" (36.8 to 73.4 cm)
Mattress Size Recommended	35 x 80" (89 x 203 cm); 35 x 82" (89 x 208 cm); 35 x 84" (89 x 213 cm)
Maximum Thickness	6.25" (15.9 cm)
Fowler Angle w/CPR Mechanism	0 to 60°
Fowler Angle w/o CPR Mechanism	0 to 65°
Knee Gatch Angle w/o Auto Contour	0 to 32°
Knee Gatch Angle w/Auto Contour	0 to 24°
Trendelenburg/Reverse Trendelenburg	+14 to -14°
Environmental Conditions	
- Transport and storage	
- Ambient Temperature	-40 to 70°C (-40 to 158°F)
- Relative humidity	10 to 100%
- Atmospheric Pressure	500 to 1060 hPa
- Operating **	
- Ambient Temperature	18.3 to 26.7°C (65 to 80°F)
- Relative humidity	20 to 80% without condensation
- Atmospheric Pressure	700 to 1060 hPa
*** Electrical Requirements - all electrical requirements meet CSA C22.2 No. 601.1-M90, UL 2601 and IEC 60601-1, 60601-2-38 specifications.	100V~, 50-60Hz, 7.5A -Two 250V, 10A Fast Acting Fuses 120V~, 50-60Hz, 4.0A (9.8A w/Auxiliary Outlet) -Two 250V, 10A Fast Acting Fuses 200V~, 50-60Hz, 3.2A -Two 250V, 6.3A Slow Blow Fuses 220V~, 50-60Hz, 2.9A -Two 250V, 6.3A Slow Blow Fuses 240V~, 50-60Hz, 2.7A -Two 250V, 6.3A Slow Blow Fuses

* Stryker pays special attention to product improvement and reserves the right to change specifications without notice.

**Operating environment recommended to ensure the scale precision.

*** The device has a 10% duty cycle.

1.2 TECHNICAL SUPPORT

For questions regarding this product, contact the following Technical Service department or your local representative:

In Canada:

Stryker Canada
1 888 233-6888
45, Innovation Drive
Hamilton, Ontario, L9H 7L8
Canada

In the United States:

Stryker Medical
1 800 327-0770
3800, East Centre Avenue
Portage, MI 49002
USA

1.3 WARNING / CAUTION / NOTE DEFINITIONS

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

**WARNING**

The personal safety of the patient or user may be involved. Disregarding this information could result in injury to the patient or user.

**CAUTION**

These instructions point out special procedures or precautions that must be followed to avoid damaging the equipment.

NOTE

Notes provide special information to make maintenance easier or important instructions clearer.

1.4 SAFETY TIPS AND GUIDELINES

Before operating the **GOBED⁺**, it is important to read and understand all information in this manual. Carefully read and strictly follow the safety guidelines listed below. It is important that all users have been trained and educated on the inherent hazards associated with the usage of electric beds.

**WARNING**

- The **GOBED⁺** is not intended for pediatric use.
- This bed is equipped with a hospital grade plug for protection against shock hazard. It must be plugged directly into a properly grounded receptacle. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Shock Hazard - Improper handling of the power cord may result in damage to the power cord and potential shock hazards. If damage has occurred to the power cord, immediately remove the bed from service, and contact the appropriate maintenance personnel. Failure to do so could result in death or serious injury.
- Serious injury can result if caution is not used when operating the bed. Operate the bed only when all people and equipment are clear of the electrical and mechanical systems.
- Always apply the brakes when a patient is on the bed (except during transport). Serious injury could result if the bed moves while a patient is getting on or off the bed. After the brake pedal is engaged, push on the bed to ensure the brakes are securely locked.
- To help reduce the number and severity of falls by patients, always leave the bed in the lowest position when the patient is unattended.
- Leave the siderails fully up and locked when the patient is unattended. When raising the siderails, be sure that you hear the "click" that signals the up and locked condition. Pull firmly on the siderail to ensure it is locked into position.
- Keep siderails in the fully raised position and the sleep surface horizontal in its lowest position when the patient is sleeping unless the patient's medical condition dictates otherwise.
- When the sleep surface sections are articulated, ensure that all patient's extremities are within the raised siderails to avoid patient injury.

- To reduce risk of injury, ensure the sleep surface is horizontal and in the lowest position with the siderails fully raised and locked when moving the bed with a patient in it.
- To avoid injury to the patient and/or user, do not attempt to move the bed laterally with the steer pedal engaged. The fifth steer wheel cannot swivel.
- Siderails, with or without their padded covers or nets, are not intended to serve as restraint devices to keep patient from exiting the bed. Siderails are designed to keep a patient from inadvertently rolling off the bed. It is the responsibility of the attending medical personnel to determine the degree of restraint necessary to ensure a patient will remain safely in bed. Failure to utilize the siderails properly could result in serious patient injury.
- When a patient's condition requires greater safety measures for his/her security, use the lockout switches in the foot board control panel to inhibit the siderail functions or remove any optional pendant control and install protective pads on the siderails.
- The CPR emergency release is for emergency use only. When activating the CPR release handle, all people and equipment must be removed from the area below and around the head, thigh and foot sections of the bed or serious personal injury or equipment damage could occur.
- Possible fire hazard exists when this bed is used with oxygen administering equipment other than nasal, mask type or half bed-length tent type. It is recommended to unplug the bed power cord from the wall when oxygen-administering equipment is used. When using a half bed-length tent type, ensure the siderails are outside the oxygen tent and oxygen tent should not extend below the mattress support level.
- The Bed Exit system (optional) is intended only to aid in the detection of a patient exiting the bed. It is not intended to replace patient monitoring protocol. The Bed Exit system signals when a patient is about to exit the bed. The addition or removal of equipment with a Bed Exit system active must be done using the "Adding or Removing Equipment with the System Active" procedure, otherwise the sensitivity of the system may be affected and the readings of the patient's movements in the bed be erroneous.
- The Bed Exit system is not designed to be used with patients weighing less than 50 lb (23 kg).
- When large fluid spills occur in the area of the circuit board, cables and motors, immediately unplug the bed power cord from the wall outlet. Remove the patient from the bed and clean up the fluid. Have maintenance completely check the bed. Fluids can have an adverse effect on operational capabilities of any electrical product. **DO NOT** put the bed back into service until it is completely dried and has been thoroughly tested for safe operation.
- Do not steam clean, hose off or ultrasonically clean the bed. Do not immerse any part of the bed. The internal electrical parts may be damaged by exposure to water. Hand wash regularly all surfaces of the bed with warm water and a mild detergent. Wipe cleaned surfaces dry to avoid build up of cleaning substance. Inspect the mattress after each use. Discontinue use if any cracks or rips are found in the mattress cover, which may allow fluid to enter the mattress. Failure to properly clean the mattress, or dispose of it if defective, may increase the risk of exposure to pathogenic substances and may bring about diseases to the patient and user.
- Preventive maintenance should be performed at least once a year to ensure all bed features are functioning properly. Ensure that any bed malfunction is promptly reported to your service personnel for immediate attention.
- Always unplug the bed power cord from the wall outlet when servicing or cleaning the bed. When working under the bed with the bed in the high position, always place blocks under the mattress support frame and apply the brakes to prevent injury in case the Bed Down switch is accidentally pressed.
- To avoid damage to the siderail mechanisms, do not move the bed using the raised siderails. Use the push/pull handles integrated to the boards to move the bed.
- Before using the optional emergency crank during a power failure, always unplug the power cord. An unexpected return of power could rotate the handle and cause injury.
- The mattress thickness should not exceed 6.25" (15.9 cm)
- Because individual beds may have different options, foot boards should not be moved from one bed to another. Mixing foot boards could result in unpredictable bed operation.

NOTE

Throughout this operations manual, the words "right" and "left" refer to the right and left sides of a patient lying face up on the bed.

1.5 WARRANTY

LIMITED WARRANTY

All Stryker products are guaranteed against material or manufacturing defects, improper operation of mechanisms, and premature wear of bed components under normal use conditions.

For questions regarding warranty, please contact the Technical Service department (see section 1.2) or your local representative.

TO OBTAIN SERVICE AND/OR PARTS

For an on-site diagnosis of a malfunction a Field Service Representatives or to order replacement parts, simply contact our Service department (see section 1.2) or your local representative. For the part ordering procedure, refer to section 1.5 of the bed maintenance manual, "To Obtain Service and/or Parts".

RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Technical Service department. An authorization number will be provided, which must be clearly printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items.

DAMAGED MERCHANDISE

Claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery notice at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Claims for any short shipment must be made within 5 days of invoice.

1.6 SYMBOLS



Warning, consult accompanying documents



Fuse rating for bed with the 100V~ and 120V~ electric systems



Fuse rating for beds with 200V~, 220V~ and 240V~ electric systems



Protective earth (ground)



Alternating Current



Type B Equipment

IPX4 Protection from liquid splash

1.7 SET-UP PROCEDURE

CHECKLIST

It is important to ensure that the bed is working properly before it is put into service. The following list will help ensure that each part of the bed is checked.

_____ Install the foot and head boards on the bed. Insert the foot board carefully so that the board and the casing connectors fit in smoothly.



WARNING

Because individual beds may have different options, foot boards should not be moved from one bed to another. Mixing foot boards could result in unpredictable bed operation.

NOTE

Boards can be fixed permanently by using two 1/4" flat washers/1/4" spring washers/1/4-20 x 1 3/4" hexagon screws for the head board and two 1/4" flat washers/1/4" spring washers/1/4-20 x 2 1/4" hexagon screws for the foot board. Simply insert a screw and the two different washers in factory installed nuts (refer to drawing 20-0130L, attached to the parts list L20-001, for the illustration of the nut (part VE20A1N -ref. no. 17)) accessible through holes located on the inner sides of the foot end and head end casings.

_____ Plug the power cord to the bed connector at the head end of the bed and into a properly grounded hospital grade wall outlet. Turn on the power switch located on the power connector.



WARNING

The bed is equipped with a hospital grade plug for protection against shock hazard. It must be plugged directly into a properly grounded receptacle. Grounding reliability can be achieved only when a hospital grade receptacle is used.

- _____ Depress the red "TOTAL BRAKE" side of the pedal at both sides of the bed fully down to apply the brakes (see page 12). Push on the bed to ensure it is immobilized. Toggle the pedal to the neutral position and ensure the brakes are released.
- _____ Depress the green "AXIAL STEER" side of the pedal at both sides of the bed fully down to put the bed in the steer mode (see page 12). Ensure the 5th steer wheel is properly engaged. Toggle the pedal to the neutral position and ensure the 5th wheel disengages.
- _____ Ensure the siderails raise, lock in the up position and lower smoothly (see page 15).
- _____ Run through each control of the foot board control panel (see page 17).
- _____ Verify the optional Weigh and Bed Exit systems for proper operation (see page 18, "Operating (Zeroing) the Scale..." and page 20 and 21, "To **activate** the Bed Exit...").
- _____ Run through each function on both head end siderails (see page 16).
- _____ Raise the bed fully up and activate the Trendelenburg function (see page 17). Ensure the head end lowers to the full down position.
- _____ Raise the bed fully up and activate the reverse Trendelenburg function (see page 17). Ensure the foot end lowers to the full down position.
- _____ Raise the Fowler fully up. Using the CPR hold-to-run handle (see page 13), lower the Fowler gradually to flat position by pulling, holding and releasing the handle several times. Ensure the Knee Gatch (if raised) also starts flattening when the Fowler is completely down. Verify that the Fowler motor resets itself when the Fowler is completely lowered. Wait 30 seconds - the time for the fowler motor to reset itself - and raise the Fowler to ensure the Fowler motor has indeed reset.
- _____ Verify the following optional equipment for proper operation: 120V auxiliary outlet, night light, two or three function pendant control, emergency crank operation, etc.

If any problems are found during bed set-up, contact our Technical Service department (see section 1.2).

1.8 CLEANING AND PREVENTATIVE MAINTENANCE

BED CLEANING AND MATTRESS CARE

**WARNING**

Always unplug the bed power cord from the wall outlet when cleaning or servicing the bed.

**CAUTION**

Do not use harsh cleaners, solvents or detergents. Do not steam clean, hose off or ultrasonically clean the bed. Do not immerse any part of the bed. The bed electrical parts may be damaged by exposure to water.

Germicidal disinfectant, used as directed, and/or Chlorine Bleach products are not considered mild detergents. These products are corrosive in nature and may cause damage to your bed if used improperly. If these types of products are used, ensure the beds are rinsed with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the beds will leave a corrosive residue on the surface of the bed, possibly causing premature corrosion of critical components. Failure to follow the above directions when using these types of cleaners may void this product warranty.

Bed Cleaning

Hand wash all surfaces of the bed with a soft cloth moistened with a solution of lukewarm water and a mild detergent.

Wipe the bed clean and dry thoroughly to avoid build up of cleaning solution.

NOTE

Do not use cleaning solutions containing a degreaser near the siderail pivots (see figure 2.2B, page 11 of the Maintenance Manual) to avoid deteriorating the properties of the grease used to ensure a smooth movement of the siderails.

Mattress Care

**WARNING**

Inspect the mattress after each use. Discontinue use if any cracks or rips are found in the mattress cover, which may allow fluid to enter the mattress. Failure to properly clean the mattress, or dispose of it if defective, may increase the risk of exposure to pathogenic substances and may bring about diseases to the patient and user.

• Inspection

Implement local policies to address regular care, maintenance, and cleaning of mattresses and covers. The cover cleaning procedure can be found below and on the mattress label.


Inspect mattress cover surface (also zip fasteners and cover inner surface if mattresses have zip fasteners) regularly for signs of damage. If the mattress cover is heavily stained or soiled, or torn, remove the mattress from service.

• Cleaning

Stains: Wash with lukewarm water using a mild detergent. Rinse with water and let dry. For tough stains use bleach diluted with ten parts of water.

PREVENTATIVE MAINTENANCE

Annual Checklist

- ___ All fasteners secure.
- ___ Inspect for excessive wear the oil-impregnated bronze shoulder spacers found at the bed hinge points. **Do not** lubricate these spacers; replace if found worn.
- ___ Engage the "TOTAL BRAKE" pedal at both sides of the bed and ensure the brakes are applied. Toggle the pedal to neutral and ensure the brakes are released.
- ___ Engage the "AXIAL STEER" pedal at both side of the bed and ensure the 5th steer wheel engages properly. Toggle the pedal to neutral and ensure the 5th wheel disengages.
- ___ Siderails move, latch and stow properly (see page 15).
- ___ All functions on the foot board working properly, including LED's (see page 17).
- ___ Calibrate the optional scale system (see section 4.6 "Scale System Calibration" of the bed Maintenance manual).
- ___ All siderail functions working properly (see page 16). Verify the nurse call function. Does the alarm reach the nurse station?
- ___ Fowler and Knee Gatch (if raised) flatten and the Fowler control motor resets itself automatically when hold-to-run CPR release handles (optional) are pulled until Fowler is flattened.
 - ___ Wait about 30 seconds - the time for the Fowler motor to reset itself - and raise the Fowler to ensure that the resetting of the motor has indeed occurred.
- ___ Verify the Fowler, Knee Gatch and Hi-lo movements to ensure that the motion interrupt switch integrated to each of the four actuators is operating properly.
- ___ Auxiliary outlet (option available only with 120V electric system beds) working properly (see page 14).
- ___ Optional night light working properly.
- ___ Foot prop rod working properly with the Knee Gatch or Auto Contour positioning activated.
- ___ No cracks or splits in the boards and siderails.
- ___ Head end bumpers tightly secured to frame and working properly.
- ___ No rips or cracks in mattress cover.
- ___ Power cord not frayed.
- ___ No cables worn or pinched.
- ___ All electrical connections tight.
- ___ All grounds secure to the frame.
- ___ All casters roll properly. Check caster for cuts, wear, etc.
- ___ Ground chain intact and in place.
- ___ Inspect and lubricate, if needed, the bed lubrication points described in section 2.2 of the Maintenance Manual.
-  ___ Measure current leakage and grounding continuity of the bed and the optional auxiliary outlet (check with our Technical Service department for the acceptable values).

NOTE

Preventive maintenance may need to be performed more frequently based on the usage level of the bed.

1.9 BED POSITION PICTOGRAMS

The following pictograms illustrate the GOBED⁺ basic positions.

Not illustrated here is the Cardiac Chair position that is obtained by combining the reverse Trendelenburg and the Auto Contour positions.




Figure 1.9

1.10 BED ILLUSTRATION




Figure 1.10

2. OPERATION GUIDE

2.1 SWITCHING ON POWER

The bed is equipped with a main power switch located at the head end of the bed where the power cord connects to the bed. Turn it on to activate the bed functions. When the bed power switch is turned off or in the event of a power failure, the settings of the lockout controls and the optional Weigh and/or Bed Exit calibration data are preserved.



WARNING

Shock Hazard - Improper handling of the power cord may result in damage to the power cord and potential shock hazards. If damage has occurred to the power cord, immediately remove the bed from service, and contact the appropriate maintenance personnel. Failure to do so could result in death or serious injury.

2.2 APPLYING THE BRAKES

The bed is equipped with a central locking system activated by a pedal (see page 11, “Bed Illustration”) located at the midpoint of the bed on both sides.



WARNING

Always apply the brakes when a patient is on the bed (except during transport) or when entering/exiting the bed. Serious injury could result if the bed moves while a patient is getting on or off the bed. After the brake pedal is engaged, push on the bed to ensure the brakes are securely locked.

Brake Pedal Operation

To **engage** the wheel brakes, fully depress the red “TOTAL BRAKE” side of the brake/steer pedal at either side of the bed.

To **disengage** the wheel brakes, toggle the pedal to the neutral position.

2.3 MOVING THE BED

The bed is equipped with a fifth steer wheel activated by a pedal (see page 11, “Bed Illustration”) located at midpoint of the bed on both sides.

The fifth steer wheel helps guiding the bed along a straight line and helps the bed pivot around corners.



WARNING

To reduce risk of injury, ensure the sleep surface is horizontal and in the lowest position with the siderails fully raised and locked when moving the bed with a patient in it.

To avoid injury to the patient and/or user, do not attempt to move the bed laterally with the steer pedal engaged. The fifth steer wheel cannot swivel.

To avoid damage to the siderail mechanisms, do not move the bed using the raised siderails. Use the push/pull handles integrated to the boards to move the bed.

Steer Pedal Operation

To **engage** the steer mode, fully depress the green “AXIAL STEER” side of the brake/steer pedal at either side of the bed.

To **remove** the steer mode, toggle the pedal to the neutral position.

2.4 FOLEY BAG HOOK USAGE

Foley bag hooks (see page 11, “Bed Illustration”) are found at three locations on both sides of the bed under the edges of the mattress support head, seat and foot sections.

NOTE

The Foley bag hooks move when the Fowler is raised or lowered. Fowler motion must be locked out when using these hooks to avoid inadvertent movement of the hooks.

2.5 FOOT PROP ROD OPERATION

A foot prop rod, integrated to the foot section, engages automatically as the Knee Gatch is raised thus maintaining the foot section in a nearly horizontal position. The prop rod enables the positioning of the sleep surface into the vascular position (see page 11, “Bed Positions”).

Since the prop rod is automatically engaged, it must be disengaged manually whenever a simple Knee Gatch raise is required. To do so, perform the following steps:

1. Before raising the Knee Gatch, reach under the foot section, lift the prop rod and hold it while raising the Knee Gatch through the foot board control panel.
2. As the Knee Gatch raises, the rod will come to pass the support brackets. Release it then and continue to raise the Knee Gatch. The foot section will then simply follow the move of the Knee Gatch without being propped.

2.6 PATIENT RESTRAINT STRAP LOCATIONS

The bed has 12 locations on the mattress support for installing patient restraint straps. Ten of them are located on the mattress support edges directly across from each other and the remaining two are located on the top edge of the head section (see page 11, “Bed Illustration”).



WARNING

Improperly adjusted restraint straps can cause serious injury to a patient. It is the responsibility of the attending medical personnel to determine proper use of restraint straps and restraint strap locations.

2.7 CPR EMERGENCY RELEASE (OPTIONAL)



WARNING

The CPR emergency release is for emergency use only. When activating the CPR release handle, all people and equipment must be removed from the area below and around the head, thigh and foot sections of the bed or serious personal injury or equipment damage could occur.

When quick access to the patient is needed, and the Fowler is raised, pull outward one of the two CPR release handles located under the upper right and left sides of the Fowler (see page 11, “Bed Illustration”). The Fowler and the Knee Gatch (if raised) will automatically flatten.

The CPR release handle can be released at any time to stop the lowering of the Fowler. But doing so will subsequently require that the Fowler be completely lowered, using the CPR release handle or the Fowler down control, to enable the Fowler motor to reset itself. **Failing to do so will prevent the Fowler from being fully raised.**

NOTE

The use of the CPR emergency release creates a situation where the course of the Fowler motor is temporarily out of sync with the actual position of the Fowler. The situation is corrected, and only then, when the Fowler is completely lowered; then the Fowler motor begins an automatic resetting process to harmonize its course with the actual (flatten) position of the Fowler. During the time the resetting process goes on, approximately 30 seconds, the Fowler controls are not available.

- **Lowering the Fowler During a Power Failure**

The CPR emergency release can also be used during a power failure to partially or completely lower the Fowler. Simply pull one of the CPR handles until the desired angle is reached. When current resume, fully lower the Fowler to enable the Fowler motor to reset itself (see above note).

If the bed is not equipped with the CPR emergency release, unplug the bed power cord from the wall outlet and use the emergency handle (optional) to lower the Fowler. A pictogram affixed on the head end casing cover illustrate clearly in which direction the handle should be turned to obtain the desired angle.

2.8 NIGHT LIGHT USAGE (OPTIONAL)

The bed may be equipped with an optional photoelectric night light to illuminate the floor area around the bed. The night light turns on as the room lights dim.

2.9 NURSE CALL USAGE (OPTIONAL)

The bed may be equipped with an optional nurse call function allowing the patient to communicate through speakerphones with the nurse station by simply pressing the nurse call red button integrated to both head end siderail inner control panels (see page 16).

The communication between the patient and the nurse station is established the moment a member of the nursing staff responds to the call signal.

Included with this option are two Ø1/4" phono jacks located on both sides of the head end casing. The phono jacks enable the use of a nurse call cord which can be placed within reach of a patient who is not in bed.

NOTE

The nurse call system will automatically send a call signal to the nurse station if the communication between the bed and the nurse station is interrupted following a power failure, the switching off of the bed or the disconnection of the nurse call communication cable.

2.10 AUXILIARY POWER OUTLET USAGE (OPTION AVAILABLE W/120V~ BED MODEL)

The North American version (120V~ electrical system) of the GOBED⁺ may be equipped with an optional 120V~ auxiliary power outlet integrating a 5A breaker. Located on the left side at the foot end of the bed (see page 11, "Bed Illustration"), this feature provides nursing staff with a convenient power source for peripheral equipments.

**WARNING**

Use only hospital grade electric equipment consuming 5A or less with the auxiliary power outlet (optional). The use of standard electric equipment may bring the current leakage to a level unacceptable for hospital equipment.

2.11 POSITIONING SIDERAILS

The siderails have one lock position, in the upper position. They can be tucked away against the mattress support when not in use.



WARNING

Leave the siderails fully up and locked when the patient is unattended. When raising the siderails, be sure that you hear the "click" that signals the up and locked condition. Pull firmly on the siderail to ensure it is locked into position.

Siderails, with or without their padded covers or nets, are not intended to serve as restraint devices to keep patient from exiting the bed. Siderails are designed to keep a patient from inadvertently rolling off the bed. It is the responsibility of the attending medical personnel to determine the degree of restraint necessary to ensure a patient will remain safely in bed. Failure to utilize the siderails properly could result in serious patient injury.

Keep siderails in the fully raised position and the sleep surface horizontal in its lowest position when the patient is sleeping unless the patient's medical condition dictates otherwise.

NOTE

The rail has to be removed from its tucked position by pulling it outward prior to raise it.

To **raise** a head siderail, grasp the rail, pull it outward and swing it upward towards the head end of the bed until it latches in the up position.

To **raise** a foot siderail, the same procedure is required as for the head siderail, however, the siderail swings to the foot end of the bed.

To **lower** a siderail, lift it up slightly, push in the release lever located under the mattress support (see page 11, "Bed Illustration") and hold siderail as it rotates down to low position. Push the siderail against the mattress support to tuck it away.

2.12 HEAD AND FOOT BOARD OPERATION

The head and foot boards slide down into two mounting sockets located at each end of the bed. Both boards can be removed or replaced easily. They may also be mounted permanently (see page 8, "Set-Up Procedures").



WARNING

Because individual beds may have different options, foot boards should not be moved from one bed to another. Mixing foot boards could result in unpredictable bed operation.

Removing Boards


Seize both ends of the board and lift it up until posts come out of the mounting sockets.

Replacing Boards

- **Head board:** Insert the board with the laminated finish side facing inside the bed.
- **Foot board:** Insert the board with the laminated finish side facing outside the bed. Be cautious when inserting the foot board to avoid damaging the board and foot end casing connectors.

2.13 HEAD SIDERAIL FUNCTION GUIDE

Outside Siderail (Patient's Right Rail)



Outside Siderail (Patient's Left Rail)





Figure 2.13A

- A:** Press to raise Fowler
- B:** Press to lower Fowler
- C:** Press to raise sleep surface

- D:** Press to lower sleep surface
- E:** Press to raise Knee Gatch
- F:** Press to lower Knee Gatch

Inside Siderail (Patient's Right Rail)



Inside Siderail (Patient's Left Rail)




Figure 2.13B

- A:** Press to raise Fowler
- B:** Press to lower Fowler
- E:** Press to activate the Nurse Call function (optional)

- C:** Press to raise Knee Gatch
- D:** Press to lower Knee Gatch




Figure 2.13C

F: Nurse call speakerphones (this panel is optional equipment and will be present on the inner side of both head siderails if the bed is equipped with the Nurse Call option, see page 14, "Nurse Call Usage").

2.14 FOOT BOARD CONTROL PANEL GUIDE




Figure 2.14

- A1:** Press to lock out Fowler controls at both head siderails. The Fowler padlock icon will come on. Press again to unlock.
- A2:** Press to lock out Knee Gatch controls at both head siderails. The Knee Gatch padlock icon will come on. Press again to unlock.
- A3:** Press to lock out the bed height adjustment (Hi-lo) at both head siderails (outer sides). The Hi-lo padlock icon will come on. Press again to unlock.

NOTE

The foot board controls are not affected by the lockout switches.

Lockout settings are automatically saved in the event of a power failure and restored following resumption of power.

B1: Press to raise Fowler

B2: Press to lower Fowler

C1: Press to raise Knee Gatch

C2: Press to lower Knee Gatch

D1: Press to raise the sleep surface

D2: Press to lower the sleep surface

E1: Press to obtain the Trendelenburg position (head down, foot up)

E2: Press to obtain the reverse Trendelenburg position (head up, foot down)

NOTE

To replace the sleep surface to horizontal position after a Trendelenburg positioning, simply use the bed height controls to either raise or lower the sleep surface to its limits. Then position the bed to the desired height.

F: Auto Contour positioning control. The Auto Contour positioning function is operational by default (padlock icon off). To deactivate it, press on the Contour switch (F) (padlock icon lights up). The Auto Contour and Knee Gatch lockout controls are interrelated, activating one activates the other.

2.15 FUNCTION LOCKOUT USAGE

Whenever a patient or visitor should be restricted from operating the siderail controls, activate the appropriate lockout switch (A1, A2, A3, figure 2.16) in the foot board control panel. The lockout switches are for the safety of the patient. Use them when appropriate.

2.16 WEIGH SYSTEM CONTROL PANEL GUIDE (OPTIONAL)




Figure 2.16

A: /ENTER key:

- Press the key once to activate the display. The Scale mode will be displayed.
- Press when instructed to press ENTER.

B/C: MENU Up/Down keys:

- Press either key repeatedly until desired mode is reached.

D: Scale Display.

Scale Menu Modes

Order of appearance of the modes when the MENU Down key is pressed repeatedly:

Mode	Display First Line
Scale	WEIGHT ANGLE
Change Equipment	CHANGE EQUIP.
Change Patient Weight	CHNG PTNT WEIGHT
Select Weight Unit (lb/kg)	UNITS
Scale Zero	SCALE ZERO

NOTE

When the bed is switched on, the scale display will show "STRYKER / SCALE V x.x" where x.x is the version of the control software.

The Weigh system does not operate when the bed is inclined more than 12° in the Trendelenburg or reverse Trendelenburg position. The LCD will display the current bed angle but not the patient's weight.

Operating (Zeroing) the Scale Before Putting a New Patient on the Bed

- Prepare the bed for a patient stay (linens, pillows, etc.) and press the key, the display will read:

WEIGHT ANGLE
XXX.X lb +/- XX.X°

- Press the MENU Up key to access the Scale Zero mode (as written on the membrane). The display will read:

SCALE ZERO / PRESS ENTER


- Press and hold the ENTER key. The display will read:
HOLD TO ZERO WT. (keep pressing the ENTER key), followed by:
RELEASE TO ZERO (release the ENTER key), followed by:
DO NOT TOUCH BED (be sure nobody touches the bed while this message appears).
- The scale system will return to the Scale mode and will display zero for the weight and the current angle value. The bed is now ready for the patient.

NOTE

The scale display turns off automatically after one minute of idle time, but the system itself remains active in the background. Pressing the key will reactivate the display in the Scale mode.


On beds equipped with the Weigh and Bed Exit systems, zeroing one system also zeroes the other. Do not zero the bed while a patient is on the bed. An inaccurate patient weight reading will result. If this should occur, use the Change Patient Weight mode to enter the correct patient weight or remove the patient from the bed and zero the bed again (see above "Operating (Zeroing) the Scale Before Putting a New Patient on the Bed").

Adding or Removing Equipment when a Patient is on the Bed

- Press the  key to activate the scale display.
- Access the Change Equipment mode using the MENU Down key, the display will read:
CHANGE EQUIP / PRESS ENTER
- Press and hold the ENTER key. The display will read:
RELEASE TO START (release the ENTER key), followed by:
DO NOT TOUCH BED (be sure nobody touches the bed while this message appears).
- The ADD/REMOVE EQUIP / PRESS ENTER message will then appear, signalling that equipment can now be added or removed from the bed.
When done, press and hold the ENTER key.
- Display will read RELEASE TO FIN. (release the ENTER key), followed by:
- DO NOT TOUCH BED (ensure the patient stands still as much as possible and no one touches the bed while this message is displayed).
- The scale will return to the Scale mode displaying the patient unchanged weight (the weight should be the same as the one displayed before the change of equipment).


Changing Manually the Numerical Value of Displayed Weight

In certain circumstances, it may be necessary to change manually the displayed patient's weight. To do so:

- Press the  key to activate the scale display.
- Access the Change Patient Weight mode using the MENU Down key. The display will read:
CHNG PTNT WEIGHT / PRESS ENTER
- Press and hold the ENTER key. The display will read:
REL. TO CHNG WT. (release the ENTER key), followed by:
DO NOT TOUCH BED (be sure no one touches the bed while this message appears).
- The display will then read ▲▼ XXX.X lb / press ENTER. Adjust the weight using the MENU Up/Down keys.
- Once done press the ENTER key.
- The scale will return to the Scale mode and will display the adjusted weight.

The Change Patient Weight mode should not be used systematically for new patients. The scale should be zeroed each time a new patient enters onto the bed (see page 18, "Operating (Zeroing) the..").

Changing the Weight Measure Unit

- Press the  key to activate the scale display. Access the Select Weight Unit (lb/kg) mode using the MENU Down key. The display will read:
UNITS / PRESS ENTER
- Press the ENTER key. The display will read:
▲▼ lb AND OR ▲▼ kg AND
PRESS ENTER PRESS ENTER
- Select the weight unit using the MENU Up or Down key and press the ENTER key. The scale will return to the Scale mode with the weight in the unit chosen.

2.17 CHAPERONE BED EXIT (OPTIONAL)




Figure 2.17

A: Press to activate/deactivate the Bed Exit function

B: Armed LED



WARNING

The Bed Exit system is intended only to aid in the detection of a patient exiting the bed. It is not intended to replace a patient monitoring protocol. The Bed Exit system signals when a patient is about to exit the bed. Adding or removing objects from the bed after having armed the Bed Exit system may cause a reduction in the sensitivity of the Bed Exit system, resulting in erroneous readings of the patient's movements in the bed.

The Bed Exit system is not designed to be used with patients weighing less than 50 lb (23 kg).

NOTE

For the bed exit alarm signal to be heard in the nurse station, the bed must be equipped with the optional nurse call function and be connected to the hospital nurse call system. For beds equipped with the Bed Exit function but not with the Nurse Call option, the alarm signal will be heard only in the patient's room.

Be sure to place the patient on the bed in its normal position before activating the bed exit.

Zeroing the Bed Exit System

The Bed Exit system must be zeroed prior to use. The following procedures will explain how to zero the Bed Exit system depending on the presence or not of the Weigh system (optional):

- **Bed with the Weigh system (optional)**

Simply zero the Weigh system through the "Operating (Zeroing) the Scale Before Putting a New Patient on the Bed" procedure described on page 18. The Bed Exit system will then be zeroed.

- **Bed without the Weigh system (optional)**


- Before putting a new patient on the bed, prepare the bed for a patient stay by adding linens and equipment to the bed.
- Press and hold the ARM/DISARM key (A) until the ARMED LED (B) start flashing, then release the key (do not touch the bed until the LED goes out completely). When the armed LED stops flashing, the system is zeroed.
- You can now place the patient in the bed. Ensure the patient lies in the bed in his usual position before activating the Bed Exit system.

Adding or Removing Equipment

If equipment is added or removed from the bed and the Bed Exit system is activated, remove the patient from the bed, add or remove the equipment and zero the system as described in "Zeroing the Bed Exit System" above. Only then can the patient reintegrate the bed.

To **activate** the Bed Exit, press the ARM/DISARM key (A). The ARMED LED (B) will come on.

To **deactivate** the Bed Exit after an alarm or simply to turn it off, press the ARM/DISARM key (A). The Armed LED (B) will turn off.

2.18 CHAPERONE BED EXIT WITH ZONE CONTROL (OPTIONAL)**Figure 2.18**

A: Press to activate/deactivate the Bed Exit function

B: Zone LEDs

C: Press to select zone

D: Zone icons (zone from left to right: less restrictive to most restrictive, see below)

NOTE

For the bed exit alarm signal to be heard in the nurse station, the bed must be equipped with the optional nurse call function and be connected to the hospital nurse call system. For beds equipped with the Bed Exit function but not with the Nurse Call option, the alarm signal will be heard only in the patient's room.

Be sure to place the patient on the bed in its normal position before activating the bed exit.

**WARNING**

The Bed Exit system is intended only to aid in the detection of a patient exiting the bed. It is not intended to replace patient monitoring protocol. The Bed Exit system signals when a patient is about to exit the bed. Adding or removing objects from the bed after having armed the Bed Exit system may cause a reduction in the sensitivity of the Bed Exit system, resulting in erroneous readings of the patient's movements in the bed. Failure to set the Chaperone zone properly could result in improper monitoring and patient injury. Verify that the proper zone is selected before leaving the patient.

The Bed Exit system is not designed to be used with patients weighing less than 50 lb (23 kg).

Zone Settings

- The first zone (left LED) is the traditional Bed Exit zone. The patient can move in the bed freely but cannot fully exit the bed or the alarm will sound.
- The second zone (middle LED) is more restrictive. It allows the patient to sit up and roll over but any attempt to exit the bed will cause the alarm to sound.
- The third zone (right LED) is the most restrictive zone. Small movements like raising an arm or lifting the shoulder off the bed will cause the alarm to sound. The third zone is used to alert staff to a change in the condition of an unconscious or paralyzed patient.

Zeroing the Bed Exit System

The Bed Exit system must be zeroed prior to use. The following procedures will explain how to zero the Bed Exit system depending on the presence or not of the Weigh system (optional):

- **Bed with the Weigh system (optional)**

Simply zero the Weigh system through the "Operating (Zeroing) the Scale Before Putting a New Patient on the Bed" procedure described on page 18. The Bed Exit system will then be zeroed.

- **Bed without the Weigh system (optional)**

- Before putting a new patient on the bed, prepare the bed for a patient stay by adding linens and equipment to the bed.
- Press and hold the ARM/DISARM key (A) until the ARMED LED (B) starts flashing, then release the key (do not touch the bed until the LED goes out completely). When the armed LED stops flashing, the system is zeroed.
- You can now place the patient in the bed. Ensure the patient lies in the bed in his usual position before activating the Bed Exit system.

Adding or Removing Equipment

If equipment is added or removed from the bed and the Bed Exit system is activated, remove the patient from the bed, add or remove the equipment and zero the system as described in "Zeroing the Bed Exit System" above. Only then can the patient reintegrate the bed.

To **activate** the Bed Exit, press the ARM/DISARM key (A). The ARMED LED (B) will come on.

To **deactivate** the Bed Exit after an alarm or simply to turn it off, press the ARM/DISARM key (A). The Armed LED (B) will turn off.

3. ACCESSORIES

Listed below are the accessories certified compliant (IEC 60601-2-38) for use with the GOBED⁺ as well as their maximum load capacity, where applicable.

- Bed extension