GoBed® II MedSurg Bed  
GoBed® II Lit MedSurg  
Model FL28C / Modèle FL28C

# stryker®

## Maintenance Manual

Guide d'entretien

For Parts or Technical Assistance:  
USA: 1-800-327-0770

CE

# Table of Contents

Symbols and Definitions . . . . .	<a href="#">1-4</a>
Warning / Caution / Note Definition . . . . .	<a href="#">1-5</a>
Introduction . . . . .	<a href="#">1-6</a>
Intended Use of Product . . . . .	<a href="#">1-6</a>
Specifications . . . . .	<a href="#">1-6</a>
Environmental Conditions . . . . .	<a href="#">1-7</a>
Product Illustration . . . . .	<a href="#">1-8</a>
Standard features . . . . .	<a href="#">1-9</a>
Optional Features . . . . .	<a href="#">1-9</a>
Accessories . . . . .	<a href="#">1-9</a>
Safety Tips and Guidelines . . . . .	<a href="#">1-10</a>
Static Discharge Precautions . . . . .	<a href="#">1-13</a>
Cleaning . . . . .	<a href="#">1-14</a>
Preventative Maintenance . . . . .	<a href="#">1-15</a>
Checklist . . . . .	<a href="#">1-15</a>
Lubrication Requirements . . . . .	<a href="#">1-16</a>
Troubleshooting Guide . . . . .	<a href="#">1-18</a>
Quick Reference Replacement Parts List . . . . .	<a href="#">1-20</a>
Electrical Service Information . . . . .	<a href="#">1-22</a>
Motor Control Board . . . . .	<a href="#">1-22</a>
Power Supply . . . . .	<a href="#">1-23</a>
Communication Interface (Optional) . . . . .	<a href="#">1-24</a>
Scale (Optional) . . . . .	<a href="#">1-25</a>
Maintenance Menu: Error Codes . . . . .	<a href="#">1-26</a>
Service Information . . . . .	<a href="#">1-27</a>
Maintenance Procedures . . . . .	<a href="#">1-27</a>
Head End Siderail Assembly Replacement . . . . .	<a href="#">1-28</a>
Foot End Siderail Assembly Replacement . . . . .	<a href="#">1-29</a>
Head End Siderail Outer Panel Membrane Replacement . . . . .	<a href="#">1-30</a>
Speaker/Siderail Inner Panel Membrane Replacement . . . . .	<a href="#">1-31</a>
Footboard Control Membrane Replacement . . . . .	<a href="#">1-33</a>
Scale System LCD Display Replacement (Optional) . . . . .	<a href="#">1-34</a>
Foot Section Replacement . . . . .	<a href="#">1-35</a>
Thigh Section Replacement . . . . .	<a href="#">1-36</a>
Seat Section Replacement . . . . .	<a href="#">1-37</a>
Fowler Replacement . . . . .	<a href="#">1-38</a>
Fowler Motor Replacement . . . . .	<a href="#">1-40</a>
Gatch Motor Replacement . . . . .	<a href="#">1-42</a>
Head End Lift Motor Replacement . . . . .	<a href="#">1-44</a>
Foot End Lift Motor Replacement . . . . .	<a href="#">1-45</a>
Control Board Replacement . . . . .	<a href="#">1-47</a>
Headwall Control Board Replacement (Optional) . . . . .	<a href="#">1-48</a>
Power Cord Connector Replacement . . . . .	<a href="#">1-49</a>
Power Cord Connector Fuse Replacement . . . . .	<a href="#">1-50</a>
Scale/Bed Exit Board/Beeper Replacement (Optional) . . . . .	<a href="#">1-51</a>

# Table of Contents

Service Information (Continued)	
Transformer Replacement (100/200/220/240 VAC Options) . . . . .	<a href="#">1-52</a>
Load Cell Replacement (Optional) . . . . .	<a href="#">1-53</a>
Scale Calibration (Optional) . . . . .	<a href="#">1-55</a>
Auto Contour Micro Switch Replacement . . . . .	<a href="#">1-57</a>
CPR Micro Switch Replacement (Optional) . . . . .	<a href="#">1-59</a>
CPR Activation Cable Replacement (Optional) . . . . .	<a href="#">1-60</a>
CPR Pneumatic Cylinder (Optional) . . . . .	<a href="#">1-61</a>
Brake/Steer Pedal Replacement . . . . .	<a href="#">1-62</a>
Steer Wheel Caster Replacement . . . . .	<a href="#">1-63</a>
Steer Wheel Assembly Replacement . . . . .	<a href="#">1-64</a>
Steer Wheel Swing Arm Assembly Replacement . . . . .	<a href="#">1-65</a>
Bed Caster Replacement . . . . .	<a href="#">1-66</a>
Brake Rod Replacement . . . . .	<a href="#">1-67</a>
Brake System . . . . .	<a href="#">1-68</a>
Fifth Wheel Assembly . . . . .	<a href="#">1-70</a>
Caster Hood Assembly . . . . .	<a href="#">1-72</a>
Base Tube Cover Assembly . . . . .	<a href="#">1-73</a>
Lift Mechanism . . . . .	<a href="#">1-74</a>
Lower Stationary Frame with Scale Option, without Network . . . . .	<a href="#">1-76</a>
Upper Stationary Frame without Scale Option . . . . .	<a href="#">1-78</a>
Upper Stationary Frame with Scale Option . . . . .	<a href="#">1-80</a>
Moving Frame Assembly . . . . .	<a href="#">1-82</a>
Auto Contour without Scale Option . . . . .	<a href="#">1-84</a>
Auto Contour with Scale Option . . . . .	<a href="#">1-86</a>
120 VAC Electrical System . . . . .	<a href="#">1-88</a>
Power Cord Options . . . . .	<a href="#">1-93</a>
Straight Plug - QDF8066 . . . . .	<a href="#">1-93</a>
90° Plug - QDF8066-90D . . . . .	<a href="#">1-93</a>
International Electrical System . . . . .	<a href="#">1-94</a>
220VAC Electrical System . . . . .	<a href="#">1-98</a>
International Power Cord Options . . . . .	<a href="#">1-99</a>
Head End Cover without Scale, Nurse Call and Pendant Port . . . . .	<a href="#">1-100</a>
Head End Cover with Nurse Call and Pendant Port, without Scale . . . . .	<a href="#">1-101</a>
Head End Cover with Pendant Port, without Scale and Nurse Call . . . . .	<a href="#">1-102</a>
Head End Cover with Nurse Call/Scale/GENIII/37 Pin Connectors . . . . .	<a href="#">1-103</a>
120 VAC Auxiliary Outlet . . . . .	<a href="#">1-104</a>
Night Light . . . . .	<a href="#">1-105</a>
CPR Mechanism . . . . .	<a href="#">1-106</a>
Foot Section with Support Rod . . . . .	<a href="#">1-110</a>
Foot End Mattress Retainer . . . . .	<a href="#">1-111</a>
Litter with Plastic Covers . . . . .	<a href="#">1-112</a>
Litter with Steel Covers . . . . .	<a href="#">1-114</a>
Head End Siderail Assembly . . . . .	<a href="#">1-116</a>
Head End Siderail Assembly, Right . . . . .	<a href="#">1-117</a>

# Table of Contents

Assembly Drawings (Continued)	
Siderail Assembly, Reverse - 28-0006	<a href="#">1-119</a>
Head End Siderail Assembly, Left	<a href="#">1-120</a>
Siderail Assembly - 28-0005	<a href="#">1-122</a>
Head End Siderail Labels without Nurse Call	<a href="#">1-123</a>
Head End Siderails with Nurse Call	<a href="#">1-124</a>
Foot End Siderail Assembly	<a href="#">1-125</a>
Foot End Siderail Assembly, Right	<a href="#">1-126</a>
Foot End Siderail Assembly, Left	<a href="#">1-128</a>
Headboard Assembly	<a href="#">1-130</a>
Footboard Assembly	<a href="#">1-131</a>
Footboard without Scale and Bed Exit	<a href="#">1-132</a>
Footboard with Scale, without Bed Exit	<a href="#">1-134</a>
Footboard with Scale and Bed Exit - Three Zone	<a href="#">1-136</a>
Bed Labeling - English	<a href="#">1-138</a>
Bed Labeling - English / French / Spanish	<a href="#">1-140</a>
Bed Labeling - English / Spanish	<a href="#">1-142</a>
Optional Accessories	<a href="#">1-144</a>
Pendant Assemblies	<a href="#">1-145</a>
Monitor Tray - FA64163	<a href="#">1-146</a>
Patient Helper - FA64148	<a href="#">1-147</a>
Emergency Crank Handle - FOHMAU	<a href="#">1-148</a>
Oxygen Bottle Holder - FA64169	<a href="#">1-149</a>
Bed Extender without Scale - FA64164	<a href="#">1-150</a>
Bed Extender with Scale - FA64173	<a href="#">1-151</a>
Removable I.V. Pole - 1/2" - FDTSH	<a href="#">1-152</a>
Removable I.V. Pole - 1" - FA61002	<a href="#">1-153</a>
Fixed Two Stage I.V. Pole, Foldable - FA64171	<a href="#">1-154</a>
Fixed Three Stage I.V. Pole, Foldable - FA64170	<a href="#">1-155</a>
Traction Sleeves	<a href="#">1-156</a>
Recycling Passport	<a href="#">1-157</a>
EMC Information	<a href="#">1-162</a>
Warranty	<a href="#">1-166</a>
Limited Warranty	<a href="#">1-166</a>
To Obtain Parts and Service	<a href="#">1-166</a>
Service Contract Coverage	<a href="#">1-166</a>
Service Contract Programs	<a href="#">1-167</a>
Return Authorization	<a href="#">1-167</a>
Damaged Merchandise	<a href="#">1-167</a>
International Warranty Clause	<a href="#">1-167</a>

# Symbols and Definitions

Warning / Caution, Refer to Accompanying Documentation

Alternating Current

Direct Current

Type B Equipment: equipment providing a particular degree of protection against electric shock, particularly regarding allowable leakage current and reliability of the protective earth connection.

Class 1 Equipment: equipment in which protection against electric shock does not rely on **basic insulation** only, but which includes an additional safety precaution in that means are provided for the connection of the **equipment** to the protective earth conductor in the fixed wiring of the installation in such a way that **accessible metal parts** cannot become live in the event of a failure of the **basic insulation**.

CE Mark – Not applicable with models with scale

**IPX4**

Protection from liquid splash

Dangerous Voltage Symbol

Protective Earth Terminal

Potential Equalization Symbol

Medical Equipment Classified by Underwriters Laboratories Inc. with respect to Electric Shock, Fire, Mechanical and Other Specified Hazards Only in Accordance with UL 60601-1, First Edition (2003) and CAN/CSA C22.2 No. 601.1-M90 with updates 1 and 2.

Safe Working Load Symbol

In accordance with **European Directive 2002/96/EC** on Waste Electrical and Electronic Equipment (**WEEE**), this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.

# Symbols and Definitions

---

## **WARNING / CAUTION / NOTE DEFINITION**

The words WARNING, CAUTION, and NOTE carry special meanings and should be carefully reviewed.

---

### **WARNING**

Alerts the reader about a situation, which if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

---

### **CAUTION**

Alerts the reader of a potentially hazardous situation, which if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

---

### **Note**

This provides special information to make maintenance easier or important instructions clearer.

# Introduction

This manual is designed to assist you with the maintenance of GoBed® II MedSurg Bed, Model FL28C. Carefully read this manual thoroughly before using the equipment or beginning maintenance on it. To ensure safe operation of this equipment, it is recommended that methods and procedures be established for educating and training staff on the safe operation of this bed.

## INTENDED USE OF PRODUCT

The GoBed® II MedSurg Bed, Model FL28C is an AC-powered adjustable hospital bed designed to position patients for procedures, therapy, and recovery in healthcare environment, transport patients between bays and procedural rooms, provide patient security, measure and display patient weight, allow patients to alert caregiver when the patient requires emergency assistance or any assistance to improve comfort levels, and monitor patient position to alert caregiver of potential fall.

Intended operators are health care professionals (nurses, nurse aids, doctors) which can use all bed operations (e.g., bed motion functions, nurse call, siderail communications, bed exit), patient and bystander, which can use bed motion functions, nurse call & siderail communications, and trained professionals for installation, service, and calibration.

The product is intended to be used in a healthcare environment, including hospitals, surgery centers, long term acute care centers, and rehabilitation centers. The intended environmental conditions are 10 to 40 degrees Celsius and 30 to 75% RH. The product is compatible with 35" x 84" surfaces, the facility nurse call system, standard Med-Surg equipment, and the facility infrastructure. Its expected life is 10 years.

Intended patient population are general hospital and acute care patients. The bed accommodates a safe working load of 500 pounds. Patients must weigh more than 50 pounds to use the bed exit system. It is not intended for pediatric use and/or patients measuring 35" or less.

## SPECIFICATIONS

	<p>Safe Working Load  <b>Note:</b> Safe Working Load indicates the sum of the patient, mattress, and accessory weight.</p>	500 lb	227 kg
<p>Scale System (Optional)                      - Accuracy of scale display*                       - Operating Angular Range</p>		<p>± 2% of the total patient weight for patient weighing greater than 100 pounds (45,3 kg)                      ± 2 lbs below 100 lb (45,3 kg)                      -12° to +12°</p>	
<p>* To meet this accuracy claim, the patient surface must be in the flat position (fowler and gatch down) and the stretcher cannot exceed 5 degrees of Trendelenburg/reverse Trendelenburg.</p>			
<p>Overall Length/Width                      - Siderails Up                      - Siderails Down</p>		<p>94 5/8" x 40" (240,34 cm x 101,6 cm)                      94 5/8" x 39" (240,34 cm x 99,06 cm)</p>	
<p>Weight with Headboard and Footboard</p>		<p>472 lb (214,1 kg)</p>	
<p>Patient Sleep Surface</p>		<p>35 x 84" (89 x 213 cm)</p>	
<p>Recommended Mattress Size                      Mattress <b>Maximum</b> Thickness</p>		<p>35 x 84" (89 x 213 cm)                      6" (15,24 cm)</p>	
<p>Minimum/Maximum Bed Height</p>		<p>14 1/2" to 29" (36,8 to 73,7 cm)</p>	
<p>Fowler Angle</p>		<p>0° to 59°</p>	

# Introduction




## SPECIFICATIONS (CONTINUED)

Gatch Angle - With Auto Contour (Optional) - Without Auto Contour	0° to 24° 0° to 32°
Trendelenburg/Reverse Trendelenburg	+14° to -14°
Electrical Requirements - all electrical requirements meet CSA C22.2 No. 601.1, UL 60601-1 and IEC 60601-1, 60601-2-38 specifications.	100Vac, 50-60Hz, 7.5A - Two 250Vac, 10A Fuses 120Vac, 50-60Hz, 4A (9.8A with 120Vac Optional Auxiliary Outlet) - Two 250V, 10A Fuses 200Vac, 50-60Hz, 3.2A -Two 250Vac, 6.3A Fuses 220Vac, 50-60Hz, 2.9A -Two 250Vac, 6.3A Fuses 240Vac, 50-60Hz, 2.7A -Two 250Vac, 6.3A Fuses
Duty Cycle	10%

**Stryker reserves the right to change specifications without notice.**

- Specifications listed are approximate and may vary slightly from unit to unit or by power supply fluctuations.

## ENVIRONMENTAL CONDITIONS

Environmental Conditions	Operation	Storage and Transportation
Ambient Temperature		
Relative Humidity (Non-Condensing)		
Atmospheric Pressure		

**Stryker reserves the right to change specifications without notice.**

- Specifications listed are approximate and may vary slightly from unit to unit.
- Operating environment recommended to ensure the scale system (optional) precision.

# Introduction

## PRODUCT ILLUSTRATION


Figure 1.9

## STANDARD FEATURES

- 14.5" (37 cm) low bed height
- Retractable fifth-wheel steering
- Retractable sleep surface
- 6" (15 cm) casters
- Trendelenburg/Reverse Trendelenburg
- Lateral and Footend Mattress Retainers
- Four independent electric motors
- Trendelenburg display on footboard (only available with scale option)
- Centrally located steer and four wheel brake mechanism
- Four foley bag hooks
- Eight I.V. pole/traction equipment sockets
- Roller bumpers
- Split Siderails
- Manual back up for head and knee control
- Integrated pump holder
- Nurse controls on footboard and siderails
- Fixed patient controls on siderails
- CPR Release
- Degree indicator for head elevation
- Patient restrain locations
- Electronic function lockout controls

## OPTIONAL FEATURES

- Chaperone center-of-gravity bed exit system
- Chaperone with Zone Control<sup>®</sup>
- In-bed scale system
- Smart TV - includes channel up/down
- Mattresses: management and prevention
- Pillow speaker interface
- Removable litter covers

## ACCESSORIES

- Monitor tray
- Patient helper
- Emergency crank handle
- Oxygen bottle holder
- Bed extender and pad (with and without scale option)
- Removable I.V. poles
- Fixed I.V. pole
- Traction sleeves and adapters

# Safety Tips and Guidelines

Before operating the GoBed® II MedSurg Bed, Model FL28C, it is important to read and understand all information in this manual. Carefully read and strictly follow the safety guidelines listed below. It is important that all users have been trained and educated on the inherent hazards associated with the usage of electric beds.

## WARNING

- Always unplug the bed power cord from the power source before moving the bed.
- This bed is not intended for pediatric use, i.e., for any patient measuring 35 inches (88,9 cm) or less.
- The maximum load capacity of the monitor tray is 40 lb (18 kg).
- The weight of the I.V. bags must not exceed 11 lb (5 kg) per hook.
- The weight of the I.V. bags must not exceed 40 lb (18 kg).
- The maximum load capacity of the oxygen bottle holder is 75 lb (34 kg).
- The maximum load capacity of the overhead trapeze is 150 lb (68 kg).
- The overhead trapeze is designed to help the patient change position in bed. The trapeze must not be used as a support to enter or exit the bed. Injuries may result due to improper use of the overhead trapeze. Medical advice may be necessary to determine if the use of an overhead trapeze is adequate for the patient's condition or treatment given.
- To prevent injury to the patient, ensure the trapeze is properly attached to the upper section.
- The overhead trapeze should not be used when in the storage position. This may result in injury to the patient or damage to the equipment.
- This bed is equipped with a hospital grade plug for protection against shock hazard. It must be plugged directly into a properly grounded receptacle. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Shock Hazard - Improper handling of the power cord may result in damage to the power cord and potential shock hazards. If damage has occurred to the power cord, immediately remove the bed from service, and contact the appropriate maintenance personnel. Failure to do so could result in serious injury or death.
- Use only hospital grade electric equipment consuming 5A or less with the auxiliary power outlet (optional). The use of standard electric equipment may bring the current leakage to an unacceptable level for hospital equipment.
- Serious injury can result if caution is not used when operating the bed. Operate the bed only when all people and equipment are clear of the electrical and mechanical systems.
- Before using the emergency crank during a power failure, always unplug the unit. An unexpected return of power could rotate the handle and cause injury.
- Ensure all the connectors stay plugged into the electronic board if it is present on the head casing cover.
- Always apply the brakes when a patient is on the bed (except during transport) or when entering/exiting the bed. Serious injury could result if the bed moves while a patient is getting on or off the bed. After the brake pedal is engaged, push on the bed to ensure the brakes are applied.
- To help reduce the number and severity of a potential fall when the patient is unattended, keep the siderails in the fully raised position and the sleep surface horizontal in its lowest position, unless its medical condition dictates otherwise. When raising the siderails, be sure that you hear the "click" that signals the locked condition. Pull firmly on the siderail to ensure it is locked into position.
- When the sleep surface sections are articulated, ensure that all the patient's limbs are within the raised siderails to avoid patient injury.
- When a patient's condition requires greater safety measures for his/her security, use the lockout controls in the foot board control panel to inhibit the siderail functions or remove any optional pendant control and install protective pads on the siderails.
- Improperly adjusted restraint straps can cause serious injury to a patient. It is the responsibility of the attending medical personnel to determine proper use of the restraint straps.
- Siderails, with or without their padded covers, are not intended to serve as restraint devices to keep patient from exiting the bed. Siderails are designed to keep a patient from inadvertently rolling off the bed. It is the responsibility of the attending medical personnel to determine the degree of restraint necessary to ensure a patient will remain safely in bed. Failure to use the siderails properly could result in serious patient injury.
- To reduce risk of injury, ensure the sleep surface is horizontal and in the lowest position with the siderails fully raised and locked when moving the bed with a patient in it.

# Safety Tips and Guidelines

## WARNING (CONTINUED)

- To avoid injury to the patient and/or user, do not attempt to move the bed laterally with the steer mode engaged. The steer wheel cannot swivel.
- The CPR emergency release is for emergency use only. When activating the CPR release handle, all people and equipment must be removed from the area below and around the head, thigh and foot sections of the bed or serious personal injury and/or equipment damage could occur.
- Possible fire hazard exists when this unit is used with oxygen administering equipment other than nasal, mask type or half bed-length tent type. Unplug the unit power cord from the wall when oxygen-administering equipment is used. When using a half bed-length tent type, ensure that the siderails are outside the oxygen tent and that the tent does not extend below the mattress support level.
- The weight system (optional) is intended to assist in the monitoring of the patient's weight variation. Under no circumstances should its readings be used as sole reference for medical treatment.
- The Bed Exit system (optional) is intended only to aid in the detection of a patient exiting the bed. It is not intended to replace patient monitoring protocol. The Bed Exit system signals when a patient is about to exit the bed. The addition or removal of equipment with a Bed Exit system armed must be done using the "Adding or Removing Equipment with the System Armed" procedure, otherwise the sensitivity of the system may be affected and the readings of the patient's movements in the bed will be erroneous.
- The surface added by the bed extender is not included in the area monitored by the Bed Exit system.
- The Bed Exit system (optional) is not designed to be used with patients weighing less than 50 lb (23 kg).
- Failure to set the zone properly could result in improper monitoring and patient injury. Verify that the proper zone is selected before leaving the patient.
- The scale and bed exit systems must be adjusted if an accessory is added when either system is in operation. Refer to the "Add / Remove Equipment" procedure.
- The bed extender must be removed before using the emergency crank (optional).
- When large fluid spills occur in the area of the circuit boards, cables and motors, immediately unplug the bed power cord from the wall outlet, remove the patient from the bed and clean up the fluid. Have maintenance completely check the bed. Fluids can have an adverse effect on operational capabilities of any electrical product. **Do not** put the bed back into service until it is completely dried and has been thoroughly tested for safe operation. Ensure, among other things, that the plastic components being used as covers for the siderail mechanism arms and the foot end casing are removed and that the parts they cover are thoroughly dried.
- Preventative maintenance should be performed at least once a year to ensure all bed features are functioning properly. Ensure that any bed malfunction is promptly reported to your service personnel for immediate attention.
- Always unplug the bed power cord from the wall outlet when servicing or cleaning the bed.
- When working under the bed with the bed in the high position, always apply the brakes and place blocks under the Hi-Low levers to prevent injury in case the bed down switch is accidentally pressed.
- Before using the optional emergency crank during a power failure, always unplug the power cord. An unexpected return of the power could rotate the handle and cause injury to the user.
- Because individual beds may have different options, foot boards should not be moved from one bed to another. Mixing foot boards could result in unpredictable bed operation.
- When servicing use only identical replacement parts provided by Stryker.
- Make sure that the ground chain is in place, intact and it's touching the floor (see Figure 1.9, [page 1-8](#)).
- Do not put anything under the bed.
- Only service technicians from Stryker or service personnel trained by Stryker should perform the procedures detailed in this maintenance manual. This is especially important for procedures related to the Scale (option) and Bed Exit (option) systems. Failure to follow this restriction can result in serious damage to the bed and/or severe injury to users or patients.
- To prevent injury when working under the bed with the bed in the high position, always place blocks under the mattress support frame and apply the brakes.

# Safety Tips and Guidelines

---

## CAUTION

- Do not use the monitor tray as a push/pull device.
  - Do not use the I.V. pole as a push/pull device.
  - Do not use an I.V. pole at the foot end of a unit equipped with this device.
  - The mattress thickness should not exceed 6 inches (15,24 cm).
  - When using a mattress thicker than 6 inches (15,24 cm) or a mattress overlay, extra caution and/or patient supervision may be required to reduce the likelihood of occurrence of a patient fall.
  - To avoid damage to the siderail mechanisms, do not move the bed using the raised siderails. Use the push/pull handles integrated to the footboard and headboard to move the bed.
  - Do not steam clean, hose off or ultrasonically clean the bed. Do not immerse any part of the bed in any kind of liquid. The internal electrical parts may be damaged by exposure to water. Hand wash regularly all surfaces of the bed with warm water and a mild detergent. Wipe cleaned surfaces dry to avoid build up of cleaning substance. Inspect the mattress after each use. Discontinue use if any cracks or rips, which may allow fluid to enter the mattress, are found in the mattress cover. Failure to properly clean the mattress, or dispose of it if defective, may increase the risk of exposure to pathogenic substances and may bring about diseases to the patient and user.
  - **SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY.** If the products suggested above are used to clean Stryker patient care equipment, measures must be taken to insure the bed is wiped with a damp cloth soaked in clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the bed will leave a corrosive residue on the surface of the bed, possibly causing premature corrosion of critical components. Failure to follow the above directions when using these types of cleaners may void this product's warranty.
- 

### Note

Throughout this Maintenance Manual, the words “right” and “left” refer to the right and left sides of a patient lying face up on the bed.

# Static Discharge Precautions

The electronic circuits in the GoBed® II MedSurg Bed, Model FL28C are completely protected from static electricity damage only while the bed is assembled. It is extremely important that all service personnel always use adequate static protection when servicing the electronic systems of the bed. Whenever you are touching wires, you should be using static protection.

## Static Protection Equipment

The necessary equipment for proper static protection is:


- 1 static wrist strap; 3M part number 2214 or equivalent,
- 1 grounding plug; 3M part number 61038 or equivalent,
- 1 test lead with a banana plug on one end and an alligator clip on the other; Smith part number N132B699 or equivalent.

## CAUTION

All electronic service parts will be shipped in static shielding bags. Do not open the bags until you have completed steps 2 and 3 of the following procedure. Do not place unprotected circuit boards on the floor. All circuit boards to be returned to Stryker Medical should be shipped in the static shielding bags the new boards were shipped in.

## Static Protection Procedure

1. Unplug the power cord from the wall receptacle.
2. Insert the grounding plug into a properly grounded hospital grade wall receptacle. Plug the banana plug of the test lead into the receptacle on the grounding plug. Connect the alligator clip on the other end of the test lead to a ground point on the bed.
3. Place the static control wrist strap on your wrist. Connect the alligator clip at the other end of the wrist strap cord to a ground point on the bed.

# Cleaning

---

---

## WARNING

Always unplug the bed power cord from the wall outlet when servicing or cleaning the bed.

---

Hand wash all surfaces of the GoBed® II MedSurg Bed, Model FL28C with warm water and mild detergent. **Dry thoroughly.** Do not steam clean or hose off the bed. Using these methods of cleaning is **not** recommended and may void this product's warranty. Do not immerse any part of the bed in any kind of liquid. Some of the internal parts of the bed are electric and may be damaged by exposure to water.

Suggested cleaners for bed surfaces:

- Quaternary Cleaners (active ingredient - ammonium chloride).
- Phenolic Cleaners (active ingredient - o-phenylphenol).
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water).

Avoid over saturation and ensure the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

---

## CAUTION

Do not steam clean, hose off or ultrasonically clean the bed. Do not immerse any part of the bed in any kind of liquid. The internal electrical parts may be damaged by exposure to water. Hand wash regularly all surfaces of the bed with warm water and a mild detergent. Wipe cleaned surfaces dry to avoid build up of cleaning substance. Inspect the mattress after each use. Discontinue use if any cracks or rips, which may allow fluid to enter the mattress, are found in the mattress cover. Failure to properly clean the mattress, or dispose of it if defective, may increase the risk of exposure to pathogenic substances and may bring about diseases to the patient and user.

---

## CAUTION

SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products suggested above are used to clean Stryker patient care equipment, measures must be taken to insure the bed is wiped with a damp cloth soaked in clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the bed will leave a corrosive residue on the surface of the bed, possibly causing premature corrosion of critical components. Failure to follow the above directions when using these types of cleaners may void this product's warranty.

---

For mattress cleaning instructions, please see the tag on the mattress, or contact the mattress manufacturer.

Clean Velcro® **after each use.** Saturate Velcro® with disinfectant and allow disinfectant to evaporate. (Appropriate disinfectant for nylon Velcro® should be determined by the hospital).

# Preventative Maintenance

Beds require an effective maintenance program. We recommend checking these items annually. Use this sheet for your records and keep on file.

## CHECKLIST

- \_\_\_\_\_ All fasteners secure (reference all assembly drawings).
- \_\_\_\_\_ Inspect for excessive wear on the oil-impregnated bronze shoulder spacers found at the bed hinge points. Do not lubricate these spacers. Replace as needed.
- \_\_\_\_\_ Check that grease is present on the components detailed in section “Lubrication Requirements”, lubricate if necessary. Lubricate at least every two years.
- \_\_\_\_\_ On both sides of the bed, depress the side of the pedal identified with a red arrow and ensure that the brakes are applied and the bed is immobilized. Toggle the pedal to neutral and ensure the brakes release. On both sides of the bed, depress the side of the pedal identified with a green arrow and ensure that the steer wheel is engaged. Toggle the pedal to neutral and ensure that the steer wheel disengages.
- \_\_\_\_\_ Siderails move, latch and stow properly.
- \_\_\_\_\_ All controls of the foot end panel working properly, including LEDs.
- \_\_\_\_\_ Calibrate the scale (optional). Refer to the “Scale Calibration Procedure”.
- \_\_\_\_\_ All siderail controls working properly.
- \_\_\_\_\_ Ensure that the nurse call (optional) alarm sounds in the nurse station.
- \_\_\_\_\_ Verify the CPR emergency release using both CPR release handles: raise the Fowler fully up and, using the CPR handle, lower the Fowler gradually to flat position by pulling, holding and releasing the handle several times. Ensure the Gatch (if raised) also starts flattening when the Fowler is completely down. Following the complete lowering of the Fowler, wait approximately 30 seconds - the time for the Fowler actuator to reset - and verify that the actuator has indeed reset itself by raising the Fowler fully up.
- \_\_\_\_\_ Verify the Fowler, Gatch and Hi-Lo movements to ensure that the motion interrupt switch integrated to the four electric motors is operating properly.
- \_\_\_\_\_ Auxiliary outlet (option available only with 120VAC beds) working properly.
- \_\_\_\_\_ Night light working properly.
- \_\_\_\_\_ No cracks in the headboard or footboard, siderails, wheel covers, steer wheel hood (optional) and plastic covers (optional) of the head and foot sections.
- \_\_\_\_\_ Head end bumpers tightly secured to frame and working properly.
- \_\_\_\_\_ No rips or cracks in mattress cover. Replace if so.
- \_\_\_\_\_ Power cord intact.
- \_\_\_\_\_ No cables pinched or worn.
- \_\_\_\_\_ All electrical connections tight.
- \_\_\_\_\_ All grounds secure to the frame.
- \_\_\_\_\_ All casters roll properly. Check caster for cuts, wear, etc.
- \_\_\_\_\_ Ground Impedance not more than 100 mΩ (milliohms).
- \_\_\_\_\_ Current leakage not more than 100 μA (microamps).
- \_\_\_\_\_ Make sure that the ground chain is in place, intact and it’s touching the floor (see Figure 1.9, [page 1-8](#)).
- \_\_\_\_\_ Check labels as specified in the Operations and Maintenance manuals for legibility, proper adherence and integrity.

Bed Serial Number:		

Completed by: \_\_\_\_\_

Date: \_\_\_\_\_

[Return To Table of Contents](#)

# Preventative Maintenance

---

## LUBRICATION REQUIREMENTS

The only components of the bed needing periodic lubrication are the four actuator screws and the clevis pins holding the head actuator to the head section lever. These components should be inspected annually and lubricated every two years.

---

### CAUTION

The bed uses oil impregnated shoulder spacers at hinge points. DO NOT lubricate these shoulder spacers. When shoulder spacers become worn, replace them.

---

## ACTUATOR SCREW LUBRICATION

### Tools Required:

- #2 Phillips Screwdriver
- Brush
- Bungee Cord (or equivalent)
- OG2 Grease (or equivalent)

### Procedure:

#### Clevis Pin

1. Raise the bed to the highest position and apply the brakes.
2. Remove the plastic head section cover (if equipped). Raise the head section to the highest position.
3. Remove the headboard.
4. Unplug the power cord from the wall receptacle.
5. Apply a coat of OG2 grease to the clevis pin (A).
6. Raise and lower the head section several times to ensure the clevis pin is properly coated.

### Procedure:

#### Head and Thigh Actuators

1. Raise the bed to the highest position and apply the brakes.
2. Remove the plastic head and foot section covers (if equipped).
3. Raise the head section to the highest position and flatten the thigh section. Lower the four siderails.
4. Unplug the power cord from the wall receptacle.
5. Lift the foot section and fold it back towards the head end of the bed. Secure it using a bungee cord.
6. Using a #2 Phillips screwdriver, remove the four screws (B) holding the cover plate to the frame and remove the plate. Remove the night light if equipped.
7. Using a #2 Phillips screwdriver, remove the two screws (C) holding the dust tube of each actuator. Push the tube back to uncover the screw threads.
8. Using a brush, apply OG2 grease to the screw threads. Ensure the grease penetrates to the bottom of the threads.
9. Replace the dust tubes and cover plate, and lower the foot section.
10. Plug the power cord into the wall receptacle.
11. Using the electronic controls, raise and lower the head and thigh section several times to spread the grease evenly.

## ACTUATOR SCREW LUBRICATION (CONTINUED)

### Procedure:

#### Hi-Lo Actuators

1. Lower the bed to the lowest position and apply the brakes. Flatten the litter.
2. Lower the four siderails.
3. Using a #2 Phillips screwdriver, remove the head and foot section plastic covers (if equipped), to reach the Hi-Lo actuators, or raise the head section to the highest position and fold the foot section back towards the head end of the bed. Secure the foot section using a bungee cord.
4. Unplug the power cord from the wall receptacle.
5. Using a brush, apply OG2 grease to the screw threads of both actuators through the openings on the sides of the dust tubes. Ensure the grease penetrates to the bottom of the threads.
6. Lower the foot section.
7. Plug the power cord into the wall receptacle.
8. Using the electronic controls, raise and lower the bed several times to spread the grease evenly.

# Troubleshooting Guide

Problem / Failure	Recommended Action
No power to bed.	<ul style="list-style-type: none"> <li>• Verify power at wall receptacle.</li> <li>• Verify the power cord connections at the wall receptacle and the bed.</li> <li>• Check the power cord for damage. Replace if necessary.</li> <li>• Check the two fuses inside the power connector.</li> </ul>
No bed up or down motion.	<ul style="list-style-type: none"> <li>• Ensure Total Lockout (padlock LED) is not enabled.</li> <li>• Check the siderail control panel cable connection to the bed receptacle (under the mattress support).</li> <li>• Verify the bed has power.</li> </ul>
No Fowler up or down motion.	<ul style="list-style-type: none"> <li>• Ensure Total Lockout (padlock LED) is not enabled.</li> <li>• Ensure Fowler Lockout (padlock LED) is not enabled.</li> <li>• Check the siderail control panel cable connection to the bed receptacle (under the mattress support).</li> <li>• Verify the bed has power.</li> </ul>
Fowler does not reach full height.	<p>This situation happens when the CPR handle is used to partly lower the Fowler. The use of the CPR mechanism for this purpose creates a situation where the Fowler motor is out of sync with the actual position of the Fowler. To correct the situation:</p> <ul style="list-style-type: none"> <li>• Completely lower the Fowler using the CPR handle or the Fowler down control to enable the Fowler motor to reset itself.</li> </ul> <p>Refer to the “Emergency CPR Release” section of the Operations Manual for more information.</p>
No Gatch up or down motion.	<ul style="list-style-type: none"> <li>• Ensure Total Lockout (padlock LED) is not enabled.</li> <li>• Ensure Gatch Lockout (padlock LED) is not enabled.</li> <li>• Check the siderail control panel cable connection to the bed receptacle (under the mattress support).</li> <li>• Verify the bed has power.</li> </ul>
No Auto Contour motion.	<ul style="list-style-type: none"> <li>• Ensure Gatch Lockout (padlock LED) is not enabled.</li> <li>• Check the two Auto Contour limit switches. Replace if needed.</li> </ul>

# Troubleshooting Guide

Problem / Failure	Recommended Action
Improper operation of the CPR positioning: Gatch does not lower and/or the Fowler actuator does not reset.	<ul style="list-style-type: none"> <li>Check the two CPR limit switches and replace if needed.</li> </ul>
The weight shown on the scale display (optional) is incorrect.	<ul style="list-style-type: none"> <li>Zero the bed without the patient on it (see the “Scale System Usage” section in the Operations Manual).</li> <li>Calibrate the scale.</li> </ul>
Scale display (optional) reads OVERLOAD or shows incorrect weight.	<ul style="list-style-type: none"> <li>The weight on the bed is greater than the scale capacity.</li> <li>If the patient weight is greater than 500 lb (227 kg), do not use the scale. If that is not the case, zero the scale without the patient on the bed.</li> <li>Ensure equipment added to the bed was done using the “Add/Remove Equipment Procedure” (see the Operations Manual).</li> <li>One or more load cell cables are not properly connected or completely disconnected from the scale control board (see Figure 13, <a href="#">page 1-51</a>) for the connecting positions of the load cell cable).</li> </ul>
Angle shown on the scale display (optional) is incorrect.	<ul style="list-style-type: none"> <li>Calibrate the scale.</li> </ul>
Scale display (optional) does not turn on when the scale is activated.	<ul style="list-style-type: none"> <li>Verify the bed has power.</li> <li>Press ENTER key twice.</li> <li>Unhook the foot board from the bed by lifting it off the connectors. Ensure the connectors match correctly and press the foot board down firmly to reconnect it to the bed.</li> <li>Check the Scale user interface cable connection to the scale control board (J6 connector, see drawing OL280005L, <a href="#">page 1-76</a>).</li> </ul>
Scale (optional) displays “No scale found”.	<ul style="list-style-type: none"> <li>Press the ENTER key. If message persists, contact our Technical Service. If message disappears, data displayed is correct and current procedure can be resumed.</li> </ul>
All lockout LEDs (padlock icons) flash or “ERROR” appears at the bottom of the LCD display.	<ul style="list-style-type: none"> <li>The scale is unable to measure the weight due to an electrical problem. Please contact our Technical Service.</li> </ul>
Bed Exit system (optional) does not operate properly.	<ul style="list-style-type: none"> <li>Verify the load cell cable connections to the scale control board (see Figure 13, <a href="#">page 1-51</a>) for the load cell cable connection positions).</li> </ul>
Bed Exit system (optional) does not work at all or alarm goes off when Bed Exit is activated.	<ul style="list-style-type: none"> <li>Check the load cell cable connections to the scale control board (see Figure 13, <a href="#">page 1-51</a>) for the load cell cable connection positions).</li> </ul>
Nurse call (optional) or Bed Exit (optional) signal does not reach the nurse desk.	<ul style="list-style-type: none"> <li>Verify the bed has power.</li> <li>Check the connecting cable.</li> <li>Check the wall outlet.</li> <li>Contact the Technical Service.</li> </ul>

# Quick Reference Replacement Parts List

## Note

The parts and accessories listed on this page are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call Stryker Customer Service USA: 1-800-327-0770, Canada: 1-888-233-6888 for availability and pricing.


Part Name	Part Number
<b>ELECTRIC/ELECTRONIC COMPONENTS</b>	
Motor Control Board	QDF23-0315
Scale Control Board (Optional)	QDF28-0265
Scale Display (Optional)	QDF28-0266
Nurse Call Control Board (Optional)	QDF20-0126
Nurse Call (Optional)/GEN III (Optional) Control Board	QDF20-0121
Micro Switch	1325P003
Push Button Micro Switch	QDF9159
Power Connector	QDF9574
Fuse Support	QDF9575
Fuse - Fast Acting 10A, 250Vac for 100/120Vac Electric Systems	QDF8078
Fuse - Slow Blow 6.3A, 250Vac for 200/220/240Vac Electric Systems	QDF8068
Power Cord with Straight North American Molded Plug	QDF8066
Power Cord with 90° North American Molded Plug (Optional)	QDF8066-90D
Foot Board/Control Board Cable	QDF28-0259
Night Light	QDF9539
Load Cell with Long Cable	QDF14-1367
Load Cell with Short Cable	QDF25-0218
Hi-Lo Actuator	28-0768
Head Section Actuator	28-0769L
Thigh Section Actuator	28-0770
<b>FOOT BOARD COMPONENTS</b>	
Round Self-Sticking Screw Cover	QDF28-0334

# Quick Reference Replacement Parts List

Part Name	Part Number
<b>SIDERAIL COMPONENTS</b>	
"Lift To Release" Self-Sticking Sticker	QDF28-0146
Oblong Self-Sticking Screw Cover	QDF28-0145
<b>BASE COMPONENTS</b>	
Foot End Casing Plastic Cover	QP28-0111-13
Plastic Base Tube Cover	QP25-0023-13
6" Caster with Locking Mechanism	R25-0388-13
Fifth Steer Wheel	RL5
Right Wheel Cover	QP28-0299-13
Left Wheel Cover	QP28-0300-13
Brake/Steer Pedal	QP28-0130-11
<b>MISCELLANEOUS</b>	
OG2 Grease	M0027
Touch-Up Paint - Opal	PD020
Cable Tie	QDF9518

# Electrical Service Information

## MOTOR CONTROL BOARD


SW1 SWITCH POSITIONS	
BED MODEL AND MAINTENANCE TEST	SW1
FL28C	1 2 3 4
MAINTENANCE TEST	ON OFF ON ON
	ON ON ON ON


## POWER SUPPLY

### 120v Power supply


### 100V-200V-220V-240V POWER SUPPLY

# Electrical Service Information


## COMMUNICATION INTERFACE (OPTIONAL)

# Electrical Service Information

SCALE (OPTIONAL)

English

[Return To Table of Contents](#)

## MAINTENANCE MENU: ERROR CODES

### Scale Error Codes

16 = Angle sensor defective.

17 = Scale chip failure. The scale chip on the scale board is faulty.

For errors 16 and 17, unplug the bed and wait 5 seconds. Then plug the bed back in. If the error persists, replace the scale board and calibrate the bed scale.

18 = Value of the foot left load cell is out of range: load cell may be disconnected, faulty or scale board is faulty.

19 = Value of the head right load cell is out of range: load cell may be disconnected, faulty or scale board is faulty.

20 = Value of the foot right load cell is out of range: load cell may be disconnected, faulty or scale board is faulty.

21 = Value of the head left load cell is out of range: load cell may be disconnected, faulty or scale board is faulty.

For errors 18 to 21, first verify load cell cables are plugged in and inspect for broken wires. If no visible problem is found, swap two load cells on the scale connectors. This will help determine whether the load cell or the scale board is faulty. If the error code changes with the load cell, the problem is with the load cell. Replace the load cell. However, if the error code remains the same, then the problem is with the scale board. Replace scale board. In both cases, the bed has to be calibrated.

## MAINTENANCE PROCEDURES

---

### **WARNING**

- Only service technicians from Stryker or service personnel trained by Stryker should perform the procedures detailed in this maintenance manual. This is especially important for procedures related to the Scale (option) and Bed Exit (option) systems. Failure to follow this restriction can result in serious damage to the bed and/or severe injury to users or patients.
  - To prevent injury when working under the bed with the bed in the high position, always place blocks under the mattress support frame and apply the brakes.
  - Always unplug the bed power cord from the wall outlet when servicing or cleaning the bed.
- 

### **Note**

Throughout this maintenance manual, the words “right” and “left” refer to the right and left sides of a patient lying face up on the bed.

# Service Information

## HEAD END SIDERAIL ASSEMBLY REPLACEMENT

### Tools Required:

- 1/2" Socket
- 3/8" Drive Ratchet
- 3/16" Allen Wrench
- Cutting Pliers


Figure 1

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. If applicable, remove the plastic litter cover (if equipped). Raise the fowler to the highest position and raise the siderail to be repaired.
3. Unplug the power cord from the wall receptacle.
4. Loosen the lock ring (A) and unplug the siderail cable.
5. Using cutting pliers, remove the cable ties holding the siderail cable to the frame.
6. Using a 3/16" Allen wrench and a 1/2" socket, remove the four locknuts/bolts (B) holding the siderail assembly to the head section and remove the assembly. Support the assembly when removing the last bolts.
7. Reverse the above steps to install the new siderail assembly.
8. Test the siderail movement and all controls of both control panels, including the nurse call (optional) and the communications package (optional), for proper operation before returning the bed to service.

## FOOT END SIDERAIL ASSEMBLY REPLACEMENT

### Tools Required:

- 1/2" Socket
- 1/2" Combination Wrench
- 3/8" Drive Ratchet
- Bungee Cord (or equivalent)


Figure 2

### Procedure:

1. Raise the bed to the highest position and apply the brakes. Raise the siderail to be repaired.
2. Raise the Gatch to the highest position and unplug the power cord from the wall receptacle.
3. If applicable, remove the plastic litter cover (optional). Lift and fold the foot section back toward the head end of the bed. Secure it using a bungee cord.
4. Using a 1/2" socket and a 1/2" combination wrench, remove the four locknuts and four bolts (A) holding the siderail assembly to the support and remove the assembly. Support the siderail assembly while removing the last bolts.
5. Reverse the above steps to install the new siderail assembly.
6. Test the siderail for proper operation before returning the bed to service.

# Service Information

## HEAD END SIDERAIL OUTER PANEL MEMBRANE REPLACEMENT

### Tools Required:

- Small Regular Screwdriver
- # 2 Phillips Screwdriver


Figure 3

### Procedure:

1. Raise the bed to the highest position and apply the brakes. Raise the siderail to be repaired.
2. Unplug the power cord from the wall receptacle.
3. Using a small regular screwdriver, lift and remove the following self-adhesive parts: the outer membrane (A), the two screws covers (C) and the "Lift to Release Siderail" Label (B) to expose the screws holding the siderail cover. Use caution when inserting the screwdriver under the labels to avoid damaging the cover.

### Note

Do not reuse the self-adhesive parts due to the coating considerably losing its effectiveness once they are removed. Replace them with new parts.

4. Using a #2 Phillips screwdriver, remove the eight screws (D) holding the cover to the siderail and remove the cover.

### WARNING

The eight screws (D) must be replaced. The Scotch-Grip coating is less effective once they have been tightened and removed thereafter.

5. Disconnect the defective membrane and remove it.
6. Reverse the above steps to install the new membrane.
7. Test all controls of the control panel for proper operation before reassembling the siderail and returning the bed to service.

## SPEAKER/SIDERAIL INNER PANEL MEMBRANE REPLACEMENT

### Tools Required:

- Small Regular Screwdriver
- #2 Phillips Screwdriver
- #3 Phillips Screwdriver
- Cutting Pliers
- Bungee Cord (or equivalent)
- 5/32" Allen Wrench

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 3 ([page 1-30](#)) for this procedure.

1. Follow steps 1 to 4 of the "Siderail Outer Panel Membrane Replacement" procedure ([page 1-30](#)).
2. Secure the siderail to be repair to the adjacent siderail using a bungee cord (see Figure 4).
3. Using cutting pliers, cut the lower cable tie (E) holding the bottom part of the grey cable (there may be two cables) to the aluminum structure. The slack in the cable will be necessary for the following steps.


Figure 4

### Note

Make sure the cable ties are inserted in the aluminum structure before reinstalling the structure.

4. Using a #2 Phillips screwdriver, remove the screw (F) holding each siderail arm cover and lift up the covers to remove them.
5. Using a 5/16" Allen wrench, remove the two Allen screws (G) holding each pommel to the siderail arms. Leave one screw loosened on each pommel until ready to remove both pommels.
6. While supporting the rail, remove the two last screws and the bungee cord.
7. Lift the siderail, pass it over the mechanism arms and lay it on the mattress support while gently pulling on the cables.
8. Using a #3 Phillips screwdriver, remove the six screws (H) holding the aluminum structure to the siderail.
9. Grasp and lift the upper part of the aluminum structure until it may be secured temporarily to one of the mechanism arms using a cable tie.

**To replace the speakerphone, proceed with steps 10 to 13.**

**To replace the inner panel membrane, proceed to steps 14 to 17.**

10. Remove the cables connected to the speakerphone. Note their position before removing them.
11. Using a #2 Phillips screwdriver, remove the four screws (I) holding the speakerphone to the siderail and remove it.
12. Reverse the preceding steps to install the new speakerphone.

### Note

Before reassembling the siderail, it is recommended the continuity is checked at the speakerphone terminals to ensure the circuit operates properly.

13. Test the nurse call and all controls of the inner and outer control panels for proper operation before reassembling the siderail and returning the bed to service.
14. Disconnect the inner membrane connector.

# Service Information

---

## **SPEAKER/SIDERAIL INNER PANEL MEMBRANE REPLACEMENT (CONTINUED)**

15. Using a small regular screwdriver, lift and remove the inner membrane. Use caution when inserting the screwdriver under the membrane to avoid scratching the plastic cover.
16. Connect and install the new membrane.
17. Test the siderail inner and outer controls for proper operation before reassembling the siderail and returning the bed to service.

## FOOTBOARD CONTROL MEMBRANE REPLACEMENT

### Tools Required:

- Small Regular Screwdriver
- #2 Phillips Screwdriver


Figure 5

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the power cord from the wall receptacle.
3. Using a small regular screwdriver, lift and remove the five round self-adhesive screw-covers (A) located on the foot board cover and the self-adhesive membrane (B) of the foot board control panel to expose the screws holding the cover to the foot board. Use caution when inserting the screwdriver under the membrane to avoid scratching the plastic cover.

### Note

Do not reuse the self-sticking parts removed since their self-adhesive coating considerably loses its efficiency once they are removed. Replace them with new parts.

4. Properly ground yourself (see the “Static Protection Procedure”, [page 1-13](#)).
5. Using a #2 Phillips screwdriver, remove the seven screws (C) holding the cover to the foot board and remove the cover after having removed the cables connected to it. Note the cable positions in order to replace them properly.
6. Remove the defective membrane.
7. Reverse the above steps to install the new membrane.

### Note

Use caution to avoid adhering the cable to the back of the membrane.

8. Test all the functions of the foot end control panel before returning the bed to service.

## SCALE SYSTEM LCD DISPLAY REPLACEMENT (OPTIONAL)

### Tools Required:

- Small Regular Screwdriver
- #2 Phillips Screwdriver

### Note

Unless otherwise indicated, refer to Figure 5 ([page 1-33](#)) for this procedure.

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the power cord from the wall receptacle.
3. Using a small regular screwdriver, lift and remove the five round self-adhesive screw-covers (B) located on the foot board cover as well as the self-sticking membrane (A) of the foot board control panel to expose the screws holding the cover to the foot board. Use caution when inserting the screwdriver under the membrane to avoid scratching the plastic cover.

### Note

Do not reuse the self-adhesive parts due to their adhesive coating considerably losing its effectiveness once they are removed. Replace them with new parts.

4. Properly ground yourself (see the “Static Protection Procedure”, [page 1-13](#)).
5. Using a #2 Phillips screwdriver, remove the seven screws (C) holding the cover to the foot board and remove the cover after having removed the cables connected to it. Note the cable positions in order to replace them properly.
6. Using a small regular screwdriver, push aside the clips (D) and remove the LCD display (E) after having disconnected its cable.
7. Reverse the above steps to install the new LCD display.
8. Test the LCD display for proper operation before replacing the foot board cover and returning the bed to service.

## FOOT SECTION REPLACEMENT

### Tools Required:

- 1/2" Combination Wrench
- #2 Phillips Screwdriver
- Bungee Cord (or equivalent)


Figure 6

### Procedure:

1. Raise the bed to the highest position and apply the brakes. Flatten all sections of the mattress support.
2. Unplug the power cord from the wall receptacle.
3. Lower the siderails.
4. If applicable, remove the foot section plastic cover (optional). Lift and fold back the foot section toward the head end of the bed and secure it to the bed using a bungee cord (or equivalent).
5. Remove the foot end mattress retainer (A) and the prop rod and install it on the new foot section. Replace the foot section to horizontal position.
6. Using a 1/2" combination wrench, remove the two locknuts, two washers, two shoulder spacers, two bolts (B) linking the foot section to the thigh section and remove the defective section.
7. Reverse the above steps to install the new foot section.
8. Test the bed for proper operation before returning it to service.

## THIGH SECTION REPLACEMENT

### Tools Required:

- Needle Nose Pliers
- 1/2" Combination Wrench
- #2 Phillips Screwdriver

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 6 ([page 1-35](#)) for this procedure.

1. Raise the bed to the highest position and apply the brakes. Flatten all sections of the mattress support.
2. Unplug the power cord from the wall receptacle.
3. Raise the siderails.
4. Using needle nose pliers, remove the rue ring cotter/washers/clevis pin (E, Figure 8, [page 1-40](#)) linking the thigh section lever arms to the thigh actuator tube.
5. Using a 1/2" combination wrench, remove the four locknuts, four washers, four shoulder spacers, and four bolts (B, C) linking the thigh section to the foot and seat sections.
6. Remove the defective thigh section and lay it upside down on a workbench.
7. Using a #2 Phillips screwdriver, remove the two screws (D) holding the micro switch activator to the thigh section. Install the micro switch activator on the new thigh section.
8. Reverse the above steps to install the new thigh section.

#### Note

Apply grease on the clevis pin and nylon washers before joining the foot section to the actuator tube.

9. Test the Gatch as well as the Auto Contour positioning (optional) before returning the bed to service.

#### Note

Check the positioning of the switch activator on the thigh section if Auto Contour (option) does not operate properly.

## SEAT SECTION REPLACEMENT

### Tools Required:

- 1/2" Long Socket
- (2) 1/2" Combination Wrenches
- 3/8" Drive Ratchet
- #2 Phillips Screwdriver

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 6 ([page 1-35](#)) for this procedure.

1. Raise the bed to the highest position and apply the brakes. Flatten the mattress support.
2. Unplug the power cord from the wall receptacle.
3. Lower the siderails.
4. Using two 1/2" combination wrenches, remove the locknut, shoulder spacers and bolt holding the mattress support lever of the CPR mechanism (N, Figure 17, [page 1-59](#)) to the seat section.
5. Using a 1/2" long socket and a 1/2" combination wrench, remove the four locknuts, four washers and four bolts (E) holding the seat section to the frame.
6. Using a 1/2" combination wrench, remove the four locknuts, four washers, four shoulder spacers and four bolts (C, F) linking the seat section to the thigh and head section and remove the defective seat section.
7. Using a #2 Phillips screwdriver, remove the three screws (G) holding the protective plate to the seat section.
8. Reverse the above steps to install the new seat section.
9. Test the bed for proper operation before returning the bed to service.

# Service Information

## FWLER REPLACEMENT

### Tools Required:

- 1/2" Socket
- (2) 1/2" Combination Wrenches
- #2 Phillips Screwdriver
- 3/16" Allen Wrench
- 7/16" Combination Wrench
- Cutting Pliers
- Strap (or equivalent)
- 3/8" Drive Ratchet

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 6, [page 1-35](#), for this procedure.

---

#### WARNING

Do not activate the CPR emergency release handles during this procedure or serious injury to the maintenance technician or damage to equipment could occur.

---

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the power cord from the wall receptacle.
3. Raise the head siderails.
4. If applicable, remove the litter plastic cover (if equipped) or raise the fowler to the highest position. Secure the head section to the bed using a strap to prevent the head section from moving during the procedure.
5. Loosen the lock rings (A, Figure 1, [page 1-28](#)) of both siderail cables and remove the siderail cables. Using cutting pliers, cut the cable ties holding the siderail cables to the frame.
6. Using a 3/16" Allen wrench and a 1/2" socket, remove the four locknuts and four bolts (B, Figure 1, [page 1-28](#)) holding each head siderail assembly to the head section and remove the siderails.
7. Using a #2 Phillips screwdriver, remove the screw holding each cable tie present on both sides of the fowler structure.
8. Using a 1/2" combination wrench, remove the two nuts, two shoulder spacers and two bolts (H) holding the fowler to the stabilizer arms. Lay the stabilizer arms down.

**If the bed is equipped with the CPR mechanism, proceed with steps 9 to 16.**

**If the bed is not equipped with the CPR mechanism, proceed to steps 17 to 21.**

9. Remove the two wires (C, Figure 17, [page 1-59](#)) from the micro switch. Note their position before removing them.
10. Using two 1/2" combination wrenches, remove the locknut, shoulder spacer and bolt (L, Figure 17, [page 1-59](#)) holding each CPR handle to the head section. To remove the cable ends from the handle slots, loosen the two locknuts holding each CPR cable to its handle using a 7/16" combination wrench.

#### Note

Insert the cable ends inside the handle slots before assembling the handle to the new head section. Then adjust the two locknuts of each cable in order to have a 1/8" play in the CPR handle movement.

## FOWLER REPLACEMENT (CONTINUED)

11. Using two 1/2" combination wrenches, remove the locknut and bolt (M, Figure 17, [page 1-59](#)) holding the pneumatic cylinder end fitting to the head section lever arms. Lay the CPR mechanism down on the frame.
12. Grasp the head section firmly. Remove the strap and lower the head section.
13. Using a 1/2" combination wrench, remove the two locknuts, two shoulder spacers and two bolts (F) holding the head section to the seat section.
14. Remove the micro switch from the defective head section by pressing its two side clips.
15. Reverse the above steps to install the new head section.
16. Test the head section and CPR mechanism for proper operation before returning the bed to service.


Figure 7

17. Using two 1/2" combination wrenches, remove the locknut and bolt (A, Figure 7) holding the lever to the head section. Lay the lever down on the frame. Keep the spacer (B, Figure 7).
18. Grasp the head section firmly. Remove the strap and lower the head section.
19. Using a 1/2" combination wrench, remove the two locknuts, two shoulder spacers and two bolts (F) holding the head section to the seat section.
20. Reverse the above steps to install the new head section.
21. Test the head section for proper operation before returning the bed to service.

# Service Information

## FOWLER MOTOR REPLACEMENT

### Tools Required:

- Needle Nose Pliers
- #2 Phillips Screwdrivers
- #3 Phillips Screwdrivers
- Pliers
- Cutting Pliers
- 5/16" Socket
- 3/8" Drive Ratchet
- Bungee Cord (or equivalent)
- Small Regular Screwdriver


Figure 8

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Lower the Fowler to the flat (0°) position.
3. Raise the Gatch to the highest position.
4. Unplug the power cord from the wall receptacle.
5. Raise the siderails.
6. Remove the headboard and footboard.
7. Lift and fold the foot section back toward the head end of the bed. Secure it to the bed using a bungee cord.
8. Using needle nose pliers, remove the rue ring cotter, washers, clevis pin (A) linking the actuator tube to the seat section lever arms.

## FWLER MOTOR REPLACEMENT (CONTINUED)

9. Using a #3 Phillips screwdriver, remove the six screws (A, Figure 11, [page 1-47](#)) holding the plastic cover to the foot end casing. Grasp both ends of the cover, move them apart to disengage the inside clips from the round accessory brackets and remove the cover.
10. Using a #2 Phillips screwdriver, remove the four screws (B, Figure 11, [page 1-47](#)) holding the protective plate to the foot end casing.
11. Properly ground yourself (see the “Static Protection Procedure”, [page 1-13](#)).
12. Using cutting pliers, cut the cable tie holding the head actuator cable to the other cables.
13. Remove the actuator cable (connector J8) from the control board.
14. Using pliers, squeeze the upper part of the strain relief bushing (B) and lift it up to remove it from its location.
15. Remove the cable from the bushing and pass the cable connector through the rear casing hole.
16. Using a 5/16” socket, remove the two screws (C) holding the retaining plate to the actuator support.
17. Remove the two pivot pins (D) holding the actuator to the support. To ease the removal of the pins, insert a small regular screwdriver into the opening at the end of the actuator and push out the pins.
18. Move the actuator toward the center of the bed to remove it from its location.
19. Reverse the above steps to install the new actuator.

---

### CAUTION

The limits of the new fowler motor must be adjusted before reconnecting the tube to the head section lever arms. An improper adjustment can damage the head section structure.

---

20. To adjust the limits of the replacement head actuator, proceed as follows:
  - A. Ensure the motor cable is connected to the control board and plug the power cord into the wall receptacle.
  - B. Position the new motor tube holes horizontally. Grasp the tube firmly to prevent it from rotating and press the “Fowler up” control for a few seconds. Then press the down control until the motor stops. This will be the lower limit of the motor.
  - C. Gently turn the tube in either direction to align the tube holes with those of the head section lever arms. Then raise the Fowler a few inches while firmly holding the tube and lower it completely.
  - D. Check the alignment of the holes. If the tube holes are not aligned with those of the lever arms, repeat steps B and C until they are. Once the holes are aligned, install the washers and the clevis pin. Slightly raise the head section manually to ease the insertion of the clevis pin.
  - E. Before installing the rue ring cotter, raise and lower the Fowler completely to ensure the Fowler section reaches the frame.
21. Test the Fowler section for proper operation before returning the bed to service.

## GATCH MOTOR REPLACEMENT

### Tools Required:

- Needle Nose Pliers
- #2 Phillips Screwdrivers
- #3 Phillips Screwdrivers
- Pliers
- Cutting Pliers
- Bungee Cord (or equivalent)
- 5/16" Socket
- 3/8" Drive Ratchet
- Small Regular Screwdriver

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 8, [page 1-40](#), for reference for this procedure.

1. Raise the bed to the highest position and apply the brakes. Flatten the litter.
2. Unplug the power cord from the wall receptacle.
3. Raise the siderails.
4. Remove the litter plastic covers (if equipped).
5. Using needle nose pliers, remove the rue ring cotter, washers and clevis pin (E) linking the motor tube to the thigh section lever arms.
6. Using a #3 Phillips screwdriver, remove the six screws (A, Figure 11, [page 1-47](#)) holding the plastic cover to the foot end casing. Grasp both ends of the cover, move them apart to disengage the inside clips from the round accessory brackets and remove the cover.
7. Using a #2 Phillips screwdriver, remove the four screws (B, Figure 11, [page 1-47](#)) holding the protective plate to the foot end casing.
8. Properly ground yourself (see the "Static Protection Procedure" section, [page 1-13](#)).
9. Using cutting pliers, cut the cable tie holding the thigh actuator cable to the other cables.
10. Remove the motor cable (connector J9) from the control board.
11. Using pliers, squeeze the upper part of the strain-relief bushing (F) and lift it up to remove it from its location.
12. Remove the cable from the bushing and pass the cable connector through the casing hole.
13. Using a 5/16" socket, remove the two screws (G) holding the retaining plate to the actuator support.
14. Remove the two pivot pins (H) holding the motor to the support. To ease the removal of the pins, insert a small regular screwdriver into the opening at the end of the motor and push out the pins.
15. Move the actuator toward the center of the bed to remove it from its location.
16. Reverse the above steps to install the new motor.

---

### CAUTION

The limits of the new gatch motor must be adjusted before reconnecting the tube to the thigh section lever arms. An improper adjustment can damage the thigh section structure.

---

## GATCH MOTOR REPLACEMENT (CONTINUED)

17. To adjust the limits of the replacement thigh motor, proceed as follows:
  - A. Ensure the motor cable is connected to the control board and plug the power cord into the wall receptacle.
  - B. Position the new actuator tube holes horizontally. Grasp the tube firmly to prevent it from rotating, and press the Gatch up control for a few seconds. Then press the down control until the motor stops. This will be the lower limit of the motor.
  - C. Gently turn the tube in either direction to align the tube holes with those of the thigh section lever arms. Then raise the Gatch a few inches while firmly holding the tube and lower it completely.
  - D. Check the alignment of the holes. If the tube holes are not aligned with those of the lever arms, repeat steps B and C until they are. Once the holes are aligned, install the washers and the clevis pin. Slightly raise the thigh section manually to ease the insertion of the clevis pin.
  - E. Before installing the rue ring cotter, raise and lower the Gatch completely to ensure the motor stops when the thigh section reaches the frame.
  
18. Test the Gatch section for proper operation before returning the bed to service.

# Service Information

## HEAD END LIFT MOTOR REPLACEMENT

### Note

In order to maintain the adjustment of the low position when replacing a lift motor, a special tool kit must be used. The kit includes alignment fixtures. To obtain this kit, contact our Customer Service department and order part number KR0154.

### Tools Required:

- Tool Kit KR0154
- Cutting Pliers
- #2 Phillips Screwdriver
- Small Regular Screwdriver
- 1/2" Combination Wrench

### Procedure:

### Note

When raising or lowering any of the litter sections, ensure that all siderails are in the highest position. Unless otherwise indicated, refer to Figure 10 for reference for this procedure.


Figure 9

1. Install the alignment fixtures (KR0154) on both sides of the brake mechanism on the base frame, under the lower fixed frame.
2. Lower the bed, until the lower fixed frame is completely supported by the U-shaped pieces of the fixture (see Figure 9).
3. Remove the headboard.
4. Raise the head section to about 45°.
5. Unplug the power cord from the wall receptacle.
6. Using cutting pliers, cut the cable ties holding the defective lift motor power cord to the bed frame. Then unplug the motor connector.
7. Using a 1/2" combination wrench, remove the two bolts (D), two washers (C) and two shoulder spacers (B) holding the molded nut support to the elevation lever arm.
8. Using a #2 Phillips screwdriver, remove the two screws (E) holding the actuator support cover (G) to the head lever arm and remove it.
9. Remove the two pivot pins (H) holding the motor (A) to the support. To ease the removal of the pins, insert a small regular screwdriver into the opening at the end of the motor and push out the pins.
10. Using a #2 Phillips screwdriver, remove the two screws (F) holding the plastic screw cover (J) to the motor. Remove the plastic screw cover and install it on the new actuator.
11. Reverse the above steps to install the new actuator.


Figure 10

### CAUTION

The limit of the new motor must be adjusted prior to reconnecting it to the lift lever. An improper adjustment can damage the lift mechanism.

## HEAD END LIFT MOTOR REPLACEMENT (CONTINUED)

12. To adjust the limits of the new motor, proceed as follows:
  - A. Ensure the motor cable is connected to the control board and plug the power cord into the wall receptacle.
  - B. Press the bed down control until the motor stops. This will be the lower limit of the motor.
  - C. Attach the molded nut support on the motor to the lift lever arms.
  - D. With the alignment fixtures still in place, raise and lower the bed completely to ensure the lower limit is properly adjusted.
13. Test the Bed Up and Bed Down functions for proper operation before returning the bed to service.

## FOOT END LIFT MOTOR REPLACEMENT

### Note

In order to maintain the adjustment of the low position when replacing a lift motor, a special tool kit must be used. The kit includes alignment fixtures. To obtain this kit, contact our Customer Service department and order part number KR0154.

### Tools Required:

- Tool Kit KR0154
- Cutting Pliers
- #2 Phillips Screwdriver
- Small Regular Screwdriver
- 1/2" Combination Wrench
- Bungee Cord (or equivalent)

### Note

When raising or lowering any of the litter sections, ensure that all siderails are in the highest position. Unless otherwise indicated, refer to Figure 10, [page 1-44](#), for reference for this procedure.

1. Install the alignment fixtures (KR0154) on both sides of the brake mechanism on the base frame, under the lower fixed frame.
2. Lower the bed, until the lower fixed frame is completely supported by the U-shaped pieces of the fixture (see Figure 9, [page 1-44](#)).
3. Remove the mattress and the plastic litter cover (if equipped) from the foot section metal frame.
4. Raise the fowler section to about 45°, and raise the thigh section to the highest position.
5. Fold the foot section back toward the head end of the bed. Secure it using a bungee cord.
6. Unplug the power cord from the wall receptacle.
7. Using a #2 Phillips screwdriver, remove the four screws holding the frame cover plate.
8. Using cutting pliers, cut the cable ties holding the defective lift motor power cord to the bed frame. Unplug the motor connector.
9. Using a 1/2" combination wrench, remove the two bolts (D), two washers (C) and two spacers (B) holding the molded nut support to the elevation lever arm.
10. Using a #2 Phillips screwdriver, remove the two screws holding the motor support cover to the lift lever arm and remove it.
11. Remove the two pivot pins (H) holding the motor (A) to the support. To ease the removal of the pins, insert a small regular screwdriver into the opening at the end of the motor and push out the pins.
12. Remove the defective motor by passing it under the bed frame.
13. Using a #2 Phillips screwdriver, remove the two screws (F) holding the plastic screw cover (J) to the motor. Remove the plastic screw cover and install it on the new actuator.

[Return To Table of Contents](#)

## FOOT END LIFT MOTOR REPLACEMENT (CONTINUED)

14. Reverse the above steps to install the new motor.
15. Plug the bed into the wall receptacle.
16. Lower the bed, and hold the bed down button until the lift motor reaches its lowest position.

### CAUTION

The limits of the new motor must be adjusted prior to reconnecting it to the lift lever. An improper adjustment can damage the lift mechanism.

---

17. To adjust the limits of the new motor, proceed as follows:
  - A. Ensure the motor cable is connected to the control board and plug the power cord into the wall receptacle.
  - B. Press the bed down control until the motor stops. This will be the lower limit of the motor.
  - C. Attach the molded nut support on the motor to the lift lever arms.
  - D. With the alignment fixtures still in place, raise and lower the bed completely to ensure the lower limit is properly adjusted.
18. Test the Bed Up and Bed Down functions for proper operation before returning the bed to service.

## CONTROL BOARD REPLACEMENT

### Tools Required:

- #2 Phillips Screwdriver
- #3 Phillips Screwdriver
- Needle Nose Pliers
- Cutting Pliers
- 3/8" Combination Wrench


Figure 11

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the bed from the wall receptacle.
3. Remove the foot board.
4. Using a #3 Phillips screwdriver, remove the six screws (A) holding the plastic cover to the foot end casing. Grasp both ends of the cover, move them apart to disengage the inside clips from the round accessory brackets and remove the cover.
5. Using a #2 Phillips screwdriver, remove the four screws (B) holding the protective plate to the foot end casing.
6. Properly ground yourself (see "Grounding Diagram" [page 1-13](#)).
7. Using cutting pliers, cut the cable ties holding together the cables.

### Note

- Carefully note the position and the grouping of the cables before clipping the cable ties.
  - Ensure Dip Switch settings of the new control board match the ones on the old board.
8. Remove the jumper from W1 on the control board and install it on the new board.
  9. Remove all cables connected to the control board. Note their location so they will be connected properly to the new board. Refer to drawing OL280120, [page 1-88](#), or OL280121, [page 1-94](#), for the connecting position of the cables on the motor control board.
  10. Using a 3/8" combination wrench, remove the nut/screw (C) holding the ground cable to the frame.
  11. Using needle nose pliers, pinch the upper part of the stand-off pins and lift the defective board up and out.
  12. Using a 3/8" wrench and a #2 Phillips screwdriver, remove the nut/screw (D) holding the ground cables to the board.
  13. Reverse the above steps to install the new control board.
  14. Test all the bed functions before returning the bed to service.

[Return To Table of Contents](#)

# Service Information

## HEADWALL CONTROL BOARD REPLACEMENT (OPTIONAL)

### Tools Required:

- #2 Phillips Screwdriver



Figure 12

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the bed power cord from the wall outlet.
3. Remove the power cord from the power connector and disengage it from the wire clip.
4. Properly ground yourself (see “Grounding Diagram” [page 1-13](#)).
5. Using a #2 Phillips screwdriver, remove the six screws (A) holding the cover to the litter. Keep the cable clips.
6. Once the screws removed, pivot the cover and lay it flat on the bottom of the litter (see opposite illustration).
7. Remove all cables connected to the control board. Note their location so they will be connected properly to the new board. Refer to drawing OL280120, [page 1-88](#), or OL280121, [page 1-94](#), for the connecting position of the cables on the power connector.
8. Turn the cover to a vertical position and. Using a #2 Phillips screwdriver, remove the four screws (B) holding the control board (C) to the head casing cover and remove the board.
9. Reverse the above steps to install the new control board. Before installing the new board, ensure the dip switch settings on the new board match the ones on the old board.
10. Test the nurse call for proper operation before returning the bed to service.

## POWER CORD CONNECTOR REPLACEMENT

### Tools Required:

- #1 Phillips Screwdriver
- #2 Phillips Screwdriver

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 12, [page 1-48](#), for reference for this procedure.

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the bed power cord from the wall outlet.
3. Remove the power cord from the power connector and disengage it from the wire clip.
4. Properly ground yourself (see “Grounding Diagram” [page 1-13](#)).
5. Using a #2 Phillips screwdriver, remove the six screws (A) holding the cover to the litter. Keep the cable clips.
6. Once the screws are removed, pivot the cover and lay it flat on the bottom of the litter (see opposite illustration).
7. Remove all cables connected to the power cord connector. Note the location of the cables in order to properly reconnect them later. Refer to drawing OL280120, [page 1-88](#), or OL280121, [page 1-94](#), for the connecting position of the cables on the power connector.
8. Turn the cover to a vertical position. Using a #1 Phillips screwdriver, remove the two screws (D) holding the power connector to the cover.
9. Unplug the connector and remove it.

### Note

Ensure the proper size fuses are used in the new power cord connector before installing it.

10. Reverse the above steps to install the new power cord connector.
11. Test all the bed functions before returning the bed to service.

## POWER CORD CONNECTOR FUSE REPLACEMENT

### Tools Required:

- Small Regular Screwdriver

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 12, [page 1-48](#), for reference for this procedure.

1. Remove the power cord from the power cord connector.
2. Using a small regular screwdriver, remove the fuse holder.
3. Remove the defective fuse and replace it with a new one.

### Note

120 VAC rated FL28C beds use fast acting, 250 Vac, 10A type fuses. For beds rated at voltages other than 120 Vac, see "Specifications", [page 1-6](#), for fuse ratings.

4. Replace the fuse holder.

### Note

Make sure to use the right fuse holder for the fuses used with the bed and make sure the fuses are correctly installed to be stable in the fuse holder.

5. Test the bed for proper operation before returning the bed to service.

## SCALE/BED EXIT BOARD/BEEPER REPLACEMENT (OPTIONAL)

### Tools Required:

- #2 Phillips Screwdriver
- 3/8" Combination Wrench
- Needle Nose Pliers


Figure 13

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. If applicable, remove the plastic litter cover (if equipped) and raise the fowler to the highest position.
3. Unplug the power cord from the wall receptacle.
4. Using a #2 Phillips screwdriver, remove the two screws (A) holding the cover plate.
5. Properly ground yourself (see “Grounding Diagram” section).  
**To replace the Bed Exit Beeper, proceed with steps 6 to 9.**  
**To replace the scale/bed exit board, proceed with steps 10 to 14**
6. Disconnect the beeper cable from J3 of the control board.
7. Lift the beeper (B) to remove it from the plate. Remove any glue residue.
8. Reverse the above steps to install the new beeper.
9. Test the Bed Exit system to make sure the beeper operates properly.
10. Remove all cables connected to the control board. Note the cable locations in order to reconnect them correctly. Refer to drawing OL280005L, [page 1-76](#)) for the connecting position of the cables on the scale/bed exit board.
11. Using a 3/8" combination wrench, remove the nut holding the ground cable to the scale board (B).
12. Using needle nose pliers, pinch the upper part of the stand-off pins and lift the defective board up and out.
13. Reverse the above steps to install the new scale control board.

### Note

Be sure to reconnect the load cell cables on the scale/bed exit board in the proper order. Refer to Figure 13 when reconnecting the load cell cables.

14. Calibrate the scale before returning the bed to service. Refer to the “Scale Calibration” procedure, [page 1-55](#).

# Service Information

## TRANSFORMER REPLACEMENT (100/200/220/240 VAC OPTIONS)

### Tools Required:

- #2 Phillips Screwdriver
- 1/2" Combination Wrench


Figure 14

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Unplug the power cord from the wall receptacle.
3. Using a #2 Phillips screwdriver, remove the six screws (A) holding the cover to the litter. Keep the cable clips. See Figure 12, [page 1-48](#).
4. Properly ground yourself (see "Grounding Diagram" section, [page 1-13](#)).
5. Using a #2 Phillips screwdriver, remove the nut, washer and bolt (A) holding the transformer to the litter. Disconnect the transformer cable and remove it.

### Note

Use caution to avoid overtightening the bolt.

6. Reverse the above steps to install the new transformer.
7. Test the bed for proper operation before returning the bed to service.

## LOAD CELL REPLACEMENT (OPTIONAL)

### Tools Required:

- Needle Nose Pliers
- Diagonal Pliers
- (2) Jack Stands (or equivalent)
- 1/4" Allen Wrench
- 1/2" Combination Wrench


Figure 15

### Procedure:

1. Plug the power cord in to the wall receptacle.
2. Raise the bed to the highest position and apply the brakes.
3. Depending on the location of the load cell to replace, proceed with the following operation:

**Head end of the bed:** Remove the head board. Remove the plastic litter cover (if equipped) of the fowler or raise the fowler to the highest position.



**Foot end of the bed:** Remove the foot board. Remove the plastic litter cover (if equipped) of the foot section or raise the gatch to the highest position.

4. Place a jack stand under each side of the frame at the head or foot end of the bed, depending on the location of the load cell to be replace (see following illustrations).
5. Using needle nose pliers, remove the two rue ring cotters, two washers and two clevis pins (A) holding the defective load cell and the one adjacent to it. It is important that the end of both load cells be freed and their cables loosened.

[Return To Table of Contents](#)

# Service Information

## LOAD CELL REPLACEMENT (OPTIONAL) (CONTINUED)

6. Using diagonal pliers, cut the two cable ties holding the defective cell cable to the frame immediately after the cell. Repeat for the adjacent cell.
7. Lower the bed until the two Allen screws holding the defective load cells to the frame are easily accessible.
8. Unplug the power cord from the wall receptacle.
9. Using a 1/4" Allen wrench and a 1/2" combination wrench, remove the two nuts and two Allen screws (B) holding the defective load cell to the frame.
10. Cut the cable ties holding the defective load cell cable the frame. Unplug the defective load cell cable from the scale/bed exit board. Remove the defective load cell.
11. Inspect the elastomer sleeves (C) to ensure that the inner circular hole is not worn. Replace if necessary. Note the orientation of the elastomer sleeves when replacing (see detail at right).
12. Reverse the above steps to install the new load cell.
13. Calibrate the scale before returning the bed to service. Refer to the "Scale Calibration" procedure, [page 1-55](#).

**Proper elastomer sleeve orientation.**

## SCALE CALIBRATION (OPTIONAL)

### VERIFYING THE SCALE ACCURACY

1. Remove all accessories, linens, mattress, etc. from the bed and zero the scale by:
  - Press  to activate the scale. The screen will display:

Weight	Angle
XXX.X KG	+/-XX.X°

- Press and hold ZERO. Display will read:  
HOLD TO ZERO WT., followed by:  
RELEASE TO ZERO (release ZERO), followed by:  
DO NOT TOUCH BED (do not touch the bed)  
Once the zeroing process is completed, the system will return to the Scale mode. The display will read zero for the weight and the current angle value.
2. Place a calibrated weight on the center of the mattress support. The weight displayed should be within the following range:  $\pm 2\%$  for a weight  $\geq 100$  lb (45.4 kg) or  $\pm 2$  lb (0.9 kg) if the weight  $< 100$  lb (45.4 kg).
  3. If the weight displayed is incorrect, remove the weight from the bed and calibrate the scale.

#### Note

If the scale calibration procedure cannot be completed once it has been started, turn the main switch off and turn it back on to restore the previous calibration data.

The scale calibration is done in three steps. Steps 2 and 3 begin automatically as soon as the preceding step is over. Once the calibration process is over, the system returns to the main maintenance menu.

In the maintenance menus, the ZERO key returns to the previous level of the maintenance menus.

The system will automatically shut off after one minute of idle time when in the main maintenance menu. Should this occur, the main maintenance menu will have to be accessed again.

#### Tools Required:

- Digital Inclinometer

#### Procedure:

1. Power off the bed using the main switch at the head end of the bed.
2. Power on the bed while simultaneously holding the Up and Down buttons on the scale system control panel. This allows access to the main maintenance menu. Display will read:  
▶ DEBOG MENU  
CALIB. MENU ▼

#### Note

The main maintenance menu can only be accessed for five seconds after the bed is powered on. The system will shut off after this delay. Repeat the process should the system time out.

# Service Information

---

## SCALE CALIBRATION (OPTIONAL) (CONTINUED)

3. Select CALIB. MENU using the Down key and press ENTER. The display will read:  
PUT BED AT 0°  
PRESS ENTER
4. Put the inclinometer on the seat section of the mattress support and position the bed to 0°. Then press ENTER. Display will read:  
DO NOT TOUCH BED
5. Do not touch the bed until the end of this step. When completed, the display will read:  
PUT BED AT +12°  
PRESS ENTER
6. Put the bed in the Trendelenburg position (head down, foot up) until the inclinometer reads +12°. Press ENTER. The display will read:  
DO NOT TOUCH BED
7. Do not touch the bed until the end of this step. When completed, the display will read:  
PUT BED AT -12°  
PRESS ENTER
8. Put the bed in the Reverse Trendelenburg position (head up, foot down) until the inclinometer reads -12°. Press ENTER. The display will read:  
DO NOT TOUCH BED
9. Once the last step is completed, the system will automatically return to the main maintenance menu.
10. To return to the normal scale display, press ZERO.
11. Before returning the bed to service, test the scale accuracy. See the “Verifying Scale Accuracy”, [page 1-55](#).

### Note

- It is possible to end the calibration process any time before it is completed by pressing the ZERO key. The preceding calibration data will be restored.
- If the bed is in the incorrect angle in step two or three, the same message will reappear on the display until the bed is placed at the angle requested on the display.

## AUTO CONTOUR MICRO SWITCH REPLACEMENT

### Tools Required:

- Bungee Cord (or equivalent)
- #1 Phillips Screwdriver


Figure 16

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Raise the thigh section. Remove the plastic litter cover (if equipped) of the fowler section. Lift and fold the foot section back toward the head end of the bed. Secure it using a bungee cord.

### 1ST AUTO CONTOUR MICRO SWITCH

This switch tells the control board that the Gatch has reached the preset angle for the Auto Contour position.

3. Remove the cables from the switch. **Note the cable locations in order to properly reconnect them to the new switch.**
4. Using the Gatch down control, slowly lower the thigh section until only light contact remains between the activator (A) attached to the knee section and the push-button switch (B). The replacement procedure will start from this position (see opposite illustration).
5. Unplug the bed power cord from the wall outlet.
6. Using a #1 Phillips screwdriver, remove the two screws (C) holding the support to the frame.
7. Press the two clips on each side of the switch to remove it from the support.
8. Reverse the above steps to install the new micro switch.
9. Test the Auto Contour positioning for proper operation before returning the bed to service.

## AUTO CONTOUR MICRO SWITCH REPLACEMENT (CONTINUED)

### 2ND AUTO CONTOUR MICRO SWITCH

This switch tells the control board to start lowering the Gatch when the Fowler is lowered while the bed is in the Auto Contour position.

1. Unplug the bed power cord from the wall receptacle.
2. Using a #1 Phillips screwdriver, remove the two screws (D) holding the support to the litter and remove the support.
3. Remove the cables from the switch. **Note the cable locations in order to properly reconnect them to the new switch.**
4. Using a #1 Phillips screwdriver, remove the two screws (E) holding the switch to the support.
5. Reverse the above steps to install the new micro switch.
6. Test the Auto Contour positioning for proper operation before returning the bed to service.

## CPR MICRO SWITCH REPLACEMENT (OPTIONAL)

### Tools Required:

- #1 Phillips Screwdriver
- #2 Phillips Screwdriver


Figure 17

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Raise the head siderails.
3. Remove the plastic litter cover (if equipped) of the fowler or raise the fowler to the highest position.
4. Unplug the power cord from the wall receptacle.

### FOWLER LEVER MICRO SWITCH

This switch informs on the state of the Fowler actuator (engaged or disengaged) when the CPR emergency release is activated.

5. Using a #1 Phillips screwdriver, remove the two screws (A) holding the micro switch to the fowler lever. Keep the mounting plate (B). Remove the defective micro switch.
6. Remove the cables from the micro switch. **Note the cable locations in order to properly reconnect them to the new switch.**
7. Reverse the above steps to install the new micro switch.
8. Test the CPR positioning before returning the bed to service.

### FOWLER STRUCTURE MICRO SWITCH

This switch signals to the motor control board on the one hand that the resetting of the Fowler actuator may begin because the Fowler is completely lowered, and on the other hand that the lowering to flat of the Gatch may also begin.

9. Remove the cable wires (C) from the micro switch. **Note the location of the wires so they will be connected correctly to the new micro switch.**
10. Press the switch clips to remove it from the support.
11. Reverse the above steps to install the new micro switch.
12. Test the CPR positioning before returning the bed to service.

## CPR ACTIVATION CABLE REPLACEMENT (OPTIONAL)

### Tools Required:

- 7/16" Combination Wrench

### Note

Unless otherwise indicated, refer to Figure 17, [page 1-59](#), for this procedure.

### Procedure:

1. Raise the bed to the highest position and apply the brakes.
2. Remove the head section plastic cover (if equipped) or raise the head section completely.
3. Remove the head board.
4. Unplug the power cord from the wall receptacle.
5. Using a 7/16" combination wrench, remove the nut, sleeve and bolt (D) holding both ends of the activation cable to the fixed lever.
6. Using a 7/16" combination wrench, remove the nut, washer (E) holding the cable jacket to the mobile lever. Remove the cable from the lever and keep the spring (F). Use caution to avoid moving the other cable from its position.
7. Unscrew the two lock nuts (G) holding the jacket of the other end of the cable to the CPR handle. Remove the wire from the handle. Once the new cable is installed, the free play of the CPR handle must not be more than 1/8".
8. Reverse the above steps to install the new cable.
9. Test the CPR positioning before returning the bed to service.

## CPR PNEUMATIC CYLINDER (OPTIONAL)

### Tools Required:

- Needle Nose Pliers
- (2) 1/2" Combination Wrenches
- 11/16" Combination Wrench
- Strap (or equivalent)

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 17, [page 1-59](#), for this procedure.

1. Raise the bed to the highest position and apply the brakes.
2. Remove the head section plastic cover (if equipped).
3. Unplug the power cord from the wall receptacle.
4. Using needle nose pliers, remove the rue ring, washer, two nylon washers and clevis pin (A, Figure 8, [page 1-40](#)) connecting the head actuator tube to the seat section lever arms.

#### Note

Apply grease to the clevis pin and the nylon washers before reconnecting the actuator tube to the head section.

5. Lift the head section completely and secure it using a strap.
6. Using two 1/2" combination wrenches, remove the nut (H) from the bolt holding the lower end of the pneumatic cylinder. Partially remove the bolt until the lower end of the cylinder is free. Keep the shoulder spacer/sleeve (I). Disengage the lower end of the cylinder from the remaining sleeve.
7. Using an 11/16" combination wrench, loosen the lock nut (J) and unscrew the cylinder completely from the end fitting (K). Note the location of the lock nut before loosening it so the new cylinder will be properly screwed into the end fitting.
8. Reverse the above steps to install the new cylinder.
9. Test the CPR mechanism for proper operation before returning the bed to service.

# Service Information

## BRAKE/STEER PEDAL REPLACEMENT

### Tools Required:

- 1/2" Combination Wrench
- Soft Hammer


Figure 18

### Procedure:

1. Raise the bed to the highest position.
2. Unplug the power cord from the wall receptacle.
3. Remove the optional steer wheel hood.

### Note

When reinstalling the steer wheel hood, ensure the colored labels affixed to the hood match the brake/steer pedal.

4. Position the brake/steer pedal to neutral.
5. Using a 1/2" combination wrench, remove the two bolts (A) holding the steer wheel support on the side of the defective pedal and lay the support down. Keep the molded spacer (B).
6. Using a soft hammer, remove the pedal from the activation lever shaft.
7. Reverse the above steps to install the new pedal.

## STEER WHEEL CASTER REPLACEMENT

### Tools Required:

- (2) 1/2" Combination Wrenches

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 18, [page 1-62](#), for this procedure.

1. Raise the bed to the highest position.
2. Unplug the power cord from the wall receptacle.
3. Remove the optional steer wheel hood.

#### Note

When reinstalling the steer wheel hood, ensure the colored labels affixed to the hood and the brake/steer pedal match.

4. Position the brake/steer pedal to neutral.
5. Using two 1/2" combination wrenches, remove the nut/bolt (C) holding the caster to the swing arm.
6. Reverse the above steps to install the new caster.

## STEER WHEEL ASSEMBLY REPLACEMENT

### Tools Required:

- 1/2" Combination Wrench
- Soft Hammer

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 18, [page 1-62](#), for this procedure.

1. Raise the bed to the highest position.
2. Unplug the power cord from the wall receptacle.
3. Remove the optional steer wheel hood.

### Note

When reinstalling the steer wheel hood, ensure the colored labels affixed to the hood and the brake/steer pedal match.

4. Position the brake/steer pedal to neutral.
5. Using a 1/2" combination wrench, remove the two nuts, two shoulder spacers, two washers, and two bolts (C1, C2, Figure 19, [page 1-66](#)) holding the brake rods to the locking levers on both sides of the mechanism. If the optional four-wheel brake system is present on the bed, also remove the fasteners holding the two other brake rods to the locking levers.
6. Using a 1/2" combination wrench, remove the two bolts (A) holding each side of the steer wheel mechanism to the base. Keep the molded spacers (B). Remove the assembly.
7. Using a soft hammer, remove the two brake/steer pedals from the activation lever shafts.
8. Keep the two locking levers (D, Figure 19, [page 1-66](#)).
9. Reverse the above steps to install the new steer wheel assembly.
10. Test the steer wheel before returning the bed to service.

## STEER WHEEL SWING ARM ASSEMBLY REPLACEMENT

### Tools Required:

- (2) 1/2" Combination Wrench
- 3/16" Allen Wrench

### Procedure:

#### Note

Unless otherwise indicated, refer to Figure 18, [page 1-62](#), for this procedure.

1. Raise the bed to the highest position.
2. Position the brake/steer pedal to neutral.
3. Remove the optional steer wheel hood.

#### Note

When reinstalling the steer wheel hood, ensure the colored labels affixed to the hood and the brake/steer pedal match.

4. Using a 1/2" combination wrench, remove the nut, shoulder spacer, washer, and bolt (C1, Figure 19, [page 1-66](#)) holding the brake rod to the locking lever on the right side of the steer wheel mechanism. If the optional four-wheel brake system is present on the bed, also remove the fasteners holding the other brake rod to the locking lever.
5. Using a 1/2" combination wrench, remove the two bolts (A) holding the right steer wheel support to the base. Lay the support down and keep the molded spacer (B).
6. Using two 1/2" combination wrenches, remove the nut, shoulder spacers and bolt (D) holding the right torsion lever (E) to the lower part of the right counter lever.
7. Using a 1/2" combination wrench and a 3/16" Allen wrench, remove the two nuts and two Allen screws (F) holding the torsion levers (E) to both ends of the torque shaft (G).

#### Note

When reassembling the steer wheel swing arm assembly, screw down the Allen screws before installing the nuts.

8. Disengage the swing arm assembly from the torsion levers and remove it.
9. Using two 1/2" combination wrenches, remove the nut/bolt (C) holding the caster to the swing arm and install the caster on the new swing arm assembly.
10. Reverse the above steps to install the new swing arm assembly.
11. Test the steer wheel before returning the bed to service.

# Service Information

## BED CASTER REPLACEMENT

### Tools Required:

- #2 Phillips Screwdriver
- Jack Stand (or equivalent)
- Needle Nose Pliers
- 1/2" Combination Wrench



Figure 19

### Procedure:

1. Raise the bed to the highest position.
2. Position the brake/steer pedal to neutral.
3. Place a jack stand under the frame where the defective caster is.
4. Using a #2 Phillips screwdriver, remove the screw holding the wheel cover and remove it.
5. Lower the bed until the caster is six inches off the ground.
6. Unplug the power cord from the wall receptacle.
7. If the caster is part of the braking mechanism, using needle nose pliers, remove the rue ring cotter, washer, and locking axle (A) linking the caster shaft to the locking lever.

### Note

The shaft of a caster that is part of the braking mechanism must be properly oriented before mounting it to the base. The opposite illustration shows the side of the shaft that should be facing the inside of the bed. It is the side showing the largest hole below the cam in the neutral position.

8. Using a 1/2" combination wrench, remove the bolt and two washers (B) holding the caster to the base.

### Note

Before tightening the bolt (B), make sure the shaft is completely inserted into the mounting socket.

9. Reverse the above steps to install the new caster.

## BRAKE ROD REPLACEMENT

### Tools Required:

- 1/2" Combination Wrench

### Procedure:

### Note

Unless otherwise indicated, refer to Figure 19, [page 1-66](#), for this procedure.

1. Raise the bed to the highest position.
2. Remove the steer wheel hood.


### Note

When reinstalling the steer wheel hood, ensure the colored labels affixed to the hood and the brake/steer pedal match.

3. Using a 1/2" combination wrench, remove the two nuts, two shoulder spacers, two washers and two bolts (C1, C2) holding each end of the brake rod to the locking levers.
4. Remove the defective rod.
5. Reverse the above steps to install the new brake rod.
6. Test the brakes before returning the bed to service.

# Brake System

OL280003 Rev M (Reference Only)



# Brake System

## Brake System - OL280003 Rev M (Reference Only)

Item	Part No.	Part Name	Qty.
1	25-0291P	Foot Connecting Rod	2
2	QPA25-0340	Spacer	8
3	25-0451	Lever Spacer	2
4	28-0165P	Head Connecting Rod	2
5	28-0340P	Brake Hammer	1
6	28-0341P	Brake Hammer Stop	1
8	QDF17-0020	Shoulder Spacer	1
9	QDF7878	Rue Ring Cotter	4
10	QPA23-0233	Locking Axle	4
11	QPA25-0411	Brake Lever	4
12	QPA25-0427	Locking Lever	2
13	RD25-0388	6" Brake Caster	4
15	VB35A1028FT	Long Neck Carriage Bolt	9
17	VE79A1O	Locknut	9
18	VV83A9G16	Pan Head Tapping Screw	1
19	VW10A10	Flat Washer	9
20	VW10C162404	Nylon Washer	4
22	M0019	OG2 Petro-Canada Grease	1

# Fifth Wheel Assembly

L28-006 Rev J (Reference Only)


# Fifth Wheel Assembly

## Fifth Wheel Assembly - L28-006 Rev J (Reference Only)

Item	Part No.	Part Name	Qty.
1	QPA25-0340	Spacer	4
2	17-0104Z	Fifth Wheel Support	1
3	17-0124P	Fifth Wheel Spring	1
4	QP28-0773	Fifth Wheel Bushing	1
5	QP25-0581	Nylon Bushing	2
6	28-0132Z	Right Fifth Wheel Support	1
7	28-0133Z	Fifth Wheel Left Support	1
10	28-0151P	Clearance Arm	1
11	28-0156P	Activation Lever	1
12	17-0082P	Torque Shaft	1
13	QE71-0496	Fifth Wheel Green Label	2
14	QE71-0511	Brake Label	2
15	QP28-0130	Pedal	2
16	QPA25-0312	Molded Spacer	2
17	QRD25-0582	Left Torsion Spring	1
18	QRD25-0585	Right Torsion Spring	1
19	RL5	Fifth Wheel	1
20	VB15A1O32	Hex Bolt	4
21	VB15A1O50	Hex Bolt	1
22	VE30A1O	Nylon Locknut	7
23	VV11A1O36	Hex Head Flat Cap Screw	2
24	VV11B1P48	Hex Head Flat Cap Screw	1
25	VVB4A1O32	Hex Head Bolt	4
26	28-0149Z	S.A. Lever	2
27	28-0150Z	Torsion Lever	4
28	QP25-0156	Fifth Wheel Cover	1
29	QE71-0904	No Oxygen Bottle Label	2

# Caster Hood Assembly

L28-020 Rev E (Reference Only)


Item	Part No.	Part Name	Qty.
1	25-0108P	Base	1
2	25-0550	Caster Hood Velcro	4
4	QP28-0299	Right Wheel Cover	2
5	QP28-0300	Left Wheel Cover	2
7	VV81A9E16-13	Flat Head Tapping Screw	4
8	VV83A9G16	Pan Head Tapping Screw	1
9	28-0469	Ground Chain	1

# Base Tube Cover Assembly

OL280090 Rev 00 (Reference Only)

English


Item	Part No.	Part Name	Qty.
1	QP25-0023	Base Tube Cover	4
2	VV23A1G24HL	Pan Head Tapping Screw	12

[Return To Table of Contents](#)

# Lift Mechanism

L25-006 Rev J (Reference Only)

# Lift Mechanism


---

## Lift Mechanism - L25-006 Rev J (Reference Only)

Item	Part No.	Part Name	Qty.
1	14662Z	Lever Clamp	4
2	25-0281P	Stabilizer	2
3	25-0402P	Hi-Lo Lever	2
4	QDF17-0020	Shoulder Spacer	4
5	QP25-0270	Molded Support Bearing	4
6	QPN-1310P009	Nylon Glide	4
7	VVB4A1O32	Rolling Thread Bolt	4
8	VB15A1O32	Hex Head Bolt	8
9	VE30A1O	Nylon Locknut	8
10	VW10A10	Flat Washer	8

# Lower Stationary Frame with Scale Option, without Network

OL280005 Rev E (Reference Only)

Orient as shown.

Connector Positions


# Lower Stationary Frame with Scale Option, without Network

## Lower Stationary Frame with Scale Option, without Network - OL280005 Rev E (Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0544	Fixed Lower Frame with Load Cells	1
2	25-0527Z	Connector Sleeve	2
3	28-0186P	Casing Cover	1
4	28-0508P	Scale Board Casing, without Network	1
5	QDF28-0265	Scale Board	1
6	QDF14-1367	Load Cell with Flexible Cable - Long	2
7	QDF17-0138	PC Board Ground Wire	1
8	QDF25-0218	Load Cell with Flexible Cable - Short	2
9	QDF28-0234	Left Harness - Scale without Network	1
10	QDF28-0259	Scale Control/Footboard Cable	1
11	QDF5096	Plastic Tie Holder	6
12	QDF7878	Rue Ring Cotter	4
13	QDF9576	PC Board Support	5
14	QDF9518	Cable Tie	19
15	QE71-0820-T	Connector Position Label	1
16	QP20-0037	Elastomer Sleeve	8
17	VE30A0G	Nylon Locknut	3
18	VE30A1O	Nylon Locknut	8
20	VG50A1244	Clevis Pin	4
21	VV10A1O44	Hex Socket Cap Bolt	8
22	VV33A0G28	Pan Head Machine Screw	2
23	VV83A9G16	Pan Head Tapping Screw	6
25	VW10A12	Washer	4
26	QDF8041	Double Plug Amp	1
27	VV33A0G16	Pan Head Machine Screw	1
28	QDF5095	Beeper	1
29	QE71-1001-T	Label - Caution	1

# Upper Stationary Frame without Scale Option

OL280029 Rev 04 (Reference Only)

Head End

Foot End


# Upper Stationary Frame without Scale Option

## Upper Stationary Frame without Scale Option - OL280029 Rev 04 (Reference Only)

Item	Part No.	Part Name	Qty.
1	QP28-0182	Protective Cap	1
2	25-0519P	Bumper Bracket	2
3	28-0475P	Head End Casing Lateral Cover	2
4	28-0210P	Upper Fixed Frame without Scale	1
5	28-0211P	Head End Casing	1
6	AAL28-0134	Head/Footboard Extrusion	2
7	AAL28-0144	Accessory Extrusion	2
8	QDF28-0235	Left Harness, without Scale/Network	1
9	QPC14-0321	Bumper Roller	6
10	QPCF1001	Crutch Tip	2
11	QPAG1801	Rod Support	2
12	VB15A1N24	Bolt	4
13	VB15A1O62	Hex Bolt	2
14	VB35A1O72-13	Carriage Bolt	8
15	VE30A0G	Nylon Locknut	2
16	VE30A1N	Nylon Locknut	4
17	VE30A1O	Nylon Locknut	10
18	VV33A0G28	Pan Head Machine Screw	2
19	VV83A9G12	Pan Head Tapping Screw	14
20	VW10A10	Flat Washer	2

# Upper Stationary Frame with Scale Option

OL280028 Rev 04 (Reference Only)


Item	Part No.	Part Name	Qty.
1	25-0519P	Bumper Bracket	2
2	28-0141P	Upper Fixed Frame	1
3	28-0475P	Head End Casing Lateral Plate	2
4	AAL28-0134	Head/Footboard Extrusion	2
5	AAL28-0144	Accessory Extrusion	2
6	QPC14-0321	Bumper Roller	6
7	QPCF1001	Crutch Tip	2
8	QPAG1801	Rod Support	2
9	VB15A1N24	Bolt	4
10	VB15A1O62	Hex Bolt	2
11	VB35A1O72-13	Carriage Bolt	8
12	VE30A1N	Nylon Locknut	4
13	VE30A1O	Nylon Locknut	10
14	VV83A9G12	Pan Head Tapping Screw	12
15	VW10A10	Flat Washer	2

# Notes


---

# Moving Frame Assembly

L28-005 Rev M (Reference Only)

# Moving Frame Assembly


Moving Frame Assembly - L28-005 Rev M (Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0420P	Frame Cover Plate	1
2	28-0107P	Mobile Frame	1
3	28-0220P	Foot End Casing Cover	1
4	28-0272Z	Anchor Plate	2
5	AAL28-0134	Head/Foot Board Extrusions	2
6	AAL28-0144	Accessory Extrusion	2
7	M0019	Petro Canada OG2 Grease	1
8	QDF7881	Outer Retaining Ring	6
10	QP23-0257	Mattress Support Glide	2
11	QP28-0111	Foot End Molded Section	1
12	QPN14-0705	Roller	6
13	QPPF1514-1	White Plastic Cap	2
14	QDF1515	Plastic Washer	2
15	VB35A1O72-13	Carriage Bolt	8
16	VE30A1O	Nylon Locknut	8
18	VV33A1N16	Pan Head Machine Screw	4
19	VV37A1N57	Truss Head Machine Screw	2
20	VV83A9G16	Pan Head Tapping Screw	14
21	QDF132X	Weather Strip Tape	1
22	QDF233X	Cushion Tape	1


# Auto Contour without Scale Option

OL280030 Rev B (Reference Only)

See OL280008

Foot End

Green

White

Orange

Yellow

Attach right wire harness here.  
See OL280008

# Auto Contour without Scale Option



---

## Auto Contour without Scale Option - OL280030 Rev B (Reference Only)

Item	Part No.	Part Name	Qty.
1	1325P003	Micro Switch	1
2	25-0465P	Micro Switch Support	1
3	25-0468Z	Micro Switch Actuator	1
4	25-0469P	Micro Switch Support	1
5	26-0067P	Micro Switch Actuator Frame	1
6	QDF9159	Limit Switch	1
7	VV83A9G16	Pan Head Tapping Screw	8
8	VV87A9A20	Truss Head Tapping Screw	2

# Auto Contour with Scale Option

OL280031 Rev 02 (Reference Only)

# Auto Contour with Scale Option


---

## Auto Contour with Scale Option - OL280031 Rev 02 (Reference Only)

Item	Part No.	Part Name	Qty.
1	1325P003	Micro Switch	1
2	25-0465P	Micro switch Support	1
3	25-0468Z	Micro switch Activator	1
4	25-0469P	Micro switch Support	1
5	25-0470P	Micro switch Activator	1
6	QDF9159	Limit Switch	1
7	VV83A9G16	Pan Head Tapping Screw	8
8	VV87A9A20	Truss Head Tapping Screw	2


# 120 VAC Electrical System

OL280120 Rev C (Reference Only)


# 120 VAC Electrical System

**Detail B**

Orient as shown.



**Detail A**



[Return To Table of Contents](#)

# 120 VAC Electrical System

Detail D

Detail C

# 120 VAC Electrical System


## 120 VAC Electrical System - OL280120 Rev C (Reference Only)

Item	Part No.	Part Name	Qty
1	28-0770	Gatch Motor	1
2	25-0250Z	Motor Pivot	8
3	25-0251Z	Head/Thigh Motor Support Cover	4
4	25-0521	Spacer	4
5	25-0527Z	Connector Sleeve	2
6	28-0768	Lift Motor	2
7	28-0769	Fowler Motor	1
8	QDF28-0334	Cap	3
9	QDF17-0138	PC Board Ground Wire	1
10	QDF9574	Power Connector	1
11	QDF25-0508	Main Grounding Wire	1
12	QDF23-0315	Control Board	1
13	QDF28-0231	Right Harness	1
14	QDF28-0259	Scale Control Board/Footboard Cable	1
16	QDF7878	Rue Ring Cotter	2
17	QDF8078	10A Fuse	2
18	QDF9518	Cable Tie	32
19	QDF9532	Wire Clip	2
20	QDF9533	Pressure Cable Clip	1
21	QDF9541	Strain Relief	3
22	QE71-0571	10A, 250 VAC Fuse Label	1
23	QE71-0572	Ground Label	1
25	QE71-1277	120 VAC CSA Label	1
26	QP23-0256	Screw Cover	4
27	QP23-0735	Control Board Stand-off Pin	8
28	QP28-0131	Liquid Retention Device	1
29	VB15A1O40	Hex Bolt	4
30	VE30A0G	Nylon Locknut	5
31	VE80A0G	K-lok Locknut	2
33	VG50A1250	Clevis Pin	2
34	VV41A1A20	Pan Head Tapping Screw	2
35	VV33A0G16	Pan Head Machine Screw	1
36	VV33A0G28	Pan Head Machine Screw	4
37	VV83A9E12	Pan Head Tapping Screw	4
38	VV83A9G12	Pan Head Tapping Screw	4
39	VV84A1H16	Hex Head Tapping Screw	10
40	VW10C122402	Nylon Washer	2
41	VW10C173602	Nylon Washer	5
42	VW20A10	Spring Washer	4
43	QDF28-0238	120 VAC Connector	1
44	QDF8042	AMP Cap Housing	1
45	QDF9575	Fuse Holder	1
46	M0019	Petro Canada OG2 Grease	1
47	QPPF1518	Double "D" Plug	1
48	QE71-0898-T	Serial Number Plate	1
49	VR11H43	Dome Head Pop Rivet	2
50	VV83A9G16	Pan Head Tapping Screw	4
51	QDF8047	PC Board Jumper	1

# Notes


---

## Straight Plug - QDF8066

For Reference Only: OL250053 Rev B


## 90° Plug - QDF8066-90D


For Reference Only: OL250055 Rev B

# International Electrical System

OL280121 Rev C (Reference Only)


[Return To Table of Contents](#)

# International Electrical System



# International Electrical System

## International Electrical System - OL280121 Rev C (Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0770	Gatch Motor	1
2	25-0250Z	Motor Pivot	8
3	25-0251Z	Head/Thigh Motor Support Cover	4
4	25-0521	Spacer	4
5	25-0527Z	Connector Sleeve	2
6	28-0768	Lift Motor	2
7	28-0769	Fowler Motor	1
8	QDF28-0334	Cap	3
9	QDF17-0138	PC Board Ground Wire	1
10	QDF9574	Power Connector	1
11	QDF25-0508	Main Grounding Wire	1
12	QDF23-0315	Control Board	1
13	QDF28-0231	Right Harness	1
14	QDF28-0259	Scale Control Board/Footboard Cable	1
16	QDF7878	Rue Ring Cotter	2
18	QDF9518	Cable Tie	32
19	QDF9532	Wire Clip	2
20	QDF9533	Pressure Cable Clip	1
21	QDF9541	Strain Relief	3
23	QE71-0572	Ground Label	1
26	QP23-0256	Screw Cover	4
27	QP23-0735	Control Board Standoff	8
28	QP28-0131	Liquid Retention Device	1
29	VB15A1O40	Hex Bolt	4
30	VE30A0G	Nylon Locknut	5
31	VE80A0G	K-lok Locknut	2
33	VG50A1250	Clevis Pin	2
34	VV41A1A20	Pan Head Tapping Screw	2
35	VV33A0G16	Pan Head Machine Screw	1
36	VV33A0G28	Pan Head Machine Screw	4
37	VV83A9E12	Pan Head Tapping Screw	4
38	VV83A9G12	Pan Head Tapping Screw	4
39	VV84A1H16	Hex Head Tapping Screw	10
40	VW10C122402	Washer	2
41	VW10C173602	Nylon Washer	5
42	VW20A10	Spring Washer	4
44	QDF8042	AMP Cap Housing	1
46	M0019	Petro Canada OG2 Grease	1
47	QPPF1518	Double "D" Plug	1
48	QE71-0898-T	Serial Number Plate	1
49	VR11H43	Pop Rivet	2
50	VV83A9AG16	Pan Head Tapping Screw	6
51	QDF8047	PC Board Jumper	1
52	28-0191P	Transformer Bracket	1
53	VE79A1O	Locknut	1
54	VB35A1O50	Carriage Bolt	1
55	QDF14-1160	Transformer	1
56	VW10A10	Washer	1










# 220VAC Electrical System

OL280107 Rev C (Reference Only)

Item	Part No.	Part Name	Qty.
1	QDF20-0083	220 VAC Connector	1
2	QDF8068	6.3A, 250 VAC Fuse	2
3	QE71-0573	Label - 6.3A, 250 VAC Fuse	1
4	QE71-1280-T	Label - CSA 220 VAC	1
5	QDF9580	Fuse Holder	1


# International Power Cord Options

[Return To Table of Contents](#)

# Head End Cover without Scale, Nurse Call and Pendant Port


OL280096 Rev D (Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0677P	Head Cover without Connector	1
2	QDF2155	Steel Cable Clamp	1
3	VV83A9G12	Pan Head Tapping Screw	6

# Head End Cover with Nurse Call and Pendant Port, without Scale

OL280101 Rev D (Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0678P	Head End Cover	1
2	28-0285P	Nurse Call Mounting Plate	1
3	QDF20-0121	GENIII Interface	1
4	QDF5097	GENIII Connector Support	1
5	QDF8047	Electronic Board Jumper	2
7	VV83A9G12	Pan Head Tapping Screw	12
8	QDF2155	Steel Cable Clamp	1

[Return To Table of Contents](#)

# Head End Cover with Pendant Port, without Scale and Nurse Call

OL280102 Rev D (Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0678P	Head End Cover	1
2	28-0285P	Nurse Call Mounting Plate	1
3	QDF20-0121	GENIII Interface	1
4	QDF5097	GENIII Connector Support	1
5	QDF8047	Electronic Board Jumper	4
7	VV83A9G12	Pan Head Tapping Screw	12
8	QDF2155	Steel Cable Clamp	1

# Head End Cover with Nurse Call/Scale/GENIII/37 Pin Connectors

OL280099 Rev E (Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0285P	Nurse Call Mounting Plate	1
2	28-0678P	Head End Cover	1
3	QDF20-0121	Nurse Call and GENIII Control Board	1
4	QDF20-0137	Bed Exit Cable	1
5	QDF2155	Steel Cable Clamp	1
6	QDF5097	GENIII Connector Support	1
7	VV83A9G12	Pan Head Tapping Screw	12

[Return To Table of Contents](#)

# 120 VAC Auxiliary Outlet

OL250029 Rev D (Reference Only)


Item	Part No.	Part Name	Qty.
1	25-0408P	Auxiliary Outlet Casing	1
2	25-0409P	Auxiliary Outlet Bracket	1
3	QDF17-0199	Auxiliary 120V Plug Wire	1
4	QDF17-0223	Circuit Breaker Wire	1
5	QDF8024	Double Socket Outlet	1
6	QDF9025	Circuit Breaker - 5A	1
7	QE71-0225-T	Circuit Breaker Label	1
8	QE71-0742-T	Auxiliary Outlet Label	1
9	VV37A1C12	Truss Head Machine Screw	1
10	VV83A9G16	Pan Head Tapping Screw	4
11	QPCK3001	Rubber Play Bumper	1

# Night Light

OL250018 Rev E (Reference Only)

## Foot End


Item	Part No.	Part Name	Qty.
1	25-0273Z	Night Light Support	1
2	QDF18825	Interlocking Outlet	1
3	QDF9539	Night Light	1
4	QE71-0712	Night Light Label	1
5	VV83A9G16	Pan Head Tapping Screw	2

[Return To Table of Contents](#)



# CPR Mechanism

OL280034 Rev T (Reference Only)

English

# CPR Mechanism

[Return To Table of Contents](#)

# CPR Mechanism

## CPR Mechanism - OL280034 Rev T (Reference Only)

Item	Part No.	Part Name	Qty.
3	20-0017Z	Activation Mechanism Socket	1
4	20-0018Z	Mobile Lever	1
5	20-0019Z	Fixed Lever	1
6	25-0362Z	Mattress Support Lever	1
7	25-0456Z	Sleeve	1
8	25-0514Z	Cylinder Inner Sleeve	1
9	25-0515Z	Cylinder Outer Sleeve	1
10	28-0242P	CPR Handle	2
11	M0008	Threadlocker (Blue)	-
12	M0019	Petro Canada OG2 Grease	1
13	QDF17-0020	Shoulder Spacer	7
14	QDF19-0354	Head Section Cable	1
15	QDF19-0815	Short Head Section Cable	1
16	QDF5090	CPR Pneumatic Cylinder	1
17	QDF9159	Cherry Switch	1
18	QDF9518	Cable Tie	3
19	QE71-0963-T	Label - CPR	2
20	QRC18579	Compression Spring	2
21	VB15A1N44	Hex Bolt	1
22	VB15A1O48	Hex Bolt	1
23	VB15A1O54	Hex Bolt	1
24	VB15A1O56	Hex Bolt	2
25	VB18A1O52	Hex Bolt	1
26	VB35A1N24	Carriage Bolt	1
27	VE20A80	Locknut	1
28	VE30A1N	Nylon Hex Locknut	2
29	VE30A1O	Nylon Hex Locknut	5
32	23-0762	Nylon Shoulder Washer	1
33	VV87A9A24	Truss Head Tapping Screw	2
34	1325P003	Micro Switch	1
35	VV60B1O12	Set Screw	1
36	17-0192	Nut	1
37	28-0402	Spacer	1
38	QDF5096	Flat Tie Holder	1


# Notes

---

# Foot Section with Support Rod

English

OL280047 Rev 01 (Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0804P	Support Rod	1
2	28-0461Z	Right Hook	1
3	28-0468Z	Left Hook	1
4	QPNI2404	Nylon Bushing with Shoulder Spacer	2
5	VV83A9G16	Pan Head Tapping Screw	2

# Foot End Mattress Retainer

OL250023 Rev 00 (Reference Only)


## Foot End Mattress Retainer - OL250023 (82"/84") Rev 00 (For Reference Only)

Item	Part No.	Part Name	Qty.
1	17-0211P	Foot End Mattress Retainer	1


[Return To Table of Contents](#)

# Litter with Plastic Covers


OL280032 Rev 03 (Reference Only)


**Detail A**

**Detail B**

**Detail C**


# Litter with Plastic Covers

## Litter with Plastic Covers - OL280032 Rev 03 (For Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0342P	Head Section Lever	2
2	28-0530P	Head Section Frame	1
3	28-0345P	Seat Section	1
4	28-0034P	Thigh Section	1
5	28-0035P	Foot Section Frame	1
6	28-0354P	Protective Plate	1
7	28-0101	Plastic Head Section	1
8	28-0103	Plastic Foot Section	1
9	QDF17-0020	Shoulder Spacer	10
10	QP14034-13	Mattress Retainer	2
11	VV11A1O44	Flat Head Hex Socket Cap Screw	4
12	VE30A1O	Nylon Hex Locknut	14
13	VV83A9G16	Pan Head Tapping Screw	3
14	VV83A9G24	Pan Head Tapping Screw	4
16	VW10A10	Flat Washer	14
18	VB35A1O32	Carriage Bolt	10

# Litter with Steel Covers

OL280007 Rev E (Reference Only)


# Litter with Steel Covers

## Litter with Steel Covers - OL280007 Rev E (For Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0342P	Head Litter Lever	2
2	28-0345P	Seat Section	1
3	28-0034P	Thigh Section	1
4	28-0354P	Protective Angle	1
5	28-0196P	Foot Section - Steel	1
6	28-0529P	Head Section - Steel	1
7	QDF17-0020	Spacer	10
8	QP14034	Mattress Retainer	2
9	VV11A1O44	Flat Head Hex Socket Cap Screw	4
10	VE30A1O	Nylon Hex Locknut	14
11	VV83A9G16	Pan Head Tapping Screw	3
12	VB35A1O32	Carriage Bolt	10
13	VW10A10	Flat Washer	14
14	QDF8086	Cap	8
16	VV83A9G24	Pan Head Tapping Screw	4

# Head End Siderail Assembly


L28-012 Rev C (Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0244	Right Head Siderail	1
2	28-0245	Left Head Siderail	1
3	QDF5096	Flat Tie Holder	4
4	QDF9518	Cable Tie	4
5	VE30A10	Nylon Hex Locknut	8
6	VV11A1048	Flat Head Hex Socket Cap Screw	8


# Head End Siderail Assembly, Right

28-0244 Rev R (Reference Only)


[Return To Table of Contents](#)

# Head End Siderail Assembly, Right

Assembly Detail on part #16

Assembly Detail on part #18

## Head End Siderail Assembly, Right - 28-0244 Rev R (Reference Only)


Item	Part No.	Part Name	Qty.
1	QPA28-0493	Siderail Sleeve	1
2	28-0006	Siderail Mechanism, Reverse	1
3	QDB28-0278Z	Siderail Latch	1
4	QDF2042	14 Position Connector	1
5	QDF2043	Cable Clamp	1
6	QDF2049	Bushing	4
7	QDF26-0108	Siderail Control Module Cable, Right	1
8	QDF28-0145	Screw Cover	2
9	QDF28-0146	Label - Siderail Latch	1
10	QDF28-0246	Label - Nurse Control, Right	1
11	QDF7826	3M Tape #950	1
12	QDF9018	Amp Pin Contact	7
13	QDF28-0491	Bushing	1
14	QDF9518	Cable Tie	4
15	QP28-0489	Siderail Handle, Right	1
16	QP26-0011	Inner Siderail Panel, Head, Right	1
17	QP28-0040	Siderail Cover	2
18	QP28-0521	Siderail, Head, Right	1
19	QPA28-0015	Siderail Sleeve	1
20	QPA28-0016	Sleeve with Lock	1
21	QPA28-0095	Siderail Frame, Head, Right	1
22	QRE25-0245	Spring	1
23	QRE28-0293	Spring	1
24	VV10B0G40-S	Hex Socket Head Cap Screw	4
25	VV10A1P20-S	Hex Socket Low Head Cap Screw	1
26	VV31A0G16	Flat Head Machine Screw	2
27	VV33A0G28-S	Pan Head Machine Screw	8
28	VVB3A9N24PF	Pan Head Screw	6
29	QDF2053	PVC Tape, Red	1
30	M0019	OG2 Grease	1

# Siderail Assembly, Reverse - 28-0006

Rev 13

English


Item	Part No.	Part Name	Qty.
1	28-0283P	Siderail Arm, Reverse	2
2	28-0284P	Siderail Base Angle Bracket	1
3	28-0014Z	Transfer Plate	1
4	QPA28-0086	Transfer Siderail Arm, Reverse	2
5	QRC28-0208	Compression Spring	2
6	28-0219Z	Sliding Guide, Reverse	2
7	QPA27-1720	Bushing	2
8	QDF28-0270	Flanged Bushing	2
9	QDF28-0268	Flanged Bushing	2
10	M0019	OG2 Grease	1
11	QDF2005-S	Ball Stud	2
12	QDF2033	Damper (not shown)	1
13	QDF7894	Bushing	2
14	VD60A1N1016-S	Shoulder Bolt	4
15	VV10A0G16-S	Hex Socket Head Cap Screw	4

[Return To Table of Contents](#)


# Head End Siderail Assembly, Left

28-0245 Rev R (Reference Only)


[Return To Table of Contents](#)

# Head End Siderail Assembly, Left

Assembly Detail on part #15

Assembly Detail on part #17


## Head End Siderail Assembly, Left - 28-0245 Rev R (Reference Only)

Item	Part No.	Part Name	Qty.
1	QPA28-0493	Siderail Sleeve	1
2	28-0005	Siderail Mechanism	1
3	QDB28-0278Z	Siderail Latch	1
4	QDF2042	14 Position Connector	1
5	QDF2043	Cable Clamp	1
6	QDF2049	Bushing	4
7	QDF26-0129	Siderail Control Module Cable, Left	1
8	QDF28-0145	Screw Cover	2
9	QDF28-0146	Label - Siderail Latch	1
10	QDF28-0247	Label - Nurse Control, Left	1
11	QDF9018	Amp Pin Contact	7
12	QDF28-0491	Bushing	1
13	QDF9518	Cable Tie	4
14	QP28-0490	Siderail Handle, Foot	1
15	QP26-0012	Inner Siderail Panel, Head, Left	1
16	QP28-0039	Siderail Cover	2
17	QP28-0522	Siderail, Head, Left	1
18	QPA28-0015	Siderail Sleeve	1
19	QPA28-0016	Sleeve with Lock	1
20	QPA28-0096	Siderail Frame, Head, Left	1
21	QRE25-0245	Spring	1
22	VV10B0G40-S	Hex Socket Head Cap Screw	4
23	VV10A1P20-S	Hex Socket Low Head Cap Screw	1
24	VV31A0G16	Flat Head Machine Screw	2
25	VV33A0G28-S	Pan Head Machine Screw	8
26	VVB3A9N24PF	Pan Head Screw	6
27	QRE28-0293	Spring	1
28	QDF7826	3M Tape #950	1
29	QDF2053	PVC Tape, Red	1
30	M0019	OG2 Grease	1

# Siderail Assembly - 28-0005

Rev 13


English

Item	Part No.	Part Name	Qty.
1	28-0282P	Siderail Arm	2
2	28-0284P	Siderail Base Angle Bracket	1
3	QPA28-0011	Transfer Siderail Arm	2
4	28-0014Z	Transfer Plate	1
5	QRC28-0208	Compression Spring	2
6	28-0218Z	Sliding Guide	2
7	QDF28-0268	Flanged Bushing	2
8	QDF28-0270	Flanged Bushing	2
9	QPA27-1720	Bushing	2
10	QDF2005-S	Ball Stud	2
11	QDF2033	Damper	1
12	QDF7894	Bushing	2
13	VD60A1N1016-S	Shoulder Bolt	4
14	VV10A0G16-S	Hex Socket Head Cap Screw	4
15	M0019	OG2 Grease	1

# Head End Siderail Labels without Nurse Call

OL280027 Rev 01 (Reference Only)


Item	Part No.	Part Name	Qty.
1	QDF28-0254	Right Patient Control Label	1
2	QDF28-0255	Left Patient Control Label	1
3	QDF7826	3M Adhesive Tape	1
4	QDF28-0335	Speaker Label	2

[Return To Table of Contents](#)

# Head End Siderails with Nurse Call

OL280026 Rev 02 (Reference Only)


Item	Part No.	Part Name	Qty.
1	QDF28-0336	Speaker Label	2
2	QDF26-0111	Speaker with Connector	2
4	QDF26-0113	Speaker/Nurse Call Cable - Siderail	2
5	QDF28-0256	Patient Control Label with Nurse Call, Right	1
6	QDF28-0257	Patient Control Label with Nurse Call, Left	1
7	QDF7826	3M Adhesive Tape	1
8	QDF9018	AMP Pin Contact	12
9	QDF9518	Cable Tie	1
10	QE71-0804	Nurse Call Isolating Label	2
11	VV23A9C12HL	Pan Head Tapping Screw	8

# Foot End Siderail Assembly

L28-002 Rev C (Reference Only)

English


Item	Part No.	Part Name	Qty.
1	28-0003	Foot Siderail, Right	1
2	28-0004	Foot Siderail, Left	1
3	28-0178P	Siderail Support	1
4	VB18A1O24	Hex Bolt	4
5	VB15A1O32	Hex Bolt	8
6	VE30A1O	Nylon Hex Locknut	12
7	VW10A10	Flat Washer	4

[Return To Table of Contents](#)

# Foot End Siderail Assembly, Right

28-0003 Rev M (Reference Only)


# Foot End Siderail Assembly, Right

## Foot End Siderail Assembly, Right - 28-0003 Rev M (Reference Only)

Item	Part No.	Part Name	Qty.
1	QDB28-0278Z	Siderail Latch	1
2	QDA28-0493	Siderail Sleeve	1
3	28-0005	Siderail Mechanism	1
4	QDF2049	Bushing	4
5	QDF28-0146	Label - Siderail Latch	1
6	QDF28-0145	Screw Cover	3
7	QDF28-0491	Bushing	1
8	QP25-0188	Inner Siderail Panel, Foot, Right	1
9	QP28-0490	Siderail Handle, Foot	1
10	QP28-0039	Siderail Cover	2
11	QP28-0519	Siderail, Foot, Right	1
12	QPA28-0015	Siderail Sleeve	1
13	QPA28-0016	Sleeve with Lock	1
14	QPA28-0097	Siderail Frame, Foot, Right	1
15	QRE25-0245	Spring	1
16	VV10B0G40-S	Hex Socket Head Cap Screw	4
17	VV10A1P20-S	Hex Socket Low Head Cap Screw	1
18	VV31A0G16	Flat Head Machine Screw	2
19	VV33A0G28-S	Pan Head Machine Screw	8
20	VVB3A9N24PF	Pan Head Screw	6
21	QRE28-0293	Spring	1
22	M0019	OG2 Grease	1

# Foot End Siderail Assembly, Left

28-0004 Rev M (Reference Only)


# Foot End Siderail Assembly, Left

## Foot End Siderail Assembly, Left - 28-0004 Rev M (Reference Only)

Item	Part No.	Part Name	Qty.
1	QDB28-0278Z	Siderail Latch	1
2	QDA28-0493	Siderail Sleeve	1
3	28-0006	Siderail Mechanism, Inverse	1
4	QDF2049	Bushing	4
5	QDF28-0145	Screw Cover	3
6	QDF28-0146	Label - Siderail Latch	1
7	QDF28-0491	Bushing	1
8	QP28-0489	Siderail Handle, Right	1
9	QP25-0189	Inner Siderail Panel, Foot, Left	1
10	QP28-0040	Siderail Cover	2
11	QP28-0520	Siderail, Foot, Left	1
12	QPA28-0015	Siderail Sleeve	1
13	QPA28-0016	Sleeve with Lock	1
14	QPA28-0098	Siderail Frame, Foot, Left	1
16	VV10B0G40-S	Hex Socket Head Cap Screw	4
17	VV10A1P20-S	Hex Socket Low Head Cap Screw	1
18	VV31A0G16	Flat Head Machine Screw	2
19	VV33A0G28-S	Pan Head Machine Screw	8
20	VVB3A9N24PF	Pan Head Screw	6
21	QRE25-0245	Spring	1
22	QRE28-0293	Spring	1
23	M0019	OG2 Grease	1

# Headboard Assembly

L28-011 Rev 01 (Reference Only)


Item	Part No.	Part Name	Qty.
1	QP21-3139	Head Board	1
2	M0005	Glue	1
3	QPC21-3855	Rubber Cap	2
4	QE14376-T	Label - Installation	1

# Footboard Assembly

L28-015 Rev C (Reference Only)

English

Item	Part No.	Part Name	Qty.
1	25-0260P	Pump Support	1
2	M0005	Glue	1
4	QDF28-0334	Screw Cover	9
5	QP28-0679	Foot Board	1
6	QP28-0676	Foot Board Cover	1
7	QPC21-3855	Rubber Cap	2
8	VVB3A9N16PF	Pan Head Machine Screw	7


[Return To Table of Contents](#)

# Footboard without Scale and Bed Exit

OP280006 Rev F (Reference Only)

Affix the cable connector to the inside of the panel using an adhesive strip (QDF7826).

# Footboard without Scale and Bed Exit

---

## Footboard without Scale and Bed Exit - OP280006 Rev F (For Reference Only)



Item	Part No.	Part Name	Qty.
1	QP28-0182	Protective Cap	1
2	28-0263P	Connector Reinforcement Plate	1
3	QDF28-0660	Foot End Control Panel Fascia	1
4	QDF26-0109	Foot Board Cable	1
5	VV83A9E24	Pan Head Tapping Screw	4
6	28-0325Z	Connector Reinforcement Wedge	1
7	QDF2053	Red Electric Tape	1
8	QDF7826	3M Tape	1

# Footboard with Scale, without Bed Exit

OP280003 Rev E (Reference Only)

Affix the cable connector to the inside of the panel using an adhesive strip (QDF7826).

**Detail A**

Orient as shown.

See Detail A

# Footboard with Scale, without Bed Exit

---

## Footboard with Scale, without Bed Exit - OP280003 Rev E (For Reference Only)


Item	Part No.	Part Name	Qty.
1	28-0263P	Connector Reinforcement Plate	1
2	QDF2053	Red PVC Tape	1
3	QDF28-0412	Foot Board Cable	1
4	QDF28-0260	Foot Board Connector	1
5	QDF28-0266	Scale Display Control Board	1
6	QDF28-0275	Foot End Control Panel Fascia	1
7	VV83A9E24	Pan Head Tapping Screw	4
8	28-0325Z	Connector Reinforcement Wedge	1
9	QDF7826	3M Tape	1
10	VV33A1N16	Pan Head Machine Screw	2

# Footboard with Scale and Bed Exit - Three Zone

OP280001 Rev E (Reference Only)

Affix the cable connector to the inside of the panel using an adhesive strip (QDF7826).

# Footboard with Scale and Bed Exit - Three Zone



---

## Footboard with Scale and Bed Exit - Three Zone - OP280001 Rev E (For Reference Only)

Item	Part No.	Part Name	Qty.
1	28-0263P	Connector Reinforcement Plate	1
2	QDF28-0412	Foot Board Cable	1
3	QDF28-0260	Foot Board Connector	1
4	QDF28-0266	Scale Display Control Board	1
5	QDF28-0274	Foot End Control Panel Fascia	1
6	VV83A9E24	Pan Head Tapping Screw	4
7	28-0325Z	Connector Reinforcement Wedge	1
8	QDF2053	Red PVC Tape	1
9	QDF7826	3M Tape	1

# Bed Labeling - English

OL280037 Rev 01 (Reference Only)



# Bed Labeling - English

## Bed Labeling - English - OL280037 Rev 01 (For Reference Only)

Item	Part No.	Part Name	Qty.
1	QE14400-F	Label - Automated Test Results	1
2	QE18545	Label - Refer to Manual	1
3	QE71-0301	Label - Lower Hd. Elevation Angle, L	1
4	QE71-0303	Label - Lower Hd. Elevation Angle, R	1
5	QE71-0346	Label - Stryker	2
6	QE71-0547	Label - Maximum Load	2
7	QE71-0699-E	Label - Serious Injury	1
8	QE71-0700-E	Label - Electric Shock Hazard	1
9	QE71-0706-E	Label - Fire Hazard	1
10	QE71-0709-E	Label - Mattress Maximum Thickness	1
11	QE71-0719-E	Label - High Hd. Elevation Angle, R	1
12	QE71-0720-E	Label - High Hd. Elevation Angle, L	1
13	QE71-0721-E	Label - Right Pedal Positions	1
14	QE71-0722-E	Label - Left Pedal Position	1
15	QE71-0875	Label - Emergency Crank Handle	1

# Bed Labeling - English / French / Spanish

OL280036 Rev C (Reference Only)

# Bed Labeling - English / French / Spanish



## Bed Labeling - English / French - OL280036 Rev C (For Reference Only)

Item	Part No.	Part Name	Qty.
1	QE14400-F	Label - Results of Automated Tests	1
2	QE18545	Label - Refer to Manual	1
3	QE71-0247-T	Label - Head Elevation Angle, Left	1
4	QE71-0301	Label - Graduated Scale, Left	1
5	QE71-0302-T	Label - Head Elevation Angle, Right	1
6	QE71-0303	Label - Graduated Scale, Right	1
7	QE71-0346	Label - STRYKER	2
8	QE71-0547	Label - Maximum Load	2
9	QE71-0737-T	Label - Serious Injury	1
10	QE71-0738-T	Label - Electric Shock Hazard	1
11	QE71-0739-T	Label - Right Pedal Positions	1
12	QE71-0740-T	Label - Left Pedal Positions	1
13	QE71-0741-T	Label - Mattress Maximum Thickness	1
14	QE71-0752-T	Label - Fire Hazard	1
15	QE71-0875	Label - Emergency Crank Handle	1

# Bed Labeling - English / Spanish

English

OL280038 Rev 01 (Reference Only)

# Bed Labeling - English / Spanish

## Bed Labeling - English / Spanish - OL280038 Rev 01 (For Reference Only)


Item	Part No.	Part Name	Qty.
1	QE14400-F	Label - Results of Automated Tests	1
2	QE18545	Label - Refer to Manual	1
3	QE71-0247-E	Label - Head Elevation Angle, Left	1
4	QE71-0301	Label - Graduated Scale, Left	1
5	QE71-0302-E	Label - Head Elevation Angle, Right	1
6	QE71-0303	Label - Graduated Scale, Right	1
7	QE71-0346	Label - STRYKER	2
8	QE71-0547	Label - Maximum Load	2
9	QE71-0737-T	Label - Serious Injury	1
10	QE71-0738-T	Label - Electric Shock Hazard	1
11	QE71-0739-T	Label - Right Pedal Positions	1
12	QE71-0740-T	Label - Left Pedal Positions	1
13	QE71-0741-T	Label - Mattress Maximum Thickness	1
14	QE71-0752-E	Label - Fire Hazard	1
15	QE71-0875	Label - Emergency Crank Handle	1

## Optional Accessories

Accessory	Part Number	Page Reference
Pendant Assemblies	FA64136 FA64137 FA64141 FA64142	<a href="#">See page 1-145</a>
Monitor Tray	FA64163	<a href="#">See page 1-146</a>
Patient Helper	FA64148L	<a href="#">See page 1-147</a>
Emergency Crank Handle	FOHMAU	<a href="#">See page 1-148</a>
Oxygen Bottle Holder	FA64169	<a href="#">See page 1-149</a>
Bed Extender without Scale	FA64164	<a href="#">See page 1-150</a>
Bed Extender with Scale	FA64173	<a href="#">See page 1-151</a>
Removable I.V. Pole, 1/2"	FDTSH	<a href="#">See page 1-152</a>
Removable I.V. Pole, 1"	FA61002-G	<a href="#">See page 1-153</a>
Fixed Two Stage I.V. Pole, Foldable	FA64171	<a href="#">See page 1-154</a>
Fixed Three Stage I.V. Pole, Foldable	FA64170	<a href="#">See page 1-155</a>
Traction Sleeve Sets	FA64188 FA64189 FA64190 FA64191	<a href="#">See page 1-156</a>

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
**Nurse Call and TV Control  
Pendant - FA64136**


**Nurse Call, Lights and TV Control  
Pendant - FA64137**

**Bed Height and Head / Foot  
Section Control Pendant - FA64141**


**Head / Foot Section Control  
Pendant - FA64142**

# Monitor Tray - FA64163


L64-085 Rev B (Reference Only)

Item	Part No.	Part Name	Qty.
1	64-0777P	Tray Support	1
2	QDF5093	Tray and Fasteners	1
3	QE14399-T	Label - Manufacturer	1
4	QE71-0533	Label - Protective Top	2
5	QE71-0578	Label - Maximum Load	4
6	QE71-0579-F	Label - Tied Up Load	1
7	QE71-0580-F	Label - Caution	1

# Patient Helper - FA64148

L64-075 Rev E (Reference Only)


Item	Part No.	Part Name	Qty.	Item	Part No.	Part Name	Qty.
1	16424C	Plunger Washer	2	21	VE40A10	Cap Nut	2
2	16243Z	Plunger	1	22	VV83A9G12	Pan Head Tapping Screw	1
3	18939Z	Plunger	1	23	VW10A10	Washer	2
4	64-0700P	Support Tube	1	24	VW10C322802	Nylon Washer	1
5	64-0725P	Trapeze Frame	1	25	VW20A10	Spring Washer	2
6	9870	Triangle Retainer	1	26	VR11H46	Pop Rivet	2
8	QDF6023	Trapeze Handle	1	27	64-0536	Chain	2
9	64-1298P	Trapeze Hook	1	28	QDF7871	Hitch Pin	2
11	QE71-0292-T	Label - Maximum Load	2	29	QE14399-T	Label - Manufacturer	2
12	QE71-0731	Label - Caution	2	30	QE71-1163-T	Label - Handle Storage	1
13	QPC14-0321	Roller Bumper	3				
15	QPPE1205	Cap	1				
16	QPPF1505	Square Tube Cap	1				
17	QR110027	Compression Spring	1				
18	QRC9793	Compression Spring	1				
19	VB15A1O60	Hex Bolt	1				
20	VE30A1O	Nylon Hex Locknut	1				

[Return To Table of Contents](#)

# Emergency Crank Handle - FOHMAU


LC-3203 Rev 03 (Reference Only)

Item	Part No.	Part Name	Qty.
1	16142C	Crank	1
2	QPN-9438	Handle	1
3	QRC16887	Compression Spring	1
4	VG10B0424	Spring Pin	1
5	VR43A9439	Shoulder Rivet	1
6	QE71-0814-F	Label - Emergency Crank	1
7	M0019	Petro Canada OG2 Grease	-
8	QE71-1313-T	Label - Manufacturer	1

# Oxygen Bottle Holder - FA64169

L64-097 Rev C (Reference Only)


Item	Part No.	Part Name	Qty.
1	64-0647	Chain	1
2	QDF64-0973	Bottle Holder	1
3	QE14399-T	Label - Manufacturer	1
4	QE71-0601	Label - Maximum Load	1
5	VR11H42	Pop Rivet	1

[Return To Table of Contents](#)

# Bed Extender without Scale - FA64164

L64-086 Rev F (Reference Only)


Item	Part No.	Part Name	Qty.	Item	Part No.	Part Name	Qty.
1	QP28-0182	Protective Cap	2	17	QP14034	Mattress Retainer	2
2	64-1337P	Foot End Mattress Retainer	1	18	QPNC0607	Nylon Shoulder Washer	2
3	25-0527Z	Connector Sleeve	2	19	VB35A1O72-13	Carriage Bolt	8
4	64-0781P	Bed Extender	1	20	VE30A0G	Nylon Hex Locknut	12
5	64-0996P	Connector Support	1	21	VE30A1O	Nylon Hex Locknut	8
6	64-0783P	Male Connector Support	1	22	VE90A08	Stationary Nut	3
7	64-0995P	Bed Extender Support	1	23	VV33A0G24	Pan Head Machine Screw	2
8	AAL28-0134	Head/Footboard Extrusions	2	24	VV33A0G28	Pan Head Machine Screw	6
9	AAL28-0144	Accessory Extrusion	2	25	VV33A1N16-S	Pan Head Machine Screw	3
10	QDF17-0139	Extension Cord	1	26	VV37A0G16	Truss Head Machine Screw	4
11	QDF2053	PVC Tape, Red	1	27	VV83A9G24	Pan Head Tapping Screw	4
12	QDF28-0259	Foot Board Cable	1	28	VW10C142402	Nylon Washer	8
13	QDF9518	Cable Tie	1				
14	QE14399-T	Label - Manufacturer	1				
15	QE71-0697-T	Label - Extension Warning	1				
16	QE71-0874-F	Label - Bed Extension	1				

[Return To Table of Contents](#)

# Bed Extender with Scale - FA64173

L64-102 Rev F (Reference Only)


Item	Part No.	Part Name	Qty.	Item	Part No.	Part Name	Qty.
1	64-1337P	Foot End Mattress Retainer	1	17	VE30A10	Nylon Hex Locknut	8
2	25-0527Z	Connector Sleeve	4	19	VV33A0G24	Pan Head Machine Screw	4
3	AAL28-0134	Head/Footboard Extrusions	2	20	VV33A0G28	Pan Head Machine Screw	4
4	AAL28-0144	Accessory Extrusion	2	21	VV37A0G16	Truss Head Machine Screw	4
5	64-0995P	Bed Extender Support	1	22	QP14034-B	Mattress Retainer	2
6	64-0781P	Bed Extender	1	23	VV83A9G24	Pan Head Tapping Screw	4
7	64-0782P	Connector Support	1	24	VW10C142402	Nylon Washer	8
8	64-0783P	Male Connector Support	1	25	QDF2053	PVC Tape, Red	1
9	QDF17-0139	Extension Cord	2	26	QE71-0873-F	Label - Bed Extension	1
10	QDF28-0259	Foot Board Cable	2	27	VE90A08	Stationary Nut	3
11	QDF9518	Cable Tie	1	28	VV33A1N16-S	Pan Head Machine Screw	3
12	QE14399-T	Label - Manufacturer	1				
13	QE71-0697-T	Label - Extension Warning	1				
14	QPNC0607	Nylon Shoulder Washer	4				
15	VB35A1O72-13	Carriage Bolt	8				
16	VE30A0G	Nylon Hex Locknut	12				

[Return To Table of Contents](#)

# Removable I.V. Pole - 1/2" - FDTSH


LC-7300 Rev E (Reference Only)

Item	Part No.	Part Name	Qty.
1	L-2113	Hook	1
2	2111A	Tightening Sleeve	1
3	2109	Compression Socket	1
4	2107A	Top Tube	1
5	2105A	Bottom Tube	1
6	QE14399-T	Manufacturer's Label	1
101	QPCE1206	Rubber cap	1
102	QE71-0246-T	Warning Label	1

# Removable I.V. Pole - 1" - FA61002

61-0014 Rev 01 (Reference Only)


## L61-0002 Rev C (For Reference Only)

Item	Part No.	Part Name	Qty.
1	61-0011	Top Rod	1
2	61-0002A	Bottom Tube	1
3	QPCE1206	Rubber cap	1
4	QE71-0246-T	Label - Warning	1
5	QE14399-T	Label - Manufacturer	1
6	VG10B0630	Spring Pin	1

[Return To Table of Contents](#)

# Fixed Two Stage I.V. Pole, Foldable - FA64171


L64-100 Rev D (Reference Only)

Item	Part No.	Part Name	Qty.
1	64-0734C	I.V. Pole Support	1
2	64-1351	I.V. Pole Stopper	1
3	QDF64-0613	Two Stage I.V. Pole	1
4	64-1215P	Storage Bracket Plate	1
5	QE14399-T	Label - Manufacturer	1
6	QE71-0246-T	Label - Warning	1
7	VB15A1O24-S	Hex Bolt	1
8	VB15A1O44	Hex Bolt	1
9	VE30A1O	Nylon Hex Locknut	3
10	VV83A9G16	Pan Head Tapping Screw	2
11	VW10A10	Washer	2
12	90-1854P	Roller Bumper Support	2
13	QDF2094	Roller Bumper	2
14	VB15A1N24	Hex Bolt	4
15	VB15A1O60	Hex Bolt	2
16	VE30A1N	Nylon Hex Locknut	4

# Fixed Three Stage I.V. Pole, Foldable - FA64170

L64-099 Rev D (Reference Only)


Item	Part No.	Part Name	Qty.
1	64-0734C	I.V. Pole Support	1
2	64-1351	I.V. Pole Stopper	1
3	QDF64-0614	Three Stage I.V. Pole	1
4	64-1215P	Storage Bracket Plate	1
5	QE14399-T	Label - Manufacturer	1
6	QE71-0246-T	Label - Warning	1
7	VB15A1O24-S	Hex Bolt	1
8	VB15A1O44	Hex Bolt	1
9	VE30A1O	Nylon Hex Locknut	3
10	VV83A9G16	Pan Head Tapping Screw	2
11	VW10A10	Washer	2
12	90-1854P	Roller Bumper Support	2
13	QDF2094	Roller Bumper	2
14	VB15A1N24	Hex Bolt	4
15	VB15A1O60	Hex Bolt	2
16	VE30A1N	Nylon Hex Locknut	4

[Return To Table of Contents](#)

# Traction Sleeves

Traction Sleeve Set Size	Part Number	Reference Part Number
4" x 1/2"	FA64188	L64-113 (Rev B)
4" x 3/4"	FA64189	L64-114 (Rev B)
8" x 1/2"	FA64190	L64-115 (Rev B)
8" x 3/4"	FA64191	L64-116 (Rev B)


Item	Part No.	Part Name	Qty.
1	64-1035C	Traction Sleeve Set - 4" x 1/2"	4
	90-1778C	Traction Sleeve Set - 4" x 3/4"	4
	64-1036C	Traction Sleeve Set - 8" x 1/2"	4
	64-1037C	Traction Sleeve Set - 8" x 3/4"	4
2	QE14399-T	Label - Manufacturer	1
3	VB15A1O24-S	Hex Bolt	4
4	VW10A10	Washer	4

# Recycling Passport

OL280005 (Reference Only)

English


Item	Recycling/Material Code	Important Information	Qty
A	QDF28-0265 Circuit Board		1
B	QDF5095 Buzzer		1

[Return To Table of Contents](#)

# Recycling Passport

OL280120 (Reference Only)


Item	Recycling/Material Code	Important Information	Qty
A	28-0769L Motor		1
B	28-0770 Motor		1
C	25-0768 Motor		2
D	QDF23-0315 Circuit Board		1

# Recycling Passport

OL280121 (Reference Only)

English


Item	Recycling/Material Code	Important Information	Qty
A	28-0769L Motor		1
B	28-0770 Motor		1
C	25-0768 Motor		2
D	QDF23-0315 Circuit Board		1
E	QDF14-1160 Transformer		1

[Return To Table of Contents](#)

# Recycling Passport

OL280097, OL280098, OL280100 (Reference Only)


Item	Recycling/Material Code	Important Information	Qty
A	QDF20-0126 Circuit Board		1

# Recycling Passport

OL280035, OL280101, OL280102 (Reference Only)

English

Item	Recycling/Material Code	Important Information	Qty
A	QDF20-0121 Circuit Board		1

[Return To Table of Contents](#)

# EMC Information

## GOBED® II MEDSURG BED, MODEL FL28C

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
The GoBed® II MedSurg Bed, Model FL28C is suitable for use in the electromagnetic environment specified below. The customer or the user of the GoBed® II MedSurg Bed, Model FL28C should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast Transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±8 kV differential mode ±2 kV common mode	±8 kV differential mode ±2 kV common mode	Main power quality is that of a typical commercial and/or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	<5%Ut (95% dip in Ut) for 0.5 cycle 40%Ut (60% dip in Ut) for 5 cycles 70%Ut (30% dip in Ut) for 25 cycles. <5% Ut (>95% dip in Ut) for 5 sec.	<5%Ut (95% dip in Ut) for 0.5 cycle 40%Ut (60% dip in Ut) for 5 cycles 70%Ut (30% dip in Ut) for 25 cycles. <5% Ut (>95% dip in Ut) for 5 sec.	Main power quality should be that of a typical commercial and/or hospital environment. If the user of the GoBed® II MedSurg Bed, Model FL28C requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.
Note: $U_T$ is the a.c. mains voltage prior to applications of the test level.			

# EMC Information

## GOBED® II MEDSURG BED, MODEL FL28C (CONTINUED)

Recommended separation distances between portable and mobile RF communications equipment and the GoBed® II MedSurg Bed, Model FL28C.			
The GoBed® II MedSurg Bed, Model FL28C is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the GoBed® II MedSurg Bed, Model FL28C can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the GoBed® II MedSurg Bed, Model FL28C as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	m		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	1.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
<b>Note 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
<b>Note 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

# EMC Information

## GOBED® II MEDSURG BED, MODEL FL28C (CONTINUED)

The GoBed® II MedSurg Bed, Model FL28C is suited for use in the electromagnetic environment specified below. The customer or the user of the GoBed® II MedSurg Bed, Model FL28C should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the GoBed® II MedSurg Bed, Model FL28C, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended Separation Distance</p> $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P}$ $d=2.3 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

**Note 1**

At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2**

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the GoBed® II MedSurg Bed, Model FL28C is used exceeds the applicable RF compliance level above, the GoBed® II MedSurg Bed, Model FL28C should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the GoBed® II MedSurg Bed, Model FL28C.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

# EMC Information

## GOBED® II MEDSURG BED, MODEL FL28C (CONTINUED)

<b>Guidance and Manufacturer's declaration - Electromagnetic Emissions</b>		
<p>The GoBed® II MedSurg Bed, Model FL28C is intended for use in an electromagnetic environment specified below. The customer or the user of the GoBed® II MedSurg Bed, Model FL28C should assure that it is used in such an environment.</p>		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment</b>
RF Emissions CISPR 11	Group 1	The GoBed® II MedSurg Bed, Model FL28C uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The GoBed® II MedSurg Bed, Model FL28C is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

# Warranty

---

## LIMITED WARRANTY

Stryker Medical Division, a division of Stryker Corporation, warrants to the original purchaser the GoBed® II MedSurg Bed, Model FL28C to be free from defects in material and workmanship for a period of One (1) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product which is, in the sole discretion of Stryker, found to be defective. If requested by Stryker, products or parts for which a warranty claim is made shall be returned prepaid to the factory. Any improper use or any alteration or repair by others in such manner as in Stryker's judgment affects the product materially and adversely shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

Stryker Medical Bed products are designed for a 10 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its Bed products will be free from structural defects for the expected 10 year life of the Bed product as long as the original purchaser owns the product.

This statement constitutes Stryker's entire warranty with respect to the aforesaid equipment. **Stryker makes no other warranty or representation, either expressed or implied, except as set forth herein. There is no warranty of merchantability and there are no warranties of fitness for any particular purpose. In no event shall Stryker be liable here under for incidental or consequential damages arising from or in any manner related to sales or use of any such equipment.**

## TO OBTAIN PARTS AND SERVICE

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative, or call Stryker Customer Service USA at 1-800-327-0770, Canada 1-888-233-6888.

## SERVICE CONTRACT COVERAGE

Stryker has developed a comprehensive program of service contract options designed to keep your equipment operating at peak performance at the same time it eliminates unexpected costs. We recommend that these programs be activated before the expiration of the new product warranty to eliminate the potential of additional equipment upgrade charges.

### A Service Contract helps to:

- Ensure equipment reliability
- Stabilize maintenance budgets
- Diminish downtime
- Establish documentation for JCAHO
- Increase product life
- Enhance trade-in value
- Address risk management and safety

# Warranty

## SERVICE CONTRACT PROGRAMS

Stryker offers the following service contract programs:

Service Agreement Options *	Premium	Complete	Standard
Annually scheduled preventative maintenance	X		X
All parts	X	X	
All labor and travel	X	X	
Unlimited emergency service calls	X	X	
Priority one contact: two hour phone response	X	X	
Most repairs completed within 3 days	X	X	
JCAHO documentation	X	X	X
On-site record of PM & emergency service	X		X
Factory-trained Stryker service technician	X	X	X
Stryker authorized parts used	X	X	X
Service during regular business hours (8–5)	X	X	X

\* Does not include maintenance due to abuse or for any disposable items. Stryker reserves the right to change options without notice.

Stryker Medical also offers personalized service contracts.  
Pricing is determined by age, location, model and condition of product.

**For more information on our service contracts,  
please call your local representative.**

## RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. **Special, modified, or discontinued, items not subject to return.**

## DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full. Claims for any short shipment must be made within thirty (30) days of invoice.

## INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

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