Pioneer®
Pressure Redistribution Surface

REF 0850

stryker®

Operations/Maintenance Manual

CE

Table of Contents

	<u>Page</u>
Symbols and Definitions	
Symbols	3
Warning / Caution / Note Definitions	4
Introduction	
Intended Use of Product	5
Product Description	5
Contraindications.....	5
Specifications	6
Contact Information	7
Summary of Safety Precautions	8-9
Installation Guide	10
Operation Guide	11-12
Cleaning and Disinfection	
Support Surface Cleaning	13
Support Surface Disinfection	14
Service Information	
Cover Replacement	15
Fire Barrier Replacement	15
Preventive Maintenance	
Checklist	16
Product Labels	17-18
Quick Reference Replacement Parts List	19
Warranty	
Limited Warranty	20
To Obtain Parts and Service	20
Return Authorization	20
Damaged Merchandise	20

Symbols and Definitions

SYMBOLS

	Warning - Consult accompanying documentation
	Caution
	Safe Working Load
	Consult Instructions for Use
	Manufacturer
	Date of Manufacture
	Model or catalog number
	Type B
	Authorized European Representative in the European Community

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

NOTE! Provides special information to make specific instructions easier to understand.

Introduction

This manual is designed to assist with the operation and maintenance of the Pioneer Non-Powered Support Surface. Carefully read this manual thoroughly before using or beginning maintenance on the support surface. To ensure safe operation of this equipment, it is recommended that methods and procedures are established for educating and training staff on the safe operation of the support surface.

INTENDED USE OF THE PRODUCT

Pioneer Support Surface is a non-powered pressure relieving stretcher pad that incorporates the combination of air and foam technologies. This support surface is intended to assist in improving human patient outcomes, focusing on equalizing pressure redistribution and enhancing comfort. Pioneer support surface is used to assist in the prevention of pressure wounds.

This product is intended to be used:

- With patients in the Emergency Department (ED).
- In acute care, general hospital care, or other locations as prescribed by a physician
- For multiple patients after appropriate cleaning/ disinfection, i.e., reusable; this product is not disposable.
- On Stryker and other stretchers (See "**Specifications**" section on **Page 6**)

This product is not intended to be used in a home health care environment.

This product is not designed to be waterproof. The excessive exposure to fluid could cause fluid intrusion through the cover into the foam core. Refer to **Page 13 and 14** for Cleaning and Disinfection guide.

The Pioneer has a safe working load maximum of 700 pounds and therapeutic weight maximum of 290 pounds. The patient must not exceed safe working load specified by the support surface, frame, and accessories. The use of this support surface for patients outside the therapeutic weight maximum should be determined by a health care professional. Patients shall meet the minimum age requirement of 2 years old.

This support surface is not intended to support a patient in a prone position, for example laying face downward. Pioneer shall be used with a mattress cover at all times.

The Pioneer is not intended to be a sterile product, nor is it intended to include a measuring function. Pioneer has an expected product life of three (3) years for the surface and three (3) years for the cover.

PRODUCT DESCRIPTION

Pioneer is a non-powered support surface focusing on equalizing pressure and enhancing comfort. It features interface-regulated air adjustment to help balance pressure redistribution and management of pressure loading in high risk zones of the body.

CONTRAINDICATIONS

None known.

Introduction

SPECIFICATIONS

	26 Inch Mattress		30 Inch Mattress	
Model Number	0850-026-000		0850-030-000	
Overall Width	26 inches	66 cm	30 inches	76 cm
Overall Length	76 inches	193 cm	76 inches	193 cm
Overall Thickness	5 inches	12.7 cm	5 inches	12.7 cm
Weight Capacity	700 lbs	318 kg	700 lbs	318 kg
Product Weight	16 lbs	7.3 kg	18 lbs	8.2 kg

NOTE! Stryker reserves the right to change specifications without notice.

Contact Information

Contact Stryker Customer Service or Technical Support at:

1-800-327-0770 or (269) 324-6500

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Summary of Safety Precautions

Carefully read and strictly adhere to the following warnings and cautions. Service of the product should only be completed by qualified personnel.

WARNING

- Inspect Pioneer support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a mattress cover becomes compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - To avoid the risk of injury or equipment damage, do not transfer the patient from one platform to another using the support surface.
 - It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance.
 - Check patient's skin regularly. Consult physician if any skin redness or breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
 - Disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
 - **Some disinfectants are corrosive in nature and may cause damage to the product if used improperly.** Only use the "Suggested Disinfectants" described on page 14 according to the manufacturer's instructions.
 - The support surface cover must be completely DRY before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material. Failure to follow recommended cleaning/disinfection procedures may void the product warranty.
 - Internal components of this support surface are not washable. If there is contamination inside the support surface, discard the support surface.
 - When cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.
 - Do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
 - Patient stability and side rail coverage may be compromised with the use of the Pioneer. When using an overlay, extra caution and/or supervision is required to help reduce the likelihood of a patient falling.
 - Leave the frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
 - The safe use of the support surface is maximized when used in conjunction with side rails; there may be an increased risk of falls when side rails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. Local policies regarding the use of side rails should be taken into account. Whether and how to use side rails is a decision that should be based on each patient's individual needs and should be made by the physician, operators, and responsible parties.
 - To reduce risk of patient or operator injury and support surface damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.
 - The risk of entrapment can develop when the support surface is placed on a stretcher frame that leave gaps of even a few inches between the support surface and the headboard, footboard and/or side rails. Using a Pioneer support surface on a stretcher frame that has a larger width than the support surface may result in sliding and patient injury (see Specifications section on page 6).
 - Extra precautions are needed for a patient at risk of a fall (i.e. agitated or confused) to help reduce the likelihood of a fall.
 - When placing a patient on the support surface, ensure the opposite side rail is raised to reduce the risk of patient fall.
-

Summary of Safety Precautions

CAUTION

- Be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover. Doing so may reduce pressure redistribution performance.
 - It is the responsibility of the caregiver team to evaluate the appropriate CPR protocol to be used with the surface.
 - The use of quaternaries containing glycol ethers and/or accelerated hydrogen peroxides may result in compromised legibility of the silk screen of the cover.
 - Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
-

NOTE! The Pioneer surface must be used with a mattress cover at all times.

Installation Guide

SUPPLIES REQUIRED:

- Rubbing or isopropyl alcohol
- Velcro® strips (mating loop Velcro® included in the package)

SECURING PIONEER MATTRESS TO STRETCHER:

1. Mate loose piece of Velcro® to the elongated piece of Velcro® that goes down the center of the mattress.
 2. Position mattress on stretcher and observe location where Velcro® will reside.
 3. Clean the area where Velcro® will reside with rubbing alcohol or isopropyl alcohol.
 4. Allow to dry.
 5. Expose adhesive backing and position mattress onto stretcher.
 6. Apply a load of 50 pounds to the area of the mattress above the Velcro® for one minute.
-

WARNING

- The risk of entrapment can develop when the support surface is placed on a stretcher frame that leave gaps of even a few inches between the support surface and the headboard, footboard and/or side rails. Using a Pioneer support surface on a stretcher frame that has a larger width than the support surface may result in sliding and patient injury (see Specifications section on **Page 6**).
-

Operation Guide

TRANSFER OF PATIENT FROM ONE PATIENT SUPPORT PLATFORM TO ANOTHER

WARNING

- Leave the frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
- Patient stability and side rail coverage may be compromised with the use of an overlay. When using an overlay, extra caution and/or supervision is required to help reduce the likelihood of a patient falling.
- The safe use of the support surface is maximized when used in conjunction with side rails; there may be an increased risk of falls when side rails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. Local policies regarding the use of side rails should be taken into account. Whether and how to use side rails is a decision that should be based on each patient's individual needs and should be made by the physician, operators, and responsible parties.
- When placing a patient on the support surface, ensure the opposite side rail is raised to reduce the risk of patient fall.
- The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard, and side rails. Do not use the support surface when such gaps are present.
- Using a Pioneer support surface on a stretcher frame that has a larger width than the support surface may result in sliding and patient injury (see Specifications section on page 6).
- Extra precautions are needed for a patient at risk of a fall (i.e., agitated or confused) to help reduce the likelihood of a fall.
- To reduce risk of patient or user injury and support surface damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.
- To avoid the risk of injury or equipment damage, do not transfer the patient from one platform to another using the support surface.

PROCEDURE

To transfer the patient from one patient support platform (for example, stretcher, gurney, operating table) to another:

1. Position one patient support platform alongside the other patient support platform as closely as possible.
2. Set the brakes to "ON" for both patient support platforms.
3. Prior to transfer of patient between fixed patient support surfaces, ensure that the surfaces are level with one another.
4. Raise the patient support platform side rail located opposite the patient transfer.
5. Transfer the patient following all applicable safety rules and institution protocols to ensure patient and operator safety.

CAUTION

- Be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover. Doing so may reduce pressure redistribution performance.
 - It is the responsibility of the caregiver team to evaluate the appropriate CPR protocol to be used with the surface.
-

Operation Guide

USING INCONTINENCE MANAGEMENT DEVICES

This support surface is not intended to manage incontinence. Therefore, it is recommended to use incontinence management devices when appropriate. Disposable diapers or incontinence pads may be used.

WARNING

- It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance.
 - Check patient's skin regularly. Consult physician if any skin redness or breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
-

Cleaning

SUPPORT SURFACE CLEANING

 **WARNING**

- To avoid the risk of equipment damage, do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
 - Inspect Pioneer support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a mattress cover becomes compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - Internal components of this support surface are not washable. If there is contamination inside the support surface, discard the support surface.
-

Stryker support surfaces are designed for long-lasting comfort and reliability. The life of the support surface can be adversely affected by an increase in frequency of usage which might include more frequent cleaning and disinfection.

1. Using a clean, soft, damp cloth, wipe down the entire support surface with a mild soap and water solution to remove foreign material.
 2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or cleaning agent.
 3. Care must be taken to thoroughly rinse and dry covers following cleaning.
 4. Disinfect as necessary with a hospital grade disinfectant AFTER cleaning has been completed. Refer to following “Support Surface Disinfection” section.
-

 **CAUTION**

To avoid the risk of equipment damage, when cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface

Disinfection

SUPPORT SURFACE DISINFECTION

WARNING

- To avoid risk of patient injury, disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
 - Some disinfectants may cause damage to the product if used improperly. If the products described below are used to disinfect the support surface, measures must be taken to ensure the entire surface is wiped with a damp cloth soaked in clean water and thoroughly dried following disinfection. The cover can be damaged when exposed to such disinfectants beyond the manufacturers' recommendations. Failure to follow these directions when using these types of disinfectants may void this product warranty.
 - To avoid the risk of equipment damage, the support surface cover must be completely dry before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.
-

SUGGESTED DISINFECTANTS

- Quaternaries are acceptable if $\text{pH} < 10$. DO NOT USE quaternaries containing sodium hydroxide (NaOH)
 - Chlorinated Bleach Solution (5.25% bleach. For daily use, dilute to 1 part bleach to 1,000 part water. For decontamination, dilute to a maximum of 1 part bleach to 100 parts water).
 - 70% Isopropyl Alcohol is acceptable as a disinfectant but pooling of the cleaning solution must be avoided and care must be taken until fully dry.
1. Ensure surface has been thoroughly cleaned and dried prior to applying disinfectants.
 2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or disinfectant.
 3. Care must be taken to thoroughly **rinse and dry** covers following disinfection.

NOTE! After servicing internal component, or replacing the Fire Sock, ensure the reinstallation of the Fire Sock encloses all internal components and the extra length of material is tucked completely under the foot end.

CAUTION

- Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
 - The use of accelerated hydrogen peroxides or quaternaries containing glycol ethers may damage the cover.
 - DO NOT USE Phenolic cleaning agents as it result in irreversible swelling of coating.
-

RECOMMENDATION: The Safety Communication lists several recommendations for inspection and maintenance including:

- Regularly check each medical bed mattress cover for any visible signs of damage or wear such as cuts, tears, cracks, pinholes, snags or stains.
- Routinely remove the medical bed mattress cover and check its inside surface. Once the mattress cover is removed, inspect the mattress for wet spots, staining, or signs of damage or wear. Check all sides and the bottom of the mattress.
- Immediately replace any medical bed mattress cover with visible signs of damage or wear to reduce the risk of infection to patients.
- DO NOT stick needles into a medical bed mattress through the mattress cover.

Service Information

COVER REPLACEMENT

Tools Required:

- None

PROCEDURE

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient's right side and then remove the foam crib assembly from the bed and set aside.
5. Discard the old cover.
6. Place the new cover, unzipped and open, on the bed with the black bottom cover on the litter and the top cover folded over the patient's right side of the bed.
7. Carefully place the foam assembly on top of the bottom part of the cover so it aligns with the cover.

NOTE! Maintain the correct "Head-Foot" orientation between the new cover and the foam assembly.

8. Fold the top cover over the top of the foam crib assembly so it aligns with the foam crib assembly.
9. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
10. Verify proper operation of the unit before returning it to service.

FIRE BARRIER REPLACEMENT

Tools Required:

- None

PROCEDURE

1. Raise the stretcher height to the full up position.
2. Lower the fowler and gatch sections to the full down positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient's right side.
5. Starting at the foot end, roll up the fire barrier on the foam crib assembly

NOTE! Work from side to side, a little at a time, to move the fire barrier to the top of the support surface.

6. Discard the old fire barrier.
7. Starting at the head end, roll the new fire barrier down and slide over the foam crib assembly.

NOTE! Align the fire barrier on the foam crib before sliding over the foam crib assembly.

8. Carefully slide the fire barrier down the foam crib assembly, working from side to side, to make sure that the fire barrier is tight on the foam crib assembly.
9. Align the foam crib assembly on top of the bottom part of the cover.

NOTE! Spread the excess fire barrier material equally below the foam crib assembly at the foot end.

10. Fold and align the top cover over the top of the foam crib assembly.
11. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
12. Verify proper operation before returning the product to service.

Preventive Maintenance

Preventive maintenance should be performed annually, at a minimum. A preventive maintenance program should be established for all Stryker Medical equipment. Preventive maintenance may need to be performed more frequently based on the usage level of the product

The following or a similar checklist should be completed before using the surface on a new patient and/or after servicing.

CHECKLIST

Product Serial Number:		OK ✓
1. Cover zipper opens and closes properly and has no visible damage		
2. No tears, rips, holes, cracks, or other openings in the mattress cover		
3. Labels are legible, and securely attached (see page 17-18)		
4. No rips or cracks in the stitching		
5. Foam and other components have not degraded, or come apart		
6. Fire barrier has no visible signs of damage; if excessive wear is observed, replacement of fire barrier is recommended		
7. Individual Air Bladders should hold air and show no cuts or leaks		
8. Bladder tubing and fittings are tightly connected and no cuts or leaks are present		
9. Velcro® strips are present and securely attached		
10. Surface is operating properly		
Complete by:		Date:

Product Labeling

All labels (shown below) are attached to the Pioneer support surface (26 inch x 76 inch):

UNDER PENALTY OF LAW THIS TAG NOT TO BE REMOVED EXCEPT BY THE CONSUMER.

ALL NEW MATERIAL
 CONSISTING OF:
POLYURETHANE FOAM

FINISHED SIZE: 5.0"x 26.0"x 76.0"	NET WT. OF FILLING MAT'L: 16 LBS.
---	--

<p>THIS ARTICLE MEETS THE FLAMABILITY REQUIREMENTS OF CALIFORNIA BUREAU OF HOME FURNISHINGS TECHNICAL BULLETIN NO 129. CARE SHOULD BE EXERCISED NEAR OPEN FLAME OR WITH BURNING CIGARETTES.</p>	<p>CONFORMS WITH FLAMMABILITY STANDARD 16 CFR 1633, 1632 REGISTRY NO. CA-35016 (CAL) CERTIFICATION IS MADE BY THE MANUFACTURER THAT THE MATERIALS IN THIS ARTICLE ARE DESCRIBED IN ACCORDANCE WITH THE LAW.</p>
---	---

stryker[®]
 3800 E. CENTRE AVE.
 PORTAGE, MI 49002-5826

CE  **MAX. 700 LBS (319 KG)**

MATTRESS 0850-026-000
 DATE OF MANUFACTURE: _____

Manufactured by:
AMF SUPPORT SURFACES, INC.
 9901 KINCAID DRIVE
 FISHERS, IND. 46037

Date of Manufacture:

Model:
Pioneer[®]

Prototype ID:
AST-IPP-AA-2676

This mattress meets the requirements of 16 CFR Part 1633 (federal flamability (open flame) standard for mattress sets) When used without a foundation.

THIS MATTRESS IS INTENDED TO BE USED WITHOUT A FOUNDATION

Label Not Required for Europe, Only US

Product Labeling

All labels (shown below) are attached to the Pioneer support surface (30 inch x 76 inch):

UNDER PENALTY OF LAW THIS TAG NOT TO BE REMOVED EXCEPT BY THE CONSUMER.

ALL NEW MATERIAL
CONSISTING OF:
POLYURETHANE FOAM

FINISHED SIZE: 5.0"x 30.0"x 76.0"	NET WT. OF FILLING MAT'L. 18 LBS.
---	--

THIS ARTICLE MEETS THE FLAMABILITY REQUIREMENTS OF CALIFORNIA BUREAU OF HOME FURNISHINGS TECHNICAL BULLETIN NO 129. CARE SHOULD BE EXERCISED NEAR OPEN FLAME OR WITH BURNING CIGARETTES.	CONFORMS WITH FLAMMABILITY STANDARD 16 CFR 1633, 1632 REGISTRY NO. CA-35016 (CAL) CERTIFICATION IS MADE BY THE MANUFACTURER THAT THE MATERIALS IN THIS ARTICLE ARE DESCRIBED IN ACCORDANCE WITH THE LAW.
--	--

stryker®
3800 E. CENTRE AVE.
PORTAGE, MI 49002-5826

CE  **MAX. 700 LBS (319 KG)**

MATTRESS 0850-030-000
DATE OF MANUFACTURE: _____

Manufactured by:
AMF SUPPORT SURFACES, INC.
9901 KINCAID DRIVE
FISHERS, IND. 46037

Date of Manufacture:

Model: **Pioneer®**

Prototype ID:
AST-IPP-AA-3076

This mattress meets the requirements of 16 CFR Part 1633 (federal flammability (open flame) standard for mattress sets) When used without a foundation.

THIS MATTRESS IS INTENDED TO BE USED WITHOUT A FOUNDATION

Label Not Required for Europe, Only US

Quick Reference Replacement Parts List

The parts and accessories listed on this page are currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call *Stryker Customer Service USA* at **1-800-327-0770** for availability and pricing.

Part Number	Part Name	Quantity
0850-026-051	26 inch Bladder Kit	1
0850-030-051	30 inch Bladder Kit	1
0850-030-052	Tubing Kit	1
0850-026-053	26 inch Cover Kit, Pioneer	1
0850-030-053	30 inch Cover Kit, Pioneer	1
0850-026-054	Fire Sock, Pioneer 26 inch	1
0850-030-054	Fire Sock, Pioneer 30 inch	1
0850-030-055	Velcro Kit, Pioneer	1

Warranty

Stryker Medical, a division of Stryker Corporation ("Stryker") warrants that its Model 0850 Pressure Redistribution Surface products will be free from defects in manufacturing and workmanship for a period of two (2) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying a product replacement, at its option, any product which is, in the sole discretion of Stryker, found to be defective. If requested by Stryker, products for which a claim is made shall be returned prepaid to the factory. Any improper use or any alterations or repair by others in such manner as in Stryker's judgment affects the product materially and adversely shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

The above noted warranty periods apply only to the original purchaser of the Pioneer and begin on the date of delivery to such original purchaser.

Warranty Exclusion and Damage Limitations

The express warranty set forth herein is the only warranty applicable to the product. **Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by Stryker.** In no event shall Stryker be liable for incidental or consequential damages.

To Obtain Parts and Service

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative or call Stryker Customer Service at **1-800-327-0770** or **(269) 324-6500**.

Return Authorization

Product cannot be returned without prior approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged Product

ICC Regulations require that claims for damaged product must be made within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

International Warranty Clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Contact your local Stryker Medical representative for additional information.

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**UNDER PENALTY OF LAW THIS TAG NOT TO
BE REMOVED EXCEPT BY THE CONSUMER.**

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FINISHED SIZE:
5.0"x 26.0"x 76.0"

NET WT. OF FILLING
MAT'L. 16 LBS.

THIS ARTICLE
MEETS THE
FLAMABILITY
REQUIREMENTS
OF CALIFORNIA
BUREAU OF HOME
FURNISHINGS
TECHNICAL
BULLETIN NO
129. CARE SHOULD
BE EXERCISED NEAR
OPEN FLAME OR
WITH BURNING
CIGARETTES.

CONFORMS WITH
FLAMMABILITY
STANDARD
16 CFR 1633, 1632
REGISTRY NO.
CA-35016 (CAL)
CERTIFICATION
IS MADE BY THE
MANUFACTURER
THAT THE MATERIALS
IN THIS ARTICLE ARE
DESCRIBED IN
ACCORDANCE WITH
THE LAW.

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Prototype ID:

AST-IPP-AA-2676

This mattress meets the requirements of
16 CFR Part 1633 (federal flammability
(open flame) standard for mattress sets)
When used without a foundation.

**THIS MATTRESS
IS INTENDED TO BE USED
WITHOUT A FOUNDATION**

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CE  **MAX. 700 LBS (319 KG)**

MATTRESS 0850-030-000
 DATE OF MANUFACTURE: _____

Manufactured by:
 AMF SUPPORT SURFACES, INC.
 9901 KINCAID DRIVE
 FISHERS, IND. 46037

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Hurtig referenceliste over reservedele

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Stryker France S.A.S.