IsoTour Pump

Operations Manual

Standard and premium

REF 2874

EN
FR

Global symbol glossary

See the Global Symbol Glossary at ifu.stryker.com for symbol definitions.

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR), disconnection point
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Catalogue number
	Serial number
	Batch code
	European medical device
	Unique device identifier
	Quantity
	Authorized representative in Switzerland
	Authorized representative in the European Community
	For US Patents see www.stryker.com/patents
	Manufacturer
	Date of manufacture
	CE mark

	UK Conformity Assessment mark
	Importer
	Mass of product
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provisions for protective earthing or reliance upon installation conditions.
	Type B applied part - 2872 IsoTour support surface
	Fuse
	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1:2005, A1:2012, and A2:2021, CAN/CSA-C22.2 No. 60601-1:14 and A2:2022 (R2022), AAMI/IEC 60601-1-8 (2006) + AMD 1 (2013) + AMD 2 (2021), and CAN/CSA C22.2 No. 60601-1-8:23.
IP21	Solids: Protection from touch by fingers and objects greater than 12 mm Liquids: Protection from the fall of vertical water drops
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Indications for use	4
Contraindications	4
Expected service life	5
Disposal/recycle	5
Specifications	5
Environmental conditions	6
Product illustration	6
Button and icon descriptions	7
Features	8
Contact information	8
Serial number location	9
Date of manufacture	9
Setup	10
Connecting the power cord	10
Connecting the hose	11
Operation	12
Powering on the pump	12
Activating CPR release	13
Powering on or off Low Air Loss (LAL)	13
Selecting TruTurn	13
Increasing or decreasing TruTurn	14
Pausing or canceling TruTurn	15
Locking or unlocking product functions	15
Powering off the pump	16
Accessories and parts	17
Troubleshooting	18
Product alarms	20
Preventive maintenance	21
Cleaning and disinfecting with SideKick	22
Cleaning	23
Disinfecting	24
EMC information	25

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always consult the operations manual for the support surface before you attach the pump.
- Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
- Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
- Always set up and use in accordance with the instructions to avoid harmful interference to other devices in the vicinity. This product radiates radio frequency energy.
- Always make sure the hose is routed down and away from the patient. Avoid accidental patient contact at the foot end of the support surface. Hose temperature may reach 48° C (118.4° F).
- Never transfer or transport a patient while TruTurn is in use.
- Always raise the bed siderails before you begin TruTurn to avoid the risk of a patient fall.
- Do not leave the patient unattended until TruTurn of the patient position is complete.
- Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.
- Always unplug the product before you clean or disinfect.
- Do not immerse the product in cleaning or disinfectant solutions.
- Do not spray or allow liquids to enter the power cord ports or the hose ports.
- Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of IsoTour, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Always allow the product to come to room temperature after exposure to extreme high or low temperatures.

- Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
 - Do not attach the power cord to or route the power cord through any moving parts of the bed frame to avoid the risk of electric shock, pinching the power cord, or interfering with bed functions.
 - Always plug the pump power cord into the bed auxiliary outlet or a grounded power receptacle. Do not use an extension cord.
 - Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
 - This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 2874 IsoTour Pump is an electric pump for use with the Stryker Model 2872 IsoTour Support Surface. The pump connects to the support surface with an air hose. Air is delivered through the integrated perforated tubing system within the support surface.

The TruTurn function is available when the pump and premium support surface are connected. The pump inflates the integrated turn bladders in the bottom cover. When the pump is turned on, the Low Air Loss (LAL) function is present for the standard or premium support surface.

Intended use

The IsoTour pump when used with the IsoTour gel support surface is also intended to help in managing microclimate via Low Air Loss technology. Premium IsoTour surfaces are TruTurn capable. TruTurn is intended to assist in pressure redistribution and patient repositioning.

Indications for use

The IsoTour pump when used with the IsoTour gel support surface is also intended to help in managing microclimate via Low Air Loss technology. Premium IsoTour surfaces are TruTurn capable. TruTurn is intended to assist in pressure redistribution and patient repositioning.

Contraindications

The TruTurn feature is contraindicated for use with patients that have unstable spine.

Note - TruTurn requires the IsoTour pump to be connected for use. The TruTurn feature is visually identified by the icons on the IsoTour pump (Figure 1).

Figure 1 – TruTurn icon

Expected service life

The IsoTour Pump has a five year expected service life under normal use conditions and with appropriate periodic maintenance.

The air hose has a two year expected service life under normal use conditions and with appropriate periodic maintenance.







Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications


Height		17.25 in.	43.8 cm
Width		9.25 in.	23.5 cm
Depth		7.25 in.	18.4 cm
Weight		16.3 lb	7.4 kg
Interface control panel		Liquid crystal display (LCD)	
Input voltage		100-240 V	
Current		0.7 amps	
Input frequency		50/60 Hz	
Noise level		<50 dba	
Mode of operation		Continuous	
Power cord	Short	3.2 ft	1 m
	Long	16.4 ft	5 m
Compliance		IEC 60601-1 3.2 Edition; RoHS Directive 2002/95/EC Reach; IEC 60601-1-8:2020, IEC 60601-1-2:2014 +AMD1:2020.	
Compatible with bed frame	IsoTour standard and premium support surface with an IsoTour pump	3002 S3, 3005 S3, InTouch, ProCeed, ProCuity, Spirit Select	

Environmental conditions

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

Product illustration


A	Graphical user interface
B	Status indicator
C	Bumper
D	Transport handle
E	Hook
F	Air filter
G	Back pump housing

H	Power inlet
I	Anti-vibration foot
J	Fuse
K	Flex grip
L	Hose port
M	Hose
N	Front pump housing

Button and icon descriptions

The buttons are located on the outside of the graphical user interface.

Icon	Name	Description
	Power button	See <i>Powering on the pump</i> (page 12) or <i>Powering off the pump</i> (page 16)
	Cancel icon	Cancels the current function, returns to the main menu
	Pause icon	Pauses the current function
	Action button, right	Allows the operator to select an action on the screen
	Action button, left	Allows the operator to select an action on the screen
	Lock or unlock icons	See <i>Locking or unlocking product functions</i> (page 15)
	Back icon	Returns to the previous screen
	Low Air Loss (LAL) icon	Low Air Loss function active. On by default. See <i>Powering on or off Low Air Loss (LAL)</i> (page 13)
	Fill default segment 2 icon	Allows the operator to turn the patient to the right or to the left
	Maximum fill icon	Maximum fill is three bars. See <i>Increasing or decreasing TruTurn</i> (page 14)
	Decrease or increase icons	Decreases or increases patient turn
	Check siderails animated icon	Reminder to the operator to check siderails before they activate TruTurn
	Check air connector alert icon	Product detects air connector error
	Kinked air hose alert icon	Product detects a bend or twist in the air hose
	Alarm off icon	Closes and resets pump for kinked hose or check connector

Features

TruTurn	Available for the premium support surface. Air bladders inflate to rotate up to 30 degrees. Helps to off-load the sacrum. Turn angle can be customized per patient tolerance and comfort.
Moisture management	The top surface cover is designed to draw excess moisture away from the patient's skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to help manage microclimate of the skin.
Standard support surface	Orange bottom cover, gel and foam core, LAL capable when connected to pump, handles, and D-rings
Premium support surface	Blue bottom cover, gel and foam core, spacer fabric for enhanced air flow, handles, D-rings, anchor straps, TruTurn and LAL capable when connected to pump

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location


Figure 2 – Serial number location

Date of manufacture

The year of manufacture is the first four digits of the serial number.

Setup

Unpack the cartons and check all items. Make sure that the product operates before you place the product into service.

WARNING

- Always consult the operations manual for the support surface before you attach the pump.
 - Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
 - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
 - Always set up and use in accordance with the instructions to avoid harmful interference to other devices in the vicinity. This product radiates radio frequency energy.
-

CAUTION

- Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
 - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
 - Do not attach the power cord to or route the power cord through any moving parts of the bed frame to avoid the risk of electric shock, pinching the power cord, or interfering with bed functions.
 - Always plug the pump power cord into the bed auxiliary outlet or a grounded power receptacle. Do not use an extension cord.
-

Connecting the power cord

To connect the power cord:

1. Bring the product to room temperature.
2. Attach the power cord to the pump.
3. Insert the power cord under the power cord retention flex grip (A) (Figure 3).


Figure 3 – Cord retention flex grip

4. Plug the pump power cord into an accessible wall outlet.
5. Place the power cord out of the path of foot traffic to avoid trip hazards.
6. Place the pump in a secure location. Allow space around the pump to connect and disconnect the hoses and power cord.
 - a. Option 1: Hang the pump on a compatible bed. Make sure that the pump is secure.

- b. Option 2: Place the pump on the floor.

Connecting the hose

WARNING - Always make sure the hose is routed down and away from the patient. Avoid accidental patient contact at the foot end of the support surface. Hose temperature may reach 48° C (118.4° F).

To connect the hose:

1. Attach the hose to the pump hose ports.
2. Open the hose port (A) on the support surface (Figure 4).


Figure 4 – Hose connection

3. Attach the hose to the support surface hose port. Make sure that the rib (B) lines up with the red line on the hose (C).

Operation

Powering on the pump

To power on the pump:

1. See *Setup* (page 10) to prepare the pump for operation.
2. Press the **Power** button once. The **Startup** screen appears (Figure 5).


Figure 5 – Startup screen - operator view

Note - The pump detects either the standard or premium support surface.

3. The pump, on startup, pumps air into the left bladder to determine which compatible mattress is connected.
4. If the pump detects that there are no turn bladders, the pump shows the standard software support surface display (Figure 6).

Figure 6 – Standard software support surface display

Note - TruTurn is available for the premium support surface only.

5. If the pump detects that there are turn bladders, the pump shows the premium software support surface display (Figure 7).

Figure 7 – Premium software support surface display

Note - TruTurn is available for the premium support surface only.

Activating CPR release

Remove the hose from the pump to deflate the support surface for CPR. Follow your hospital protocol for CPR.

Powering on or off Low Air Loss (LAL)

LAL is on by default. You can turn this function on or off.

To turn LAL on or off:

1. Press the **Action** button next to the active (green) **LAL** icon to turn off this function (Figure 8).

Figure 8 – LAL on

2. Press the **Action** button next to the inactive (white) **LAL** icon to turn on this function (Figure 9).

Figure 9 – LAL off

Selecting TruTurn

WARNING

- Never transfer or transport a patient while TruTurn is in use.
 - Always raise the bed siderails before you begin TruTurn to avoid the risk of a patient fall.
 - Do not leave the patient unattended until TruTurn of the patient position is complete.
-

To select TruTurn:

1. Read the IsoTour Gel Support Surface Operations Manual before you start TruTurn.
2. Make sure that the head of bed is at zero degrees before you start TruTurn.
3. Press the **Action** button next to the patient left or patient right **TruTurn** icon (Figure 10).

Figure 10 – TruTurn

4. Raise the siderails to the highest height position (Figure 11).


Figure 11 – Check siderail

5. Press the **Action** button next to the time to set the turn duration.

- 02:00:00 = 2 hours
- 00:30:00 = 30 minutes
- 00:00:15 = 15 seconds

Note

- Before you transfer or transport the patient, or complete other patient handling activity, cancel the TruTurn feature. You can either cancel the TruTurn feature, turn off the pump, or pull the CPR release.
- The support surface will turn up to 30 degrees.
- The countdown for the time selected will start after the pump has completed the fill.

Increasing or decreasing TruTurn

After TruTurn has reached the default angle fill, you can increase or decrease the turn.

To increase the turn, press and hold the **Action** button next to the **Increase** arrow.

Figure 12 – Increase or decrease turn

To decrease the turn, press and hold the **Action** button next to the **Decrease** arrow.

Note - The **Increase** and **Decrease** arrows illuminate green when active.

Pausing or canceling TruTurn

You can only pause TruTurn when the pump fills the support surface with air. You can cancel TruTurn at any time.

To pause TruTurn, press the **Action** button next to the **Pause** icon.


Figure 13 – TruTurn fill

To cancel TruTurn, press the **Action** button next to the **Cancel** icon. The turn bladder deflates and returns the surface to a flat position.


Figure 14 – Active TruTurn

Locking or unlocking product functions

The **Lock** icon is orange and closed when active.

To lock or unlock product functions:

1. Press and hold the **Unlock** icon for one second to turn on the lock (Figure 15).

Figure 15 – Unlocked

2. Press and hold the **Lock** icon for two seconds to turn off the lock (Figure 16).

Figure 16 – Locked

Note - When the lock is active, you cannot access the other functions.

Powering off the pump

Press and hold the **Power** button for two seconds to power off the pump.

Note

- Allow the support surface to deflate before you unplug the power cord from the wall outlet.
- For optional deflation, remove the CPR plug from the support surface.


Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Name	Part number
Air filter	2874-007-026
IsoTour air hose	2874-001-200
Power cord, type B, 1 meter	2874-007-001
Power cord, type B, 5 meter	2874-007-002
Power cord, type E/F, 1 meter	2874-007-003
Power cord, type E/F, 5 meter	2874-007-004
Power cord, type B, 220V, 1 meter	2874-007-005
Power cord, type B, 220V, 5 meter	2874-007-006
Power cord, type J, 1 meter	2874-007-007
Power cord, type J, 5 meter	2874-007-008
Power cord, type I, 1 meter	2874-007-009
Power cord, type I, 5 meter	2874-007-010
Power cord, type G, 1 meter	2874-007-011
Power cord, type G, 5 meter	2874-007-012
Power cord, type L, 1 meter	2874-007-013
Power cord, type L, 5 meter	2874-007-014
Power cord, type N, 1 meter	2874-007-015
Power cord, type N, 5 meter	2874-007-016
Power cord, type D, 1 meter	2874-007-037
Power cord, type D, 5 meter	2874-007-038
Power cord, type M, 1 meter	2874-007-041
Power cord, type M, 5 meter	2874-007-042

Troubleshooting

Problem	Screen	Cause	Recommended action
Kinked air hoses		Air hose is bent or an obstruction in the hose may cause air flow to be compromised TruTurn was attempted but you have a standard model (bottom cover is orange)	If you have a premium model: 1. Make sure that the hoses are straight and air flow is not obstructed. 2. Press the Action button next to the Alarm off icon to reset the pump. Note - If you have a standard model, contact sales support for the TruTurn option.
Missed connection with air hoses		Air hoses are not connected to the pump or the support surface	1. Make sure that the hoses are seated all the way onto the hose ports on the pump or the support surface. 2. Press the Action button next to the Alarm off icon to reset the pump.
Power loss, product does not turn on		Power cord not seated, power cord unplugged from wall outlet, possible internal damage, button is stuck	1. Make sure that the power cord is plugged into the product and a wall outlet. 2. Secure the power cord to the flex grip. See <i>Connecting the power cord</i> (page 10). 3. Replace the LCD screen.
Power loss while TruTurn is in use		Unintentional or intentional power loss	Pull the hose from the pump or the support surface to activate CPR. The bladder deflates and returns the support surface to a flat position.
Button not responsive		Lock function active, button is stuck	1. Press and hold the Action button next to the Locked icon for two seconds to turn off the lock. 2. Replace the LCD screen. 3. Power cycle the pump.

Problem	Screen	Cause	Recommended action
TruTurn not shown on screen		Kinked hose during setup or IsoTour standard model is connected	If you have a premium model: 1. Make sure that the hoses are not kinked. 2. Verify that you have a premium model (bottom cover is black). Note - If you have a standard model, contact sales support for the TruTurn option.
TruTurn shown on screen with standard model		Graphical user interface shows TruTurn but you have a standard model (bottom cover is orange)	1. Turn off the pump to reset the pump features. 2. Make sure that the air hose is connected before you turn the pump on.

Product alarms

The product alarms for this product are display only. There are no audible alarms.

Icon	Priority alarm	Name	Delay	Therapy interrupted
	Low	Check air connector	Less than 45 seconds	Yes
	Low	Kinked air hose	Less than 5 seconds	Yes

Note

- Verify that the air connector alarm is active when **TruTurn** is active.
- Alarms display on the interface control panel of the pump.
- Alarms reset if the pump loses power.
- Press the **Action** button next to the **Alarm off** icon to reset the alarm.

Preventive maintenance

WARNING - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.

Remove product from service before you perform preventive maintenance. At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Note - Consult your local regulations to dispose of electronic equipment.

Inspect the following items:

- _____ All fasteners are secure
- _____ Pump housing or components (hoses, power cords, or case) are free of cracks, holes, or damage
- _____ Bed frame pump hooks are free of damage
- _____ Pump and attached connectors or hoses are free of air leaks
- _____ Graphical user interface are free of cracks or damage
- _____ Fuse
- _____ Air filter (replace each year)
- _____ All functions on graphical user interface operate
- _____ Run diagnostic test (all tests must pass with OK)

Product serial number:
Completed by:
Date:

Cleaning and disinfecting with SideKick

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred 6" x 10" wipes (2060-000-001) include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface with SideKick wipe:

To clean:

1. Wipe down the external product surface with a fresh, clean wipe to remove all visible soils.
2. Repeat as necessary until the external product surface is visibly clean.
3. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Note - Use as many wipes as necessary.

To disinfect:

1. Clean first.
2. Wipe down the external product surface with a fresh, clean wipe until wet.
3. Allow the external product surface to remain wet for two minutes at room temperature.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Cleaning

WARNING

- Always unplug the product before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports or the hose ports.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
-

CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
-

To clean the external side of the pump:

1. Unplug the power cord from the wall outlet.
2. Wipe the outside of the product with a clean, soft cloth moistened with a mild soap and water solution to remove foreign material.
3. Follow specified contact time in accordance with the manufacturer's instructions for use.
4. Using a clean, soft cloth, wipe the outside of the product with water.
5. Dry the external surface with a clean, dry cloth.

Disinfecting

WARNING

- Always unplug the product before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports or the hose ports.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
-

CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
-

Suggested disinfectants:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Accelerated Hydrogen Peroxide (AHP) 10,000 ppm hydrogen peroxide plus additional additives $\leq 0.13\%$ peroxyacetic acid
- Potassium peroxymonosulfate
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients to avoid the risk of cross-contamination and infection.

To disinfect the pump:

1. Unplug the power cord from the wall outlet.
2. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes. Do not soak the product.

Note - Follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.

3. Rinse the product with a clean, damp cloth to remove excess disinfectant.
4. Dry the product with a clean, dry cloth to remove excess liquid or disinfectant.
5. Allow the product to dry before you return to service.

EMC information

WARNING

- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of IsoTour, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

CAUTION

- This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

Note

- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- This equipment is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Guidance and manufacturer's declaration - electromagnetic emissions		
The 2874 IsoTour pump is intended for use in the electromagnetic environment specified below. The customer or the user of the 2874 IsoTour pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The 2874 IsoTour pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The 2874 IsoTour pump is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	Note - The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communication equipment and the 2874 IsoTour pump			
The 2874 IsoTour pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 2874 IsoTour pump can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 2874 IsoTour pump as recommended below, according to the maximum output power of the communications equipment.			
Band (MHz)	Service	Maximum Power (W)	Minimum Separation Distance (m)

Recommended separation distances between portable and mobile RF communication equipment and the 2874 IsoTour pump			
380-390	TETRA 400	1.8	0.3
430-470	GMRS 460; FRS 460	2.0	0.3
704-787	LTE Band 13, 17	0.2	0.3
800-960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	2.0	0.3
1,700-1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	2.0	0.3
2,400-2,570	Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE Band 7	2.0	0.3
5,100-5,800	WLAN 802.11 a/n	0.2	0.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The 2874 IsoTour pump is suitable for use in the electromagnetic environment specified below. The customer or the user of the 2874 IsoTour pump should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/ output lines	± 2 kV for power supply lines ± 1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV lines to lines	± 0.5 kV, ± 1 kV lines to lines	Main power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	0%U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%U _T for 1 cycle 70%U _T (30% dip in U _T) for 25/30 cycles 0% U _T for 250/300 cycles	0%U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%U _T for 1 cycle 70%U _T (30% dip in U _T) for 25/30 cycles 0% U _T for 250/300 cycles	Main power quality should be that of a typical commercial or hospital environment. If the user of the 2874 IsoTour pump requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note - U _T is the a.c. mains voltage before applications of the test level.			

Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 2874 IsoTour pump, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance $D=(1.2) (\sqrt{P})$ 80 MHz to 800 MHz $D=(2.3) (\sqrt{P})$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Proximity Magnetic Fields IEC 61000-4-39	65 A/m 134.2 kHz	65 A/m	RFID readers and similar generators of magnetic fields should not be operated closer than 50 mm from the electronics or cables of the IsoTour pump.
	7.5 A/m 13.56 MHz	7.5 A/m	
Note • At 80 MHz and 800 MHz, the higher frequency range applies. • These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. • The ISM (Industrial, Scientific, and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.			

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 2874 IsoTour pump is used exceeds the applicable RF compliance level above, the 2874 IsoTour pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 2874 IsoTour pump.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 Vrms.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **InTouch, IsoTour, ProCeed, ProCuity, S3, SideKick, Spirit Select, Stryker, TruTurn**. All other trademarks are trademarks of their respective owners or holders.