Bariatric Non-Powered Support Surface for use with MV3® Bariatric Bed

Operations Manual

REF 2885

CE

EN
FR

Symbols

	Refer to instruction manual/booklet
	Operating instructions / Consult instructions for use
	General warning
	Caution
	Catalogue number
	Serial number
	Authorized representative in the European Community
	CE mark
	European medical device
	For US Patents see www.stryker.com/patents
	Manufacturer
	Safe working load
	Type B applied part
	Mass of product
	Wash by hand
	Do not tumble dry
	Do not dry clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
	Keep dry

	Do not stack more than 6 high
	This side up
	Fragile
	Do not use sharp objects to open the package.
	Arrow, patient left
	Arrow, patient right
	Bolster replacement

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Clinical benefits	4
Expected life	5
Disposal/recycle	5
Contraindications	5
Specifications	5
Product illustration	6
Contact information	7
Serial number location	7
Date of manufacture	7
Setup	8
Setup the support surface	8
Setup the side bolster	8
Operation	10
Storing the bolster sleeve	10
Transferring a patient from one patient support platform to another	11
Cleaning and disinfecting with wipes	11
Care and maintenance	12
Preventive maintenance	13
Accessories and parts	13
Cover replacement	14

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always place the product with the branding at the head end of the frame.
- Always have two people move the product and use the four handles on the back cover.
- Do not place foreign objects inside the cover.
- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
- Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on non-compatible bed frames. Make sure the deck bed matches the width of the support surface.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
- Always keep your fingers away from the path of the zipper when you open or close the zipper on this product.
- Do not use the support surface as a transfer device.
- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
- Do not wash the internal components of this support surface. Discard the support surface if a contamination is found inside.
- Do not immerse the product in cleaning or disinfectant solutions.
- Do not allow liquid to pool on the product.

- Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these cleaning instructions may void your warranty.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the support surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Always expand or retract the litter width before you place a patient on the product.
 - Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ether >3% as these chemicals may damage the support surface cover.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
 - Always wear personal protective equipment to avoid possible skin irritation from contact with the fire barrier.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Bariatric non-powered support surface is for bariatric patients. The product uses three layers of foam with a viscoelastic top layer. Two 6 inch (15.2 cm) wide by 7 inch (17.8 cm) high foam bolsters can be attached to either side of the support surface to fit an expandable bed frame of the **MV3**® bariatric bed. An expandable cover with sealed seams allows the surface to expand with the attached foam bolsters. Articulation cuts provide a break for head of bed elevation, knee Gatch, and foot end articulation points. The bottom cover is compiled of nonskid material. The cover has a sewn in zipper that allows for access to the foam cores with reinforced handles. The main support surface and side bolsters are fitted with fire barriers.

Intended use

The intent of the Bariatric support surface is to assist in the improvement of human patient outcomes by redistributing the pressure and maximizing the patient surface area. This support surface assists in the prevention and treatment of pressure injuries or pressure ulcers (all stages, Unstageable injury, and Deep Tissue injury). Stryker recommends that you implement this product in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional. The dimensions of this product are intended to match the expandable litter deck of the **MV3** bariatric bed frame.

The product is intended for use by patients in an acute care setting. This may include emergency department, Pre-Op, Transport, Endoscopy, GI, critical care, step down, progressive care, med/surg, sub-acute care, and post anesthesia care unit (PACU), OR, or other locations as prescribed.

Indications for use of this product include:

- Human skin
- Patients with a BMI of 30 kg/m² or greater and weigh 250 lb (113.4 kg) or more

This product is not intended to be used with the behavioral health patient populations or pediatric patient populations.

This product is not intended to be sterile, include a measuring function, or be used in a home health care environment.

Clinical benefits

Assists in the prevention and treatment of all pressure ulcers or pressure injuries

Expected life

The main foam core and foam bolster cores of this product have 3 year expected life under normal use, conditions, and with appropriate periodic maintenance.

The main cover and bolster covers of this product have a 3 year expected life under normal use, conditions, and with appropriate periodic maintenance.

Disposal/recycle







Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Contraindications

None known.

Specifications

Safe working load 	1102 lb	500 kg
Note - The patient and accessories must not exceed the safe working load of the Bariatric Non-Powered support surface and the MV3 bed frame.		
Bariatric patient population	BMI > 30 kg/m ²	
Patient weight	≥ 250 lb	≥ 113.4 kg
Minimum patient weight	250 lb	113.4 kg
Model	288505550001	
Length	79.5 in.	201.9 cm
Width	45 in.	114.3 cm
Thickness	7 in.	17.8 cm
Product weight	60.0 lb	27.2 kg
Bolster width, patient right and patient left	6 in.	15.2 cm
Bolster high, patient right and patient left	7 in.	17.8 cm
Top cover material for mattress and bolster	END406 by Dartex ®	
Mattress material	Gel viscoelastic foam	21.5%
	Polyurethane foam	78.5%
Product compliance with fire barrier	16CFR1632, 16CFR1633, CAN Method 27.7- CAN 2-4.2 M77, BS EN 597-1, BS EN 597-2, BS 6807, CAL TB129, Boston BFD IX-11	

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

Product illustration


Figure 1 – Front view


Figure 2 – Back cover

A	Patient right, bolster	C	Patient left, bolster
B	Core	D	Handles, back cover

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Date of manufacture

The year of manufacture is the first four digits of the serial number.

Setup

Setup the support surface

WARNING

- Always place the product with the branding at the head end of the frame.
 - Always have two people move the product and use the four handles on the back cover.
 - Do not place foreign objects inside the cover.
 - Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
 - Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
 - Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on non-compatible bed frames. Make sure the deck bed matches the width of the support surface.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
-

CAUTION

- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the support surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
-

Note

- A minimum of two people is recommended to move this product.
- When you lift this product, use all four handles on the back cover.


To setup the support surface:

1. Make sure that the support surface fits on the bed frame inside the mattress retainers.
2. Make sure that the arrows on the support surface are placed at the foot end of the bed frame.
3. Place linens on the support surface per hospital protocols.


Setup the side bolster

WARNING - Always keep your fingers away from the path of the zipper when you open or close the zipper on this product.


1. Raise the bed height to the full up position and lower the siderails.


2. Use the white pull tab (indicated by the arrow on the cover) to unzip the left or right bolster. Start at the foot end of the support surface and stop at the head end.

3. Place the left or right bolster where indicated.

4. Zip the bolster starting at the head end down to the foot end.


Operation

Storing the bolster sleeve


WARNING - Always keep your fingers away from the path of the zipper when you open or close the zipper on this product.

CAUTION - Always expand or retract the litter width before you place a patient on the product.

1. Raise the bed height to the full up position and lower the siderails.


2. Use the white pull tab (indicated by the arrow on the cover) to unzip the left or right bolster. Start at the foot end of the support surface and stop at the head end.

3. Remove the left or right bolster.

4. Tuck in the bolster sleeve and zip to close starting at the head end down to the foot end.

Transferring a patient from one patient support platform to another

WARNING

- Do not use the support surface as a transfer device.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
 - Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
 - Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
-

CAUTION - Always expand or retract the litter width before you place a patient on the product.

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes to **ON** for both patient support platforms.
3. Adjust the patient support platform heights so that they are level with one another.
4. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.

Cleaning and disinfecting with wipes

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred wipes (2060-000-001 6" x 10" or 2060-000-002 9" x 12") include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface:

1. To clean, wipe external surfaces with a fresh, clean wipe to remove all visible soils. Repeat as necessary until the product is clean.

Note

- Use as many wipes as necessary.
 - Complete step 1 before you disinfect.
2. To disinfect, wipe external surfaces with a fresh, clean wipe until wet. Allow the external surface to remain wet for two minutes at room temperature.
 3. Allow the product to dry before you return it to service.

Care and maintenance

WARNING

- Do not wash the internal components of this support surface. Discard the support surface if a contamination is found inside.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not allow liquid to pool on the product.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these cleaning instructions may void your warranty.
-

CAUTION

- Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ether >3% as these chemicals may damage the support surface cover.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
-

The support surface cover is resistant to the following chemical solutions:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (10000 ppm)
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients, to avoid the risk of cross-contamination and infection.

Preventive maintenance

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage.

Remove product from service before performing preventive maintenance.

Note - Perform support surface care before inspection, if applicable.

Inspect the following items:

- _____ Zipper and covers (top and bottom) are free of tears, cuts, holes, or other openings
- _____ Internal components for signs of staining from fluid ingress or contamination by fully unzipping the covers
- _____ Labels for legibility, proper adherence, and integrity
- _____ Handles are free of rips or cracks and stitching is intact
- _____ Foam has not degraded or come apart
- _____ Fire barrier for rips, cracks, or other visible signs of damage

Product serial number:	
Completed by:	Date:

Accessories and parts

These parts are available for purchase. Call Stryker Customer Service: 1-800-327-0770 for availability and pricing.

Name	Number
Mattress cover assembly	288505550003
Bolster assembly	288505550005

Cover replacement

CAUTION - Always wear personal protective equipment to avoid possible skin irritation from contact with the fire barrier.

Tools required:

- None

Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and Gatch sections to the full down positions.
3. Unzip the cover.
4. Fold the top of the cover to the patient's right side and then remove the foam crib assembly from the bed and set aside.
5. Remove and discard the cover. Leave the fire barrier in place.
6. Place the replacement cover, unzipped and open, on the bed with the black bottom cover on the litter and the top cover folded over the patient's right side of the bed.
7. Place the foam crib assembly on top of the bottom part of the cover to make sure that the foam crib is aligned with the cover.
8. Fold the top cover over the top of the foam crib assembly to make sure that the top cover is aligned with the foam crib assembly.
9. Zip the cover to close.
10. Verify proper operation before you return the product to service.

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