IsoAir® System

Operations Manual

REF 2941

EN
DA
DE
ES
FR
IT
NL
NO
PL
PT
SV

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR)
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Type B applied part, support surface
	Catalogue number
	Serial number
	Batch code
	Manufacturer
	Importer
	European medical device
	Authorized representative in the European Community
	CE mark
	Mass of product
	Safe working load
	Wash by hand
	Do not tumble dry



	Do not dry clean
	Do not iron
	Allow to air dry
	Chlorinated bleach
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provisions for protective earthing or reliance upon installation conditions.
	Fuse
	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1:2012 and CAN/CSA-C22.2 No. 60601-1:16. IEC 60601-1-8:2012
IP21	Solids: Protection from touch by fingers and objects greater than 12 mm Liquids: Protection from the fall of vertical water drops
	For US Patents see www.stryker.com/patents
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Clinical benefits	4
Contraindications	4
Expected service life	4
Disposal/recycle	5
Specifications	5
IsoAir 84 in. (213.3 cm) system	5
Environmental conditions	6
Product illustration - support surface	6
Product illustration - pump	8
Features	8
Button and icon descriptions	9
Contact information	11
Serial number location	11
Setup	12
Setup the support surface	12
Setup the pump	14
Setup with ProCuity bed	16
Setup with InTouch bed	17
Setup with 3002 S3 bed or 3005 S3 bed	18
Setup with GoBed II bed	19
Setup with Spirit Plus bed or Spirit Select bed	20
Connecting the CPR connector and AST cable	21
Applying linens	21
Operation	22
Powering on the pump	22
Transferring a patient from one patient support platform to another	22
Cardiopulmonary resuscitation (CPR)	23
Starting or canceling Max Inflate	23
Starting the static mode without AST	25
Starting the ALP mode without AST	26
Starting or stopping the AST mode	28
Turning the LAL on or off	28
Locking or unlocking the screen	29
Powering off the pump	29
Storing the support surface	29
Storing the pump	30
Troubleshooting	31
Product alarms	33
Accessories and parts	34
Preventive maintenance	35
Cleaning	36
Disinfecting	37
Top cover replacement	38
Fire barrier replacement	38
EMC information	39

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always secure the support surface to the bed frame to prevent the shifting of the product while in use.
- Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
- Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
- Always plug the power cord into the bed auxiliary outlet or a grounded power receptacle. Do not use an extension cord.
- Do not touch patients and the AST socket on the pump at the same time.
- Always secure the pump to the bed frame to prevent the shifting of the product while in use.
- Always check the product before use to make sure there are no broken parts. If the product is damaged call technical support.
- Always make sure that you do not exceed the safe working load of the bed frame before you mount the pump.
- Always connect the AST cable to the AST socket on the pump to avoid the risk of electrical shock.
- Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
- Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
- This product radiates radio frequency energy. If not setup and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity.
- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Do not use this product with an overlay. Patient stability and siderail coverage may be compromised.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Do not use the support surface as a transfer device.
- Do not use the support surface handles to lift or move the support surface with a patient on board.
- Do not stick needles or other sharp objects into a support surface through the support surface cover. Holes may allow fluids to enter the inside (inner core) of the support surface, product damage, or product malfunction.

- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
 - Always make sure that the opposite siderail is raised when you transfer a patient onto the support surface to reduce the risk of patient fall.
 - Always check the CPR connector before use, if visibly soiled, replace. Direct skin contact with soiled material may increase the risk of infection.
 - Always inflate the product and unplug the pump from the mains before you begin cleaning or disinfecting.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports, the hose ports, or the AST socket.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - The **IsoAir** system may interfere with patient ECG measurements. You may need to disconnect the AST cable, turn off the **IsoAir** system, or move the patient to a non-powered system.
 - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
 - Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of the **IsoAir** system, that includes cables specified by the manufacturer.
 - Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.
 - This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating **IsoAir** system or shielding the location.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not tie the support surface to the bed frame with the bottom center D-ring straps to avoid the risk of equipment damage. The bottom center-D ring straps are for storage only.
 - Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Do not place hoses or power cords in walkways to avoid trip hazards.
 - Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
 - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
 - Always be aware that an X-ray taken on this product may cause visible image artifacts. An image artifact may interfere with patient diagnoses.
 - Always store the product within the specified environmental condition values.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not overexpose the product to higher concentration chemical solutions as these may degrade the product.
 - Always follow the manufacturing instructions, failure may affect the useful life of the product.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 2941, **IsoAir**® system consists of a support surface and pump. The support surface is an air powered support surface that consists of air cells and bladders that run across the surface to support the patient and side bolsters. The surface is outfitted with a pump that provides the pressure source to inflate and deflate the air cells. The pump, connected to the air cells through flexible hoses, monitors and automatically adjusts the air in the air cells dependent on the type of therapy selected. This product includes alternating pressure (ALP), low air loss (LAL), max inflate, active sensor technology (AST), and CPR deflate. Comfort control is provided through five levels of alternating low pressure and static pressure.

Intended use

The **IsoAir** system helps in the prevention and treatment of pressure injuries or pressure ulcers (all stages, unstageable injury, and deep tissue injury). Stryker recommends this product to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional. This product provides pressure redistribution and is intended to assist in management of the microclimate of the patient skin.

The **IsoAir** system is intended for use with adult, non-behavioral patients. This product is for use with a bed frame and a top sheet. The system is intended for use in an acute care and long-term care units, including general hospital use, critical care, step down care, Med/Surg, PACU, and progressive care settings.

This product is not intended for use in home healthcare environments. This product is not sterile.

Clinical benefits

Assists in the prevention and treatment of all pressure ulcers or pressure injuries

Contraindications

Unstable spine

Expected service life

The **IsoAir** support surface has a 3 year expected service life under normal use conditions and with appropriate periodic maintenance.

The **IsoAir** pump has a 5 year expected service life under normal use conditions and with appropriate periodic maintenance.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

IsoAir 84 in. (213.3 cm) system

Number	Plug type	Power cord 1M	Power cord 5M
2941-000-004	E/F	2874-007-003	2874-007-004
2941-000-007	G	2874-007-011	2874-007-012
2941-000-010	I	2874-007-009	2874-007-010
2941-000-013	L	2874-007-013	2874-007-014
2941-000-015	N	2874-007-015	2874-007-016
2941-000-017	K	2941-700-008	2941-700-009
2941-000-019	J	2874-007-007	2874-007-008
2941-000-021	B	2874-007-001	2874-007-002

 <p>Note - Safe working load indicates the sum of the patient, mattress, and accessory weight.</p>	500 lb	226 kg
Support surface		
Model	2941-001-100	
Length	84.0 in.	213.3 cm
Width	35.0 in.	88.9 cm
Thickness	8.0 in.	20.3 cm
Weight	25 lb	11.3 kg
Cover material	Dartex® Endurance cover	
Compatible with bed frame	ProCuity™ , 3002 S3® , 3005 S3® , InTouch® , Spirit Plus® , Spirit Select® , GoBed® II	

Dartex is a registered trademark of the Dartex Coatings, Ltd.







Pump		
Model	2941-001-001	
Height	17.25 in.	43.8 cm
Width	9.25 in.	23.5 cm

Pump			
Depth		7.25 in.	18.4 cm
Weight		17 lb	7.7 kg
Interface control panel		LCD	
Input voltage		100-240 V	
Electric current		1.2 A	
Input frequency		50/60 Hz	
Noise level		<50 dBa	
Mode of operation		Continuous	
Power cord	Short	3.2 ft	1 m
	Long	16.4 ft	5 m
Compliance		IEC 60601-1 3.1 Edition; RoHS Directive 2015/863/EU, Annex II to directive 2011/65EU; REACH (EC) 1907/2006 EU; IEC 60601-1-8:2012, IEC 60601-1-2:2014	

Note - The power cord for this product has a ground conductor that is only used as a functional earth.


Stryker reserves the right to change specifications without notice.

Environmental conditions

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Product illustration - support surface


CAUTION - Do not tie the support surface to the bed frame with the bottom center D-ring straps to avoid the risk of equipment damage. The bottom center-D ring straps are for storage only.

A	Top cover
B	Strap handle
C	Law label location
D	AST cable

E	CPR connector
F	Tie down strap
G	D-ring (see Applying linens)
H	D-ring strap, bottom center
I	Bottom cover

Product illustration - pump

A	Graphical user interface
B	Status indicator
C	Bumper
D	Transport handle
E	Hook
F	Power inlet

G	Anti-vibration foot
H	HEPA filter
I	Flex grip
J	Fuse
K	Hose port
L	AST port connector

Features

Shear management

The product moves with the patient to assist in the reduction of shear forces at the interface between the surface and the patient.

Low Air Loss (LAL)

The product provides air flow to help manage microclimate of the skin.

Alternating low pressure (ALP)

The support surface is constructed with rows of air cells connected to two different air lines (C). Odd row air cells (A) are connected to the A air line and even row cells (B) are connected to the B air line (Figure 1).

Active sensor technology (AST)

When AST is on, sensors in the sacral region (D) of the support surface help control patient immersion (Figure 1).

Side bolster

The left and right side bolsters (E) provide exterior support to the patient. This helps to secure the patient within the support surface (Figure 1).




Figure 1 – Air cells

A	A air cell	D	AST air cell
B	B air cell	E	Left and right side bolster
C	Air line, two line manifold	F	Bottom foam assembly

Button and icon descriptions

Sample screens

Icon	Name	Description
	Power button	See <i>Powering on the pump</i> (page 22) or <i>Powering off the pump</i> (page 29)
	Action button, right	Allows the operator to select an action on the screen

Icon	Name	Description
	Action button, left	Allows the operator to select an action on the screen
	Lock or unlock	See <i>Locking or unlocking the screen</i> (page 29)
	Home	Returns to the home screen
	Low Air Loss (LAL) on	See <i>Turning the LAL on or off</i> (page 28)
	Low Air Loss (LAL) off	
	Static mode on	Adjusts comfort settings when used without AST. in this mode, all cells are set to the same pressure
	Static mode off	Static mode off when you choose to cancel this therapy
	Decrease or increase	Adjusts pressure in the air cells
	ALP on	Provides therapy with pressure differences between alternating air cells
	ALP off	ALP mode off when you choose to cancel this therapy
	Mode menu	Mode menu screen
	Max Inflate full	Inflates the support surface to a maximum pressure
	Max Inflate filling	Shows the surface filling to max inflate
	Max Inflate mode	Max Inflate off, you choose to cancel max inflate or the max inflate timer has expired
	AST on	Optimal air pressure recalculates and adjusts as needed every 3 minutes
	AST off	AST mode off when you choose to cancel this therapy
	Disconnected air hose alarm	Detects air hose connection error
	Kinked air hose alarm	Detects a bend or twist in the air hose
	AST malfunction alarm	Detects the AST malfunction error

	Remove from use alarm	Detects error
	Alarm off	Cancels alarm

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.



Stryker Medical
 3800 E. Centre Avenue
 Portage, MI 49002
 USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Setup

Setup the support surface

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

WARNING

- Always secure the support surface to the bed frame to prevent the shifting of the product while in use.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
-

CAUTION

- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
-

To setup the support surface:

1. Allow the product to reach room temperature.
2. Lower the siderails to the lowest position of the bed frame.
3. Raise the Fowler to 30 degrees.
4. Remove the headboard and footboard from the bed frame.
5. Loosen the straps that secure the support surface.
6. Unroll and place the support surface with the hose at the foot end of the litter deck. Place the bottom cover face down.
7. Start on one side of the support surface and litter deck. Locate the head end and foot end retaining slot locations on the litter deck. See the setup for the appropriate bed.
 - a. *Setup with ProCuity bed (page 16)*
 - b. *Setup with InTouch bed (page 17)*
 - c. *Setup with 3002 S3 bed or 3005 S3 bed (page 18)*
 - d. *Setup with GoBed II bed (page 19)*
 - e. *Setup with Spirit Plus bed or Spirit Select bed (page 20)*

Note

- This product has ten tie down straps. A minimum of four should be used. See the appropriate bed setup for the number of straps and tie down locations.
 - Do not attach the tie down straps to the bed siderails.
8. Thread the tie down strap (E) through the Fowler hole (D) (Figure 2).


Figure 2 – Insert tie down straps

9. Put the tie down strap (E) through the D-rings (F) (Figure 3).


Figure 3 – Tie down strap

10. Loop the tie down strap between the D-rings (G) (Figure 3). Pull the tie down strap until secure (Figure 4).


Figure 4 – Pull tie down strap

11. Insert the loose end of the tie down strap behind the Fowler.
12. Repeat steps 8 to 10 for the other tie down strap.
13. Check that all the tie down straps are secure.
14. Replace the headboard and footboard that were removed in step 4.
15. Place linens on the support surface per hospital protocols.

Setup the pump

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

WARNING

- Always plug the power cord into the bed auxiliary outlet or a grounded power receptacle. Do not use an extension cord.
 - Do not touch patients and the AST socket on the pump at the same time.
 - Always secure the pump to the bed frame to prevent the shifting of the product while in use.
 - Always check the product before use to make sure there are no broken parts. If the product is damaged call technical support.
 - Always make sure that you do not exceed the safe working load of the bed frame before you mount the pump.
 - Always connect the AST cable to the AST socket on the pump to avoid the risk of electrical shock.
 - Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
 - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
 - This product radiates radio frequency energy. If not setup and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity.
-

CAUTION

- Do not place hoses or power cords in walkways to avoid trip hazards.
 - Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
 - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
-

To setup the pump:

1. Allow the product to reach room temperature.
2. Attach the power cord to the pump.
3. Insert the power cord under the power cord retention flex grip (A) (Figure 5).


Figure 5 – Cord retention flex grip

4. Plug the pump power cord into an accessible hospital grade outlet.
5. Place the pump power cord out of the path of foot traffic to avoid trip hazards.
6. Place the pump on the bed.
 - a. *Setup with ProCuity bed* (page 16)

- b. *Setup with InTouch bed* (page 17)
 - c. *Setup with 3002 S3 bed or 3005 S3 bed* (page 18)
 - d. *Setup with GoBed II bed* (page 19)
 - e. *Setup with Spirit Plus bed or Spirit Select bed* (page 20)
7. Route the hose on the patient left side of the footboard.

Note

- Allow space around the pump to connect and disconnect the CPR connector and power cord.
- Allow space on the integrated pump rack or footboard handle for other devices.

Setup with ProCuity bed


Figure 6 – Retaining slot locations for tie down straps

Figure 7 – Pump location and hose routing

Setup with InTouch bed


Figure 8 – Retaining slot locations for tie down straps


Figure 9 – Pump location and hose routing

Setup with 3002 S3 bed or 3005 S3 bed


Figure 10 – Retaining slot locations for tie down straps

Figure 11 – Pump location and hose routing

Setup with GoBed II bed


Figure 12 – Retaining slot locations for tie down straps


Figure 13 – Pump location and hose routing

Setup with Spirit Plus bed or Spirit Select bed


Figure 14 – Retaining slot locations for tie down straps


Figure 15 – Pump location on the accessory hanger and hose routing


Figure 16 – Pump location on footboard handle and hose routing

Connecting the CPR connector and AST cable

To connect the CPR connector and AST cable:

1. Connect the CPR connector to the pump (C) (Figure 17).

Figure 17 – Connect CPR connector and AST cable

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Applying linens

To apply linens:

1. Part the D-rings (C) (Figure 18).
2. Thread the four linen corners through the D-rings (C) attached to the bottom cover to secure the linens (A) to the support surface (B).


Figure 18 – Apply linens

Note - This product is designed for use with non-fitted sheets. For best results, do not pull the linens too tight to avoid the hammock effect.

Operation

WARNING

- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Do not use this product with an overlay. Patient stability and siderail coverage may be compromised.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.

CAUTION - Always be aware that an X-ray taken on this product may cause visible image artifacts. An image artifact may interfere with patient diagnoses.

Powering on the pump

To power on the pump:

1. Follow *Setup the pump* (page 14) to prepare the pump for operation.
2. Press and hold the **Power** button until **IsoAir** appears on the screen (Figure 19).


Figure 19 – Startup screen

Figure 20 – Home screen, operator view

Note - The ALP, AST, and the LAL modes that are powered on are shown in green on the home screen (Figure 20).

Transferring a patient from one patient support platform to another

WARNING

- Do not use the support surface as a transfer device.
- Do not use the support surface handles to lift or move the support surface with a patient on board.
- Do not stick needles or other sharp objects into a support surface through the support surface cover. Holes may allow fluids to enter the inside (inner core) of the support surface, product damage, or product malfunction.
- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always make sure that the opposite siderail is raised when you transfer a patient onto the support surface to reduce the risk of patient fall.

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Press and wait for **Max Inflate** to fully inflate. See *Starting or canceling Max Inflate* (page 23).
2. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
3. Set the brakes to **ON** for both patient support platforms.
4. Adjust the patient support platform heights so that they are level with one another.
5. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.

Cardiopulmonary resuscitation (CPR)

CPR compressions may begin as soon as you remove the CPR connector.

Note - It is recommended that you remove the CPR connector. Follow your hospital protocols to complete CPR.

1. Depress the tabs on both sides of the CPR connector. Pull straight out to disconnect the CPR from the pump (Figure 21).


Figure 21 – CPR connector

2. Remove the AST cable to prevent damage to the cable (recommended).

Starting or canceling Max Inflate

To start or cancel Max Inflate:

1. On the home screen, press the **Action** button next to the **Mode** icon (Figure 22).

Figure 22 – Home screen

2. To start maximum inflation, press the **Action** button next to the **Max Inflate** icon (Figure 23).

Figure 23 – Max Inflate

Note - The **Max Inflate** icon blinks until the support surface inflates to the maximum pressure setting.


Figure 24 – Max Inflate blinking

Note - The countdown is displayed and pressure remains for 15 minutes.


Figure 25 – Max Inflate timer

Note

- The support surface inflates to the maximum pressure setting.
 - The manual adjustments of the air cells are disabled.
3. To cancel Max Inflate, press the **Action** button next to the **Cancel** icon or wait until the time expires (Figure 26).

Figure 26 – Cancel Max Inflate

Note - No other functions are available until Max Inflate time expires or is canceled.

Starting the static mode without AST

The static mode keeps the A/B air cells at constant pressure and allows for manual adjustment of air cell pressure.

To start the static mode:

1. On the home screen, press the **Action** button next to the **Mode** icon (Figure 27).


Figure 27 – Mode

2. To start the static mode, press the **Action** button next to **Static Mode** icon (Figure 28).


Figure 28 – Static mode off

3. Press the **Action** button next to the **up-arrow** or **down arrow** to increase or decrease the static pressure level from 1 to 5 (Figure 29).


Figure 29 – Static mode on

Note - The adjustment arrows illuminate green when you press the **Action** button.

Static mode	Pressure
1	10 mmHg
2	15 mmHg
3 (default)	20 mmHg
4	25 mmHg
5	30 mmHg

Starting the ALP mode without AST

ALP mode alternately inflates and deflates the A/B air cells to relieve pressure and allows for manual adjustment of air cell pressure.

To start the ALP mode:

1. On the home screen, press the **Action** button next to the **Mode** icon (Figure 30).


Figure 30 – Mode

2. To start ALP mode, press the **Action** button next to the **ALP mode** icon (Figure 31).

Figure 31 – ALP mode off

Note - The support surface alternates between zone A and zone B in approximately six minute intervals.

3. Press the **Action** button next to the **ALP mode** icon (Figure 32).

Figure 32 – ALP mode on

4. Press the **Action** button next to the **up-arrow** or **down arrow** to increase or decrease the ALP level from 1 to 5 (Figure 33).


Figure 33 – ALP mode level 3

Note - The adjustment arrows illuminate green when you press the **Action** button.

ALP mode	Setting (A/B or B/A)
1	13/7 mmHg
2	18/12 mmHg
3 (default)	23/17 mmHg
4	28/22 mmHg
5	33/27 mmHg

Starting or stopping the AST mode

The AST mode is on by default. The AST mode automatically adjusts air cell pressures to control patient immersion. The AST mode can be on when static mode or ALP mode is on. You can start or stop the AST mode from either the home screen or the mode screen.

1. To start the AST mode, press the **Action** button next to the **AST mode** icon (Figure 34).


Figure 34 – AST mode off

2. To stop the AST mode, press the **Action** button next to the **AST mode** icon (Figure 35).

Figure 35 – AST mode on

Turning the LAL on or off

LAL is on by default. You can turn LAL on or off from either the home screen or the mode screen.

1. To turn on LAL, press the **Action** button next to the **LAL off** icon (Figure 36).


Figure 36 – LAL off

2. To turn off LAL, press the **Action** button next to the **LAL on** icon (Figure 37).


Figure 37 – LAL on

Locking or unlocking the screen

You are unable to change modes when the screen is locked. The alarms continue to function with the lock on.

1. To lock the screen, press and hold the **Action** button next to the **Unlocked** icon for approximately one second (Figure 38).

Figure 38 – Unlocked screen

2. To unlock the screen, press and hold the **Action** button next to the **Locked** icon for two seconds (Figure 39).

Figure 39 – Locked screen

Powering off the pump

Press the **Power** button and hold for two seconds to power off the pump.

Note - The pump will power up and default to the previous setting if you unplug and re-plug the power cord.

Storing the support surface

CAUTION - Always store the product within the specified environmental condition values.

To store the support surface:

1. Start Max Inflate, see *Starting or canceling Max Inflate* (page 23).
2. Clean and dry the product.
3. Power off the pump.
4. Press the tab in on the CPR connector to disconnect.
5. Remove the CPR connector from the pump.
6. Deflate the support surface.
7. Roll the support surface from the head end to the foot end of the bed.

Note - Make sure that you insert the AST cable and hose when you roll the support surface.

8. Use the straps on the bottom center to tie the rolled support surface.

Note - Use the foot end handles to lift and move the support surface for storage.

9. Follow your hospital protocols to store the product.


Storing the pump


CAUTION - Always store the product within the specified environmental condition values.

To store the pump:

1. Power off the pump.
2. Unplug the power cord from the hospital grade outlet.
3. Clean and dry the product.
4. Follow your hospital protocols to store the product.

Troubleshooting

Problem	Screen	Cause	Recommended action
AST malfunction alarm		AST cable is not plugged in	<ol style="list-style-type: none"> 1. Make sure that the AST cable is connected. 2. Press the Action button next to the Alarm off icon.  3. You may continue to use the product. Turn off the AST mode within three seconds of canceling the alarm. 4. If the condition persists, contact a qualified service technician.
		AST cable is broken	Contact a qualified service technician.
		AST circuit inside the support surface is broken or shorted	Contact a qualified service technician.
Kinked air hose alarm		Bend or some obstruction in the air hose that may cause a compromise in the air flow	<ol style="list-style-type: none"> 1. Make sure that the air hose is straight and air flow is not obstructed. 2. Make sure that the support surface is secure to the bed. See <i>Setup the support surface</i>. 3. Route the pump and hose. See <i>Setup the support surface</i>. 4. Press the Action button next to the Alarm off icon.  5. If the condition persists, contact a qualified service technician.

Problem	Screen	Cause	Recommended action
Disconnected air hose alarm		Air hose is disconnected from the pump or there is a leak in the system	<ol style="list-style-type: none"> 1. Make sure that the air hoses are seated all the way into the air hose ports on the pump. 2. Press the Action button next to the Alarm off icon.  3. If the condition persists, contact a qualified service technician.
Remove from use alarm		Software or calibration data corruption error	Contact a qualified service technician.
Power loss, product does not turn on		Power cord not seated, power cord unplugged from outlet	<ol style="list-style-type: none"> 1. Make sure that the power cord is plugged into the product and the outlet. 2. Secure the power cord to the flex grip. See <i>Setup the pump</i>. 3. If the condition persists, contact a qualified service technician.
Button not responsive		Lock may be on	<ol style="list-style-type: none"> 1. Press the Action button next to the Lock on icon to unlock the screen. 2. If the condition persists, contact a qualified service technician.

Product alarms

The product alarms for this product are visual only. There are no audible alarms.

Icon	Priority alarm	Name	Delay	Therapy interrupted
	Low	Disconnected air hose alarm	Less than 16 minutes	Yes
	Low	Kinked air hose alarm	Less than 5 seconds	Yes
	Low	AST malfunction alarm	Less than 5 seconds	Yes
	Low	Remove from use alarm	Less than 5 seconds	Yes

Note

- Alarms display on the interface control panel of the pump.
- Alarms reset if there is power loss to the pump.
- Press the **Action** button next to the **Alarm off** icon  to turn off active alarms and reset the alarm system.

Accessories and parts

These parts are currently available for purchase. Call Stryker Customer Service at 1-800-327-0770 for availability and pricing.

Name	Number
IsoAir support surface, 84 in. (213.3 cm), END406	2941-001-100
IsoAir pump	2941-001-001
IsoAir, top cover, 84 in. (213.3 cm) END406	2941-700-006
IsoAir, fire barrier, 84 in. (213.3 cm)	2941-700-007
Power cord, type B, 1 meter	2874-007-001
Power cord, type B, 5 meter	2874-007-002
Power cord, type E/F, 1 meter	2874-007-003
Power cord, type E/F, 5 meter	2874-007-004
Power cord, type G, 1 meter	2874-007-011
Power cord, type G, 5 meter	2874-007-012
Power cord, type I, 1 meter	2874-007-009
Power cord, type I, 5 meter	2874-007-010
Power cord, type L, 1 meter	2874-007-013
Power cord, type L, 5 meter	2874-007-014
Power cord, type N, 1 meter	2874-007-015
Power cord, type N, 5 meter	2874-007-016
Power cord, type K, 1 meter	2941-700-008
Power cord, type K, 5 meter	2941-700-009
Power cord, type J, 1 meter	2874-007-007
Power cord, type J, 5 meter	2874-007-008
Fuse	2874-007-021

Preventive maintenance

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product use. Service only by qualified service technicians.

CAUTION - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note

- Remove product from service before performing preventive maintenance.
- Clean and disinfect the exterior of the support surface before inspection, if applicable.

- _____ Zipper and covers are free of tears, cuts, holes, or other openings
- _____ Unzip the cover to inspect the inside surface and core for signs of staining due to fluid ingress or contamination
- _____ Air cells are free of excessive wear, cracks, tears, or other damage
- _____ Fire barrier cover for excessive wear
- _____ All connectors are free of damage
- _____ AST cable free from visible cracks or damage
- _____ Pump housing or components (hoses, power cords, or case) for cracks, holes, or not damaged
- _____ Pump hooks that hang the pump on the bed frame are not damaged
- _____ No air leaks from the pump or the attached connectors or hoses
- _____ Graphical user interface screen is not cracked or damaged
- _____ HEPA filter (replace each year)
- _____ Fuse
- _____ AST function
- _____ Max inflate function
- _____ CPR release
- _____ Run diagnostic test (All OK)

Note - Replace worn or damaged components as necessary.

Pump serial number:
Completed by:
Date:

Support surface serial number:
Completed by:
Date:

Cleaning

WARNING

- Always check the CPR connector before use, if visibly soiled, replace. Direct skin contact with soiled material may increase the risk of infection.
 - Always inflate the product and unplug the pump from the mains before you begin cleaning or disinfecting.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports, the hose ports, or the AST socket.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not overexpose the product to higher concentration chemical solutions as these may degrade the product.
 - Always follow the manufacturing instructions, failure may affect the useful life of the product.
-

To clean the support surface and pump:

1. Inflate the support surface with the pump.
2. Unplug the power cord from the hospital grade outlet.

Note - Leave the hose connected to the support surface and pump to retain the air.

3. Complete steps 4 - 6 for the support surface and repeat for the pump.
4. Wash the outside of the product with a clean, soft cloth, moistened with a mild soap and water solution, to remove foreign material.
5. Rinse the outside of the product using a clean, soft cloth, and water.
6. Dry the product with a clean, dry cloth.
7. Insert the power cord under the power cord retention flex grip.
8. Plug the pump power cord into a hospital grade outlet.

Disinfecting

WARNING

- Always check the CPR connector before use, if visibly soiled, replace. Direct skin contact with soiled material may increase the risk of infection.
 - Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inflate the product and unplug the pump from the mains before you begin cleaning or disinfecting.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not overexpose the product to higher concentration chemical solutions as these may degrade the product.
 - Always follow the manufacturing instructions, failure may affect the useful life of the product.
-

Suggested disinfectants:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Potassium peroxymonosulfate
- Accelerated hydrogen peroxide
- 70% isopropyl alcohol

To disinfect the support surface and pump:

1. Inflate the support surface with the pump.
2. Unplug the power cord from the hospital grade outlet.

Note - Leave the hose connected to the support surface and pump to retain the air.

3. Clean and dry the support system and pump before you apply disinfectants.
4. Complete steps 5 - 8 for the support surface and repeat for the pump.
5. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes. Do not soak the product.

Note - Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.

6. Rinse the product with a clean, damp cloth to remove excess disinfectant.
7. Dry the product with a clean, dry cloth to remove excess liquid or disinfectant.
8. Allow the product to dry before you return to service.
9. Insert the power cord under the power cord retention flex grip.
10. Plug the pump power cord into a hospital grade outlet.

Top cover replacement

Tools required:

- None

Procedure:

1. Raise the bed height to the highest position.
2. Lower the Fowler, Gatch sections, and the siderails to the full down positions.
3. Unzip the cover. The zipper starts in the middle of the foot end of the support surface.

Note - Use caution not to damage the fire barrier.

4. Remove and discard the top cover per hospital protocols.
5. Reverse to install.
6. Verify proper operation before you return the product to service.

Fire barrier replacement

Tools required:

- None

Procedure:

1. Raise the bed height to the highest position.
2. Place the litter surface in a flat, horizontal position.
3. Unzip and remove the support surface top cover.
4. Unclip the head end and foot end center belt retainers from the bottom cover.
5. Remove the CPR connector to allow the support surface to deflate.
6. Unzip the fire barrier.
7. Move the fire barrier over the side bolsters and air cells toward the head end to remove the fire barrier.
8. Reverse steps to reinstall.
9. Run the pump diagnostic test. All tests must pass.
10. Verify proper operation before you return the product to service.

EMC information

WARNING

- The **IsoAir** system may interfere with patient ECG measurements. You may need to disconnect the AST cable, turn off the **IsoAir** system, or move the patient to a non-powered system.
- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of the **IsoAir** system, that includes cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.
- This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating **IsoAir** system or shielding the location.

Note

- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- This equipment is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
The IsoAir system is suitable for use in the electromagnetic environment specified below. The customer or the user of the IsoAir system should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+8 kV contact ±15 kV air	+8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/ burst IEC 61000-4-4	+2 kV for power supply lines ±1 kV for input/ output lines	+2 kV for power supply lines ±1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+1 kV lines to lines +2 kV lines to earth	+1 kV lines to lines +2 kV lines to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	0% U_T for 0.5 cycle at 0°, 45, 90, 135, 180, 225, 270, and 315 0% U_T for 1 cycle 70% U_T (30% dip in U_T) for 25 cycles 0% U_T for 250 cycles	0% U_T for 0.5 cycle at 0°, 45, 90, 135, 180, 225, 270, and 315 0% U_T for 1 cycle 70% U_T (30% dip in U_T) for 25 cycles 0% U_T for 250 cycles	Main power quality should be that of a typical commercial or hospital environment. If the user of the IsoAir system requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note - U_T is the AC mains voltage before applications of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and IsoAir System.			
The IsoAir system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the IsoAir system can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IsoAir system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=[3.5/\sqrt{1}] \cdot \sqrt{P}$	80 MHz to 800 MHz $d=[3.5/\sqrt{E1}] \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d=[7/\sqrt{E1}] \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

IsoAir system is suited for use in the electromagnetic environment specified below. The customer or the user of IsoAir system should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000- 4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms</p> <p>6 Vrms in ISM bands 150 kHz to 80 MHz</p> <p>10 V/m 80 MHz to 2.7 GHz</p>	<p>3 V</p> <p>6 V in ISM bands</p> <p>10 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of IsoAir, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended Separation Distance</p> $d=(0.35)\sqrt{P}$ <p>80 MHz to 800 MHz</p> $d=(0.70)\sqrt{P}$ <p>800 MHz to 2.7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note - At 80 MHz and 800 MHz, the higher frequency range applies.

Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **IsoAir** is used exceeds the applicable RF compliance level above, **IsoAir** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **IsoAir**.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

Guidance and Manufacturer's declaration - Electromagnetic Emissions

IsoAir is intended for use in an electromagnetic environment specified below. The customer or the user of **IsoAir** should make sure that it is used in such an environment.




Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	IsoAir uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	IsoAir is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

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