IsoAir® support surface replacement kit instructions

If you ordered Model 2941 IsoAir System, disregard these instructions and refer to the Model 2941 IsoAir System Operations Manual.

Use these instructions if you are replacing Model 2940 IsoAir support surface with Model 2941-001-100 IsoAir support surface.

Tools required

- None

Procedure

1. Loosen the tie down straps that secure the Model 2940 support surface to the bed frame.
2. Remove and discard the Model 2940 support surface.
3. Unpack the 2941-001-100 support surface and discard the packaging.
4. Setup the 2941-001-100 support surface follow the Model 2940 AO-SM-70-INT Operations Manual.

Note

- For instructions about how to use the Model 2941 IsoAir Support Surface with the Model 2940 IsoAir Pump, see the AO-SM-70-INT Operations Manual. The content in this instruction sheet is a supplement to your Model 2940 IsoAir Operations Manual.
- Do not service the 2941-001-100 support surface with the Model 2940 Operations/Maintenance/Service Manuals.
- Do not launder the 2941-001-100 support surface.
- Use this instruction sheet for cleaning and disinfecting instructions of your 2941-001-100 support surface.
- To view your Operations or Maintenance Manual online, see <https://techweb.stryker.com/>.

For available service parts, contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.


WARNING

- Always secure the support surface to the bed frame to prevent the shifting of the product while in use.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
 - Always check the CPR connector before use. Replace if visibly soiled. Direct skin contact with soiled material may increase the risk of infection.
 - Always inflate the product and unplug the pump from the mains before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluids may cause premature degradation of components and may cause unpredictable safety and performance of this product.
 - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this instruction sheet.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
- Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
- Do not tie the support surface to the bed frame with the bottom center D-ring straps to avoid the risk of equipment damage. The bottom center-D ring straps are for storage only.
- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
- Do not overexpose the product to higher concentration chemical solutions to avoid the risk of product degradation.
- Always follow the manufacturing instructions, failure may affect the useful life of the product.

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR)
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Type B applied part, support surface - 2941-001-100
	Catalogue number
	Serial number
	Batch code
	Manufacturer
	European medical device
	Authorized representative in the European Community

	CE mark
	Mass of product
	Safe working load
	Wash by hand
	Do not tumble dry
	Do not dry clean
	Do not iron
	Allow to air dry
	Chlorinated bleach
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	For US Patents see www.stryker.com/patents
	Keep dry
	Stacking limit by number
	Do not use sharp objects to open the package

Contraindications

Unstable spine

Expected service life

The IsoAir support surface has a 3 year expected service life under normal use conditions and with appropriate periodic maintenance.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications


 Note - Safe working load indicates the sum of the patient, mattress, and accessory weight.	500 lb	226 kg
Model	2941-001-100	
Length	84.0 in.	213.3 cm
Width	35.0 in.	88.9 cm
Thickness	8.0 in.	20.3 cm
Weight	25 lb	11.3 kg
Cover material	Dartex® Endurance cover	
Compatible with bed frame	ProCuity™ , 3002 S3® , 3005 S3® , InTouch® , Spirit Plus® , Spirit Select® , GoBed® II	

Dartex is a registered trademark of the Dartex Coatings, Ltd.

Stryker reserves the right to change specifications without notice.

Product illustration - support surface

CAUTION - Do not tie the support surface to the bed frame with the bottom center D-ring straps to avoid the risk of equipment damage. The bottom center-D ring straps are for storage only.


A	Top cover
B	Strap handle
C	Law label location
D	AST cable

E	CPR connector
F	Tie down strap
G	D-ring (see Applying linens)
H	D-ring strap, bottom center
I	Bottom cover

Serial number and date of manufacture

You will find the serial number and date of manufacture on the law label (A) on the foot end of the product (B).

Setup the support surface

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

WARNING

- Always secure the support surface to the bed frame to prevent the shifting of the product while in use.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
-

CAUTION

- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Do not tie the support surface to the bed frame with the bottom center D-ring straps to avoid the risk of equipment damage. The bottom center-D ring straps are for storage only.
-

To setup the support surface:

1. Allow the product to reach room temperature.
2. Lower the siderails to the lowest position of the bed frame.
3. Raise the Fowler to 30 degrees.
4. Remove the headboard and footboard from the bed frame.
5. Loosen the straps that secure the support surface.
6. Unroll and place the support surface with the hose at the foot end of the litter deck. Place the bottom cover face down.
7. Start on one side of the support surface and litter deck. Locate the head end and foot end retaining slot locations on the litter deck. See the setup for the appropriate bed.
 - a. *Setup Model 2940 pump with ProCuity bed (page 9)*
 - b. *Setup Model 2940 pump with InTouch bed (page 10)*
 - c. *Setup Model 2940 pump with 3002 S3 bed or 3005 S3 bed (page 11)*
 - d. *Setup Model 2940 pump with GoBed II bed (page 12)*
 - e. *Setup Model 2940 pump with Spirit Plus bed or Spirit Select bed (page 13)*

Note

- This product has ten tie down straps. A minimum of four should be used. See the appropriate bed setup for the number of straps and tie down locations.
 - Do not attach the tie down straps to the bed siderails.
8. Thread the tie down strap (E) through the Fowler hole (D) (Figure 1).


Figure 1 – Insert tie down straps

9. Put the tie down strap (E) through the D-rings (F), (Figure 2).


Figure 2 – tie down strap

10. Loop the tie down strap between the D-rings (G) (Figure 2). Pull the tie down strap until secure (Figure 3).


Figure 3 – Pull tie down strap

11. Insert the loose end of the tie down strap behind the Fowler.
12. Repeat steps 8 to 10 for the other tie down strap.
13. Check that all the tie down straps are secure.
14. Replace the headboard and footboard that were removed in step 4.
15. Place linens on the support surface per hospital protocols.

Setup Model 2940 pump with ProCuity bed


Figure 4 – Retaining slot locations for tie down straps

Figure 5 – Pump location and hose routing

Setup Model 2940 pump with InTouch bed


Figure 6 – Retaining slot locations for tie down straps

Figure 7 – Pump location and hose routing

Setup Model 2940 pump with 3002 S3 bed or 3005 S3 bed


Figure 8 – Retaining slot locations for tie down straps

Figure 9 – Pump location and hose routing

Setup Model 2940 pump with GoBed II bed


Figure 10 – Retaining slot locations for tie down straps

Figure 11 – Pump location and hose routing

Setup Model 2940 pump with Spirit Plus bed or Spirit Select bed


Figure 12 – Retaining slot locations for tie down straps


Figure 13 – Pump location on the accessory hanger and hose routing


Figure 14 – Pump location on footboard handle and hose routing

Applying linens

To apply linens:

1. Part the D-rings (C) (Figure 15).
2. Thread the four linen corners through the D-rings (C) attached to the bottom cover to secure the linens (A) to the support surface (B).


Figure 15 – Apply linens

Note - This product is designed for use with non-fitted sheets. For best results, do not pull the linens too tight to avoid the hammock effect.

Cleaning

WARNING

- Always check the CPR connector before use. Replace if visibly soiled. Direct skin contact with soiled material may increase the risk of infection.
 - Always inflate the product and unplug the pump from the mains before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluids may cause premature degradation of components and may cause unpredictable safety and performance of this product.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not overexpose the product to higher concentration chemical solutions to avoid the risk of product degradation.
 - Always follow the manufacturing instructions, failure may affect the useful life of the product.
-

To clean the support surface covers after each patient use, follow these steps in order:

1. Inflate the support surface with the pump.
2. Disconnect the mains power from the receptacle.
Note - Leave the hose connected to the support surface to retain the air.
3. Moisten a clean, soft cloth with a mild soap and water solution.
4. Wash the outside of the product with a clean, soft cloth, moistened with a mild soap and water solution, to remove foreign material.
5. Rinse the outside of the product using a clean, soft cloth and water.
6. Dry the product with a clean, dry cloth.

Disinfecting

WARNING

- Always check the CPR connector before use. Replace if visibly soiled. Direct skin contact with soiled material may increase the risk of infection.
 - Always inflate the product and unplug the pump from the mains before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluids may cause premature degradation of components and may cause unpredictable safety and performance of this product.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not overexpose the product to higher concentration chemical solutions to avoid the risk of product degradation.
 - Always follow the manufacturing instructions, failure may affect the useful life of the product.
-

Suggested disinfectants:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Potassium peroxymonosulfate
- Accelerated hydrogen peroxide
- 70% isopropyl alcohol

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Inflate the support surface with the pump.
2. Disconnect the mains power from the receptacle.

Note - Leave the hose connected to the support surface to retain the air.

3. Clean and dry the support surface covers before disinfectants are applied.
4. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes. Do not soak the support surface.

Note - Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.

5. Rinse the product with a clean, damp cloth to remove excess disinfectant.
6. Dry the product with a clean, dry cloth to remove any excess liquid or disinfectant.
7. Allow the product to dry completely before you return to service.

Troubleshooting

CAUTION - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Service only by qualified personnel.

Problem	Cause	Recommended action
Pressure alarm indicator blinking on the Model 2940 IsoAir pump	<ol style="list-style-type: none">1. Bad CPR connection2. Bad hose connection at a fitting3. A hole in an air cell or air hose	<ol style="list-style-type: none">1. Inspect the CPR connection 1a. Listen for leaks around the CPR connection For 2. and 3., Inspect all air hoses, fittings, and air cells

Note - The A and B air circuits on the Model 2940 IsoAir pump are reverse from the A and B air cells on the 2941 IsoAir support surface (Figure 16).


Figure 16 – 2941 IsoAir support surface A and B air cells

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