## Isolibrum® Support Surface

### Operations Manual

Integrated with InTouch® FL27 (2131 / 2141 / 2151 / 2152) version 4.0 or higher

REF 2971

REF 2972

# Symbols

	Refer to instruction manual/booklet
	Operating instructions/Consult instructions for use
	General warning
	Caution
	Warning; electricity
	Alternating current
	Catalogue number
	Serial number
	Manufacturer
	Date of manufacture
	Safe working load
	Mass of equipment
	Protective earth terminal
<b>IPX4</b>	Protection from liquid splash
	Defibrillation proof type B applied part
	Medical Equipment recognized by UL LLC with respect to electric shock, fire, and mechanical hazards in accordance with ANSI/AAMI ES60601-1:2005/(R)2012, A1:2012, C1:2009/(R)2012, A2:2010/(R)2012, and CAN/CSA-C22.2 No. 60601-1 (2008)
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Wash by hand
	Do not tumble dry

	Do not dry-clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
<b>US Patents</b>	For US Patents see <a href="http://www.stryker.com/patents">www.stryker.com/patents</a>
	Pressure redistribution
	Moisture management, Low Air Loss
	Lateral rotation
	Shear
	Turn assist

# Table of Contents

Symbols .....	1
Warning/Caution/Note Definition .....	2
Summary of safety precautions .....	2
Introduction .....	4
Product description .....	4
Intended use .....	4
Expected service life .....	4
Contraindications .....	4
Specifications .....	5
Contact information .....	6
Serial number location .....	6
Date of manufacture .....	8
Product illustration .....	8
Isolibrium features .....	8
Setup .....	10
Operation .....	13
Applying the linens .....	13
Positioning a patient on the support surface .....	13
Managing incontinence and drainage .....	14
Activating and resetting CPR .....	15
Relocating the support surface .....	15
Transferring of patient from one patient support platform to another .....	16
Storing the product .....	16
Cleaning .....	17
Disinfecting .....	18
Preventive maintenance .....	19
Cover replacement, top .....	20
Accessories and parts .....	22
EMC Information .....	23

# Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

## **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

## **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

**Note** - Provides special information to make maintenance easier or important instructions clearer.

## Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

### **WARNING**

- Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
- Do not move the support surface without assistance.
- Do not operate the support surface without the **Isolibrium** top and bottom covers and linens.
- Always position the support surface on the **InTouch** bed frame to that the support surface does not contact the headboard.
- Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the **Isolibrium** support surface with **InTouch** bed frames.
- Do not use fitted sheets with this support surface.
- Do not attach the power cord to any moving parts of the **InTouch** bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
- Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
- Do not use extension cords with **Isolibrium**. It is only intended to be powered by **InTouch** with the power cord supplied.
- Always use the power source provided by the **InTouch** foot end auxiliary power outlet with **Isolibrium**.
- Always make sure that the operator has access to the CPR straps.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Always center the patient on the support surface. Align the patient's head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
- Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
- Always raise all of the **InTouch** bed siderails before starting Turn Assist or Lateral Rotation functions.
- Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
- Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.
- Always monitor the patient condition at regular intervals for patient safety.
- Do not use the support surface as a transfer device to avoid the risk of patient injury.

- Do not exceed the safe working load of the **InTouch** bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
  - Always unplug the support surface power cord before cleaning or disinfecting to avoid the risk of shock.
  - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
  - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
  - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
  - Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
  - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
  - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
  - Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.
  - Do not service or perform maintenance while the product is in use.
  - This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating **Isolibrium** or shielding the location.
- 

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
  - Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
  - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
  - Do not iron, dry-clean, or tumble dry the support surface covers.
  - Do not power wash the support surface.
  - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
  - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
  - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
  - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
-

# Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

---

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- 

## Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

## Product description

**Isolibrium**®, Stryker Model 2971/2972, is a powered support surface with features that provide Pressure Redistribution, Low Air Loss, and Lateral Rotation. **Isolibrium** is for use with the **InTouch**® (revisions 4.0 or higher software) bed frame. **Isolibrium** and **InTouch** are an integrated bed and surface system that caregivers use to assist in patient care.

## Intended use

The **Isolibrium** support surface assists in the prevention and treatment of all pressure injury stages (including stages 1, 2, 3, 4, Unstageable and Deep Tissue Pressure Injury). Before implementing **Isolibrium**, a healthcare professional should complete a clinical evaluation of risk factors and skin assessments. The **Isolibrium** patient-specific immersion to assist in managing microclimate of the patient skin, patient repositioning, early mobility, and the prevention of pulmonary complications.

The **InTouch** and **Isolibrium** system is for use by patients in an acute care setting. The acute care settings may include critical care, step down, progressive care, medical / surgical, sub-acute care, and post anesthesia care unit (PACU), or other locations as prescribed. The patient, accessories, and any extra loading must not exceed the safe working load of **InTouch** or **Isolibrium**.

The operators for the system are healthcare professionals (nurses, nurse aids, and doctors) who can use all system operations (for example: bed motion functions, nurse call, siderail communications, bed exit, or therapy options). Patients and bystanders can use the bed motion functions, nurse call, and siderail communications but are restricted from the use of the powered support surface functions. Trained professionals install, service, and calibrate the **InTouch** and **Isolibrium** system.

## Expected service life

**Isolibrium** support surface has a five year expected service life under normal use, conditions, and with appropriate periodic maintenance.

**Isolibrium** covers (top and bottom) have a two year expected service life under normal use, and conditions.

## Contraindications

Stryker promotes the clinical assessment of each patient and appropriate usage by the operator.

The **InTouch** and **Isolibrium** system is not intended:

- for use with patients 35 inches or less in height
- for use with patients that weigh less than 50 pounds
- for use with a full oxygen tent
- to support more than one individual at a time
- for use with a patient less than two years old
- in a home healthcare environment
- as a sterile product
- for use in the presence of flammable anesthetic mixture with air rich in oxygen or nitrous oxide

**Note** - When **Isolibrium** is on the **InTouch** bed, a half-oxygen tent within the sidrails is allowable. Do not enclose any part of the support surface foot box with the half-oxygen tent.

**Isolibrium** is not intended for use on patients with the following:

- unstable fractures
- unstable spinal cord injuries

**Isolibrium** Lateral Rotation and Turn Assist features are for use only after you complete a clinical assessment. This is to make sure that no adverse reaction would result in patients with the following:







- spinal cord injury
- skeletal traction
- significant hemoptysis
- bleeding disorders
- rib fractures
- for whom the techniques cause increased dyspnea or wheezing
- hemodynamic instable
- increased intracranial pressure
- skeletal traction

## Specifications

Model	2971/2972 <b>Isolibrium</b>	
Length	84 in.	213,4 cm
Width	35 in.	88,9 cm
Thickness	8 - 10.5 in.	20.3 - 26.7 cm
Weight	90 lb	40.8 kg
Therapeutic weight	50 lb to 350 lb	22.7 kg to 158.7 kg
  Safe working load <b>Isolibrium</b> (Sum of the patient and accessory weight)	460 lb	208.6 kg

Electrical requirements	ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012 CAN/CSA-C22.2 No. 60601-1 (2008)	120VAC 50/60Hz, 2.5A
Power cord	3 foot, 16 AWG cord with hospital grade plug for use with the <b>InTouch</b> foot end auxiliary power outlet	
Low Air Loss flow rate	100 LPM	
Leakage current	300 uA Maximum	
Noise level	50 dBA maximum	

**Stryker** reserves the right to change specifications without notice.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

See the **InTouch** manual for bed specifications.

## Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

## Serial number location


The serial number (A) for your support surface is located on the patient right side of the foot box on the specification label. Unzip the support surface top cover (see *Cover replacement, top* (page 20)).

## Date of manufacture

The year of manufacture is the first four digits of the serial number.

## Product illustration

A	Head end
B	Foot end
C	D-rings
D	CPR straps
E	Power cord

F	Integration cable
G	Handles
H	Top support surface cover
J	Bottom support surface cover
K	Retainer straps (D-rings not shown)

To view the power cord label, the operator must position the cord by rotating it until the label is legible from a distance of up to 1m.

## Isolibrium features

Physical Concepts Pressure Redistribution

The ability of a support surface to distribute load over the contact areas of the human body.

Pressure Redistribution setpoint

Allows the operator to input patient weight range resulting in automatic change of zoned sensor technology allowing immersion and envelopment.

Lateral Rotation Therapy	A feature of a support surface that provides rotation about a longitudinal axis as characterized by degree of patient turn, duration, and frequency.
<b>BackSmart®</b> Turn Assist	A feature of a support surface that provides a single position turn about a longitudinal axis.
Max Inflate	Allows maximum firmness inflation for use when boosting a patient and on lateral transfer, Promotes patient ease with dangling, egress, and ingress.
Moisture management	The top surface cover draws excess moisture away from the patients skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to assist in the management of (microclimate) the skin.
<b>HeelGel®<sup>1</sup></b> Technology	A feature of the support surface to address pressure redistribution in the heel area.
<b>Patient Specific Immersion™</b>	Assists in managing microclimate of the patient skin, patient repositioning, early mobility, and prevention of pulmonary complications.

<sup>1</sup>**Intelli-Gel®** hollow column configuration licensed by EdiZONE LLC of Alpine, UT

# Setup

The product must reach the ambient operating temperature before installation or operation to avoid the risk of permanent damage to this product. Upon initial installation, reset the CPR.

---

## WARNING


- Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
- Do not move the support surface without assistance.
- Do not operate the support surface without the **Isolibrium** top and bottom covers and linens.
- Always position the support surface on the **InTouch** bed frame to that the support surface does not contact the headboard.
- Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the **Isolibrium** support surface with **InTouch** bed frames.
- Do not use fitted sheets with this support surface.
- Do not attach the power cord to any moving parts of the **InTouch** bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
- Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
- Do not use extension cords with **Isolibrium**. It is only intended to be powered by **InTouch** with the power cord supplied.
- Always use the power source provided by the **InTouch** foot end auxiliary power outlet with **Isolibrium**.

---

**Note** - Do not plug into the power source until the assembly is complete.

**Prerequisite:** Minimum of two operators are required to install the support surface.


1. Put the **InTouch** bed into the flat position before installing the support surface.
2. If powered, unplug the **InTouch** from the wall outlet and turn the battery switch to off.
3. Remove the **InTouch** headboard and footboard from the bed.
4. Position an operator at each end of the support surface and lift onto the **InTouch** bed frame.
5. Thread the retainer straps (A), down through the **InTouch** Fowler retainer holes (B) (Figure 1).

**Figure 1 – Insert retainer straps**

**Note** - The reverse side of the support surface is shown for reference. Make sure that the support surface is installed with the top cover facing up.

6. Put the retainer strap through both D-Rings (C), both sides.
7. Thread the retainer strap through the bottom D-Ring, both sides.
8. Pull the retainer strap until the loop around the Fowler is tight, both sides.
9. Insert the loose end of the retainer straps behind the Fowler.
10. Reinstall the **InTouch** headboard and footboard.
11. Position the support surface against the footboard. Allow a small air gap at the headboard.
12. Connect the integration cable from the support surface foot end to the **InTouch** cable adapter (B) (Figure 2).
13. Turn the locking collars clockwise to secure the connection.
14. Place the linens on the support surface as loose as possible to avoid creating a **hammocking** effect and increasing interface pressures.
15. Plug the power cord for the support surface into the **InTouch** foot end auxiliary power outlet (A) (Figure 2).

**Figure 2 – InTouch foot end auxiliary power outlet and mattress connector**

16. Make sure that the power cord is routed to the side of the bed and not out the foot end of the **InTouch**.
17. Plug in the **InTouch** power cord into a properly grounded hospital grade three-prong power outlet.
18. Calibrate the bed following the procedures in the **InTouch** maintenance manual.

**Note**

- **Isolibrium** is only operational when the power cord for the **InTouch** bed is in a 120V wall outlet.
- When the support surface is plugged into the **InTouch** bed, the **InTouch** footboard will respond. See the **InTouch** manual for more information.

# Operation

## Applying the linens

---

### WARNING

- Always make sure that the operator has access to the CPR straps.
  - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
  - Do not use fitted sheets with this support surface.
- 

To apply the linens:

1. Apply the linens using the “D” rings for the flat sheet.
2. Thread the four linen corners through the “D” rings (C) attached to the bottom cover to secure the linens (A) to the support surface (B).


Figure 3 – Apply linens

**Note** - Do not pull the linens tight. Keep the linens loose and as smooth as possible on top of the support surface to effectively use the Turn Assist and Lateral Rotation functions.

## Positioning a patient on the support surface

---

### WARNING

- Always center the patient on the support surface. Align the patient’s head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
  - Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
  - Always raise all of the **InTouch** bed siderails before starting Turn Assist or Lateral Rotation functions.
  - Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
  - Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.
- 

### CAUTION

- Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
  - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
-

To position the patient:

1. Max Inflate the support surface.
2. Position the patient in the center of the support surface, align the patient's head toward the head board (Figure 4).
3. Check the patient frequently during Lateral Rotation for proper positioning and support surface inflation (Figure 5 and Figure 6).


Figure 4 – Center the patient


Figure 5 – Correct patient positioning


Figure 6 – Incorrect patient positioning

## Managing incontinence and drainage

---

**WARNING** - Always monitor the patient condition at regular intervals for patient safety.

---

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

## Activating and resetting CPR

When you activate for CPR, the system deflates and the InTouch control display will change.

---

**WARNING** - Always make sure that the operator has access to the CPR straps.

---

To activate for CPR, pull the red CPR strap straight out from the support surface (Figure 7). The CPR straps are located on the left and right side of the foot end of the support surface. You can activate either one or both CPR straps.

Figure 7 – CPR activation

Figure 8 – CPR reset

To reset for CPR, pull down on red CPR strap (Figure 8). The support surface will begin to reinflate.

## Relocating the support surface

Relocating the support surface is for the purpose of cleaning or moving to another InTouch bed frame.

---

### WARNING

- Do not move the support surface without assistance.
  - Do not use the support surface as a transfer device to avoid the risk of patient injury.
- 

**Prerequisite:** Minimum of two operators are required to relocate the support surface and setup this product.

1. Make sure that there is not a patient on the support surface.
2. Remove the linens.
3. Unplug the power cord for the support surface from the foot end auxiliary power outlet.
4. Unplug the support surface connector.
5. Remove the headboard and footboard
6. Position an operator at each end of the support surface and remove the support surface and relocate.

**Note** - Use all four handles when relocating the support surface.

**Note** - See (*Setup* (page 10)) of the support surface.

7. Reinstall the InTouch headboard and footboard.

# Transferring of patient from one patient support platform to another

---

## WARNING

- Do not use the support surface as a transfer device to avoid the risk of patient injury.
  - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
  - Do not exceed the safe working load of the **InTouch** bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- 

To transfer the patient from one patient support surface to another:

**Prerequisite:** Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes to on for both patient support platforms.
3. Activate maximum inflation of the support surface until the inflation level is reached. See the **InTouch** Starting and Stopping Max Inflate section of the manual for use instructions.
4. Adjust the patient support platform heights so that they are level with one another.
5. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.
6. Turn off Max Inflate.

## Storing the product

Before storing the product, deflate the support surface by activating CPR.

# Cleaning

---

## WARNING

- Always unplug the support surface power cord before cleaning or disinfecting to avoid the risk of shock.
  - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
  - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
  - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
  - Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
  - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
- 

## CAUTION

- Do not iron, dry-clean, or tumble dry the support surface covers.
  - Do not power wash the support surface.
  - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
- 

**Prerequisite:** Minimum of two operators are required to clean the bottom support surface cover.

Always follow hospital protocol for cleaning and disinfecting.

To clean the support surface covers between patient uses, follow these steps in order:

1. If powered, unplug the **InTouch** from the wall outlet and turn the battery switch to off.
2. Unplug the support surface power cord.
3. Disconnect the support surface integration cable.
4. Using a clean, soft, damp cloth, wipe the support surface covers with a mild soap and water solution to remove foreign material.
5. Wipe the support surface covers with a clean, dry cloth to remove any excess liquid or cleaning agent.

**Note** - To clean the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.

6. Thoroughly rinse and dry the support surface covers.

# Disinfecting

---

## WARNING

- Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
  - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
  - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
  - Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.
- 

## CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
  - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
  - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
- 

**Prerequisite:** Minimum of two operators are required to disinfect the support surface.

Suggested disinfectants:

- Quaternaries
- Phenolic disinfectants
- Chlorinated bleach solution (5.25% bleach diluted 1 part bleach to 100 parts water)
- Chlorinated bleach solution for 2972 Endurance cover is 10,000ppm
- 70% Isopropyl Alcohol

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Unplug the InTouch power cord.
2. Unplug the support surface.
3. Disconnect the integration cable.
4. Thoroughly clean and dry the support surface covers (see *Cleaning* (page 17)) before disinfectants are applied.
5. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes (do not soak the support surface).

### Note

- Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.
  - To disinfect the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.
6. Rinse the support surface covers with a clean, damp cloth to remove excess disinfectant.
  7. Wipe the support surface covers with a clean, dry cloth to remove any excess liquid or disinfectant.
  8. Allow the support surface covers to dry completely before returning to service.

# Preventive maintenance

**WARNING** - Do not service or perform maintenance while the product is in use.

---

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage. Service only by qualified personnel.

Remove product from service before you perform preventive maintenance inspection.

**Note** - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

\_\_\_\_\_ Zipper, snaps, and covers (top, bottom, and fire barrier) are free of tears, cuts, holes, or any other damage

**Note** - If excessive wear is observed on the cover or the fire barrier, it is strongly recommended to replace the cover.

\_\_\_\_\_ Support surface cover labels are legible, have proper adherence, and free of damage.

\_\_\_\_\_ Handles are free of rips or cracks

\_\_\_\_\_ Power cord, integration cable, and plug are free of damage

\_\_\_\_\_ \*\*Perform a functional test of the Isolibrium functions

\_\_\_\_\_ \*\*Perform system diagnostics and confirm there are no errors

\_\_\_\_\_ \*Internal components for signs of staining from fluid ingress or contamination by fully unzipping the cover

\_\_\_\_\_ Pods are not leaking or cracking

\_\_\_\_\_ Hose connections to the manifold and sensor tubes (both ends) are seated fully

\_\_\_\_\_ Foam is free from large tears or large gouges

\_\_\_\_\_ Turn bladders are free of excessive wear and function properly

\_\_\_\_\_ Low Air Loss fans rotate freely and are free of debris and dust

\_\_\_\_\_ Left and Right CPR releases function properly

\_\_\_\_\_ Both foot box cooling fans rotate freely and are free of debris and dust

\_\_\_\_\_ Pump spring mount has visible clearance from the foot box bottom

\_\_\_\_\_ Listen for abnormal pump vibration and noise level

\_\_\_\_\_ Clean the foot box intake metal screen

\_\_\_\_\_ Clean the mesh on the bottom cover

\_\_\_\_\_ Clean the fan filter screen on the foot box


\_\_\_\_\_ Standoffs (feet of the foot box) are present

\_\_\_\_\_ Cable tie is present on the power cord

\_\_\_\_\_ Current leakage not more than 200  $\mu$ A

\_\_\_\_\_ Ground impedance not more than 100 m $\Omega$  (Figure 9)

**Note** - The ground lug can be accessed through the slit in the cover for the CPR, patient left.

**Figure 9 – Ground lug (patient left side)**

**Note**

- \*A zipper pull tool or equivalent is required to access the internal components of the support surface.
- \*\*You cannot perform these preventive maintenance checks without the **InTouch** footboard. See the **InTouch** integrated Models 2131 and 2141, software version 4.0 or higher manual.
- The zipper overlaps at the end and may appear to be misaligned.

Replace the following items every two years:

- Cover assembly, top and bottom

Product serial number:
Completed by:
Date:

**Cover replacement, top**


**Tools required:**

- Zipper pull tool or equivalent

**Procedure:**

1. Apply the **InTouch** brakes.
2. Unplug the support surface power cord.
3. Remove the headboard and footboards.
4. Lower all the siderails.
5. Unsnap the two secure snaps (C) at the foot end corners (Figure 10).
6. Unsnap the four corner retainers (B) (Figure 10).
7. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
8. Unzip the support surface top cover from the bottom using the zipper pull tool or equivalent (A) (Figure 10).

**Note** - The zipper starts on the patient right side near the foot end.

**Figure 10 – Top cover**

9. Remove and discard the top cover.

10. Reverse to install the new cover.

**Note**

- After installation, remove the zipper pull tool or equivalent from the zipper.
- Connect the corner retainers by placing the strap below each D-ring and snaps to the cover. Repeat for the other three corners.
- Make sure that the watershed covers the zipper.

11. Verify proper operation before returning the product to service.

## Accessories and parts

These parts are currently available for purchase. Call Stryker Customer Service: 1-800-327-0770 for availability and pricing.

Part name	Part number
Board, main	2971-022-140
Board, sensor	2971-021-033
Cover, bottom	2971-700-008
Cover, top	2971-021-004
Cover, top (Alternative <b>Dartex</b> )	2972-021-004
Fan, foot box	2971-022-188
Fan, LAL	2971-021-189
Integration cable	2971-022-133
Kit, bottom cover	2972-021-018
Kit, power cord	2971-700-003
Kit, pump assembly	2971-700-001
Kit, resonator assembly	2971-700-002
Kit, valve manifold	2971-700-007
Pods assembly	2971-021-003
Power supply	2971-022-171
Kit, solenoid valve, SLI, BLI, BLD, BRI, BRD, SRD, SRI, SLD	2971-700-004
Kit, solenoid valve, TRI, TRD, EX, TRI, TU	2971-700-005
Kit, solenoid valve, IN	2971-700-006

# EMC Information

## Guidance and Manufacturer's declaration - Electromagnetic Immunity

**Isolibrium** is suitable for use in the electromagnetic environment specified below. The customer or the user of **Isolibrium** should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/ burst IEC 61000-4-4 *	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5 *	+1 kV lines to lines +2 kV lines to earth	+1 kV lines to lines +2 kV lines to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11 *	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Main power quality should be that of a typical commercial and/or hospital environment. If the user of <b>Isolibrium</b> requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.

**Note** -  $U_T$  is the a.c. mains voltage before applications of the test level.

**Recommended separation distances between portable and mobile RF communications equipment and Isolibrium.**

**Isolibrium** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of **Isolibrium** can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and **Isolibrium** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <b>W</b>	Separation distance according to frequency of transmitter <b>m</b>		
	150 kHz to 80 MHz $d=[3.5/\sqrt{f}] \sqrt{P}$	80 MHz to 800 MHz $d=[3.5/\sqrt{f}] \sqrt{P}$	800 MHz to 2.5 GHz $d=[7/\sqrt{f}] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note** - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note** - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Isolibrium** is suited for use in the electromagnetic environment specified below. The customer or the user of **Isolibrium** should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
---------------	----------------------	------------------	--

<p>Conducted RF IEC 61000-4-6 *</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of <b>Isolibrium</b>, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended Separation Distance</p> $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ <p>80 MHz to 800 MHz</p> $d=2.3\sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	----------------------------	---

**Note** - At 80 MHz and 800 MHz, the higher frequency range applies.

**Note** - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **Isolibrium** is used exceeds the applicable RF compliance level above, **Isolibrium** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **Isolibrium**.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

**Guidance and Manufacturer's declaration - Electromagnetic Emissions**

**Isolibrium** is intended for use in an electromagnetic environment specified below. The customer or the user of **Isolibrium** should make sure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	<b>Isolibrium</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	<b>Isolibrium</b> is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

**WARNING** - This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating **Isolibrium** or shielding the location.

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