Isolibrium Support Surface

Operations Manual

REF 2972

Integrated with InTouch FL27
(2131/2141/2151/2152)

Version 4.0 or higher

Global symbol glossary

See the Global Symbol Glossary at ifu.stryker.com for symbol definitions.

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Warning; electricity
	Alternating current
	Catalogue number
	Serial number
	Unique device identifier
	Quantity
	Manufacturer
	Date of manufacture
	Safe working load
	Mass of equipment
	Protective earth terminal
IPX4	Protection from liquid splash
	Defibrillation proof type B applied part

	Medical Equipment recognized by UL LLC with respect to electric shock, fire, and mechanical hazards in accordance with ANSI/AAMI ES60601-1:2005/(R)2012, A1:2012, C1:2009/(R)2012, A2:2010/(R)2012, and CAN/CSA-C22.2 No. 60601-1:14.
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Wash by hand
	Do not tumble dry
	Do not dry-clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
	For US Patents see www.stryker.com/patents
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package
	Center of gravity/weight distribution
	Two person lift
	Pressure redistribution
	Moisture management, Low Air Loss

	Lateral rotation
	Shear
	Turn assist

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Indications for use	4
Contraindications	4
Expected service life	4
Specifications	4
European REACH	5
Product illustration	6
Isolibrium features	6
Contact information	7
Serial number location	7
Setup	8
Operation	11
Applying the linens	11
Positioning a patient on the support surface	11
Managing incontinence and drainage	12
Activating and resetting CPR	13
Relocating the support surface	13
Transferring a patient from one patient support platform to another	14
Storing the product	14
Accessories and parts	15
Cleaning	16
Disinfecting	17
Preventive maintenance	18
EMC Information	20

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
- Do not move the support surface without assistance.
- Do not operate the support surface without the Isolibrium top and bottom covers and linens.
- Always position the support surface on the InTouch bed frame so that the support surface does not contact the headboard.
- Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the Isolibrium support surface with InTouch bed frames.
- Do not use fitted sheets with this support surface.
- Do not attach the power cord to any moving parts of the InTouch bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
- Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
- Do not use extension cords with Isolibrium. It is only intended to be powered by InTouch with the power cord supplied.
- Always use the power source provided by the InTouch foot end auxiliary power outlet with Isolibrium.
- Always make sure that the operator has access to the CPR straps.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Always center the patient on the support surface. Align the patient's head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
- Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
- Always raise all of the bed siderails before you begin Turn Assist or Lateral Rotation.
- Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
- Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.
- Always monitor the patient condition at regular intervals for patient safety.
- Do not use the support surface as a transfer device to avoid the risk of patient injury.
- Do not exceed the safe working load of the InTouch bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always unplug the support surface power cord before cleaning or disinfecting to avoid the risk of shock.
- Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
- Do not allow fluid to pool on the support surface. Fluid may cause degradation of components and may cause unpredictable safety and performance of this product.
- Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
- Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
- Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
- Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
- Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.

- Do not service or perform maintenance while the product is in use.
 - This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating Isolibrium or shielding the location.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
 - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
 - Do not iron, dry clean, or tumble dry the support surface covers.
 - Do not power wash the support surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
 - This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

Stryker Model 2972 Isolibrium is a powered support surface with features that provide Pressure Redistribution, Low Air Loss, and Lateral Rotation. Isolibrium is for use with the InTouch (revisions 4.0 or higher software) bed frame. Isolibrium and InTouch are an integrated bed and surface system that caregivers use to assist in patient care.

Indications for use

The Stryker Model 2972 Isolibrium support surface is intended to assist in the prevention and treatment of all pressure injury or pressure ulcers (includes all stages, unstageable and deep tissue pressure injury). As a recommendation, healthcare professional should complete a clinical evaluation of risk factors and skin assessments. The intended users is both healthcare HCPs (nurses, nurse aids, and medical doctors), transporters, and human patients. This product is for use by human patients in a health delivery organization (HDO).

The Isolibrium patient-specific immersion is intended to assist in the management of microclimate of the patient's skin, patient repositioning, and the support in prevention of early mobility and pulmonary complications.

The Isolibrium powered support surface is not intended for use in or with:

- Oxygen rich environments
- Sterile environments
- Home care or long-term care facility settings
- Behavioral health populations

Contraindications

Isolibrium is contraindicated for use on patients with:

- Unstable fractures
- Unstable spinal cord injuries

Expected service life

Isolibrium support surface has a five year expected service life under normal use and conditions and with appropriate periodic maintenance.






Isolibrium covers (top and bottom) have a two year expected service life under normal use and conditions.

Specifications

  Safe working load Note - Safe working load indicates the sum of the patient and accessory weight.	460 lb	208.6 kg
Therapeutic weight	50 lb to 350 lb	22.7 kg to 158.7 kg
Length	84 in.	213.4 cm
Width	35 in.	88.9 cm
Thickness	8 - 10.5 in.	20.3 - 26.7 cm
Weight	90 lb	40.8 kg

Electrical requirements	ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012 CAN/CSA-C22.2 No. 60601-1:14	120VAC 50/60Hz, 2.5A
Power cord	3 foot, 16 AWG cord with hospital grade plug for use with the InTouch foot end auxiliary power outlet	
Low Air Loss flow rate	100 LPM	
Leakage current	300 uA Maximum	
Noise level	50 dBA maximum	

Stryker reserves the right to change specifications without notice.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		


See the InTouch Critical Care Bed Operations Manual for bed specifications.

European REACH

In accordance with the European REACH regulation and other environmental regulatory requirements, the components that contain declarable substances are listed.

Description	Number	Substance of very high concern (SVHC) chemical name
Sensor assembly	2971-021-033	1,3,5-Tris(oxiran-2-ylmethyl)-1,3,5-triazinane-2,4,6-trione (TGIC)
		2-Ethoxyethyl acetate
		Cadmium
		Diboron-trioxide
		Lead-monoxide
Main board PCB assembly	2971-022-140	Cadmium
		Diboron-trioxide
		Lead-monoxide
Footbox bottom weldment	2971-022-009	Tris(2,3-epoxypropyl) isocyanurate
Footbox top	2971-022-200	Tris(2,3-epoxypropyl) isocyanurate
Pump assembly	297300220010	Lead-steel
		Lead-brass
		Lead-aluminum

Product illustration

A	Head end
B	Foot end
C	D-rings
D	CPR straps
E	Power cord

F	Integration cable
G	Handles
H	Top support surface cover
J	Bottom support surface cover
K	Retainer straps (D-rings not shown)

To view the power cord label, the operator must position the cord by rotating it until the label is legible from a distance of up to 1m.

Isolibrium features

Physical Concepts Pressure Redistribution	The ability of a support surface to distribute load over the contact areas of the human body.
Pressure Redistribution setpoint	Allows the operator to input patient weight range resulting in automatic change of zoned sensor technology allowing immersion and envelopment.
Lateral Rotation Therapy	A feature of a support surface that provides rotation about a longitudinal axis as characterized by degree of patient turn, duration, and frequency.
BackSmart Turn Assist	A feature of a support surface that provides a single position turn about a longitudinal axis.
Max Inflate	Allows maximum firmness inflation for use when boosting a patient and on lateral transfer, Promotes patient ease with dangling, egress, and ingress.
Moisture management	The top surface cover draws excess moisture away from the patients skin and support surface interface.

Low Air Loss (LAL)

Provides air flow to assist in the management of (microclimate) the skin.

HeelGel¹ Technology

A feature of the support surface to address pressure redistribution in the heel area.

Patient Specific Immersion

Assists in managing microclimate of the patient skin, patient repositioning, early mobility, and prevention of pulmonary complications.

¹Intelli-Gel® hollow column configuration licensed by EdiZONE LLC of Alpine, UT

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

The serial number (A) for your support surface is located on the patient right side of the foot box on the specification label. Unzip the support surface top cover to access the specification label. The serial number (B) is for the support surface cover.

Setup

The product must reach the ambient operating temperature before installation or operation to avoid the risk of permanent damage to this product. Upon initial installation, reset the CPR.

WARNING

- Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
 - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
 - Do not move the support surface without assistance.
 - Do not operate the support surface without the Isolibrium top and bottom covers and linens.
 - Always position the support surface on the InTouch bed frame so that the support surface does not contact the headboard.
 - Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the Isolibrium support surface with InTouch bed frames.
 - Do not use fitted sheets with this support surface.
 - Do not attach the power cord to any moving parts of the InTouch bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
 - Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
 - Do not use extension cords with Isolibrium. It is only intended to be powered by InTouch with the power cord supplied.
 - Always use the power source provided by the InTouch foot end auxiliary power outlet with Isolibrium.
-

Note - Do not plug into the power source until the assembly is complete.

Prerequisite: A minimum of two operators is required to install the support surface.

1. Put the InTouch bed into the flat position before installing the support surface.
2. If powered, unplug the InTouch from the wall outlet and turn the battery switch to off.
3. Remove the InTouch headboard and footboard from the bed.
4. Position an operator at each end of the support surface and lift onto the InTouch bed frame.
5. Thread the retainer straps (A), down through the InTouch Fowler retainer holes (B) (Figure 1).


Figure 1 – Insert retainer straps

Note - The reverse side of the support surface is shown for reference. Make sure that the support surface is installed with the top cover facing up.

6. Put the retainer strap through both D-Rings (C), both sides.
7. Thread the retainer strap through the bottom D-Ring, both sides.
8. Pull the retainer strap until the loop around the Fowler is tight, both sides.
9. Insert the loose end of the retainer straps behind the Fowler.
10. Reinstall the InTouch headboard and footboard.
11. Position the support surface against the footboard. Allow a small air gap at the headboard.
12. Connect the integration cable from the support surface foot end to the InTouch cable adapter (B) (Figure 2).
13. Turn the locking collars clockwise to secure the connection.
14. Place the linens on the support surface as loose as possible to avoid creating a hammocking effect and increasing interface pressures.
15. Plug the power cord for the support surface into the InTouch foot end auxiliary power outlet (A) (Figure 2).


Figure 2 – InTouch foot end auxiliary power outlet and mattress connector

16. Make sure that the power cord is routed to the side of the bed and not out the foot end of the InTouch.

17. Plug in the InTouch power cord into a properly grounded hospital grade three-prong power outlet.

18. Calibrate the bed following the procedures in the InTouch maintenance manual.

Note

- Isolibrium is only operational when the power cord for the InTouch bed is in a 120V wall outlet.
- When the support surface is plugged into the InTouch bed, the InTouch footboard will respond. See the InTouch Critical Care Bed Operations Manual for more information.

Operation

Applying the linens

WARNING

- Always make sure that the operator has access to the CPR straps.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not use fitted sheets with this support surface.
-

To apply the linens:

1. Apply the linens using the “D” rings for the flat sheet.
2. Thread the four linen corners through the “D” rings (C) attached to the bottom cover to secure the linens (A) to the support surface (B).


Figure 3 – Apply linens

Note - Do not pull the linens tight. Keep the linens loose and as smooth as possible on top of the support surface to effectively use the Turn Assist and Lateral Rotation functions.

Positioning a patient on the support surface

WARNING

- Always center the patient on the support surface. Align the patient's head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
 - Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
 - Always raise all of the bed siderails before you begin Turn Assist or Lateral Rotation.
 - Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
 - Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.
-

CAUTION

- Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
 - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
-

To position the patient:

1. Max Inflate the support surface.
2. Position the patient in the center of the support surface, align the patient's head toward the head board (Figure 4).
3. Check the patient frequently during Lateral Rotation for proper positioning and support surface inflation (Figure 5 and Figure 6).


Figure 4 – Center the patient


Figure 5 – Correct patient positioning


Figure 6 – Incorrect patient positioning

Managing incontinence and drainage

WARNING - Always monitor the patient condition at regular intervals for patient safety.

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

Activating and resetting CPR

When you activate for CPR, the system deflates and the InTouch control display changes.

WARNING - Always make sure that the operator has access to the CPR straps.

To activate for CPR, pull the red CPR strap straight out from the support surface (Figure 7). The CPR straps are located on the left and right side of the foot end of the support surface. You can activate either one or both CPR straps.

Figure 7 – CPR activation

Figure 8 – CPR reset

To reset for CPR, pull down on red CPR strap (Figure 8). The support surface reinflates.

Relocating the support surface

Relocating the support surface is for the purpose of cleaning or moving to another InTouch bed frame.

WARNING

- Do not move the support surface without assistance.
 - Do not use the support surface as a transfer device to avoid the risk of patient injury.
-

Prerequisite: A minimum of two operators is required to relocate the support surface and set up this product.

1. Make sure that there is not a patient on the support surface.
2. Remove the linens.
3. Unplug the power cord for the support surface from the foot end auxiliary power outlet.
4. Unplug the support surface connector.
5. Remove the headboard and footboard
6. Position an operator at each end of the support surface and remove the support surface and relocate.

Note - Use all four handles when relocating the support surface.

Note - See (*Setup* (page 8)) of the support surface.

7. Reinstall the InTouch headboard and footboard.

Transferring a patient from one patient support platform to another

WARNING

- Do not use the support surface as a transfer device to avoid the risk of patient injury.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not exceed the safe working load of the InTouch bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
-

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes to on for both patient support platforms.
3. Activate maximum inflation of the support surface until the inflation level is reached. See the InTouch Critical Care Bed Operations Manual, *Starting and stopping Max Inflate*.
4. Adjust the patient support platform heights so that they are level with one another.
5. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.
6. Turn off Max Inflate.

Storing the product

Before storing the product, deflate the support surface by activating CPR.

Accessories and parts

These parts are currently available for purchase. Call Stryker Customer Service: 1-800-327-0770 for availability and pricing.

Name	Number
Board, main	2971-022-140
Board, sensor	2971-021-033
Cover, bottom assembly	2972-021-017
Cover, top assembly	2972-021-004
Fan, foot box	2971-022-188
Fan, LAL	2971-021-189
Integration cable	2971-022-133
Kit, power cord	2971-700-003
Kit, pump assembly	297307000001
Kit, resonator assembly	2971-700-002
Kit, solenoid valve, IN	2971-700-006
Kit, solenoid valve, SLI, BLI, BLD, BRI, BRD, SRD, SRI, SLD	2971-700-004
Kit, solenoid valve, TRI, TRD, EX, TRI, TU	2971-700-005
Kit, valve manifold	2971-700-007
Pods assembly	297300210003
Power supply	2971-022-171

Cleaning

WARNING

- Always unplug the support surface power cord before cleaning or disinfecting to avoid the risk of shock.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface. Fluid may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
-

CAUTION

- Do not iron, dry clean, or tumble dry the support surface covers.
 - Do not power wash the support surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
-

Prerequisite: A minimum of two operators is required to clean the bottom support surface cover.

Always follow hospital protocol for cleaning and disinfecting.

To clean the support surface covers between patient uses, follow these steps in order:

1. If powered, unplug the InTouch from the wall outlet and turn the battery switch to off.
2. Unplug the support surface power cord.
3. Disconnect the support surface integration cable.
4. Using a clean, soft, damp cloth, wipe the support surface covers with a mild soap and water solution to remove foreign material.
5. Wipe the support surface covers with a clean, dry cloth to remove any excess fluid or cleaning agent.

Note - To clean the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.

6. Thoroughly rinse and dry the support surface covers.

Disinfecting

WARNING

- Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluid may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
-

Prerequisite: A minimum of two operators is required to disinfect the support surface.

Suggested disinfectants:

- Quaternaries
- Phenolic disinfectants
- Chlorinated bleach solution (5.25% bleach diluted 1 part bleach to 100 parts water)
- Chlorinated bleach solution for 2972 Endurance cover is 10,000ppm
- 70% Isopropyl Alcohol

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Unplug the InTouch power cord.
2. Unplug the support surface.
3. Disconnect the integration cable.
4. Thoroughly clean and dry the support surface covers (see *Cleaning* (page 16)) before disinfectants are applied.
5. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes (do not soak the support surface).

Note

- Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.
 - To disinfect the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.
6. Rinse the support surface covers with a clean, damp cloth to remove excess disinfectant.
 7. Wipe the support surface covers with a clean, dry cloth to remove any excess fluid or disinfectant.
 8. Allow the support surface covers to dry completely before returning to service.

Preventive maintenance

WARNING - Do not service or perform maintenance while the product is in use.

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage. Service only by qualified personnel.

Remove product from service before you perform preventive maintenance inspection.

Note - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

_____ Zipper, snaps, and covers (top, bottom, and fire barrier) are free of tears, cuts, holes, or any other damage

Note - If excessive wear is observed on the cover or the fire barrier, it is strongly recommended to replace the cover.

_____ Support surface cover labels are legible, have proper adherence, and are free of damage

_____ Handles are free of rips or cracks

_____ Power cord, integration cable, and plug are free of damage

_____ **Perform a functional test of the Isolibrium functions

_____ **Perform system diagnostics and confirm there are no errors

_____ *Internal components for signs of staining from fluid ingress or contamination by fully unzipping the cover

_____ Pods are not leaking or cracking

_____ Hose connections to the manifold and sensor tubes (both ends) are seated fully

_____ Foam is free from large tears or large gouges

_____ Turn bladders are free of excessive wear and function properly

_____ Low Air Loss fans rotate freely and are free of debris and dust

_____ Left and Right CPR releases function properly

_____ Both foot box cooling fans rotate freely and are free of debris and dust

_____ Pump spring mount has visible clearance from the foot box bottom

_____ Listen for abnormal pump vibration and noise level

_____ Clean the foot box intake metal screen

_____ Clean the mesh on the bottom cover

_____ Clean the fan filter screen on the foot box

_____ Standoffs (feet of the foot box) are present

_____ Cable tie is present on the power cord

_____ Current leakage not more than 200 μ A

_____ Ground impedance not more than 100 m Ω (Figure 9)

Note - The ground lug can be accessed through the slit in the cover for the CPR, patient left.


Figure 9 – Ground lug (patient left side)

Note

- *A zipper pull tool or equivalent is required to access the internal components of the support surface.
- **You cannot perform these preventative maintenance checks without the user interface of an integrated bed's footboard. See the bed's maintenance manual to access the service menu on the footboard.
- The zipper overlaps at the end and may appear to be misaligned.

Replace the following items every two years:

- Cover assembly, top and bottom

Product serial number:
Completed by:
Date:

EMC Information

WARNING - This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating Isolibrium or shielding the location.

CAUTION - This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
Isolibrium is suitable for use in the electromagnetic environment specified below. The customer or the user of Isolibrium should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/ burst IEC 61000-4-4 *	± 2 kV for power supply lines ± 1 kV for input/ output lines	± 2 kV for power supply lines ± 1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5 *	± 0.5 kV, ± 1 kV lines to lines ± 0.5 kV, ± 1 kV, ± 2 kV lines to earth	± 0.5 kV, ± 1 kV lines to lines ± 0.5 kV, ± 1 kV, ± 2 kV lines to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11 *	$<5\%$ U_T ($>95\%$ dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles $<5\%$ U_T ($>95\%$ dip in U_T) for 5 sec	$<5\%$ U_T ($>95\%$ dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles $<5\%$ U_T ($>95\%$ dip in U_T) for 5 sec	Main power quality should be that of a typical commercial and/or hospital environment. If the user of Isolibrium requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.
Note - U_T is the a.c. mains voltage before applications of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and Isolibrium.
Isolibrium is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Isolibrium can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Isolibrium as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=[3.5/\sqrt{P}]$	80 MHz to 800 MHz $d=[3.5/\sqrt{P}]$	800 MHz to 2.5 GHz $d=[7/\sqrt{P}]$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Isolibrium is suited for use in the electromagnetic environment specified below. The customer or the user of Isolibrium should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6 * Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of Isolibrium, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended Separation Distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

Note - At 80 MHz and 800 MHz, the higher frequency range applies.

Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Isolibrium is used exceeds the applicable RF compliance level above, Isolibrium should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Isolibrium.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

Guidance and Manufacturer's declaration - Electromagnetic Emissions

Isolibrium is intended for use in an electromagnetic environment specified below. The customer or the user of Isolibrium should make sure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	Isolibrium uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	Isolibrium is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **BackSmart, HeelGel, InTouch, Isolibrium, Stryker**. All other trademarks are trademarks of their respective owners or holders.