

REF NUMBER  
2060-000-001

Manufactured for:  
Stryker Medical  
3800 E. Centre Ave.  
Portage, MI 49002  
Toll-free: 1-800-STRYKER  
EPA REG. NO. 70144-2-92728  
EPA EST. 74058-CA-001,  
74058-OH-001

stryker

# SideKick®

DISINFECTING WIPES

TUBERCULOCIDAL - BACTERICIDAL  
FUNGICIDAL - VIRUCIDAL\*



ACTIVE INGREDIENTS:		
n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride	.....	0.154%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride	.....	0.154%
Isopropanol	.....	21.000%
OTHER INGREDIENTS:	.....	78.692%
TOTAL:	.....	100.000%

KEEP OUT OF REACH OF CHILDREN

**CAUTION** - SEE SIDE PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS  
Net Contents **100 WIPES 6x10 INCHES**  
Net Weight 1 lb 4 oz (567.3g)

These wipes are disposable, single-use towelette wipes pre-saturated with disinfectant solution. This product is intended for use on pre-cleaned, hard, non-porous surfaces in areas where the disinfection of surfaces is required.

This product can be used in the following locations: ambulances, bathrooms, EMS units, extended care facilities, field hospitals, GI laboratories, hospitals, healthcare facilities, human transport vehicles, long-term care facilities, medical offices, prisons, radiology centers, surgical<sup>^</sup> centers and areas/facilities concerned with the hazards of cross-contamination from infectious microorganisms. This product may be used on the following surfaces: patient care equipment surfaces, EMS equipment/work surfaces, infant care equipment (bassinet, crib, infant incubator surfaces), laboratory equipment/surfaces, OR surfaces, nursing station surfaces, patient bed surfaces, patient stretchers, patient transport surfaces, stair chairs, wheel chairs and devices/surfaces potentially contaminated with disease-causing microorganisms and use sites made of glass, metal, Plexiglas®, plastics (i.e., polycarbonate, polypropylene, polyvinylchloride and polystyrene), clear plastic, baked-on painted surfaces, stainless steel, brass, chrome, copper, acrylic, non-porous vinyl, powder painted coated surfaces, Formica®, glazed tile, laminated surfaces and other similar surfaces.

#### Directions For Use

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. **DO NOT USE THIS PRODUCT AS A BABY OR SKIN WIPE, OR ON CONTACT LENSES!** Pre-clean all surfaces to be treated with this product. Apply this product by wiping to thoroughly wet hard, non-porous surfaces.

#### BACTERICIDAL

This product kills the following bacteria in 2 minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

*Mycobacterium bovis*, BCG (tuberculosis) (TB)  
*Escherichia coli*-with ESBL (E. coli)  
Carbapenem-Resistant *Klebsiella pneumoniae* (CRKP)  
Methicillin-Resistant *Staphylococcus aureus* (MRSA)  
*Salmonella enterica*  
Vancomycin-Resistant *Enterococcus faecium* (VRE)

*Acinetobacter baumannii*  
*Enterobacter aerogenes*  
*Listeria monocytogenes*  
*Pseudomonas aeruginosa*  
*Serratia marcescens*  
*Staphylococcus aureus*

#### VIRUCIDAL\*

This product kills the following viruses\* in 2 minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

Avian Influenza A Reassortant Virus (H3N2)  
Hepatitis B Virus\* (HBV)  
Herpes Simplex II Virus (HSV 2)  
Rhinovirus Type 14 (RV-14)

Rotavirus Strain WA  
Hepatitis C Virus\* (HCV)  
Influenza A Virus (H1N1)  
Influenza A Virus Strain AZ/Hong Kong

SARS-Related Coronavirus 2 SARS-CoV-2  
Human Immunodeficiency Virus\* (HIV-1) (associated with AIDS)

For disinfection of Avian Influenza A, use this product only for surfaces that are not conducive to treatment by immersion or excess liquid (e.g., cages, equipment tables, lab benches, faucet handles, door knobs and similar "high touch" surfaces). Food and drink dispensing and preparation surfaces must be rinsed with potable water prior to use.

#### FUNGICIDAL

This product kills the following fungi in 2 minutes on pre-cleaned hard, non-porous surfaces at room temperature when used as directed:

*Trichophyton interdigitale* (ringworm)

#### Dispensing Directions

To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation.

#### Cleaning Directions

For surface cleaning: To clean surfaces, use a fresh wipe on surfaces. Repeat as necessary until surfaces are visibly clean.

#### Disinfection Directions

For use as a hard, non-porous non-food contact surface disinfectant/virucide/fungicide: Pre-clean hard, non-porous surfaces to be disinfected with a clean wipe. Use a fresh, clean wipe to disinfect hard, non-porous surfaces until thoroughly wet. Allow surfaces to remain wet for 2 minutes at room temperature.

#### \*HIV-1, HBV and HCV Precautions

This product kills HIV-1, HBV, and HCV on pre-cleaned hard, non-porous surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special instructions for cleaning against HIV-1, HBV, and HCV on surfaces/objects soiled with blood/body fluids.**

**Personal Protection:** Wear appropriate barrier protection including gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow hard, non-porous surfaces to remain wet for 2 minutes to kill HIV-1, HBV, and HCV. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

<sup>^</sup> This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean critical and semi-critical medical devices prior to sterilization or high-level disinfection.

### STORAGE AND DISPOSAL

Do not contaminate food or feed by storage or disposal.

**Storage:** Store in a cool, well-ventilated area. **Product Disposal:** Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. **Container Disposal:** Non-refillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard container in trash.

#### PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

**CAUTION:** Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wear appropriate personal protective equipment based on safe work practices. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

#### Physical or Chemical Hazards

Do not use or store near heat or open flame.

#### FIRST AID

**If Inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. **If on Skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **If in Eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact the National Pesticides Information Center (NPIC) at 1-800-858-7378 for non-emergency information concerning this product.



(01) 37613327405256

M60043F-85468B-1223