InTouch® Critical Care Bed
InTouch® Intensivpflegebett
InTouch®-bed voor kritieke zorg
Cama de cuidados intensivos InTouch®

REF FL27 (2151/2152/2153/2154/2155/2156) InTouch CE

Version 3

stryker®

Maintenance Manual

Wartungshandbuch
Onderhoudshandleiding
Manual de manutenção

For parts or technical assistance call:
Ersatzteile oder technische Unterstützung:
Voor onderdelen of technische bijstand belt u:
Para encomendar peças ou solicitar assistência técnica:
USA/États-Unis/EE. UU./USA/USA/VS/EUA/USA/USA/USA: 1-800-327-0770

CE 0086

International Addresses

EUROPE HEADQUARTERS

Stryker SA
Cite-Centre, Grand-Rue 90
P.O. Box 1568
1820 Montreux, Switzerland
Phone: 011-41-21-966-12-01
Fax: 011-41-21-966-12-00

EASTERN EUROPE

Stryker SA - Export Business
Grand-Rue 90
P.O. Box 1567
1820 Montreux, Switzerland
Phone: 011-41-21-966-14-00
Fax: 011-41-21-966-14-01

AUSTRALIA

Stryker Australia
8 Herbert Street
St Leonards NSW 2065
Australia
Phone: 61-2-9467-1000
Fax: 61-2-9467-1010

AUSTRIA

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SOUTH AFRICA

Stryker Osteonics PTY. LTD.
3 Susan Street
Strydom Park
Johannesburg, 2194 - South Africa
Mailing Address:
P.O. Box 48039
2129 Roosevelt Park
Johannesburg, South Africa
Phone: 011-27-11-791-4644
Fax: 011-27-11-791-4696

SPAIN

Stryker Iberia SL
c/Sepulveda n17
28108 Alcobendas
Madrid
Spain
Phone: 011-34-91-728-35-00
Phone: 011-34-91-358-20-44
Fax: 011-34-91-358-07-48

SWEDEN

Stryker AB/Scandinavia
Box 50425
SE-204 14 Malmo
Sweden
Phone: 011-46-40-691-81-00
Fax: 011-46-40-691-81-91

Table of Contents

English

Symbols and Definitions	1-3
Warning / Caution / Note Definition.	1-4
Introduction	1-5
Intended Use	1-5
Brief Product Description and Intended Use Statement	1-5
Product Illustration	1-6
Specifications	1-7
Summary of Safety Precautions	1-9
Optional Zoom® Drive System (Model 2152/2154/2156 Only)	1-11
Cleaning.	1-12
Bed Cleaning.	1-12
Mattress Care	1-12
Preventative Maintenance.	1-13
Checklist	1-13
Quick Reference Replacement Parts List	1-14
Maintenance Menu Guide.	1-16
Accessing Configuration Screen.	1-16
Configuration Screen	1-17
Troubleshooting Guide	1-32
Scale Troubleshooting	1-39
Maintenance Error Messages	1-40
Error Handling	1-40
Error Messages	1-40
CPU/Power Board - QDF75-0440	1-43
Fuse Specification	1-44
Bed Electrical Diagram	1-45
Service Information.	1-62
Bed lift Actuator (Head) Removal and Replacement - (Base)	1-62
Bed lift Actuator (Foot) Removal and Replacement - (Base)	1-63
Fowler Actuator Removal and Replacement - (Litter)	1-64
Gatch Actuator Removal and Replacement - (Litter)	1-65
Foot Actuator Removal and Replacement - (Litter)	1-66
Zoom® Drive Actuator Removal and Replacement (2152 Model Only) - (Base)	1-67
CPU / Power Board Removal and Replacement - (Litter)	1-68
Load Cell (Head End) Removal and Replacement - (Litter)	1-69
Load Cell (Foot End) Removal and Replacement - (Litter)	1-70
Display Removal and Replacement - (Footboard)	1-71
Brake Control Board Removal and Replacement - (Footboard)	1-72
Function Selection/LED Board Removal and Replacement - (Footboard)	1-73
Touch Screen Removal and Replacement - (Footboard)	1-74
Battery Removal and Replacement - (Footboard)	1-75
Brake / Neutral / Drive Potentiometer Removal and Replacement - QDF27-2024	1-76
Battery Removal and Replacement - (Litter)	1-77
Fowler Angle Sensor Removal and Replacement - (Litter)	1-78
Gatch Angle Sensor Removal and Replacement - (Litter)	1-79

Table of Contents

English

Service Information (Continued)	
Foot Angle Sensor Removal and Replacement - (Litter)	1-80
Base Angle Sensor Removal and Replacement - (Base)	1-81
Trend Angle Sensor Removal and Replacement - (Litter)	1-82
Headwall Communication Board Removal and Replacement - (Litter)	1-83
EMC Information.	1-84
InTouch® Critical Care Bed	1-84
Warranty	1-88
Limited Warranty	1-88
To Obtain Parts and Service	1-88
Return Authorization.	1-88
Damaged Merchandise	1-88

Symbols and Definitions

English

	Warning/Caution, consult accompanying documentation		
	Safe Working Load		
~	Alternating Current		
	Fuse Rating for Beds with 230V~ Electric System		
	<p>Type B Equipment: equipment providing a particular degree of protection against electric shock, particularly regarding allowable leakage current and reliability of the protective earth connection.</p> <p>Class 1 Equipment: equipment in which protection against electric shock does not rely on BASIC INSULATION only, but which includes an additional safety precaution in that means are provided for the connection of the EQUIPMENT to the protective earth conductor in the fixed wiring of the installation in such a way that ACCESSIBLE METAL PARTS cannot become live in the event of a failure of the BASIC INSULATION.</p> <p>Mode of Operation: Continuous</p>		
IPX4	Protection from liquid splash		
	Dangerous Voltage		
	Protective Earth Terminal		
	Potential Equalization		
	Medical Equipment approved by the Canadian Standards Association with Respect to Electric Shock, Fire, Mechanical and Other Specified Hazards.		
	In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.		
	Model Number		Measuring Device
	Manufacturer		Accuracy Class
	Batch Code		Instructions for Use

[Return To Table of Contents](#)

Symbols and Definitions

WARNING / CAUTION / NOTE DEFINITION

The words WARNING, CAUTION, and NOTE carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation, which if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation, which if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note

This provides special information to make maintenance easier or important instructions clearer.

INTENDED USE

This manual is designed to assist you with the maintenance of Stryker *InTouch*® Critical Care Bed Model FL27 (2151/2152, 2153/2154, 2155/2156). Carefully read this manual thoroughly before using the equipment or beginning maintenance on it. To ensure safe operation of this equipment, it is recommended that methods and procedures be established for educating and training staff on the safe maintenance of this bed.

This Maintenance Manual is an integral part of the bed and should be included if the bed is sold or transferred.

BRIEF PRODUCT DESCRIPTION AND INTENDED USE STATEMENT

INTENDED MEDICAL INDICATIONS: In-touch is an AC-powered adjustable hospital bed designed to position patients for procedures, therapy, and recovery in healthcare environment, transport patients between bays and procedural rooms, provide patient security, measure and display patient weight, allow patients to alert caregiver when the patient requires emergency assistance or any assistance to improve comfort levels, and monitor patient position to alert caregiver of a deliberate exit or potential fall. The bed has thirty-nine prerecorded clinical phrases in 24 languages and offers music therapy.

INTENDED PATIENT POPULATION: The intended patient population is acute-care human patients. The safe working load (i.e. the sum of the patient, the mattress, and accessory weight) for InTouch is 550lbs (250kg). This bed is not intended for use with patients 35 inches or less and/or patients that weigh less than 50 pounds.


INTENDED PART OF THE BODY: InTouch is intended to support a patient. It is intended to be used with a sleep surface. The Frame can come in contact with human skin but a patient should never be on the frame without a support surface.

INTENDED USER PROFILE: InTouch is intended for use in Acute Care settings. These settings may include critical care, step down, progressive care, med/surg, sub-acute care, and post anesthesia care unit (PACU) or other locations as prescribed. Intended operators are health care professionals (nurses, nurse aids, doctors) which can use all bed operations (e.g., bed motion functions, nurse call, siderail communications, bed exit, therapy options), patient and bystander, which can use bed motion functions, nurse call & siderail communications, and trained professionals for installation, service, and calibration.

INTENDED CONDITION OF USE AND EXPECTED LIFE: Stryker Medical Bed products are designed for a 10 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its Bed products will be free from structural defects for the expected 10 year life of the Bed product as long as the original purchaser owns the products.

The product is intended to be used in a healthcare environment, including hospitals, surgery centers, long term acute care centers, and rehabilitation centers. The intended environmental conditions are 10 to 40 degrees Celsius and 30 to 75% RH. The product is compatible with 35" x 84" surfaces, the facility nurse call system, standard Med-Surg equipment, and the facility infrastructure. InTouch is intended to be used with a 6"-8.5" sleep surface; asleep surface or overlay greater than 6" that offers therapeutic value may be used with, patient supervision.

PRODUCT ILLUSTRATION

Introduction

English

SPECIFICATIONS

	Safe Working Load		
	Note: Safe Working Load indicates the sum of the patient, mattress and accessory weight.		
		249 kg	550 lbs
Bed Weight		340.2 kg	
Overall Bed Length		228.6 cm	90"
Overall Bed Width	Siderails Up	106.7 cm	42"
	Siderails Down	102.9 cm	40"
Base	Under Bed Clearance	12.7 cm	5"
Litter	Patient Surface		
	• Width	88.9 cm	35"
	• Length	213.4 cm	84"
	• Length (with Optional Bed Extender)	228.6 cm	90"
	Seat Section		
	• Depth	47 cm	18.5"
	Foot Section		
	• Length	73.7 cm	29"
	• Angle	0° to 50°	0° to 50°
	Fowler Section		
	• Length	91.4 cm	36"
	• Width	86.4 cm - 88.9 cm	34" to 35"
	• Angle	0° to 70°	0° to 70°
	Gatch Section		
	• Length	45.7 cm	18"
	• Width	86.4 cm - 88.9 cm	34" to 35"
	• Angle	0° to 15°	0° to 15°
	Cardiac Chair		
	• Standard Cardiac Chair Position	Head: 65°, Seat: 17°, Foot: 30°, Trend: 3°	
	• Enhanced Cardiac Chair Position	Head: 70°, Seat: 19°, Foot: 47°, Trend: 3°	
Fowler	Length	92.7 cm	36.5"
Bed Lift System	Height (high) to top of litter	83.8 cm	33"
	Height (low) to top of litter	40.6 cm	16"
	Trendelenburg/Reverse Trendelenburg	12°	
	Bed Lift Time	35 seconds maximum from lowest to highest position	
Scale System	Maximum Load	249 kg	
	Minimum Load	20 kg	
	e (Error of Accuracy)	2 kg	
	T (Tare Value)	60 kg	
CPR System	Speed to level bed from any position		
	• Fowler	15 seconds	
	• Foot/Seat	60 seconds	
Drive System	Speed (Optional Zoom® Drive)		
	• Forward	4.8 km/h	2.98 mph
	• Backwards	2.88 km/h	1.79 mph

[Return To Table of Contents](#)

Introduction

English

SPECIFICATIONS (CONTINUED)

Mattress	Recommended Mattress Size With Bed Extender DM64197 (XPRT™) (35" x 6" x 5.5")	88,9 cm x 213,4 cm x 15,2 cm	35" x 84" x 6"
		88,9 cm x 228,6 cm x 13,9 cm	35" x 90" x 5.5"
	Recommended Air Mattress Size With Bed Extender DM64197 (XPRT™) (35" x 6" x 5.5")	88,9 cm x 213,4 cm x 15,2 cm - 21,6 cm	35" x 84" x 6" - 8.5"
		88,9 cm x 228,6 cm x 13,9 cm	35" x 90" x 5.5"
Maximum Current Consumption	Without Optional Auxiliary Outlet (120VAC Only)	4.8 Amps	
	With Optional Auxiliary Outlet(s) (120VAC Only)	9.8 Amps	
Environmental Conditions	Storage • Ambient Temperature • Relative Humidity • Atmospheric Pressure	-40 °C to 70 °C 10 to 100% 500 to 1060 hPa	-40 °F to 158 °F
	Operating • Ambient Temperature • Relative Humidity • Atmospheric Pressure • Scale Accuracy	10 °C to 40 °C 30 to 75% 700 to 1060 hPa 18 °C to 25 °C	50 °F to 104 °F 64 °F to 77 °F
Electrical Requirements	All electrical requirements meet CSA C22.2 No. 601.1, UL 60601-1 and IEC 60601-1.60601-2-38 specifications.	230V~, 50-60 Hz, 4.0A (with Auxiliary Outlet) - Two 250V, 10A Fuses	
Battery	12V, 17.6 Ahr, Sealed Lead-Acid Battery (Part Number QDF9188) 3.0V 220mAh Lithium Battery, Size 20 mm - Varta Int. CR2032 Coin Cell PC Battery Holder, Size 20 mm - MDP Int. BA2032		

Stryker reserves the right to change specifications without notice.

Specifications listed are approximate and may vary slightly from unit to unit or by power supply fluctuations.

Summary of Safety Precautions

Before operating the bed, it is important to read and understand all information in this manual. Carefully read and strictly follow the safety guidelines listed below.

Note

To ensure safe operations of the bed, methods and procedures must be established for educating and training hospital staff on the intrinsic risks associated with the usage of electric beds.

WARNING

- To prevent permanent damage to this unit, the unit must reach room temperature prior to conducting any setup and/or unit operations.
- The Weigh system is intended to assist in the monitoring of the patient's weight variation. Under no circumstances should its reading be used as sole reference for medical treatment.
- For an accurate measurement please ensure that the patient is completely supported by the surface of the bed. Any weight under 20 kg may not be accurate.
- Preventative maintenance should be performed at least once a year to ensure all bed features are functioning properly.
- This bed is not intended for pediatric use.
- This bed is equipped with a hospital grade plug for protection against shock hazard. It must be plugged directly into a properly grounded receptacle. Grounding reliability can be achieved only when a hospital grade receptacle is used.
- Shock Hazard - Improper handling of the power cord may result in damage to the power cord and potential shock hazards. If damage has occurred to the power cord, immediately remove the bed from service and contact the appropriate maintenance personnel. Failure to do so could result in serious injury or death.
- Serious injury can result if caution is not used when operating the bed. Operate the bed only when all people and equipment are clear of the electrical and mechanical systems.
- Always apply the brakes when a patient is on the bed or entering/exiting the bed. Serious injury could result if the bed moves while a patient is getting on or off the bed. Once the brake pedal is engaged, push on the bed to ensure the brakes are securely applied. Do the same test when using the electrical brake.
- Do not apply the brakes when the bed is in motion.
- Keep siderails in the fully raised position and the sleep surface horizontal in its lowest position when the patient is unattended, unless his/her medical condition dictates otherwise. When raising the siderails, be sure that you hear the "click" that signals the locked condition. Pull firmly on the siderail to ensure it is locked into position.
- When the sleep surface sections are articulated, ensure that all the patient's limbs are within the raised siderails to avoid patient injury.
- When a patient's condition requires greater safety measures for his/her security, use the lockout controls in the footboard control panel to inhibit the siderail functions or remove any optional pendant control and install protective pads on the siderails.
- Siderails, with or without their padded covers, are not intended to serve as restraint devices to keep patient from exiting the bed. Siderails are designed to keep a patient from inadvertently rolling off the bed. It is the responsibility of the attending medical personnel to determine the degree of restraint necessary to ensure a patient will remain safely in bed. Failure to use the siderails properly could result in serious patient injury.
- To reduce risk of injury, ensure the sleep surface is horizontal and in the lowest position with the siderails fully raised and locked when moving the bed with a patient in it.
- To avoid injury to the patient and/or user, do not attempt to move the bed laterally with the steer mode engaged. The steer wheel cannot swivel.
- The CPR emergency pedal is for emergency use only. When activating the CPR pedal, all people and equipment must be removed from the area below and around the head, thigh and foot sections of the bed or serious personal injury and/or equipment damage could occur.
- The manual backup brake is for emergency use only. It should not be used for any other situation than an emergency or it might get overused when the time comes to use it for an emergency.
- Possible fire hazard exists when this bed is used with oxygen administering equipment other than nasal, mask type or half bed-length tent type. Unplug the bed power cord from the wall when oxygen administering equipment is used. When using a half bed-length tent type, ensure that the siderails are outside the oxygen tent and that the tent does not extend below the mattress support level.

Summary of Safety Precautions

WARNING (CONTINUED)

- The Bed Exit system is intended only to aid in the detection of a patient exiting the bed. It is not intended to replace patient monitoring protocol. The Bed Exit system signals when a patient is about to exit the bed.
- The Bed Exit system is not designed to be used with patients weighing less than 50 lbs (23 kg).
- Do not steam clean, hose off or ultrasonically clean the bed. Do not immerse any part of the bed. The internal electrical parts may be damaged by exposure to water. Hand wash regularly all surfaces of the bed with warm water and a mild detergent. Wipe cleaned surfaces dry to avoid build-up of cleaning substance. Inspect the mattress after each use. Discontinue use if any cracks or rips, which may allow fluid to enter the mattress, are found in the mattress cover. Failure to properly clean the mattress, or dispose of it if defective, may increase the risk of exposure to pathogenic substances and may bring about diseases to the patient and user.
- Always unplug the bed power cord from the wall socket when servicing or cleaning the bed. When working under the bed with the bed in the high position, always apply the brakes and place blocks under the Bed Lift levers to prevent injury in case the Bed Down switch is accidentally pressed.

CAUTION

- Preventative maintenance should be performed at least once a year to ensure all bed features are functioning properly. Ensure that any bed malfunction is promptly reported to your service personnel for immediate attention.
- When large fluid spills occur in the area of the circuit boards, cables and motors, immediately unplug the bed power cord from the wall socket, remove the patient from the bed and clean up the fluid. Have maintenance completely check the bed. Fluids can have an adverse effect on operational capabilities of any electrical product. DO NOT put the bed back into service until it is completely dried and has been thoroughly tested for safe operation. Ensure, among other things, that the plastic components being used as covers for the siderail mechanism arms and the foot end casing are removed and that the parts they cover are thoroughly dried.
- To avoid damage to the siderail mechanisms, do not move the bed using the raised siderails. Move the bed using the handles integrated to the boards.
- Because individual beds may have different options, footboards should not be moved from one bed to another. Mixing footboards could result in unpredictable bed operation.
- Leave the frame in the lowest position when the patient is unattended to help reduce the number and severity of a potential fall.
- The mattress thickness should not exceed 6 inches (15.24 cm).
- When using a mattress thicker than 6 inches or mattress overlay, extra caution and/or patient supervision may be required to reduce the likelihood of occurrence of a patient fall.
- When servicing use only identical replacement parts provided by Stryker.
- For the batteries to maintain sufficient charge, you must plug the bed into a wall outlet (regulated AC power source) when not in use to maximize bed performance while on battery power.
- If batteries are observed to have corrosion at the terminals, display cracking, have expanded or are bulging at the sides, or no longer can maintain a full charge, they should be replaced immediately.
- When replacing the batteries, use only Stryker authorized batteries. Use of non-Stryker batteries may lead to unpredictable system performance.
- Upon a Battery Low alarm (Battery Low LED on Footboard and/or audible beep), stop using the Zoom® drive and recharge the batteries immediately. Ignoring the Battery Low alarms may cause your batteries to degrade quicker than normal and may decrease battery life.

Note

- Throughout this Maintenance Manual, the words “right” and “left” refer to the right and left sides of a patient lying on his/her back on the bed.
- The addition of accessories affects the motion of the bed.
- The iBed option is only used as an informational feature and is not supposed to replace normal patient monitoring protocol.

Summary of Safety Precautions

English

OPTIONAL ZOOM® DRIVE SYSTEM (MODEL 2152/2154/2156 ONLY)

In addition to the previous warnings and cautions, all of the following warnings and cautions apply to units equipped with the Zoom®.

WARNING

- The 2152/2154/2156 *InTouch*® Bed is intended for use by trained hospital personnel only. Failure to properly train personnel could result in injury.
 - USE CAUTION while maneuvering the unit with the drive wheel activated. Always ensure there are no obstacles near the unit while the drive wheel is activated. Injury to the patient, user or bystanders or damage to the frame or surrounding equipment could occur if the unit collides with an obstacle.
 - Use caution when transporting the unit down halls, through doors, in and out of elevators, etc. Damage to the siderails or other parts of the unit could occur if the unit comes in contact with walls or door frames.
 - Put the drive wheel in the neutral position and release the brake before pushing the unit manually. For 2152/2154/2156 models; push the Brake Off button to disengage drive wheel (Zoom®) before pushing the unit manually. Do not attempt to push the unit manually with the drive wheel engaged. The unit will be difficult to push and injury could result.
 - If unanticipated motion occurs, unplug the power cord from the wall socket, push the battery power on/off switch to the “OFF” (0) position (the LED will not be illuminated), actuate the drive wheel pedal to the neutral position and call maintenance.
 - In the event the batteries become discharged while using the Zoom® drive, press the “N” (Brake Off) button to place the drive wheel in Neutral and proceed to push the bed manually. Recharge the batteries prior to using Zoom® again. Continuing to use the Zoom® drive could lead to battery damage and the drive wheel getting stuck in the down position.
 - The power save mode is activated after one hour on battery power with no motion release switch activation. Functions including Bed Exit, Scale and Motion will cease to operate when the unit enters the power save mode. Injury to the patient could occur if proper patient monitoring protocol is not observed.
 - Always unplug the power cord and push the battery power on/off switch to the “OFF” (0) position before service or cleaning. When working under the frame, always support the litter frame to prevent injury in case the Bed Down switch is accidentally activated.
 - Battery posts, terminals and related accessories contain lead and lead compounds, chemicals known to the State of California to cause cancer and birth defects or other reproductive harm. Wash hands after handling.
 - Do not modify the *InTouch*® Model 2152/2154/2156 bed. Modifying the unit can cause unpredictable operation resulting in injury to the patient or user. Modifying the unit will also void this warranty.
-

Cleaning

BED CLEANING

Hand wash all surfaces of the bed with warm water and mild detergent. Wipe the bed clean and DRY THOROUGHLY to avoid build up of cleaning solution. Do not steam clean or hose off the bed. Do not immerse any part of the bed. Some of the internal parts of the bed are electric and may be damaged by exposure to water.

Suggested cleaners for bed surfaces:

- Quaternary Cleaners (active ingredient - ammonium chloride).
- Phenolic Cleaners (active ingredient - o-phenylphenol).
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water).

Avoid over saturation and ensure the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

CAUTION

Some cleaning products are corrosive in nature and may cause damage to the product if used improperly. If the products described above are used to clean Stryker patient care equipment, measures must be taken to ensure the beds are wiped with a damp cloth soaked in clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the beds will leave a corrosive residue on the surface of the bed, possibly causing premature corrosion of critical components. Failure to follow the above directions when using these types of cleaners may void this product's warranty.

MATTRESS CARE

Inspect the mattress after each use. Discontinue use if any cracks or rips are found in the mattress cover as these openings may allow fluid to enter the mattress.

Preventative Maintenance

English

Beds require an effective maintenance program, we recommend checking these items annually. Use this sheet for your records. Keep on file.

CHECKLIST

- _____ All fasteners secure (reference all assembly prints).
- _____ Engage brake pedal and push on the bed to ensure all casters lock securely.
- _____ "Brake Not Set" LED on the footboard and head end siderails when brakes are not engaged.
- _____ Locking steer caster engages and disengages properly (2151/2153/2155 models only).
- _____ Siderails move, latch and stow properly.
- _____ CPR release working properly.
- _____ I.V. pole working properly (if on bed).
- _____ No cracks or splits in head or footboards.
- _____ No rips or cracks in mattress cover.
- _____ All functions on head end siderails working properly (including LEDs).
- _____ All functions on footboard working properly (including LEDs).
- _____ Scale and bed exit system working properly.
- _____ Night light working properly.
- _____ Power cord not frayed.
- _____ No cables worn or pinched.
- _____ All electrical connections tight.
- _____ All grounds secure to the frame.
- _____ Ground impedance not more than 100 milliohms.
- _____ Current leakage not more than 300 microamps.
- _____ Engage drive wheel and ensure it is operating properly (Zoom® option - 2152/2154/2156 model only).
- _____ Motion release switches working properly (Zoom® option - 2152/2154/2156 model only).
- _____ Confirm head end Zoom® handle functionality (2152/2154/2156 model only).
- _____ Replace batteries every two years.
- _____ Replace batteries immediately if batteries are observed to have corrosion at the terminals, display cracking, have expanded or are bulging at the sides, or no longer can maintain a full charge. Use only QDF9188 for battery replacement.
- _____ Ensure ground chains are clean, intact, and have at least two links touching the floor.
- _____ Check fowler angle for accuracy 0° - 70°.
- _____ Check gatch angle for accuracy 0° - 15°.
- _____ Check foot section angle for accuracy 0° - 50°.
- _____ Siderail switches working properly (iBed Awareness option).
- _____ iBed Awareness Light Bar LEDs on footboard and siderails working properly (iBed Awareness option).
- _____ Inspect footboard control labeling for signs of degradation.
- _____ Ensure calibration of the touch screen is accurate.
- _____ Ensure calibration of the bed is accurate.

Bed Serial Number:

_____	_____
_____	_____
_____	_____
_____	_____

Completed by: _____

Date: _____

[Return To Table of Contents](#)

Quick Reference Replacement Parts List

English

The parts and accessories listed on this page are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call Stryker Customer Service USA: 1-800-327-0770 for availability and pricing.

Part Name	Part Number
Electrical Components	
Actuator, Brake	QDF27-1227
Actuator, Foot	QDF27-1216
Actuator, Fowler	QDF27-1214
Actuator, Gatch	QDF27-1215
Actuator, Lift, Foot End	QDF27-1251
Actuator, Lift, Head End	QDF27-1252
Angle Sensor, Foot, Fowler, Base, Gatch, Lift	27-2477
Batteries (Replace both at same time)	QDF9188
Board, CPU/Power	QDF75-0440
Board, Drive (2152 Zoom® Only)	QDF27-1430
Board, Headwall	QDF75-0600
Board, Brake Control	QDF27-1097
Board, Footboard, Function Section/LEDs	QDF75-0010
Board, Siderail, Outside (Bed Motion)	QDF27-1099
Board, Touch Screen, Foot	QDF2125
Fuse, 8Amp Ceramic	QDF2120
Fuse, 10Amp, Main Power	QDF8078
Fuse, 25Amp Cartridge	QDF2119
Load Cell	QDF27-1372
Motor, Drive (2152 Zoom® Only)	QDF27-1445
Motor, Drive Actuator (2152 Zoom® Only)	27-2593
Position Sensor, Brake/Brake Off/Drive	QDF27-2024
Power Cord	QDF8066
Speaker, Right Siderail	QDF27-2216
Speaker, Left Siderail	QDF27-2175
Touch Screen, Footboard	QDF27-2193
Transformer	QDF27-2658
Other Components	
Caster (2152 Model - all four casters) (2151 Model - both head end casters)	RD27-1970
Caster, Steer (2151 Model Only - both foot end casters)	RD27-1971
Headboard Assembly	27-2583K
InTouch 3.0 Footboard Assembly without iBed, without iAudio	27-2760K
InTouch 3.0 Footboard Assembly with iBed, without iAudio	27-2761K
InTouch 3.0 Footboard Assembly with iBed, with iAudio	27-2762K

[Return To Table of Contents](#)

Maintenance Menu Guide

The Maintenance Menu is accessed through the Touch Screen and contains additional features of the product. This menu provides an interface to the user and/or service personnel in order to provide the ability to control and access maintenance features.

ACCESSING CONFIGURATION SCREEN

WARNING

Please ensure patient is not in the bed prior to starting bed calibration. In calibration mode, the software does not control the interferences between the mechanical parts of the bed. Mechanical damage could occur without supervision. Only qualified personnel should perform the calibration.

Note

Verify that the bed is on a level surface which does not have any slopes or inclines prior to entering into the calibration mode.

1. Unseat and reseat the footboard and wait until the main control screen is displayed.
2. Push and hold the Main Menu button located in the upper right corner of the footboard control panel (see Figure 1 below). Continue pressing on the Main Menu button while executing steps 3-5.
3. Push and hold the **Brake** button for 5 seconds and then release.
4. Push the **HOB 30°+** button once and release.
5. Push the **Drive** button once and release.
6. Release the Main Menu button and you will be taken into the Configuration Screen shown in Figure 2.

Figure 1: Footboard Control Panel


Figure 2: Configuration Screen

CONFIGURATION SCREEN

The following items A, B, and C are configuration buttons available and displayed on the Configuration screen (Figure 3).

- A. Bed Calibration
- B. Full Diagnostic (see page 1-28)
- C. Touch Screen Calibration (see page 1-29)

The following items D and E are configuration options available but are **not displayed** on the Configuration screen. To access these configuration options, follow the procedures identified on page listed below.

- D. Bed Options Configuration (see page 1-30)
- E. Serial Number Configuration (see page 1-31)


Figure 3: Configuration Screen Menu Items

A. Bed Calibration

CAUTION

Prior to calibration, connect the bed to an A/C outlet.

Note: During calibration, if the backlight shuts off touch the screen to continue.

WARNING

Refer to Figure 4 for WARNING details.

1. To start the bed calibration, press the Bed Calibration button on the Configuration Screen (refer to Figure 3).
2. After pressing the Bed Calibration button, Figure 4 will be displayed. Press "Next" to continue with calibration.


Figure 4: Bed Calibration - WARNING

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

- After pressing "Next" to continue calibration, Figure 5 will appear.

Figure 5: Bed Calibration - Step #1 of 6

- Place the Litter surface to a flat position by pressing simultaneously the **Foot Up**, **Fowler Down** and **Gatch Down** buttons as shown in Figure 6.


Figure 6: HE Siderail Control Panel - (Right Side Shown)

- Lower the bed to the lowest position. Press the **Bed Height Down (Head Lift Down)** and **Reverse Trendelenburg (Foot Lift Down)** buttons simultaneously as shown in Figure 7.


Figure 7: HE Siderail Control Panel - (Right Side Shown)

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

6. Press the  button when done. The “Do Not Touch Bed” screen will appear as shown in Figure 8.


Figure 8: Do Not Touch Bed Screen

7. When step 1 of the calibration procedure has completed, step 2 of the calibration procedure will begin and Figure 9 will be displayed as shown below.

Figure 9: Bed Calibration - Step 2

Maintenance Menu Guide

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

8. Raise the bed height to 20 inches measuring from the top of the seat section to the floor. Press the **Bed Height Up (Head Lift Up)** and **Trendelenburg (Foot Lift Up)** buttons as shown in Figure 10.


Figure 10: Head Siderail Control Panel
(Right Side Shown)

9. Standing on either side of the bed at the fowler section, raise the foot siderail and then position the digital level/inclinometer on the bottom of the litter below the mattress retainer (refer to Figure 11). Using the digital level/inclinometer, verify the bed is level (0.0 ± 0.2).

Note: Cycle power on the digital level/inclinometer prior to placing it on the bottom of the litter frame and do not zero/calibrate the digital level/inclinometer.

Placement
on Litter

Figure 11: Inclinometer placement

10. Press the **Next** button when done. The "Do Not Touch Bed" screen will appear as shown in Figure 12.


Figure 12: Do Not Touch Bed Screen

A. Bed Calibration (Continued)

- Bed Calibration

stryker

Step 3 of 6

Put Bed:

- At +12 degrees Trendelenburg
+/- 0.1 degree

- Litter surface flat

Back

Next

Controls

Scale

Bed Exit

Surfaces

Options

Lockouts

Bed

**Trendelenburg
(Foot Lift Up)**

The diagram shows a patient in the Trendelenburg position (head down, feet up). The patient is lying on their back with their feet raised and supported by a wedge. The patient's head is resting on a pillow. The diagram is labeled 'Trendelenburg (Foot Lift Up)'.

1-21

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

13. Press the  button when done.

14. The “Do Not Touch Bed” screen will appear as shown in Figure 15.


Figure 15: Do Not Touch Bed Screen

15. When step 3 of the calibration procedure has completed, step 4 of the calibration procedure will begin and Figure 16 will be displayed as shown below.

Figure 16: Bed Calibration - Step 4

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

16. Level the bed back out to zero degrees by pushing the **Reverse Trendelenburg** button until the litter is level while referencing the inclinometer. Note: Confirm inclinometer reads zero degrees.
17. Place the bed at -12 degrees Reverse Trendelenburg by pressing the **Bed Height Up (Head Lift Up)** button as shown in Figure 17. Verify -12 degrees +/- 0.1 with the inclinometer you previously placed on the bottom of the litter frame in step 9.


Figure 17: Head Siderail Control Panel - (Right Side Shown)

18. Press the  button when done.

19. The "Do Not Touch Bed" screen will appear as shown in Figure 18.


Figure 18: Do Not Touch Bed Screen

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

20. When step 4 of the calibration procedure has completed, step 5 of the calibration procedure will begin and Figure 19 will appear on the screen as shown below.

Figure 19: Bed Calibration - Step 5

21. Place the bed at the highest height by pressing simultaneously the **Trendelenburg (Foot Lift Up)** button and **Bed Height Up (Head Lift Up)** button as shown in Figure 20.


Figure 20: Head Siderail Control Panel - (Right Side Shown)

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

22. Press the **Fowler Up** button until the Fowler reaches the highest height, next press the **Gatch Up** button until the Gatch reaches its highest height, lastly press the **Foot Up** button until the foot section reaches its highest height (refer to Figure 21).
The foot section should be at a flat position.


Figure 21: Head Siderail Control Panel - (Right Side Shown)

23. Press the  button when done.

24. The “Do Not Touch Bed” screen will appear as shown in Figure 22.


Figure 22: Do Not Touch Bed Screen

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

25. When step 5 of the calibration procedure has completed, step 6 of the calibration procedure will begin and Figure 23 will be displayed as shown below.

Figure 23: Bed Calibration - Step 6

26. Place the bed at the highest height by pressing simultaneously the **Trendelenburg (Foot Lift Up)** button and **Bed Height Up (Head Lift Up)** button (see Figure 24).
27. Place the Fowler and Gatch section of the bed at the highest height by first pressing the **Fowler Up** button then the **Gatch Up** button (See Figure 24).
28. Place the foot at the lowest position by pushing the **Foot Down** button until limit is met (see Figure 24) . **CAUTION: The Gatch must be raised to the highest height prior to running the Foot Down or damage could occur.**


Figure 24: Head Siderail Control Panel - (Right Side Shown)

CONFIGURATION SCREEN (CONTINUED)

A. Bed Calibration (Continued)

29. When the calibration procedure is completed, Figure 25 will be displayed as shown below.
30. Press the Close button to exit the Calibration Procedure Menu.
31. Level the litter flat using the CPR pedal. Verify all display readings are zero degrees while the bed is at 33" (84 cm) (highest height). Lower the bed to 16" (41 cm) (lowest height) and verify the display readings maintain a constant reading of zero degrees with no fluctuations. When the bed reaches its lowest height at 16", verify the display readings are still at zero degrees.

Note: If readings are not all zero, you will need to recalibrate the bed by repeating steps 1-30.

Figure 25: Calibration Complete

Maintenance Menu Guide

CONFIGURATION SCREEN (CONTINUED)

B. Full Diagnostic

- From the Configuration Screen, press the Full Diagnostic button. Figure 26 will appear.
- The following menu items may be selected by pressing their button.
 - BOARDS** (Figure 27a)
Provides information on the switch boards and the touch screen's software version.
 - ERROR CODES** (Figure 27b)
Provides information on errors which the CPU board has identified.
 - INPUT STATES** (Figure 27c)
Provides information on the status of all switches and jumpers on the bed.
 - MOTORS** (Figure 27d)
Provides information on what a motor is doing when a function button is pushed.
Note: This requires assistance to press the buttons on the head siderails or at the head end control.
 - BUTTONS PRESSED** (Figure 27e)
Provides information on when the CPU has detected a button being pressed.
Note: This requires assistance to press the buttons on the head siderails or at the head end control.
 - SIGNAL VALUES** (Figure 27f)
Provides information on CPU voltages, load cell values, and angle sensor values.


Figure 26: Full Diagnostic Screen


Figure 27a: Boards Connected

Figure 27b: Error Codes

Figure 27c: Input States


Figure 27d: Motors


Figure 27e: Buttons Pressed


Figure 27f: Signal Values

CONFIGURATION SCREEN (CONTINUED)

C. Touch Screen Calibration

Ensure patient is not present on the bed before performing the Touch Screen Calibration.

1. From the Configuration Screen, press the Touch Screen Calibration button, Figure 28 will appear.
NOTE: If the touch screen will not respond, you will need to push the HOB 30°+, Brake and Drive buttons on the footboard at the same time (this will enable you to get directly into the touch screen calibration).
CAUTION: When pushing the buttons above the fowler, brakes, or drive actuators may activate.
2. Carefully press and briefly hold a stylus or your finger on the center of the target (cross-hair) shown on the screen. Repeat as the target moves around the screen. There will be four different locations to press: upper center, lower center, center left and center right.
3. When the last target has been touched, the Figure 29 will be displayed to inform you the 'New calibration settings have been measured'.
4. Tap the screen to register saved data.


Figure 28: Touch Screen Calibration: Start


Figure 29: Touch Screen Calibration: Completed

5. If you do not want to save the new data, wait for 30 seconds to cancel saved data and keep the current setting.

NOTE: The TOUCH SCREEN CALIBRATION SCREEN is a validation of the touch screen's calibration. When the user presses anywhere on the screen, the coordinates x and y are displayed on the lower part of the screen.

CONFIGURATION SCREEN (CONTINUED)

D. Bed Options Configuration

Access the configuration screen as shown on [page 1-16](#). Once in the configuration screen you will need to follow the steps below to access the options screen or the serial number screen.

To access the Bed Options screen:

1. Press and hold the **HOB 30** button and the **BRAKE** button at the same time for 5 seconds then release both buttons. The message screen shown in Figure 30a will be displayed. Press the **Close** button.
2. Once the **Close** button is pressed in step 1, the Bed Options screen will be displayed as shown in Figure 30b.


Figure 30a: European Community Message


Figure 30b: Bed Options Screen

3. On the Bed Options screen, select all of the options that apply to the bed configuration then press the **Continue** button. The 'European Community' option must be selected if you have an InTouch CE bed.

Note: In Figure 30b, the options selected are Zoom, iBed and iAudio, iDoc and European Community. Options selected will turn green.

4. After pressing **Continue** in step 3, a confirmation screen as shown in Figure 31 will be displayed verifying the options you selected.
 - a. If the options on the confirmation screen do not match the options on the bed, press the **Back** button.
 - b. If the options on the confirmation screen match the options on the bed, press the **OK** button.
5. After pressing the **OK** button, the options will be saved. A power cycle screen (Figure 34) will be displayed stating to cycle power on the bed.
6. Cycle the power to the bed by first turning the battery disconnect switch to **OFF** (O) then unplug the power cord from the wall outlet. Plug the power cord back into the wall outlet and then turn the battery disconnect switch back **ON** (I).
7. Test bed functionality prior to returning the bed into service.


Figure 31: Bed Options Confirmation Screen

CONFIGURATION SCREEN (CONTINUED)

E. Serial Number Configuration

To access the Serial Number screen:

1. Press and hold the **HOB 30** button and the **Vascular Position** button at the same time for 25 seconds then release both buttons. The Bed Configuration - Serial Number screen (Figure 32) will be displayed.


Figure 32: Bed Configuration Serial Number Screen

2. Enter or confirm the serial number of the bed (Figure 32).
 - a. **To enter the serial number:**
 - i. Press the pencil button located to the right of the serial number field. The Edit screen will be displayed.
 - ii. Enter the serial number of the bed in the serial number field then press the **OK** button. You will be returned to the serial number main screen as shown in Figure 32
 - iii. Press the **Continue** button and proceed to step 3.
 - b. **To confirm the serial number:**
 - i. Review the serial number displayed in the serial number field.
 - ii. If serial number is correct, press the **Back** button and proceed to step 6
 - iii. If the serial number is incorrect, press the pencil icon button to open the Edit screen then enter the correct serial number.
 - iv. Press the **OK** button and you will be returned to the serial number main screen.
 - v. Press the **Continue** button and proceed to step 3.
3. After pressing the **Continue** button, the serial number will be saved and the serial number Confirmation screen (Figure 33) will be displayed. Press the **OK** button.
4. After pressing the **OK** button a power cycle screen (Figure 34) will be displayed stating to cycle power on the bed.
5. To cycle power on the bed, turn the battery disconnect switch to **OFF** (O) then unplug the power cord from the wall outlet. Plug the power cord back into the wall outlet then turn the battery disconnect switch back **ON** (I).
6. Test bed functionality prior to returning the bed into service.


Figure 33: Serial Number Confirmation Screen


Figure 34: Cycle Power Screen

Troubleshooting Guide

Note

See “CPU/Power Boards” section for an outline of bed PCB’s and voltage test points.

Problem / Failure		Recommended Action
No Power to Bed. (On wall voltage 120VAC)		<div>1. Verify the bed is plugged into a functional wall outlet. A. Check your country voltage option at wall outlet.<div>I. If your country voltage option is present, go to step 2. II. If your country voltage option is not present, contact hospital maintenance staff and try another outlet.</div></div> <div>2. Verify the bed main power fuses are good, located in drawer where power cord plugs into the bed. A. Check for continuity of each 10A fuse.<div>I. If each fuse (see CV Table) is good, go to step 3. II. If either fuse (see CV Table) does not have continuity, replace the fuse.</div></div> <div>3. Verify there is power at the transformer connection (J11) on the CPU/Power board at the foot end. A. Check for 24VAC at J11 between the blue and red wires.<div>I. If 24VAC is present, go to step 4. II. If 24VDC is not present, check the 25 Amp fuse in the fuse holder on the red wire from the transformer. If bad, replace the fuse. III. If 24VDC is not present, check the power cable quick connection going to the transformer for (see CV Table) VAC. If no voltage, follow the cable and repair or replace the damaged component. IV. If (see CV Table) VAC is present, replace the transformer assembly.</div> B. Check for 30VAC at J11 between the yellow and orange wires for the transformer.<div>I. If 30VAC is present, go to step 5. II. If 30VAC is not present, check the power cable quick connection going to the transformer for (see CV Table) VAC. If no voltage, follow the cable and repair or replace the damaged component. III. If (see CV Table) VAC is present, replace the transformer assembly.</div> C. Check fuse F1 on the CPU/Power board.<div>I. If fuse is good, replace CPU/Power board. II. If fuse does not have continuity, replace the fuse (Littelfuse 215008.P).</div></div>
COUNTRY VOLTAGE (CV) TABLE		
VOLTAGE	FUSE	
100V	10A	
110V	10A	
120V	10A	
200V	10A	
220V	10A	
230V	10A	
240V	10A	

Troubleshooting Guide

English

Problem / Failure		Recommended Action
No Bed Up Motion.	FOOT	<ol style="list-style-type: none"> Put the bed into the Bed Calibration menu. <ol style="list-style-type: none"> Using one of the head siderails, push the trend button. <ol style="list-style-type: none"> If the foot lift motor runs up, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> If recalibration does not work, replace the litter angle sensor and recalibrate. If the foot lift motor does not run up, check for 24VDC at connector J6 while pressing the trend button. Black lead to pin 1, red lead to pin 2. <ol style="list-style-type: none"> If voltage is present, replace the motor. If voltage is not present, replace the CPU/Power board.
	HEAD	<ol style="list-style-type: none"> Put the bed into the Bed Calibration menu. <ol style="list-style-type: none"> Using one of the head siderails, push the bed up button. <ol style="list-style-type: none"> If the head lift motor runs up, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> If recalibration does not work, replace the litter angle sensor and recalibrate. If the head lift motor does not run up, check for 24VDC at connector J4 while pressing the trend button. Black lead to pin 1, red lead to pin 2. <ol style="list-style-type: none"> If voltage is present, replace the motor. If voltage is not present, replace the CPU/Power board.
No Bed Down Motion.	FOOT	<ol style="list-style-type: none"> Put the bed into the Bed Calibration menu. <ol style="list-style-type: none"> Using one of the head siderails, push the reverse trend button. <ol style="list-style-type: none"> If the foot lift motor runs down, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> If recalibration does not work, replace the litter angle sensor and recalibrate. If the foot lift motor does not run down, check for 24VDC at connector J6 while pressing the reverse trend button. Red lead to pin 1, black lead to pin 2. <ol style="list-style-type: none"> If voltage is present, replace the motor. If voltage is not present, replace the CPU/Power board.
	HEAD	<ol style="list-style-type: none"> Put the bed into the Bed Calibration menu. <ol style="list-style-type: none"> Using one of the head siderails, push the bed down button. <ol style="list-style-type: none"> If the head lift motor runs down, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> If recalibration does not work, replace the litter angle sensor and recalibrate. If the head lift motor does not run down, check for 24VDC at connector J4 while pressing the trend button. Red lead to pin 1, black lead to pin 2. <ol style="list-style-type: none"> If voltage is present, replace the motor. If voltage is not present, replace the CPU/Power board.

[Return To Table of Contents](#)

Troubleshooting Guide

English

Problem / Failure	Recommended Action
No Fowler Up Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the fowler up button.<ol style="list-style-type: none">I. If the fowler motor runs up, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the fowler angle sensor and recalibrate.II. If the fowler motor does not run up, check for 24VDC at connector J5 while pressing the fowler up button. Black lead to pin 1, red lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.
No Fowler Down Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the fowler down button.<ol style="list-style-type: none">I. If the fowler motor runs down, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the fowler angle sensor and recalibrate.II. If the fowler motor does not run down, check for 24VDC at connector J5 while pressing the fowler down button. Red lead to pin 1, black lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.

Troubleshooting Guide

English

Problem / Failure	Recommended Action
No Gatch Up Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the gatch up button.<ol style="list-style-type: none">I. If the gatch motor runs up, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the gatch angle sensor and recalibrate.II. If the gatch motor does not run up, check for 24VDC at connector J3 while pressing the gatch up button. Red lead to pin 1, black lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.
No Gatch Down Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the gatch down button.<ol style="list-style-type: none">I. If the gatch motor runs down, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the gatch angle sensor and recalibrate.II. If the gatch motor does not run down, check for 24VDC at connector J3 while pressing the gatch down button. Black lead to pin 1, red lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.

[Return To Table of Contents](#)

Troubleshooting Guide

English

Problem / Failure	Recommended Action
No Foot Up Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the foot up button.<ol style="list-style-type: none">I. If the foot motor runs up, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the foot angle sensor and recalibrate.II. If the foot motor does not run up, check for 24VDC at connector J1 while pressing the foot up button. Red lead to pin 1, black lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.
No Foot Down Motion.	<ol style="list-style-type: none">1. Put the bed into the Bed Calibration menu.<ol style="list-style-type: none">A. Using one of the head siderails, push the foot down button.<ol style="list-style-type: none">I. If the foot motor runs down, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17).<ol style="list-style-type: none">a. If recalibration does not work, replace the foot angle sensor and recalibrate.II. If the foot motor does not run down, check for 24VDC at connector J1 while pressing the foot down button. Black lead to pin 1, red lead to pin 2.<ol style="list-style-type: none">a. If voltage is present, replace the motor.b. If voltage is not present, replace the CPU/Power board.

Troubleshooting Guide

English

Problem / Failure	Recommended Action
No Trendelenburg Motion.	<ol style="list-style-type: none"> 1. Check the touch screen trend angle display for accuracy of the level of the litter. <ol style="list-style-type: none"> I. If not accurate, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> a. If recalibration does not work, replace the trend angle sensor and recalibrate. b. If replacement of the trend angle sensor did not resolve the problem, replace the CPU/Power board.
No Reverse Trendelenburg Motion.	<ol style="list-style-type: none"> 1. Check the touch screen trend angle display for accuracy for the accuracy of the level of the litter. <ol style="list-style-type: none"> I. If not accurate, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> a. If recalibration does not work, replace the trend angle sensor and recalibrate. b. If replacement of the trend angle sensor did not resolve the problem, replace the CPU/Power board.
No Cardiac Chair Motion.	<ol style="list-style-type: none"> 1. Check the touch screen fowler and foot section angle display for accuracy for the accuracy of the angle of the fowler and foot section. <ol style="list-style-type: none"> I. If not accurate, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> a. If recalibration does not work, replace the fowler or foot section angle sensor depending on the one which was not accurate, then recalibrate. b. If replacement of the fowler or foot section angle sensor did not resolve the problem, replace the CPU/Power board.
HOB 30°	<ol style="list-style-type: none"> 1. Check the touch screen fowler angle display for accuracy of the angle of the fowler. <ol style="list-style-type: none"> I. If not accurate, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> a. If recalibration does not work, replace the fowler angle sensor and recalibrate. b. If replacement of the fowler angle sensor did not resolve the problem, replace the CPU/Power board.
No Vascular Motion	<ol style="list-style-type: none"> 1. Check the touch screen fowler, gatch, foot, and trend/rev. trend display for accuracy of the angle of all. <ol style="list-style-type: none"> I. If not accurate, recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). <ol style="list-style-type: none"> a. If recalibration does not work, replace the angle sensor of the section that is not accurate and then recalibrate the bed (refer to the Bed Calibration procedures located on page 1-17). b. If replacement of the fowler angle sensor did not resolve the problem, replace the CPU/Power board.

[Return To Table of Contents](#)

Troubleshooting Guide

English

Problem / Failure	Recommended Action
No Electric Brake Motion.	<ol style="list-style-type: none">1. Verify the Brake Not Set LED is flashing and the Brake Set LED is OFF.<ol style="list-style-type: none">I. If the Brake Set LED is ON, check manual brake position.<ol style="list-style-type: none">a. If manual brake pedal is in the brake position, the bed should not move and is okay; go to step 2.b. If manual brake is not in the brake position, check the switch on the patients left side in the middle below the base hood.2. If the brake motor does not run down when the brake button is pushed, check for 24VDC at connector J7 while pressing the brake button. Black lead to pin 1, red lead to pin 2.

SCALE TROUBLESHOOTING

When the Scale System is unable to correctly weigh the patient weight due to a problem with the electronics, Figure 40 will be displayed. It also appears when there is a problem with the Trendelenburg angle sensor; thus the value for the weight and the angle cannot be displayed.

When the weight exceeds 550 lbs., Figure 41 will be displayed. If the weight is less than 2 lbs., the screen will display "0 lb."


Figure 40


Figure 41

When the Trendelenburg or Reverse Trendelenburg angle is above 12° or below -12° , Figure 42 will be displayed.

The display in Figure 43 will appear as long as the correct weight has not been taken and also when the patient moves too much for the weight to be taken properly.


Figure 42


Figure 43

Maintenance Error Messages

ERROR HANDLING

- There are two different CAN networks; each of the networks is divided into two connectors.
- A “safe” error without command can be obtained when a bad connection is carried out on the Power Sensor Board when the control panel is defective (there is no message sent on the network when a button is pressed) or when the network is defective (short-circuit between the signals, open circuit or defective network circuit).
- A command error without “safe” is obtained when a bad connection is done on the Power Sensor board, when the control panel is defective (a message is sent over the network but not the safe signal) or when the safe signal is defective (short-circuit panel or open circuit or safe circuit on the Power Sensor board is broken).
- These errors can be present when a button is pressed, or at any time, and will thus cause the Call Maintenance LED to light up.
- There is a LED on the Power Sensor Board (DS2 SAFE) which is active when there is at least one safe signal which is active. There are also four LEDs on the Power Sensor board, which shows the activity of the network. Every time a message is sent on the network, the LEDs will flash. If the network is defective, the LEDs will remain on or remain off. If there is nothing connected to the network, the LEDs will remain on or off.
- A safe error can be obtained if a command button is pressed without having the Power Sensor Board receiving a safe signal or if a safe signal is received by the Power Sensor board without having any command button pressed.
- The control panel for the weighing scale of the Power Sensor board can also send errors to the micro-controller which will display them. An error of the weighing scale control panel will light on the Call Maintenance LED.

ERROR MESSAGES

Touch Screen Error Messages Name	Definition
Atd Invalid Values	Digital to analog converter is damaged (Replace the Control Board).
Brake pot bad range	Verify if potentiometer is still in place.
Brake unable to elec ctl	Verify wiring to Brake motor and limit switch for manual engage.
Brake pot Disconnected or short	Verify wiring to potentiometer and replace potentiometer. Ensure potentiometer is still in place.
Brake motor time out	Verify wiring to brake motor.
Calibration Error	Previous calibration step performed incorrectly, redo calibration procedure.
Cmd WO safe from nurse SRR	Right siderail outside board has a network communication error (Check network connections, dip-switch configuration of board or if a button is stuck).
Cmd WO safe from nurse SRL	Left siderail outside board has a network communication error (Check network connections, dip-switch configuration of board or if a button is stuck).
Cmd WO safe from nurse SRH	Head end board has a network communication error (Check network connections).
Cmd WO safe from pat pend R	Optional pendant control has a network communication error (Check network connections).
Cmd WO safe from pat pend L	Optional pendant control has a network communication error (Check network connections).
Cmd WO safe from pat pend H	Optional pendant control has a network communication error (Check network connections).
Cmd WO safe from TS	Touch screen has a network communication error (Check network connections).
Cmd WO safe from room	Communication board has a network communication error (Check network connections).

Maintenance Error Messages

English

ERROR MESSAGES (CONTINUED)

Touch Screen Error Messages Name	Definition
GPIO Failure Init	Initialization of the PCA9555 (GPIO expansion chip) failed (Replace control board).
GPIO Failure Read	Reading from the PCA9555 (GPIO expansion chip) failed (Replace control board).
GPIO Failure Write	Writing to the PCA9555 (GPIO expansion chip) failed (Replace control board).
Limit switch head side rail right	Verify wiring to limit switch and replace limit switch.
Limit switch head side rail left	Verify wiring to limit switch and replace limit switch.
Limit switch foot side rail right	Verify wiring to limit switch and replace limit switch.
Limit switch foot side rail left	Verify wiring to limit switch and replace limit switch.
Foot Right Load cell over range	Foot Right load cell or cabling is damaged (Replace load cell).
Head Right Load cell over range	Head Right load cell or cabling is damaged (Replace load cell).
Foot Left Load cell over range	Foot Left load cell or cabling is damaged (Replace load cell).
Head Left Load cell over range	Head Left load cell or cabling is damaged (Replace load cell).
Motor Brake Overheat	Brake motor has ran too long (Leave motor stationary for 54 minutes).
Motor Brake Overload	Brake motor is drawing too many AMPS (remove restriction or replace motor).
Motor Foot Overheat	Foot motor has ran too long (Leave motor stationary for 54 minutes).
Motor Foot Overload	Foot motor is drawing too many AMPS (remove restriction or replace motor).
Motor Gatch Overheat	Gatch motor has ran too long (Leave motor stationary for 54 minutes).
Motor Gatch Overload	Gatch motor is drawing too many AMPS (remove restriction or replace motor).
Motor Head Overheat	Fowler motor has ran too long (Leave motor stationary for 54 minutes).
Motor Head Overload	Fowler motor is drawing too many AMPS (remove restriction or replace motor).
Motor HL Foot Overheat	Foot lift motor has ran too long (Leave motor stationary for 54 minutes).
Motor HL Foot Overload	Foot lift motor is drawing too many AMPS (remove restriction or replace motor).
Motor HL Head Overheat	Head lift motor has ran too long (Leave motor stationary for 54 minutes).
Motor HL Head Overload	Head lift motor is drawing too many AMPS (remove restriction or replace motor).
Motor Zoom® Overheat	Zoom® motor has ran too long (Leave motor stationary for 54 minutes).
Motor Zoom® Overload	Zoom® motor is drawing too many AMPS (remove restriction or replace motor).
No Error	Angle sensor failure or calibration issue.
One Motor Drive Short	Short on drive motor detected.
Safe WO Cmd from A1	Network A1 has a wiring issue (Check cabling).
Safe WO Cmd from A2	Network A2 has a wiring issue (Check cabling).
Safe WO Cmd from B1	Network B1 has a wiring issue (Check cabling).
Safe WO Cmd from B2	Network B2 has a wiring issue (Check cabling).
Scale ADC Error	Scale chip not calibrated. DC Control Board needs to be replaced.
Scale Chip Failure	Control board is bad (Replace Control Board).
Tilt base over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.


[Return To Table of Contents](#)

Maintenance Error Messages

English

ERROR MESSAGES (CONTINUED)

Touch Screen Error Messages Name	Definition
Tilt Error Base	Angle sensor is damaged (Replace sensor).
Tilt Error Foot	Angle sensor is damaged (Replace sensor).
Tilt Error Gatch	Angle sensor is damaged (Replace sensor).
Tilt Error Head	Angle sensor is damaged (Replace sensor).
Tilt Error HiLo Foot	Angle sensor is damaged (Replace sensor).
Tilt Error Trend	Angle sensor is damaged (Replace sensor).
Tilt foot over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.
Tilt gatch over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.
Tilt head over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.
Tilt HiLo foot over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.
Tilt Trend over range	Verify if tilt sensor is still in place. Inspect for a damaged or improperly assembled tilt sensor or recalibrate bed.
Zoom® time out for switch	Verify wiring to Zoom® motor.

CPU/Power Board - QDF75-0440

English

FUSE SPECIFICATION


Location	Description	Amp
F1	30VAC from transformer for the battery charger and the 5VDC / 12VDC supply	10 Amp - 215008.P
F2	24 VDC from batteries for power supply	40 Amp - 142.6185.5402
F3	24 VAC from transformer for power supply	40 Amp - 142.6185.5402

Cable Location	Voltage	Positive Lead	Negative Lead	Description
J11	24-31VAC	Blue	Red	Stepped down voltage from transformer for power supply
J11	30-39VAC	Yellow	Orange	Stepped down voltage from transformer for battery charger and power supply
J12	24-29VDC	Pin 1 - Red	Pin 2 - Black	Battery/Charger
J1	24-28VDC	Pin 1 - Black	Pin 2 - Brown	Foot Actuator Up
J1	24-28VDC	Pin 2 - Brown	Pin 1 - Blue	Foot Actuator Down
J3	24-28VDC	Pin 1 - White	Pin 2 - Black	Gatch Actuator Up
J3	24-28VDC	Pin 2 - Black	Pin 1 - White	Gatch Actuator Down
J6	24-28VDC	Pin 2 - White	Pin 1 - Black	Bed Lift-Foot Up
J6	24-28VDC	Pin 1 - Black	Pin 2 - White	Bed Lift-Foot Down
J4	24-28VDC	Pin 2 - White	Pin 1 - Black	Bed Lift-Head Up
J4	24-28VDC	Pin 1 - Black	Pin 2 - White	Bed Lift-Head Down
J5	24-28VDC	Pin 2 - White	Pin 1 - Black	Fowler Actuator Up
J5	24-28VDC	Pin 1 - Black	Pin 2 - White	Fowler Actuator Down
J7	24-28VDC 24-28VDC 24-28VDC	Pin 1 - Black Pin 1 - Black Pin 1 - Black	Pin 2 - White Pin 2 - White Pin 2 - White	Brake/Brake Off/Drive Actuator • Brake • Brake Off • Drive
J8	Bed Unplugged 25VDC 25VDC	Pin 2 - White Pin 1 - Black	Pin 1 - Black Pin 2 - White	Zoom® Drive Actuator (Model 2152/2154/2156 Only)


Bed Electrical Diagram

English

For Reference Only: Part Number 27-1136 Rev P


[Return To Table of Contents](#)


Bed Electrical Diagram


English


Bed Electrical Diagram

English


Bed Electrical Diagram


Bed Electrical Diagram

English


Bed Electrical Diagram


English

Bed Electrical Diagram

English


English


Bed Electrical Diagram

English

Bed Electrical Diagram


English

[Return To Table of Contents](#)


Bed Electrical Diagram

English

Bed Electrical Diagram

English

[Return To Table of Contents](#)

Bed Electrical Diagram

English

Bed Electrical Diagram - 27-1136 Rev P (reference only)

Item	Part No.	Part Name	Qty.
1	QDF5095	Sound Alarm	1
2	QDF75-0440	DC Power Control	1
3	QDF27-1185	Number 3 Harness	1
4	QDF27-2229	Footboard Brake Cable	1
5	QDF27-2181	Number 2 Harness	1
6	QDF27-2182	Base Structure Extension	1
7	QDF27-1381	Battery Switch Cable	2
8	QDF9188	12V, 17.9 Ah Battery	2
9	QDF27-1208	Limit Switch Siderail Cable	2
10	QDF27-1430	CSI 1109 Board	1
11	QDF27-1097	Brake / Neutral / Drive Board	3
12	QDF27-1382	Speaker	5
13	QDF27-1429	Nurse Call Board	2
14	QDF75-0230	CSI 1157 Board	1
15	QDF9136	Load Cell Handle	1
16	QDF9130	Right Handle	1
17	QDF9131	Left Handle	1
18	QDF8066	Power Cord	2
19	QDF2087	Rocker Switch	1
20	QDF8078	10A (100V, 120V) Fuse	2
21	QDF8068	6.3A (200V, 220V, 240V) Fuse	2
22	QDF27-1607	I.V. Pole and Bed Extender Cable	1
23	QDF27-1682	Nurse Call Board Cable	2
24	QDF8024	120V Auxiliary Outlet	1
25	QDF27-2024	Linear Position Sensor	1
26	QDF9025	Breaker	1
27	QDF27-2214	Footboard Control Board Cable	1
28	QDF27-2038	Toroidal Transformer	1
29	QDF27-1646	Battery Cable	2
30	QDF9004	Micro-Switch	2
31	QDF27-2232	Footboard Cable	1
32	QDF27-1834	Footboard LBS Cable	2
33	QDF27-1841	100V Adapter	1
34	QDF27-1159	SID CAN Board	1
35	QDF27-1372	Load Cell	4
36	QDF21-1151	Exterior Side Head COM CAN-MOD Board	3
37	QDF21-2895	12" Network Cable	4
38	QDF21-1160	Voice Control CAN Board	1
39	QDF75-0600	Without GEN III/37 BRO CAB-CONN Board	1
40	QDF27-1521	Siderail Limit Switch	2
41	QDF27-2213	Number 1 Harness Wire	1
42	QDF27-2230	Touch Screen Cable	2
43	QDF75-0290	Touch Screen	1
44	QDF27-2049	Multivoltage Transformer	1
45	QDF75-0140	Angle Sensor	6
46	QDF27-1156	Cable for Brake Board	2
47	QDF27-2025	12V 3 LEDs Night Light	1
48	QDF27-2432	Cable for Serial Interface Connector	1
49	QDF27-2212	Siderail "Y" Wire	2

Bed Electrical Diagram

English

Bed Electrical Diagram - 27-1136 Rev P (reference only) (Continued)


Item	Part No.	Part Name	Qty.
50	QDF21-4109	20 PINS Auxiliary Cable Port	1
51	QDF9186	Stryker Pendant Port	1
52	QDF2060	Smart TV Board	1
53	QDF21-2901	SIP-N-PUFF Cable	2
54	QDF27-1099	Fowler Side COM CAN-MOD Board	3
55	QDF75-0010	Touch Screen CAN-MENU Board	1
56	QDF9573	XPRT Mattress Outlet	-
57	QDF27-1842	200 V Adaptor	1
58	QDF27-1843	220 V Adaptor	1
59	QDF27-1102	Patient Pendant Control	2
60	---	Patient Pendant Connector	1
61	QDF27-1525	Patient Pendant Cable	2
62	QDF9574	Schurter PEM TYPE KEC #4303.0001	1
63	QDF9575	Schurter PEM Fuse Holder #4303.2001	1
64	QDF27-1840	240 V Adaptor	1
65	QDF27-2253	LBS Cable	1
66	QDF27-1562	Local Bed Status Board	2
67	QDF9571	Medical Filter	1
68	QDF27-1524	Filter and Receptacle Connector	1
69	QDF27-1681	Interior Siderail Wire	2
70	QDF27-1526	Speaker	4
71	---	Patient Pendant Cable (Optional)	1
72	QDF27-2231	USB Cable	1
73	QDF27-1976	Auxiliary Outlet Cable + Mattress	1
74	QDF27-2542	Zoom® Interface Board Harness	1
75	QDF27-0290	Zoom® Interface Board	1

Service Information

BED LIFT ACTUATOR (HEAD) REMOVAL AND REPLACEMENT - (BASE)

Tools Required:

- Needle-Nose Pliers
- Diagonal Pliers
- Jack (if needed)
- 2 x 4 (20" recommended)

Step #8:

Locate Jack Stand between floor and underside of Foot Frame.

Floor

Procedure:

1. Plug bed into wall outlet.
2. Raise the head end siderails to the full up position and set the brakes using one of the brake control locations or by using the manual brake pedal.
3. Remove head end base cover by pulling up and out.
4. Remove center base cover by pulling up and out.
5. Using diagonal pliers, cut the zip tie securing the actuator cable to the base frame.
6. Using Needle-nose pliers, remove the rue clips from the clevis pins securing the actuator.

CAUTION

Do not remove the clevis pins.

7. If the bed will lower down, run it all the way down so the litter is supported by the base litter stop.
8. If the bed will not lower all the way down; using a jack and a 2 x 4, take tension off of the actuator clevis pins by placing the jack and 2 x 4 between the bottom of the right side of the litter and the floor and jack up just enough to take the litter weight off of the base frame (reference picture).

CAUTION

Make sure that prior to jacking, the 2 x 4 is perpendicular to the floor.

9. Using the bed up/down controls, tap the up or down button to remove tension on the clevis pins and remove clevis pins.
10. Unplug the cable quick connect and remove the actuator.
11. Reverse the steps to install the new actuator.
12. Test all bed functionality prior to putting bed back in service.

BED LIFT ACTUATOR (FOOT) REMOVAL AND REPLACEMENT - (BASE)

Tools Required:

- Needle-Nose Pliers
- Diagonal Pliers

Procedure:

1. Plug bed into wall outlet.
2. Raise head end siderails to full up position and ensure brakes are activated.
3. Remove head end base cover by pulling up and out.
4. Remove center base cover by pulling up and out.
5. Using diagonal pliers, cut the zip tie securing the actuator cable to the base frame.
6. Using Needle-nose pliers, remove the rue clips from the clevis pins securing the actuator.

CAUTION

Do not remove the clevis pins.

7. Lower the bed all of the way down so the litter is supported by the base litter stops.
8. Using the bed up/down controls, tap the up or down button to remove tension on the clevis pins and remove clevis pins.
9. Unplug the cable quick connect and remove the actuator.
10. Reverse the steps to install the new actuator.
11. Test all bed functionality prior to putting bed back in service.

Service Information


FOWLER ACTUATOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Diagonal Pliers
- Regular Screwdriver
- Phillips Screwdriver #2
- 3/8" Combination Wrench
- Needle-nose pliers

Procedure:

1. Plug bed into wall outlet.
2. Raise bed to the full up position.
3. Raise the head end siderails to the full up position and set the brakes using one of the brake control locations or by using the manual brake pedal.
4. Working from the right side of the bed, step on the CPR pedal to lower the fowler all the way down.
5. Raise the gatch to the full up position.
6. Remove the three center base hoods (head, center, foot), then remove the right base hood.
7. Using Needle-nose pliers, remove the two rue clips from the clevis pins holding the actuator.
8. Using diagonal pliers, cut the zip ties securing the CPR cable to the actuator and the litter frame. Also cut the zip ties securing the actuator cable to the wire harness.
9. Unplug the actuator from the quick connect.
10. Holding the actuator with one hand remove the clevis pins from the actuator and carefully lower the actuator down.
11. Using a Phillips screwdriver and a 3/8" combination wrench, remove the CPR cable retaining screw from the right CPR pedal.
12. Using a regular screwdriver, loosen the CPR cable retainer on the actuator just enough so it will move.
13. Slide the gray cable retainer lock towards the CPR cable retainer and unclip the CPR cable from the red actuator release activator.
14. Remove actuator.
15. Reverse these steps to reinstall.
16. When re-adjusting the cable retainer lock set, set it to 1.00 ± 0.06 as shown (right). When the actuator is back in place, use needle-nose pliers to reattach the spring to the center of the clevis pin.
17. Test all bed functionality prior to putting bed back into service.

GATCH ACTUATOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Needle-Nose Pliers

Procedure:

1. Plug bed into wall outlet.
2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
3. Remove the mattress and set aside.
4. Raise bed to the full up position.
5. If the actuator will run, lower the gatch down to take the tension off of the actuator mounting pins. If the actuator will not run, support the gatch section to take tension off of the actuator mounting pins.
6. Unplug the actuator cable from the quick connect near the actuator.
7. Using Needle-nose pliers, remove the two rue clips from the clevis pins holding the actuator in place.
8. Holding the actuator with one hand, remove the clevis pins securing the actuator to the bed then remove the actuator.
9. Reverse these steps to reinstall.


Note: If new actuator shaft does not line up with mounting holes, mount the base of the actuator then plug the actuator in and run it electrically in or out until lined up.

10. Test all bed functionality prior to putting bed back in service.

FOOT ACTUATOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Needle-Nose pliers
- Diagonal Pliers
- Phillips Screwdriver #2
- Jack Stand

Step #4:

Locate Jack Stand between floor and underside of Foot Frame.

Floor

Procedure:

1. Plug bed into wall outlet.
2. Set the brakes by using one of the brake control locations or the manual brake pedal.
3. Remove the mattress and set aside.
4. Using a jack stand, support the foot section by lowering the bed height down enough to take the tension off the actuator mounting pins (reference drawing above).
5. Using diagonal pliers, cut the zip ties securing the actuator cable to the bed.
6. Remove the mattress assembly or fold back to expose the foot section.
7. Using a Phillips screwdriver, remove the two screws (reference drawing above, item A and B) securing the electrical cover from the foot section and remove the cover.

Note: Use caution as the cover is large and heavy.

8. Unplug the actuator from J1 and feed cable down to actuator.
9. Using needle-nose pliers, remove the two ratchet clips from the clevis pins holding the actuator in place.
10. Holding the actuator with one hand, remove the clevis pins securing the actuator to the bed then remove the actuator.
11. Reverse these steps to reinstall.

Note: If new motor shaft does not line up with mounting holes, mount the base of the actuator then plug the actuator in and run it electrically in or out until lined up.

12. Test all bed functionality prior to putting bed back in service.

Zoom® DRIVE ACTUATOR REMOVAL AND REPLACEMENT (2152 MODEL ONLY) - (BASE)

Tools Required:

- Needle-Nose Pliers
- Pry Bar or 2 x 4
- Diagonal Pliers

Procedure:

1. Plug bed into wall outlet.
2. Position the bed height in the middle range.
3. Remove the center base cover then remove the Zoom® drive actuator cover and set aside.
4. Using Needle-nose pliers, remove the two rue clips from the clevis pins securing the actuator to the base frame.
5. Using diagonal pliers, cut the zip ties securing the actuator cable and then unplug the actuator from the quick connect.
6. Using a pry bar or a 2 x 4, pry upward on the Zoom® drive actuator frame and push the actuator clevis pins out.
7. Remove actuator.
8. Reverse steps to install new actuator.
9. Test all bed functionality prior to putting bed back into service.

Service Information

CPU / POWER BOARD REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:



- Phillips Screwdriver #2
- ESD System
- 3/8" Nutdriver

Procedure:

1. Plug bed into wall outlet.
2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
3. Raise the bed to the full upright position.
4. Remove the mattress assembly or fold back to expose the foot section.
5. Using a Phillips screwdriver, remove the two screws (reference drawing, item A and B) securing the electrical cover and remove the cover.

Note: Use caution as the cover is large and heavy.

6. Unplug the bed from the wall outlet and turn battery switch OFF. The bed should now have no power.
7. Using an ESD system, properly ground yourself.
8. Unplug all cable connections from the CPU/Power board.
9. Using a 3/8" nutdriver and a Phillips screwdriver, unscrew the four screws (C) securing the metal CPU/Power board mounting plate to the electrical tray (reference drawing above).
10. Remove CPU/Power board and discard.
11. Reverse steps 9-4 to install the new CPU/Power board (QDF75-0440).
12. Plug bed back into wall outlet and turn the battery switch back ON.
13. Configure bed options (refer to the Bed Options Configuration procedure on [page 1-30](#)).
14. Configure the bed serial number (refer to the Serial Number Configuration procedure on [page 1-31](#)).
15. Recalibrate the bed (refer to Bed Calibration procedure on [page 1-17](#)).
16. Test all bed functionality prior to putting the bed back into service.

LOAD CELL (HEAD END) REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- Needle-Nose Pliers
- 3/8" Drive Ratchet
- 1/2" Socket
- 1/2" Combination Wrench


Figure 1


Figure 2

Procedure: (procedure the same for left and right sides)

1. Plug bed into wall outlet.
 2. Raise bed and fowler to full up position, gatch down and foot up.
 3. Raise head end siderails.
 4. Using a Phillips screwdriver, remove the two screws securing the wire cover to the head frame (Figure 1).
 5. Using a Phillips screwdriver, remove the six screws from the head frame cover, lift up on the cover and lay it on the litter (Figure 2)
 6. Unscrew the load cell connection and feed the connection back towards the cell.
 7. Using Needle-nose pliers, remove the rue clip and the flat washer.
 8. Using one hand, grab the upper frame cross bar to remove tension on the litter then push the clevis pin out.
 9. Using a ratchet with a 1/2" socket and a 1/2 combination wrench, remove the two bolts securing the load cell then remove the load cell.
 10. Reverse procedures to install the new load cell.
- CAUTION: Ensure cables do not get pinched when securing the head end fame cover together.**
11. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
 12. Test all bed functionality prior to putting bed back into service.

LOAD CELL (FOOT END) REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- Needle-Nose Pliers
- 1/2" Combination Wrench
- 3/8" Ratchet
- 1/2" Shallow Well Socket
- Diagonal Pliers
- Jack
- 2 x 4 (20" recommended)


Procedure: (procedure the same for left and right sides)

1. Plug bed into wall outlet.
2. Raise bed to around 24" (reference the footboard display) and gatch to full up position.
3. Remove the mattress assembly.
4. Using a 1/2" combination wrench, remove the four bolts securing the gatch section cover and remove by lifting up on the foot end and then pushing backwards on the cover.
5. Raise the foot end siderails up.
6. Using needle-nose pliers, remove the ratchet clip and flat washer.
7. Remove the foot end base cover.
8. Using a Phillips screwdriver, remove the four screws securing the foley bag bracket.
9. Using a jack and 2 x 4, take the tension off of the load cells and remove the clevis pins from both foot end load cells.
10. Unscrew load cell cable connector.
11. Using diagonal pliers, cut the two cable ties securing the main wire harness and the siderail cable located just behind the load cell.
12. Using a ratchet and 1/2" socket, remove the two nuts and bolts securing the load cell to the litter frame.
13. Remove the load cell by pulling back and upward towards the middle of the bed while feeding the cable through the litter frame.
14. Reverse procedures to install the new load cell.
CAUTION: When installing new load cell do not pinch the cable with the load cell.
15. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
16. Test all bed functionality prior to putting bed back into service.

DISPLAY REMOVAL AND REPLACEMENT - (FOOTBOARD)

Tools Required:

- Phillips Screwdriver #2
- Stubby Phillips Screwdriver #2
- Small Regular Screwdriver
- Utility Knife
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and level.
3. Unplug the bed and turn the battery disconnect "OFF" (0).
4. Using a small regular screwdriver or a utility knife, remove the rear footboard label.


NOTE: You will need to replace this label: QDF27-2756.

5. Using a Phillips screwdriver, remove the six screws from the top of the control panel which the label was covering.
6. Raise the display up and using a stubby Phillips screwdriver, remove the three screws on the bottom.
7. Remove the display housing and unplug the two cables from the function control board and set the display assembly on a bench to work on.
8. Using a Phillips screwdriver, remove the four screws securing the display to the display housing.
9. Unplug the ribbon cable going to the touch screen by pulling outward on the black locking tab and then pulling out on the ribbon cable.
10. Unplug the two cables going to the display, making note of their location and orientation.
11. Carefully work your way around the support bracket display, pulling upward evenly; then remove the display.
12. Reverse the procedures to install the new display.
13. Recalibrate the touch screen (refer to Touch Screen Calibration procedures on [page 1-29](#)).
14. Test all bed functionality prior to putting the bed back into service.

BRAKE CONTROL BOARD REMOVAL AND REPLACEMENT - (FOOTBOARD)

Tools Required:

- Phillips Screwdriver #2
- Stubby Phillips Screwdriver #2
- Small Regular Screwdriver
- Utility Knife
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and level.
3. Unplug the bed and turn the battery disconnect "OFF" (0).
4. Using a small regular screwdriver or a utility knife, remove the rear footboard label.


NOTE: You will need to replace this label: QDF27-2756.

5. Using a Phillips screwdriver, remove the six screws from the top of the control panel which the label was covering.
6. Raise the display up and using a stubby Phillips screwdriver, remove the three screws on the bottom.
7. Remove the display housing and unplug the two cables from the function control board and set the display assembly on a bench to work on.
8. Using a Phillips screwdriver, remove the six screws securing the brake control board.
9. Unplug the cable going to the brake control board.
10. Reverse the procedures to install the new brake control board.
11. Test all bed functionality prior to putting the bed back into service.

FUNCTION SELECTION/LED BOARD REMOVAL AND REPLACEMENT - (FOOTBOARD)

Tools Required:

- Phillips Screwdriver #2
- Stubby Phillips Screwdriver #2
- Small Regular Screwdriver
- Utility Knife
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and level.
3. Unplug the bed and turn the battery disconnect "OFF" (0).
4. Using a small regular screwdriver or a utility knife, remove the rear footboard label.

NOTE: You will need to replace this label: QDF27-2756.


5. Using a Phillips screwdriver, remove the six screws from the top of the control panel which the label was covering.
6. Raise the display up and using a stubby Phillips screwdriver, remove the three screws on the bottom.
7. Remove the display housing and unplug the two cables from the function control board and set the display assembly on a bench to work on.
8. Using a Phillips screwdriver, remove the four screws securing the function selection/LED board.
9. Unplug the cable going to the function selection/LED board.
10. Reverse the procedures to install the new function selection/LED board.
11. Test all bed functionality prior to putting the bed back into service.

Service Information

TOUCH SCREEN REMOVAL AND REPLACEMENT - (FOOTBOARD)

Tools Required:

- Phillips Screwdriver #2
- Stubby Phillips Screwdriver #2
- Small Regular Screwdriver
- Utility Knife
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and level.
3. Unplug the bed and turn the battery disconnect "OFF" (0).
4. Using a small regular screwdriver or a utility knife, remove the rear footboard label.

NOTE: You will need to replace this label: QDF27-2756.

5. Using a Phillips screwdriver, remove the six screws from the top of the control panel which the label was covering.
6. Raise the display up and using a stubby Phillips screwdriver, remove the three screws on the bottom.
7. Remove the display housing and unplug the two cables from the function control board and set the display assembly on a bench to work on.
8. Using a Phillips screwdriver, remove the four screws securing the display to the display housing.
9. Unplug the ribbon cable going to the touch screen by pulling outward on the black locking tab and then pulling out on the ribbon cable.
10. Unplug the three cables going to the display, making note of their location and orientation.
11. Carefully work your way around the support bracket display, pulling upward evenly; then remove the display.
12. Using a small regular screwdriver or utility knife, remove the foot control label.

NOTE: You will need to replace the following labels:

* **Without /Bed Awareness** - QDF27-2188-ENG; * **With /Bed Awareness** - QDF27-2188-ENG

13. From the back of the display housing push out on the touch screen and remove.
14. Reverse the procedures to install the new touch screen.
15. Recalibrate the touch screen (refer to Touch Screen Calibration procedures on [page 1-29](#)).
16. Test all bed functionality prior to putting the bed back into service.

[Return To Table of Contents](#)

BATTERY REMOVAL AND REPLACEMENT - (FOOTBOARD)

Tools Required:

- Phillips Screwdriver #2
- Stubby Phillips Screwdriver #2
- Small Regular Screwdriver
- Utility Knife
- ESD System

Procedure:


1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and level.
3. Unplug the bed and turn the battery disconnect switch to the "OFF" (0) position.
4. Using a small regular screwdriver or a utility knife, remove the rear footboard label.

NOTE: You will need to replace this label: QDF27-2756.

5. Using a Phillips screwdriver remove the six screws from the top of the control panel which the label was covering.
6. Raise the display up and using a stubby Phillips screwdriver, remove the three screws on the bottom.
7. Remove the display housing and unplug the two cables from the function control board and set the display assembly on a bench to work on.
8. Remove the battery and replace with new battery (see figure below).
9. Reverse the above steps to reassemble.
10. Referencing the operation manual complete the Time / Date Setup procedure.


Battery

BRAKE / NEUTRAL / DRIVE POTENTIOMETER REMOVAL AND REPLACEMENT - QDF27-2024

Tools Required:

- Diagonal Pliers

Procedure:

1. Plug bed into wall outlet.
2. Raise the bed lift and patient's left siderails all the way up.
3. Remove the base center covers (foot, center, and head).
4. Remove the patient's left base frame cover.
5. Cut the two wire ties securing the potentiometer to the potentiometer bracket; then cut the two wire ties securing the potentiometer wires to the frame.
6. Remove the electrical tape holding the quick connection together.
7. Reverse the above procedures for installation of the new potentiometer.

NOTE: Make sure to install new wire ties in the same location.

8. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
9. Test all bed functionality prior to putting the bed back into service.

BATTERY REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- 5/16" Nutdriver
- 5/16" Combination Wrench


Figure 2


Figure 1

Procedure:

Note

Batteries should always be replaced in pairs.

1. Plug bed into wall outlet.
2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
3. Raise the bed to the full up position.
4. Remove the mattress assembly or fold back to expose the foot section.
5. Using a Phillips screwdriver, remove the two screws securing the electrical cover and remove the cover (reference Figure 1 above, item A and B).

Note: Use caution as the cover is large and heavy.

6. Unplug the bed from the wall outlet and turn battery switch "OFF" (0). The bed should now have no power.
7. Unplug J12 connector from the CPU/Power board.
8. Lift up on both batteries and stand them upright (reference Figure 2 above).
9. Using a 5/16" nutdriver and a 5/16" combination wrench, unbolt the wires from the battery posts.
10. Reverse the steps to install new batteries.

Note: New battery information should be facing the right side when installing and upright when installed.


11. Test all bed functionality while plugged into wall outlet and when unplugged prior to putting back into service.

Service Information

FOWLER ANGLE SENSOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Stubby Phillips Screwdriver #2
- ESD System


Procedure:

1. Plug bed into wall outlet.
 2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
 3. Raise the bed to the full up position and raise the patient left siderails (head end and foot end).
 4. If fowler motor will run, raise fowler up to approximately 20 degrees. If the fowler motor will not run, put the bed into the calibration mode step one as shown on [page 1-17](#). Run the Fowler up to approximately 20 degrees.
 5. Working from the left side of the bed, use a stubby Phillips screwdriver and remove the two screws securing the fowler angle sensor to the bottom of the fowler frame.
 6. Using an ESD system, properly ground yourself.
 7. Unclip the three clips holding the board cover on.
 8. Unplug the cable from the board.
 9. Reverse the steps to install the new fowler angle sensor.
- Note: Do not over-tighten the two screws.**
10. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
 11. Test all bed functionality prior to putting the bed back into service.

GATCH ANGLE SENSOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Stubby Phillips Screwdriver #2
- ESD System

Procedure:

1. Plug bed into wall outlet.
 2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
 3. Raise the bed to the full up position.
 4. If gatch motor will run, raise gatch up to approximately 20 degrees. If the gatch motor will not run, then put the bed into the calibration mode step one as shown on [page 1-17](#).
 5. Working from the bottom left side of the litter, under the gatch section, use a stubby Phillips screwdriver and remove the two screws securing the gatch angle sensor to the bottom of the gatch frame.
 6. Using an ESD system, properly ground yourself.
 7. Unclip the three clips holding the board cover on.
 8. Unplug the cable from the board.
 9. Reverse the steps to install new gatch angle sensor.
- Note: Do not over tighten the two screws.**
10. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
 11. Test all bed functionality prior to putting the bed back into service.

Service Information

FOOT ANGLE SENSOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- ESD System


Figure 2



Figure 1


Procedure:

1. Plug the bed into a wall outlet.
 2. Raise the bed to the full up position and level the litter surface.
 3. Remove or fold back the mattress, exposing the foot section.
 4. Using a Phillips screwdriver, remove the two screws securing the electrical cover and then remove the cover (refer to Figure 1 above, item A and B).
 5. Using a Phillips screwdriver, remove the two screws securing the angle sensor to the foot frame (refer to Figure 2 above). **Note: Use caution as the cover is large and heavy.**
 6. Using an ESD system, properly ground yourself.
 7. Unclip the three clips holding the board cover on.
 8. Unplug the cable from the board.
 9. Reverse the steps to install the new foot angle sensor.
- Note: Do not over-tighten the two screws.**
10. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
 11. Test all bed functionality prior to putting the bed back into service.

BASE ANGLE SENSOR REMOVAL AND REPLACEMENT - (BASE)

Tools Required:

- Stubby Phillips Screwdriver
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Set the brakes using one of the brake control locations or by using the manual brake pedal.
3. Raise the bed to the full up position and raise the patient right siderails (foot end and head end).
4. Remove the head end, enter, foot end, and right base cover.
5. Using a stubby Phillips screwdriver, remove the two screws securing the base angle sensor to the base frame (reference Figure above).
6. Using an ESD system, properly ground yourself.
7. Unclip the three clips holding the board cover on.
8. Unplug the cable from the board.
9. Reverse the steps to install new base angle sensor.


Note: Do not over-tighten the two screws.

10. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
11. Test all bed functionality prior to putting the bed back into service.

TREND ANGLE SENSOR REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position and raise the patient right siderail.
3. Working from the patient's right side, use a Phillips screwdriver and remove the two screws securing the trend angle sensor to the litter frame (refer to Figure above),
4. Using an ESD system, properly ground yourself.
5. Unclip the three clips holding the board cover on.
6. Unplug the cable from the board.
7. Reverse the steps to install new trend angle sensor.


Note: Do not over-tighten the two screws.

8. Recalibrate the bed (refer to Bed Calibration procedures on [page 1-17](#)).
9. Test all bed functionality prior to putting the bed back into service.

HEADWALL COMMUNICATION BOARD REMOVAL AND REPLACEMENT - (LITTER)

Tools Required:

- Phillips Screwdriver #2
- ESD System

Procedure:

1. Plug the bed into a wall outlet.
2. Raise the bed to the full up position.
3. Remove the head board and any accessories located on the head end.
4. Using a Phillips screwdriver, remove the six screws securing the head frame cover (reference figure above).
5. Using a Phillips screwdriver, remove the three screws securing the headwall board to the head frame (reference figure above).
6. Using an ESD system, unplug the seven cables (make note of their location; they are not all labeled). Compare the old board and new board and mirror all of the dip-switch settings from the old board to the new board. Verify if any jumpers need to be moved.
7. Reverse procedures to install new headwall communication board.
8. Test all bed functionality prior to putting the bed back into service.

INTOUCH® CRITICAL CARE BED

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
The InTouch® Critical Care Bed is suitable for use in the electromagnetic environment specified below. The customer or the user of the InTouch® Critical Care Bed should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	+2 kV for power supply lines +1 kV for input/output lines	+2 kV for power supply lines +1 kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 8 kV differential mode ± 2 kV common mode	± 8 kV differential mode ± 2 kV common mode	Main power quality is that of typical commercial and/or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	<5% U_t (>95% dip in U_t) for 0,5 cycle 40% U_t (60% dip in U_t) for 5 cycles 70% U_t (30% dip in U_t) for 25 cycles. <5% U_t (>95% dip in U_t) for 5 sec.	<5% U_t (>95% dip in U_t) for 0,5 cycle 40% U_t (60% dip in U_t) for 5 cycles 70% U_t (30% dip in U_t) for 25 cycles. <5% U_t (>95% dip in U_t) for 5 sec.	Main power quality should be that of a typical commercial and/or hospital environment. If the user of the InTouch® Critical Care Bed requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.
Note: U_T is the a.c. mains voltage prior to applications of the test level.			

INTOUCH® CRITICAL CARE BED - (CONTINUED)

Recommended separation distances between portable and mobile RF communications equipment and the InTouch® Critical Care Bed.			
The InTouch® Critical Care Bed is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the InTouch® Critical Care Bed can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the InTouch® Critical Care Bed as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d=1,2\sqrt{P}$	$d=1,2\sqrt{P}$	$d=2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

INTOUCH® CRITICAL CARE BED - (CONTINUED)

The InTouch® Critical Care Bed is suited for use in the electromagnetic environment specified below. The customer or the user of the InTouch® Critical Care Bed should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the InTouch® Critical Care Bed, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p style="text-align: center;">Recommended Separation Distance</p> $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ <p style="text-align: center;">80 MHz to 800 MHz</p> $d=2,3\sqrt{P}$ <p style="text-align: center;">800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

Note 1

At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the InTouch® Critical Care Bed is used exceeds the applicable RF compliance level above, the InTouch® Critical Care Bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the InTouch® Critical Care Bed.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

INTOUCH® CRITICAL CARE BED - (CONTINUED)

Guidance and Manufacturer's declaration - Electromagnetic Emissions		
The InTouch® Critical Care Bed is intended for use in an electromagnetic environment specified below. The customer or the user of the InTouch® Critical Care Bed should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	The InTouch® Critical Care Bed uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The InTouch® Critical Care Bed is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Warranty

LIMITED WARRANTY

Stryker Medical Division, a division of Stryker Corporation, warrants to the original purchaser the *InTouch®* Critical Care Bed, Model FL27 (2151/2152, 2153/2154, 2155/2156) be free from defects in material and workmanship for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product which is, in the sole discretion of Stryker, found to be defective. If requested by Stryker, products or parts for which a warranty claim is made shall be returned prepaid to the factory. Any improper use or any alteration or repair by others in such manner as in Stryker's judgment affects the product materially and adversely shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

Stryker Medical Bed products are designed for a 10 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its Bed products will be free from structural defects for the expected 10 year life of the Bed product as long as the original purchaser owns the product.

This statement constitutes Stryker's entire warranty with respect to the aforesaid equipment. **Stryker makes no other warranty or representation, either expressed or implied, except as set forth herein. There is no warranty of merchantability and there are no warranties of fitness for any particular purpose. In no event shall Stryker be liable here under for incidental or consequential damages arising from or in any manner related to sales or use of any such equipment.**

TO OBTAIN PARTS AND SERVICE

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative, or call Stryker Customer Service USA at 1-800-327-0770, option 2, Canada 1-888-233-6888.

RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. **Special, modified, or discontinued, items not subject to return.**

DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full. Claims for any short shipment must be made within thirty (30) days of invoice.

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



Parts
Ersatzteile
Onderdelen
Peças

Use 2 Ty-Raps to attach the case of the potentiometer to the steel support. (See Drawing 27-1660)

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
To assemble the casters, direct the smile (chrome part) towards the line of the bed.

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


The angle sensor orientation is very important. Wire should exit towards the bottom.

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

Parts
Ersatzteile
Onderdelen
Peças

A

To assemble casters, direct the smile (chrome part) toward the left side of the bed.

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

Parts
Ersatzteile
Onderdelen
Peças

**A**


To assemble casters, direct the smile (chromed part) toward the left side of the bed.

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

Parts
 Ersatzteile
 Onderdelen
 Peças

Note

To avoid any risk of pinching when assembling these cables:

- the bed must be in high position
- the Gatch must be at 15 degrees
- the Foot Litter must be at 45 degrees

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




















Parts
Ersatzteile
Onderdelen
Peças

Parts
Ersatzteile
Onderdelen
Peças



Sensor orientation is important!
The wire must leave downwards.

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











Parts
Ersatzteile
Onderdelen
Peças

The sensor orientation is important!
(The wire must leave downwards.)

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







Parts
Ersatzteile
Onderdelen
Peças


CABLE CONNECTIONS QDF27-2248

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-2248	J17	To	QDF27-1960	J17 DC Card
QDF27-2248	J17A	To	Box Hole	(27-1200)


Parts
Ersatzteile
Onderdelen
Peças

Mattress Connector
Matratzenanschluss
Matrasconnector
Conector do colchão


Mattress Connector - Groove towards inside
Matratzenanschluss – Kerbe nach innen
Matrasconnector – groef naar binnen
Conector do colchão – ranhura virada para o interior

1	QDF27-2248	1
2	QDF9518	2
3	VV83A9G16	3
4	VV37A1C12	2
5	QDF9025	1
6	QDF9573	1
7	27-1967P	1


8	QE71-0964-tri	1 (English/Français/Español)
	QE71-0964-ger	1 (Deutsch) or
	QE71-0964-ita	1 (Italiano)
	QE71-0964-dut	1 (Nederlands)
	QE71-0964-por	1 (Português)
	QE71-0964-pol	1 (Polski)
9	QE71-1348-tri	1 (English/Français/Español)
	QE71-1348-ger	1 (Deutsch) or
	QE71-1348-ita	1 (Italiano)
	QE71-1348-dut	1 (Nederlands)
	QE71-1348-por	1 (Português)
	QE71-1348-pol	1 (Polski)

[Return To Table of Contents](#)

Parts
Ersatzteile
Onderdelen
Peças

1	27-2687	1
2	27-2695	1

Parts
Ersatzteile
Onderdelen
Peças


Note: To avoid all risks of cable jamming during assembly ensure the bed is in the highest position, the head litter is at 70 degrees, the knee gatch litter is at 15 degrees and the foot litter is at 45 degrees.

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

QDF27-2025#12 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-2025	Light	To	QDF75-0440	J20 DC Board

QDF27-2213 #2 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-2213	J22	To	QDF75-0440	J22 DC Board
QDF27-2213	J21	To	QDF75-0440	J21 DC Board
QDF27-2213	J20A	To	QDF75-0440	J20A DC Board
QDF27-2213	J20B	To	QDF75-0440	J20B DC Board
QDF27-2213	J20C	To	QDF75-0440	J20C DC Board
QDF27-2213	J19	To	QDF75-0440	J19 DC Board
QDF27-2213	J29	To	QDF75-0440	J29 DC Board
QDF27-2213	J15A	To	QDF75-0440	J15A DC Board
QDF27-2213	J15B	To	QDF75-0440	J15B DC Board
QDF27-2213	J15C	To	QDF75-0440	J15C DC Board
QDF27-2213	J15D	To	QDF75-0440	J15D DC Board
QDF27-2213	J204	To	QDF27-2038	Transfo J204
QDF27-2213	J27	To	QDF75-0440	J27 DC Board
QDF27-2213	J3	To	QDF75-0440	J3 DC Board
QDF27-2213	J5	To	QDF75-0440	J5 DC Board
QDF27-2213	J113	To	QDF75-0140	Foot Sensor (27-2688)
QDF27-2213	J202	To	---	Gatch Ground
QDF27-2213	J114	To	QDF75-0140	Trend Sensor
QDF27-2213	J153	To	QDF27-1372	Right Foot Cell
QDF27-2213	J103	To	27-2598	Removable Control
QDF27-2213	J172	To	QDF27-1208	Right Foot Siderail
QDF27-2213	J102	To	QDF27-2212	Right Head Siderail
QDF27-2213	J152	To	QDF27-1372	Right Head Cell
QDF27-2213	J202A	To	---	Foot Ground
QDF27-2213	J144	To	QDF27-1215	Gatch Motor
QDF27-2213	J142	To	QDF27-1214	Fowler Motor
QDF27-2213	J112	To	QDF75-0140	Fowler Sensor
QDF27-2213	J118	To	QDF75-0140	Gatch Sensor
QDF27-2213	J155	To	QDF27-1372	Left Foot Cell
QDF27-2213	J171	To	QDF27-1208	Left Foot Siderail
QDF27-2213	J105	To	QDF27-2212	Left Head Siderail
QDF27-2213	J150	To	QDF27-1372	Left Head Cell
QDF27-2213	J23	To	QDF75-0440	J23 DC Board
QDF27-2213	J26	To	QDF75-0440	J26 DC Board
QDF27-2213	J200A Black	To	QDF9571	Line "L"
QDF27-2213	J200B White	To	QDF9571	Line "N"
QDF27-2213	J200C Green/Yellow	To	Ground	Head Screw

Parts
Ersatzteile
Onderdelen
Peças

QDF27-1607 #3 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-1607	Con. MTA 5 Pos.	To	QDF75-0440	J25 DC Board
QDF27-1607	Contact	To	QDF27-1213	Contact 7 or 8
QDF27-1607	Contact	To	QDF27-1213	Contact 7 or 8
QDF27-1607	Terminal	To	QDF9535	Right Limit Switch No/Com
QDF27-1607	Terminal	To	QDF9535	Right Limit Switch Com/No
QDF27-1607	Terminal	To	QDF9535	Right Limit Switch No/Com
QDF27-1607	Terminal	To	QDF9535	Right Limit Switch Com/No

QDF5095 #5 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF5095	Con. MTA	To	QDF75-0440	J13 DC Board
QDF5095	Buzzer	To	Electrical Box	See Position on Drawing

QDF27-1524 #7 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-1524	Black 90° Terminal	To	QDF9571	Line "L"
QDF27-1524	White 90° Terminal	To	QDF9571	Line "N"
QDF27-1524	Green 90° Terminal	To	QDF9571	Line "E"
QDF27-1524	Black Right Terminal	To	QDF9574	"L" (Power Supply)
QDF27-1524	White Right Terminal	To	QDF9574	"N" (Power Supply)
QDF27-1524	Green Right Terminal	To	QDF9574	"E" (Power Supply)
QDF27-1524	Cross Grommet	To	Head Ground	Screw
QDF27-1524	Cross Grommet	To	Head Ground	Screw

QDF27-2214 #10 CONNECTION

Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-2214	Con. MTA 6 Pos.	To	QDF75-0440	J18 DC Board
QDF27-2214	Con. Metrimate	To	27-1547	Panel (27-1547)

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
























DIP SWITCH SETUP ON QDF27-1099			
SW13	OFF	ON	ON
SW14	OFF	OFF	-

CABLES CONNECTION				
Cable No.	Connector No.	To	Cable No.	Connector No.
QDF27-1681	MTA – 3 Positions	To	QDF27-1099	J5
QDF27-1156	MTA – 6 Positions	To	QDF27-1099	J3
QDF27-1156	MTA – 6 Positions	To	QDF27-1097	J1
QDF27-1784	MTA – 3 Positions	To	QDF27-2212	MTA-2 Positions
QDF27-2212	Mate N Lock 12 Pos.	To	QDF27-2213	J105 (L27-038)
QDF27-2212	MTA – 6 Positions	To	QDF27-1099	J6
QDF27-2212	MTA – 3 Positions	To	See OL270164, OL270175	
QDF27-2212	MTA – 2 Positions	To		
QDF27-2212	Connector 2 Pins	To		

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


















Position on Wire 27-1834	Position on Connector ④	Position on Wire 27-1208	Signal
1	1	1	LED 1
2	2	2	LED 2
3	3	3	+12V


1	QDF27-2030-eng	1 (English) or	3	QDF27-1834	2
	QDF27-2030-bil	1 (English/Français) or	4	QDF2080	2
	QDF27-2030-spa	1 (Español) or	5	QDF27-1562	4
	QDF27-2030-ger	1 (Deutsch) or	6	QP27-1831	2
	QDF27-2030-dut	1 (Nederlands) or	7	VV23A9C12HL	4
	QDF27-2030-por	1 (Português) or			
2	QDF27-2031-eng	1 (English) or			
	QDF27-2031-bil	1 (English/Français) or			
	QDF27-2031-spa	1 (Español) or			
	QDF27-2031-ger	1 (Deutsch) or			
	QDF27-2031-dut	1 (Nederlands) or			
	QDF27-2031-por	1 (Português) or			

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
















































Manufactured for
UNITED STATES
Stryker Medical
3800 E. Centre Ave.,
Portage, Michigan USA
49002

EC	REP
----	-----

Stryker France S.A.S.
ZAC – Avenue de Satolas Green
69881 MEYZIEU Cedex
France

stryker[®]