

# LIFEPAK<sup>®</sup>15

## Performance Inspection Procedure (PIP) Checklist



Model # \_\_\_\_\_ Department/Location \_\_\_\_\_  
 Serial # \_\_\_\_\_ Performed By \_\_\_\_\_  
 Type of PIP \_\_\_\_\_ Post-Repair  Annual  Date \_\_\_\_\_

### Manual Mode Access

1. Manual Mode Access
  - a. Record customer-selected MANUAL ACCESS configuration

### Exterior Physical Inspection

#### 2. Exterior physical inspection

	Pass	Fail	NA	Comments
a. Device exterior damage (general)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Check device for loose/rattling hardware	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
c. Check for damaged or missing rubber feet	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
d. Inspect battery pins as specified in the Service Manual	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
e. Check if battery pins were replaced during this servicing event	<input type="checkbox"/>	Battery Pins Replaced		_____
f. Inspect therapy cable pins and connector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
g. Confirm spring button on therapy connector is functional	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
h. Inspect device connectors for damage	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
i. Inspect keypads and overlays for damage	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
j. Check device accessories for condition and expiration dates	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
k. Inspect carrying case and carrying strap for damage	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____

#### 3. Device Setup

- a. Insert two fully charged Li-ion batteries into the device
- b. Install a roll of 100-mm printer paper
- c. Connect therapy cable or standard paddles to the device  **Completed** \_\_\_\_\_

#### 4. Power On/Self –Test

- a. All items are conforming    \_\_\_\_\_

#### 5. Auxiliary Power Switching (if Auxiliary Power Connector is installed)

- a. Battery icons appear but neither is highlighted.    \_\_\_\_\_

# LIFEPAK<sup>®</sup>15

## Performance Inspection Procedure (PIP) Checklist



6. Power Source Management				
a. Confirm battery status indicator switching	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
7. User Test and Date/Time Verification				
a. Confirm device passes User Test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Confirm Time and Date are correct	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
<b>Miscellaneous Function</b>				
8. Temperature Calibration Check Test (if Temp option is installed)	<b>Pass</b>	<b>Fail</b>	<b>NA</b>	<b>Comments</b>
a. Confirm Temperature Cal Check complete	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
9. CO2 Tests (if CO2 option is installed)				
a. Confirm change in vacuum reading is less than 15 mmHg	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Record CO2 concentration reading is 5.0% ±0.5%	Measured Value _____		<input type="checkbox"/>	_____
10. NIBP Tests (if NIBP option is installed)				
a. Confirm LEAKAGE TEST OK message	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Confirm 50 mmHg readings agree within ±20 mmHg	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
c. Confirm 150 mmHg readings agree within ±20 mmHg	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
d. Confirm the overpressure switch activates at 290 ±20 mmHg	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
PPV 6SHHG3ULQ7HWW				
D Confirm printer test strip and CHECK PRINTER message	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
PPV 6SHHG3ULQ7HWW				
D Confirm printer 12.5 mm/s test strip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
.H \SDGHVW				
D Confirm all control text boxes are highlighted and TEST COMPLETE message appears	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
14. Audio Test				
a. Confirm voice messages and tones are clear and not distorted	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
15. Invasive Blood Pressure Verification (if IP option is installed)				
a. Confirm P1 pressure channel zero	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Record P1 pressure reading of 250 ±8 mmHg	Measured Value _____		<input type="checkbox"/>	_____
c. Record P1 pressure reading of 100 ±5 mmHg	Measured Value _____		<input type="checkbox"/>	_____
d. Record P1 pressure reading of 20 ±3 mmHg	Measured Value _____		<input type="checkbox"/>	_____

# LIFEPAK<sup>®</sup>15

## Performance Inspection Procedure (PIP) Checklist



- |  |  |                          |       |
|--|--|--------------------------|-------|
| e. Record P1 pressure reading of -20 ±3 mmHg | Measured Value _____   | <input type="checkbox"/> | _____ |
| f. Confirm P2 pressure channel zero          | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |                          | _____ |
| g. Record P2 pressure reading of 250 ±8 mmHg | Measured Value _____   | <input type="checkbox"/> | _____ |
| h. Record P2 pressure reading of 100 ±5 mmHg | Measured Value _____   | <input type="checkbox"/> | _____ |
| i. Record P2 pressure reading of 20 ±3 mmHg  | Measured Value _____   | <input type="checkbox"/> | _____ |
| j. Record P2 pressure reading of -20 ±3 mmHg | Measured Value _____   | <input type="checkbox"/> | _____ |

16. SpO2/SpCO/SpMet Tests	Pass	Fail	NA	Comments
a. Confirm SpO2 reading is between 50% and 100% (if SpO2 is installed)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Confirm SpCO reading is between 0% and 40% (if SpCO is installed)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
c. Confirm SpMet reading is between 0% and 15% (if SpMet is installed)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____

### 17. Record Operating Data (Optional)

#### Total Shocks:

		Fault Messages	_____
		Power Cycle Count	_____
360J Shocks	<input type="text"/>	Pacing Count	_____
	_____	Shock Count	_____
225-325J Shocks	<input type="text"/>	Power On Time	_____
	_____	Printer On Time	_____
0-200J Shocks	<input type="text"/>	SPO2 Operating Time (if installed)	_____
	_____	CO2 Operating Time (if installed)	_____
		NIBP Inflation Cycles (if installed)	_____

### ECG Performance Testing

18. ECG Tests (12-lead, 3-lead or 5-wire ECG tests)	Pass	Fail	NA	Comments
a. Confirm LEADS-OFF screen messages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. Record Lead I gain (tolerance 25 to 31 mm)	Measured Value _____		<input type="checkbox"/>	_____
c. Record Lead II gain (tolerance 36 to 44 mm)	Measured Value _____		<input type="checkbox"/>	_____
d. Record Lead V1/C gain (tolerance 36 to 44 mm) (5-wire, 12-lead)	Measured Value _____		<input type="checkbox"/>	_____

e. Record Lead V2 gain (tolerance 36 to 44 mm) (12-lead)	Measured Value _____	<input type="checkbox"/>	_____
f. Record Lead V3 gain (tolerance 36 to 44 mm) (12-lead)	Measured Value _____	<input type="checkbox"/>	_____
g. Record Lead V4 gain (tolerance 36 to 44 mm) (12-lead)	Measured Value _____	<input type="checkbox"/>	_____
h. Record Lead V5 gain (tolerance 36 to 44 mm) (12-lead)	Measured Value _____	<input type="checkbox"/>	_____
i. Record Lead V6 gain (tolerance 36 to 44 mm) (12-lead)	Measured Value _____	<input type="checkbox"/>	_____

## 19. ECG Analog Output (optional, perform as required)

a. Record signal amplitude (tolerance 0.90 to 1.10 Vp-p)	Measured Value _____	<input type="checkbox"/>	_____
--	----------------------	--------------------------	-------

## Defibrillator/Pacing Testing

## 20. Delivered Energy Test

	Pass	Fail	NA	Comments
a. 10 J – Record delivered energy (tolerance 9.1 to 10.9 J)	Measured Value _____		<input type="checkbox"/>	_____
b. 200 J – Record delivered energy (tolerance 186.0 to 214.0 J)	Measured Value _____		<input type="checkbox"/>	_____
c. 360 J – Record delivered energy (tolerance 334.9 to 384.9 J)	Measured Value _____		<input type="checkbox"/>	_____

## 21. Charge Time to 360J Test

a. Confirm device charges to 360 J in less than 10 seconds	Measured Value _____		<input type="checkbox"/>	_____
--	----------------------	--	--------------------------	-------

## 22. Synchronous Cardio version Test

a. Record Sync delay (maximum 60ms)	Measured Value _____		<input type="checkbox"/>	_____
-------------------------------------	----------------------	--	--------------------------	-------

## 23. Therapy ECG Characteristics

a. Record ECG paddle lead gain (tolerance 1mV = 36 to 44 mm)	Measured Value _____		<input type="checkbox"/>	_____
b. Fast-Restore baseline in 0.5 seconds	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
c. Fast-Restore amplitude restored is >50% within 3 seconds	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
d. Positive R-wave test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____

## 24. Standard Paddles User Test (N/A for QUIK-COMBO-only device)

a. Confirm device passes test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
-------------------------------	--------------------------	--------------------------	--------------------------	-------

## 25. Pacer Option Characteristics

a. Confirm leads-off detection	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	_____
b. 10 mA– Record current (tolerance 5 to 15 mA)	Measured Value _____		<input type="checkbox"/>	_____
c. 100 mA – Record current (tolerance 91 to 109 mA)	Measured Value _____		<input type="checkbox"/>	_____
d. 200 mA – Record current (tolerance 181 to 219 mA)	Measured Value _____		<input type="checkbox"/>	_____
e. Record pulse width (tolerance 19.2 to 20.8 ms)	Measured Value _____		<input type="checkbox"/>	_____

**26. Patient Impedance Test**

- |  |                          |                          |                          |       |
|--|--------------------------|--------------------------|--------------------------|-------|
| a. Verify the PADDLES LEADS OFF message is not visible (50 ohms)   | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | _____ |
| b. Verify the device displays PADDLES LEADS OFF message (370 ohms) | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | _____ |
| c. Verify the PADDLES LEADS OFF message is not visible (238 ohms)  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | _____ |

**Data Management**

- |   |                          |                          |                          |       |
|---|--------------------------|--------------------------|--------------------------|-------|
| <b>27. Bluetooth Wireless Technology</b> (if Bluetooth option is installed) | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | _____ |
| a. Verify Bluetooth Pairing Successful                                      |                          |                          |                          |       |

**Leakage Current Test**

**28. Leakage Test Battery Powered**

- |  | <b>Pass</b>          | <b>Fail</b> | <b>NA</b>                | <b>Comments</b> |
|--|----------------------|-------------|--------------------------|-----------------|
| a. ECG Direct Applied Part at 120 or 240 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 45 µA)       | Measured Value _____ |             | <input type="checkbox"/> | _____           |
| b. Therapy Direct Applied Part at 120 or 240 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 2625 µA) | Measured Value _____ |             | <input type="checkbox"/> | _____           |
| c. SpO2 Direct Applied Part at 120 or 240 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 2625 µA)    | Measured Value _____ |             | <input type="checkbox"/> | _____           |

**29. Leakage Test AC Powered Device at 120VAC** (If Aux power is installed)

- |  |                      |  |                          |       |
|--|----------------------|--|--------------------------|-------|
| a. Direct Equipment Leakage at 120 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Open Earth</b> , (15 µA - 270 µA) | Measured Value _____ |  | <input type="checkbox"/> | _____ |
| b. ECG Direct Applied Part at 120 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 45 µA)        | Measured Value _____ |  | <input type="checkbox"/> | _____ |
| c. Therapy Direct Applied Part at 120 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 2625 µA)  | Measured Value _____ |  | <input type="checkbox"/> | _____ |
| d. SpO2 Direct Applied Part at 120 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 2625 µA)     | Measured Value _____ |  | <input type="checkbox"/> | _____ |

**30. Leakage Test AC Powered Device at 240 VAC** (if Aux power is installed)

- |  |                      |  |                          |       |
|--|----------------------|--|--------------------------|-------|
| a. Direct Equipment Leakage at 240 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Open Earth</b> , (15 µA - 450 µA) | Measured Value _____ |  | <input type="checkbox"/> | _____ |
| b. ECG Direct Applied Part at 240 VAC<br>Polarity <b>NC/RM</b> , Condition <b>Normal</b> , (5 µA - 45 µA)        | Measured Value _____ |  | <input type="checkbox"/> | _____ |

**LIFEPAK®15**  
**Performance Inspection Procedure (PIP) Checklist**



- c. Therapy Direct Applied Part at 240 VAC  
Polarity **NC/RM**, Condition **Normal**, (5  $\mu$ A - 2625  $\mu$ A) Measured Value \_\_\_\_\_  \_\_\_\_\_
- d. SpO2 Direct Applied Part at 240 VAC  
Polarity **NC/RM**, Condition **Normal**, (5  $\mu$ A - 2625  $\mu$ A) Measured Value \_\_\_\_\_  \_\_\_\_\_

---

**31. LIFEPAK 15 Maintenance Instruction**

**Pass**

- a. Maintenance prompt disabled or reset  **Completed** \_\_\_\_\_

**Comments:**