

**POWER PRO TL**  
**Ambulance Cot**  
**Model 6550 RUGGED®**

# stryker®

## Operations/Maintenance Manual



For Parts or Technical Assistance:  
USA: 1-800-327-0770 (option 2)  
Canada: 1-888-233-6888

CE



# Table of Contents

---

Introduction . . . . .	6
Specifications . . . . .	6
WARNING / CAUTION / NOTE DEFINITION. . . . .	7
Symbols and Definitions . . . . .	9
Warranty . . . . .	10
DeWALT® Product Warranty . . . . .	10
Stryker EMS Return Policy. . . . .	11
Return Authorization. . . . .	11
Damaged Merchandise . . . . .	11
International Warranty Clause. . . . .	11
Patent Information . . . . .	11
Summary of Safety Precautions . . . . .	12
Component Identification . . . . .	16
Product Inspection . . . . .	17
General Inspection. . . . .	17
Setup Procedures. . . . .	19
Cot Fastener Installation . . . . .	20
Cot Features . . . . .	22
Using Patient Restraint Straps . . . . .	22
Using The Patient Restraint Belt Extension (Optional Equipment). . . . .	24
Pedi-Mate™ Infant Restraint System (Optional Equipment) Attachment Instructions . . . . .	25
Operating the Wheel Lock(s) . . . . .	27
Adjusting the Wheel Locking Force. . . . .	28
Adjusting the Footrest and Knee Gatch. . . . .	29
Operating the Backrest . . . . .	30
Operating the Siderails. . . . .	30
Proper Battery Management: DeWALT® Battery Management and Charger Care . . . . .	31
Proper Battery Management: DeWALT® Battery Frequently Asked Questions . . . . .	32
Battery Power Indicator . . . . .	33
Hour Usage Meter . . . . .	34
Installing the Backrest Storage Pouch (Optional Equipment) . . . . .	35
Using the Equipment Hook (Optional Equipment). . . . .	36
Operating the Rigid Push Bars and Storage Pouch (Optional Equipment). . . . .	37
Operating the 2-Stage I.V. Pole (Optional Equipment) . . . . .	38
Operating the 3-Stage I.V. Pole (Optional Equipment) . . . . .	39
Cot Operation. . . . .	40
Operating Guidelines . . . . .	40
Proper Lifting Techniques. . . . .	40
Raising and Lowering the Cot (unoccupied) . . . . .	41
Raising and Lowering the Cot (occupied) . . . . .	41
Transporting a Patient . . . . .	41
Transferring the Patient to the Cot . . . . .	42
Using the <b>▲RUGGED▲</b> Transfer Flat (Optional Equipment) . . . . .	42

# Table of Contents

---

Ambulance Cot Motion . . . . .	42
Tail Lift Loading and Unloading . . . . .	43
Ramp Loading and Unloading . . . . .	44
Engaging/Releasing the Cot into the Fastener . . . . .	45
High Speed Retract/Extend . . . . .	46
Caster Steer Lock . . . . .	46
Manual Raising and Lowering Of The Cot . . . . .	47
Adjusting Cot Height . . . . .	48
Using Additional Assistance . . . . .	48
Operating the Retractable Head/Foot Sections . . . . .	49
Battery Operation . . . . .	50
Hydraulic Sub-Assembly (6550-001-030) Access Instructions . . . . .	51
Cleaning . . . . .	52
Washing Procedure . . . . .	52
Washing Limitations . . . . .	52
Removal of Iodine Compounds . . . . .	53
Preventative Maintenance . . . . .	54
Regular Inspection and Adjustments . . . . .	55
Maintenance Record . . . . .	57
Training Record . . . . .	58
Cot Assembly . . . . .	59
Base Assembly . . . . .	71
Telescoping Foot End . . . . .	79
Telescoping Head Section . . . . .	81
Fowler Assembly . . . . .	83
Gatch Assembly . . . . .	84
Telescoping Gatch Assembly . . . . .	89
Hitch Assembly, Foot End . . . . .	90
Hitch Assembly, Head End . . . . .	94
Mounted Hydraulics Assembly . . . . .	97
Hydraulics Sub-Assembly . . . . .	99
Outer Rail, Patient Right . . . . .	101
Outer Rail, Patient Left . . . . .	102
Inner Lift Tube Assembly . . . . .	103
Outer Lift Tube Assembly . . . . .	104
Inner Lift Tube Assembly, Litter Pivot, Patient Right . . . . .	105
Inner Lift Tube Assembly, Litter Pivot, Patient Left . . . . .	106
Caster, Steer Lock, Wheel Lock Assembly . . . . .	107
Corner Handle Assembly . . . . .	108
2-Stage I.V. Assembly, Patient Left - 6550-211-000 . . . . .	111
2-Stage I.V. Pole Assembly, Patient Left - 6500-001-042 . . . . .	112
2-Stage I.V. Assembly, Dual - 6550-212-000 . . . . .	113
3-Stage I.V. Assembly, Patient Right - 6550-215-000 . . . . .	114



# Table of Contents



---







3-Stage I.V. Pole Assembly, Patient Right - 6500-001-043	115
3-Stage I.V. Assembly, Patient Left - 6550-216-000	116
3-Stage I.V. Pole Assembly, Patient Left - 6500-001-044	117
3-Stage I.V. Assembly, Dual - 6550-217-000	118
No I.V. Pole Assembly Option - 6550-218-000	119
Optional Push Bar Assembly - 6550-040-000	120
Trouble Shooting Guide	121
POWER PRO Electronics and Hydraulics Locator	121
POWER PRO Hydraulic Manifold Components Locator	122
POWER PRO Electrical System Block Diagram	123
POWER PRO Troubleshooting Guide Table of Contents	125
POWER PRO Troubleshooting Guide	126
Recycling Passport.	130
Battery Pack: 6500-700-006.	130
Switch Housing Assembly: 6550-001-036	131
Foot End Assembly: 6550-001-022.	132
Electronics Assembly: 6550-001-014	133
Hydraulic Power Unit: 6550-001-030.	134
Main Cable: 6550-001-172	135
Quick Reference Replacement Part List	136
EMC Information.	137
POWER PRO TL	137

# Introduction

This manual is designed to assist you with the operation and maintenance of the **RUGGED** POWER PRO TL Ambulance Cot. Read it thoroughly before using the equipment or beginning any maintenance on it.

## SPECIFICATIONS

Maximum Cot Load Capacity 	50 stones	317.5 kg	700 pounds
Maximum Unassisted Lift Capacity 	50 stones	317.5 kg	700 pounds
Backrest Articulation/Shock Position	0° to 73°		
Overall Length/Minimum Length/Width	194.3 cm / 139.7cm / 58.4 cm	76.5" / 55" / 23"	
Height <sup>1</sup>	Adjustable from:	48.2 cm to 96.5 cm	19.5" to 38"
Weight <sup>2</sup>	10.35 stones	67 kg	145 pounds
Caster Diameter/Width	15.2 cm / 5.1 cm	6" / 2"	
Minimum Operators Required for Loading/ Unloading an Occupied Cot	2		
Minimum Operators Required for Loading/ Unloading an Unoccupied Cot	1		
Recommended Fastener Systems	Model 6385, 6386, 6387		
Double Wheel Lock / Four Wheel Lock	2 Wheel Lock Standard, 4 Optional		
Hydraulic Oil	Stryker Part Number 6500-001-293		
Battery System	DeWALT® 24 Volt NiCd Battery		
- Battery	Stryker Part Number 6500-700-006		
- Charger	230 Volt - Stryker Part Number 6500-071-000		
Standards (Cots and Chargers)	IEC 60601-1 CAN/CSA-C22.2 No. 601.1-M90 UL 60601-1 IEC 60601-1-2:2001 BS EN:1789		

Environmental Conditions	Operation	Storage and Transportation
Temperature		
Relative Humidity		
Atmospheric Pressure		

<sup>1</sup> Height measured from bottom of mattress at seat section to ground level.

<sup>2</sup> Cot is weighed with 1 battery and without mattress and restraints.

Stryker reserves the right to change specifications without notice.

The POWER PRO TL is designed to be compatible with competitive cot fastener systems.

DeWALT® is a registered trademark of Black & Decker Inc.

Patents pending.

# Introduction

---

## **WARNING / CAUTION / NOTE DEFINITION**

The words WARNING, CAUTION and NOTE carry special meanings and should be carefully reviewed.

---

### **WARNING**

Alerts the reader about a situation, which if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

---

### **CAUTION**

Alerts the reader of a potentially hazardous situation, which if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

### **NOTE**

This provides special information to make maintenance easier or important instructions clearer.


# Introduction

---

Contact Stryker Customer Service or Technical Support at: +44 (0) 1635 262431.

Stryker UK Limited  
Hambridge Rd.  
Newbury, Berkshire  
RG14 5EG, England

**Please have the serial number of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.**



**Figure 1 - Cot Serial Number & Location**

# Symbols and Definitions

---



Consult accompanying documentation



Dangerous Voltage Symbol



Pinch Point



Extend



Retract



Type B Equipment: equipment providing a particular degree of protection against electric shock, particularly regarding allowable leakage current and reliability of the protective earth connection.

Internally Powered Equipment: equipment able to operate from an internal (removable) electrical power source.

Mode of Operation: 10% (33 sec. On / 5 Min. Off)

**IPX6**

Protection from powerful jets of water.



Safe Working Load Symbol



The DeWALT® Battery System is a registered trademark of Black & Decker Inc.



In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.



Medical Equipment Classified by Underwriters Laboratories Inc. with Respect to Electric Shock, Fire, Mechanical and Other Specified Hazards Only in Accordance with UL 60601-1, First Edition (2003) and CAN/CSA C22.2 No. 601.1-M90 with updates 1 and 2.



BS EN 1789:2000 Certified.

# Warranty

**Stryker EMS**, a division of the Stryker Corporation, offers two distinct warranty options in the United States:

**One (1) year parts and labor.** Under this option, Stryker EMS warrants to the original purchaser that its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective.

**Two (2) year parts.** Under this option, Stryker EMS warrants to the original purchaser that non-expendable components of its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of two (2) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts for, or replacing, at its option, any product which is, in the sole discretion of Stryker, found to be defective. Expendable components, i.e. mattresses, restraints, IV poles, storage nets, storage pouches, O2 straps, and other soft goods, have a one (1) year limited warranty with this option.

Under either warranty option, Stryker Medical EMS products are designed for a 7 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its EMS products will be free from structural defects for the expected 7 year life of the EMS product as long as the original purchaser owns the product.

If Stryker requests, products or parts for which an original purchaser makes a warranty claim, the purchaser shall return the product or part prepaid freight to Stryker's factory.

Any improper use or alteration or repair by unauthorized service providers in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

This statement constitutes Stryker EMS's entire warranty with respect to the aforesaid equipment. STRYKER MAKES NO OTHER WARRANTY OR REPRESENTATION EITHER EXPRESSED OR IMPLIED, EXCEPT AS SET FORTH HERIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL STRYKER BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

## **DeWALT® PRODUCT WARRANTY**

Any DeWALT® product purchased from Stryker EMS is covered for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective due to faulty materials or workmanship.

# Warranty

## STRYKER EMS RETURN POLICY

Cots, Stair Chairs, Evacuation Chairs, Cot Fasteners and Aftermarket Accessories may be returned up to 180 days of receipt if they meet the following guidelines:

### Prior to 30 Days

- 30 day money back guarantee in effect
- Stryker EMS is responsible for all charges
- Returns will not be approved on modified items

### Prior to 90 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 10% restocking fee

### Prior to 180 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 25% restocking fee

## RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items.

**SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.**

## DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Claims for any short shipment must be made within thirty (30) days of invoice.

## INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

## PATENT INFORMATION

**▲RUGGED▲** products are manufactured under the following patents:

United States	5,575,026	6,276,010	6,648,343	6,908,133
	5,537,700	6,125,485	6,735,794	

Other Patents Pending

# Summary of Safety Precautions

---

The following is a list of safety precautions that must be observed when operating or servicing this unit. The precautions are repeated throughout the manual, where applicable. Carefully read this list before using or servicing the unit.

---

## WARNING

- Do not modify the **▲RUGGED▲** ambulance cot or any components of the cot, including the hydraulic unit. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product will also void its warranty.
- It is the responsibility of the cot operator to ensure the ambulance cot being used in the Stryker cot fastener system meets the installation specifications listed on [page 20](#). Injury may result if a non-compatible ambulance cot is used in the Stryker cot fastener system.
- To avoid accidental release of the Pedi-Mate™, and possible injury to the infant, ensure that the restraint buckle is located away from obstructions on the ambulance cot or accessories.
- Never apply the optional wheel lock(s) while a patient is on the ambulance cot. Tipping could occur if the ambulance cot is moved while wheel locks are applied, resulting in injury to the patient or operator and/or damage to the product.
- Never install or use wheel locks on an ambulance cot with excessively worn wheels. Installing or using wheel locks on wheels with less than a 6" (16 cm) diameter could compromise the holding ability of the wheel lock, resulting in injury to the patient or operator and/or damage to the product or other equipment.
- To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
- Do not remove the battery when the ambulance cot is activated.
- Avoid contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.
- Improper usage of the **▲RUGGED▲** ambulance cot can cause injury to the patient or operator. Operate the ambulance cot only as described in this manual.
- Entanglement in powered ambulance cot mechanisms can cause serious injury. Operate the ambulance cot only when all persons are clear of the mechanisms.
- Practice changing height positions and loading and unloading the ambulance cot until operation of the product is fully understood. Improper use can cause injury.
- Do not allow untrained assistants to assist in the operation of the ambulance cot. Untrained technicians/assistants can cause injury to the patient or themselves.
- Do not ride on the base of the **▲RUGGED▲** ambulance cot. Damage to the product could occur, resulting in injury to the patient or operator.
- Grasping the **▲RUGGED▲** ambulance cot improperly can cause injury. Keep hands, fingers and feet away from moving parts.
- Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. **Transporting the cot in a lowered position, head or foot end first, will minimize the potential of a cot tip.**
- Always use all restraint straps to secure the patient to the ambulance cot. An unrestrained patient may fall from the ambulance cot and be injured.
- Do not attach restraints to the base or cross tubes, improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.
- Never leave a patient unattended on the ambulance cot or injury could result. Hold the ambulance cot securely while a patient is on the cot.
- Siderails are not intended to serve as a patient restraint device. Refer to pages 24 and 25 for proper restraint strap usage. Failure to utilize the siderails properly could result in patient injury.
- High obstacles such as curbing, steps or rough terrain can cause the ambulance cot to tip, possibly causing injury to the patient or operator. **Transporting the cot in lower positions can reduce the potential of a cot tip.** If possible, obtain additional assistance ([see page 48](#) for a reference chart) or take an alternate route. Once the weight is off the ground, the operator(s) must support the load of the patient, ambulance cot and any accessories. Failure to support the load properly may cause injury to the patient or operator.
- Explosion Hazard - Do not use cot in presence of flammable anesthetics.
- Do not attempt to operate the ambulance cot when loaded into a cot fastener. Damage to the product or injury to the patient or operator may occur.



# Summary of Safety Precautions

---

## WARNING (CONTINUED)

---

- Use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the ambulance cot.
- Failure to properly clean or dispose of contaminated mattress or other ambulance cot components will increase the risk of blood borne pathogens and may cause injury to the patient or operator.
- Hydraulically raising or lowering the cot may temporarily affect electronic patient monitoring equipment. For best results, patient monitoring should be conducted when the cot is idle.
- **SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY.** If the products described above are used to clean Stryker patient care equipment, measures must be taken to insure the cots are wiped with clean water and thoroughly dried following cleaning. Failure to properly rinse and dry the cots will leave a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.
- Escaping fluid under pressure can penetrate the skin causing serious injury. Avoid the hazard by relieving pressure before disconnecting hydraulic or other lines. Tighten all connections before applying pressure. If an accident occurs, see a doctor immediately. Any fluid injected into the skin must be surgically removed within a few hours or gangrene may result. Doctors unfamiliar with this type of injury should reference a knowledgeable medical source.
- Ensure the push bar is properly locked into its mount to prevent the bar from coming out during use and possibly injuring the operator or patient.
- Do not lift the cot with the push bars. Lifting the cot by the push bars may result in failure of the push bar lock mechanism and may injure the operator or patient.
- When not in use, the push bars should be stored in the base storage pouch provided with the push bar option.
- A lift gate stop which is not properly functioning can result in injury to the patient or operator; ensure that the cot is unable to roll back off the lift before using the lift with a cot and patient. Verify the lift gate stop is maintained and functioning properly and ensure that the cot is secure at all times when on the tail lift.
- It is important to ensure that the push bar is properly locked into its mount to prevent the bar from coming out during use and possibly injuring the operator or patient.
- Medical electrical equipment requires special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided on [page 137](#) to prevent equipment malfunction.
- Portable and mobile RF communication equipment can affect Medical Electrical Equipment.

## CAUTIONS

---

- Ensure that restraints are not entangled in the base frame when raising and lowering the cot.
- Wheel locks are only intended to help prevent the ambulance cot from rolling while unattended. Wheel locks may not provide sufficient resistance on all surfaces or under loads.
- Only use the battery and charger as specified.
- The POWER PRO TL is not for use with an AC adaptor.
- Ensure that the battery is charged prior to placing into service. An uncharged or depleted battery may cause poor ambulance cot performance.
- A preventative maintenance program should be established for all Stryker EMS equipment. Preventative maintenance may need to be performed more frequently based on the usage level of the product. Close attention should be given to safety features including, but not limited to:
  - Hydraulic power mechanism
  - All electrical controls return to off or neutral position when released
  - For additional maintenance information, refer to the preventative maintenance section.
- Do not store items under the ambulance cot mattress. Storing items under the mattress can interfere with the operation of the ambulance cot.
- The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 1.4 stones (20 pounds), (9 kg).
- To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 1.4 stones (20 pounds), (9 kg).
- The weight of the I.V. bags or equipment must not exceed 2.9 stones (40 pounds), (18 kg) for usage on the I.V. poles (if equipped).

# Summary of Safety Precautions

---

## CAUTIONS (CONTINUED)

---

- Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.
- Loading, unloading or changing the position of a loaded ambulance cot requires a minimum of **two trained operators**. The operator(s) must be able to lift the total weight of the patient, cot and any other items on the cot (if additional assistance is needed, see Using Additional Assistance under the Cot Operations section).
- Remove the battery if the cot is not going to be used for an extended period of time (over 24 hours).
- Remove the battery before washing the cot.
- Do not steam clean or ultrasonically clean this ambulance cot. Maximum water temperature should not exceed 180 °F/82 °C. Maximum air dry temperature (cart washers) is 240 °F/115 °C. Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 61 centimeters from the unit. Towel dry all casters and interface points. Failure to comply with these instructions may invalidate any/all warranties.
- Improper maintenance can cause injury or damage to the product. Maintain the ambulance cot as described in this manual. Use only Stryker approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty.
- Failure to use authorized parts, lubricants, etc. could cause damage to the ambulance cot and will void the warranty of the product.
- Hydraulic hose connections and lines can fail due to physical damage, kinks, age, and exposure. Check hoses and lines regularly to avoid damage to the cot. Check and tighten loose connections.
- Hydraulic fluid connections can loosen due to physical damage and vibration. Check and tighten loose connections.
- When charging batteries in an ambulance vehicle, locate the charger either in the forward cab or an enclosed compartment (i.e. cabinet).
- The unit must be in its fully down position in order to properly engage the cot fastener (not included). In order to ensure the cot is in the lowest position press the retract (-) button until the cot stops moving downward. This should be done prior to attempting to engage the cot into the cot fastener (not included). **In order to reduce risk of damage to the cot or cot fastener, do not attempt to activate the cot height activation while it is engaged in the cot fastener.**
- When not in use, the push bars should be stored in the base storage pouch provided with the push bar option.
- To avoid malfunction, the model 6550 **▲RUGGED▲** POWER PRO TL ambulance cot should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the model 6550 **▲RUGGED▲** POWER PRO TL ambulance cot should be observed to verify normal operation in the configuration in which it will be used.
- Possible fire hazard when used with oxygen administering equipment of other than the nasal, mask, or 1/2 bed length tent type. Oxygen tent should not extend below mattress support level. An oxygen rich environment is defined as an environment where the oxygen concentration levels are greater than 25% for ambient pressures up to 110 kPa or the partial pressure of oxygen is greater than 27.5 kPa at ambient pressures exceeding 110kPa.

## NOTES

- This manual should be considered a permanent part of the ambulance cot and should remain with the product even if the cot is subsequently sold.
- These are general instructions for installation of the Pedi-Mate™. Safe and proper use of the Pedi-Mate™ is solely at the discretion of the user. Stryker recommends all users be trained on the proper use of the Pedi-Mate™ before using it in an actual situation.
- The operators must lift the cot slightly off the ground to use the manual extend or retract while a patient is on the cot.
- Keep your spare battery on the charger at all times. Batteries slowly lose power when not on the charger.
- Failure to follow the cleaning directions when using the specified types of cleaners may void this product's warranty.
- Activation of the manual release may cause the ambulance cot to drop slowly if less than 9.2 stones (130 pounds), (59 kg) are on the cot.
- When operating the manual release, avoid rapid lifting or lowering of the base or movement may appear sluggish; lift with a slow constant motion.

# Summary of Safety Precautions




Figure 2 - Potential Pinch Points

# Component Identification




Figure 3 - Cot Components

# Product Inspection

## GENERAL INSPECTION

The condition of the ambulance cot is the responsibility of the owner. It is important that the ambulance cot is working properly before the product is put into service. Have a qualified service person use the following list and the operation instructions to check the ambulance cot before the product is put into service.

Unpack the cartons and check all items for proper operation.

Item	Routine	Page(s)
Battery	Unpack batteries and charger Charge battery according to DeWALT® instructions	P/N 6500-001-206

### CAUTION

When charging batteries in an ambulance vehicle, locate the charger either in the forward cab or an enclosed compartment (i.e. cabinet).

Prior to checking the features and condition of the cot, the battery must be charged until the red LED lights continuously to ensure a full charge. Refer to the DeWALT® instructions for further charging information.

Once the battery is fully charged, inspect the ambulance cot for the following points:

Item	Routine	Page(s)
Battery	Charge spare battery (if necessary) according to DeWalt® instructions.	P/N 6500-001-206
	Install battery into foot-end enclosure, battery indicator operates.	50
	Ensure the battery remains firmly secured.	50
	Release and remove battery from foot-end enclosure.	50
	Reinstall battery into foot-end enclosure.	50
Hydraulics	Inspect motor mount, all fasteners secure	131
	Check cylinder attachments at both ends, all fasteners secure.	121
	Inspect main cable, all connections secure.	121
	Inspect hoses and cylinder seal for leaks.	121
Electronic Controls	Check battery indicator, charged.	33
	Extend cot to raised position.	48
	Check high speed retract.	46
	Extend cot to full height, no drift.	48
Manual Release	Verify the manual release lever functions properly, adjust accordingly.	47, 129
	With the cot empty, check the raise/lower function.	47
	With the cot loaded with a minimum of 9.2 stones (130 lbs), (59 kg), check the raise/lower function.	47
Litter	All fasteners secure (reference all assembly drawings).	
	All welds intact, not cracked or broken.	
	No bent, broken, or damaged components	
	Verify siderails operate and latch properly.	30
	Verify backrest cylinder operates properly through range of motion.	30
	Verify the leg rest operates properly.	29
	Install body restraints. Restraints intact and operating properly.	22-23
	No rips or tears in mattress cover.	

# Product Inspection

Item	Routine	Page(s)
Head/Foot Section	All fasteners secure (reference all assembly drawings).	
	No bent or broken tubing or sheet metal.	
	Verify the head/foot section extends and retracts properly.	49
Base	All fasteners secure (reference all assembly drawings).	
	All welds intact, not cracked or broken.	
	No bent, broken, or damaged components.	
Wheels and Tires	No debris in wheels.	
	All wheels secure, rolling and swiveling properly.	
	Operate wheel locks; wheel secure when engaged, rolls freely when disengaged.	27-28
Accessories	Verify I.V. Pole (if equipped) operates properly.	38-39
	Verify Pedi-Mate™ restraint package (if equipped) intact.	25-26
	Verify equipment hook (if equipped) is installed properly.	36
	Verify pocketed backrest storage pouch (if equipped) is installed properly.	35
	Verify 36" restraint extender (if equipped) is included.	23
	Verify the Bariatric transfer flat (if equipped) is included.	42

# Setup Procedures

---

The patient compartment of the vehicle in which the ambulance cot will be used must have:

- A level floor large enough for the folded cot.
- Stryker 6385 cot fastener (not included).

When necessary, modify the vehicle to fit the cot. Do not modify the cot.

---

 **WARNING**

- Do not modify the **ARUGGED** POWER PRO TL ambulance cot or any components of the cot, including the hydraulic unit. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product will also void its warranty.
- Refer to DeWALT® manual (Stryker part number 6500-001-206) for battery and charger operation.

**NOTE**

- This manual should be considered a permanent part of the ambulance cot and should remain with the product even if the cot is subsequently sold.
- Stryker continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your ambulance cot and this manual.
- If you have any questions, please contact your Stryker Customer Service or Technical Support representative.

# Cot Fastener Installation

---

The **▲RUGGED▲** POWER PRO TL Cot Fastener system, model 6385, 6386, and 6387 are designed to be compatible only with cots which conform to the Cot Fastener Installation specifications listed on the following page.

The **▲RUGGED▲** POWER PRO TL Cot Fastener system, model 6385, 6386 and 6387 are not designed for any purpose other than to restrict movement of an ambulance cot being transported in the patient compartment of an ambulance under normal conditions. Usage of this product in any other way becomes the complete responsibility of the owner/user. Caution must be used at all times during placement of the cot into the ambulance.

---

 **WARNING**

It is the responsibility of the cot operator to ensure that the cot being used in the **▲RUGGED▲** POWER PRO TL Cot Fastener system, model 6385, 6386, or 6387 meets the installation specifications listed on the following page. Injury may result if a non-compatible cot is used in the **▲RUGGED▲** POWER PRO TL Cot Fastener system, model 6385, 6386 or 6387.

**NOTE**

For more detailed instruction and operation instructions for the **▲RUGGED▲** POWER PRO TL Cot Fastener system, model 6385, 6386 or 6387, refer to the 6385 **▲RUGGED▲** Cot Fastener Installation/Operation Manual.



# Cot Fastener Installation




Figure 6 - Installation Specifications

# Cot Features

## USING PATIENT RESTRAINT STRAPS

### WARNING

Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.

1. Secure the patient on the cot with **all** the restraint straps.
2. Buckle the restraints across the patient's chest/shoulders, waist and legs (Figure 14).
3. Keep the restraint straps buckled when the cot is not being used with a patient to avoid damage to the buckles and straps.

### WARNING

Do not attach restraints to the base or cross-tubes, improper restraint attachment could result in damage to the cot further resulting in injury to the patient or operator.




Figure 14 - Safety Restraints




Figure 15 - Head Section Restraints




Figure 16 - Foot Section Restraints

To attach the patient restraint straps to the cot litter frame, feed the loop end of each strap through the litter frame at each of the designated locations in Figure 15 and 16 and then feed the buckle or latch plate back through the loop end of the strap. The arrows indicate alternate attachment areas.

When attaching the restraint straps to the cot, remember the attachment points should provide both strong anchorage and proper restraint position while not interfering with equipment and accessories.

# Cot Features

## USING PATIENT RESTRAINT STRAPS (CONTINUED)




Figure 17 - Buckling the Safety Restraints




Figure 18 - Lengthening the Safety Restraint

### CAUTION

Ensure that restraints are not entangled in the base frame when raising and lowering the cot.




Figure 19 - Shortening the Safety Restraint

When the cot is put into service, open the restraints and place them at either side of the cot until the patient is positioned on the cot mattress. Lengthen the restraints, buckle them around the patient and shorten them until the required tightness is achieved.

#### To open the restraint:

1. Press the red button (A) on the front of the buckle “receiver”.
2. This releases the buckle latch plate (B) which can then be pulled out of the receiver (Fig. 17).

#### To close the restraint:

1. Push the latch plate into the receiver until a “click” is heard.
2. When fastening the chest restraint ensure the latch plate passes through both links (C) on the shoulder strap (Fig. 17).

#### To lengthen the restraint:

1. Grasp the buckle latch plate, turn it at an angle to the webbing, then pull it out (Fig. 18).
2. A hemmed tab at the end of the webbing prevents the latch plate from coming off the strap.

#### To shorten the restraint:

1. Grasp the hemmed tab, pull the webbing back through the latch plate until the required tightness is achieved (Fig. 19).

#### Note


- Whenever a restraint is buckled on a patient, verify the latch plate is fully engaged and any extra webbing is not tangled in the cot or hanging loose.
- Inspection of the restraints should be done **at least** once a month (more frequently if used heavily). Inspection should include checking for a bent or broken receiver or latch plate, torn or frayed webbing, etc. Any restraint showing wear or not operating properly **must** be replaced immediately.

# Cot Features

---

## USING THE PATIENT RESTRAINT BELT EXTENSION (OPTIONAL EQUIPMENT)

Use the restraint belt extension (Stryker part number 6082-160-050) for extra length when buckling the lap belt around large patients.



**Figure 20 - Attaching the Restraint Belt Extension**

# Cot Features

---

## PEDI-MATE™ INFANT RESTRAINT SYSTEM (OPTIONAL EQUIPMENT) ATTACHMENT INSTRUCTIONS

Refer to the Pedi-Mate™ users manual for the manufacturer's recommendations for the use, operation and care of the Pedi-Mate™ Infant Restraint System.

Securing the Pedi-Mate™ to the cot:

1. Remove any restraints already attached to the cot.
2. Raise the cot backrest to the full upright position.
3. Position the Pedi-Mate™ pad flat on the backrest with the black backrest straps out (see Figure 21).



Figure 21 - Positioning the Pedi-Mate™

4. Wrap the straps around the backrest and insert the ends of the straps through the brackets. Securely fasten the buckle (see Figure 22).




Figure 22 - Fastening the Pedi-Mate™ Buckle

---

### **WARNING**

To avoid accidental release of the Pedi-Mate™, and possible injury to the infant, ensure the buckle is located away from obstructions on the cot or accessories.

# Cot Features

---

## PEDI-MATE™ INFANT RESTRAINT SYSTEM ATTACHMENT INSTRUCTIONS (CONTINUED)

5. Pull firmly on the end of the adjustable backrest strap and tighten it securely.
6. Insert the mainframe straps between the cot frame and the mattress. To ensure the release button is toward the foot end of the cot, insert the buckle behind the litter crossbrace and bring it up in front of the crossbrace. Secure the buckle around the crossbrace, leaving a little slack in the strap for final adjustment (see Figure 23).



Figure 23 - Securing the Safety Restraints on a **RUGGED** Cot

---

### **WARNING**

To avoid accidental release of the Pedi-Mate™, and possible injury to the infant, ensure the buckle is located away from obstructions on the cot.

7. Verify all the straps are snug and fastened securely (see Figure 24).



Figure 24- Pedi-Mate™ Strapped to a **RUGGED** Cot

### **NOTE**

- These are general instructions for installation of the Pedi-Mate™. Safe and proper use of the Pedi-Mate™ is solely at the discretion of the user. Stryker recommends all users be trained on the proper use of the Pedi-Mate™ before using it in an actual situation.
- Retain these instructions for future reference. Include them with the product in the event of transfer to new users.
- Pedi-Mate™ is a trademark of Ferno-Washington Inc.

# Cot Features

---

## OPERATING THE WHEEL LOCK(S)




Figure 25 - **RUGGED** Wheel Lock

The Power Pro TL includes corner to corner wheel locks as a standard feature with an option to specify locks for all four wheels.

1. To secure the wheel lock(s), press down on the base of the locking pedal (A) until it rests firmly against the surface of the wheel.
2. To release the wheel lock(s), press the upper face of the pedal until it returns to its upright position. The upper portion of the pedal will rest against the caster frame when the wheel lock is released.

---

### **WARNING**

- Never apply the wheel lock(s) while a patient is on the cot. Tipping could occur if the cot is moved while a wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
- Never leave a patient unattended on the cot or injury could result. Hold the cot securely while a patient is on the cot.
- Never install or use a wheel lock on a cot with excessively worn wheels. Installing or using a wheel lock on a wheel with less than a 6" diameter could compromise the holding ability of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the cot or other equipment.

---

### **CAUTION**

Wheel lock(s) are only intended to help prevent the cot from rolling while unattended. A wheel lock may not provide sufficient resistance on all surfaces or under loads.

# Cot Features

## ADJUSTING THE WHEEL LOCKING FORCE




1. To adjust the wheel locking force, remove the socket screw from the center of the lock pedal. The wheel lock is initially assembled with the pedal set at the minimum locking force. The marker on the pedal (item A) is aligned with the marker on the octagonal sleeve (item B).
2. Remove the sleeve (B). Rotate the sleeve counterclockwise to increase the pedal locking force and clockwise to decrease the locking force. Insert the sleeve into the pedal. Reinstall the socket screw.
3. Test the pedal locking force and verify the pedal holds properly before returning the cot to service.



# Cot Features

---

## ADJUSTING THE FOOTREST AND KNEE GATCH



**Figure 28 - Gatch Raised**

The footrest is adjustable to allow for elevation of the patient's legs.

**To raise the Footrest**, lift the foot rest frame (A) as high as possible until it locks in place. Be sure the foot rest frame is securely locked before releasing your grip.

**To lower the Footrest**, lift the foot rest frame (A) and, while holding the frame, push up on the red release handle (C). Slowly lower the footrest to the flat position.

**To raise the Knee Gatch**, lift either of the red lifting loops (B) until the knee gatch is in its fully raised position then slowly lower the knee gatch to allow the support bracket to engage in the locking mechanism. Check to be sure the lock is fully engaged before releasing the lifting loop.

**To lower the Knee Gatch**, lift either of the red lifting loops (B) to relieve pressure on the locking mechanism and while holding the loop, push on the red release handle (C) until the bracket disengages. Lower the knee gatch to the flat position.

# Cot Features

---

## OPERATING THE BACKREST

To raise or lower the backrest, squeeze the red release handle (A) and move the backrest into the desired position. Release the handle at any point, and the backrest will lock into position.

**Note:** Lowering the backrest without a patient may require slightly more pressure.

## OPERATING THE SIDERAILS

To raise, swing the siderail up toward the head end of the cot until it locks in place. When a patient is on the cot, always keep the siderails in the raised position unless the patient is being transferred.


To lower, squeeze the red release handle (B) to release the siderail latch. Guide the siderail down toward the foot end until flat.

---

### WARNING

Siderails are not intended to serve as a patient restraint device. Refer to [pages 22](#) and [23](#) for proper restraint strap usage. Failure to utilize the siderails properly could result in patient injury.

---



**Figure 29 - Backrest Elevated and Siderails Raised**

# Cot Features

## PROPER BATTERY MANAGEMENT: DeWALT® BATTERY MANAGEMENT AND CHARGER CARE

The 24 Volt DeWALT® battery pack is designed for power cordless tools throughout the world. Care and maintenance procedures for the DeWALT® battery pack may differ slightly between the power tool industry and the emergency service provider industry. This document is intended to assist with proper battery pack and charger management in order to provide longer life to the rechargeable battery.

### The following facts apply to the DeWALT® 24 Volt Battery system:

All batteries have a life span. This life span is dependent on proper care and maintenance. When a battery is placed in a charger or the charger is turned on-and-off with the battery in it, this is considered a charging cycle. A typical battery has between 800-900 recharging cycles. Keeping a battery in a 220-240V charger on a truck that does not have a continuous power supply to the charger will start a recharge cycle every time the power is switched off and on. For example, vehicles that are connected to a shoreline, unplugging and plugging the truck back in will count as a recharge cycle.

Having a fully charged battery in a charger, while unplugging and plugging in the charger may have damaging effects on the battery that will reduce the battery life. (This also includes inverter and/or generator powered patient compartments that can be switched on and off.)

Moisture and debris in the battery may reduce the battery life. A wet battery should not be recharged until completely dry.

Water on the charger may cause an electrical shock. Keep moisture away from all battery chargers.


### Stryker recommends the following procedures for charging the DeWALT® 24 Volt Battery:

When charging the battery in the vehicle, it is important to provide continuous-uninterrupted power to the battery charger and battery. Connect a 12V charger to a hot lead directly connected to the truck battery system (e.g. cigarette lighter adapter). Do not modify the 12V DeWALT® battery charger in any way. A full charge will take place in approximately one (1) hour, but when the battery is not in use, it is best to leave the battery in the charger for a minimum of (8) hours to take advantage of the “Tune-Up” mode which automatically conditions the battery. The 12V DeWALT® charger draws 2.8 Amps in full charge mode and is protected to not remove power from the truck battery system if the onboard charger fails.

When charging the battery at the station, connect a 110V charger into a station outlet and rotate batteries through the charger. A full charge will take place in approximately one (1) hour, but when the battery is not in use, it is best to leave batteries in the charger to take advantage of the “Tune-Up” mode which automatically conditions the batteries.

If a 110V charger is kept on a truck, do not leave the battery in the charger. Use this option as an “On-Scene” recharge system, not as a battery up-keep system.

**Like any piece of rescue equipment, proper care and maintenance will provide the rescuer the best tools to perform the task at hand. Contact Stryker Technical Support with any additional questions at 1-800-784-4336.**



# Cot Features

---

## PROPER BATTERY MANAGEMENT: DeWALT® BATTERY FREQUENTLY ASKED QUESTIONS

### Is it better for DeWALT® batteries to be completely discharged before charging?

When using the DeWALT® battery to power the POWER PRO ambulance cot, it is beneficial to run the cot until the battery power indicator is flashing red. When the indicator is flashing red, the battery will power the cot through approximately (3) additional service calls. It is recommended that the battery be changed at this time.

### What is Memory, and do DeWALT® batteries have it?

Memory is one of many conditions which causes a loss of runtime. Memory is created from repetitive light use in the exact same application (i.e. Cordless Phones, Video Cameras, Electric Shavers, etc.) The DeWALT® batteries rarely see light use or the exact same loads, due to variability from the user, the accessory size, as well as the product. The same variability which causes different runtimes prevents our cells from developing memory. The POWER PRO is considered a high-drain application. Memory typically develops in lower-drain rate applications, such as cordless phones, laptops, etc... because the rate at which the battery is draining is continuously the same. The POWER PRO will draw a higher current and have sporadic drain rates minimizing the opportunity for the battery to develop a memory.

### Does it hurt DeWALT® batteries to leave them in the charger?

No. The DeWALT® chargers have a maintenance mode which allows batteries to remain in the charger, maintaining a fully charged battery until the user is ready to use it. If the batteries are stored outside of the charger they will discharge naturally, 15-20% in the first 24 hours, 7-10% the next day, and about 1% every day thereafter. NiCd batteries lose the bulk of the capacity when outside of the charger in the first 3 days. In fact, it is better for the battery to leave it in the charger to be sure it goes through equalization and maintenance modes which takes (8) hours minimum.

### What can I do to improve the runtime of my battery?

If no permanent damage has been done to your battery, you may be able to improve its runtime. The correct procedure for charging your batteries is as follows:

1. Discharge the battery under normal use until the battery indicator light is flashing red or you feel a loss of power.
2. Remove the battery.
3. Let the battery sit out of the charger until the battery is at room temperature.
4. Place the battery in the charger for a minimum of 8 hours to allow for a full charge on each individual cell. If there is no difference in runtime, there is either permanent damage or the battery has reached the end of its usable life. In either case, the battery should be replaced.

### Does the outside temperature affect batteries? How?

Yes. If the batteries are too hot or too cold, the batteries will not take a full charge. The optimum **charging** temperature range is between 4°C and 40°C. Charging batteries outside this range may result in a permanent loss of runtime. The maximum **operating** temperature range is between 6°C and 115°C with an optimum range between -10°C and 60°C. Operating batteries outside of this range may result in a permanent loss of runtime. When batteries are being charged and discharged, a chemical reaction is taking place, and if it is too hot or cold the chemical reaction is disturbed causing a loss of runtime.

### Can the DeWALT® charger be used with a generator?

Yes. All DeWALT® chargers, excluding the DW9106, have been designed to handle the variations in voltage and current delivered by generators.

### What should be done with batteries once they have gone bad?

RECYCLE THEM. DeWALT® is an active participant with RBRC (Recharge Battery Recycling Corporation), the organization which is the international leader in the collection, transportation, and recycling of NiCd cells. Old batteries should be disposed of at a DeWALT® Service Center. For more information call 1-800-8-BATTERY OR 1-800-8-228-8379.

DeWALT® is a registered trademark of Black & Decker Inc.

# Cot Features

---

## BATTERY POWER INDICATOR




Figure 30 - Extend and Retract Button

To check the battery power level, depress **lightly** on the minus (–) button (A) to activate the battery power indicator light (B). The battery power indicator is located on the button at the foot end, represented by the red or green light.

- The indicator lights GREEN when the battery has adequate charge.
- The indicator flashes RED when the battery needs to be changed and recharged.

---

### WARNING

- To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
- Do not remove the battery when the ambulance cot is activated.
- Avoid contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.

---


### CAUTION

- Only use the battery and charger as specified.
- The POWER PRO TL ambulance cot is not for use with an AC adapter.
- When charging batteries in an ambulance vehicle, locate the charger either in the forward cab or an enclosed compartment (i.e. cabinet).
- Ensure that the battery is fully charged prior to placing into service. An uncharged or depleted battery may cause poor ambulance cot performance.
- Refer to the DeWALT® manual (Stryker part number 6500-001-206) for battery and charger information.

# Cot Features

---

## HOUR USAGE METER



**Figure 31 - Hour Usage Meter**

The cot has an hour usage meter at the foot-end control enclosure that indicates the amount of time (HHH.H hours) that the cot has been in use under power, and supports scheduling of preventative maintenance activities found on [page 54](#) and [56](#).

---


 **CAUTION**

- A preventative maintenance program should be established for all Stryker EMS equipment. Preventative maintenance may need to be performed more frequently based on the usage level of the product.
- Close attention should be given to safety features including, but not limited to:
  - Hydraulic power mechanism
  - All electrical controls return to off or neutral position when released
- For additional maintenance information, refer to the preventative maintenance section (pages 48-51).

# Cot Features

---

## INSTALLING THE BACKREST STORAGE POUCH (OPTIONAL EQUIPMENT)



**Figure 33 - Backrest Storage Pouch**

1. Install the optional backrest storage pouch using the Velcro® straps. Insert each strap through the corresponding holes in the backrest skin and mount the pouch flat against the backrest.

---


**⚠ CAUTION**

- Do not store items under the ambulance cot mattress. Storing items under the mattress can interfere with the operation of the ambulance cot.
- The weight of the equipment in the pocketed backrest storage pouch (if equipped) must not exceed 1.4 stones (20 pounds), (9 kg).
- Make sure the pouch does not interfere with the operation of the retractable head section.

# Cot Features

---

## USING THE EQUIPMENT HOOK (OPTIONAL EQUIPMENT)



**Figure 34 - Equipment Hook**

The equipment hook is used to hang additional accessories or equipment such as defibrillators or monitors.

---

**⚠ CAUTION**

To avoid damage to the equipment hook, the weight of the accessories or equipment must not exceed 1.4 stones (20 pounds), (9 kg).



# Cot Features

---

## OPERATING THE RIGID PUSH BARS AND STORAGE POUCH (OPTIONAL EQUIPMENT)

The **optional** rigid push bars are valuable tools for manoeuvring the cot at any height position.

### Installing the Rigid Push Bar

1. Insert the bars into each of the four corner sockets and push down until each locking mechanism is fully engaged (refer to Figure 35a).

---

#### **WARNING**

- It is important to ensure that the push bar is securely attached before mounting the cot to prevent injury to operation or patient.
- Do not lift the cot with the push bars. Lifting the cot by the push bars may result in failure of the push bar lock mechanism and may injure the operator or patient.



Figure 35a - Installed Push Bar

### Removing the Rigid Push Bar

1. Press the red release button and lift out of its socket.
2. When not in use, the push bars should be stowed in the storage pouch located on the base frame. (refer to Figure 35b).

---

#### **CAUTION**

When not in use, the push bars should be stowed in the base storage pouch located on the base frame.



Figure 35b - Store Push Bar

### Attaching the Storage Pouch

1. To attach the storage pouch, connect the corresponding buckles to secure the pouch.
2. Check to be sure the pouch, and contents do not interfere with the cot operation before raising or lowering the cot, or attempting transport.

# Cot Features

---

## OPERATING THE 2-STAGE I.V. POLE (OPTIONAL EQUIPMENT)




Figure 36 - 2-Stage I.V. Pole Storage Position




Figure 37 - 2-Stage I.V. Pole

1. Lift and pivot the pole from the storage position and push down until it is locked into receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the telescoping portion in place.
4. Hang the I.V. bags on the I.V. hook (D).

---

**⚠ CAUTION**

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 2.8 stones, (40 pounds), (18 kg).

5. Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
6. Lift up and pivot the pole down into the storage position.

# Cot Features

## OPERATING THE 3-STAGE I.V. POLE (OPTIONAL EQUIPMENT)

1. Lift and pivot the pole from the storage position and push down until it is locked into receptacle (A).
2. To raise the height of the pole, turn the lock actuator (B) counterclockwise and pull up on the bottom telescoping portion (C) of the pole to raise it to the desired height.
3. Turn the lock actuator (B) clockwise to lock the bottom telescoping portion in place.
4. For a higher I.V. pole, pull up on section (D) until the spring clip (E) engages.
5. Hang I.V. bags on the I.V. hook (F).

### CAUTION

To avoid damage to the I.V. pole, the weight of the I.V. bags or equipment must not exceed 2.8 stones, (40 pounds), (18 kg).

6. To lower the I.V. pole, push in on the spring clip (E) and slide section (D) down into section (C). Turn the lock actuator (B) counterclockwise and slide section (C) into the bottom tube.
7. Lift up and pivot the I.V. pole down into the storage position.




Figure 38 - 3-Stage I.V. Pole Storage Position




Figure 39 - 3-Stage I.V. Pole

# Cot Operation

---

## OPERATING GUIDELINES

- Use the cot only as described in this manual.
- Read all labels and instructions on the cot before using the cot.
- Use a minimum of two (2) operators to manipulate the cot while a patient is on the cot. If additional assistance is needed, see the reference chart on page 48.
- Do not adjust, roll or load the cot into a vehicle without advising the patient. Stay with the patient and control the cot at all times.
- The ambulance cot can be transported in any position. Stryker recommends transporting the patient in as low a position as is comfortable for the operators to maneuver the cot.
- Only use the wheel lock(s) during patient transfer or without a patient on the ambulance cot.
- Always use the restraint straps.
- Use properly trained helpers when necessary to control the cot and patient.

---

## WARNING

- Improper usage of the **▲RUGGED▲** POWER PRO TL ambulance cot can cause injury to the patient or operator. Operate the ambulance cot only as described in this manual.
- Entanglement in powered ambulance cot mechanisms can cause serious injury. Operate the ambulance cot only when all persons are clear of the mechanisms.
- Practice changing height positions and loading the ambulance cot until operation of the product is fully understood. Improper use can cause injury.
- Do not allow untrained assistants to assist in the operation of the ambulance cot. Untrained technicians/assistants can cause injury to the patient or themselves.
- Do not ride on the base of the **▲RUGGED▲** POWER PRO TL ambulance cot. Damage to the product could occur, resulting in injury to the patient or operator.
- Transporting the cot sideways can cause the cot to tip, resulting in possible damage to the product and/or injury to the patient or operator. Transporting the cot in a lowered position, head or foot end first, will minimize the potential of a cot tip.
- Grasping the **▲RUGGED▲** POWER PRO TL ambulance cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the ambulance cot.

---

## CAUTION

Before operating the cot, clear any obstacles that may interfere and cause injury to the operator or patient.

## PROPER LIFTING TECHNIQUES

- When lifting the ambulance cot and patient, there are five basic guidelines to help you avoid injury:
- Keep your hands close to your body.
- Keep your back straight.
- Coordinate your movements with your partner and lift with your legs.
- Avoid twisting.
- Always operate the ambulance cot as described in this manual.

# Cot Operation

---

## RAISING AND LOWERING THE COT (UNOCCUPIED)

**To raise the cot**, simply press and hold the PLUS (+) button and the cot will rise to its highest position in about 6 seconds.

**To lower the cot**, press and hold the MINUS (-) button and the cot will move from its highest position to the lowest position in about 6 seconds.

## RAISING AND LOWERING THE COT (OCCUPIED)

---

### WARNING

---

A minimum of two operators are required to raise or lower an occupied cot.

**To raise the cot**, both operators should grasp the cot frame to lend stability as the foot end operator presses the PLUS (+) button on the control module.

**To lower the cot**, maintain a firm grip on the cot frame, as the foot end operator presses the MINUS (-) button on the control module. When the cot reaches the desired height, release the button and the cot will stop.

## TRANSPORTING A PATIENT

---

### WARNING

---

- A minimum of two operators are required to transport a patient on the cot.
- In addition to the head and foot end, larger patients may require additional operators.

The Power Pro TL cot can transport a patient in any height position, however, lower heights lend to improved balance and stability, so Stryker recommends keeping the cot in the lowest position comfortable for the cot operators.

During transport, approach doorways and other low obstacles squarely and lift each set of wheels over the obstacle separately. High obstacles should be avoided. Stryker recommends finding an alternate route or seeking additional help with transport.

# Cot Operation

---

## TRANSFERRING THE PATIENT TO THE COT

1. Roll the cot to the patient.
2. Place the cot beside the patient and raise or lower the cot to the level of the patient.
3. Lower the siderails closest to the patient and open the patient restraint straps.
4. Engage wheel locks to prevent cot from moving during patient transfer.
5. Transfer the patient to the cot using accepted EMS procedures.
6. Use all the restraint straps to secure the patient to the cot (refer to pages 22 and 23 for usage instructions).
7. Raise the siderails and adjust the backrest and leg rest as necessary.
8. Disengage wheel locks for transport.

---

### WARNING

- Always use all restraint straps to secure the patient on the cot. An unrestrained patient may fall from the cot and be injured.
- Never leave a patient unattended on the ambulance cot or injury could result. Hold the ambulance cot securely while a patient is on the product.
- Never apply the optional wheel lock(s) while a patient is on the cot. Tipping could occur if the cot is moved while the wheel lock is applied, resulting in injury to the patient or operator and/or damage to the cot.
- Siderails are not intended to serve as a patient restraint device. Refer to pages 22 and 23 for proper restraint strap usage. Failure to utilize the siderails properly could result in patient injury.
- Hydraulically raising or lowering the cot may temporarily affect electronic patient monitoring equipment. For best results, patient monitoring should be conducted when the cot is idle.

## USING THE TRANSFER FLAT (OPTIONAL EQUIPMENT)

When transferring large patients, use of the  Transfer Flat (Stryker part number 6083-001-200) is recommended.

## AMBULANCE COT MOTION

- Make sure all the restraint straps are securely buckled around the patient (refer to Using Restraint Straps under the Cot Operations section).
- The cot can be in any position for rolling.
- When rolling the cot with a patient on it, position an operator at the foot end and one at the head end of the cot **at all times**.
- During transport, approach door sills and/or other low obstacles squarely and lift each set of wheels over the obstacle separately.

---

### WARNING

- High obstacles such as curbing, steps or rough terrain can cause the ambulance cot to tip, possibly causing injury to the patient or operator.
- Transporting the cot in lower positions can reduce the potential of a cot tip. If possible, obtain additional assistance (refer to page 48) or take an alternate route.



# Cot Operation

## TAIL LIFT LOADING AND UNLOADING

Loading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**.

### WARNING

A lift gate stop which is not properly functioning can result in injury to the patient or operator; ensure that the cot is unable to roll back off the lift before using the lift with a cot and patient. Verify the lift gate stop is maintained and functioning properly and ensure that the cot is secure at all times when on the tail lift.

1. Ensure that the patient is secured at all times while on the cot.
2. Move the cot to its lowest position for greater stability by pressing the retract or minus button.
3. Push the cot forward onto the tail lift, head end first. Ensure that the wheels are in the proper location on the tail lift to allow the safety stop to properly rotate in place and prevent the cot from roll backwards (refer to Figure 40).
4. Before raising the tail lift, verify that there is adequate distance between the cot and rear of ambulance, and that there is nothing hanging down from the cot.
5. As one operator raises the tail lift, the second operator should maintain a firm grip on the cot frame to provide greater stability.
6. Both operators should carefully guide the cot into the patient compartment and into the cot fastener (refer to Figure 41). Cot Fastener not included.



Figure 40 - Load Cot



Figure 41 - Push Cot in Patient Compartment

### CAUTION

- The unit must be in its fully down position in order to properly engage the cot fastener (not included). In order to ensure the cot is in the lowest position press the retract (-) button until the cot stops moving downward. This should be done prior to attempting to engage the cot into the cot fastener (not included).
- **In order to reduce risk of damage to the cot or cot fastener, do not attempt to activate the cot height activation while it is engaged in the cot fastener.**

### To Unload a Cot:

1. Ensure the safety gate stop is in position to prevent the cot from rolling off of the lift.
2. Raise the tail lift to the ambulance floor level and disengage the cot from the cot fastener system.
3. Both operators should firmly grasp the cot frame and roll the cot into the tail lift, being sure that the cot is fully removed from the patient compartment and the cot wheels are in the proper position to allow the lift to lower freely.
4. Lower the tail lift to the ground and check to be sure that it is fully lowered and stopped before disengaging the tail lift safety gate and allowing the cot to be rolled off the tail lift.

# Cot Operation

## RAMP LOADING AND UNLOADING

Loading an occupied cot into a vehicle requires a minimum of **two (2) trained operators**.

1. Ensure that the patient is secured at all times while on the cot.
2. Before pushing up the ramp, lower the cot to its lowest position by pushing the retract (-) button for maximum stability. Also verify that there is nothing hanging down from the cot such as blankets or straps (refer to Figure 42).
3. Using the optional push bars, both operators should push/pull the cot up the ramp, head end first. The operators should take care to guide the cot up the center of the ramp (refer to Figure 43).
4. Both operators should then push the cot into the patient compartment, until the cot engages the cot fastener (not included).



Figure 42 - Load Cot

### CAUTION

- The unit must be in its fully down position in order to properly engage the cot fastener (not included).
- In order to ensure the cot is in the lowest position press the retract (-) button until the cot stops moving downward. This should be done prior to attempting to engage the cot into the cot fastener (not included).
- In order to reduce risk of damage to the cot or cot fastener, do not attempt to activate the cot height activation while it is engaged in the cot fastener.



Figure 43 - Push Cot in Patient Compartment



# Cot Operation

## ENGAGING/RELEASING THE COT INTO THE FASTENER

### Engaging the Cot into the Fastener:

1. Before attempting to engage the cot into the fastener ensure that the pedal is fully depressed (refer to Figure 46).
2. Guide the cot into both the head and foot end portions of the fastener until the locking mechanism is engaged. (refer to Figure 47).
3. Check to be sure the cot is securely fastened at both the head and foot end before releasing your grip on the cot.




Figure 46 - Pedal fully depressed



Figure 47 - Cot in Fastener

### Removing the Cot from the Fastener:

1. Press down firmly on the foot pedal until the locking mechanism disengages (refer to Figure 48).
2. Roll the cot out of the patient compartment (refer to Figure 49).




Figure 48 - Press on foot end




Figure 49 - Remove cot

# Cot Operation

## HIGH SPEED RETRACT/EXTEND

- The ambulance cot is equipped with a high-speed retract mode to expedite lifting the cot over obstacles.
- The undercarriage **rapidly** retracts towards its uppermost position once the weight of the ambulance cot and patient is off the wheels. Press the retract (-) button to actuate the control switch.
- The undercarriage **rapidly** extends towards its lowermost position once the weight of the ambulance cot and patient is off the wheels. Press the extend (+) button to actuate the control switch and extend the base.

## WARNING

- Whenever the weight of the ambulance cot and patient is off the wheels, the ambulance cot will **automatically** enter the high speed retract mode if the retract (-) button is pressed.
- Once the weight is off the ground, the operator(s) must support the load of the patient, ambulance cot and any accessories. Failure to support the load properly may cause injury to the patient or operator.

## CASTER STEER LOCK

The caster steer lock feature locks the head end caster to improve control in steering the cot and unlocks the caster to allow free swiveling for greater mobility.

## CAUTION

Engaging the steer lock and attempting to push the cot sideways may affect the cot stability.

### To engage the caster steer lock:

1. Press the front portion of the steer lock pedal (refer to Figure 50).
2. Push the cot towards the head end.
3. The caster will swivel and align in the forward position and the steer lock will engage.

### To disengage the caster steer lock:

1. Press the rear portion of the steer lock pedal (refer to Figure 51).
2. After the pedal is in the unlocked position, push the cot in any direction.
3. The casters will swivel freely.



Figure 50 - Engage Caster Steer Lock



Figure 51 - Disengage Caster Steer Lock

## NOTE

When engaging the cot into the fastener, it is advised to disengage the caster steer lock feature to allow the head end casters to swivel freely and allow the head end of the cot to more easily align with the fastener.

# Cot Operation

---

## MANUAL RAISING AND LOWERING OF THE COT




Figure 52 - Manual Release Handle

The Power Pro TL is equipped with a manual override function to allow two operator manual raising and lowering of the cot with or without a patient.

The **red** manual release lever (A) is located along the patient left side of the lift bar at the foot end of the cot.

### To raise or lower the cot with the manual release:

Both Operators - lift the cot during the raise/lower operation to support the weight of the cot.

With the ambulance cot supported by an operator at each end, the foot end operator pulls the release handle toward the lift bar as both operators lift the cot slightly to take pressure off of the cot base. While the release handle is pulled, raise or lower the ambulance cot to the desired position and then release the handle to lock ambulance cot into position.

### NOTE

- The operators must lift the cot weight slightly off the wheels to use the manual extend or retract while a patient is on the cot.

---

### CAUTION

Activation of the manual release may cause the ambulance cot to drop slowly if less than 9.2 stones (130 pounds), (59 kg) are on the cot.

# Cot Operation

## ADJUSTING COT HEIGHT

Changing the height of the cot while a patient is on the cot requires a minimum of **two (2) operators**, positioned at each end of the ambulance cot.





**Operator 1** - Grasp the ambulance cot frame at the foot end. Actuate the control switch, depress either the (+) or (-) button depending on desired travel direction, and allow the litter to raise/lower to the desired position.

**Operator 2** - Maintain a firm grip on the outer rail until the ambulance cot is in position.

### WARNING

Grasping the **RUGGED** ambulance cot improperly can cause injury. Keep hands, fingers and feet away from moving parts. To avoid injury, use extreme caution when placing your hands and feet near the base tubes while raising and lowering the ambulance cot.

## USING ADDITIONAL ASSISTANCE

	Changing Levels	Rolling
<p>Two Operators</p> <p>Two Helpers</p>		
<p>Two Operators</p> <p>Four Helpers</p>		

# Cot Operation

## OPERATING THE RETRACTABLE HEAD/FOOT SECTIONS

The head section telescopes from an extended position suitable for loading the ambulance cot into an emergency vehicle to a second position retracted within the litter frame.

### To Extend the Head or Foot Section:

1. Squeeze the red release levers (A) on each side of the frame as you pull the head or foot section out of its retracted position.
2. Check to be sure the head or foot section is securely locked in position.




Figure 61 - Head and Foot Section Extended

### To Retract the Head or Foot Section:

1. Squeeze the red release levers (A) on each side of the frame as you push the head or foot section out of its extended position.
2. Check to be sure the head or foot section is securely locked in position.




Figure 62 - Head and Foot Section Retracted

---

 **CAUTION**

Never operate the cot with the head or foot section in a partially extended and unlocked position.



# Cot Operation

---

## BATTERY OPERATION




Figure 64 - Battery Removal & Replacement

The ambulance cot is supplied with two removable 24 Volt DEWALT® batteries as the power source. To insert the battery, align it with the battery mounting bracket slot and push the battery into the enclosure until it locks into place.

To remove the battery, locate the **RED** battery release (A) along the patient right side of the foot end control enclosure. Push the battery release button to release the latch. Slide the released battery out of the enclosure to the right (not shown).

To reinstall the battery, align the tabs in the battery enclosure and push the battery into the enclosure until the latch clicks into place. Once the battery is installed, a green indicator light located adjacent to the foot end controls will confirm the battery has adequate charge. If the indicator flashes red, the battery needs to be changed and recharged.

### NOTE

Keep your spare battery on the charger at all times. Batteries slowly lose power when not on the charger.

---

### WARNING

- To avoid risk of electric shock, never attempt to open the battery pack for any reason. If the battery pack case is cracked or damaged, do not insert it into the charger. Return damaged battery packs to a service center for recycling.
- Do not remove the battery when the ambulance cot is activated.
- Avoid contact with a wet battery or battery enclosure. Contact may cause injury to the patient or operator.

---

### CAUTION

Remove the battery if the cot is not going to be used for an extended period of time (over 24 hours).


# Hydraulic Sub-Assembly (6550-001-030) Access Instructions

## Tools Required:

- T27 Torx Wrench

## Procedure:

1. Using a T27 Torx wrench, remove the six button head cap screws (p/n 0004-592-000) from the six well nuts (p/n 0055-100-074) that hold the mid-section (p/n 6550-001-111) to the motor mounts (p/n 6500-001-294) and (p/n 6500-001-194).
2. Remove the mid-section skin (p/n 6550-001-111) and lay to the side out of the way along with the six button head cap screws (p/n 0004-592-000). The six well nuts will stay with the two motor mounts (p/n 6500-001-294) and (p/n 6500-001-194).
3. Remove the four button head cap screws (p/n 0004-596-000) and the four nylock hex nuts (p/n 0016-102-000) from the two straight “T” pivots (p/n 6100-003-125) that hold the gatch assembly (p/n 6550-001-019) to the litter cross brace (p/n 6500-001-196).
4. Flip the gatch assembly (p/n 6550-001-019) toward the foot end of the cot until it rests on the telescoping foot end (p/n 6550-001-015).
5. Reverse the above procedures to reassemble.



# Cleaning

---

The **RUGGED** POWER PRO TL ambulance cot is designed to be power washable. The unit may show some signs of oxidation or discoloration from continuous washing, however, no degradation of the cot's performance characteristics or functionality will occur due to power washing as long as the proper procedures are followed.

Thoroughly clean the cot once a month. Clean Velcro AFTER EACH USE. Saturate Velcro with disinfectant and allow disinfectant to evaporate. (Appropriate disinfectant for nylon Velcro should be determined by the service.)

## WASHING PROCEDURE

- **Remove the battery.** The battery and charger are **not** immersible or power washable.
- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- The preferred method Stryker Medical recommends for power washing the POWER PRO TL ambulance cot is with the standard hospital surgical cart washer or hand held wand unit.

## WASHING LIMITATIONS

---

### WARNING

Use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the cot.

---

### CAUTION

- DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT
- Maximum water temperature should not exceed 82°C /180°F.
- Maximum air dry temperature (cart washers) is 115°C /240°F.
- Maximum water pressure should not exceed 130.5 bar/1500 psi/ If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 61 cm (24 inches) from the unit.
- Towel dry all casters and interface points.
- Failure to comply with these instructions may invalidate any/all warranties.
- Remove the battery before washing the cot.



# Cleaning

---

In general, when used in those concentrations recommended by the manufacturer, either phenolic type or quaternary type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may result.

Suggested cleaners for the **RUGGED** POWER PRO TL cot surfaces:

- Quaternary Cleaners (active ingredient - ammonium chloride)
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)

Avoid oversaturation and ensure the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

---

## WARNING

- SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products described above are used to clean Stryker patient care equipment, measures must be taken to insure the cots are wiped with a cloth soaked in clean water and thoroughly dried following cleaning.
- Failure to properly rinse and dry the cots will leave a corrosive residue on the surface of the cots, possibly causing premature corrosion of critical components.

## NOTE

Failure to follow the above directions when using these types of cleaners may void this product's warranty.

## REMOVAL OF IODINE COMPOUNDS

Use a solution of 1/2 tablespoons Sodium Thiosulfate in a pint of warm water to clean the stained area. Clean as soon as possible after staining occurs. If stains are not immediately removed, allow solution to soak or stand on the surface. Rinse surfaces which have been exposed to the solution in clear water before returning unit to service.

---

## WARNING

Failure to properly clean or dispose of contaminated mattress or cot components will increase the risk of exposure to blood borne pathogens and may cause injury to the patient or the operator.

# Preventative Maintenance

---

The **RUGGED** POWER PRO TL ambulance cot requires regular maintenance. Establish and follow a maintenance schedule and keep records of maintenance activity (refer to page 57 for a form).

---

## WARNING

- Do not modify the **RUGGED** POWER PRO TL ambulance cot or any components of the cot, including the hydraulic unit. Modifying the product can cause unpredictable operation resulting in injury to the patient or operator. Modifying the product will also void its warranty.
- Escaping fluid under pressure can penetrate the skin causing serious injury. Avoid the hazard by relieving pressure before disconnecting hydraulic or other lines. Tighten all connections before applying pressure.
- If an accident occurs, see a doctor immediately. Any fluid injected into the skin must be surgically removed within a few hours or gangrene may result. Doctors unfamiliar with this type of injury should reference a knowledgeable medical source.

When using maintenance products, follow the directions of the manufacturer and reference all material safety data sheets.

---

## CAUTION

- Improper maintenance can cause injury or damage to the product. Maintain the ambulance cot as described in this manual. Use only Stryker approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty.
- Failure to use authorized parts, lubricants, etc. could cause damage to the ambulance cot and will void the warranty of the product.
- Hydraulic lines, hoses, and connections can fail or loosen due to physical damage, kinks, age, and exposure. Check hoses and lines regularly to avoid damage to the cot. Check and tighten loose connections.
- Do not tip the ambulance cot on end and actuate the product as this will allow air to enter the hydraulic system.

# Preventative Maintenance

## REGULAR INSPECTION AND ADJUSTMENTS

### Maintenance Intervals

The following schedule is intended as a general guide to maintenance. Bear in mind that such factors as weather, terrain, geographical location, and individual usage will alter the required maintenance schedule. If you are unsure as to how to perform these checks please contact your Stryker service technician. If you are in doubt as to what intervals to follow in maintaining your product, consult your Stryker service technician.

Item	Routine	Every (whichever comes first)			
		1 Month or 2 hours	3 Months or 6 hours	6 Months or 12 hours	12 Months or 24 hours
Cylinder	Verify the ambulance cot and fastener fit and function properly.				X
	All fasteners are secure (reference all assembly drawings). Verify the cylinder is adjusted so the lock nut is tight and the cot stops moving when it hits the dead stops.		X		
Hydraulics	Inspect for and verify that there are no hydraulic fluid (red) leaks; inspect the fittings and tighten as necessary. Wipe away any residual fluids.		X		
	Inspect motor mount and verify that all fasteners are secure.		X		
	Verify that there are no hydraulic fluid leaks. Wipe any fluid drops.		X		
	Inspect the reservoir and verify that there are no leaks.		X		
	Inspect hoses and fittings for damage or wear; replace as necessary.			X	
	Verify the hydraulic velocity fuse - Place a weight of approximately 3.6 stones (50 pounds), (22.7 kg) on the cot, raise the cot, lift the cot with two operators, pull the manual release, rapidly set the cot down, verify that the cot does not drop.				X
Electronic Controls	Verify high speed retract is working.			X	
Switches	Verify there is no damage or wear to switch.			X	
	Verify switch operates correctly - Replace if necessary.			X	
Cables/Wires	Verify there is no damage or pinching of wiring harness, cables or lines.		X		
	Check routing(s) and connection(s), verify there are no hanging wires.	X			
	Verify there are no damaged connectors, replace if necessary.		X		
Manual Release	Verify that the manual release functions properly.	X			
	Verify the release handle returns to the stowed position.				X
	Verify the base extends/retracts smoothly when the manual release is engaged.		X		
	With 7.1 stones (100 lbs) or more on the cot, verify the cot does not lower when the manual release is pulled.		X		

# Preventative Maintenance


Item	Routine	Every (whichever comes first)			
		1 Month or 2 hours	3 Months or 6 hours	6 Months or 12 hours	12 Months or 24 hours
Litter	Inspect the ambulance cot frame/litter.	X			
	Verify all welds intact, not cracked or broken.				X
	Verify no bent, broken or damaged components.			X	
	Verify all fasteners secure (reference all assembly drawings).		X		
	Verify warning labels present, legible (reference assembly drawings).				X
	Verify the siderails operate and latch properly.		X		
	Verify the backrest cylinder operates properly.		X		
Mattress	Adjust pneumatic cylinder for full range of motion, if required.		X		
	Verify the footrest operates properly.			X	
Restraints	Verify no cracks or tears on cot mattress.			X	
	Inspect patient restraints for proper function and no excessive wear. Check for bent or broken receiver or latch plate, torn or frayed webbing. <b>Note:</b> Restraints showing wear or not operating properly must be replaced immediately.	X			
Base	Inspect the ambulance cot frame/base.	X			
	Verify all welds intact, not cracked or broken.				X
	Verify no bent, broken, or damaged components.			X	
	Verify all fasteners secure.		X		
Wheels	Verify wheels are free of debris.			X	
	Verify tyres in good condition.				X
	Verify all wheels secure, rolling and swiveling properly.	X			
	Check and adjust wheel locks as necessary.				X
X-Frame	Verify smooth operation of X-frame.		X		
Head Section	Verify all fasteners secure.		X		
	Verify no bent, broken, or damaged components.			X	
	Verify the head and foot sections extend and lock properly.		X		
	Verify the grip bar has no excessive damage.			X	
Accessories	Verify the IV pole (optional) operates properly.		X		
	Verify the restraint extender (optional) operates properly.		X		
	Verify the push bars (optional) lock and unlock properly.		X		






# Cot Assembly

For Reference Only: 6550-001-010



# Cot Assembly


For Reference Only: 6550-001-010






# Cot Assembly

For Reference Only: 6550-001-010




# Cot Assembly

For Reference Only: 6550-001-010





# Cot Assembly

For Reference Only: 6550-001-010




# Cot Assembly

For Reference Only: 6550-001-010



# Cot Assembly




For Reference Only: 6550-001-010

# Cot Assembly



For Reference Only: 6550-001-010

## CABLE ROUTING





# Cot Assembly

For Reference Only: 6550-001-010



# Cot Assembly

For Reference Only: 6550-001-010





# Cot Assembly

## Cot Assembly - Components - 6550-001-010 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-585-000	Button Head Cap Screw	10
B	0004-589-000	Button Head Cap Screw	12
C	0004-592-000	Button Head Cap Screw	6
D	0004-593-000	Button Head Cap Screw	4
E	0004-596-000	Button Head Cap Screw	8
F	0004-850-000	Button Head Cap Screw	2
G	0004-852-000	Socket Head Cap Screw	8
H	0011-004-000	Washer	1
J	0011-013-000	Washer	1
K	0014-040-000	Washer	2
L	0016-028-000	Nylock Hex Nut	2
M	0016-035-000	Nylock Hex Nut	1
N	0016-049-000	Nylock Hex Nut	2
P	0016-102-000	Nylock Hex Nut	8
R	0025-079-000	Rivet	25
T	0025-133-000	Rivet	9
U	0028-181-000	Truarc Ring	1
V	0038-574-000	Crest-to-Crest Spring	2
W	0946-001-155	Bumper	5
Y	1010-031-077	Gas Spring	1
AA	2030-009-901	WEEE Label	1
AC	6060-090-004	'RUGGED', Small Label	1
AD	6080-090-101	Brake Warning Label	5
AE	6082-026-010	Siderail Assembly	2
AF	6082-090-043	Stryker Label	2
AG	6100-003-125	Straight "T" Pivot	4
AH	6500-001-111	Mid Section Skin	1
AJ	6500-001-118	Siderail Nut	6
AK	6500-001-123	Hall Effects Slider	2
AL	6500-001-128	Plastic Extrusion Spacer	2
AM	6500-001-157	Flange Bearing	2
AN	6500-001-168	Rod Attachment Pin	1
AP	6500-001-170	Fowler Cylinder Pin	1
AR	6500-001-171	Cylinder Mount Cross Tube	1
AT	6500-001-191	Fowler Cylinder Spacer	2
AU	6500-001-195	Motor Mount Casting	2
AV	6500-001-196	Litter Cross Brace	1
AW	6500-001-225	"D" Washer	2
AY	6500-001-258	DeWalt Label	1
BA	6550-001-012	Base Assembly	1
BB	6550-001-015	Telescoping Foot End	1
BC	6550-001-018	Fowler Assembly	1
BD	6550-001-019	Gatch Assembly	1
BE	6550-001-020	Telescoping Head Section	1
BF	6550-001-022	Hitch, FE	1
BG	6550-001-031	Mounted Hydraulics	1
BH	6550-001-032	Outer Rail, Patient Right	1
BJ	6550-001-033	Outer Rail, Patient Left	1
BK	6550-001-090	Litter Hitch Interface, HE	1
BL	6550-001-091	Litter Cross Tube, FE	1

## Cot Assembly - Components - 6550-001-010 (reference only)


# Cot Assembly

---

<b>Item</b>	<b>Part No.</b>	<b>Part Name</b>	<b>Qty.</b>
BM	6550-001-099	Fowler Release Handle	1
BN	6550-001-102	Backrest Skin	1
BP	6550-001-170	Outer Rail Bumper	2
BR	6550-001-197	Velcro Strap	2
BT	6550-001-233	Power Pro TL Label	2
BU	6550-001-234	Power Pro TL Specification Label	1
BV	6550-001-259	Weight Capacity Label	2
BW	0059-211-000	11" Nylon Cable Tie	3
BY	0004-634-000	Button Head Cap Screw	1
CA	6550-001-201	O <sup>2</sup> Caution Label	1



# Base Assembly

For Reference Only: 6550-001-012





# Base Assembly

For Reference Only: 6550-001-012




# Base Assembly

For Reference Only: 6550-001-012




# Base Assembly

For Reference Only: 6550-001-012



# Base Assembly

For Reference Only: 6550-001-012







# Base Assembly

## Base Assembly - Components - 6550-001-012 (reference only)

Item	Part No.	Part Name	Qty.
A	0003-205-000	Hex Head Cap Screw	4
B	0003-388-000	Hex Head Cap Screw	2
C	0004-319-000	Flat Head/ Hex Socket Bolt	4
D	0004-438-000	Socket Head Cap Screw	2
E	0004-376-000	Button Head Cap Screw	2
F	0004-587-000	Button Head Cap Screw	12
G	0004-634-000	Button Head Cap Screw	1
H	0004-848-000	Button Head Cap Screw	6
J	0025-142-000	Large Flange Open End Pop Rivet	2
K	0011-062-000	Washer	4
L	0014-002-000	Washer	4
M	0014-040-000	Washer	2
N	0014-115-000	Washer	2
P	0015-051-000	Square Nut	12
R	0016-002-000	Fiberlock Hex Nut	4
T	0016-049-000	Nylock Hex Nut	2
U	0016-060-000	Toplock Hex Jam Nut	4
V	0016-102-000	Nylock Hex Nut	2
W	0021-180-000	Set Screw	2
Y	0025-133-000	Rivet	2
AA	0037-083-000	Tube Plug	4
AB	0037-243-000	Tube Plug	2
AC	0037-244-000	Hole Plug	2
AD	0038-382-000	Compression Spring	1
AE	0081-244-000	Flange Bearing	4
AF	6060-002-010	6" Molded Wheel Assembly	4
AG	6090-001-009	Caster Nut	4
AH	6100-003-022	Brake and Caster Assembly	1
AJ	6100-003-136	Steer Lock Pin	1
AK	6100-203-135	Pin Sleeve-Steer Lock	1
AL	6100-203-136	Steer Lock Plastic Bushing	1
AM	6500-001-056	Frame, Base, Tube, Inner	1
AN	6500-001-090	Head End Cross Tube	1
AP	6500-001-145	Caster Mount Spacer	4
AR	6500-001-162	Flange Bearing	1
AT	6500-001-165	Cylinder Mount, Bottom Pivot	14
AU	6500-001-166	Flange Bearing	14
AV	6500-001-178	Spacer, Base, Outer	2
AW	6500-001-182	Cross Tube, Stiffener Bar	1
AY	6500-001-183	Spacer, Base, Outer	2
BA	6500-001-225	"D" Washer	2
BB	6500-001-226	Bearing, Pivot, Base, Tube	4
BC	6500-001-227	Post, Pivot, Base, Tube	2
BD	6500-001-228	Inner Lift Tube Sleeve	4
BE	6500-001-229	Base Tube, Foot	2
BF	6550-001-203	Spacer, Solid, FE	2
BG	6500-001-308	Base Strap, Patient Right	1
BH	6500-001-309	Base Strap, Patient Left	1
BJ	6500-001-310	Base Strap Clamp	2

# Base Assembly


---

## Base Assembly - Components - 6550-001-012 (reference only)

Item	Part No.	Part Name	Qty.
BK	6500-001-341	Post, Pivot, Base, Tube	2
BL	6550-001-013	Outer Base Tube Assembly	2
BM	6550-001-016	Wheel Lock Caster, Steer Lock	1
BN	6550-001-021	Hitch, HE	1
BP	6550-001-023	<a href="#">Inner Lift Tube Assembly (Page 103)</a>	2
BR	6550-001-028	<a href="#">Outer Lift Tube Assembly (Page 104)</a>	2
BT	6550-001-034	Inner Lift Tube Ass'y, Litter Pivot, PR	1
BU	6550-001-035	Inner Lift Tube Ass'y, Litter Pivot, PL	1
BV	6550-001-085	Spacer, Base, Outer	1
BW	6550-001-103	Steer Lock Pedal	1
BY	6550-001-104	Steer Lock Cross Piece	1
CA	6550-001-105	Steer Lock Spring	1
CB	6550-001-107	Steer Lock Pivot	1
CC	6550-001-204	Spacer, Solid, HE	2
CD	6550-001-171	Rectangular Cross Tube, FE	1
CE	6550-001-179	Steer Lock Cover	1
CF	6550-001-200	Support Link	2


# Telescoping Foot End

For Reference Only: 6550-001-015



# Telescoping Foot End

For Reference Only: 6550-001-015




**Telescoping Foot End Assembly - Components - 6550-001-015 (reference only)**

Item	Part No.	Part Name	Qty.
A	0004-848-000	Button Head Cap Screw	2
B	0004-849-000	Socket Head Cap Screw	2
C	0011-065-000	Washer	2
D	0025-079-000	Rivet	2
E	0025-126-000	Rivet	2
F	0025-185-000	Rivet	8
G	0055-100-075	Riv Nut	2
H	0721-032-049	Bumper Roller	2
J	6500-001-026	Head Section Lock Assembly	2
K	6500-001-087	Cap Bearing	2
L	6500-001-088	Internal Bearing	4
M	6550-001-025	Frame, FE, Bonded	1
N	6550-001-156	Section Release Trigger, HE/FE	2
P	6550-001-162	Link Release, FE	2
R	6550-001-163	Link Guide	2
T	6550-001-169	Bumper Pivot	2


# Telescoping Head Section

For Reference Only: 6550-001-020



# Telescoping Head Section

For Reference Only: 6550-001-020




Telescoping Head Section Assembly - Components - 6550-001-020 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-848-000	Button Head Cap Screw	2
B	0004-849-000	Socket Head Cap Screw	2
C	0011-065-000	Washer	2
D	0025-079-000	Rivet	2
E	0025-126-000	Rivet	2
F	0025-185-000	Rivet	8
G	0055-100-075	Riv Nut	2
H	0721-032-049	Bumper Roller	2
J	6500-001-026	Head Section Lock Assembly	2
K	6500-001-087	Cap Bearing	2
L	6500-001-088	Internal Bearing	4
M	6550-001-024	Frame, HE, Bonded	1
N	6550-001-156	Section Release Trigger, HE/FE	2
P	6550-001-161	Link Release, HE	2
R	6550-001-163	Link Guide	2
T	6550-001-169	Bumper Pivot	2

# Fowler Assembly

For Reference Only: 6550-001-018




## Fowler Assembly - Components - 6550-001-018 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-597-000	Button Head Cap Screw	1
B	0015-050-000	Hex Nut	1
C	0016-028-000	Fiberlock Hex Nut	1
D	0021-119-000	Set Screw	1
E	0021-138-000	Set Screw	1
F	0028-076-000	External Retaining Ring	2
G	0946-035-025	Liner	2
H	6060-032-038	Gas Spring Yoke	1
J	6060-032-040	Fowler Lift Pivot - handle - semi	1
K	6550-001-051	Fowler Weldment	1

# Gatch Assembly


For Reference Only: 6550-001-019






# Gatch Assembly

For Reference Only: 6550-001-019





# Gatch Assembly

For Reference Only: 6550-001-019



# Gatch Assembly

For Reference Only: 6550-001-019



# Gatch Assembly


---

## Gatch Assembly - Components - 6550-001-019 (reference only)

Item	Part No.	Part Name	Qty.
A	6550-001-186	Gatch Pivot Pin	2
B	0004-853-000	Button Head Cap Screw	4
C	0011-448-000	Plain Washer	4
D	0014-020-000	Washer	2
E	0016-028-000	Fiberlock Hex Nut	2
F	0023-162-000	Screw	8
G	0025-079-000	Rivet	34
H	0025-132-000	Rivet	2
J	0081-255-000	Split Bearing	2
K	6100-031-096	Trend Lift Strap	2
L	6100-031-108	Pivot, Knee, Litter	4
N	6550-001-017	Gatch Telescoping Assembly	1
P	6550-001-053	Gatch Pivot Weldment	2
R	6550-001-057	Gatch Lock Tube Weldment	1
T	6550-001-110	Foot Section Skin	1
U	6550-001-111	Thigh Section Tube	2
V	6550-001-112	Thigh Section Skin	1
W	6550-001-116	U-Tube, Telescoping Foot Section	1
Y	6550-001-124	Front Gatch Release	1
AA	6550-001-125	Back Gatch Release	1
AB	6550-001-126	Lever Front Gatch Release	1
AC	6550-001-129	Cross Tube, Pivot, Gatch	1
AD	6550-001-131	End Cap, Bearing, Gatch	1
AE	6550-001-193	Gatch Handle Guard, PL	1
AF	6550-001-196	Gatch Handle Guard, PR	1

# Telescoping Gatch Assembly

For Reference Only: 6550-001-017




## Telescoping Gatch Assembly - Components - 6550-001-017 (reference only)

Item	Part Name	Qty.
A	Rivet	1
B	Rivet	1
C	Semi-Tubular Rivet	4
D	Head Section Lock Assembly	1
E	Internal Bearing	2
F	Gatch Link	1
G	Inner Gatch Tube	1


# Hitch Assembly, Foot End

For Reference Only: 6550-001-022




# Hitch Assembly, Foot End

For Reference Only: 6550-001-022



# Hitch Assembly, Foot End

For Reference Only: 6550-001-022





# Hitch Assembly, Foot End


---

## Hitch Assembly, Foot End - Components - 6550-001-022 (reference only)

Item	Part No.	Part Name	Qty.
A	0003-391-000	Hex Head Cap Screw	2
B	0023-162-000	Screw	15
C	0025-133-000	Rivet	2
D	0038-572-000	Compression Spring	2
E	0946-001-155	Bumper	2
F	6500-001-138	Battery Release Button	1
G	6500-001-139	Battery Release Lock	1
H	6500-001-216	Terminal Block, Reterminated	1
J	6550-001-014	Electronics Assembly	1
K	6550-001-029	Switch Assembly	1
L	6550-001-060	Tube Weldment, FE Hitch	1
M	6550-001-092	Board Enclosure	1
N	6550-001-141	Sleeve, Bar, FE Hitch	1
P	6550-001-146	Bar, FE Hitch	1
R	6550-001-167	Wear Guard, FE Hitch	2
T	6550-001-172	Cable Assembly	1
U	6550-001-173	Cover, FE Housing	1
V	6550-001-178	Bottom Plate, Foot End Enclosure	1


# Hitch Assembly, Head End

For Reference Only: 6550-001-021



# Hitch Assembly, Head End

For Reference Only: 6550-001-021



# Hitch Assembly, Head End


---

## Hitch Assembly, Head End - Components - 6550-001-021 (reference only)


Item	Part No.	Part Name	Qty.
A	0003-391-000	Hex Head Cap Screw	2
B	0004-376-000	Button Head Cap Screw	2
C	0004-634-000	Button Head Cap Screw	8
D	0016-028-000	Fiberlock Hex Nut	2
E	0025-133-000	Rivet	2
F	0026-385-000	Clevis Pin	1
G	0038-588-000	Extension Spring	1
H	0055-100-075	Riv Nut	4
J	0056-016-000	Bumper	2
K	0081-012-002	Flange Bearing	2
L	6550-001-054	Crash Connect Center Weldment	1
M	6550-001-055	Case Weldment, HE Hitch	1
N	6550-001-059	Crash Connect Center Wldmnt, PL	1
P	6550-001-063	Crash Connect Center Wldmnt, PR	1
R	6550-001-083	Label, Pinch Point, HE Hitch	1
T	6550-001-087	Wear Strip, HE Hitch	1
U	6550-001-134	Bar, HE Hitch	1
V	6550-001-140	Sleeve, Bar, HE Hitch	1
W	6550-001-174	Cover, HE Hitch	2
Y	6550-001-175	Guide, HE Hitch	2

# Mounted Hydraulics Assembly

For Reference Only: 6550-001-031



# Mounted Hydraulics Assembly




## Mounted Hydraulics Assembly - Components - 6550-001-031 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-577-000	Button Head Cap Screw	2
B	0004-589-000	Button Head Cap Screw	6
C	0008-030-000	Socket Head Shoulder Bolt	2
D	0015-052-000	Hex Jam Nut	1
E	0016-002-000	Fiberlock Hex Nut	2
F	0016-102-000	Nylock Hex Nut	2
G	0055-100-074	Well Nut	6
H	6500-001-105	Cross Tube, Litter Support	1
J	6500-001-164	Pivot, Cylinder Mount, Top	1
K	6500-001-169	Rod End, Cylinder	1
L	6500-001-194	Motor Mount	1
M	6500-001-212	Cross Bar, Motor Mount	1
N	6500-001-249	Plastic Extrusion, Spacer	1
P	6500-001-250	Plastic Extrusion, Spacer	2
R	6500-001-294	Motor Mount	1
T	6550-001-030	Hydraulics Assembly	1

# Hydraulics Sub-Assembly

For Reference Only: 6550-001-030



# Hydraulics Sub-Assembly

---


## Hydraulics Sub-Assembly - Components - 6550-001-030 (reference only)

Item	Part No.	Part Name	Qty.
A	6500-001-214	Manifold Assembly	1
B	6550-001-213	Cylinder	1
C	6500-001-211	Hydraulic Hose, Rod Side,	1
D	6500-001-210	Hydraulic Hose, Cap Side	1
E	6500-001-282	Hose Fitting, Rod Side, Manifold	1
F	6500-001-283	Hose Fitting, Cap Side, Manifold	1
G	6500-001-296	Hose Fitting, Cap Side, Cylinder	1
H	6500-001-297	Hose Fitting, Rod Side, Cylinder	1
J	6500-001-290	Pressure Switch	1
AJ	6500-001-286	A Valve	1
AM	6500-001-295	Motor	1
AN	6500-001-289	Manual Valve, Non-Locking	1
AP	6500-001-284	Solenoid, A Valve	1
AR	6500-001-285	Solenoid, B Valve	1
AV	6500-001-288	Manual Valve, Locking	1
BA	6500-001-287	B Valve	1



# Outer Rail, Patient Right

For Reference Only: 6550-001-032




**Outer Rail Assembly, Patient Right - Components - 6550-001-032 (reference only)**

Item	Part No.	Part Name	Qty.
A	0004-589-000	Button Head Cap Screw	1
B	0004-591-000	Socket Head Cap Screw	12
C	0004-6120-000	Button Head Cap Screw	2
D	0004-613-000	Socket Head Cap Screw	2
E	0016-102-000	Nylock Hex Nut	2
F	0055-100-076	Riv Nut	5
G	6550-001-029	Sensor Housing, Empty	1
H	6550-001-102	Bracket, Base/Litter Interface	2
J	6550-001-104	Support Bracket, Litter, Outside	2
K	6550-001-106	Support Bracket, Litter, Inside	2
L	6550-001-116	Siderail Bracket	3
M	6550-001-117	Siderail Clamp	2
N	6550-001-064	Outer Rail, PR	1
P	6550-001-096	Clamp, Outer Rail	1
R	6550-001-166	Dead Stop	1

# Outer Rail, Patient Left

For Reference Only: 6550-001-033




**Outer Rail Assembly, Patient Left - Components - 6550-001-032 (reference only)**

Item	Part No.	Part Name	Qty.
A	0004-589-000	Button Head Cap Screw	1
B	0004-591-000	Socket Head Cap Screw	12
C	0004-6120-000	Button Head Cap Screw	2
D	0004-613-000	Socket Head Cap Screw	2
E	0016-102-000	Nylock Hex Nut	2
F	0055-100-076	Riv Nut	5
G	6500-001-029	Sensor Housing, Empty	1
H	6500-001-102	Bracket, Base/Litter Interface	2
J	6500-001-104	Support Bracket, Litter, Outside	2
K	6500-001-106	Support Bracket, Litter, Inside	2
L	6500-001-116	Siderail Bracket	3
M	6500-001-117	Siderail Clamp	2
N	6550-001-065	Outer Rail, PL	1
P	6550-001-096	Clamp, Outer Rail	1
R	6550-001-166	Dead Stop	1
T	0055-100-075	Riv Nut	1

# Inner Lift Tube Assembly

For Reference Only: 6550-001-023




## Inner Lift Tube Assembly - Components - 6550-001-023 (reference only)

Item	Part No.	Part Name	Qty.
A	0014-115-000	Washer	4
B	0025-079-000	Rivet	4
C	6500-001-186	Upper Frame Tube Bearing	1
D	6500-001-187	Lower Frame Tube Bearing	1
E	6550-001-066	Inner Lift Tube Wldmnt, Base Pivot	1

# Outer Lift Tube Assembly

For Reference Only: 6550-001-028




## Outer Lift Tube Assembly - Components - 6550-001-028 (reference only)

Item	Part No.	Part Name	Qty.
A	0014-115-000	Washer	4
B	0025-079-000	Rivet	4
C	6500-001-186	Upper Frame Tube Bearing	1
D	6500-001-187	Lower Frame Tube Bearing	1
E	6500-001-311	Slide Bearing	1
F	6550-001-056	Outer Lift Tube Wldmnt, Base Pivot	1

# Inner Lift Tube Assembly, Litter Pivot, Patient Right

For Reference Only: 6550-001-034




## Inner Lift Tube Assembly, Litter Pivot, PR - Components - 6550-001-034 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-634-000	Button Head Cap Screw	4
B	0055-100-075	Riv Nut	4
C	6500-001-125	Dead Stop, Base	1
D	6550-001-058	Inner Lift Tube Wldmnt, Litter Pivot	1

# Inner Lift Tube Assembly, Litter Pivot, Patient Left

For Reference Only: 6550-001-035




## Inner Lift Tube Assembly, Litter Pivot, PL - Components - 6550-001-035 (reference only)


Item	Part No.	Part Name	Qty.
A	0004-634-000	Button Head Cap Screw	4
B	0055-100-075	Riv Nut	4
C	6500-001-125	Dead Stop, Base	1
D	6550-001-058	Inner Lift Tube Wldmnt, Litter Pivot	1

# Caster, Steer Lock, Wheel Lock Assembly

For Reference Only: 6550-001-016



**Note:** Item H can be rotated incrementally in the counter clockwise direction to increase locking force.




## Caster, Steer Lock, Wheel Lock Assembly - Components - 6550-001-016 (reference only)

Item	Part No.	Part Name	Qty.
A	0004-098-000	Button Head Cap Screw	1
B	0011-454-000	Washer	2
C	0011-456-000	Washer	1
D	0016-118-000	Hex Lock Nut	1
E	0025-079-000	Rivet	2
F	6080-100-032	Test Model	1
G	6080-200-030	Adjustable Caster Lock Pedal	1
H	6080-200-041	Octagonal Sleeve - Adj. Caster Lock	1
J	6550-001-050	Caster Horn, FE	1

# Corner Handle Assembly

For Reference Only: 6550-001-026




## Corner Handle Assembly - Components - 6550-001-026 (reference only)

Item	Part No.	Part Name	Qty.
A	0026-387-000	Slotted Spring Pin	1
B	0038-589-000	Compression Spring	1
C	6510-001-119	Handle Ball	1
D	6550-001-067	Handle Weldment	1
E	6550-001-100	Push Bar Lock Button	1



# 2-Stage I.V. Assembly, Patient Right - 6550-210-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 2-Stage I.V. Pole Assembly, PR - Components - 6550-210-000


Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-041	<a href="#">I.V. Pole Assembly, 2 Stage, PR (pg. 110)</a>	1
C	6550-001-117	I.V. Pole Spacer Plate	1
D	6550-001-118	I.V. Clip Spacer Plate	1

## 2-Stage I.V. Pole Assembly, Patient Right - 6500-001-041



Item	Part No.	Part Name	Qty.
A	0025-133-000	Rivet	2
B	6070-090-105	Caution Label	1
C	6070-210-040	Pole Assembly	1
D	6070-210-045	Sleeve	1
E	6070-210-046	Pivot	1
F	6070-210-049	Ring	1
G	6100-115-050	Socket Weldment	1
H	6500-001-253	Label	1

# 2-Stage I.V. Assembly, Patient Left - 6550-211-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 2-Stage I.V. Pole Assembly, PL - Components - 6550-211-000


Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-042	<a href="#">I.V. Pole Assembly, 2 Stage, PR (pg. 112)</a>	1
C	6550-001-117	I.V. Pole Spacer Plate	1
D	6550-001-118	I.V. Clip Spacer Plate	1

# 2-Stage I.V. Pole Assembly, Patient Left - 6500-001-042



Item	Part No.	Part Name	Qty
A	0025-133-000	Rivet	2
B	6070-090-105	Caution Label	1
C	6070-210-040	Pole Assembly	1
D	6070-210-045	Sleeve	1
E	6070-210-046	Pivot	1
F	6070-210-049	Ring	1
G	6100-115-050	Socket Weldment	1
H	6500-001-254	Label	1

# 2-Stage I.V. Assembly, Dual - 6550-212-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 2-Stage I.V. Pole Assembly, Dual - Components - 6550-212-000

Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-041	I.V. Pole Assembly, 2 Stage, PR (pg. 110)	1
C	6500-001-042	I.V. Pole Assembly, 2 Stage, PL (pg. 112)	1

# 3-Stage I.V. Assembly, Patient Right - 6550-215-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 3-Stage I.V. Pole Assembly, PR - Components - 6550-215-000


Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-043	I.V. Pole Assembly, 3 Stage, PR (pg. 115)	1
C	6550-001-117	I.V. Pole Spacer Plate	1
D	6550-001-118	I.V. Clip Spacer Plate	1

# 3-Stage I.V. Pole Assembly, Patient Right - 6500-001-043



Item	Part No.	Part Name	Qty.
A	0025-133-000	Rivet	2
B	6070-090-105	Caution Label	1
C	6070-210-045	Sleeve	1
D	6070-210-046	Pivot	1
E	6070-210-049	Ring	1
F	6070-215-040	Pole Assembly	1
G	6100-115-050	Socket Weldment	1
H	6500-001-255	Label	1

# 3-Stage I.V. Assembly, Patient Left - 6550-216-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 3-Stage I.V. Pole Assembly, PL - Components - 6550-216-000

Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-044	I.V. Pole Assembly, 3 Stage, PL (pg. 117)	1
C	6550-001-117	I.V. Pole Spacer Plate	1
D	6550-001-118	I.V. Clip Spacer Plate	1




# 3-Stage I.V. Pole Assembly, Patient Left - 6500-001-044



Item	Part No.	Part Name	Qty.
A	0025-133-000	Rivet	2
B	6070-090-105	Caution Label	1
C	6070-210-045	Sleeve	1
D	6070-210-046	Pivot	1
E	6070-210-049	Ring	1
F	6070-215-040	Pole Assembly	1
G	6100-115-050	Socket Weldment	1
H	6500-001-256	Label	1

# 3-Stage I.V. Assembly, Dual - 6550-217-000




**NOTE: 0004-589-000 shown in Cot Assembly**

## 3-Stage I.V. Pole Assembly, Dual - Components - 6550-217-000

Item	Part No.	Part Name	Qty.
A	6100-115-060	I.V. Pole Clip	1
B	6500-001-043	I.V. Pole Assembly, 3 Stage, PR (pg. 115)	1
C	6550-001-044	I.V. Pole Assembly, 3 Stage, PL (pg. 117)	1

# No I.V. Pole Assembly Option - 6550-218-000

---




**NOTE: 0004-589-000 shown in Cot Assembly**

## No I.V. Pole Assembly Option - Components - 6550-218-000

Item	Part No.	Part Name	Qty.
A	6550-001-117	I.V. Pole Spacer Plate	1
B	6550-001-118	I.V. Clip Spacer Plate	1

# Optional Push Bar Assembly - 6550-040-000




## Push Bar Assembly Option - Components - 6550-040-000

Item	Part No.	Part Name	Qty.
A	0021-179-000	Set Screw	4
B	6550-001-026	<a href="#">Corner Assembly Handle (pg. 108)</a>	4
C	6550-001-108	Push Bar Sleeve	4
D	6550-001-199	Push Bar Storage Pouch	1

# Trouble Shooting Guide


## POWER PRO ELECTRONICS AND HYDRAULICS LOCATOR

Note: Some Components removed for clarity




# Trouble Shooting Guide

## POWER PRO HYDRAULIC MANIFOLD COMPONENTS LOCATOR



**DETAIL B**


**Wiring Schematics**  
Main Cable 8-Pin Connector



# Trouble Shooting Guide

## POWER PRO ELECTRICAL SYSTEM BLOCK DIAGRAM

Lower and Retract Functions




# Trouble Shooting Guide

---

## POWER PRO ELECTRICAL SYSTEM BLOCK DIAGRAM

Lift and Extend Functions





# Trouble Shooting Guide

## POWER PRO TROUBLESHOOTING GUIDE TABLE OF CONTENTS

Check for proper operation after each step. When the problem is fixed, return the cot to service. If assistance is needed at any time during troubleshooting, please contact a service technician at +44 (0) 1635 262431.

Intended Function	Problem	See Page, Section
Hold height position	Cot drifts down	126, 3
	Leaking hydraulic fluid	126, 4
Lift litter (with or without patient)	Does not lift (patient weight of 50 stones (700 lbs), (318 kg) or less)	126, 1
	Lifts slowly	126, 1
	Does not lift smoothly	127, 5
	Excessive motor noise	127, 6
Lower litter (with or without patient, powered)	Does not lower	127, 7
	Lowers too quickly	127, 8
	Does not lower smoothly	127, 8
	Cot is noisy while lowering	127, 5
Extend base (powered)	Base does not extend	126, 1
	Base extends slowly	126, 3
	Base doesn't extend smoothly	127, 5
	Base compresses significantly	127, 5 and 9
	Excessive motor noise	127, 6
Retracts base when loading (powered)	Base doesn't retract	128, 10
	Base retracts slowly	128, 10
	Base doesn't retract smoothly	128, 11
	Base doesn't stay retracted (drifts down)	127, 9
	Excessive motor noise	127, 6
Extend base (manual)	Does not allow extend	128, 12
	Extends too slowly	128, 12
Retract base (manual)	Does not allow retract	128, 13
Lower litter (with patient, manual)	Does not lower	128, 14
	Does not lower smoothly	128, 14
Lower litter (without patient, manual)	Does not lower	128 15
	Does not lower smoothly	128 15

# Trouble Shooting Guide

## POWER PRO TROUBLESHOOTING GUIDE

1. Does the battery indicator light when the button is pressed?

a.	<b>NO:</b> Replace battery with fully charged battery. If the indicator still does not light, the problem is with the battery terminals, or the control board. Contact a service technician at +44 (0) 1635 262431.
b.	<b>YES,</b> the indicator flashes red at a constant rate. Replace battery with fully charged battery.
c.	<b>YES,</b> the indicator flashes red with long and short flashes. This indicates an error in the electrical system. Check all electrical connections. If this doesn't eliminate the error code, there is a connection broken at some point in the electrical system. Please contact a service technician at +44 (0) 1635 262431.
d.	<b>YES,</b> the indicator lights solid green. Continue to number 2.

2. Does the motor turn on when the (+) button is pressed and held?

a.	<b>NO:</b> Replace the switch. If this doesn't work, the problem could be a faulty A-valve solenoid, A-valve, control board, or motor. Replace these components in that order, trying the function after each.
b.	<b>YES:</b> Does the cot drive up, and then drift down? i. <b>NO:</b> The problem is a faulty A-valve solenoid, A-valve, or motor. Replace these components in that order, trying the function after each. ii. <b>YES:</b> Continue to number 3.

3. The problem is most likely with the cap-side manual valve or the B-valve.

a.	Use the manual release lever to raise and lower the cot a couple of times.
b.	Push the cap-side manual valve in by tapping the stem with a punch and hammer.
c.	Loosen the manual release cable until it will no longer pull the valves. If this fixes the drift, readjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable).
d.	Replace the B-valve.
e.	Replace the cap-side manual valve.
f.	If none of the above fixes the drift, the manifold may need to be replaced. Please contact a service technician at +44 (0) 1635 262431.

4. For leaking hydraulic fluid, wipe the area clean so the source of the leak can be found, then tighten or replace the leaking component.

# Trouble Shooting Guide

---

5. There is most likely air in the hydraulics.

a.	Run the cot up and down several times.
b.	There may not be enough hydraulic fluid in the system. Refer to page 129 for instructions on filling the reservoir.
c.	If the problem persists, go to section 10 b.

6. Excessive motor noise is caused by a faulty motor, misaligned coupling, or pump. The latter two are internal to the hydraulics manifold. Please contact a service technician at +44 (0) 1635 262431.
7. Replace the switch. If this doesn't work, the problem could be a faulty B-valve solenoid, B-valve, cylinder, or control board. Replace these components in that order, trying the function after each.
8. Make sure all the casters are on the ground and that the load wheels and foot end handles are not being supported. Does the motor turn on?

a.	<b>NO:</b> See page 126, number 4 and number 5 above.
b.	<b>YES:</b> The problem could be a short in the cable to the pressure switch, a bad pressure switch or a sticking hydraulic cylinder. Inspect and replace these components as necessary.

9. The problem is most likely with the rod-side manual valve or the A-valve.

a.	Use the manual release lever to raise and lower the cot a couple of times.
b.	Push the rod-side manual valve in by tapping the stem with a punch and hammer.
c.	Loosen the manual release cable until it will no longer pull the valves. If this fixes the drift, readjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable).
d.	Replace the A-valve.
e.	Replace the rod-side manual valve.

# Trouble Shooting Guide

---

10. Press the (+) button to extend the legs briefly, then try the retract operation again.

a.	Use the manual release lever to raise and lower the cot a couple of times.
b.	Push the rod-side manual valve in by tapping the stem with a punch and hammer.
c.	Loosen the manual release cable until it will no longer pull the valves. If this fixes the drift, readjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable).
d.	Replace the rod-side manual valve.
e.	Replace the cylinder.
f.	The problem may be with the manifold assembly. Please contact a service technician at +44 (0) 1635 262431.

11. Run the cot up and down several times to remove any air in the system. If the problem persists, the problem could be the pressure switch, main cable, board, or cylinder. Replace these components in that order, trying the function after each.

12. Adjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable). If this does not fix the problem, replace the rod-side manual valve.

13. Adjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable). If this does not fix the problem, replace the cap-side manual valve.

14. Make sure the weight of the patient is being supported by the operators while actuating the manual release lever.

a.	Manually raise the cot litter slightly and try to lower again.
b.	Adjust the manual release cable until proper function is achieved (refer to page 129 on how to adjust the cable).
c.	Replace the cap-side manual valve.
d.	If the problem persists, replace the cylinder.

15. The operator may have to lift slightly on the foot end while lowering an empty cot. If this does not fix the problem, see number 14.

# Trouble Shooting Guide

---

## Adjusting the Manual Release Cable

Tools:

- 8mm wrench
  - 10mm wrench
1. Support the litter so no weight is on the base.
  2. Ensure the cable is intact.
  3. Using a 10mm wrench, loosen the cable lock nut.
  4. Using a 8mm hex wrench, adjust the tension on the manual release cable.
  5. Tighten cable lock nut.

## Filling the Reservoir - Use only Mobil Mercon Synthetic Blend Oil (6500-001-293)

**Note:** Any time you work with the hydraulics you may lose some oil.

Tools:


- 3/16" Allen wrench
1. Raise the cot to the full up position.
  2. Ensure the fill port is horizontal and lined up with the hole in the motor mount.
  3. Remove the port plug using a 3/16" allen wrench.
  4. Fill the reservoir up to the bottom of the fill port.
  5. Replace the plug and run the cot up and down a few times.

# Recycling Passport

BATTERY PACK: 6500-700-006



View of Battery



Item	Recycling/Material Code	Important Information	Qty
A	Battery Pack	NiCd	2

**NOTE**


Reference DeWALT® Manual (Part Number 6500-001-206) for battery recycling information.

# Recycling Passport

SWITCH HOUSING ASSEMBLY: 6550-001-036



Exploded view of Switch Housing Assembly




Item	Recycling/Material Code	Important Information	Qty
A	Printed Circuit Board		2

# Recycling Passport

**FOOT END ASSEMBLY: 6550-001-022**



Exploded view of Foot End Assembly



Item	Recycling/Material Code	Important Information	Qty
A	Printed Circuit Board	Contains Liquid Crystal Display	1




# Recycling Passport

**ELECTRONICS ASSEMBLY: 6550-001-014**



Exploded view of Electronics Control Board




Item	Recycling/Material Code	Important Information	Qty
A	External Electrical Cable		1
B	External Electrical Cable		1
C	External Electrical Cable		1
D	External Electrical Cable		1
E	External Electrical Cable		1
F	External Electrical Cable		1

# Recycling Passport

**HYDRAULIC POWER UNIT: 6550-001-030**



Exploded view of Hydraulics Assembly




Item	Recycling/Material Code	Important Information	Qty
A	Motor	Contains Automatic Transmission Fluid*	1
B	External Electrical Cable		1
C	External Electrical Cable		1
D	External Electrical Cable		1
E	External Electrical Cable		1

\* Mobil Mercon V Synthetic Blend or equivalent

# Recycling Passport

MAIN CABLE: 6550-001-172



Item	Recycling/Material Code	Important Information	Qty
A	External Electrical Cable		1

# Quick Reference Replacement Part List

## Note

The parts and accessories listed on this page are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call Stryker Customer Service at +44 (0) 1635 262431 for availability and pricing.

Part Name	Part Number
Backrest Pouch Option	6500-130-000
DC Battery Charger, 230V, International	6500-071-000
DC Battery Charger 12V/24V, In-Ambulance	6500-072-000
Equipment Hook	6500-700-003
Gas Cylinder, Fowler	1010-031-077
Hydraulic Oil	6500-001-293
I.V. Pole, 2-Stage, Right	6500-210-000
I.V. Pole, 2-Stage, Left	6500-211-000
I.V. Pole, 2-Stage, Dual	6500-212-000
I.V. Pole, 3-Stage, Right	6500-215-000
I.V. Pole, 3-Stage, Left	6500-216-000
I.V. Pole, 3-Stage, Dual	6500-217-000
Label, I.V. Pole Caution	6070-090-005
Label, Damage Warning	6080-090-009
Manual Installation / Operation, Cot Fastener	6385-009-001
Mattress, Bolster	6550-001-084
NiCd Battery Pack	6500-700-006
Restraint Belt Extension	6082-160-050
Restraint, Chest	6060-260-046
Restraint, Lap Belt (2 Used per Unit)	6060-160-044
Restraint Package	6082-260-010
Restraint, Shoulder Harness	6060-260-045
Restraint Strap Plastic Cap (short)	6082-160-051
Restraint Strap Plastic Cap (long)	6082-160-055
Touch-up Paint (Yellow)	6060-199-010
Touch-up Paint (Black)	6060-199-011
Wheel Lock	6082-200-010

# EMC Information

## POWER PRO TL

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
The model 6550 <b>RUGGED</b> Ambulance Cot is suitable for use in the electromagnetic environment specified below. The customer or the user of the model 6550 <b>RUGGED</b> Ambulance Cot should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	$\pm 6$ kV contact $\pm 8$ kV air	$\pm 6$ kV contact $\pm 8$ kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Applies to: <ul style="list-style-type: none"> <li>Cot,</li> <li>DeWalt AC/DC Charger (6500-070-000),</li> <li>DeWalt AC/DC Charger (6500-071-000),</li> <li>DeWalt DC/DC Charger (6500-072-000).</li> </ul>
Electrostatic fast Transient/burst IEC61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment. Applies to: <ul style="list-style-type: none"> <li>DeWalt AC/DC Charger (6500-070-000),</li> <li>DeWalt AC/DC Charger (6500-071-000).</li> </ul>
Surge IEC 61000-4-5	$\pm 8$ kV differential mode $\pm 2$ kV common mode	$\pm 8$ kV differential mode $\pm 2$ kV common mode	Main power quality is that of a typical commercial and/or hospital environment. Applies to: <ul style="list-style-type: none"> <li>DeWalt AC/DC Charger (6500-070-000),</li> <li>DeWalt AC/DC Charger (6500-071-000).</li> </ul>
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	$< 5\%U_t$ (95% dip $U_t$ ) for 0,5 cycle $40\%U_t$ (60% dip in $U_t$ ) for 5 cycles $70\%U_t$ (30% dip in $U_t$ ) for 25 cycles. $< 5\% U_t$ ( $> 95\%$ dip in $U_t$ ) for 5 sec.	$< 5\%U_t$ (95% dip $U_t$ ) for 0,5 cycle $40\%U_t$ (60% dip in $U_t$ ) for 5 cycles $70\%U_t$ (30% dip in $U_t$ ) for 25 cycles. $< 5\% U_t$ ( $> 95\%$ dip in $U_t$ ) for 5 sec.	Main power quality should be that of a typical commercial and/or hospital environment. If the user of the DeWalt AC/DC Charger (6500-070-000) or (6500-071-000) requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment. Applies to: <ul style="list-style-type: none"> <li>Cot,</li> <li>DeWalt AC/DC Charger (6500-070-000),</li> <li>DeWalt AC/DC Charger (6500-071-000),</li> <li>DeWalt DC/DC Charger (6500-072-000).</li> </ul>
Note: $U_T$ is the a.c. mains voltage prior to applications of the test level.			

# EMC Information


## POWER PRO TL (CONTINUED)

Recommended separation distances between portable and mobile RF communications equipment and the POWER PRO TL.			
The model 6550 <b>RUGGED</b> Ambulance Cot is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model 6550 <b>RUGGED</b> Ambulance Cot can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model 6550 <b>RUGGED</b> Ambulance Cot as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	m		
	150 kHz to 80 MHz $d=1,2 \sqrt{P}$	80 MHz to 800 MHz $d=1,2 \sqrt{P}$	8000 MHz to 2,5 GHz $d=2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
<b>NOTE 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
<b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

# EMC Information

## POWER PRO TL (CONTINUED)

The model 6550 **▲RUGGED▲** Ambulance Cot is suited for use in the electromagnetic environment specified below. The customer or the user of the model 6550 **▲RUGGED▲** Ambulance Cot should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 6100-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the model 6550 <b>▲RUGGED▲</b> Ambulance Cot, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.  Recommended Separation Distance  $d=1,2 \sqrt{P}$ Applies to Chargers: <ul style="list-style-type: none"> <li>• DeWalt AC/DC Charger (6500-070-000),</li> <li>• DeWalt AC/DC Charger (6500-071-000).</li> </ul>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d=1,2 \sqrt{P}$ 80 MHz to 800 MHz  Applies to Cot and Chargers: <ul style="list-style-type: none"> <li>• DeWalt AC/DC Charger (6500-070-000),</li> <li>• DeWalt AC/DC Charger (6500-071-000),</li> <li>• DeWalt DC/DC Charger (6500-072-000).</li> </ul> $d=2,3 \sqrt{P}$ 800 MHz to 2.5 GHz where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:  <div style="text-align: center;">  </div> Applies to Cot and Chargers: <ul style="list-style-type: none"> <li>• DeWalt AC/DC Charger (6500-070-000),</li> <li>• DeWalt AC/DC Charger (6500-071-000),</li> <li>• DeWalt DC/DC Charger (6500-072-000).</li> </ul>

# EMC Information

## POWER PRO TL (CONTINUED)

The model 6550 **▲RUGGED▲** Ambulance Cot is suited for use in the electromagnetic environment specified below. The customer or the user of the model 6550 **▲RUGGED▲** Ambulance Cot should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
---------------	----------------------	------------------	--

(Continued from previous page)

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model 6550 **▲RUGGED▲** Ambulance Cot is used exceeds the applicable RF compliance level above, the model 6550 **▲RUGGED▲** Ambulance Cot should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model 6550 **▲RUGGED▲** Ambulance Cot.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.



# EMC Information

## POWER PRO TL (CONTINUED)

<b>Guidance and Manufacturer's declaration - Electromagnetic Emissions</b>		
<p>The model 6550 <b>▲RUGGED▲</b> Ambulance Cot is intended for use in an electromagnetic environment specified below. The customer or the user of the model 6550 <b>▲RUGGED▲</b> Ambulance Cot should assure that it is used in such an environment.</p>		
Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	The model 6550 <b>▲RUGGED▲</b> Ambulance Cot uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Cot: Class A	The model 6550 <b>▲RUGGED▲</b> Ambulance Cot is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
	DeWalt AC/DC Charger (6500-070-000): Class B DeWalt AC/DC Charger (6500-071-000): Class B	The model 6550 <b>▲RUGGED▲</b> Ambulance Cot is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
	DeWalt DC/DC Charger (6500-072-000): Class B	
Harmonic Emissions IEC 61000-3-2	Cot: N/A	
	DeWalt AC/DC Charger (6500-070-000): Class A DeWalt AC/DC Charger (6500-071-000): Class A	
	DeWalt DC/DC Charger (6500-072-000): N/A	
Voltage Fluctuations Flicker Emissions IEC 6100-3-3	Cot: N/A	
	DeWalt AC/DC Charger (6500-070-000): Complies DeWalt AC/DC Charger (6500-071-000): Complies	
	DeWalt DC/DC Charger (6500-072-000): N/A	



STRYKER UK LIMITED  
Hambridge Rd.  
Newbury, Berkshire  
RG14 5EG, England

UNITED STATES  
Stryker Medical  
3800 E. Centre Ave.,  
Portage, Michigan USA  
49002

EC	REP
----	-----

**European Representative**  
Stryker France  
ZAC Satolas Green Pusignan  
Av. De Satolas Green  
69881 MEYZIEU Cedex  
France

**stryker**<sup>®</sup>