Stair-PRO®
Stair Chair/Chaise-civière/Treppenrollstuhl
Trapstoel/Sedia portantina/Silla para escaleras
Cadeira para escada/Trappstol/Trappestol/Porrastuoli
Αμαξίδιο για σκάλες/Krzeseło do przenoszenia po schodach

REF Model/Modèle/Modell/Modello/Modelo/Malli/Μοντέλο

6250/6251/6252

stryker®

Operations/Maintenance Manual

6250

6252

6251

Manuel d'utilisation et d'entretien
Bedienungs- und Wartungshandbuch
Gebruiks-/onderhoudshandleiding
Manuale d'uso e manutenzione
Manual de uso y mantenimiento
Manual de funcionamento/manutenção
Betjenings-/vedligeholdelsesmanual
Användar-/underhållshandbok
Käyttö- ja huolto-ohjekirja
Εγχειρίδιο λειτουργίας/συντήρησης
Instrukcja obsługi / konserwacji

International Addresses

EUROPE HEADQUARTERS

Stryker SA

Cite-Centre, Grand-Rue 90
P.O. Box 1568
1820 Montreux, Switzerland
Phone: 011-41-21-966-12-01
Fax: 011-41-21-966-12-00

EASTERN EUROPE

Stryker SA - Export Business

Grand-Rue 90
P.O. Box 1567
1820 Montreux, Switzerland
Phone: 011-41-21-966-14-00
Fax: 011-41-21-966-14-01

AUSTRIA

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Table of Contents

Symbols and Definitions	1-4
Symbols	1-4
Warning / Caution / Note Definition.	1-4
Introduction	1-5
Specifications	1-5
Component Identification	1-6
Warranty	1-7
Stryker EMS Return Policy.	1-8
Return Authorization.	1-8
Damaged Merchandise	1-8
International Warranty Clause.	1-8
Patent Information	1-8
Summary of Safety Precautions	1-9
Setup Procedures.	1-10
Operation Guide	1-11
Operating Guidelines	1-11
Unfolding the Chair	1-12
Folding the Chair	1-13
Transferring the Patient to the Stair-PRO® Stair Chair	1-14
Using Restraint Straps	1-15
Proper Lifting Techniques.	1-18
Transporting the Patient on Flat Surfaces.	1-19
Transporting the Patient Down Stairs	1-20
Transporting the Patient Up Stairs	1-23
Operating the Wheel Locks	1-24
Adjusting the Wheel Locking Force.	1-25
Operating the Optional Locking Rear Lift Handles.	1-26
Using the Optional Head Support (Model 6252 Only)	1-27
Using Additional Assistance	1-28
Cleaning.	1-29
Washing Procedure	1-29
Washing Limitations	1-29
Preventative Maintenance.	1-31
Checklist	1-31
Maintenance Record.	1-32
Training Record	1-33
Quick Replacement Part List.	1-34
Service Information.	1-35
Vinyl Seat Replacement	1-35
Track Belt Replacement (For Model 6252)	1-36
Track Belt Reconditioning (For Model 6252)	1-39
Rear Wheel Replacement.	1-41
Upper Control Handle Cable Replacement (For Model 6252).	1-42

Symbols and Definitions

English

SYMBOLS

Warning/Caution: consult accompanying documentation

Safe Working Load

WARNING / CAUTION / NOTE DEFINITION

The words WARNING, CAUTION and NOTE carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation, which if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation, which if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

NOTE

This provides special information to make maintenance easier or important instructions clearer.

Introduction

INTRODUCTION

This manual is designed to assist you with the operation and maintenance of the Model 6250, 6251, and 6252 Stair-PRO® Stair Chairs. Read it thoroughly before using the equipment or beginning any maintenance on it.

SPECIFICATIONS

Model	6250		6251		6252	
Maximum Load ³ 	500 lbs.	228 kg.	500 lbs.	228 kg.	500 lbs.	228 kg.
Height	36"	91 cm.	36"	91 cm.	37.5"	95.25 cm.
Width	20.5"	52 cm.	20.5"	52 cm.	20.5"	52 cm.
Depth ²	28"	71 cm.	28"	71 cm.	28"	71 cm.
Folded Depth	8"	20 cm.	8"	20 cm.	8"	20 cm.
Weight	20 lbs.	9 kg.	23 lbs.	10 kg.	31.5 lbs.	14 kg.

¹ Dimensions are measured from the outermost edges of the main frame. Specifications are rounded. Conversions are calculated before rounding.

² Depth dimensions are measured with extendable handles retracted.




³ Maximum load capacity is total weight distributed in accordance to basic human anatomy. Operators must consider the weight of the patient, equipment, and accessories when determining the total load on the product.

Stryker reserves the right to change specifications without notice.

Introduction

English

COMPONENT IDENTIFICATION

[Return To Table of Contents](#)

Stryker EMS, a division of the Stryker Corporation, offers two distinct warranty options in the United States:

One (1) year parts and labor. Under this option, Stryker EMS warrants to the original purchaser that its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of one (1) year after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts and labor for, or replacing, at its option, any product that is, in the sole discretion of Stryker, found to be defective.

Two (2) year parts. Under this option, Stryker EMS warrants to the original purchaser that non-expendable components of its products should be free from manufacturing non-conformances that affect product performance and customer satisfaction for a period of two (2) years after date of delivery. Stryker's obligation under this warranty is expressly limited to supplying replacement parts for, or replacing, at its option, any product which is, in the sole discretion of Stryker, found to be defective. Expendable components, i.e. mattresses, restraints, IV poles, storage nets, storage pouches, O2 straps, and other soft goods, have a one (1) year limited warranty with this option.

Under either warranty option, Stryker EMS products are designed for a 7 year expected service life under normal use, conditions, and with appropriate periodic maintenance as described in the maintenance manual for each device. Stryker warrants to the original purchaser that the welds on its EMS products will be free from structural defects for the expected 7 year life of the EMS product as long as the original purchaser owns the product. Original purchasers will also obtain a three (3) year limited parts warranty for the X-frame components of the MX-PRO R3 stretcher provided they also purchase X-frame guards at the time of the original purchase and the guards are installed on the MX-PRO before it is put into service.

If Stryker requests, products or parts for which an original purchaser makes a warranty claim, the purchaser shall return the product or part prepaid freight to Stryker's factory.

Any improper use or alteration or repair by unauthorized service providers in such a manner as in Stryker's judgment affects the product materially and adversely, shall void this warranty. Any repair of Stryker products using parts not provided or authorized by Stryker shall void this warranty. No employee or representative of Stryker is authorized to change this warranty in any way.

This statement constitutes Stryker EMS's entire warranty with respect to the aforesaid equipment. STRYKER MAKES NO OTHER WARRANTY OR REPRESENTATION EITHER EXPRESSED OR IMPLIED, EXCEPT AS SET FORTH HERIN. THERE IS NO WARRANTY OF MERCHANTABILITY AND THERE ARE NO WARRANTIES OF FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL STRYKER BE LIABLE HEREUNDER FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN ANY MANNER RELATED TO SALES OR USE OF ANY SUCH EQUIPMENT.

Warranty

English

STRYKER EMS RETURN POLICY

Cots, Stair Chairs, Evacuation Chairs, Cot Fasteners and Aftermarket Accessories may be returned up to 180 days of receipt if they meet the following guidelines:

Prior to 30 Days

- 30 day money back guarantee in effect
- Stryker EMS is responsible for all charges
- Returns will not be approved on modified items

Prior to 90 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 10% restocking fee

Prior to 180 Days

- Product must be **unused, undamaged** and in the original packaging
- Customer is responsible for a 25% restocking fee

RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items.

SPECIAL, MODIFIED, OR DISCONTINUED ITEMS NOT SUBJECT TO RETURN.

DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **DO NOT ACCEPT DAMAGED SHIPMENTS UNLESS SUCH DAMAGE IS NOTED ON THE DELIVERY RECEIPT AT THE TIME OF RECEIPT.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full.

Claims for any short shipment must be made within thirty (30) days of invoice.

INTERNATIONAL WARRANTY CLAUSE

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for additional information.

PATENT INFORMATION

Stryker products are covered by one or more of the following patents:

United States	5,575,026	6,276,010	6,648,343	6,908,133	6,796,757
	5,537,700	6,125,485	6,735,794	7,100,224	7,398,571
	D527,103				

Other Patents Pending

Summary of Safety Precautions

The following is a list of safety precautions that must be observed when operating or servicing this unit. The precautions are repeated throughout the manual, where applicable. Carefully read this list before using or servicing the unit.

WARNING

- Do not modify the Stair-PRO® Stair Chair. Modifying the chair can cause unpredictable operation resulting in injury to the patient or operator. Modifying the chair will also void its warranty.
 - Do not allow untrained helpers to assist in the operation of the Stair-PRO® Stair Chair. Untrained technicians/helpers can cause injury to the patient or themselves.
 - An unlocked chair can fold during use, causing injury to the patient or operator. Always make sure the chair is locked in the unfolded position before use.
 - To avoid injury, always verify the lift handles are locked in place before using them to lift the chair.
 - Always use all restraint straps to secure the patient on the chair. An unrestrained patient may fall from the chair and be injured.
 - Do not push the Model 6252 with the upper control handle in the fully extended position. Pushing the chair with the handle in the fully extended position may cause the chair to tip when obstacles are encountered.
 - Never leave a patient unattended on the chair or injury could result. Hold the chair securely while a patient is on the chair.
 - The Stair-PRO® Stair Chair is not recommended for use with suspected cervical, spinal, or fracture injuries.
 - To avoid injury, transporting the patient on stairs requires a minimum of two operators. If more people are required to safely control the chair, use the chart on [page 1-28](#) as reference for proper positioning of the helpers.
 - Only use the wheel locks during patient transfer or without a patient on the chair. Tipping could occur if the chair is moved while wheel locks are applied, resulting in injury to the patient or operator and/or damage to the chair.
 - Never use a wheel lock on a chair with excessively worn wheels. Using a wheel lock on a wheel with less than a 5" diameter (Model 6252) or a 6" diameter (Models 6250/6251) could compromise the holding ability of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the chair or other equipment.
 - To avoid injury, always verify the Stair-TREAD™ system on the Model 6252 is securely locked in place before transporting the patient.
 - Water, ice and debris on the stairs can affect operator footing and proper operation of the Stair-TREAD™ system. To avoid injury, clear the path or consider an alternate route.
 - Condensation, water, ice and/or debris on the Stair-TREAD™ system can cause unpredictable performance, resulting in a sudden change in the weight the operators must support. To avoid injury, and to aid proper operation of the Stair-TREAD™ system, ensure the belts are clean and dry before transporting the patient.
 - Never lubricate the Stair-TREAD™ system. Lubrication on the system can cause inconsistent operation possibly resulting in injury to the patient or operator.
 - To avoid injury to the operators and/or the patient, operators should never attempt to transport patient loads greater than what they can safely lift.
 - Use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the chair.
-

CAUTION

- Improper usage of the Stair-PRO® Stair Chair can cause injury to the patient or operator. Operate the Stair-PRO® Stair Chair only as described in this manual.
 - Improper maintenance can cause injury or damage to the unit. Maintain the Stair-PRO® Stair Chair as described in this manual. Use only Stryker approved parts and maintenance procedures. Using unapproved parts and procedures could cause unpredictable operation and/or injury and will void the product warranty.
 - Casters are not suitable for all surfaces. Caution should be used at all times.
 - Wheel locks are only intended to help prevent the empty chair from rolling while unattended and to aid in patient transfer. The wheel lock may not provide sufficient resistance on all surfaces or under loads.
 - Release the red track release bar before clicking the Stair-TREAD™ system into the locked position. Failure to follow this procedure could result in the track failing to lock. Always verify the Stair-TREAD™ system is locked by trying to fold it before descending stairs.
-

[Return To Table of Contents](#)

Setup Procedures

Unpack the cartons and check all items for proper operation. It is important that the Stair-PRO® Stair Chair is working properly before it is put into service. Have a qualified service person use the following list and the operation instructions to check the chair before it is put into service.

- All fasteners secure (reference all assembly drawings).
- All welds intact, not cracked or broken.
- No bent or broken tubing or sheet metal.
- No debris in wheels.
- All wheels secure and rolling properly.
- Chair unfolds and locks properly.
- No rips or cracks in seat or backrest.
- Patient restraints intact and working properly.
- Wheel locks operating properly.
- Foot end lift handles extend and lock properly.
- Head end lift handles fold and unfold.
- Front casters secure, rolling and swiveling properly (if equipped).
- Upper control handle extends and locks in all positions (if equipped).
- Stair-TREAD™ system unfolds and locks (if equipped).
- Track belts roll properly (if equipped).
- Optional accessories intact and operating properly.

The storage compartment on the vehicle in which the Stair-PRO® Stair Chair will be transported must be large enough to accommodate the folded dimensions of the chair (see [page 1-5](#)).

When necessary, modify the vehicle to fit the chair. Do not modify the chair.

 **WARNING**

Do not modify the Stair-PRO® Stair Chair. Modifying the chair can cause unpredictable operation resulting in injury to the patient or operator. Modifying the chair will also void its warranty.


OPERATING GUIDELINES

- Use the Stair-PRO® Stair Chair only as described in this manual.
- Read all labels and instructions on the chair before using the chair.
- When a patient is on the chair, use a minimum of two operators to manipulate the chair on stairs. If more people are required to safely control the chair, use the chart on [page 1-28](#) as reference for proper positioning of the helpers.
- Do not roll the chair, ascend, or descend stairs without advising the patient. Stay with the patient and control the chair at all times.
- Only use the wheel locks during patient transfer or without a patient on the chair.
- Always use the restraint straps when a patient is on the chair.
- Use properly trained helpers when necessary to control the chair and patient.

WARNING

- Always use all restraint straps to secure the patient on the chair. An unrestrained patient may fall from the chair and be injured.
 - Only use the wheel locks during patient transfer or without a patient on the chair. Tipping could occur if the chair is moved while the wheel locks are applied, resulting in injury to the patient or operator and/or damage to the chair.
 - Condensation, water, ice and/or debris on the Model 6252 Stair-TREAD™ system can cause unpredictable performance, resulting in a sudden change in the weight the operators must support. To avoid injury, and to aid proper operation of the Stair-TREAD™ system, ensure the belts are clean and dry before transporting the patient.
 - To avoid injury to the operators and/or the patient, operators should never attempt to transport patient loads greater than what they can safely lift.
-

UNFOLDING THE CHAIR

**Figure 1a - Unfolding the Chair
(6250 Shown)**


To unfold the chair (refer to Figure 1a above):

1. Stand behind the chair.
2. Apply the wheel locks (if desired).
3. Pull the backrest and the extension handle apart. The lock mechanism will automatically engage when the chair is completely unfolded.
4. Verify the lock is engaged by pulling up on the seat. If the lock is properly engaged, the chair will not fold.

⚠ WARNING

An unlocked chair can fold during use, causing injury to the patient or operator. Always make sure the chair is locked in the unfolded position before use.

FOLDING THE CHAIR

**Figure 1b - Folding the Chair
(6250 Shown)**

To fold the chair (refer to Figure 1b above):

1. Apply the wheel locks (if desired). Buckle the restraint straps and fold them neatly to prevent them from interfering with proper folding of the chair. If the chair is equipped with the optional head support strap, secure it behind the chair frame.
2. Stand at the side of the chair.
3. Pull up on the red lock bar at the rear of the chair.
4. Tip the chair forward.
5. Fold the seat up to the backrest until the front legs lock in the clips on the bottom of the seat tube.

Note: On models 6251 and 6252, rotate the front casters so they do not interfere with folding the chair.

⚠ WARNING

An unlocked chair can fold during use, causing injury to the patient or operator. Always make sure the chair is locked in the unfolded position before use.

TRANSFERRING THE PATIENT TO THE STAIR-PRO® STAIR CHAIR

To transfer the patient to the chair:

1. Place the chair beside the patient.
2. Apply the wheel locks to prevent the chair from moving.
3. Open the restraint straps.
4. Transfer the patient to the chair using accepted EMS procedures.
5. Use all the restraints to secure the patient on the chair ([page 1-15](#)).
6. Disengage the wheel locks before transporting.

WARNING

The Stair-PRO® is not recommended for use with suspected cervical, spinal, or fracture injuries.

USING RESTRAINT STRAPS


Figure 2
Cross-Chest Configuration




Figure 3
Chest and Lap Strap Configuration

Vinyl Seat
Top View

ABS Seat
Top View

Figure 4
Lap Strap
Attachment Points

WARNING

Always use all restraint straps to secure the patient on the chair. An unrestrained patient may fall from the chair and be injured.

The restraints can be used in two configurations.

For the cross-chest configuration (Figure 2): Buckle the chest restraints in an “X” pattern across the patient’s chest.

For the chest and lap configuration (Figure 3): Buckle one restraint across the patient’s chest and the other across the patient’s lap.

For either configuration: Match up the color coded restraints and buckle the ankle restraint across the patient’s legs.

To avoid damage to the buckles and straps, keep the restraint straps buckled when the chair is not being used with a patient.

When attaching the restraint straps to the chair, remember the attachment points must provide strong anchorage and proper restraint position while not interfering with equipment and accessories.

[Return To Table of Contents](#)

USING RESTRAINT STRAPS - CONTINUED

To attach the chest straps in the cross-chest configuration:

1. Wrap each strap around the chair frame, insert the end through the loop on the end of the strap and pull it tight.
2. Pull the strap across the patient's chest, lengthening the strap as necessary.
3. Buckle the strap.
4. Pull the loose end of the strap to tighten it securely around the patient.
5. Repeat for the second strap and the lap belt.

Figure 5 - Insert the end through the loop

Figure 6 - Pull the strap tight

Figure 7 - Lengthen strap as necessary

Figure 8 - Buckle strap

Figure 9 - Tighten strap securely

USING RESTRAINT STRAPS - CONTINUED

To attach the ankle strap:

1. Remove the three plastic pieces from the strap.
2. Thread the strap through one of the plastic “D” rings and slide the ring down to the end of the strap up to the white label.

3. Loop the strap around the front leg of the chair and pull the loose end of the strap through the ring. Pull until the ring rests against the foot rest tube.

4. Thread the strap through the other “D” ring.
5. Loop the strap around the other front leg of the chair and pull the loose end of the strap through the loop. Pull until the “D” ring rests against the foot rest tube.

6. Attach the final “male” end clip on the strap.


Figure 10-18 - Attaching the leg strap

USING RESTRAINT STRAPS - CONTINUED

- To lengthen the restraint, grasp the buckle, turn it at an angle to the webbing and pull it out. A hemmed tab at the end of the webbing prevents the buckle from coming off the strap. To shorten the restraint, grasp the hemmed tab and pull the webbing back through the buckle until the required tightness is achieved.
- When the chair is put into service, open the restraints and place them at either side of the chair until the patient is positioned on the seat. Lengthen the restraint, buckle it around the patient and shorten it until the required tightness to properly secure the patient is achieved.
- To open the restraint, press the tabs on the side of the buckle to release the buckle and pull the tang out of the receiver. To close the restraint, push the tang into the receiver until a "click" is heard.
- Whenever a restraint is buckled on a patient, the attendant should verify the tang is locked and the extra webbing is not tangled in the chair or hanging loose.
- Inspection of the restraints should be done at least once a month (more frequently if used heavily). Inspection should include checking for a bent or broken receiver or tang, torn or frayed webbing, etc. Any restraint showing wear or not operating properly must be replaced immediately.

PROPER LIFTING TECHNIQUES

When lifting the Stair-PRO® Stair Chair and patient, remember these five basic guidelines:

- Keep your hands close to your body.
- Keep your back straight.
- Coordinate your movements with your partner and lift with your legs.
- Avoid twisting.
- Always operate the Stair-PRO® Stair Chair as described in this manual.

TRANSPORTING THE PATIENT ON FLAT SURFACES

Figure 19 - Transporting the Patient on a Model 6250

Figure 20 - Transporting the Patient on a Model 6252

CAUTION

Casters are not suitable for all surfaces. Caution should be used at all times.

Model 6250

To roll the Model 6250 Stair-PRO® over flat surfaces, use either the locking rear lift handles (if equipped), or the backrest tube to tip the chair back. Once the chair is tipped back, use either the head end lift handles or the backrest tube to roll and guide the chair. Lift the chair over and around obstructions with the head end and foot end lift handles.

Model 6251

To roll the Model 6251 Stair-PRO® over flat surfaces, push and guide the chair from the rear of the patient, using either the head end lift handles or the backrest tube. Lift the chair over and around obstructions with the head end and foot end lift handles.

Model 6252

To roll the Model 6252 Stair-PRO® over flat surfaces, follow the same method as above for the Model 6251. In addition to the head end lift handles, the extendable upper control handle can be used in any position to roll and guide the chair. Extend the upper control handle by pulling the red release cable with one hand, and pulling up on the control handle with the other. Release the cable to lock the handle in either the intermediate or fully extended position. Lift the chair over and around obstructions with the head end and foot end lift handles.

WARNING

Do not push the Model 6252 with the upper control handle in the fully extended position. Pushing the chair with the handle in the fully extended position may cause the chair to tip when obstacles are encountered.

TRANSPORTING THE PATIENT DOWN STAIRS


Figure 21 - Foot end lift handles


Figure 22 - Head end lift handles

WARNING

To avoid injury, transporting the patient on stairs requires a minimum of two operators. If more people are required to safely control the chair, use the chart on [page 1-28](#) as reference for proper positioning of the helpers. To avoid injury, always verify the lift handles are locked in place before using them to lift the chair.

Models 6250 and 6251

1. Foot end operator - descend in front of the Stair-PRO® Stair Chair.
2. Head end operator - approach the stairs squarely.
3. Foot end operator - push each red foot end lift handle button and pull out each handle until it stops. Release the button and verify the handle is locked securely in position.
4. Head end operator - unfold the head end lift handles and face down the stairs.
5. Foot end operator - either face backward for improved patient monitoring or forward for an improved view of the stairs and easier maneuvering around obstacles.

Note: Although the front operator can face either direction while lifting and carrying, any applicable protocols for carrying chairs should be followed.


Figure 23 and 24 - Transporting down stairs

6. Both operators - simultaneously lift the chair, using the head and foot end lift handles and following the proper lifting techniques (see [page 1-18](#)). Carry the chair slowly down the stairs, avoiding any obstructions.

TRANSPORTING THE PATIENT DOWN STAIRS - CONTINUED

Model 6252

WARNING

- To avoid injury, always verify the Stair-TREAD™ system on the Model 6252 is locked in place before transporting the patient.
- To avoid injury, transporting the patient on stairs requires a minimum of two operators. If more people are required to safely control the chair, use the chart on [page 1-28](#) as reference for proper positioning of the helpers.
- To avoid injury, always verify the lift handles are locked in place before using them to lift the chair.

1. Roll the chair to the stairs and align it squarely with the edge of the first step.
2. Foot end operator - Extend the foot end lift handles by pushing the red release buttons and pulling the handles out until they stop. Release the buttons and verify the handles are locked.
3. Head end operator - Use one hand to pull the red upper control handle release cable while using the other hand to pull up and fully extend the handle. Release the cable and verify the handle is locked on both sides in the fully extended position.
4. Head end operator - Squeeze the red track release bar against the black cross tube. Relax your grip on the release bar and forcefully pull the Stair-TREAD™ system to the fully extended position until both sides lock securely. Always verify both sides of the Stair-TREAD™ system are locked by trying to fold it back up.

CAUTION

Release the red track release bar before clicking the Stair-TREAD™ system into the locked position. Failure to follow this procedure could result in the track failing to lock. Always verify the Stair-TREAD™ system is locked by trying to fold it before descending stairs.

5. Operators face each other while descending the stairs.
6. Head end operator - Tilt the chair back just far enough to allow the Stair-TREAD™ system to contact the floor.


Figure 25 - Foot end lift handles


Figure 26 - Upper control handle release cable


Figure 27 - Red track release bar


Figure 28 - Tilt the chair

TRANSPORTING THE PATIENT DOWN STAIRS - CONTINUED

Model 6252 - Continued

7. Both operators - Maintaining the angle, guide the Stair-PRO® Stair Chair over the edge of the stairs, allowing the Stair-TREAD™ system to engage the first step.
8. Both operators - Glide down the stairs until the treads are level across the edges of two or three steps.
9. Head end operator - Apply slight downward pressure on the extendable upper control handle while the foot end operator applies slight upward pressure on the foot end lift handles to keep the chair from rocking forward as it glides down the stairs.


Figure 29 - Track engaging first step

WARNING

- The Stair-TREAD™ system may not work the same on all stair surfaces and in all environmental conditions. Based on conditions, varying amounts of resistance may be encountered. Avoid getting dirt or other obstructions inside the tracks. Water, ice and/or debris on the stairs can affect operator footing and proper operation of the Stair-TREAD™ system. To avoid injury, clear the path or consider an alternate route. Condensation, water, ice and/or debris on the Stair-TREAD™ system can cause unpredictable performance, resulting in a sudden change in the weight the operators must support.
- To avoid injury, and to aid proper operation of the Stair-TREAD™ system, ensure the belts are clean and dry before transporting the patient.
- To avoid injury to the operators and/or the patient, operators should never attempt to transport patient loads greater than what they can safely lift.


Figure 30 - Transporting down the stairs

10. Foot end operator - when the track reaches the last step, release the front handles. Head end operator - allow the chair to tip forward until all four wheels are on the ground. Roll the chair as described on [page 1-19](#).
11. To fold the Stair-TREAD™ system, pull the red track release bar toward the black cross bar and fold the track up toward the chair. Verify the Stair-TREAD™ system is locked in place.


Figure 31 - "Resting position"

CAUTION

Release the red track release bar before clicking the Stair-TREAD™ system into the locked position. Failure to follow this procedure could result in the track failing to lock. Always verify the Stair-TREAD™ system is locked by trying to fold it before descending stairs.

If, while descending the stairs, either operator needs to pause or rest, tilt the chair forward just enough to allow the rear wheels to rest in the crook of the stair. To continue down the stairs from the resting position, the head end operator exerts slight downward pressure on the upper control handle while the foot end operator provides slight upward pressure to tilt the chair back and engage the Stair-TREAD™ system.


Figure 32 - Bottom of stairs

[Return To Table of Contents](#)

TRANSPORTING THE PATIENT UP STAIRS

WARNING

- To avoid injury, transporting the patient on stairs requires a minimum of two operators. If more people are required to safely control the chair, use the chart on [page 1-28](#) as reference for proper positioning of the helpers.
- To avoid injury, always verify the lift handles are locked in place before using them to lift the chair.

All Models

1. Roll the chair to the bottom of the stairs with the patient's back to the stairs.
2. Foot end operator - extend the foot end lift handles by pushing the red buttons and pulling the handles until they stop. Release the button and verify the handle is locked.
3. Head end operator - unfold the head end lift handles.
4. The foot end operator faces up the stairs. The head end operator may either face backward for improved patient monitoring or forward for an improved view of the stairs and easier maneuvering around obstacles.

Note: Although the head end operator can face either direction while carrying, any applicable protocols for carrying chairs should be followed.

5. Both operators - simultaneously lift the chair, using the head and foot end lift handles and following proper lifting techniques (see [page 1-18](#)). Carry the chair slowly up the stairs, avoiding any obstructions.


Figure 33 - Foot end lift handles


Figure 34 - Head end lift handles

Figure 35 - Transporting up the stairs



Figure 36 - Transporting up the stairs

Operation Guide

English

OPERATING THE WHEEL LOCKS

**Figure 37 - Wheel Locks
(6250 Shown)**


To activate the wheel locks:

1. Press down on the pedals (A) until they stop.
2. To release the wheel locks, depress the upper face of the pedal with your foot or lift up with your toe under the pedal. The upper portion of the pedal will rest against the chair frame when the wheel lock is released.

WARNING

- Only use wheel locks during patient transfer or without a patient on the chair. Tipping could occur if the chair is moved while the wheel locks are applied, resulting in injury to the patient or operator and/or damage to the chair.
- Wheel locks are only intended to help prevent the empty chair from rolling while unattended, and to aid in patient transfer. A wheel lock may not provide sufficient resistance on all surfaces or under loads.
- Never leave a patient unattended on the chair or injury could result. Hold the chair securely while a patient is on the chair.
- Never use a wheel lock on a chair with excessively worn wheels. Using a wheel lock on a wheel with less than a 5" diameter (Model 6252) or a 6" diameter (Models 6250/6251) could compromise the holding ability of the wheel lock, possibly resulting in injury to the patient or operator and/or damage to the chair or other equipment.

ADJUSTING THE WHEEL LOCKING FORCE

To adjust the wheel locking force:

1. Remove the screw from the center of the lock pedal.
2. Remove the octagonal sleeve (B) from the pedal (A).
3. Rotate the octagonal sleeve counterclockwise to increase the pedal locking force or clockwise to decrease the locking force.
4. Insert the octagonal sleeve (B) into the pedal (A).
5. Reinstall the screw into the center of the lock pedal.
6. Test the pedal locking force and verify that it holds properly before returning the chair to service.

Note: If, after adjustment, the pedal still does not hold the wheel properly, replace the wheel.

Operation Guide

English

OPERATING THE OPTIONAL LOCKING REAR LIFT HANDLES

WARNING

To avoid injury, always verify the lift handles are locked in place before using them to tip the chair back.


Figure 39 - Rotate the handle up



Figure 40 - Pull the trigger

To use the locking rear lift handles: Rotate the handles up until they lock into place. Verify the handles are locked securely into place before tipping the chair back.

To lower the handles: Lift up on the handle, pull the red trigger toward you with your thumb and fold the handle down against the chair frame.

USING THE OPTIONAL HEAD SUPPORT (MODEL 6252 ONLY)

Figures 41-43 - Attaching and using the head support

WARNING

The Stair-PRO® Stair Chair is not recommended for use with suspected cervical, spinal, or fracture injuries.

Before using the optional head support, the upper control handle must be extended. First, pull the red upper control handle release cable with one hand. Then, pull up on the handle with the other hand. Release the cable and verify the handle is securely locked into one of the two available positions.

To attach the optional head support to the extendable upper control handle:

1. Wrap the loose ends of the strap around the vertical portions of the handle, then feed them through the plastic loops.
2. Pull tight, and secure the strap to itself.
3. Adjust the height by loosening the strap, moving it to the desired location, and tightening it again.

To support the patient's head, position at the base of the head. For unconscious, or semiconscious patients, secure the head using the other two parts of the strap. Wrap around the patient's head, and overlap the straps to the desired tightness to secure. When not in use, these straps can be wrapped back around the handle and attached to the back of the support.


Operation Guide

English

USING ADDITIONAL ASSISTANCE

⚠ WARNING

To avoid injury, transporting the patient on stairs requires a minimum of two operators. If more people are required to safely control the chair, use the chart below as reference for proper positioning of the helpers.

[Return To Table of Contents](#)

The Model 6250, 6251 and 6252 Stair Chairs are designed to be power-washable. The unit may show some signs of oxidation or discoloration from continuous washing, however, no degradation of the chair's performance characteristics or functionality will occur due to power washing as long as the proper procedures are followed.

WASHING PROCEDURE

- Follow the cleaning solution manufacturer's dilution recommendations exactly.
- The preferred method Stryker Medical recommends for power washing stair chairs is with the standard hospital surgical cart washer or hand held wand unit.

WASHING LIMITATIONS

WARNING

Use any appropriate personal safety equipment (goggles, respirator, etc.) to avoid the risk of inhaling contagion. Use of power washing equipment can aerate contamination collected during the use of the chair.

- DO NOT STEAM CLEAN OR ULTRASONICALLY CLEAN THE UNIT.
 - Maximum water temperature should not exceed 180°F/82°C. Maximum air dry temperature (cart washers) is 240°F/115°C.
 - Maximum water pressure should not exceed 1500 psi/130.5 bar. If a hand held wand is being used to wash the unit, the pressure nozzle must be kept a minimum of 24 inches/61 centimeters from the unit. Failure to comply with these instructions may invalidate any and or all warranties.
-

If a foreign material gets between the Stair-TREAD™ system belt and track frame on the Model 6252, the track frame must be cleaned.

1. Loosen the track belts and remove the rear wheels (see [page 1-36](#) for instructions).
2. Clean the track frame completely with rubbing alcohol.
3. Use water at high pressure (see above) to rinse the belts. Be sure to clean both the inside and outside belt surfaces.
4. Allow the belts to completely dry before reassembly.
5. Reassemble the track belts and the rear wheels (see [page 1-36](#) for instructions).
6. Following the appropriate cautions and warnings, test the performance of the chair using a simulated patient weight while descending a flight of stairs.
7. If performance does not return to the original condition, the belts may need to be reconditioned or replaced.

Cleaning

In general, when used in those concentrations recommended by the manufacturer, either phenolic type or quaternary (excluding Virex® TB) type disinfectants can be used. Iodophor type disinfectants are not recommended for use because staining may result.

Suggested cleaners for the 6250, 6251, and 6252 surfaces:

- Quaternary Cleaners (active ingredient - ammonium chloride)
- Phenolic Cleaners (active ingredient - o-phenylphenol)
- Chlorinated Bleach Solution (5.25% - less than 1 part bleach to 100 parts water)

Avoid over saturation and ensure the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

WARNING

- SOME CLEANING PRODUCTS ARE CORROSIVE IN NATURE AND MAY CAUSE DAMAGE TO THE PRODUCT IF USED IMPROPERLY. If the products described above are used to clean Stryker patient care equipment, measures must be taken to insure the chairs are wiped with a cloth soaked in clean water and thoroughly dried following cleaning.
 - Failure to properly rinse and dry the chairs will leave a corrosive residue on the surface of the chairs, possibly causing premature corrosion of critical components.
-

Note: Failure to follow the above directions when using these types of cleaners may void this product's warranty.

Preventative Maintenance

Operation	Schedule	Procedure
Cleaning and Disinfecting	Each Use.	See page 1-30
Inspection	For 1-25 calls per month , inspect chair every 6 months For 26-200 calls per month , inspect chair every 3 months For 201+ calls per month , inspect chair monthly	See below for checklist
Stair-TREAD™ system belt reconditioning (6252)	After usage on approximately 500 flights of stairs or if the Stair-TREAD™ system performance is sluggish	See page 1-39

Note: Keep up to date maintenance records using the Maintenance Record form on [page 1-32](#).

Parts, Service or Technical Assistance

Contact Stryker Customer Service at 1-800-327-0770 or
 Stryker Medical
 3800 E. Centre Ave. Portage, MI 49002
 ATTN: Customer Service

CHECKLIST

- _____ All fasteners secure (reference all assembly drawings)
- _____ All welds intact, not cracked or broken
- _____ No bent or broken tubing or sheet metal
- _____ No debris in wheels
- _____ All wheels secure and rolling properly
- _____ Front casters secure, rolling and swiveling properly (if equipped)
- _____ Wheel locks hold wheels securely when on and clear the wheels when off
- _____ Chair unfolds and locks properly
- _____ No rips or cracks in seat or back rest
- _____ Restraint straps intact and working properly
- _____ Foot end carrying handles extend and lock properly
- _____ Head end carrying handles fold and unfold
- _____ Upper control handle extends and locks in all positions (if equipped)
- _____ Stair-TREAD™ system mechanism unfolds and locks properly (if equipped)
- _____ Stair-TREAD™ system belt rolls properly (if equipped)
- _____ Stair-TREAD™ system belt inner cords not showing - replace if necessary (if equipped)
- _____ Stair-TREAD™ system performs as desired - recondition belts if necessary (if equipped) (see [page 1-39](#))
- _____ No lubricants present on the Stair-TREAD™ system belts or the track frame surfaces
- _____ Upper release handle cable not worn or frayed - replace, if necessary (Model 6252)
- _____ Optional accessories intact and operating properly

Serial Number:		

Completed by: _____ Date: _____

Maintenance Record

English

Date	Maintenance Operation Performed	By	Hours

Training Record

English

Trainee Name	Training Date		Training Method
	Basic Training	Refresher Update	Owner's Manual, In-Service, Formal Class, Etc.

[Return To Table of Contents](#)

Quick Replacement Part List

English

The parts and accessories listed on this page are all currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call Stryker Customer Service (800)-327-0770 (Option 2), for availability and pricing.

Part Name	Part Number
Back Rest, Molded ABS	6250-001-116
Belt, Track	6252-001-085
Cable, Upper Handle Release	6252-001-016
Caster, Front (6251/6252)	6251-001-083
Handle Grip	6250-001-089
Restraint, Ankle	6250-001-127
Restraint, Chest, Black	6250-001-126
Restraint, Chest, Green	6250-001-125
Restraint Set, Polypropelene	6250-160-000
Restraint Set, Vinyl	6250-001-019
Seat, Molded ABS	6250-021-000
Seat, Vinyl	6250-020-000
Touch-Up Paint (Yellow)	6060-199-010
Touch-Up Paint (Black)	6060-199-011
Wheel, Rear, 5" (6252)	6252-001-114
Wheel, Rear, 6" (6250/6251)	6060-002-010
Wheel Lock Pedal	6080-200-030

VINYL SEAT REPLACEMENT

Tools Required:

- None

Figure 44 - Seat section of vinyl seat

Figure 45 - Back section of vinyl seat

Procedure:

1. Place the seat section on the frame, wrap it around the seat tube and fit the slots in the seat section over the clips on the underside of the seat tube. Be sure the “shiny” side of the seat is on the inside.
2. Tip the chair on its back to allow access to the bottom of the seat tube.
3. Insert both ends of the front strap through the loops in the seat section. Thread the end of the strap through the buckle and pull it tight.
4. Thread the end of the back strap through the buckle and pull it tight.
5. Tip the chair back up on its wheels.
6. Wrap the back section around the back tube.
7. Insert both ends of the top strap through the loops in the back section. Thread the end of the strap through the buckle and pull it tight.
8. Thread the end of the bottom strap through the buckle and pull it tight.

TRACK BELT REPLACEMENT (FOR MODEL 6252)

Tools Required:

- (2) 7/16" Wrenches
- 3/16" Hex Key
- 1/2" Wrench
- 5/32" Hex Key

Procedure:

1. Using a 7/16" wrench and a 3/16" hex key, remove the socket head cap screw, nut and washer from the track support bracket **on only one side** of the chair.

Figure 46 - Remove the screw, nut and washer

2. Turn the chair upside down so it is resting on the seat frame and back frame as shown in Figure 47.

Figure 47 - Turn the chair upside down

3. Using a 1/2" wrench and a 3/16" hex key, remove the button head cap screw, nut, spacers and wheel from the frame on the same side as the parts removed in step one.

Figure 48 - Remove the wheel assembly

TRACK BELT REPLACEMENT (FOR MODEL 6252) - CONTINUED

4. Using a 5/32" hex key, remove the two button head cap screws from the track spacer pivots and remove the pivots.

Figure 49 - Remove the track spacer pivots

5. Using two 7/16" wrenches, remove the hex head cap screw, nut and washers from both head end track rollers and remove the rollers.

Figure 50 - Remove the rollers

6. Remove the belts by threading them through the wheel support plate on the side where the parts were removed.


Figure 51 - Remove the worn belts

TRACK BELT REPLACEMENT (FOR MODEL 6252) - CONTINUED

7. Thread the new belts through the wheel support plate.

Figure 52 - Attach the new track belt

8. Replace the track rollers, washers, hex head cap screws and nuts.
9. Use two 7/16" wrenches to apply the desired tension to each belt and tighten securely.

Figure 53 - Use wrenches to apply tension to the belts

Note: To check the tension, pull up on the center of the belt until it is taut. The gap between the belt and the track frame should measure between 3/8" and 1" as shown in Figure 54.

10. Spin the belts to verify they roll freely.


Figure 54 - Proper belt tension

TRACK BELT RECONDITIONING (FOR MODEL 6252)

Tools Required:

- 7/16" Wrenches
- Permanent Marker
- 50 Grit Sandpaper
- Sanding Block

Procedure:

1. Extend the upper control handle, open the Stair-TREAD™ system tracks and tip the chair forward until it rests on the handle and the seat section as shown in Figure 55.

Figure 55 - Tip the chair forward

2. Using two 7/16" wrenches, remove the hex head bolt, hex nut, two washers, spacer and wheel from both sides of the track frame. Remove the wheel and spacer down in the direction of the floor.


Figure 56 - Remove the wheel spacer in the direction of the floor

3. Turn both belts over so the small internal teeth are on the outside. Put the belts on the track frame for support while sanding them.

Note: New track belts must be installed if the inner cords are exposed on the wear surfaces.

4. For a start/end point reference while sanding, use a permanent marker to color a tooth on the belt.

Figure 57 - Mark a belt tooth for reference

TRACK BELT RECONDITIONING (FOR MODEL 6252) - CONTINUED

- Using a sanding block with 50 Grit sandpaper, sand both wear surfaces on the outside edges of each belt. The purpose is to roughen the surface of the belt just enough to remove the shiny spots caused by normal wear. Do not over sand the belts.

 **WARNING**

Do not sand the track teeth. Deformation of the teeth can cause unpredictable chair performance resulting in injury to the operators and/or patient.

Figure 58 - Sand the wear surfaces

- Remove all sanding debris from the belts and the track frame.
- Turn both belts back over so the large teeth are on the outside.
- Replace the hex head bolt, hex nut, two washers, spacer and wheel on both sides of the track frame.


Figure 59 - Belt sanding surfaces

- Use two 7/16" wrenches to apply the desired tension to each belt and tighten securely.

Figure 60 - Use wrenches to apply tension to the belt

Note: To check the tension, pull up on the center of the belt until it is taut. The gap between the belt and the track frame should measure between 3/8" and 1" as shown in Figure 61.

- Spin the belts to verify they roll freely.


Figure 61 - Proper belt tension

REAR WHEEL REPLACEMENT

Tools Required:

- 1/2" Wrench
- 3/16" Hex Key

Procedure:

1. Turn the chair upside down so it is resting on the seat frame and back frame as shown in Figure 62.

Figure 62 - Turn the chair upside down

2. Using a 1/2" wrench and a 3/16" hex key, remove the button head cap screw, nut, spacers and wheel from the frame. Discard the cap screw and nut.
3. Use the new fasteners provided to attach the new wheel to the frame.
4. Repeat for the other wheel, if necessary.

Figure 63 - Remove the wheel assembly

Service Information

English

UPPER CONTROL HANDLE CABLE REPLACEMENT (FOR MODEL 6252)

Tools Required:

- 5/8" Wrench

Procedure:

1. Pull the cable and raise the extendable upper control handle to the fully extended position.
2. Using a 5/8" wrench, remove one side of the cable.
3. Place the washer on the end of the new cable.
4. Using a 5/8" wrench, attach one side of the new cable to the chair frame in the hole where the old cable was removed.
5. Repeat steps 2-4 for the other end of the cable.

Figure 64 - Raise the extendable handle


Figure 65 - Remove one side of the cable

Figure 66 - Place the washer on the new cable

Figure 67 - Attach one side of the new cable

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