# Floor plate installation instructions

This procedure shows you how to install the floor plate with a wooden floor. These instructions assume the vehicle has a wood floor with a thickness of at least 9/16 in. (1.4 cm). There are many different vehicle configurations and you should plan your installation accordingly.

## Floor plate assembly kit (6390-700-001, 639007000100, or 639207000001) or Universal floor plate assembly kit (6390-028-000) contents

(2) Flat head cap screw (0001-194-000)	(2) Set screw (0021-197-000)	(1) Rubber grommet (0037-247-000)
(1) Floor plate attachment bracket (6390-101-108)	(2) Floor plate end cap (639000010111)	(1) Universal safety hook (639000010148)
(1) Wire protection loom (6390-001-153)	(1) Drain tube (6390-001-170)	(1) Floor plate drain tube (6390-001-183)
(6) P-style clamp (6390-001-202)	(1) Floor plate (6390-001-107)	(1) Anchor-to-vehicle cable (6390-001-135) - for 6390-700-001 kit only <b>OR</b> (1) Anchor-to-vehicle cable (639000010135) - for 639007000100 or 639207000001 kits only

(1) Floor plate cap (6390-001-109) - not included with the 6390-028-000 kit

## Additional parts used (not included)

Five sets of 3/8-16 UNC flat head cap screws (ASTM-F835 or SAE grade 8), flat washers, and lock nuts	Washers, 3/8" ID (shims to correct floor unevenness)	1/4" screws (for wire protection loom)
--	--	--

## Tools required:

Pencil	Saw	5/32" hex wrench
Router	1/2" to 1" router bit	Shop vacuum
Drill	1/2" drill bit	Silicone sealant
Tape measure	7/16" drill bit	Rubber hammer
1/2" drive 7/32" hex	9/16" drill bit	Safety glasses
Face mask	1/4" hex wrench	Extension cord
Multimeter	1/2" drive torque wrench (ft-lb) > 40 ft-lb (54.3 N-m)	

## Procedure:

**Note** - Five sets of 3/8-16 UNC flat head cap screws (ASTM-F835 or SAE grade 8), flat washers, and lock nuts are required, but are not included with your installation kit. Examine your vehicle support structure to determine the best length bolt hardware for your floor plate installation. Use these bolts to secure the floor plate as described in step 17.

1. Identify the back-to-front mounting location relative to the door sill. Start at the location where the door sill first drops off (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)). Measure 2-3/4 in. (7 cm) inwards. Use a pencil to mark the location along the width of the vehicle patient compartment.

**Note** - You can install the floor plate closer to the rear loading door if your door sill meets these requirements:

- A maximum vertical door sill height (distance between the vehicle deck and the door sill) of 1 in. (2.54 cm)
  - A minimum door sill depth (beyond the floor plate) of 5 in. (12.7 cm)
  - The door sill is clear of all obstructions in the area where the cot load wheels travel
2. Using a 5/32" hex wrench, install the two set screws (0021-197-000) to attach the floor plate end caps (639000010111) on the rear end and front end of the floor plate (6390-001-107).
  3. Position the floor plate the desired distance from the vehicle patient compartment walls (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)).
    - Align the edge of the installed end cap with the marked location (determined in step 1).
    - Use a pencil to mark the center location of the floor plate on the vehicle floor at the rear end and front end of the end caps.
    - Measure from these marks to the vehicle wall to make sure that they are parallel.

**Note** - Install the optional wheel guide (6390-027-000) if the cot center line is 17-1/2" or less from the vehicle wall.

---

**WARNING** - Always secure the item that you are cutting. Be aware of the area around your cutting location. Always wear appropriate eye protection while you operate a saw.

---

**Note** - If the floor plate does not fit lengthwise in the vehicle patient compartment, use a saw to shorten the floor plate (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)). Before cutting the floor plate, remove the front end floor plate end cap. Replace the cap after cutting the floor plate.

4. Hold the floor plate in place. Use the floor plate as a template, trace the inset edge of the floor plate and end caps onto the vehicle patient compartment floor, see Figure 1.

**Note** - If you are installing multiple floor plates, you can make a template to expedite installation.

**Figure 1 – Trace the floor plate**

5. Remove the floor plate and inspect the trace. Follow this trace when cutting the floor plate pocket for the floor plate.

**Note**

- Slightly oversize the cutout to make sure that the floor plate fits into the floor plate pocket. The flange of the floor plate will cover any minor widening of the floor plate pocket.
- It may be necessary to notch your sill plate, so your floor plate sits flat on the floor.

Steps 6 through step 8 explain the process for creating the floor plate pocket in the vehicle floor to accept the floor plate.

Here are a few alternative installations:

Floor thickness	Action required
Less than 7/16 in. (1.1 cm)	Build up your floor to a thickness of at least 7/16 in. (1.1 cm) for the entire width of the ambulance floor. The cot load wheels must make full contact with the floor during loading and unloading.

<p><b>7/16 in. (1.1 cm) to 9/16 in. (1.4 cm)</b></p>	<p>Route (or cut) down to the vehicle support structure. When you install the floor plate, there will be a gap of not more than 1/8 in. (0.32 cm) between the bottom of the flange of the floor plate and the floor. Fill this gap with silicone sealant.</p>
<p><b>Greater than 9/16 in. (1.4 cm) or with a foam structure</b></p>	<p>Route (or cut) the floor plate pocket to the depth of the solid vehicle support structure and then use spacers (for example, 3/8 in. (0.95 cm) washers – not provided) to build the vehicle support structure up to the bottom of the floor plate. For example, if your vehicle floor is 3/4 in. (1.9 cm), you need to install 3/16 in. (0.5 cm) spacers or washers to support the floor plate at the five mounting holes. This has the following advantages:</p> <ul style="list-style-type: none"> <li>• Easier identification of the patient vehicle compartment support structure</li> <li>• Less cutting of the patient vehicle compartment</li> <li>• May eliminate the need to router an additional floor plate pocket for the electrical rubber grommet (see step 14).</li> </ul>

**WARNING** - Always wear safety glasses and a face mask while operating a router, during installation.

6. Using a router and 1/2" to 1" router bit, router out the floor by following the lines marked on the floor in step 4.

**Note** - Set the router bit depth to 9/16 in. (1.4 cm). Use a firm stance when using the router to avoid damaging the floor.

7. Set the floor plate into the floor plate pocket to check the fit. Make any adjustments, as necessary.

**Note** - The flange of the floor plate will cover any minor widening of the floor plate pocket.

8. Identify and mark at least one of the five mounting holes in each of the five (5) five-hole patterns that align with the vehicle support structure. If none of the five mounting holes in any one hole pattern lines up with a support structure, you need to create a support structure for that location, such as a plate or bracket.

**Note**

- These mounting holes will be drilled out to 7/16 in. (1.1 cm) diameter.
- Each bolt and support structure must be able to withstand 2,750 lbf (12,230 N) in upward tension and 600 lbf (2,669 N) in shear in all horizontal directions.
- Consider any restrictions due to clearances with any vehicle component.

**WARNING** - Always be aware of items around and below the electrical inlet during floor plate installation. Consult the vehicle manufacturer before installing. Make sure that you do not damage or interfere with brake lines, oxygen lines, fuel lines, fuel tank, or electrical wiring of the vehicle.

**Note** - The electrical inlet hole and drain tube assembly hole do not need to align with the support structure.

9. Identify and mark the electrical inlet hole (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)). This hole will be drilled to 1/2 in. (1.3 cm) diameter.

10. Identify and mark the fifth hole from the foot end for the drain tube assembly (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)). This hole will be drilled to 9/16 in. (1.4 cm) diameter.

**Note** - If this hole location will not work for your drain tube installation, determine the fastener assembly mounting position before selecting an unused drain tube hole. Make sure that the drain does not interfere with the fastener attachment mechanism.

11. Remove the floor plate from the floor plate pocket.

12. Using a drill and appropriate drill bit, drill the following marked holes through the floor plate:

Hole	Quantity	Diameter	Drill bit
Mounting	5	7/16 in. (1.1 cm)	7/16 in.
Electrical inlet	1	1/2 in. (1.3 cm)	1/2 in.
Drain tube	1	9/16 in. (1.4 cm)	9/16 in.

13. Install the floor plate into the floor plate pocket. Transfer holes onto the floor plate pocket in the vehicle patient compartment floor. Remove the floor plate. Drill the corresponding holes into the vehicle patient compartment floor.
14. Using a drill or router, cut a 1 in. (2.54 cm) diameter hole to a depth of 1/4 in. (.635 cm) below the floor plate pocket for the electrical rubber grommet.
15. Insert the supplied electrical rubber grommet (0037-247-000) into the floor plate at the electrical input location.
16. Install the appropriate length 3/8-16 UNC flat head cap screws, ASTM-F835 or SAE grade 8 with a flat washer and lock nut.
 

**Note** - This hardware is required, but not included with your installation kit, because vehicle support structures vary.
17. Using a 1/2" drive torque wrench with a 1/2" drive 7/32" hex, tighten the floor plate tie-down bolts to a minimum torque of 40 ft-lb (54.3 N-m) in at least a two-step tightening process.
18. After you tighten the bolts, make sure that the floor plate is flat.
  - a. If the floor plate is not flat, remove the bolts and reinstall using washers.
  - b. If the floor plate is flat, remove the bolts, apply silicone sealant to the underside of the floor plate flange, and press the floor plate firmly into the floor plate pocket. Reinstall the bolts.

---

**WARNING** - Always make sure that all gaps to the exterior of the vehicle are sealed to prevent exhaust fumes from entering the vehicle patient compartment.

---

**Note** - You can use a silicone adhesive to secure the washers to the floor

19. Attach the drain tube (6390-001-170) to the floor plate drain tube (6390-001-183) and insert it into the floor plate drain tube assembly hole.
 

**Note**

  - You can cut the drain tube to meet the needs of your installation.
  - Alternatively, you can install the drain tube from the underside of the vehicle.
20. Using a 1/4" hex wrench, turn the drain a quarter turn to lock the drain in place.
21. Route the drain tube under the vehicle away from the exhaust system ensuring that the open end of the tube is lower than the drain for proper drainage. Secure the drain tube to the vehicle using appropriate tie downs (not supplied).
22. If you are not immediately installing a cot fastener, install the floor plate cap (6390-001-109) over the floor plate.

## Installing the safety hook

---

**WARNING** - Failure to install the safety hook can cause injury to the patient or operator. Install and use the safety hook as described in this manual.

---

### Tools required:

- 5/32" hex wrench

### Procedure:

1. Insert a floor plate attachment bracket (6390-101-108) at the foot end of the vehicle patient compartment to install the safety hook (639000010148). Line the notches with the tabs to insert, then slide the bracket into the floor plate assembly.


Figure 2 – Insert floor plate attachment bracket


Figure 3 – Secure safety hook

- Using a 5/32" hex wrench, install two flat head cap screws (0001-194-000) to secure the safety hook to the floor plate assembly.

## Connecting the electrical circuit

Follow these requirements and recommendations to connect the 12V DC electrical circuit between the vehicle and your product (Model 6390 Power-LOAD and Model 6392 Performance-LOAD).

### WARNING

- Always make sure that the vehicle circuit can deliver 12.8V-15.6V DC under a maximum load of 10A.
- Do not connect to a 24V DC vehicle circuit.
- Always connect **Power-LOAD** to a 15 A fuse/resettable breaker that is installed in the vehicle to prevent electrical hazards. Stryker recommends that you install the fuse/resettable breaker in the fuse box.
- Always be aware of items around and below the anchor-to-vehicle cable during floor plate installation.

**Note** - If you do not follow these requirements and recommendations, **Power-LOAD** battery failure may result.


### Wiring installation

#### Preferred method:

Wire the product to the vehicle battery bank (Figure 4).


#### In rare cases:

If your vehicle sits without idling or connecting to the shoreline for longer than one month, wire the product to the ignition or shore circuit (Figure 5) to prevent vehicle battery failure. For example, a show vehicle.

A	Vehicle/module battery bank	E	Shoreline in
B	Shoreline	F	Cot fastener
C	Modular/box disconnect	G	Continuous 12V power
D	15A fuse/resettable breaker	H	Extra module/box functions

Figure 4 – Connecting to the battery bank

A	Vehicle/module battery bank	F	Cot fastener
B	Shoreline	G	Continuous 12V power
C	Modular/box disconnect	H	Extra module/box functions
D	15A fuse/resettable breaker	J	Vehicle ignition
E	Shoreline in		

Figure 5 – Connecting to the vehicle ignition

To prevent vehicle system or modular battery depletion, the product is equipped with safeguards to protect its power source, such as:

- Charging will not start if the  $V_{source} < 12.7V$

- Charging will be limited to a two-hour time limit if  $V_{source} < 12.2V$ . After the two-hour time limit, charging will not start until  $V_{source} > 12.5V$ .
  - Charging will cease if  $V_{source} < 11.7V$
  - Anchor electronics will enter Sleep-Mode if  $V_{source} < 11.1V$
  - The maximum drain rate of the circuit is 10A.
1. Identify the point of connection between the supplied anchor-to-vehicle cable (6390-001-135 or 639000010135) and the vehicle electrical system (Electrical requirements: 12.8V-15.6V DC, 15A fuse/resettable breaker, 2 conductor 10 AWG cable).
  2. Install the anchor-to-vehicle cable through the electrical rubber grommet in the electrical inlet. Leave enough slack (approximately 6 in. (15 cm)), so you can connect the spade-connector end of the anchor-to-vehicle cable to the anchor cable.
  3. Route the anchor-to-vehicle cable back to the vehicle connection point. Twelve feet of wire protection loom (6390-001-153) and six p-clamps (6390-001-202) are provided.
  4. Attach the 15A fuse/resettable breaker to the circuit.
  5. Using a multimeter, check the vehicle voltage at the spade-connector end of the anchor-to-vehicle cable (6390-001-135 or 639000010135) to make sure that the voltage is 12.8V-15.6V DC when the vehicle is off, idle, and while revving the engine.
 

**Note** - Before you check for voltage, secure any additional interconnects (turn on the vehicle ignition, modular connection, or dedicated switch), if applicable.
  6. Connect the anchor-to-vehicle cable to the vehicle electrical system.

## Floor plate configuration

Select your configuration (*Floor plate with single standard floor mount cot fastener* (page 8) or *Floor plate with single standard wall mount cot fastener* (page 9)). All measurements are shown in inches.

**Note** - Five sets of 3/8-16 UNC flat head cap screws (ASTM-F835 or SAE grade 8), flat washers, and lock nuts are required, but are not included with your installation kit. Examine your vehicle support structure to determine the best length bolt hardware for your floor plate installation. Use bolts to secure the floor plate as described in step 17.

1. The foot end of the floor plate should be 2-3/4 in. (7 cm) from the location where the door sill first drops off.
 

**Note** - You can install the floor plate closer to the rear loading door if your door sill meets these requirements:

  - A maximum vertical door sill height (distance between the vehicle deck and the door sill) of 1 in. (2.54 cm)
  - A minimum door sill depth (beyond the floor plate) of 5 in. (12.7 cm)
  - The door sill is clear of all obstructions in the area where the cot load wheels travel
2. You can cut the floor plate, if necessary (minimum 89.50 in. (227.33 cm) for **Power-LOAD**).
3. Router depth is 9/16 in. (1.4 cm) for the floor plate pocket. Floor plate pocket width is approximately 2-5/8 in. (6.5 cm).
4. Floor mount plate - 28-7/8 in. (73.3 cm) / wall mount plate - 25-1/4 in. (64.1 cm).
 


**Note**

  - Dimension 4 should match the dimension above when you install the anchor on the floor plate in the rear most location (toward cot foot end). If you install the anchor forward of the rear most location, dimension 4 will increase by the same amount.
  - For a pre-existing rail clamp plate, if dimension 4 measures shorter than the above values, move the rail clamp plate to the appropriate location.
5. Choose one hole to drill through  $\varnothing$  1/2 in. (1.3 cm) for electrical input.
6. Drill the fifth hole from the foot end through  $\varnothing$  9/16 in. (1.4 cm) for the drain tube assembly.
7. Select any one bolt location from each of the five bolt patterns.
8. Measure this dimension from the sill edge to the floor plate pocket.

**Note**

- Use 3/8-16” UNC flat head cap screws, ASTM-F835 or SAE grade 8 bolts. Each bolt should have a flat washer and lock nut.
- Install each bolt with a minimum of 40 ft-lb (54.3 N-m) torque.
- Each bolt and support structure must be able to withstand 2,750 lbf (12,230 N) in upward tension and 600 lbf (2,669 N) in shear in all horizontal directions.
- Alternatively, you can install a tap plate. Secure the 3/8-16” flat head cap screws directly into the tap plate (tap plate must meet the strength requirements). Use an appropriate thread locker to prevent the screws from loosening.

**Floor plate with single standard floor mount cot fastener**

**Figure 6 – Floor plate with single standard floor mount cot fastener (optional)**

A	Sill depth
B	Sill edge
C	Sill
D	Sill height
E	Foot end (rear)
F	Floor plate (mini rail clamp option)
G	Wall of vehicle
H	Cot perimeter
I	Cot centerline/ floor plate centerline

J	Universal floor plate
K	Head end (front)

## Floor plate with single standard wall mount cot fastener


Figure 7 – Floor plate with single standard wall mount cot fastener (optional)

A	Sill depth
B	Sill edge
C	Sill
D	Sill height
E	Foot end (rear)
F	Floor plate (mini rail clamp option)
G	Wall of vehicle

H	Cot perimeter
I	Cot centerline/ floor plate centerline
J	Universal floor plate
K	Head end (front)
L	Wall plate (mini rail clamp option)
M	Wall mounting bracket
N	Floor of vehicle

# Installation checklist

Follow this checklist after you install your floor plate.

- \_\_\_\_\_ Make sure that you do not have any unused components after installation. Your floor plate does not ship with any extra components. If you have any unused components after installation, call Stryker service.
- \_\_\_\_\_ Make sure that all floor plate tie-down bolts are set to a minimum torque of 40 ft-lb (54.3N-m) in at least a two step tightening process.
- \_\_\_\_\_ Visually check that the floor plate sits flat on the floor deck.
- \_\_\_\_\_ Make sure that the drain tube was installed under the vehicle and away from the exhaust system so the open end of the tube is lower than the drain for proper drainage.
- \_\_\_\_\_ Make sure that the safety hook is fastened to the floor plate assembly and is secure.
- \_\_\_\_\_ Confirm that the vehicle circuit can deliver 12.8V-15.6V DC under a maximum load of 10A (for **Power-LOAD** applications).

Product serial number:			
Installed by:		Date:	
Installed by:		Date:	

**Note** - Maintain a copy of this record for at least seven years.

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