0086

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<u>Section</u>	<u>Description</u>	<u>Page</u>
1.0	Warranty.....	2
2.0	Symbols.....	3
3.0	Indications for Use.....	3
4.0	Contraindications.....	3
5.0	Safety Precautions.....	4
6.0	Description.....	5
7.0	Features	5
8.0	Control Panel.....	6
9.0	Setting Up the System.....	7
10.0	Instructions for Use	8
11.0	CPR (Cardio Pulmonary Resuscitation).....	9
12.0	System Cleaning, Care and Storage.....	9
13.0	Troubleshooting.....	10
14.0	Specifications, Control Unit.....	10

Before you begin...

 **Important**

Before using the Air Express® system, please read and understand this manual and all safety precautions prior to each application.

Only qualified medical service personnel should attempt to repair this device. A *Service Manual* is available for use by qualified personnel. For assistance contact your local dealer. If additional assistance is needed, contact Gaymar's technical service department.

Customers outside USA:

Direct +1 (716) 662-2551
Fax: +1 (716) 662-0730

Customers inside USA:

Toll Free 1 (800) 828-7341
Fax: 1 (800) 993-7890

1.0 Warranty

The Control Unit is warranted free of defects in material and workmanship for a period of two (2) years.

The Mattress and Overlay are warranted free of defects in material and workmanship each for a period of one (1) year.

The Control Unit, Mattress and Overlay are warranted under the terms and conditions of the Gaymar warranty in place at the time of purchase. A copy of the warranty is available upon request. Gaymar disclaims all implied warranties including, but not limited to, the implied warranties of merchantability and of fitness for a particular purpose.

Control Units may be returned to the factory for servicing. Contact your local dealer to arrange servicing with the manufacturer.

Please contact Gaymar Customer Service if you have warranty questions.

Customers inside USA: 1 (716) 662-2551 Option 1

Toll Free (Customers inside USA): 1 (800) 828-7341

Customers outside USA: +1 (716) 662-8636

2.0 Symbols

 Attention, consult accompanying documents

 Type BF equipment

 Dangerous voltage

 Functional earth

 Class II Equipment (Double Insulated)

 Float/Static mode

 Alternating Low Pressure mode

 Loss of Pressure

 Audible/Visible Alarm On

 Audible/Visible Alarm Off

 Disposal by Separate Collection

 Latex free

3.0 Indications for Use

This device is for pressure redistribution to assist in the prevention and treatment of pressure ulcers.

4.0 Contraindications

Air support treatment is not recommended for patients with unstable spinal cords.

5.0 Safety Precautions

Review the following SAFETY PRECAUTIONS prior to using the Air Express system.

DANGER

- Risk of electric shock. Refer servicing to qualified service personnel.
- Explosion hazard. There is an explosion risk if used in the presence of flammable anesthetics. An explosion risk also exists when used with certain oxygen delivery equipment. The use of a nasal mask or ½ length oxygen tent should be used to deliver oxygen in the presence of this equipment. A full length oxygen tent that extends past the top surface of the mattress/overlay should NOT be used.

WARNING

- Disinfect the Air Express system between patient installations. Failure to disinfect may risk cross-contamination and infection.
- Check patient at least every 8 hours or once per shift, to assure proper system operation.
- Disconnect the AC power from the wall outlet before attempting to clean the Control Unit. Do not heat or steam autoclave any component of the system.

CAUTION

- Some medical conditions may not respond to treatment of this type. Patient's skin condition should be inspected regularly. Consult physician if any redness or skin breakdown occurs.
- For grounding reliability, plug only into a properly grounded outlet.
- Make certain all mattress/overlay straps are secured to the bed frame to prevent mattress/overlay from sliding and causing patient injury.
- The Air Express mattresses and overlays are not intended to be AND DO NOT FUNCTION AS a patient fall safety device. SIDE RAILS MUST BE USED WITH THE AIR EXPRESS MATTRESS AND OVERLAYS TO PREVENT PATIENT FALLS, unless determined unnecessary based on the facility protocol or the patient's medical needs as determined by the facility. IN SUCH CASES THE USE OF OTHER SUITABLE PATIENT SAFETY MEASURES ARE RECOMMENDED. When using side rails, the use of side rail pads is strongly recommended. Failure to use a side rail pad may result in patient injury.
- Use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the support surface will reduce the effectiveness of the system.
- Do not pull linens tightly over mattress/overlay. Tight sheets can cause "hammocking" and reduce effectiveness of treatment.
- Portable and RF communications equipment can affect medical electrical equipment.
- During mattress set-up, insure that all mattress air cells are inflated and that air cells are secured to the mattress base with the Velcro® type straps located on each side of the air cell. Air cells not properly secured or inflated may result in compromised therapy.
- Ensure that the Air Express mattress or overlay properly fits the bed frame on which it is being placed. Using the Air Express mattress/overlay on larger or smaller than the stated width of the bed frame may result in the mattress/overlay sliding and patient injury.

6.0 Description

The Air Express system is a portable air flotation system with alternating low pressure, designed to provide pressure redistribution to patients suffering from, or at risk of developing pressure ulcers. The system, consisting of Control Unit, Air Flotation Mattress and top cover or Overlay and top cover, is designed to provide pressure redistribution and patient comfort.

The Control Unit inflates and maintains the M280 Overlay or M1000 Mattress to user-selected settings, and provides alternating low pressure treatment if desired. Additionally, the C1000DF Control Unit provides a user selectable float/static mode.

Overlay

The M280 Overlay contains a nylon cell which provides patient support, and is enclosed in a low-shear, vapor permeable nylon cover.

Mattress

The M1000 Mattress consists of twenty (20) transverse air cells constructed of low shear nylon. Each of the cells is approximately 20cm (8") high when inflated. The air cells are positioned above a 5cm (2") convoluted foam base.

Covering the entire mattress assembly is a low-shear, vapor permeable, nylon top cover. The top cover surface creates a membrane impermeable to liquids though still permeable to water vapor. Spun bonded polyester fiberfill is quilted onto the base of the mattress top cover. When the patient moves, the top cover will move along with the patient, minimizing shearing effects. Additionally, the quilting acts to provide a diffusion layer by which the water vapor can be drawn away from the patient.

7.0 Features

7.1 C1000DF Control Unit

- Comfort control adjustment for customized treatment
- User selectable mode for Alternating Low Pressure Mode or Float/Static Mode
- CPR quick release for rapid deflation
- Durable plastic housing
- Power fail alarm (audible and visual)
- Loss of pressure alarm (audible and visual)
- Alarm ON/OFF switch
- Filter

7.2 M280 Overlay




- Single internal cell designed for maximum support and pressure redistribution
- Low shear top cover to aid in prevention of skin breakdown
- Dimensions - 86cm W x 203cm L x 13cm H (34" W x 80" L x 5" H) when fully inflated

7.3 M1000 Air Flotation Mattress

- Horizontal, tubular cell construction for optimal comfort and pressure redistribution
- Individual air cells designed for maximum support and pressure redistribution
- High density foam base for additional protection
- Low shear top cover to aid in prevention of skin breakdown
- Dimensions - 86cm W x 203cm L x 25cm H (34" W x 80" L x 10" H) when fully inflated


- Optional Foot Stop/ZIP™ cell to minimize heel interface pressure while maintaining appropriate positioning of feet. Contact your dealer for ordering information
- Removable cells for greater flexibility and easy maintenance

8.0 Control Panel

- ① **POWER ON/OFF SWITCH** on side of unit
- ② **LIGHTED MODE SELECT SWITCH** toggles between operating modes (Float/Static or Alternating Low Pressure)
 - **FLOAT/STATIC MODE (OOOOO)**—All cells within mattress/overlay will be uniformly inflated when this function is selected. Mode Select Switch light is OFF when this mode is activated.
 - **ALTERNATING LOW PRESSURE MODE (OOOO)**—Adjacent cells within mattress/overlay will be alternately inflated and deflated in this mode. Mode Select Switch light is ON when this mode is activated.
- ③ **COMFORT CONTROL DIAL** adjusts for patient comfort
 - **SOFT** ()—Minimum inflation level of the mattress/overlay
 - **FIRM** ()—Maximum inflation level of the mattress/overlay
- ④ **LOSS OF PRESSURE WARNING LIGHT (∞∞)**
- ⑤ **POWER ON LIGHT (I)**
- ⑥ **LIGHTED ALARM SWITCH** will activate ( Alarm Switch ON) or deactivate ( Alarm Switch OFF) all audible and visual alarms.
 - If power cord is disconnected with unit **ON** and alarm switch activated (), an **Alarm** will sound and flash intermittently. To avoid inadvertent Power Fail Alarm, turn the unit off before disconnecting power.
 - If there is a loss of mattress/overlay pressure with the unit **ON** and alarm switch activated (), an **Alarm** will sound and flash intermittently. In addition, the Loss of Pressure Warning Light () will be illuminated.

CAUTION

- The Audible Low Pressure Alarm is not active during initial mattress or overlay inflation. The Audible Low Pressure Alarm is active after approximately 50 minutes has elapsed from the time the unit has been turned on.
- Always leave the Alarm Switch in the **ON** position to activate the audible and visual Loss of Pressure Alarm and the Loss of Power Alarm.

C1000DF Control Panel

9.0 Setting Up the System

M280 Overlay

1. Place the overlay on the existing bed mattress with the connecting hoses at the foot of the bed.
2. Secure overlay to bed by looping the nylon straps on the overlay around the corners of the bed mattress. Adjust straps leaving enough slack so the corners of the mattress are not compressed when the bed frame is in a horizontal position. Make sure the side rails can still move up and down without interference.

M1000 Mattress

1. Place the mattress on the bed frame with the connecting hoses at the foot of the bed.
2. Secure mattress to bed frame by looping the nylon straps on mattress through the bed frame and then back through the D-rings on mattress. Make sure the side rails can still move up and down without interference.
3. Lower frame to its horizontal position and tighten straps.
4. Cover the mattress with the low-shear top cover and attach to the mattress by looping the Velcro® straps on the top cover through the D-rings on the mattress and tighten straps. Top cover should be loose to prevent "hammocking".

Control Unit / Linens / Patient

1. Suspend Control Unit on footboard of bed using the mounting brackets on back of unit. Control Unit can also be placed on the floor in an upright position.
2. Connect the dual quick-connect couplings of the mattress or overlay hoses to the mating connectors located on the side of the Control Unit.
3. Check that the hoses are secure by pulling gently. Ensure the hoses are not kinked or tucked under the mattress or overlay.
4. Plug the AC power cord into a properly grounded outlet.
5. Before turning on the Control Unit, insure that the Alarm ON/OFF Switch is activated by positioning the switch in the Alarm ON (⏻ – lighted) position. Activate the Mode Select Switch to either the Alternating Mode or Float/Static Mode.
6. Turn on the Control Unit. Turn Comfort Control Dial completely clockwise for quickest inflation. The mattress will begin to inflate. This will take approximately 50 minutes for complete inflation. After 50 minutes, turn the Comfort Control Dial completely counter-clockwise. Wait an additional 10 minutes, then perform the Hand Check procedure as outlined in Section 10.0.
7. Apply hospital linens and/or incontinence pad over the top of the mattress or overlay. Linens should be loose to prevent "hammocking".

CAUTION: Use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the support surface will reduce the effectiveness of the system.

8. Place patient on mattress and position along center of bed. Raise the bed side rails.

10.0 Instructions for Use

10.1 With M280 Overlay

Alternating Low Pressure Mode


In Alternating Low Pressure Mode, every other longitudinal section within the overlay will alternately increase and decrease pressure every 6 minutes. To select alternating low pressure, toggle the Mode Select Switch on the C1000DF Control Unit toward the Alternating Low Pressure symbol (OOOO). Adjust the inflation level for patient comfort using the Comfort Control Dial.

Float/Static Mode

In Float/Static Mode the overlay will be uniformly inflated. To select the float/static function, toggle the Mode Select Switch on the C1000DF Control Unit toward the Float/Static symbol (OOOOO). Adjust the inflation level for patient comfort using the Comfort Control Dial.

Hand Check for Overlay

Perform a hand check to ensure that the overlay is properly inflated. Place a flat hand, palm upward, under the overlay directly underneath the patient's pelvic region. If you can not feel the patient's body on your flat hand, no adjustment is needed. If you can feel the patient's body resting on your hand, increase the inflation level to the next highest setting, wait 10 minutes and perform the hand check again.

WARNING

Check patient at least every 8 hours or once per shift, to assure proper system inflation.

10.2 With M1000 Mattress

Alternating Low Pressure Mode


In Alternating Mode, every other cell within the mattress will alternately increase and decrease pressure every 6 minutes. To select alternating low pressure, toggle the Mode Select Switch on the C1000DF Control Unit toward the Alternating Low Pressure symbol (OOOO). Adjust the inflation level for patient comfort using the Comfort Control Dial.

Float/Static Mode

In Float/Static Mode each cell of the mattress will be uniformly inflated. To select the float/static function, toggle the Mode Select Switch on the C1000DF Control Unit toward the Float/Static symbol (OOOOO). Adjust the inflation level for patient comfort using the Comfort Control Dial.

Hand Check for Mattress

Perform a hand check to ensure that the mattress is properly inflated. Detach the center mattress cover strap and lift the mattress cover to access the air cells. Using a vertical hand, insert your hand between the air cells directly underneath the patient's pelvic region. Determine if four fingers worth of clearance is present between the mattress base and the patient. If clearance of four fingers exists, no adjustment is needed. If you can feel the patient's body resting on your hand, increase the inflation level to the next highest setting, wait 10 minutes and perform the hand check again.

WARNING

Check patient at least every 8 hours or once per shift, to assure proper system inflation.

11.0 CPR (Cardio Pulmonary Resuscitation)

To deflate mattress or overlay for CPR:

1. Disconnect the dual quick-connect mattress or overlay hoses from the Control Unit.
2. Begin CPR immediately. Do not wait for mattress or overlay to completely deflate.

12.0 System Cleaning, Care and Storage

12.1 Cleaning

WARNING

Disconnect the AC power cord from the wall outlet before attempting to clean the Control Unit. Do not heat or steam autoclave any component of the system.

Do not submerge the Control Unit.

1. To clean, use soap, water and a clean cloth to wipe down the Control Unit, power cord, hoses and mattress or overlay. Do not use abrasive cleaners on the mattress. Wipe dry with a clean, dry cloth. Note: Blood and other body fluids must be thoroughly cleaned from all surfaces before applying disinfectants.
2. Gaymar surfaces are designed for reliability and long-lasting comfort. In order to maintain these surfaces, they should be cleaned regularly using warm water and a hospital grade disinfectant. Gaymar recommends Steris® Wescodyne® for both cleaning and disinfecting.
3. Wipe down the mattress or overlay with a clean, dry cloth to remove any excess disinfectant.
4. Top Cover of mattress or overlay may also be laundered between patient uses or as required. Fill the washing machine with warm water (70 -140°F or 21 - 60°C). Add one cup of laundry detergent. Place no more than four top covers in a single extra large load capacity washing machine. When wash cycle is complete, remove promptly from machine and ensure all excess water is drained from load. Place dryer on LOWEST heat setting, or AIR FLUFF if available until dry. Verify top cover is completely dry before placing under patient.
5. If individual air cells of mattress become soiled, clean and disinfect as described above or simply replace air cell with a clean replacement. A plug is available from Gaymar to prevent water from entering air cells during laundering (P/N 30287).

12.2 Storage

Control Unit

1. Check the AC power cord and plug for abrasions or excessive wear.
2. Coil power cord up and place both cord and Control Unit in plastic bag for storage.

Mattress/Overlay

1. Check the air manifold of the mattress/overlay for kinks or breaks. Replace if necessary.
2. To deflate, disconnect hose assembly from Control Unit.
3. Coil the hose assembly loosely and roll the deflated mattress or overlay starting at the head end. When rolling up the mattress, use the straps for containment. Place in plastic bag for storage.

13.0 Troubleshooting

Symptom:	Corrective Action
Mattress or overlay does not inflate or is too soft.	Make certain the air hoses are not kinked, cracked, or split. Verify the hose couplings are fully engaged to Control Unit. Check that the air cells are not damaged. Verify that the Control Unit is ON and that air is coming from either connector outlet. Additionally on the mattress, verify that all 20 air cells are connected to the air manifold within the mattress.
Loss of Pressure Alarm is sounding.	
Audible Power Fail Alarm is sounding.	Power has been disconnected from Control Unit while in an operational mode. To restore power, reconnect the cord.

14.0 Specifications, Control Unit

Specification	C1000DF
Enclosure Dimensions	31cmW x 23cmH x 15cmD (12"W x 9"H x 6"D)
Weight	2.4kg (5.3 lb)
Power Cord	Detachable 4.6m (15'), 18 AWG minimum, with ground wire
Overcurrent Protection	Primary 1.0A, 250V, T
Input	230 VAC, 50Hz <37W
Operating Ambient Temperature Range	16 to 32°C (60 to 90°F)
Classification	Class II grounded equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
	Type BF equipment 
	Double insulated 
	IEC, EN 60601-1
	IPX0, enclosed equipment without protection against ingress of water.
	Continuous operation
Electromagnetic Compatibility	Meets EN60601-1-2:2001. (CISPR 11 Classified as Class B, Group 1 ISM equipment)

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©2007. Gaymar Industries, Inc.
Gaymar and Air Express are trademarks of Gaymar Industries, Inc.
STERIS and Wescodyne are registered trademarks of STERIS Corporation
Velcro is a registered trademark of Velcro Industries B.V., Inc.
Gaymar Industries, Inc.
Stryker France S.A.S.