Auto Sure-Float™

Floatation Pressure Redistribution with
Alternating Low Pressure and Active
Sensor Technology®

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REF M1001S/M1002S-seriens madrass

Operatörshandbok [SVENSKA]

* (See inside cover for applicable serial numbers)



Table Of Contents

<u>Section</u>	<u>Description</u>	<u>Page</u>
1.0	Warranty.....	2
2.0	Symbols.....	2
3.0	Indications for Use.....	2
4.0	Contraindications.....	2
5.0	Safety Precautions.....	3
6.0	Description.....	4
7.0	Features	5
8.0	Operator Control Panel.....	6 - 7
9.0	Alarm Conditions.....	8
10.0	Setting Up the System.....	8 - 9
11.0	Instructions for Use	9 - 11
12.0	CPR (Cardio Pulmonary Resuscitation).....	11
13.0	System Cleaning, Care and Storage.....	11 - 12
14.0	Troubleshooting.....	12 - 13
15.0	Specifications, Control Unit.....	14

Before you begin...

 **Important**

This manual version is applicable for the C1000MS (Configuration A) and C1000MES (Configuration B). The products can be identified by the serial number bar code on the back of each device. See serial number examples below for location and identification.

Configuration identified here

Serial Number Examples

Before using the Auto Sure-Float™ Mattress System, please read and understand this manual and all safety precautions prior to each application.

Only qualified medical service personnel should attempt to repair this device. A Service Manual is available for use by qualified personnel. For assistance contact your local dealer. If additional assistance is needed, contact Gaymar's technical service department.

1.0 Warranty

The Control Unit is warranted free of defects in material and workmanship for a period of two (2) years.

The Air Flotation Mattress is warranted free of defects in material and workmanship each for a period of one (1) year.

The Control Unit and Mattress are warranted under the terms and conditions of the Gaymar warranty in place at the time of purchase. A copy of the warranty is available upon request. Gaymar disclaims all implied warranties including, but not limited to, the implied warranties of merchantability and of fitness for a particular purpose.

Control Units may be returned to the factory for servicing. Contact your local dealer to arrange servicing with the manufacturer. Please contact Gaymar Customer Service if you have warranty questions:

Customers inside USA: 1 (716) 662-2551 Option 1
Toll Free (Customers inside USA): 1 (800) 828-7341
Customers outside USA: +1 (716) 662-8636

2.0 Symbols

Attention, Consult Accompanying Document.

Type BF Equipment

Dangerous Voltage

Protective Earth

Latex Free

Disposal by Separate Collection

Manual Mode

Auto Mode

Float / Static Mode

Alternating Low Pressure Mode

3.0 Indications for Use

This device is intended to provide pressure redistribution to aid in the prevention and treatment of pressure ulcers.

4.0 Contraindications

Air support treatment is not recommended for patients with unstable spinal cords.

5.0 Safety Precautions

Review the following **SAFETY PRECAUTIONS** prior to using the Auto Sure-Float™ System.

DANGER

- Risk of electric shock. Refer servicing to qualified service personnel.
- Explosion hazard. There is an explosion risk if used in the presence of flammable anesthetics. An explosion risk also exists when used with certain oxygen delivery equipment. The use of a nasal mask or ½ length oxygen tent should be used to deliver oxygen in the presence of this equipment. A full length oxygen tent that extends past the top surface of the mattress should NOT be used.

WARNING

- Disinfect the Auto Sure-Float™ system between patient installations. Failure to disinfect may risk cross-contamination and infection.
- Check patient at least every 8 hours or once per shift, to assure proper system operation.
- Disconnect the AC power from the wall outlet before attempting to clean the Control Unit. Do not heat or steam autoclave any component of the system.

CAUTION

- Some medical conditions may not respond to treatment of this type. Patient's skin condition should be inspected regularly. Consult physician if any redness or skin breakdown occurs.
- For grounding reliability, plug only into a properly grounded outlet.
- Make certain all mattress straps are secured to the bed frame to prevent mattress from sliding and causing patient injury.
- The mattresses are not intended to be AND DO NOT FUNCTION AS a patient fall safety device. SIDE RAILS MUST BE USED WITH THE MATTRESS TO PREVENT PATIENT FALLS, unless determined unnecessary based on the facility protocol or the patient's medical needs as determined by the facility. IN SUCH CASES THE USE OF OTHER SUITABLE PATIENT SAFETY MEASURES ARE RECOMMENDED.
- Use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the support surface will reduce the effectiveness of the system.
- Do not pull linens tightly over mattress. Tight sheets can cause "hammocking" and reduce effectiveness of treatment.
- Portable and RF communications equipment can affect medical electrical equipment.
- During mattress set-up, insure that all mattress air cells are inflated and that air cells are secured to the mattress base with the Velcro® type straps located on each side of the air cell. Air cells not properly secured or inflated may result in compromised therapy.
- Insure that the mattress properly fits the bed frame on which it is being placed. Using the mattress on larger or smaller than the stated width of the bed frame may result in the mattress sliding and patient injury.

6.0 Description

The Auto Sure-Float™ System is a portable air floatation system with on-demand alternating low pressure treatment, designed to provide benefit to patients suffering from, or at risk of, developing pressure ulcers. The system, consisting of a Control Unit, Air Floatation Mattress and Top Cover, is designed to provide pressure redistribution and patient comfort.

Control Unit (C1000MS)

The Auto Sure-Float™ Control Unit (C1000MS) inflates and maintains air floatation mattress (M1001S & M1002S) to user-selected settings, and can also provide on demand alternating low pressure treatment if desired. Featuring Active Sensor Technology®, the system automatically adjusts the mattress to the optimal fill level ensuring complete floatation pressure redistribution. Additionally, the Control Unit provides a user selectable Float/Static mode.

Air Floatation Mattress (M1001S & M1002S)

The Auto Sure-Float™ Mattress Series consists of twenty (20) transverse air cells constructed of low shear nylon. Each of the cells is approximately 20cm (8") high when inflated. The air cells are positioned above a 5cm (2") convoluted foam base (M1001S Series Only). The mattress utilizes 4 specialized air cells in the center region of the mattress which contain active sensors for detecting patient floatation level.

Covering the entire mattress assembly is a low-shear, vapor permeable, nylon top cover. The top cover surface creates a membrane impermeable to liquids though still permeable to water vapor. Spun bonded polyester fiberfill is quilted onto the base of the mattress top cover. When the patient moves, the top cover will move along with the patient, minimizing shearing effects. Additionally, the quilting acts to provide a diffusion layer by which the water vapor can be drawn away from the patient.

7.0 Features

Control Unit (C1000MS)

- Max Inflate for easy patient transfer.
- Soft/Firm comfort control adjustment for customized treatment.
- User selectable mode for Alternating Low Pressure mode or Float/Static mode
- Automatic patient floatation control via Active Sensor Technology® to ensure optimal pressure redistribution.
- CPR quick release providing for rapid deflation.
- Durable plastic housing.
- Mounting brackets for placement on the foot board of bed.
- Illuminated indicators on the control panel for easy viewing.
- Control panel Lockout mode to prevent inadvertent patient setting changes.
- Power fail alarm (audible and visual)
- Active Sensor disconnected in Auto mode alarm (audible and visual)
- Comfort level not achieved alarm (audible and visual)
- Alarm Silence feature to temporarily suspend audible alarm.

Air Floatation Mattress (M1001S & M1002S)

- Horizontal, tubular cell construction for optimal comfort and pressure redistribution.
- Individual air cells designed for maximum support and pressure redistribution.
- High density foam base for additional protection (M1001S Series Only).
- Low shear top cover to aid in prevention of skin breakdown.
- Dimensions - Base Model 34" W x 80" L x 10" H (86cm W x 203cm L x 25cm H) when fully inflated. Other sizes available. Contact your dealer for ordering information.
- Contains (4) air cells with active sensors for automatic patient floatation control.
- Optional Foot Stop/ZIP™ cell to minimize heel interface pressure while maintaining appropriate positioning of feet. Contact your dealer for ordering information.
- Removable cells for easy maintenance

8.0 Operator Control Panel

Figure 1

- a) **POWER ON/STANDBY BUTTON** switches the Control Unit between On and Standby modes.
- b) **POWER STANDBY MODE INDICATOR** indicates AC power is present to the Control Unit.
- c) **MODE SELECT BUTTON** Press repeatedly to toggle through operating modes, Auto-Float/Static, Auto-Alternating Low Pressure, Manual-Float/Static and Manual-Alternating Low Pressure.
- d) **AUTO MODE INDICATOR** (Sensor cord from M1001S/M1002S series mattress must be connected to control unit.) Mattress will automatically inflate to optimal level based on information received from active sensors in mattress. The indicator adjacent to the symbol will illuminate while Auto mode is active.

Note: Connecting the mattress sensor cord while the unit is not in Auto mode will cause the Control Unit to enter Auto mode immediately.

If the sensor cord from mattress becomes disconnected from control unit while in Auto mode, an audible alarm will activate and the Auto mode indicator will flash.
- e) **MANUAL MODE INDICATOR** Mattress inflation level user-determined by Comfort Adjust (-, +) setting. The indicator adjacent to the symbol will illuminate while Manual mode is active.
- f) **FLOAT/STATIC MODE INDICATOR** all cells within mattress will be uniformly inflated when the indicator adjacent to the symbol is illuminated.

- g) **ALTERNATING LOW PRESSURE MODE INDICATOR** Adjacent cells within mattress will be alternately inflated and deflated in this mode. The indicator adjacent to the symbol will illuminate while Alternating Low Pressure mode is active. Note: Alternating Low Pressure mode is temporarily suspended in Max Inflate mode.
- h) **SOFT (-) BUTTON** decreases the inflation level of the bed while Control Unit is in Manual mode only. Pressing once will decrease the setting by one comfort level.
- i) **COMFORT CONTROL INDICATOR** displays the current inflation level set point.
Note: An audible alarm will activate and the comfort level indicator(s) will flash if a set point not achieved alarm is activated.
- j) **FIRM (+) BUTTON** increases inflation level of the bed while in Manual mode only. Pressing once will increase the setting by one comfort level setting (~2.5mmHg).
- k) **MAX INFLATE BUTTON AND INDICATOR** selects/deselects the Max Inflate mode (maximum mattress firmness). The indicator adjacent to the button will illuminate while Max Inflate mode is active. After 30 minutes, Max Inflate mode will cancel.
- l) **LOCKOUT BUTTON AND INDICATOR** activates or cancels Lockout mode. Press and hold for three seconds to activate or cancel mode. The indicator adjacent to the button will illuminate while Lockout mode is active. In Lockout mode, only the Max Inflate button is operational. All other buttons are temporarily inactive.
- m) **ALARM SILENCE BUTTON** when this button is momentarily pressed it silences indefinitely any existing audible alarm.
- n) **ALARM INDICATOR** flashes indefinitely when an alarm exists. If no alarm conditions are present this indicator will extinguish.
- o) **POWER FAIL INDICATOR** flashes synchronously with an audible alarm when the power supplied from line cord is interrupted while the Control Unit is operating. To avoid an inadvertent power fail alarm, place the unit in Standby mode before disconnecting power from Control Unit.


9.0 Alarm Conditions - Unit will alarm for any of the following conditions

- If device is unplugged while in the On mode, an alarm will sound for a limited period of time. This alarm can be identified by the unit power fail indicator (Letter O in Figure 1) flashing synchronously with the audible alarm. This alarm cannot be silenced by the Alarm Silence button. To correct problem, re-attach power and put unit into Standby mode, then unplug.
- If Active Sensor cord is unplugged while in the Auto mode, an alarm will sound indefinitely. This alarm can be identified by the Auto mode indicator (Letter D in Figure 1) toggling with the alarm indicator and audible alarm. This alarm can be silenced by the Alarm Silence button (Letter M in Figure 1) but the alarm indicators will continue to flash until the problem is corrected. To correct problem, place unit in Manual mode or re-attach Active Sensor cord.
- If device has been running for 30 minutes and proper inflation level cannot be attained (filling and/or venting continuously for more than 30 minutes), an alarm will sound indefinitely. This alarm can be identified by the Comfort Control Indicator(s) (Letter I in Figure 1) toggling with the alarm indicator and audible alarm. This alarm can be silenced by the Alarm Silence button (Letter M in Figure 1) but the alarm indicators will continue to flash until the problem is corrected. To correct problem, check for disconnected or damaged air cells or disconnected hose.

10.0 Setting Up the System

- Place the mattress on the bed frame with the connecting hoses at the foot of the bed.
- Before securing the mattress straps, raise the head and foot ends of the bed frame. Raise the side rails. Loosely, secure mattress to bed frame by looping the nylon straps on mattress through the bed frame and then back through the D-rings on mattress. Make sure the side rails can still move up and down without interference.
- Lower frame to its horizontal position and tighten straps.
- Zip the mattress with the low-shear quilted cover and verify the cover is loose fitting to prevent “hammocking”.
- Suspend Control Unit on footboard of bed using the mounting brackets on back of unit.
- Connect the dual quick-connect couplings of the mattress hoses to the mating connectors located on the side of the Control Unit.
- Check that the hoses are secure by pulling gently. Ensure the hoses are not kinked or tucked under the mattress.
- Connect the plug end of the mattress sensor cord into the jack on the side of the Control Unit. Plug the AC power cord into a properly grounded outlet. The Control Unit will enter the last mode it was in when power was disconnected.
- Press the On/Standby button on the lower right of the control panel. The mattress will begin to inflate. This could take approximately 30 to 60 minutes for complete inflation. After inflation is complete perform the Hand Check.

- Hand Check, to perform a hand check to ensure that the mattress is properly inflated. Unzip and lift the mattress cover to access the air cells. Using a vertical hand, insert your hand between the air cells directly underneath the patient's pelvic region. Determine if four fingers worth of clearance is present between the mattress base and the patient. If clearance of four fingers exists, no adjustment is needed. If you can feel the patient's body resting on your hand, increase the inflation level to the next highest setting, wait 10 minutes and perform the hand check again.
- Apply hospital linens and/or incontinence pad over the top of the mattress. Linens should be loose to prevent "hammocking". CAUTION: Use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the support surface will reduce the effectiveness of the system.
- Place patient on the mattress and position along center of bed. Raise the bed side rails.

WARNING

Check patient at least every 8 hours or once per shift, to assure proper system inflation.

11.0 Instructions for Use

Automatic Mode

- When the Mattress System is in the Automatic Mode the mattress will automatically inflate to optimal levels based on information received from active sensors air cells in the mattress and in conjunction with the Control Unit.
- The Automatic Mode of operation is automatically engaged at the press of the 'On' button (when the M1001S or M1002S series Air Mattress Sensor Cord is plugged into the Control Unit). The indicator adjacent to the Auto symbol will illuminate indicating the Auto mode is active.
- Connecting the Air Mattress Sensor Cord to the Control Unit while the system is not in Auto mode will engage the Auto mode automatically and illuminate the indicator adjacent to the Auto symbol.

Theory

- While in Automatic mode, the Control Unit continually repeats a deflation/inflation cycle. The inflation level of the mattress is lowered until the lowest allowable floatation height is sensed (via feedback from the active sensors in the mattress). It is then slowly increased until the appropriate patient floatation height is sensed.

- If the patient immerses into the mattress for any reason, (such as elevation of the head section) the Control Unit will automatically increase the inflation level until the appropriate floatation height is again realized.
- Disconnecting the mattress sensor cord from the Control Unit while in Auto mode will cause an audible alarm to sound and indicators adjacent to the Alarm Silence and Auto symbols to flash. This alarm can be indefinitely silenced by pressing the Alarm Silence button. The alarm will be reset if the mattress sensor cord is reconnected to the jack on the right side of the unit, or by pressing the Mode Select button to disengage the Auto mode.
- To prevent accidental changes in settings, press and hold the Lockout button for approximately 3 seconds.

Automatic Float / Static Mode

- During the Automatic Float/Static Mode all of the air cells of the mattress will be uniformly inflated.
- To select Float/Static, press the Mode Select button until the indicators adjacent to both the Auto and Float/Static symbols are illuminated.
- The Automatic Float/Static Mode requires no user input to determine the appropriate inflation level for the patient. The active sensors in the mattress located under the patient's pelvic region, in conjunction with the Control Unit automatically determine the optimal inflation level for the patient to ensure complete floatation.

Automatic Mode - Alternating Low Pressure

- In Auto mode the mattress will automatically inflate to optimal level based on information received from active sensors in mattress. In Alternating Low Pressure mode, every other longitudinal section within the mattress will alternately increase and decrease pressure every 8½ minutes.
- If alternating low pressure therapy is also desired, press the Mode Select button until the indicators adjacent to both the Auto and Alternating Low Pressure symbols are illuminated.
- In Auto – Alternating Low Pressure mode, alternate cells within the mattress are inflated and deflated every 8½ minutes. At the end of each 4¼ minute cycle the system will evaluate the floatation level, (determined by the automatic sensors within the mattress) and make comfort level adjustments as necessary before resuming the Alternating Low Pressure operation.

Manual Mode – Alternating Low Pressure

- In Alternating Low Pressure mode, every other latitudinal section within the mattress will alternately increase and decrease pressure every 8½ minutes.
- If alternating low pressure redistribution is desired, press the Mode Select button until the indicators adjacent to both the Manual and Alternating Low Pressure symbols are illuminated. If

Alternating Low Pressure is selected in combination with Manual mode, every other cell within the mattress will alternately increase and decrease pressure every 8½ minutes.

Transferring Patients In and Out of Bed

For patient transfer press the Max Inflate button. This will maximally inflate the mattress to provide a firm and stable support surface (and suspend Alternating Low Pressure and Auto modes). Pressing the Max Inflate button again will cancel this mode and return the mattress to previously chosen settings. While in the Max Inflate mode, the indicator adjacent to the Max Inflate button will illuminate to remind the caregiver that the system is in this state. After approximately 30 minutes, the Control Unit will automatically cancel the Max Inflate mode and return to the previously chosen settings.

12.0 CPR (Cardio Pulmonary Resuscitation)

- Disconnect the dual quick-connect mattress hoses from the Control Unit.
- Begin CPR immediately. Do not wait for mattress to completely deflate.

13.0 System Cleaning, Care and Storage

Cleaning

WARNING
<ul style="list-style-type: none">• Disconnect the AC power cord from the wall outlet before attempting to clean the Control Unit. Do not heat or steam autoclave any component of the system.• Do not submerge control unit.

- To clean, dampen a clean cloth with soap and water and wipe down the Control Unit, power cord, hoses and mattress. Do not use abrasive cleaners on the mattress. Wipe dry with a clean, dry cloth. **Note: Blood and other body fluids must be thoroughly cleaned from all surfaces before applying disinfectants.**
- Apply hospital grade disinfectant to the external surfaces of the Control Unit, hoses and mattress. Follow manufacturer’s instructions. The solution contact time is what makes disinfection effective. **Note: Gaymar surfaces are designed for long-lasting comfort and reliability. In order to maintain these surfaces, they should be cleaned regularly using warm water and Gaymar recommended STERIS® Wescodyne® for both cleaning and disinfecting. If this solution is not available in your area, please contact Gaymar Technical Services for a listing of other approved cleaning and disinfecting agents.**
- Wipe down the mattress with a clean, dry cloth to remove any excess disinfectant.
- Top Cover of mattress may also be laundered between patient uses or as required. Fill the washing machine with warm water (70 -176°F or 21 - 80°C). Add one cup of laundry detergent.

Place no more than four top covers in a single extra large load capacity washing machine. When wash cycle is complete, remove promptly from machine and ensure all excess water is drained from load. Place dryer on LOWEST heat setting, or AIR FLUFF if available until dry. Verify top cover is completely dry before placing under patient.

- If individual air cells of mattress become soiled, clean and disinfect as described above or simply replace air cell with a clean replacement. A plug is available from Gaymar to prevent water from entering air cells during laundering (P/N 30287).
- The Plexus mattress cover should be inspected for tears, punctures, and excessive wear each time the cover is cleaned. Compromised covers should be discarded and replaced with a new top cover and the internal components of the mattress must be cleaned and disinfected using warm water and a hospital grade disinfectant before being placed back in service.

WARNING
Do not launder the Claret colored Active Sensor Technology® air cells of the M1001S /M1002S series mattress. Damage may result. If one or more of the Active Sensor Technology® air cells become damaged, they can only be replaced as a complete set.

Control Unit Storage

- Check the AC power cord and plug for abrasions or excessive wear.
- Coil power cord up and place both cord and Control Unit in plastic bag for storage.


Air Floatation Mattress Storage

- Check the air manifold of the mattress for kinks or breaks. Replace if necessary.
- To deflate, disconnect hose assembly from Control Unit.
- Coil the hose assembly loosely and roll the deflated mattress starting at the head end. When rolling up the mattress, use the straps for containment. Place in plastic bag for storage.


14.0 Troubleshooting

Symptom:	Corrective Action
Mattress does not inflate or is too soft.	<ul style="list-style-type: none"> • Verify the air hoses are not kinked, cracked, or split. • Verify the hose couplings are fully engaged to the Control Unit. • Verify the air cells are not damaged. • Verify that the Control Unit is ON and that air is coming from either connector outlet. • Verify that all 20 air cells are connected to the air manifold within the mattress.
Other.	See the following Basic Operator Trouble Shooting flow chart.

Basic Operator Trouble Shooting

Unit Is Alarming

15.0 Specifications, Control Unit

Specification	C1000MS
Enclosure Dimensions	10" W x 12" H x 6" D (26cm W x 31cm H x 16cm D)
Weight	< 10 lbs (4.6 kg)
Power Cord	Detachable 14 ft (4.25m) , #18 AWG minimum, with ground wire
Over Current Protection	Primary 2.5A, 250V, T
Input	115/230VAC, 60/50Hz, 115V < 0.44A, 230V < 0.33A
Operating Ambient Temperature Range	60 to 90°F (16 to 32°C)
Classification	Class I grounded equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
	Type BF equipment 
	UL2601-1, CSA C22.2 NO. 601.1 EN 60601-1
	IPX0, enclosed equipment without protection against ingress of water.
	Continuous operation
Electromagnetic Compatibility	Meets EN60601-1-2:2001 (CISPR 11 Classified as Class B, Group 1 ISM equipment)

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©2009. Gaymar Industries, Inc.
STERIS and Wescodyneare registered trademarks of STERIS Corporation
Active Sensor Technology is a registered trademark of Gaymar Industries, Inc.
Velcro is a registered trademark of Velcro Industries B.B., Inc.
Gaymar and Auto Sure-Float are trademarks of Gaymar Industries, Inc.
Gaymar Industries, Inc.