**REF** P102 Plexus<sup>®</sup> Jr. Alternating Pressure System

**REF** C102 Plexus<sup>®</sup> Jr. APM Control Unit

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
English  
 Gaymar P102  
 Alternating Pressure System  
 Instructions for Use  
 Figure 1

**CONTRAINDICATIONS**

Air support therapy is not recommended when spinal stability is a concern.

**INDICATIONS**

This device (P102) is intended to provide pressure redistribution to aid in the prevention and treatment of pressure ulcers. A system consists of the Plexus® Jr. control unit (C102) combined with Plexus® overlay M100. This overlay contains two sets of air cells that alternately inflate and deflate. The air cells redistribute the patient's weight over the surface and effectively reduce tissue interface pressure.

	<b>CAUTION</b>
<ul style="list-style-type: none"> <li>• Check patient's skin regularly. Consult physician if any redness or skin break occurs. Serious injury could result if the patient's skin condition is left untreated.</li> <li>• Do not place the control unit in the patient's bed, in contact with the patient, or under sheets or other coverings. This could cause serious injury or could affect control unit performance.</li> <li>• Do not use in the presence of a flammable anesthetic mixture or with oxygen (O<sub>2</sub>) or nitrous oxide (N<sub>2</sub>O).</li> <li>• Verify bed overlay is in good condition.</li> <li>• Verify bed side rails are compatible with bed frame and existing mattress.</li> <li>• Use with appropriate top sheet and minimize layers of bedding between patient and overlay.</li> </ul>	

**CONTROL PANEL**


1. Power ON/OFF Switch on side of unit.
2. COMFORT CONTROL DIAL adjusts for patient comfort.

- Soft ( )—Minimum inflation level of the overlay
- Firm ( )—Maximum inflation level of the overlay.

**INSTRUCTIONS**

1. Place control unit on flat surface or suspend control unit on end of bed using attached hooks. See Figure 1. Position overlay on existing mattress. Tuck end flaps of overlay under ends of bed mattress.
2. Connect the hose assembly between the overlay and the control unit.
3. Plug the power cord into a properly grounded outlet. Adjust pressure control knob to highest setting for quick inflation and turn control unit on using green on/off switch. Unit will take approximately 40 minutes to inflate the overlay.
4. Position patient on overlay. Adjust pressure control knob for patient comfort.
5. Perform Hand Check.  
 A Hand Check must be performed every 8 hours to verify proper operation.  
 To perform the Hand Check:  
 With the patient on his or her back, slide hand, flat and palm up, between the overlay and the mattress. Hand should be directly under the air cell that is under the patient's buttocks (or other bony area). See Figure 2.


Figure 2

Wait for full inflation of the air cell directly above hand. If the patient's body is not in direct contact with hand, the system is operating correctly. If during full inflation of air cell, the patient's body is in direct contact with flat hand, the system is not operating properly. Adjust the pressure control to a higher setting. Wait 10 minutes and repeat the Hand Check. If the Hand Check fails, check that the hoses are not kinked or pinched. If repeated Hand Check fails and hoses are not kinked, contact Dealer or Gaymar for further instruction.

	<b>WARNING</b>
<b>Deflate before CPR or CPR will be ineffective.</b>	
To deflate overlay for CPR: Disconnect the hoses from the control unit. See Figure 3. The overlay will deflate in approximately 20 seconds. Proceed with CPR procedures.	

Figure 3

**CLEANING**

The control unit housing, tubing and overlay should be cleaned between patients.

1. To clean, use soap, water and a clean cloth to wipe down the Control Unit, power cord, hoses and overlay. Do not use abrasive cleaners on the overlay. Wipe dry with a clean, dry cloth. Note: Blood and other body fluids must be thoroughly cleaned from all surfaces before applying disinfectants.
2. Apply disinfectant to the external surfaces of the Control Unit, hoses and overlay. Gaymar recommends these disinfecting agents: Phenolic Disinfectant Cleaner or a 10% Chlorinated Bleach Solution. Allow contact time according to manufacturers instructions. The solution contact time is what makes disinfection effective.
3. Wipe down the overlay with a clean, dry cloth to remove any excess disinfectant.
4. Top sheets of overlay may also be laundered between patient uses or as required to maintain good patient hygiene. Fill the washing machine with warm water (21 - 60 °C [70 - 140 °F]). Add one cup of laundry detergent. Place no more than four top sheets in a single extra large load capacity washing machine. When wash cycle is complete, remove promptly from machine and ensure all excess water is drained from load. Place dryer on LOWEST heat setting, or AIR FLUFF if available until dry. Verify top sheet is completely dry before placing under patient.

**TROUBLE SHOOTING**

Problem: Air cells fail to inflate.

Solution: Make certain the air hose is not kinked, cracked or split. Verify the power switch is illuminated, signifying the control unit has power. Verify the air hoses are fully inserted with a positive connection.

Problem: Loss of power.

Solution: Make certain the power receptacle has power. The automatic reset circuit breaker may have tripped. Wait five (5) minutes and check power again.

## TECHNICAL DATA

### C102

Dimensions: 23cm x 10cm x 8cm (9" x 4" x 3")

Weight: 1.22kg (2.7 lb.)

Power Cord: Hospital grade, 457cm (15 ft.)

Voltage: 230 VAC

Current: .05A

Wattage: 10 Watts

Over-current Protection: Fused, Rated 1A 250/125 VAC

Frequency: 50 Hz

Operating Temperature: 15°C to 32°C (60°F to 90°F)

Storage Temperature: -15°C to 50°C (5°F to 122°F)

Classification:

EN 60601-1, EN60601-1-2

IPX0, enclosed equipment without protection against ingress of water

Not suitable for use in the presence of a flammable anesthetic mixture or with oxygen (O<sub>2</sub>) or nitrous oxide (N<sub>2</sub>O).

Continuous operation

The control unit is warranted for one (1) year.

### M100

Dimensions: 9cm x 86cm x 203cm (3.5" to 6" x 34" x 80")

Operating Temperature: 15°C to 32°C (60°F to 90°F)

Storage Temperature: -15°C to 50°C (-5°F to 122°F)

Maximum Weight: 148kg (300 lb.) There is no minimum weight for use.

NOTE: If the Plexus® Jr. APM (P102) support surface is stored at low temperatures, it is important that it be brought to room temperature before unfolding. This prevents accidental cracking and splitting of the polyvinyl material.

## SYMBOLS

Caution, consult accompanying documents

Dangerous voltage

Type BF equipment

Consult Operating Instructions

Double Insulation

Disposal: Contact local distributor who will take the necessary steps according to your national market.

## CUSTOMER SERVICE

If you have purchased this control unit through a local dealer and are experiencing problems, please contact that local dealer for assistance. For additional assistance, contact our Customer Service Department at:

Gaymar Industries, Inc.  
10 Centre Drive  
Orchard Park, NY 14127-2295 USA  
Telephone: +1 716 662-8636  
FAX: +1 716 662-0730

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©2009. Gaymar Industries, Inc.
Plexus is a registered trademark of Gaymar Industries, Inc.
Stryker France S.A.S.