
















ProCeed 8002 Support Surface

Operations/Maintenance Manual

| | |
|-----|------------|
| REF | 8002-0-100 |
| REF | 8002-0-101 |
| REF | 8002-0-102 |
| REF | 8002-0-103 |
| REF | 8002-0-104 |
| REF | 8002-0-105 |
| REF | 8002-0-106 |
| REF | 8002-0-107 |
| REF | 8002-0-108 |
| REF | 8002-0-109 |



Symbols

| | |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Refer to instruction manual/booklet |
|  | Consult instructions for use |
|  | General warning |
|  | <p>Conforms to BS 7177:2008 for medium hazard</p> <p>Note - The BS 7177 label shown on the product is blue.</p> <p>REF 8002-0-104</p> <p>REF 8002-0-105</p> |
|  | Caution |
|  | Catalogue number |
|  | Serial number |
|  | Model number |
|  | Manufacturer |
|  | Date of manufacture |
|  | CE mark |
|  | UK Conformity Assessment Mark |
|  | Importer |
|  | Unique device identifier |
|  | Authorized representative in the European Community |









| | |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
| CH REP | Authorized representative in Switzerland |
| MD | European medical device |
|  | Maximum patient weight |
|  | Type B applied part |
|  | Wash by hand |
|  | Do not tumble dry |
|  | Do not dry clean |
|  | Do not iron |
|  | Air dry |
|  | Chlorinated bleach |
| US Patents | For US Patents see www.stryker.com/patents |

Table of Contents

| | |
|-----------------------------------------------------------------|----|
| Warning/Caution/Note Definition | 2 |
| Summary of safety precautions | 2 |
| Introduction | 3 |
| Product description | 3 |
| Intended use | 3 |
| Indications for use | 3 |
| Intended users | 3 |
| Clinical benefits | 3 |
| Contraindications | 4 |
| Expected service life | 4 |
| Disposal/recycle | 4 |
| Specifications | 4 |
| Contact information | 5 |
| Serial number location | 6 |
| Setup | 7 |
| Operation | 8 |
| Protecting the support surface cover from X-ray cassettes | 8 |
| Managing incontinence and drainage | 8 |
| Preventive maintenance | 9 |
| Cleaning | 10 |
| Disinfecting | 11 |

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not load the support surface above the safe working load of 225 kg.
 - Always use the support surface with the compatible bed frame as indicated in the specifications section of this manual.
 - Always secure the support surface between the support surface retainers on the bed frame.
 - Do not keep a concentrated load on the support surface for an extended period.
 - Always use the support surface in appropriate environmental conditions.
 - Do not use this product without the support surface cover.
 - Do not allow sharp edges from the X-ray cassette to come into contact with the support surface cover. If the cover of the support surface is damaged, remove the support surface from use and replace the product to prevent cross-contamination.
 - Follow your hospital protocols and complete preventive maintenance each time you clean the support surface cover. If the cover of the support surface is compromised, remove the support surface from use and replace the product to prevent cross-contamination.
 - Always wipe the product with clean water and dry after you clean or disinfect the product. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not use cleaning agents and disinfectants with aggressive chemicals as they will reduce the expected service life of the support surface cover.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service.

Product description

The Stryker **ProCeed 8002** support surface is a non-powered support surface that is used in combination with a hospital bed. This product is available in standard and **Visco** foam options and contains a flame barrier option.

Intended use

The Stryker **ProCeed 8002** support surface is intended to assist in the redistribution of pressure specific to the individual needs of the patient when implemented as a part of a thorough pressure ulcer prevention and risk-based patient management program. The product is intended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional.

Indications for use

The Stryker **ProCeed 8002** support surface is indicated for use with human patients with existing or at risk of incurring pressure ulcers in a healthcare facility.

Intended users

Operators of this product include healthcare professionals (such as nurses, nurses aides, and doctors).

Clinical benefits

The Stryker **ProCeed 8002** support surface enables the healthcare professional to perform pressure ulcer and patient management programs in accordance with the intended use and current standards of care.

Contraindications

None known.

Expected service life

The Stryker ProCeed 8002 support surface has a one year expected life under normal use conditions and with appropriate periodic maintenance.


Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

CAUTION - Do not load the support surface above the safe working load of 225 kg.

See the support surface label for applicable flammability standards.







| | | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|-----------------------------------------------------|---------------------------------------------------------------------|
| Safe working load  | 225 kg | | |
| Note - Do not exceed the maximum patient weight of the bed frame. See the bed frame operations/maintenance manual for product specifications. | | | |
| Standard support surface | 8002-0-100 | 8002-0-102 | 8002-0-104 |
| Length | 200 cm | 200 cm | 200 cm |
| Width | 87 cm | 87 cm | 87 cm |
| Thickness | 14 cm | 14 cm | 14 cm |
| Weight | 8.5 kg | 9.1 kg | 9.1 kg |
| Flammability standards | N/A | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175 | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175, BS7177 capable |
| Visco elastic support surface | 8002-0-101 | 8002-0-103 | 8002-0-105 |
| Length | 200 cm | 200 cm | 200 cm |
| Width | 87 cm | 87 cm | 87 cm |
| Thickness | 14 cm | 14 cm | 14 cm |
| Weight | 11.0 kg | 11.6 kg | 11.6 kg |
| Flammability standards | N/A | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175 | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175, BS7177 capable |

| Standard foam bed extender support surface | 8002-0-106 | 8002-0-107 | 8002-0-108 | 8002-0-109 |
|--------------------------------------------|------------|-----------------------------------------------------|------------|-----------------------------------------------------|
| Length | 33 cm | 33 cm | 33 cm | 33 cm |
| Width | 71 cm | 71 cm | 71 cm | 71 cm |
| Thickness | 14 cm | 14 cm | 20 cm | 20 cm |
| Weight | 1.3 kg | 1.5 kg | 1.6 kg | 1.8 kg |
| Flammability standards | N/A | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175 | N/A | BS EN 597-1, BS EN 597-2, BS EN 6807:2006, UNI 9175 |

| | | |
|---------------------|------------------------------|---------------------------------------------------|
| Compatible products | ProCeed 8002 support surface | ProCeed 8002 foam bed extender support surface |
| | ProCeed Hospital Bed | ProCeed Hospital Bed with the bed extender option |

Stryker reserves the right to change specifications without notice.

Specifications listed are approximate and may vary slightly from product to product.

| Environmental conditions | Operation | Storage and transportation |
|------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Ambient temperature |  |  |
| Relative humidity (non-condensing) |  |  |
| Atmospheric pressure |  |  |

Contact information

Contact Stryker Customer Service or Technical Support at: +1 800-327-0770.

Stryker Medical International

Kayseri Serbest Bölge Şubesi

2. Cad. No:17 38070

Kayseri, Turkey

Email: infosmi@stryker.com

Phone: + 90 (352) 321 43 00 (pbx)

Fax: + 90 (352) 321 43 03

Web: www.stryker.com

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling your Stryker Customer Service. Include the serial number in all written communication.

Serial number location

The serial number is located inside the support surface cover. Unzip the support surface cover to read the product label.




Figure 1 – Serial number location

Setup

CAUTION

- Always use the support surface with the compatible bed frame as indicated in the specifications section of this manual.
 - Always secure the support surface between the support surface retainers on the bed frame.
-

To set up the support surface:

1. Place the support surface on the compatible bed frame so that the Stryker logo is facing up and is at the head end of the bed.
2. Align the support surface cover so the watershed cover barrier points down (Figure 2).



Figure 2 – Watershed cover barrier

3. Make sure that the support surface is between the support surface retainers on each side of the bed frame (Figure 3).



Figure 3 – Support surface retainer

4. Pull the watershed cover barrier over the zipper.
5. Place linens on the support surface before patient use following your hospital protocols.

To set up the bed extender for the bed extender support surface, see the **ProCeed Hospital Bed Operations/Maintenance Manual**.

1. Center the bed extender support surface with the support surface on the bed frame.

Operation

CAUTION

- Do not keep a concentrated load on the support surface for an extended period.
 - Always use the support surface in appropriate environmental conditions.
-

Protecting the support surface cover from X-ray cassettes

CAUTION

- Do not use this product without the support surface cover.
 - Do not allow sharp edges from the X-ray cassette to come into contact with the support surface cover. If the cover of the support surface is damaged, remove the support surface from use and replace the product to prevent cross-contamination.
-

To protect the support surface from an X-ray cassette:

1. If you use the support surface on a bed with the X-ray cassette holder option, use the bed's X-ray cassette holder. See the bed's operations/maintenance manual for instructions.
2. If you use the support surface on a bed without the X-ray cassette holder option, place the X-ray cassette into a pillowcase or other device to cover the sharp edges.
3. Complete the X-ray cassette procedure following hospital protocols.
4. Remove the X-ray cassette from the support surface.

Managing incontinence and drainage

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

Preventive maintenance

Remove the product from use before you perform preventive maintenance.

CAUTION - Follow your hospital protocols and complete preventive maintenance each time you clean the support surface cover. If the cover of the support surface is compromised, remove the support surface from use and replace the product to prevent cross-contamination.

Note - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

- _____ Zipper and support surface cover (top and bottom) for tears, cuts, holes, or other openings
- _____ Internal components for signs of staining from fluid ingress or contamination by fully unzipping the support surface cover
- _____ Foam and other components for degradation or other damage

| |
|------------------------|
| Product serial number: |
| Completed by: |
| Date: |

Cleaning

CAUTION

- Always wipe the product with clean water and dry after you clean or disinfect the product. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not use cleaning agents and disinfectants with aggressive chemicals as they will reduce the expected service life of the support surface cover.
-

Always follow hospital protocols for cleaning and disinfecting.

To clean the support surface cover:

1. Using a clean, soft, damp cloth, wipe the support surface cover with a mild soap and water solution to remove foreign material.
2. Using a clean, soft cloth and water, wipe the support surface cover to remove excess cleaning agent.
3. Using a clean, soft, dry cloth, wipe the support surface cover to remove excess liquid or cleaning agent.
4. Allow the support surface cover to dry before you return to service.

Disinfecting

CAUTION

- Always wipe the product with clean water and dry after you clean or disinfect the product. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not use cleaning agents and disinfectants with aggressive chemicals as they will reduce the expected service life of the support surface cover.
-

Suggested disinfectants:

- Quaternaries without glycol ethers
- Phenolic cleaners (active ingredient - o-phenylphenol)
- Chlorinated bleach solution (5.25% - less than 1 part bleach to 100 parts water)
- 70% Isopropyl alcohol

Always follow hospital protocols for cleaning and disinfecting.

To disinfect the support surface cover:

1. Clean and dry the support surface cover before you disinfect. See *Cleaning* (page 10).
2. Using pre-soaked wipes or a damp cloth, apply the recommended disinfectant solution to the support surface cover. Do not soak the support surface.

Note - Make sure that you follow the manufacturer's instructions for appropriate contact time and rinse requirements.

3. Using a clean, dry cloth, wipe the support surface cover to remove excess liquid or disinfectant.
4. Allow the support surface cover to dry before you return to service.



Stryker EMEA Supply Chain Services B. V.
Frans Maasweg 2
Venlo 5928 SB, The Netherlands



Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **ProCeed, Stryker**. All other trademarks are trademarks of their respective owners or holders.



Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA