RELEASED
CR-F0120
Effective Date:
9/14/15

ProForm™
Non-Powered Support Surface
Model 2710

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ROMANIA

Stryker Osteonics Romania S.R.L.
19, Leonida Str.
District 2
7000 Bucharest
Romania
Phone: 40-2-12-12-11-22
Fax: 40-2-12-12-11-33

SOUTH AFRICA

Stryker Osteonics PTY. LTD.
3 Susan Street
Strydom Park
Johannesburg, 2194 - South Africa
Mailing Address:
P.O. Box 48039
2129 Roosevelt Park
Johannesburg, South Africa
Phone: 27-11-791-4644
Fax: 27-11-791-4696

SPAIN

Stryker Iberia SL
c/Sepulveda n17
28108 Alcobendas
Madrid
Spain
Phone: 34-91-728-35-00
Phone: 34-91-358-20-44
Fax: 34-91-358-07-48

SWEDEN


Stryker AB/Scandinavia
Box 50425
SE-204 14 Malmo
Sweden
Phone: 46-40-691-81-00
Fax: 46-40-691-81-91

Table of Contents

	<u>page</u>
Symbols and Definitions	1-2
Symbols.....	1-2
Warning / Caution / Note.....	1-3
Introduction	1-4
Intended Use of the Product	1-4
Product Description.....	1-4
Specifications	1-5
Contact Information.....	1-6
Serial Number Location.....	1-6
Summary of Safety Precautions.....	1-7
Operation Guide.....	1-9
Transfer of Patient from One Patient Support Platform to Another	1-9
Using Incontinence Management Devices.....	1-10
Cleaning	1-11
Cleaning Support Surface.....	1-11
Disinfecting Support Surface	1-12
Service Information.....	1-13
Cover Replacement	1-13
Preventive Maintenance	1-14
Checklist.....	1-14
Checking for Compression Set	1-15
Product Labeling.....	1-16
Quick Reference Replacement Parts List.....	1-18
Mattress Assembly with Cover	1-19
Mattress Assembly with Cover & Fire Barrier.....	1-20
Warranty.....	1-21
Limited Warranty	1-21
To Obtain Parts and Service.....	1-21
Return Authorization	1-21
Damaged Merchandise	1-21

Symbols and Definitions

SYMBOLS

	Warning/Caution - Consult accompanying documentation
	Safe Working Load
	Do Not Iron
	Damp Wipe Only
	Chlorinated Bleach
	Do Not Tumble Dry
	Do Not Dry Clean
	Allow to Completely Air Dry
	Consult Instructions for Use
	Manufacturer
	Product Weight

WARNING / CAUTION / NOTE

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

NOTE!

Provides special information to make specific instructions easier to understand.

Introduction

This manual is designed to assist with the operation and maintenance of the ProForm Non-Powered Support Surface. Carefully read this manual thoroughly before using or beginning maintenance on the support surface. To ensure safe operation of this equipment, it is recommended that methods and procedures are established for educating and training staff on the safe operation of the support surface.

INTENDED USE OF THE PRODUCT

ProForm is a non-powered support surface intended to assist in improving human patient outcomes, focusing on equalizing pressure redistribution and enhancing comfort. ProForm support surface is used to assist in the prevention and treatment of all categories/stages of pressure ulcers (including stages I, II, III, IV, Unstageable, and Deep Tissue Injury) and is recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a health care professional.

This product is intended to be used:

- With patients in a general hospital environment and who are at risk for developing pressure ulcers, as well as those who require therapy for pre-existing pressure ulcers
- In acute care, general hospital care, or other locations as prescribed by a physician
- On Stryker and other flat-deck frames (See "**Specifications**" section on page 7)
- For multiple patients after appropriate cleaning/disinfection, i.e., reusable; This product is not disposable

This product is not intended to be used in a home health care environment.

The safe working load for ProForm is 500 pounds; the patient must not exceed safe working load specified by the support surface, frame, and accessories. The use of this support surface for patients outside the therapeutic weight range of 50 lbs to 350 lbs should be determined by a health care professional. Patients shall meet the minimum age requirement of 2 years old.

This support surface is not intended to support a patient in a prone position, i.e., laying face downward. ProForm shall be used with a mattress cover at all times.

The support surface is not intended to be a sterile product nor is it intended to include a measuring function. ProForm has an expected product life of three (3) years for the surface and three (3) years for the cover.

PRODUCT DESCRIPTION

ProForm is a non-powered support surface focusing on equalizing pressure and enhancing comfort. It has been designed to include the following features:

Anatomic Convolute – Support surface layer that has been biometrically engineered to address comfort and support of the patient using a convoluted foam layer.

Calcaneus Curve™ Section (CAL CRV) – Design feature that provides pressure relief in the heel section.

Introduction

SPECIFICATIONS

Non-Powered Support Surface		
Dimensions	35 x 76 x 6 (inches) 35 x 80 x 6 (inches) 35 x 84 x 6 (inches)	88,9 x 193,0 x 15,2 (cm) 88,9 x 203,2 x 15,2 (cm) 88,9 x 213,4 x 15,2 (cm)
Safe Working Load	500 lbs.	226,8 kg.
Therapeutic Weight Range	50 – 350 lbs.	22,7 – 158,8 kg.
Product Compliance	Without Fire Barrier	With Fire Barrier
Flammability Compliance Reference	16 CFR 1632	16 CFR1632 16 CFR 1633 CAL TB 129 CGSB/CAN 2-4.2 Method 27.7-M77 BOSTON IX - 11

Note: See page 0-18 for all available Replacement Parts List.

Stryker reserves the right to change specifications without notice.

Introduction

CONTACT INFORMATION

Contact Stryker Customer Service or Technical Support at: **(800) 327-0770** or **(269) 324-6500**.

Stryker Medical

3800 E. Centre Avenue
Portage, MI 49002
USA

Please have the serial number of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

SERIAL NUMBER LOCATION

Standard Model: The serial number (**A**) is located inside the cover, on the foam core, near the foot right corner of the support surface as shown in **Figure 1**. To reference the serial number, unzip the cover about one foot to access the foam crib. Serial number can be read through the fire barrier if model contains a fire barrier.

Psych Model: This model does not have a removable cover and the serial number is located on the labels attached to the cover at foot right corner of the support surface.


Figure 1

SERIAL NUMBER FORMAT (15 DIGITS):

Serial Number Example: 051227102600001

0	5	1	2	2	7	1	0	2	6	0	0	0	0	1
M	M	Y	Y	S	S	S	S	X	X	N	N	N	N	N

Manufacture Date (MM/YY): May 2012

Base Model (S): 2710 (ProForm)

Model Number (X): 26

Size: 84" x 35"

Features: CAL CRV, Stretch Cover, Fire Barrier

Sequential Number (N): 00001

Year Legend (Y)	
2012	12
2013	13
2014	14
2015	15
2016	16

Base Model Legend (S)
2710

Sequential # Legend (N)
00001 - 99999

Month Legend (M)	
January	01
February	02
March	03
April	04
May	05
June	06
July	07
August	08
September	09
October	10
November	11
December	12

Model Number Legend (X)		
01	13	25
02	14	26
03	15	27
04	16	28
05	17	29
06	18	30
07	19	31
08	20	32
09	21	33
10	22	34
11	23	35
12	24	36

Summary of Safety Precautions

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Inspect ProForm support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a mattress cover becomes compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
- **To avoid the risk of injury or equipment damage, do not use handles located on the bottom of the support surface (see page 0-19 & 0-20) to lift or move the support surface with a patient on it.**
- To avoid the risk of injury or equipment damage, do not transfer the patient from one bed to another using the support surface with a patient on it.
- A minimum of two operators is required when transferring a patient. Operators need to be positioned so that they can control patient positioning.
- If the space between the two patient support platforms is greater than three (3) inches, use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another. Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. **Do not use the transfer bridge to support patient load.**
- It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance.
- Check patient's skin regularly. Consult physician if any skin redness or breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
- **Some disinfectants are corrosive in nature and may cause damage to the product if used improperly.** Only use the "Suggested Disinfectants" described on page 0-12 according to the manufacturer's instructions.
- The support surface cover must be completely DRY before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.
- Failure to follow recommended cleaning/disinfection procedures may void the product warranty.
- Internal components of this support surface are not washable. If there is contamination inside the support surface, discard the support surface.
- When cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.
- Do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
- Patient stability and siderail coverage may be compromised with the use of an overlay. When using an overlay, extra caution and/or supervision is required to help reduce the likelihood of a patient falling.
- Leave the frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
- The safe use of the support surface is maximized when used in conjunction with siderails; there may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Local policies regarding the use of siderails should be taken into account. Whether and how to use siderails is a decision that should be based on each patient's individual needs and should be made by the physician, operators, and responsible parties.
- To reduce risk of patient or operator injury and support surface damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.
- The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard and siderails. Do not use the support surface when such gaps are present.
- Extra precautions are needed for a patient at risk of a fall (i.e. agitated or confused) to help reduce the likelihood of a fall.
- Using a ProForm support surface on a bed frame of a larger or smaller width than the stated ProForm width may result in the support surface sliding and patient injury.
- When placing a patient on the support surface, ensure the opposite siderail is raised to reduce the risk of patient fall.

Summary of Safety Precautions

CAUTION

- Be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
- Do not put overlays or accessories inside the cover. Doing so may reduce pressure redistribution performance.
- It is the responsibility of the caregiver team to evaluate the appropriate CPR protocol to be used with the surface.
- The use of quaternaries containing glycol ethers and/or accelerated hydrogen peroxides may result in compromised legibility of the silk screen of the cover.

NOTE! The ProForm surface must be used with a mattress cover at all times.

TRANSFER OF PATIENT FROM ONE PATIENT SUPPORT PLATFORM TO ANOTHER

WARNING

- **To avoid risk of injury or equipment damage, do not transfer a patient from one bed to another using the handles located on the bottom of the ProForm support surface (see page 0-19 & 0-20).** This could result in injury to the patient and/or damage to the product. If the space between the two patient support platforms is greater than three (3) inches, use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another. Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. Do not use the transfer bridge to support patient load.
- Leave the frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
- Patient stability and siderail coverage may be compromised with the use of an overlay. When using an overlay, extra caution and/or supervision is required to help reduce the likelihood of a patient falling.
- The safe use of the support surface is maximized when used in conjunction with siderails; there may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Local policies regarding the use of siderails should be taken into account. Whether and how to use siderails is a decision that should be based on each patient's individual needs and should be made by the physician, operators, and responsible parties.
- When placing a patient on the support surface, ensure the opposite siderail is raised to reduce the risk of patient fall.
- The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard, and siderails. Do not use the support surface when such gaps are present.
- Using a ProForm support surface on a bed frame of a larger or smaller width than the stated ProForm width may result in the support surface sliding and patient injury.
- Extra precautions are needed for a patient at risk of a fall (i.e., agitated or confused) to help reduce the likelihood of a fall.
- To reduce risk of patient or user injury and support surface damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.

To transfer the patient from one patient support platform (for example, bed, stretcher, gurney, operating table) to another:

1. Position one patient support platform alongside the other patient support platform as closely as possible.
2. Set the brakes to "ON" for both patient support platforms.
3. Prior to transfer of patient between fixed patient support surfaces, ensure that the surfaces are level with one another.
4. Raise the patient support platform siderail located opposite the patient transfer.
5. Transfer the patient following all applicable safety rules and institution protocols to ensure patient and operator safety.

CAUTION

- Be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.

USING INCONTINENCE MANAGEMENT DEVICES

This support surface is not intended to manage incontinence. Therefore, it is recommended to use incontinence management devices when appropriate. Disposable diapers or incontinence pads may be used.

WARNING

- It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance.
- Check patient's skin regularly. Consult physician if any skin redness or breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.

CLEANING SUPPORT SURFACE

WARNING

- To avoid the risk of equipment damage, do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
 - Do not iron, dry-clean, or tumble dry the support surface, as this will cause malfunction and damage the product.
 - Do not stick needles into a medical bed mattress through the mattress cover. Holes will allow body fluids from entering the inside (inner core) of the mattress and could cause cross-contamination.
 - To avoid the risk of patient injury, covers (top and bottom) should be inspected for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a support surface cover becomes compromised, the support surface should be removed from service immediately and replaced to prevent cross-contamination.
-

Stryker support surfaces are designed for long-lasting comfort and reliability. The life of the support surface can be adversely affected by an increase in frequency of usage which might include more frequent cleaning and disinfection.

1. Using a clean, soft, damp cloth, wipe down the entire support surface with a mild soap and water solution to remove foreign material.
 2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or cleaning agent.
 3. Care must be taken to thoroughly rinse and dry covers following cleaning.
 4. Disinfect as necessary with a hospital grade disinfectant AFTER cleaning has been completed. Refer to "Disinfecting Support Surface" on page 1-12.
-

CAUTION

To avoid the risk of equipment damage, when cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.

Cleaning

DISINFECTING SUPPORT SURFACE

WARNING

- To avoid the risk of patient injury, disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
- Some disinfectants may cause damage to the product if used improperly. If the products described below are used to disinfect the support surface, measures must be taken to ensure the entire surface is wiped with a damp cloth soaked in clean water and thoroughly dried following disinfection. The cover can be damaged when exposed to such disinfectants beyond the manufacturers' recommendations. Failure to follow these directions when using these types of disinfectants may void this product warranty.
- To avoid the risk of equipment damage, the support surface cover must be completely dry before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.

Suggested Disinfectants:

- Quaternaries are acceptable if ph < (less than) 10. DO NOT USE quaternaries containing sodium hydroxide (NaOH)
- Chlorinated Bleach Solution (5.25% bleach. For daily use, dilute to 1 part bleach to 1,000 part water. For decontamination, dilute to a maximum of 1 part bleach to 100 parts water).
- 70% Isopropyl Alcohol is acceptable as a disinfectant but pooling of the cleaning solution must be avoided and care must be taken until fully dry.

1. Ensure surface has been thoroughly cleaned and dried prior to applying disinfectants.
2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or disinfectant.
3. Care must be taken to thoroughly **rinse and dry** covers following disinfection.

CAUTION

- Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
- The use of accelerated hydrogen peroxides or quaternaries containing glycol ethers may damage the cover.
- DO NOT USE Phenolic cleaning agents as it result in irreversible swelling of coating.

RECOMMENDATION: The Safety Communication lists several recommendations for inspection and maintenance including:

- Regularly check each medical bed mattress cover for any visible signs of damage or wear such as cuts, tears, cracks, pinholes, snags or stains.
- Routinely remove the medical bed mattress cover and check its inside surface. Once the mattress cover is removed, inspect the mattress for wet spots, staining, or signs of damage or wear. Check all sides and the bottom of the mattress.
- Immediately replace any medical bed mattress cover with visible signs of damage or wear to reduce the risk of infection to patients.
- DO NOT stick needles into a medical bed mattress through the mattress cover.

COVER REPLACEMENT

Tools Required:

- None

Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient's right side and then remove the foam crib assembly from the bed and set aside.
5. Discard the old cover.
6. Place the new cover, unzipped and open, on the bed with the black bottom cover on the litter and the top cover folded over the patient's right side of the bed.
7. Carefully place the foam assembly on top of the bottom part of the cover so it is square with the cover.

NOTE! Maintain the correct "Head-Foot" orientation between the new cover and the foam assembly.

8. Fold the top cover over the top of the foam crib assembly so it is square with the foam crib assembly.
9. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
10. Verify proper operation of the unit before returning it to service.

Preventive Maintenance

Preventive maintenance should be performed annually, at a minimum. A preventive maintenance program should be established for all Stryker Medical equipment. Preventive maintenance may need to be performed more frequently based on the usage level of the product.

CHECKLIST

- _____ Cover zipper opens and closes properly and has no visible damage
- _____ No tears, rips, holes, cracks, or other openings in the mattress cover
- _____ No rips or cracks in the optional fire barrier (if present)
- _____ Check labels (see page 0-16 – 0-17) for legibility, proper adherence and integrity
- _____ No rips or cracks in the handles or stitching
- _____ Support surface cover straps are intact and are not damaged
- _____ Foam and other components have not degraded or come apart
- _____ Fire barrier (if present) has no visible sign of damage. If excessive wear is observed, replacement of mattress is recommended.
- _____ Compression set of the foam is less than 1.25" (3.2 cm) (see page 0-15)
- _____ Verify proper operation of the unit before returning it to service

Product Serial Number:		

Complete by: _____ **Date:** _____

CHECKING FOR COMPRESSION SET

Tools Required:

- 35" + Straight Edge
- Tape Measure

Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down position.

NOTE! Before taking the compression set measurement, ensure that there has not been a patient on the support surface for a minimum of two (2) hours.

3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient's right side.
5. If the support surface has the optional fire barrier, starting at the foot end, roll up the fire barrier on the foam crib assembly.

NOTE! Work from side to side a little at a time to get the fire barrier past the seat section of the foam crib.

6. Place a straight edge (at least 35" (88.9 cm) long) across the foam crib from left to right at the sacral region (patient seat section) of the support surface as shown in **Figure 2**.
7. Using a tape measure, measure the maximum depth (lowest point) from the bottom of the straight edge to the top of the foam as shown in **Figure 3**.

NOTE! Do not push down on the tape measure when taking the measurement; the tape measure should just lightly touch the top of the foam.

Figure 2


Figure 3

8. Document the measurement, serial number and date the measurement was taken. This information will be required for warranty information, if applicable.

NOTE! If the measurement is more than 1.25" (3.2 cm), the compression set is out of tolerance and the support surface should be replaced. If the support surface is still in warranty, call Customer Service (see "Contact Information" on page 0-6) with the information recorded above for a replacement support surface.

9. Reverse steps 1 - 5 to reinstall.
10. Verify proper operation of the unit before returning it to service.

Product Labeling

All labels (shown below) are attached to the ProForm support surface:

<p style="text-align: center;">①</p> <p style="text-align: center; font-size: small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">UNDER PENALTY OF LAW THIS TAG NOT TO BE REMOVED EXCEPT BY THE CONSUMER</p> <hr/> <p style="text-align: center;">ALL NEW MATERIAL Consisting of</p> <p style="text-align: center;">POLYURETHANE FOAM....100%</p> <hr/> <p>REG. NO. 36016</p> <p style="font-size: x-small;">Certification is made by the manufacturer that the materials in this article are described in accordance with law.</p> <hr/> <p style="text-align: center;">Distributed by: STRYKER MEDICAL 3800 E. Centre Avenue Portage, MI 49002 USA</p> <hr/> <p style="font-size: x-small;">MODEL/MODELE: 2710-PROFORM MATTRESS OPTIONS/OPTIONS DU MATELAS</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; font-size: x-small;">Finished Size:</th> <th style="text-align: left; font-size: x-small;">Net Wt.</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> 35"W x 76"L x 6"H</td> <td>21 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 80"L x 6"H</td> <td>22 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 84"L x 6"H</td> <td>23 lb.</td> </tr> </tbody> </table> <hr/> <p><input type="checkbox"/> WITH FLAME BARRIER/ AVEC PROTECTION ANTI-FEU</p> <p><input type="checkbox"/> WITHOUT FLAME BARRIER/ SANS PROTECTION ANTI-FEU</p> <p><input type="checkbox"/> SERVICE PART/PIÈCE DE SERVICE</p> <p>Date of Delivery:..... MM/YYYY</p> <hr/> <p style="font-size: x-small; text-align: center;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Finished Size:	Net Wt.	<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.	<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.	<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.	<p style="text-align: center;">②</p> <p style="text-align: center; font-size: x-small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">Not to be removed until delivered to the consumer</p> <p style="text-align: center;">Ne pas enlever avant la livraison au consommateur</p> <hr/> <p style="text-align: center;">This article contains NEW MATERIAL ONLY</p> <p style="text-align: center;">Cet article contient des MATERIAUX NEUFS SEULEMENT</p> <hr/> <p style="font-size: x-small; text-align: center;">This label is affixed in compliance with Provincial Law</p> <p style="font-size: x-small; text-align: center;">Cette étiquette est apposée conformément à la loi de la province</p> <hr/> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">Made by</td> <td style="width: 50%;">Fabriqué par</td> </tr> <tr> <td>Reg. No.</td> <td>No. de permis:</td> </tr> <tr> <td style="text-align: center;">05T-00532911</td> <td></td> </tr> </table> <hr/> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; font-size: x-small;">Content</th> <th style="text-align: left; font-size: x-small;">Contenu</th> </tr> </thead> <tbody> <tr> <td>Polyurethane Foam 100%</td> <td>Mousse de polyuréthane</td> </tr> </tbody> </table> <hr/> <p style="font-size: x-small; text-align: center;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Made by	Fabriqué par	Reg. No.	No. de permis:	05T-00532911		Content	Contenu	Polyurethane Foam 100%	Mousse de polyuréthane	<p style="text-align: center;">③</p> <p style="text-align: center; font-size: x-small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">Manufactured By: Tridien Medical 9901 Kincaid Drive, Suite 300 Fishers, IN 46037 USA</p> <hr/> <p>Date of Manufacture:..... MM/YYYY</p> <hr/> <p style="font-size: x-small;">Model: ProForm</p> <p style="font-size: x-small;">Prototype ID:PM12-001-02</p> <hr/> <p style="font-size: x-small;">This mattress meets the requirements of 16 CFR Part 1633 (federal Flammability (open flame) standard for mattress sets) when used without a foundation.</p> <div style="border: 1px solid black; padding: 10px; text-align: center; margin: 10px 0;"> <p>THIS MATTRESS IS INTENDED TO BE USED WITHOUT A FOUNDATION</p> </div> <p style="font-size: x-small; text-align: center;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>
Finished Size:	Net Wt.																			
<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.																			
<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.																			
<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.																			
Made by	Fabriqué par																			
Reg. No.	No. de permis:																			
05T-00532911																				
Content	Contenu																			
Polyurethane Foam 100%	Mousse de polyuréthane																			

Note: Product Labels 1 to 3 (shown above) are required by US or Canada.

All labels (shown below) are attached to the ProForm support surface:

④

STITCHING TO BE ABOVE THIS LINE

MODEL
2710 PROFORM

VALID WITHOUT FLAME BARRIER ONLY

CONFORMS WITH FLAMMABILITY
STANDARD 16 CFR 1632
CARE SHOULD BE EXERCISED NEAR
OPEN FLAME OR WITH BURNING
CIGARETTES

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VALID WITH FLAME BARRIER ONLY

CONFORMS WITH FLAMMABILITY
STANDARD US-16 CFR 1633, 1632
BOSTON -BFD IX-11, CANADA-METHOD
27.7-1979 of CAN 2-4.2 M77, CALIFORNIA-TB 129
WHEN USED WITHOUT A FOUNDATION

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SERIAL NUMBER IS LOCATED INSIDE THE
MATTRESS AT THE FOOT END

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WARNING

AVERTISSEMENT

PATIENT STABILITY AND SIDERAILS COVERAGE
MAY BE COMPROMISED WITH THE USE OF AN
OVERLAY

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350 lbs/
158.8 kg

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WARNING

DO NOT TRANSFER PATIENT FROM ONE
BED TO ANOTHER USING THE SUPPORT
SURFACE WITH A PATIENT ON IT

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
Made in U.S.A. / Fabriqué aux Etats-Unis

Quick Reference Replacement Parts List


The parts and accessories listed on this page are currently available for purchase. Some of the parts identified on the assembly drawing parts in this manual may not be individually available for purchase. Please call *Stryker Customer Service USA* at **1-800-327-0770** for availability and pricing.

Part Name	Size (Inches)	Part Number	Part Description	
Psych Mattress Assembly	35 x 76	2710-000-033	Dartex [®] W/ Fire Barrier W/ CAL CRV W/ Box	
		2710-000-034	Dartex W/ Fire Barrier W/ Box	
	35 x 80	2710-000-035	Dartex W/ Fire Barrier W/ CAL CRV W/ Box	
		2710-000-036	Dartex W/ Fire Barrier W/ Box	
		2710-000-037	Dartex W/ Fire Barrier	
	35 x 84	2710-000-038	Dartex W/ Fire Barrier	
		2710-000-039	Dartex W/ Fire Barrier W/ Box	
	Mattress Assembly	35 x 80	2710-000-001	Nylon/Vinyl W/ Fire Barrier W/ CAL CRV
			2710-000-002	Nylon/Vinyl W/ Fire Barrier W/ CAL CRV W/ Box
2710-000-005			Nylon/Vinyl W/ Fire Barrier	
2710-000-006			Nylon/Vinyl W/ Fire Barrier W/ Box	
2710-000-003			Nylon/Vinyl W/O Fire Barrier W/ CAL CRV	
2710-000-004			Nylon/Vinyl W/O Fire Barrier W/ CAL CRV W/ Box	
2710-000-007			Nylon/Vinyl W/O Fire Barrier	
2710-000-008			Nylon/Vinyl W/O Fire Barrier W/ Box	
2710-000-009			Dartex W/ Fire Barrier W/ CAL CRV	
2710-000-010			Dartex W/ Fire Barrier W/ CAL CRV W/ Box	
2710-000-013			Dartex W/ Fire Barrier	
2710-000-014			Dartex W/ Fire Barrier W/ Box	
2710-000-011			Dartex W/O Fire Barrier W/ CAL CRV	
2710-000-012			Dartex W/O Fire Barrier W/ CAL CRV W/ Box	
2710-000-015			Dartex W/O Fire Barrier	
2710-000-016			Dartex W/O Fire Barrier W/Box	
35 x 84		2710-000-017	Nylon/Vinyl W/ Fire Barrier W/ CAL CRV	
		2710-000-018	Nylon/Vinyl W/ Fire Barrier W/ CAL CRV W/Box	
		2710-000-021	Nylon/Vinyl W/ Fire Barrier	
		2710-000-022	Nylon/Vinyl W/ Fire Barrier W/ Box	
		2710-000-019	Nylon/Vinyl W/O Fire Barrier W/ CAL CRV	
		2710-000-020	Nylon/Vinyl W/O Fire Barrier W/ CAL CRV W/ Box	
		2710-000-023	Nylon/Vinyl W/O Fire Barrier	
		2710-000-024	Nylon/Vinyl W/O Fire Barrier W/ Box	
		2710-000-025	Dartex W/ Fire Barrier W/ CAL CRV	
		2710-000-026	Dartex W/ Fire Barrier W/ CAL CRV W/ Box	
		2710-000-029	Dartex W/ Fire Barrier	
		2710-000-030	Dartex W/ Fire Barrier W/ Box	
		2710-000-027	Dartex W/O Fire Barrier W/ CAL CRV	
		2710-000-028	Dartex W/O Fire Barrier W/ CAL CRV W/ Box	
		2710-000-031	Dartex W/O Fire Barrier	
		2710-000-032	Dartex W/O Fire Barrier W/Box	
Cover Assembly	35 x 80	2710-000-100	Nylon/Vinyl W/ CAL CRV	
		2710-000-101	Nylon/Vinyl	
		2710-000-102	Dartex W/ CAL CRV	
		2710-000-103	Dartex	
	35 x 84	2710-000-104	Nylon/Vinyl W/ CAL CRV	
		2710-000-105	Nylon/Vinyl	
		2710-000-106	Dartex W/ CAL CRV	
		2710-000-107	Dartex	

Mattress Assembly with Cover


Top View

Cover Assembly Bottom View

<u>ITEM</u>	<u>PART NAME</u>	<u>QTY.</u>
A	Mattress Cover Assembly	1
B	Foam Assembly	1

Mattress Assembly with Cover & Fire Barrier

Top View

Cover Assembly Bottom View

<u>ITEM</u>	<u>PART NAME</u>	<u>QTY.</u>
A	Mattress Cover Assembly	1
B	Foam Assembly	1
C	Fire Barrier Sleeve	1

LIMITED WARRANTY

Stryker Medical's ProForm Support Surface is designed for an expected service life as listed below under normal use, conditions, and with appropriate periodic maintenance as described in this manual.

- **Foam Assembly:** Three (3) years
- **Mattress Cover Assembly:** Three (3) years

The foam components of a support surface will naturally compress over time. For warranty purposes, only compression sets that are greater than 1.25" (3.2 cm), and verified by Stryker will be covered (see "Preventive Maintenance" on page 0-15).

This statement constitutes Stryker's entire warranty with respect to the aforesaid equipment. **Stryker makes no other warranty or representation, either expressed or implied, except as set forth herein. There is no warranty of merchantability and there are no warranties of fitness for any particular purpose. In no event shall Stryker be liable here under for incidental or consequential damages arising from or in any manner related to sales or use of any such equipment.**

TO OBTAIN PARTS AND SERVICE

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative or call Stryker Customer Service USA at **1-800-327-0770**.

RETURN AUTHORIZATION

Merchandise cannot be returned without approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned merchandise. Stryker reserves the right to charge shipping and restocking fees on returned items. **Special, modified, or discontinued items not subject to return.**

DAMAGED MERCHANDISE

ICC Regulations require that claims for damaged merchandise must be made with the carrier within fifteen (15) days of receipt of merchandise. **Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt.** Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claim will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the merchandise, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full. Claims for any short shipment must be made within thirty (30) days of invoice.

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
AVERTISSEMENT/ATTENTION/REMARQUE

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Français

①	②	③																								
<p>STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p>UNDER PENALTY OF LAW THIS TAG NOT TO BE REMOVED EXCEPT BY THE CONSUMER</p> <hr/> <p>ALL NEW MATERIAL Consisting of</p> <p>POLYURETHANE FOAM....100%</p> <hr/> <p>REG. NO. 36016</p> <p>Certification is made by the manufacturer that the materials in this article are described in accordance with law.</p> <hr/> <p>Distributed by: STRYKER MEDICAL 3800 E. Centre Avenue Portage, MI 49002 USA</p> <hr/> <p>MODEL/MODELE: 2710-PROFORM MATTRESS OPTIONS/OPTIONS DU MATELAS</p> <table border="0"> <thead> <tr> <th style="text-align: left;">Finished Size:</th> <th style="text-align: left;">Net Wt.</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> 35"W x 76"L x 6"H</td> <td>21 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 80"L x 6"H</td> <td>22 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 84"L x 6"H</td> <td>23 lb.</td> </tr> <tr> <td colspan="2"><input type="checkbox"/> WITH FLAME BARRIER/ AVEC PROTECTION ANTI-FEU</td> </tr> <tr> <td colspan="2"><input type="checkbox"/> WITHOUT FLAME BARRIER/ SANS PROTECTION ANTI-FEU</td> </tr> <tr> <td colspan="2"><input type="checkbox"/> SERVICE PART/PIÈCE DE SERVICE</td> </tr> </tbody> </table> <p>Date of Delivery:..... MM/YYYY</p> <hr/> <p>Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Finished Size:	Net Wt.	<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.	<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.	<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.	<input type="checkbox"/> WITH FLAME BARRIER/ AVEC PROTECTION ANTI-FEU		<input type="checkbox"/> WITHOUT FLAME BARRIER/ SANS PROTECTION ANTI-FEU		<input type="checkbox"/> SERVICE PART/PIÈCE DE SERVICE		<p>STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p>Not to be removed until delivered to the consumer Ne pas enlever avant la livraison au consommateur</p> <hr/> <p>This article contains NEW MATERIAL ONLY</p> <p>Cet article contient des MATERIAUX NEUFS SEULEMENT</p> <hr/> <p>This label is affixed in compliance with Provincial Law Cette étiquette est apposée conformément à la loi de la province</p> <hr/> <table border="0"> <tr> <td>Made by</td> <td>Fabriqué par</td> </tr> <tr> <td>Reg. No.</td> <td>No. de permis:</td> </tr> <tr> <td colspan="2" style="text-align: center;">05T-00532911</td> </tr> </table> <hr/> <table border="0"> <tr> <td style="text-align: center;">Content</td> <td style="text-align: center;">Contenu</td> </tr> <tr> <td>Polyurethane Foam 100%</td> <td>Mousse de polyuréthane</td> </tr> </table> <hr/> <p>Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Made by	Fabriqué par	Reg. No.	No. de permis:	05T-00532911		Content	Contenu	Polyurethane Foam 100%	Mousse de polyuréthane	<p>STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p>Manufactured By: Tridien Medical 9901 Kincaid Drive, Suite 300 Fishers, IN 46037 USA</p> <p>Date of Manufacture:..... MM/YYYY</p> <hr/> <p>Model: ProForm Prototype ID:PM12-001-02</p> <hr/> <p>This mattress meets the requirements of 16 CFR Part 1633 (federal Flammability (open flame) standard for mattress sets) when used without a foundation.</p> <div style="border: 1px solid black; padding: 5px; text-align: center;"> <p>THIS MATTRESS IS INTENDED TO BE USED WITHOUT A FOUNDATION</p> </div> <hr/> <p>Made in U.S.A. / Fabriqué aux Etats-Unis</p>
Finished Size:	Net Wt.																									
<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.																									
<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.																									
<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.																									
<input type="checkbox"/> WITH FLAME BARRIER/ AVEC PROTECTION ANTI-FEU																										
<input type="checkbox"/> WITHOUT FLAME BARRIER/ SANS PROTECTION ANTI-FEU																										
<input type="checkbox"/> SERVICE PART/PIÈCE DE SERVICE																										
Made by	Fabriqué par																									
Reg. No.	No. de permis:																									
05T-00532911																										
Content	Contenu																									
Polyurethane Foam 100%	Mousse de polyuréthane																									

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①

STITCHING TO BE ABOVE THIS LINE

**UNDER PENALTY OF LAW
THIS TAG NOT TO BE REMOVED
EXCEPT BY THE CONSUMER**

ALL NEW MATERIAL

Consisting of

POLYURETHANE FOAM....100%

REG. NO. 36016

Certification is made by the manufacturer that the materials in this article are described in accordance with law.

Distributed by:
STRYKER MEDICAL
3800 E. Centre Avenue
Portage, MI 49002 USA

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②

STITCHING TO BE ABOVE THIS LINE

Not to be removed until
delivered to the consumer

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This article contains
**NEW MATERIAL
ONLY**

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This label is affixed in
compliance with Provincial Law

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③

STITCHING TO BE ABOVE THIS LINE

Manufactured By:
Tridien Medical
9901 Kincaid Drive, Suite 300
Fishers, IN 46037 USA

Date of Manufacture:.....
MM/YYYY

Model: ProForm

Prototype ID:PM12-001-02

This mattress meets the requirements of 16
CFR Part 1633 (federal Flammability (open
flame) standard for mattress sets) when used
without a foundation.

**THIS MATTRESS IS
INTENDED TO BE USED
WITHOUT A FOUNDATION**

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Português

<p style="text-align: center;">①</p> <p style="text-align: center; font-size: small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">UNDER PENALTY OF LAW THIS TAG NOT TO BE REMOVED EXCEPT BY THE CONSUMER</p> <hr/> <p style="text-align: center;">ALL NEW MATERIAL Consisting of</p> <p style="text-align: center;">POLYURETHANE FOAM....100%</p> <hr/> <p>REG. NO. 36016</p> <p style="font-size: x-small;">Certification is made by the manufacturer that the materials in this article are described in accordance with law.</p> <hr/> <p style="text-align: center;">Distributed by: STRYKER MEDICAL 3800 E. Centre Avenue Portage, MI 49002 USA</p> <hr/> <p style="font-size: x-small;">MODEL/MODELE: 2710-PROFORM MATTRESS OPTIONS/OPTIONS DU MATELAS</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; font-size: x-small;">Finished Size:</th> <th style="text-align: left; font-size: x-small;">Net Wt.</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> 35"W x 76"L x 6"H</td> <td>21 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 80"L x 6"H</td> <td>22 lb.</td> </tr> <tr> <td><input type="checkbox"/> 35"W x 84"L x 6"H</td> <td>23 lb.</td> </tr> </tbody> </table> <hr/> <p><input type="checkbox"/> WITH FLAME BARRIER/ AVEC PROTECTION ANTI-FEU</p> <p><input type="checkbox"/> WITHOUT FLAME BARRIER/ SANS PROTECTION ANTI-FEU</p> <p><input type="checkbox"/> SERVICE PART/PIÉCE DE SERVICE</p> <p>Date of Delivery:..... MM/YYYY</p> <hr/> <p style="text-align: center; font-size: x-small;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Finished Size:	Net Wt.	<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.	<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.	<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.	<p style="text-align: center;">②</p> <p style="text-align: center; font-size: x-small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">Not to be removed until delivered to the consumer</p> <p style="text-align: center;">Ne pas enlever avant la livraison au consommateur</p> <hr/> <p style="text-align: center;">This article contains NEW MATERIAL ONLY</p> <p style="text-align: center;">Cet article contient des MATERIAUX NEUFS SEULEMENT</p> <hr/> <p style="font-size: x-small;">This label is affixed in compliance with Provincial Law</p> <p style="font-size: x-small;">Cette étiquette est apposée conformément à la loi de la province</p> <hr/> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">Made by Reg. No.</td> <td style="width: 50%;">Fabriqué par No. de permis:</td> </tr> <tr> <td style="text-align: center;">05T-00532911</td> <td></td> </tr> </table> <hr/> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; font-size: x-small;">Content</th> <th style="text-align: left; font-size: x-small;">Contenu</th> </tr> </thead> <tbody> <tr> <td style="font-size: x-small;">Polyurethane Foam 100%</td> <td style="font-size: x-small;">Mousse de polyuréthane</td> </tr> </tbody> </table> <hr/> <p style="text-align: center; font-size: x-small;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>	Made by Reg. No.	Fabriqué par No. de permis:	05T-00532911		Content	Contenu	Polyurethane Foam 100%	Mousse de polyuréthane	<p style="text-align: center;">③</p> <p style="text-align: center; font-size: x-small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">Manufactured By: Tridien Medical 9901 Kincaid Drive, Suite 300 Fishers, IN 46037 USA</p> <hr/> <p>Date of Manufacture:..... MM/YYYY</p> <hr/> <p style="font-size: x-small;">Model: ProForm Prototype ID:PM12-001-02</p> <hr/> <p style="font-size: x-small;">This mattress meets the requirements of 16 CFR Part 1633 (federal Flammability (open flame) standard for mattress sets) when used without a foundation.</p> <div style="border: 1px solid black; padding: 10px; text-align: center; margin: 10px 0;"> <p>THIS MATTRESS IS INTENDED TO BE USED WITHOUT A FOUNDATION</p> </div> <hr/> <p style="text-align: center; font-size: x-small;">Made in U.S.A. / Fabriqué aux Etats-Unis</p>
Finished Size:	Net Wt.																	
<input type="checkbox"/> 35"W x 76"L x 6"H	21 lb.																	
<input type="checkbox"/> 35"W x 80"L x 6"H	22 lb.																	
<input type="checkbox"/> 35"W x 84"L x 6"H	23 lb.																	
Made by Reg. No.	Fabriqué par No. de permis:																	
05T-00532911																		
Content	Contenu																	
Polyurethane Foam 100%	Mousse de polyuréthane																	

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<p style="text-align: center;">④</p> <p style="text-align: center; font-size: small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <p style="text-align: center;">型号 2710 PROFORM</p> <hr/> <p>VALID WITHOUT FLAME BARRIER ONLY CONFORMS WITH FLAMMABILITY STANDARD 16 CFR 1632 CARE SHOULD BE EXERCISED NEAR OPEN FLAME OR WITH BURNING CIGARETTES</p> <p style="text-align: center;">检测协议号：4912</p> <p style="text-align: center;">仅对无明火屏障有效 CONFORMS 易燃性标准 16 CFR 1632 靠近明火或燃烧的烟头时需小心。</p> <hr/> <p>VALID WITH FLAME BARRIER ONLY CONFORMS WITH FLAMMABILITY STANDARD US-16 CFR 1633, 1632 BOSTON -BFD IX-11, CANADA-METHOD 27.7-1979 of CAN 2-4.2 M77, CALIFORNIA-TB 129 WHEN USED WITHOUT A FOUNDATION</p> <p style="text-align: center;">检测协议号：7892, 4852, 2641</p> <p style="text-align: center;">仅对无明火屏障有效 CONFORMS 易燃性标准 US-16 CFR 1632, 1632 BOSTON -BFD IX-11, CANADA-METHOD 27.7-1979 of CAN 2-4.2 M77, CALIFORNIA-TB 129 在无基础情时使用</p> <hr/> <p style="text-align: center;">SERIAL NUMBER IS LOCATED INSIDE THE MATTRESS AT THE FOOT END</p> <p style="text-align: center;">序列号位于垫套的内部的底座端</p> <hr/> <p style="text-align: center;">Made in U.S.A. / 美国制造</p>	<p style="text-align: center;">⑤</p> <p style="text-align: center; font-size: small;">STITCHING TO BE ABOVE THIS LINE</p> <hr style="border-top: 1px dashed black;"/> <div style="display: flex; align-items: center;"> <div> <p>WARNING</p> <p>警告</p> <p>PATIENT STABILITY AND SIDERAILS COVERAGE MAY BE COMPROMISED WITH THE USE OF AN OVERLAY</p> <p>使用覆层可能使患者的稳定性和护栏覆盖范 围受到损害。</p> </div> </div> <hr/> <p style="text-align: center;">CARE LABEL 注意标签</p> <div style="display: flex; justify-content: space-around; align-items: center;"> </div> <div style="display: flex; justify-content: space-around; align-items: center; margin-top: 10px;"> </div> <hr/> <div style="display: flex; align-items: center; justify-content: center;"> <div> <p>350 lbs/ 158.8 kg</p> </div> </div> <hr/> <p style="text-align: center;">APPROVED BY / 批准者</p> <hr/> <div style="display: flex; align-items: center;"> <div> <p>WARNING</p> <p>DO NOT TRANSFER PATIENT FROM ONE BED TO ANOTHER USING THE SUPPORT SURFACE WITH A PATIENT ON IT</p> <p>警告</p> <p>切勿在存在患者的情况下， 使用支撑面将患者移至其他病床。</p> </div> </div> <hr/> <p style="text-align: center;">Made in U.S.A. / 美国制造</p>
---	---

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