ProForm® non-powered support surface

Operations Manual

REF 2815

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Catalogue number
	Serial number
	Manufacturer
	Date of manufacture
	Unique device identifier
	Safe working load
	Wash by hand
	Do not tumble dry
	Do not dry-clean
	Do not iron
	Allow to air dry
	Chlorinated bleach
	For US Patents see www.stryker.com/patents
	Keep dry
	Stack limit by number
	This side up

	Fragile
	Do not use sharp objects to open the package
	Mass of equipment

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Contraindications	5
Expected service life - standard	5
Expected life - behavioral health	5
Specifications	5
Standard	5
Behavioral health	6
Environmental conditions	6
Contact information	6
Serial number location	7
Setup	8
Operation	9
Transferring a patient from one patient support platform to another	9
Preventive maintenance	10
Care and maintenance	10
Accessories and Parts	12
Cover replacement - standard model	12

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Never use the standard model support surface in the behavioral health use environment. The zipper allows access to the internal components which could be used for harm.
- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
- Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on non-compatible bed frames.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Do not use the support surface as a transfer device.
- Do not use the support surface handles to lift or move the support surface with a patient on board.
- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. If the space between the two patient support platforms is greater than 3 in. (7.6 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another.
- Always raise the opposite siderail when you place a patient on the support surface to reduce the risk of patient fall.
- Do not wash the internal components of this support surface. Discard the support surface if a contamination is found inside.
- Do not immerse the product in cleaning or disinfectant solutions.
- Do not allow liquid to pool on the product.

- Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these cleaning instructions may void your warranty.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ether >3% as these chemicals may damage the support surface cover.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
 - Always take caution when you replace the support surface cover. The internal fire barrier contains fiberglass fibers. Dust particles from the fibers may cause skin irritation.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 2815 **ProForm**® is a 6" thick, non-powered foam support surface for the general patient population. There are articulation cuts on the mattress core. The mattress cover material options include nylon material or END406 by **Dartex**® fabric. All configurations of the product include a heel slope at the foot end of the support surface. The dimensions of this product are compatible with Stryker bed frames. See the compatible bed frame list in the specifications section of this manual.

Intended use

ProForm uses foam to help with pressure redistribution and immersion.

This product is a non-powered support surface that assists in the prevention and treatment of pressure injuries or pressure ulcers (all stages, Unstageable injury, and Deep Tissue injury). We recommend you implement this product in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional.

The standard Model (281505550001, 281505550003, 281505550005, 281505550007) is intended for the general hospital environment, specifically the acute care healthcare facilities.

The behavioral health Model (281505550009, 281505550010, 281505550011) are intended for use in behavioral care facilities.

ProForm shall be used with a support surface cover at all times. The support surface cover can interact with all external skin. A top sheet should be used with this product.

The standard general hospital environment models and the behavioral health models are not intended for:

- Patients that exceed the safe working load
- Bariatric patients
- Home health environment setting
- Not a sterile product
- Does not contain a measuring function

Stryker promotes the clinical assessment of each patient and appropriate usage by the operator.

Contraindications

None known.

Expected service life - standard

The **ProForm** support surface core and covers have a 3 year expected service life under normal use, conditions, and with appropriate periodic maintenance.

Expected life - behavioral health

The **ProForm** support surface has a 3 year expected life under normal use, conditions, and with appropriate periodic maintenance. There are no serviceable parts for this product.

Specifications

Safe working load 	500 lb	226.8 kg
Minimum patient weight	50 lb	22.7 kg

Note - The patient must not exceed safe working load specified by the support surface.

Standard

Model general patient	281505550001 or 281505550020 bagged		281505550003	
Length	84 in.	213.4 cm	80 in.	203.2 cm
Width	35 in.	88.9 cm	35 in.	88.9 cm
Thickness	6 in.	15.2 cm	6 in.	15.2 cm
Product weight	30 lb	13.6 kg	29 lb	13.2 kg
Top cover material	END406 by Dartex®			
Model general patient	281505550005		281505550007	
Length	84 in.	213.4 cm	80 in.	203.2 cm
Width	35 in.	88.9 cm	35 in.	88.9 cm
Thickness	6 in.	15.2 cm	6 in.	15.2 cm
Product weight	30 lb	13.6 kg	29 lb	13.2 kg
Top cover material	Nylon			

Mattress material	Polyurethane foam
Product compliance with fire barrier	16CFR1632, 16CFR1633, CGSB CAN 2-4.2 Method 27.7-M77, CAL TB129, BFD IX-11
Compatible frames	3002 S3® , 3005 S3® , GoBed® II , Spirit Select® with high quad 4 zone siderails

Behavioral health






Model	281505550010, 281505550011		281505550009	
Length	80 in.	203.2 cm	84 in.	213.4 cm
Width	35 in.	88.9 cm	35 in.	76 cm
Thickness	6 in.	15.2 cm	6 in.	7.6 cm
Product weight	29 lb	13.2 kg	30 lb	13.6 cm
Top cover material	END406 by Dartex®			

Mattress material	Polyurethane foam
Product compliance with fire barrier	16CFR1632, 16CFR1633, CGSB CAN 2-4.2 Method 27.7-M77, CAL TB129, BFD IX-11
Compatible frames	3002 S3® , 3005 S3® , GoBed® II , Spirit Select® with high quad 4 zone siderails, Spirit® Behavioral Health Platform Bed

Note - The Behavioral health covers do not have zippers. The product is not serviceable.

Environmental conditions

WARNING - Never use the standard model support surface in the behavioral health use environment. The zipper allows access to the internal components which could be used for harm.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure	Not applicable	

Stryker reserves the right to change specifications without notice.

Contact information


Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
 3800 E. Centre Avenue
 Portage, MI 49002
 USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Setup

WARNING

- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
 - Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
 - Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on non-compatible bed frames.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
-

CAUTION

- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
-

To setup the support surface:

1. Make sure that the support surface fits the bed frame.
2. Place the support surface with the artwork face up at the head end of the bed frame.
3. Place the support surface between the mattress retainers on the bed frame.
4. Place linens on the support surface per hospital protocols.

Operation

Transferring a patient from one patient support platform to another

WARNING

- Do not use the support surface as a transfer device.
 - Do not use the support surface handles to lift or move the support surface with a patient on board.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
 - Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. If the space between the two patient support platforms is greater than 3 in. (7.6 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another.
 - Always raise the opposite siderail when you place a patient on the support surface to reduce the risk of patient fall.
-

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Apply the brakes for both patient support platforms.
3. Adjust the patient support platform heights so that they are level with one another.
4. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.

Preventive maintenance

Remove the product from service before you perform the preventive maintenance inspection. Check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Note - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items for the standard models:

- _____ Zipper and covers (top and bottom) for tears, cuts, holes, or other openings
- _____ Internal components for signs of staining from fluid ingress or contamination by fully unzipping the covers
- _____ Labels for legibility, proper adherence, and integrity
- _____ Handles for rips, cracks, stitching, or other visible signs of damage
- _____ Foam has not degraded or come apart
- _____ Fire barrier for rips, cracks, or other visible signs of damage

Inspect the following items for the behavioral health models:

- _____ Covers (top and bottom) for tears, cuts, holes, or other openings
- _____ Labels for legibility, proper adherence, and integrity

Product serial number:
Completed by:
Date:

Care and maintenance

WARNING

- Do not wash the internal components of this support surface. Discard the support surface if a contamination is found inside.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not allow liquid to pool on the product.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these cleaning instructions may void your warranty.
-

CAUTION

- Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ether >3% as these chemicals may damage the support surface cover.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
-

The support surface cover is resistant to the following chemical solutions:

- Quaternaries (active ingredient - ammonium chloride) that contain \leq 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (10000 ppm) for use with END406 mattress covers
- Chlorinated bleach solution (1000 ppm) for use with Nylon mattress covers
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients to avoid the risk of cross-contamination and infection.

Note - Chlorinated bleach solution may cause discoloration of the Nylon mattress covers.

Accessories and Parts

These parts are available for purchase for the standard models, but subject to change. Call Stryker Customer Service: 1-800-327-0770 for available parts and prices.

Name	Number
Cover assembly, heel slope, END406 80 in.	281505550014
Cover assembly, heel slope, END406 84 in.	281505550012
Cover assembly, heel slope, Nylon 80 in.	281505550018
Cover assembly, heel slope, Nylon 84 in.	281505550016

Cover replacement - standard model

CAUTION - Always take caution when you replace the support surface cover. The internal fire barrier contains fiberglass fibers. Dust particles from the fibers may cause skin irritation.

Tools required:

- None

Procedure:


1. Raise the bed to the highest height position.
2. Lower the Fowler and Gatch sections to the lowest height positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient right side, and then remove the foam crib assembly from the bed.
5. Remove and discard the cover.
6. Place the replacement cover, unzipped and open, on the bed with the black bottom cover on the litter. Fold the top cover over the patient right side of the bed.
7. Place the foam crib assembly on top of the bottom part of the cover. Make sure that the foam crib aligns with the cover.
8. Fold the top cover over the top of the foam crib assembly. Make sure that the top cover aligns with the foam crib assembly.
9. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
10. Verify proper operation before you return the product to service.

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