IsoFlex LAL support surface

Operations Manual

REF 2860

Global symbol glossary

See the Global Symbol Glossary at ifu.stryker.com for symbol definitions.

Symbols

	Consult instructions for use
	General warning
	Caution
	Catalogue number
	Serial number
	Batch code
	For US Patents see www.stryker.com/patents
	Manufacturer
	Date of manufacture
	Unique device identifier
	Safe working load
	Mass of product
	Class II with the powered low air loss (LAL)
	Wash by hand
	Do not tumble dry
	Do not dry clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
	Keep dry


	Stacking limit by number
	Do not use sharp object to cut open the package

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Contraindications	4
Expected service life	4
Specifications	5
Product illustration	6
Product features	6
Contact information	7
Serial number location	7
Date of manufacture	7
Setup	8
Operation	9
Applying the linens	9
Transferring a patient from one patient support platform to another	9
Managing incontinence and drainage	10
Connecting and disconnecting the pump air hose option	10
Accessories and parts	11
Preventive maintenance	12
Checking the tolerance of the gel compression set	12
Top 360 cover replacement	13
Bottom 360 cover replacement	14
Fire barrier replacement	15
Cap with tether replacement	16
Cleaning and disinfecting with SideKick	17
Cleaning	18
Disinfecting	19

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
- Always leave the bed frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails and there may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Always use extra caution when reading radiology images taken of a patient on this support surface as internal components can cause artifacts and distort readings.
- Always install the **IsoFlex LAL** support surface on Stryker bed frames. See compatible frames in the specification table. This is to avoid the risk of a safety hazard including but not limited to bodily injury.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches. between the support surface and the headboard, footboard, and siderails.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Always evaluate the appropriate CPR protocol to use with this product.
- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
- Do not use the support surface as a transfer device.
- Do not use the support surface handles to lift or move the support surface with a patient on board.
- Do not transfer patient from one bed to another using the support surface with a patient on it.
- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. If the space between the two patient support platforms is greater than 3 in. (7,62 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another.

- Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
 - Always monitor the patient condition at regular intervals for patient safety.
 - Do not immerse the support surface in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If compromised immediately remove the support surface covers from service and replace the covers to prevent cross-contamination.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning instructions may void your warranty.
 - Do not iron, dry clean, or tumble dry the support surface covers.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface. Drying the product helps to prevent the performance of the product from being impaired.
 - Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Do not allow fluid to seep into the zipper area or watershed cover barrier when cleaning the underside of the support surface. Fluid allowed to come in contact with the zipper may leak into the support surface.
 - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover and reduce the legibility of the artwork.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The non-powered Stryker Model 2860 IsoFlex LAL support surface assists in improving human patient outcomes by equalizing pressure and enhancing comfort. IsoFlex LAL is available in different dimensions see the specifications table. IsoFlex LAL has a breathable cover to assist in moisture management of the skin. There is an optional low air loss (LAL) pump.

Intended use

The IsoFlex LAL support surface assists in the prevention and treatment of all pressure injuries (including stages 1,2,3,4, Unstageable injury, and Deep Tissue Injury) and is recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a health care professional.

This support surface is for use with human patients with existing or at risk of incurring pressure injuries. The patient must not exceed the safe working load as specified by the support surface. A healthcare professional should determine the use of this support surface by patients outside of the therapeutic weight range. The minimum patient age requirement for this support surface is two years.

IsoFlex LAL shall be used with a support surface cover at all times. The support surface cover can interact with all external skin.

This support surface is for use by patients in an acute care setting. This may include critical care, step down, progressive care, med/surg, subacute care, and post anesthesia care unit (PACU), or other locations as prescribed by a physician. Operators of this support surface include healthcare professionals (such as nurses, nurse aids, or doctors).

IsoFlex LAL support surface is not intended for use in a home health environment, for use as a sterile product, for use with a patient less than two years of age, and does not include a measuring function.

Stryker promotes the clinical assessment of each patient and appropriate usage by the operator.

Contraindications

None known.

Expected service life

The IsoFlex LAL support surface has a ten year expected service life under normal use, conditions, and with appropriate periodic maintenance.







The IsoFlex LAL covers have a three year expected service life under normal use, conditions, and with appropriate periodic maintenance.

Specifications


  Safe working load Note - The patient must not exceed the safe working load specified by the support surface.	500 lb	226.8 kg
Therapeutic weight	50 lb - 350 lb	22.7 kg - 158.7 kg
Length	84 in.	213.4 cm
Width	35 in.	88.9 cm
Thickness	6 in.	15.2 cm
Weight	57 lb	25.9 kg
Compatible frames	Epic II, GoBed II, InTouch, ProCuity, 3002 S3, 3005 S3	
	Bed frames that support 84 in. x 35 in. (213.4 cm x 88.9 cm) support surfaces should be assessed by the healthcare facility or end user. An assessment is to be completed if the support surface is used with other flat decks not listed. The combination of the support surface installed on a flat deck bed frame must meet the applicable IEC bed standards.	
Compatible pump	AIR+ Pump Model 2863	

Stryker reserves the right to change specifications without notice.

For more information about the AIR+ Pump, see the Model 2863 Operations Manual.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Product illustration

A	Head end	D	Strap handles (located at the head end and the foot end)
B	Foot end	E	D-ring (located on each corner of the bottom cover)
C	Support surface	F	Pump air attachment

Product features

1	Pressure redistribution	Redistributes pressure in the sacral region with its CoreGel technology by allowing immersion and envelopment to occur
2	Shear management	Moves freely with the patient to assist in the reduction of shear forces at the interface between the surface and the patient with ShearGel support layer
3	Moisture management	Provides air flow with the LAL system to assist managing the heat and humidity of the patient's skin

This product is made with Intelli-Gel®* hollow column configuration and Duragel™ elastomeric material.

*Intelli-Gel® is a registered trademark of EdiZONE, LLC of Alpine, UT

Duragel™ is a trademark of EdiZONE, LLC of Alpine, UT

Contact information


Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Date of manufacture

The date of manufacture is the first four digits of the serial number.

MMYY (MM = month and YY = year)

Setup

WARNING

- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
- Always leave the bed frame in the lowest position when the patient is unattended to help reduce the number and severity of falls.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails and there may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Always use extra caution when reading radiology images taken of a patient on this support surface as internal components can cause artifacts and distort readings.
- Always install the **IsoFlex LAL** support surface on Stryker bed frames. See compatible frames in the specification table. This is to avoid the risk of a safety hazard including but not limited to bodily injury.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches. between the support surface and the headboard, footboard, and siderails.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Always evaluate the appropriate CPR protocol to use with this product.
- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.

CAUTION

- Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
-

Prerequisite: A minimum of two operators is required for this task.

To install the support surface:

1. Make sure that the support surface properly fits on the bed frame where you are placing the product.
2. Make sure that you install the dedicated sloped heel section at the foot end of the bed frame.
3. Place the linens on the support surface per hospital protocol.

Operation

Applying the linens

To apply the linens:

1. Apply the linens using the “D” rings for a flat sheet.
2. Thread the four linen corners through the “D” rings (C) on the bottom cover to secure the linens (A) to the support surface (B).


Figure 1 – Apply linens

Transferring a patient from one patient support platform to another

WARNING

- Always check patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
 - Do not use the support surface as a transfer device.
 - Do not use the support surface handles to lift or move the support surface with a patient on board.
 - Do not transfer patient from one bed to another using the support surface with a patient on it.
 - Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
 - Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient. If the space between the two patient support platforms is greater than 3 in. (7,62 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another.
 - Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
-

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one patient support platform to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes for both patient support platforms.
3. Adjust the patient support platform heights so that they are level with one another.
4. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.

Note - Do not use the IsoFlex LAL support surface to transfer the patient to other surfaces.

5. Center the patient on the support surface (Figure 2).


Figure 2 – Center patient

Managing incontinence and drainage

WARNING - Always monitor the patient condition at regular intervals for patient safety.

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

Connecting and disconnecting the pump air hose option

Connect the AIR+ Pump air hose to provide Low Air Loss (LAL) to the support surface.

To connect the pump air hose to the support surface:

1. Pull the air attachment box (B) outside of the support surface (A) (Figure 3).


Figure 3 – Air attachment

2. Grasp and turn the air attachment cap to the left to unscrew the air attachment cap from the body valve (C).
3. Connect the pump air hose. Turn the air hose cap to the right to connect the pump air hose to the body valve.

To disconnect the pump air hose from the support surface:

4. Disconnect the pump air hose. Turn the air hose cap to the left to disconnect the pump air hose from the body valve.
5. Reinstall the air attachment cap removed in step 2.

Note - Always close the air attachment cap when the pump is not connected.

6. Push the air attachment box inside the support surface.

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

84" support surface	Part number
Kit, 360 cover assembly	2860-700-007
Fire barrier sleeve	2860-030-427
Cable tie fitting	7600-001-405
Cap screw with tether	2860-030-425
Kit, cap screw with tether (includes 2860-030-425 quantity of 10)	2860-700-001
Kit, O-Ring fitting, quantity 10	2860-700-011
O-ring fitting	2860-030-430
Reduction coupler fitting	7600-001-404
AIR+ pump	2863-000-000

Preventive maintenance

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage.

Remove product from service before performing preventive maintenance.

Note - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

_____ Zipper, covers (top and bottom), and fire barriers are free of tears, cuts, holes, or other openings

Note - If excessive wear is observed, replacement is recommended.

_____ Internal components for signs of staining from fluid ingress or contamination by fully unzipping the covers

_____ Labels for legibility, proper adherence, and integrity

_____ Handles are free of rips or cracks and stitching is intact

_____ Foam and other components have not degraded or come apart

_____ Compression set of the gel is less than 1.25" (3.2 cm) (see *Checking the tolerance of the gel compression set* (page 12))

_____ Optional LAL connections for cracks, disconnection, or other visible signs of damage

Note

- The fire barrier will contain a hole for the LAL tubing. This is acceptable.
- Foam discoloration due to oxidation is normal. This does not affect the performance or integrity of the support surface. The discoloration may be greater in areas of adhesive application due to chemical interaction.
- During the molding process, the gel forms a grid pattern. The cosmetic appearance of the corners may look like an issue, but are not defects in the product.

Product serial number:
Completed by:
Date:

Checking the tolerance of the gel compression set

Tools required:

- 35" + straight edge
- Tape measure

Procedure:

1. Raise the bed height to the highest height position.
2. Lower the Fowler backrest and Gatch sections to the lowest height positions.

Note - Before you take measurements for the compression set, make sure that there has not been a patient on the support surface for a minimum of one hour.

3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient right side.

Note - Work from side to side, a little at a time, to move the fire barrier to the top of the support surface.

5. Place a straight edge (at least 35" (88.9 cm) long) across the gel from left to right at the sacral region (patient seat section) of the support surface (Figure 4).

- Using a tape measure, measure the maximum depth (lowest point) from the bottom of the straight edge to the top of the foam (Figure 5).

Note - Do not push down on the tape measure when you take the measurement. The tape measure should just lightly touch the top of the gel grid.


Figure 4 – Position straight edge


Figure 5 – Measure maximum depth

- Document the measurement, serial number, and date the measurement was taken. This information will be required for warranty information, if applicable.

Note

- If the measurement is greater than 1.25" (3.2 cm), the compression set is out of tolerance and you should replace the support surface. See the warranty page for replacement requirements. Call Customer Service at 1-800-327-0770 with the information recorded above for a replacement support surface.
 - The support surface internal components will naturally compress over time. Should a body indentation or compression set measurement greater than 1 1/4" (3.2 cm), as measured by an authorized Stryker representative, be realized within the warranty period noted above, Stryker will provide a replacement product. Any normal body indentations or compression set of less than 1 1/4" (3.2 cm) will not be replaced.
 - Any damage to the foam assembly or gel that results due to usage of a cover assembly beyond its warranty period of three years, or is a result of abnormal wear and tear that may include cleaning processes that are inconsistent with those recommended in this Operations/Maintenance manual, shall invalidate the warranty on the support surface at Stryker's sole discretion.
- Reverse steps to reinstall.
 - Verify proper operation before you return the product to service.

Top 360 cover replacement

Tools required:

- None

Procedure:

- Raise the product to the highest height position.
- Lower the Fowler backrest and Gatch sections to the lowest height positions.
- Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.

Note - Use caution not to damage the fire barrier.

- Remove the top cover.
- Clean or discard the top 360 cover.
 - Clean the top 360 cover. See *Cleaning* (page 18).
 - Discard the top 360 cover per hospital protocols. See *Preventive maintenance* (page 12).
- Reverse the steps to install the supplied top 360 cover.

7. Verify proper operation before you return the product to service.

Bottom 360 cover replacement

Tools required:

- Diagonal pliers
- Zip tie gun

Procedure:

1. Raise the bed height to the full up position.
2. Lower the Fowler and Gatch sections to the full down positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.
4. Remove and set aside the top cover.
5. Using the diagonal pliers, carefully cut the cable tie (C) closest to the foot end of the support surface cover that secures the LAL hose and reduction coupler fitting (A) to the valve body (B) (Figure 6).

Note - Use caution not to damage the fire barrier by cutting or snagging with the diagonal pliers.


Figure 6 – Reduction coupler fitting

6. Remove the reduction coupler fitting from the valve body.
7. Remove the foam crib assembly from the bed. Save the foam crib assembly.
8. Discard the bottom cover per hospital protocol.
9. Place the supplied cover on the bed with the black bottom cover on the litter.
10. Place the foam crib assembly on top of the black bottom cover. Align the foam crib to the cover.
11. Install the reduction coupler fitting by pushing firmly until the fitting is fully seated into the valve body.

Note - Make sure that the flame barrier is tight around the hose assembly.

12. Using a zip tie gun, center the supplied cable tie (E) around the valve body (D) and the reduction coupler fitting (Figure 7).


Figure 7 – Install the cable tie to the valve body

13. Place the top cover over the top of the foam crib assembly. Make sure that the top cover aligns with the foam crib assembly.
14. Zip the cover to close. The zipper starts in the middle of the patient right side of the support surface.
15. Verify proper operation before you return the product to service.

Fire barrier replacement

Tools Required:

- Diagonal pliers
- Utility knife
- Zip tie gun

Procedure:

1. Raise the bed height to the highest height position.
2. Lower the Fowler backrest and Gatch sections to the lowest height positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.
4. Fold the top of the cover to the patient right side.
5. Using the diagonal pliers, carefully cut the cable tie (C) closest to the foot end of the support surface cover that secures the LAL hose and reduction coupler fitting (A) to the valve body (C) (Figure 6).

Note - Use caution not to damage the fire barrier by cutting or snagging it with the diagonal pliers.

6. Remove the reduction coupler fitting from the valve body.
7. Starting at the foot end, roll up the fire barrier on the foam crib assembly.

Note - Work from side to side, a little at a time, to move the fire barrier to the top of the support surface.

8. Discard the fire barrier.
9. Starting at the head end, roll the new fire barrier down and slide the fire barrier over the foam crib assembly.

Note - Align the fire barrier on the foam crib before sliding it over the foam crib assembly.

10. Carefully slide the fire barrier down the foam crib assembly, working from side to side, to make sure that the fire barrier is tight on the foam crib assembly.
11. Align the foam crib assembly on top of the bottom part of the cover.

Note - Spread the excess fire barrier material equally below the foam crib assembly at the foot end.

12. Locate the LAL tubing and, using a utility knife, cut a hole through the fire barrier to make sure that you do not obstruct the connection point for the hose.
13. Install the reduction coupler fitting. Push the reduction coupler fitting until the coupler is fully seated into the valve body.

Note - Make sure that the flame barrier is tight around the hose assembly.

14. Using a zip tie gun, center the supplied cable tie (E) around the valve body (D) and the reduction coupler fitting (Figure 7).
15. Fold and align the top cover over the top of the foam crib assembly.
16. Zip the cover to close.
17. Verify proper operation before you return the product to service.

Cap with tether replacement

Tools required:

- None

Procedure:

1. Disconnect the air pump hose from the patient right side foot end of the support surface, if attached.
2. Grasp the body valve (C) and pull outward so the air attachment box (B) extends outside the support surface cover (A) (Figure 8).
3. Push the tether (E) (Figure 9) over the body valve (C) (Figure 8). Turn the tether (E) in a clockwise motion until it is past the threads on the body valve (C).
4. Make sure that the cap (D) (Figure 9) threads onto the body valve (C) (Figure 8).
5. Push the air attachment box (B) back inside the support surface.


Figure 8 – Air attachment box and body valve


Figure 9 – Cap and tether

Cleaning and disinfecting with SideKick

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred 6" x 10" wipes (2060-000-001) include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface with SideKick wipe:

To clean:

1. Wipe down the external product surface with a fresh, clean wipe to remove all visible soils.
2. Repeat as necessary until the external product surface is visibly clean.
3. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Note - Use as many wipes as necessary.

To disinfect:

1. Clean first.
2. Wipe down the external product surface with a fresh, clean wipe until wet.
3. Allow the external product surface to remain wet for two minutes at room temperature.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Cleaning

WARNING

- Do not immerse the support surface in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If compromised immediately remove the support surface covers from service and replace the covers to prevent cross-contamination.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning instructions may void your warranty.
 - Do not iron, dry clean, or tumble dry the support surface covers.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface. Drying the product helps to prevent the performance of the product from being impaired.
-

CAUTION

- Do not allow fluid to seep into the zipper area or watershed cover barrier when cleaning the underside of the support surface. Fluid allowed to come in contact with the zipper may leak into the support surface.
-

Note - If the optional pump is installed, unplug the product before cleaning.

Always follow hospital protocol for cleaning and disinfecting.

To clean the support surface covers between patient uses, follow these steps in order:

1. Using a clean, soft, damp cloth, wipe the support surface covers with a mild soap and water solution to remove foreign material.
2. Wipe the support surface covers with a clean, dry cloth to remove any excess fluid or cleaning agent.
3. Dry thoroughly.

Disinfecting

WARNING

- Always disinfect the support surface between patients to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning instructions may void your warranty.
 - Do not immerse the support surface in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface. Drying the product helps to prevent the performance of the product from being impaired.
-

CAUTION

- Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover and reduce the legibility of the artwork.
-

Suggested disinfectants:

- Quaternaries
- Phenolic disinfectants
- Chlorinated bleach solution (5.25% bleach diluted 1 part bleach to 100 parts water)
- 70% isopropyl alcohol

Note - If the optional pump is installed, unplug the product before cleaning.

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Thoroughly clean and dry the support surface covers before you apply disinfectants (See *Cleaning* (page 18)).
2. Apply recommended disinfectant solution by spray or pre-soaked wipes. Do not soak the support surface.

Note - Make sure that you follow the disinfectant's instructions for appropriate contact time and rinsing requirements.

3. Wipe the support surface covers with a clean, dry cloth to remove any excess fluid or disinfectant.
4. Allow the support surface covers to dry completely before returning to service.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **AIR+**, **CoreGel**, **Epic**, **GoBed**, **InTouch**, **IsoFlex LAL**, **ProCuity**, **S3**, **ShearGel**, **SideKick**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.