IsoGel AIR™


Support Surface/Surface de soutien/Patientenunterlage  
Superficie de apoyo/Superfície de suporte/Superficie  
antidecubito/Steunoppervlak/Støtteoverflade/  
Stödyta/-tukialusta/Podkład/Surface de soutien/  
Superficie de suporte/Støttemadrass/支持面

REF 2860

# stryker®

## Operations/Maintenance Manual

Manuel d'utilisation et d'entretien  
Bedienungs- und Wartungshandbuch  
Manual de uso y mantenimiento  
Manual de funcionamento/manutenção  
Manuale d'uso e manutenzione  
Gebruiks-/onderhoudshandleiding  
Betjenings-/vedligeholdelsesmanual  
Användar-/underhållshandbok  
Käyttö- ja huolto-ohjekirja  
Instrukcja obsługi/konserwacji  
Manuel d'utilisation/d'entretien  
Manual de funcionamento/manutenção  
Bruks-/vedlikeholdshåndbok  
取扱・保守説明書

CE

# International Addresses

## EUROPE HEADQUARTERS

**Stryker SA**  
Cite-Centre, Grand-Rue 90  
P.O. Box 1568  
1820 Montreux, Switzerland  
Phone: 41-21-966-12-01  
Fax: 41-21-966-12-00

## EASTERN EUROPE

**Stryker SA - Export Business**  
Grand-Rue 90  
P.O. Box 1567  
1820 Montreux, Switzerland  
Phone: 41-21-966-14-00  
Fax: 41-21-966-14-01

## AUSTRALIA

**Stryker Australia**  
8 Herbert Street  
St Leonards NSW 2065  
Australia  
Phone: 61-2-9467-1000  
Fax: 61-2-9467-1010

## AUSTRIA

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## ROMANIA

**Stryker Osteonics Romania S.R.L.**  
19, Leonida Str.  
District 2  
7000 Bucharest  
Romania  
Phone: 40-2-12-12-11-22  
Fax: 40-2-12-12-11-33

## SOUTH AFRICA

**Stryker Osteonics PTY. LTD.**  
3 Susan Street  
Strydom Park  
Johannesburg, 2194 - South Africa  
Mailing Address:  
P.O. Box 48039  
2129 Roosevelt Park  
Johannesburg, South Africa  
Phone: 27-11-791-4644  
Fax: 27-11-791-4696

## SPAIN

**Stryker Iberia SL**  
c/Sepulveda n17  
28108 Alcobendas  
Madrid  
Spain  
Phone: 34-91-728-35-00  
Phone: 34-91-358-20-44  
Fax: 34-91-358-07-48

## SWEDEN

**Stryker AB/Scandinavia**  
Box 50425  
SE-204 14 Malmo  
Sweden  
Phone: 46-40-691-81-00  
Fax: 46-40-691-81-91

# Table of Contents

---

English

Symbols . . . . .	<a href="#">1-2</a>
Warning/Caution/Note Definition . . . . .	<a href="#">1-3</a>
Introduction . . . . .	<a href="#">1-4</a>
Intended Use of Product . . . . .	<a href="#">1-4</a>
Product Description . . . . .	<a href="#">1-4</a>
Product Illustration . . . . .	<a href="#">1-5</a>
Product Features . . . . .	<a href="#">1-6</a>
Specifications . . . . .	<a href="#">1-7</a>
Contact Information . . . . .	<a href="#">1-8</a>
Product Serial Number Location/Identification . . . . .	<a href="#">1-8</a>
Summary of Safety Precautions . . . . .	<a href="#">1-9</a>
Setup Procedures . . . . .	<a href="#">1-11</a>
Installing the Support Surface . . . . .	<a href="#">1-11</a>
Relocating Support Surface . . . . .	<a href="#">1-12</a>
Transferring of Patient from One Patient Support Platform to Another . . . . .	<a href="#">1-12</a>
Managing Incontinence/Drainage . . . . .	<a href="#">1-13</a>
Cleaning . . . . .	<a href="#">1-14</a>
Cleaning Support Surface . . . . .	<a href="#">1-14</a>
Disinfecting Support Surface . . . . .	<a href="#">1-15</a>
Cleaning Pump and Hoses . . . . .	<a href="#">1-15</a>
Service Information . . . . .	<a href="#">1-16</a>
Replacing Mattress Cover . . . . .	<a href="#">1-16</a>
Replacing Fire Barrier . . . . .	<a href="#">1-17</a>
Preventative Maintenance . . . . .	<a href="#">1-18</a>
Checking the Body Indentation of the Gel . . . . .	<a href="#">1-18</a>
Checklist . . . . .	<a href="#">1-19</a>
Product Labels . . . . .	<a href="#">1-20</a>
Quick Reference Replacement Parts List . . . . .	<a href="#">1-21</a>
Warranty . . . . .	<a href="#">1-22</a>
To Obtain Parts and Service . . . . .	<a href="#">1-23</a>
Return Authorization . . . . .	<a href="#">1-23</a>
Damaged Product . . . . .	<a href="#">1-23</a>

# Symbols

English

	Warning - Consult accompanying documentation
	Caution
	Safe Working Load
	Do Not Iron
	Damp Wipe Only
	Chlorinated Bleach
	Do Not Tumble Dry
	Do Not Dry Clean
	Allow to Completely Air Dry
	Consult Instructions for Use
	Manufacturer
	Product Weight
	Class II with the powered low air loss (LAL) accessory.
<b>IPX4</b>	Protection from liquid splash

# Symbols and Definitions

---

## WARNING/CAUTION/NOTE DEFINITION

The words **WARNING**, **CAUTION** and **NOTE** carry special meanings and should be carefully reviewed.

---

### **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

---

### **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

---

### **NOTE**

Provides special information to make maintenance easier or important instructions clearer.

# Introduction

---

English

This manual is designed to assist with the operation and maintenance of the Model 2860 **IsoGel AIR™** support surface. Carefully read this manual thoroughly before using or beginning maintenance on the support surface. To ensure safe operation of this equipment, it is recommended that methods and procedures are established for educating and training staff on the safe operation of the support surface.

## INTENDED USE OF PRODUCT

### INTENDED MEDICAL INDICATION

**IsoGel AIR** support surface is intended to assist in improving human patient outcomes, focusing on redistributing pressure, and enhancing comfort. **IsoGel AIR** is an appropriate support surface intended to assist in improving human patient outcomes, focusing on equalizing pressure redistribution and enhancing comfort. **IsoGel AIR** is an appropriate support surface used to assist in the prevention and treatment of all categories/stages of pressure ulcers (including stages I, II, III, IV, Unstageable, and Deep Tissue Injury) and is recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a health care professional. **IsoGel AIR** will offer 35"x 84" and 35" x 80" support surface options. The 35" x 84" option is intended for use on Stryker InTouch, GoBed II, Secure III, and Epic II and other flat deck 35" x 84" frames. The 35" x 80" option is intended for use on the Stryker Florence bed, Rose bed and other flat deck 35" x 80". **IsoGel AIR** will offer a Low Air Loss (LAL) feature and a breathable cover to assist in microclimate management of the skin.

### INTENDED PATIENT POPULATION

This surface is intended to be used with patients at risk for developing pressure ulcers, those who require therapy for pre-existing pressure ulcers, or any other patient under a healthcare professional's care, in a general hospital environment or acute care environment. The safe working load for IsoGel AIR is 500 lbs (226,8 kg), the patient must not exceed safe working load specified by the support surface, frame and accessories. The use of this support surface by patients outside the therapeutic weight range of 50 lbs to 350 lbs (22,7 kg to 158,7 kg) should be determined by a health care professional. The minimum patient age requirement for this surface is two (2) years old.

### INTENDED PART OF THE BODY

The surface is not intended to support a patient in a prone position. **IsoGel AIR** support surface must be used with a support surface cover at all times.

### INTENDED USER PROFILE

The device will be used in acute care, general hospital care, or alternate care. This device is not intended to be used in a home environment; patients will be under the care of a physician. This surface is for multiple patient uses and is not a disposable device.

### INTENDED CONDITION OF USE AND EXPECTED LIFE

The device is not intended to be a sterile product nor is it intended to include a measuring function. Please refer to the specification table for environmental requirements. **IsoGel AIR** product life is 10 years for the product and 3 years for the cover.


## PRODUCT DESCRIPTION

The **IsoGel AIR** support surface is intended to assist in the prevention and treatment of pressure ulcers caused by pressure and shear.

# Introduction

## PRODUCT ILLUSTRATION

English


1	Support Surface Cover	Equilibrium cover by Dartex® wicks away moisture from the patient.
2	Strap Handles	Strap handles are located at the head end and foot end of the support surface. Strap handles are intended for operator to move the unoccupied surface from one bed frame to another.
3	Pump Air Attachment	Allows the connection to the optional air pump.

[Return To Table of Contents](#)

# Introduction

## PRODUCT FEATURES

English

<b>1</b>	Pressure Redistribution	Redistributes pressure in the sacral region with its CoreGel™ Technology allowing immersion and envelopment to occur.
<b>2</b>	Shear Management	ShearGel™ moves freely with the patient to assist in the reduction of shear forces at the interface between the surface and the patient.
<b>3</b>	Microclimate Management	Low Air Loss system provides air flow to assist in the management of the patient's skin by managing the microclimate of heat and humidity.

# Introduction

English

## SPECIFICATIONS

Support Surface with Sloped Heel		
Mattress Cover Material	Equilibrium cover by Dartex®	
Model 2860-000-006	84" x 35" x 6"	213,4 cm x 88,9 cm x 15,2 cm
Model 2860-000-008	80" x 35" x 6"	203,2 cm x 88,9 cm x 15,2 cm
Weight Model 2860-000-006	72 lb without Stryker Air pump	32,7 kg without Stryker Air pump
Weight Model 2860-000-008	70 lb without Stryker Air pump	31,8 kg without Stryker Air pump
Safe Working Load	500 lb	226,8 kg
Therapeutic Weight	50 lb - 350 lb	22,7 kg - 158,7 kg
Stryker Air pump (2861-000-002) - Optional (2861-000-001 when purchased with IsoGel AIR)		
Weight	20 lbs	9,1 kg
<b>Note: For more information on Stryker Air pump, please refer to the manufacturers instructions for use.</b>		
Environmental Conditions for IsoGel AIR Support Surface	Operation	Storage and Transportation
Ambient Temperature		
Relative Humidity (Non-Condensing)		
Atmospheric Pressure		
Product Compliance		
Classification	16 CFR 1632 16 CFR 1633 CAL TB 129 Boston BFD IX-11 Canada - Method 27.7-1979 of 2-4.2 M77	

**Stryker reserves the right to change specifications without notice.**

**Dartex® is a registered trademark of Dartex Coatings, Ltd.**

[Return To Table of Contents](#)

# Summary of Safety Precautions

---

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

English

## WARNING

- To avoid the risk of patient injury or equipment damage, do not X-ray the patient while on Isogel AIR.
  - To avoid the risk of patient injury, covers (top and bottom) should be inspected for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a support surface cover becomes compromised, the support surface should be removed from service immediately and replaced to prevent cross-contamination.
  - To avoid the risk of patient injury, do not transfer the patient from one bed to another using the support surface.
  - To avoid the risk of patient and operator injury, a minimum of two (2) operators is required when transferring a patient. Operators need to be positioned so that they can control patient positioning.
  - To avoid the risk of patient injury, do not use the transfer bridge to support patient load. If the space between the two patient support platforms is greater than 3 in (7,62 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another. Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient.
  - It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance. Consult physician if erythema or skin breakdown occurs.
  - Some disinfectants may cause damage to the product if used improperly. If the products described in "Cleaning" on [page 1-15](#) are used to disinfect the support surface, measures must be taken to ensure the entire surface is wiped with a damp cloth soaked in clean water and thoroughly dried following disinfection. The cover can be damaged when exposed to such disinfectants beyond the manufacturers' recommendations. Failure to follow these directions when using these types of disinfectants may void this product warranty.
  - To avoid the risk of equipment damage, do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
  - To avoid the risk of patient injury, disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
  - To avoid the risk of equipment damage, the support surface cover must be completely dry before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.
  - To avoid the risk of patient injury, when using an overlay, extra caution and/or supervision is required. Patient stability and siderail coverage may be compromised with the use of an overlay.
  - To avoid the risk of patient injury, leave the frame in the lowest position when the patient is unattended to help reduce falls.
  - Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. The safe use of the support surface is maximized when used in conjunction with siderails; there may be an increased risk of falls when siderails are not present. Local policies regarding the use of siderails should be taken into account. Whether and how to use siderails is a decision that should be based on each patient's individual needs and should be made by the physician, operators, and responsible parties.
  - To reduce risk of patient or operator injury or equipment damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.
  - The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard, and siderails. The support surface is not to be used when such gaps are present.
  - Extra precautions are needed for a patient at risk of a fall (i.e. agitated or confused) to help reduce the likelihood of a fall.
  - To avoid the risk of patient injury, do not use the unit on a bed frame of a larger or smaller width than the stated width as this may cause the support surface to slide.
  - To avoid the risk of patient injury, ensure the opposite siderail is raised when placing a patient on the support surface.
  - To avoid the risk of patient injury, the operators should evaluate the appropriate CPR protocol to be used with the unit.
- 

[Return To Table of Contents](#)

# Summary of Safety Precautions

---

English

## WARNING (CONTINUED)

- To avoid the risk of equipment damage, be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, excessive heat generated by the equipment, or sharp edges on the equipment.
  - To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
  - Do not iron, dry-clean, or tumble dry the support surface cover, as this will cause malfunction and damage to the product.
  - Do not stick needles in to a medical bed mattress through the mattress cover. Holes will allow body fluids from entering the inside (inner core) of the mattress and could cause cross-contamination.
  - Consult Stryker Air pump instruction for use manual BEFORE connecting the system to the **IsoGel AIR** support surface.
- 

## CAUTION

- To avoid the risk of equipment damage, when cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.
  - The use of quaternaries containing glycol ethers and/or accelerated hydrogen peroxides may damage the cover.
  - To avoid the risk of equipment damage, do not put overlays or accessories inside the cover. Doing so may reduce pressure redistribution performance.
  - Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
-

# Setup Procedures

---

## INSTALLING THE SUPPORT SURFACE

---

### **WARNING**

- The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard, and siderails. The support surface is not to be used when such gaps are present.
  - To avoid the risk of patient injury, do not use the unit on a bed frame of a larger or smaller width than the stated width as this may cause the support surface to slide.
- 

Prerequisite: Per OSHA, two (2) operators are required to setup the support surface.

1. Ensure that the support surface properly fits the bed frame on which it is being placed.
  2. Ensure the dedicated sloped heel section is installed at the foot end of the bed frame.
  3. Place the linens on the support surface as loose as possible to avoid creating a “hammocking” effect and increasing interface pressures.
- 

### **CAUTION**

To avoid the risk of equipment damage, do not put overlays or accessories inside the cover. Doing so may reduce pressure redistribution performance.

---

# Operation Guide

---

English

## RELOCATING SUPPORT SURFACE

---

### WARNING

- Using a **IsoGel AIR** support surface on a bed frame of a larger or smaller width than the stated **IsoGel AIR** width may result in the support surface sliding and patient injury.
  - The risk of entrapment can develop when the support surface is placed on bed frames that leave gaps of even a few inches between the support surface and the headboard, footboard, and siderails. The support surface is not to be used when such gaps are present.
- 

Prerequisite: Two (2) operators will be required for this task.

1. Position an operator at the foot end and one (1) at the head end of the unoccupied support surface.
2. Each operator will lift the support surface by using the two (2) handles located at the foot end and head end.

**Note:** Operators shall at all times use both handles when relocating the **IsoGel AIR** support surface.

3. Place support surface on bed frame following the “Installing the Support Surface” on [page 1-11](#).

## TRANSFERRING OF PATIENT FROM ONE PATIENT SUPPORT PLATFORM TO ANOTHER

---

### WARNING

- To avoid the risk of patient injury or equipment damage, do not X-ray the patient while on IsoGel AIR.
  - To avoid the risk of patient injury, do not use the transfer bridge to support patient load. If the space between the two patient support platforms is greater than 3 in (7,62 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another. Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient.
  - To avoid the risk of patient injury, leave the frame in the lowest position when the patient is unattended to help reduce falls.
  - Patient stability and siderail coverage may be compromised with the use of an overlay. When using an overlay, extra caution and/or supervision is required to help reduce the likelihood of a patient falling.
  - The safe use of the support surface is maximized when used in conjunction with siderails; there may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Local policies regarding the use of siderails should be taken into account. Whether and how to use siderails is a decision that should be based on each patient’s individual needs and should be made by the physician, operators, and responsible parties.
  - To avoid the risk of patient injury, ensure the opposite siderail is raised when placing a patient on the support surface.
  - Extra precautions are needed for a patient at risk of a fall (i.e., agitated or confused) to help reduce the likelihood of a fall.
  - To reduce risk of patient or user injury or equipment damage, do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface.
-

# Operation Guide

## TRANSFERRING OF PATIENT FROM ONE PATIENT SUPPORT PLATFORM TO ANOTHER (CONTINUED)

### WARNING

- To avoid the risk of patient injury, do not transfer the patient from one bed to another using the support surface.
- To avoid the risk of patient injury, do not use the transfer bridge to support patient load. If the space between the two patient support platforms is greater than 3 in (7,62 cm), use the transfer bridge to fill the gap. The transfer bridge is meant to ease transfer of a patient from one patient support platform to another. Always make sure that the patient support platforms and their respective transfer gaps are adequate to support the patient.

To transfer the patient from one patient support platform (for example, bed, stretcher, gurney, operating table) to another:

Prerequisite: Two (2) operators will be required for this task.

1. Center the patient on the support surface.
2. Position one patient support platform alongside the other patient support platform as closely as possible.
3. Set the brakes to "ON" for both patient support platforms. Only use **IsoGel AIR** for patient transfers between fixed patient support surfaces that are level with one another.
4. Raise the patient support platform siderail located opposite the patient transfer.
5. Transfer the patient following all applicable safety rules and institution protocols to ensure patient and operator safety. Do not use IsoGel AIR support surface to transfer the patient to other surfaces.

### CAUTION

- To avoid the risk of patient injury, the operators should evaluate the appropriate CPR protocol to be used with the unit.

## MANAGING INCONTINENCE/DRAINAGE

This support surface is not intended to manage incontinence. Therefore, it is recommended to use incontinence management devices when appropriate. Disposable diapers or incontinence pads may be used. Ensure appropriate skin care is provided following each episode.

### WARNING

- To avoid the risk of equipment damage, be aware of devices or equipment placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, excessive heat generated by the equipment, or sharp edges on the equipment.
- It is the responsibility of the operator to monitor the patient and the patient's skin condition at regular intervals, per hospital protocols, to ensure patient safety and proper support surface performance. Consult a physician if erythema or skin breakdown occurs.

# Cleaning

---

## CLEANING SUPPORT SURFACE

---

### WARNING

- To avoid the risk of equipment damage, do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
  - Do not iron, dry-clean, or tumble dry the support surface cover, as this will cause malfunction and damage to the product.
  - Do not stick needles in to a medical bed mattress through the mattress cover. Holes will allow body fluids from entering the inside (inner core) of the mattress and could cause cross-contamination.
  - To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
  - To avoid the risk of patient injury, covers (top and bottom) should be inspected for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a support surface cover becomes compromised, the support surface should be removed from service immediately and replaced to prevent cross-contamination.
- 

Stryker support surfaces are designed for long-lasting comfort and reliability. The life of the support surface can be adversely affected by an increase in frequency of usage which might include more frequent cleaning and disinfection.

Note: If pump is installed, unplug from unit prior to cleaning.

1. Using a clean, soft, damp cloth, wipe down the entire support surface with a mild soap and water solution to remove foreign material.
  2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or cleaning agent.
  3. Care must be taken to thoroughly **rinse and dry** covers following cleaning.
  4. Disinfect as necessary with a hospital grade disinfectant AFTER cleaning has been completed. Refer to “Disinfecting Support Surface” on [page 1-15](#).
- 

### CAUTION

To avoid the risk of equipment damage, when cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.

---

# Cleaning

---

## DISINFECTING SUPPORT SURFACE

---

### **WARNING**

- To avoid the risk of patient injury, disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
  - Some disinfectants may cause damage to the product if used improperly. If the products described below are used to disinfect the support surface, measures must be taken to ensure the entire surface is wiped with a damp cloth soaked in clean water and thoroughly dried following disinfection. The cover can be damaged when exposed to such disinfectants beyond the manufacturers' recommendations. Failure to follow these directions when using these types of disinfectants may void this product warranty.
  - To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
  - To avoid the risk of equipment damage, the support surface cover must be completely dry before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.
- 

Suggested Disinfectants:

- Quaternaries
- Phenolic Disinfectant
- Chlorinated Bleach Solution (5.25% bleach diluted 1 part bleach to 100 parts water)
- 70% Isopropyl Alcohol

**Note:** If pump is installed, unplug from unit prior to disinfection.

1. Ensure surface has been thoroughly cleaned and dried prior to applying disinfectants.
  2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or disinfectant.
  3. Care must be taken to thoroughly **rinse and dry** covers following disinfection.
- 

### **CAUTION**

- Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
  - The use of quaternaries containing glycol ethers or accelerated hydrogen peroxides may damage the cover.
- 

## CLEANING PUMP AND HOSES

---

### **WARNING**

- To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
  - To avoid the risk of cross-contaminations, disinfect Stryker Air pump and hose between patient installations, utilizing standard hospital protocol and disinfectants and manufacturers recommendations for cleaning.
- 

1. Clean and disinfect the hose, pump, and power cord utilizing standard hospital protocol and disinfectants after referring to Stryker Air pump manufacturer's operations manual. Do not autoclave.
2. Refer to Stryker Air pump manufacturer's operations manual for storage recommendations.

# Service Information

English

## REPLACING MATTRESS COVER

### Tools Required:

- Diagonal Pliers
- Zip Tie Gun

### Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Using diagonal pliers, carefully cut the cable tie closest to the foot end of the mattress cover securing the low air loss hose and reduction coupler fitting to the valve body (See Figure 3).

**Note:** Use caution not to damage the flame barrier.

5. Remove the reduction coupler fitting from the valve body.
6. Fold the top of the cover to the patient's right side and then remove the foam crib assembly from the bed and set aside.
7. Discard the old cover per hospital protocol.
8. Install the new cover, unzipped and open, on the bed with the black bottom cover on the litter and the top cover folded over the patient's right side of the bed.


Figure 3


Figure 4

9. Reinstall the foam crib assembly by placing it on top of the bottom part of the cover so it is square with the cover.
10. Install the reduction coupler fitting by pushing firmly until it is fully seated into the valve body (See Figure 4).

**Note:** Verify the flame barrier is tight around the hose assembly.

11. Using a zip tie gun, install the supplied cable tie around the valve body, centered on the reduction coupler fitting (See Figure 4).
12. Fold the top cover over the top of the foam crib assembly so it is square with the foam crib.
13. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
14. Verify proper operation of the unit before returning it to service.

## REPLACING FIRE BARRIER

### Tools Required:

- Diagonal Pliers
- Utility Knife
- Zip Tie Gun

### Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down positions.
3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient's right side.
5. Using diagonal pliers, carefully cut the cable tie closest to the foot end of the mattress cover securing the low air loss hose and reduction coupler fitting to the valve body (See Figure 3 on page 20).

**Note:** Use caution not to damage the flame barrier.

6. Remove the reduction coupler fitting from the valve body.
7. Starting at the foot end, roll up the fire barrier on the foam crib assembly.

**Note:** Work from side to side a little at a time to get the fire barrier to the top.

8. Discard the old fire barrier, per hospital protocol.
9. Starting at the head end, roll the new fire barrier up and slide it over the foam crib assembly.

**Note:** Square the fire barrier on the foam crib before sliding it too far over the foam crib assembly.

10. Carefully slide the fire barrier down the foam crib assembly, working side to side to make sure the fire barrier is tight on the foam crib assembly.
11. Carefully adjust the foam crib assembly on top of the bottom part of the cover so it is square with the cover.
12. Locate the LAL tubing and, using a utility knife, cut a hole through the fire barrier to ensure the connection point for the hose is not obstructed.

**Note:** Spread the excess fire barrier material equally below the foam crib assembly at the foot end.

13. Install the reduction coupler fitting by pushing firmly until it is fully seated into the valve body (See Figure 4 on page 20).

**Note:** Verify the flame barrier is tight around the hose assembly.

14. Using a zip tie gun, install the supplied cable tie around the valve body, centered on the reduction coupler fitting (See Figure 4 on page 20).
15. Fold the top cover over the top of the foam crib assembly to ensure that it is square with the foam crib assembly.
16. Zip the cover to close. Start at the head end patient right corner and stop at the foot end patient right corner.
17. Verify proper operation of the unit before returning it to service.

# Preventative Maintenance

English

## CHECKING THE BODY INDENTATION OF THE GEL

### Tools Required:

- 35" + Straight Edge
- Tape Measure

### Procedure:

1. Raise the bed height to the full up position.
2. Lower the fowler and gatch sections to the full down position.

**Note:** Before taking the compression set measurement, ensure that there has not been a patient on the support surface for a minimum of one hour.

3. Unzip the cover. Start at the foot end patient right corner of the support surface and stop at the head end patient right corner.
4. Fold the top of the cover to the patient right side.

**Note:** Work from side to side a little at a time to get the fire barrier past the seat section of the foam crib.

5. Place a straight edge (at least 35" (88,9 cm) long) across the foam from left to right at the sacral region (patient seat section) of the support surface as shown in Figure 5.
6. Using a tape measure, measure the maximum depth (lowest point) from the bottom of the straight edge to the top of the gel as shown in Figure 6.

**Note:** Do not push down on the tape measure when taking the measurement; the tape measure should just lightly touch the top of the gel grid.


Figure 5


Figure 6

7. Document the measurement, serial number, and date the measurement was taken. This information will be required for warranty information, if applicable.

**Note:** If the measurement is greater than 1 1/4" (3.2 cm), the compression set is out of tolerance and the support surface should be replaced. Refer to the warranty page for replacement requirements. Call Customer Service (see "Contact Information" on [page 1-8](#)) with the information recorded above for a replacement support surface.

8. Reverse steps to reinstall.
9. Verify proper operation of the unit before returning it to service.

[Return To Table of Contents](#)

# Preventative Maintenance

Preventative maintenance should be performed annually at a minimum. A preventative maintenance program should be established for all Stryker Medical equipment. Preventative maintenance may need to be performed more frequently based on the usage level of the product.

## CHECKLIST

- \_\_\_\_\_ Verify the cover zipper opens and closes properly and has no visible damage
- \_\_\_\_\_ Inspect to ensure there are no tears, rips, holes, cracks, or other openings in the cover (top and bottom)  
**Note:** If excessive wear is observed on the cover, it is strongly recommended to replace the cover.
- \_\_\_\_\_ Inspect to ensure that there are no rips or excessive wear in the fire barrier
- \_\_\_\_\_ Check labels as specified in the operation/maintenance manual for legibility, proper adherence, and integrity
- \_\_\_\_\_ Inspect to ensure that there are no rips or cracks in the handles or stitching
- \_\_\_\_\_ Inspect optional LAL connections to ensure they are not cracked, disconnected or inoperable

**Note:**

- The fire barrier will contain a hole for the LAL tubing which is acceptable.
- Foam discoloration is normal due to oxidation, but does not affect the performance or integrity of the surface. The discoloration may be greater in areas of adhesive application due to chemical interaction.
- The Gel is molded into a grid pattern. The gel when molded may lead to cosmetic issues in the corners of the grid but are not defects in the product.










**Support Surface Serial No.**

Completed by: \_\_\_\_\_ Date: \_\_\_\_\_

# Product Labels

All labels (shown below) are attached to the IsoGel AIR™ support surface.

English

<p><b>LOT</b></p>										
<p><b>English</b></p>	<p><b>MODEL 2860 - IsoGel AIR™</b></p>									
<p>This product made with Intelli-Gel®* polymeric columned material and Gelastic™* polymeric columned material.</p> <p>*EdiZONE, LLC of Alpine, UT</p> <p>This product is covered by one of more of the following US patents: 7,964,664; 7,076,822; 6,026,527; 5,749,111 and other US and/or foreign patents and pending patent applications.</p>										
<p><b>CONFORMS WITH FLAMMABILITY STANDARDS WHEN USED WITHOUT A FOUNDATION:</b></p> <ul style="list-style-type: none"> <li>- USA-16 CFR 1632</li> <li>- USA-16 CFR 1633</li> <li>- BOSTON-BFD IX-11</li> <li>- CANADA-METHOD 27.7-1979 of CAN 2-4.2 M77</li> <li>- CALIFORNIA-TB 129</li> <li>- EUROPE BS-7177: 2008 for medium hazard</li> <li>- ITALY UNI-9175</li> </ul> <p>DHF-1316-31 JUL 2012</p>										
<p><b>REF</b></p>	<p><b>2860</b></p> <p>SERIAL NUMBER IS LOCATED INSIDE THE MATTRESS AT THE FOOT END</p> <p>PER OSHA REQUIREMENTS, TWO (2) OPERATORS ARE REQUIRED TO SETUP THIS SUPPORT SURFACE</p>									
<p><b>CE</b>  <b>Stryker Medical</b> 3800 E. Centre Avenue Portage, MI 49002 USA <b>MADE IN USA</b></p>										
<table border="1"> <tr> <td data-bbox="852 594 1112 772"> <p>84 in/ 213.4 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p> </td> <td data-bbox="1161 594 1258 741">  72 lbs/ 32.7 kg             </td> </tr> <tr> <td data-bbox="852 793 1112 972"> <p>80 in/ 203.2 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p> </td> <td data-bbox="1161 793 1258 940">  70 lbs/ 31.8 kg             </td> </tr> <tr> <td colspan="2" data-bbox="990 982 1177 1056">  500 lbs/ 226.8 kg             </td> </tr> </table>			<p>84 in/ 213.4 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p>	 72 lbs/ 32.7 kg	<p>80 in/ 203.2 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p>	 70 lbs/ 31.8 kg	 500 lbs/ 226.8 kg			
<p>84 in/ 213.4 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p>	 72 lbs/ 32.7 kg									
<p>80 in/ 203.2 cm</p>  <p>35 in/ 88.9 cm</p> <p>6 in/ 15.2 cm</p>	 70 lbs/ 31.8 kg									
 500 lbs/ 226.8 kg										
<p>APPROVED BY:</p> <p><b>MATTRESS OPTIONS</b></p> <p><input type="checkbox"/> 84 in LENGTH: 84 in x 35 in x 6 in (213.4 cm x 88.9 cm x 15.2 cm)</p> <p><input type="checkbox"/> 80 in LENGTH: 80 in x 35 in x 6 in (203.2 cm x 88.9 cm x 15.2 cm)</p> <p><input type="checkbox"/> SERVICE PART - 84 in LENGTH: 84 in x 35 in x 6 in (213.4 cm x 88.9 cm x 15.2 cm)</p> <p><input type="checkbox"/> SERVICE PART - 80 in LENGTH: 80 in x 35 in x 6 in (203.2 cm x 88.9 cm x 15.2 cm)</p>										
<p><b>WARNING</b></p> <p>PATIENT STABILITY AND SIDERAILS COVERAGE MAY BE COMPROMISED WITH THE USE OF AN OVERLAY</p> <p>DO NOT TRANSFER PATIENT FROM ONE BED TO ANOTHER USING THE SUPPORT SURFACE WITH A PATIENT ON IT</p>										
<p><b>CARE LABEL</b></p> <table border="0"> <tr> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> </tr> </table>										
										
										
<p>2860-030-905 C.0</p>										

# Quick Reference Replacement Parts List

---

These parts are currently available for purchase. Some of the parts depicted in the assembly drawing may not be available for purchase. Please call Stryker Customer Service at 1-800-327-0770 for availability and pricing.

English

<b>Part Names for 84" Support Surface</b>	<b>Part Number</b>
Cover Assembly 84"	2860-700-130
<b>Part Names for 80" Support Surface</b>	<b>Part Number</b>
Cover Assembly 80"	2860-700-135
<b>Common Components</b>	
Fire Barrier Sleeve	2850-035-001
Cap Screw	2860-030-422
Reduction Coupler Fitting	7600-001-404

[Return To Table of Contents](#)

# Warranty

---

English

Stryker Medical, a division of Stryker Corporation (“Stryker”), warrants that its **IsoGel AIR™** product will be free from defects in material and workmanship. This Stryker warranty covers only the following items of the Stryker **IsoGel AIR™** product (each known individually as a “Part” and collectively as the “Product” or “**IsoGel AIR™** product”) during normal use\* as follows:

- **Mattress (Foam Assembly and Gel) Warranty Period: 10 years**
  - *Note:* The Mattress will naturally compress over time. Should a body indentation or compression set measurement greater than 1 1/4" (3,2 cm), as measured by an authorized Stryker representative, be realized within the warranty period noted above, Stryker will provide a replacement Product. Any normal body indentations or compression set of less than 1 1/4" (3,2 cm) will not be replaced.
  - *Note:* Any damage to the Foam Assembly or Gel which results due to usage of a Cover Assembly beyond its warranty period of 3 years, or is a result of abnormal wear and tear which may include cleaning processes which are inconsistent with those recommended in this Operations/Maintenance Manual, shall invalidate the warranty on the Mattress at Stryker’s sole discretion.
- **Cover Assembly Warranty Period: 3 years**
- **Fire Barrier Sleeve Warranty Period: 3 years**

The above noted warranty periods apply only to the original purchaser of the **IsoGel AIR™** product and begin on the date of delivery to such original purchaser.

If Stryker determines, in its sole discretion, that one or more Parts is defective within the above noted warranty periods, then Stryker may, at its option, either repair or replace the **IsoGel AIR™** product or Part.

In addition, if requested by Stryker, the Part of the **IsoGel AIR™** product subject to a warranty claim shall be returned prepaid to Stryker as noted under the Return Authorization section below. No employee or representative of Stryker is authorized to change the warranty on the **IsoGel AIR™** product in any way.

The warranty set forth above does not include or cover the following:

- Abnormal wear and tear on the Product, or wear which indicates that the Product was not properly maintained in accordance with this Operations/Maintenance Manual, or which Product has been subject to unusual stress; or
- Product that has been misused, modified, refurbished or repaired without the prior written consent of Stryker; or
- Damage or Product failure due to causes beyond Stryker’s control, including but not limited to, abuse, theft, fire, flood, wind, lightning, freezing, clogging of mattress pores due to tobacco smoke, unusual atmosphere conditions, or material degradation due to exposure to moisture; or
- Damage which is determined to have resulted through the use of the Product for patient transfer or transport; or
- Product which serial numbers or other identification marks have been removed or destroyed.

\*\*“Normal use” is defined as use of the Product in typical or normal use settings in a hospital or medical facility under normal conditions. Damage to the Product which arises from abnormal use, which may include but is not limited to, damage to the Product that may be caused by needle punctures, burns, chemicals, negligent use or improper care or improper cleaning (proper cleaning to help sustain the life of the Product is as detailed in this Operations/Maintenance Manual) or staining resulting from such abnormal use(s) are exempt from the above note warranty coverage.

# Warranty

---

English

THE EXPRESS WARRANTY SET FORTH HEREIN IS THE ONLY WARRANTY APPLICABLE TO THE STRYKER **IsoGel AIR™** PRODUCT. ANY AND ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED BY STRYKER. IN NO EVENT SHALL STRYKER BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

## **TO OBTAIN PARTS AND SERVICE**

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative or call Stryker Customer Service USA at 1-800-327-0770.

## **RETURN AUTHORIZATION**

Product cannot be returned without prior approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

## **DAMAGED PRODUCT**

ICC Regulations require that claims for damaged product must be made with within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the Product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

[Return To Table of Contents](#)

## Stryker IsoGel AIR™ Cleaning Instructions

*Always follow hospital protocol for cleaning mattresses.*

### WARNING

- To avoid the risk of equipment damage, do not immerse support surface in cleaning or disinfectant solutions. Do not allow liquid to pool on the support surface.
- Do not iron, dry-clean, or tumble dry the support surface cover, as this will cause malfunction and damage to the product.
- Do not stick needles in to a medical bed mattress through the mattress cover. Holes will allow body fluids from entering the inside (inner core) of the mattress and could cause cross-contamination.
- To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
- To avoid the risk of patient injury, covers (top and bottom) should be inspected for tears, punctures, excessive wear, and misaligned zippers each time the covers are cleaned. If a support surface cover becomes compromised, the support surface should be removed from service immediately and replaced to prevent cross-contamination.

Stryker support surfaces are designed for long-lasting comfort and reliability. The life of the support surface can be adversely affected by an increase in frequency of usage which might include more frequent cleaning and disinfection.

1. Using a clean, soft, damp cloth, wipe down the entire support surface with a mild soap and water solution to remove foreign material.
2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or cleaning agent.
3. Care must be taken to thoroughly **rinse and dry** covers following cleaning.
4. Disinfect as necessary with a hospital grade disinfectant AFTER cleaning has been completed. Refer to “Disinfecting Support Surface” on [page 1-15](#).

### CAUTION

To avoid the risk of equipment damage, when cleaning the underside of the support surface, ensure that no liquid is allowed to seep into the zipper area and watershed cover barrier; fluids allowed to come in contact with the zipper may leak into the support surface.

## Stryker IsoGel AIR™ Disinfecting Instructions

*Always follow hospital protocol for cleaning mattresses.*

### WARNING

- To avoid the risk of patient injury, disinfect the support surface between patients. Failure to do so could result in cross-contamination and infection.
- Some disinfectants may cause damage to the product if used improperly. If the products described below are used to disinfect the support surface, measures must be taken to ensure the entire surface is wiped with a damp cloth soaked in clean water and thoroughly dried following disinfection. The cover can be damaged when exposed to such disinfectants beyond the manufacturers' recommendations. Failure to follow these directions when using these types of disinfectants may void this product warranty.
- To avoid the risk of shock, unplug unit prior to cleaning or disinfecting.
- To avoid the risk of equipment damage, the support surface cover must be completely dry before storage or adding linens. Failure to remove excess disinfectant could cause degradation of the cover material.

### Suggested Disinfectants:

- Phenolic Disinfectant
  - Chlorinated Bleach Solution (5.25% bleach diluted 1 part bleach to 100 parts water)
  - 70% Isopropyl Alcohol
1. Ensure surface has been thoroughly cleaned and dried prior to applying disinfectants.
  2. Wipe down the support surface with a clean, dry cloth to remove any excess liquid or disinfectant.
  3. Care must be taken to thoroughly **rinse and dry** covers following disinfection.

### CAUTION

- Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the cover fabric.
- The use of quaternaries containing glycol ethers and/or accelerated hydrogen peroxides may damage the cover.

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Stryker France S.A.S.