## AIR+™ Pump


Operations Manual


REF 2863

# Symbols

	Refer to instruction manual/booklet
	Operating instructions/Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR)
	Catalogue number
	Serial number
	Manufacturer
	Mass of product
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provisions for protective earthing or reliance upon installation conditions.
	Type B applied part – <b>SPR Plus</b> ® 2790-100-000 and <b>IsoFlex LAL</b> ® Support Surface 2860
	Fuse
	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1:2012 and CAN/CSA-C22.2 No. 60601-1:14.
<b>IP21</b>	Solids: Protection from touch by fingers and objects greater than 12 mm Liquids: Protection from the fall of vertical water drops
	For US Patents see <a href="http://www.stryker.com/patents">www.stryker.com/patents</a>

	Keep dry
	Do not stack more than 6 high
	This side up
	Fragile
	Do not use sharp objects to open the package

# Table of Contents

Warning/Caution/Note Definition .....	2
Summary of safety precautions .....	2
Introduction .....	3
Product description .....	3
Intended use .....	3
Expected service life .....	3
Contraindications .....	3
Specifications .....	4
Environmental conditions .....	4
Product illustration .....	5
Button and icon descriptions .....	5
Features .....	6
Contact information .....	6
Serial number .....	7
Date of manufacture .....	7
Setup .....	8
Connecting the power cord .....	8
Connecting the hose .....	8
Operation .....	9
Powering ON the pump .....	9
Increasing or decreasing pressure for SPR Plus .....	9
Activating CPR release .....	9
Turn the lock ON or OFF .....	9
Powering OFF the pump .....	10
Troubleshooting .....	11
Preventive maintenance .....	12
Product alarms .....	13
Cleaning and disinfecting with wipes .....	14
Cleaning .....	15
Disinfecting .....	16
EMC information .....	17

# Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

## **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

## **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

**Note** - Provides special information to make maintenance easier or important instructions clearer.

## Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

### **WARNING**

- Consult the operations manual for the support surface before you attach the pump.
- Always connect this product to a hospital grade receptacle. Failure to provide proper grounding could result in electrical shock.
- Never plug the power cord for this product into an extension cord.
- Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
- Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.
- Always unplug the product before you begin cleaning or disinfecting.
- Do not immerse the product in cleaning or disinfectant solutions.
- Do not spray or allow liquids to enter the power cord ports or the hose ports.
- Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **AIR+**, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

---

### **CAUTION**

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
  - Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
  - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
  - Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
-

# Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

---

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- 

## Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

## Product description

The Stryker Model 2863 **AIR+**<sup>™</sup> Pump is an electric pump for use with Model 2860 **IsoFlex LAL**<sup>®</sup> support surface or Model 2790-100-000 **SPR Plus**<sup>®</sup> overlay. This pump provides Low Air Loss (LAL) for **IsoFlex LAL** support surface and continuous air to the **SPR Plus** overlay. LAL provides air flow to help manage microclimate of the skin. The hooks on the back of this product allows you to attach the pump to the footboard of the compatible bed frames.

## Intended use

The **AIR+** pump when used with the **IsoFlex LAL** support surface or the **SPR Plus** overlay assists in the prevention and treatment of pressure injuries or pressure ulcers (all stages, Unstageable injury, and Deep tissue injury). We recommend you implement this product in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional.

This pump is for use in an acute care setting. This may include critical care, step down, progressive care, med/surg, sub-acute care, and post anesthesia care unit (PACU), or other locations as prescribed by a physician. Operators of this pump include healthcare professionals (such as caregivers, nurses, physicians, patient transporters, and other healthcare providers).

This product is not intended to be sterile, does not include a measuring function, and is not for use in a home healthcare environment.

## Expected service life

The **AIR+** pump has a 5 year expected service life under normal use, conditions, and with appropriate periodic maintenance.







## Contraindications

**AIR+** pump when used with **SPR Plus** overlay, air support therapy is not recommended when spinal stability is a concern.

## Specifications

Height	17.25 in.	43.8 cm	
Width	9.25 in.	23.5 cm	
Depth	7.25 in.	18.4 cm	
Weight	13.2 lb	6.0 kg	
Interface control panel	LED		
Input voltage	120 V		
Electric current	0.6 A		
Input frequency	60 Hz		
Peak noise level	<50 dba		
Mode of operation	Continuous		
Power cord	Short	3.2 ft	1 m
	Long	16.4 ft	5 m
Compliance	IEC 60601-1 3.1 Edition; RoHS Directive 2002/95/EC Reach; IEC 60601-1-8:2012, IEC 60601-1-2:2014; CAN/CSA C22.2 No. 601.1		
Compatible with bed frame	<b>3002 S3™, 3005 S3™, InTouch®, Spirit Select™</b>		


## Environmental conditions

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

## Product illustration

A	User interface	F	Power inlet
B	Bumpers	G	Anti-vibration foot
C	Transport handle	H	Hose port
D	Hooks	I	Back pump housing
E	HEPA filter	J	Front pump housing

## Button and icon descriptions

Icon	Name	Description
	Power button on, green	Press once to power ON
	Power button off	Press and hold 2 seconds to power OFF
	Lock off icon	Lock is off and allows you to change settings
	Lock on icon	Lock on locks the screen and does not allow you to change settings

Icon	Name	Description
	Decrease or increase firmness buttons	Press down to decrease or press up to increase firmness for <b>SPR Plus</b>
18, 21, 24, 27, 30	Firmness level icons	5 levels to set firmness in mmHg. Icons are green when active.
	Pressure alarm on icon	Pressure alarm on when pressure is too low or too high
	Pressure alarm off icon	Pressure alarm off
	<b>SPR Plus</b> button	Activates <b>SPR Plus</b> and allows selection of the firmness levels when selected
	<b>IsoFlex LAL</b> button	Activates Low Air Loss for <b>IsoFlex</b> when selected

## Features

Low Air Loss (LAL)

Provides air flow to help manage microclimate of the skin.

## Contact information


Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

## Serial number

## Date of manufacture

The year of manufacture is the first four digits of the serial number.

# Setup

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

---

## WARNING

- Consult the operations manual for the support surface before you attach the pump.
  - Always connect this product to a hospital grade receptacle. Failure to provide proper grounding could result in electrical shock.
  - Never plug the power cord for this product into an extension cord.
  - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
- 

## CAUTION

- Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
  - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
- 
1. Bring the product to room temperature.
  2. Place the pump in a secure location. Allow space around the pump to connect and disconnect the hoses and power cord. Optional placements:
    - a. Option 1: Hang the pump on a compatible bed. Make sure that the pump is secure.
    - b. Option 2: Place the pump on the floor. Make sure that the pump is away from foot traffic to avoid the risk of a trip hazard.

## Connecting the power cord

1. Attach the power cord to the pump.
2. Plug the power cord into an accessible receptacle.

## Connecting the hose

1. Attach the hose to the pump hose port.
2. Attach the hose to the support surface hose port.

# Operation

## Powering ON the pump

1. See *Setup* (page 8) to prepare the pump for operation.
2. Press the power button once to turn on the product. (Figure 1).

Figure 1 – Status view, operator view

**Note** - The pump will power ON with either **SPR Plus** or **IsoFlex LAL** based on last use.

## Increasing or decreasing pressure for SPR Plus

You can increase or decrease the pressure by 3 mmHg increments when the **SPR Plus** mode is selected.

To increase pressure, press the up arrow.

Figure 2 – Increase or decrease

To decrease pressure, press the down arrow.

## Activating CPR release

Unscrew the hose from the **SPR Plus** overlay to deflate the support surface for CPR.

**Note** - You do not need to remove the hose from **IsoFlex LAL** support surface to perform CPR.

## Turn the lock ON or OFF

To turn the ON or OFF.

1. Press and hold the lock button for 1 second to turn ON the lock.

**Figure 3 – Lock ON**

**Note** - Lock ON locks the screen and does not allow you to change settings.

2. Press and hold the lock button for 2 seconds to turn OFF the lock.

**Figure 4 – Lock OFF**

## Powering OFF the pump

Press the power button and hold for 2 seconds to power OFF the pump.

# Troubleshooting

Problem	Screen	Cause	Recommended action
Power loss, product does not turn ON		Power cord not seated, power cord unplugged from outlet, or possible internal damage	Make sure that the power cord is plugged into the product and the outlet.
Button not responsive		Lock is ON	Press and hold the lock button for 2 seconds to turn lock OFF.
Alarm icon is ON		Pressure alert	Check air hose for kink. Check the hose for connection. If not resolved, call service.

# Preventive maintenance

**WARNING** - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.

## Note

- At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage.
- Remove product from service before you perform preventive maintenance.
- Consult your local regulations to dispose of electronic equipment.

Inspect the following items:

- All fasteners are secure
- Pump housing or components (hose, power cords, or case) for cracks, holes, or damaged
- Hooks that hang the pump on the bed frame are not damaged
- Check for secure hose connections
- User interface is not cracked or damaged
- HEPA filter (replace each year)

Product serial number:
Completed by:
Date:

# Product alarms

The product alarm for this product is display only, no audible alarms.

Icon	Priority alarm	Name	Delay	Therapy interrupted
	Low	Pressure alarm	15 minutes	Yes

## Note

- Alarm displays on the interface control panel of the pump.
- Alarm resets if there is power loss to the pump.

## Cleaning and disinfecting with wipes

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred wipes (2060-000-001 6" x 10" or 2060-000-002 9" x 12") include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

**Note** - For safety information, read the product label.

To clean or disinfect the external product surface:

1. To clean, wipe external surfaces with a fresh, clean wipe to remove all visible soils. Repeat as necessary until the product is clean.

**Note**

- Use as many wipes as necessary.
  - Complete step 1 before you disinfect.
2. To disinfect, wipe external surfaces with a fresh, clean wipe until wet. Allow the external surface to remain wet for two minutes at room temperature.
  3. Allow the product to dry before you return it to service.

# Cleaning

---

## WARNING

- Always unplug the product before you clean or disinfect.
  - Do not immerse the product in cleaning or disinfectant solutions.
  - Do not spray or allow liquids to enter the power cord ports or the hose ports.
  - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- 

## CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
- 

Recommended cleaning agents include:

- A clean damp cloth and mild soapy water

To clean the external side of the pump:

1. Unplug the controller.
2. Wipe the outside of the product with a clean, soft cloth moistened with a mild soap and water solution to remove foreign material.
3. Follow specified contact time in accordance with the manufacturer's instructions for use.
4. Using a clean, soft cloth, wipe the outside of the product with water.
5. Dry the external surface with a clean, dry cloth.

# Disinfecting

---

## WARNING

- Always unplug the product before you clean or disinfect.
  - Do not immerse the product in cleaning or disinfectant solutions.
  - Do not spray or allow liquids to enter the power cord ports or the hose ports.
  - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- 

## CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
- 

Recommended chemical solutions:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Accelerated hydrogen peroxide
- Potassium peroxymonosulfate
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients to avoid the risk of cross-contamination and infection.

**Note** - Make sure that you unplug the product before you clean or disinfect.

## EMC information

**WARNING** - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

### Note

- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- This equipment is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Guidance and manufacturer's declaration - electromagnetic emissions		
The 2863 <b>SPR+</b> pump is intended for use in the electromagnetic environment specified below. The customer or the user of the 2863 <b>SPR+</b> pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The 2863 <b>AIR+</b> pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The 2863 <b>AIR+</b> pump is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

### WARNING

- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **AIR+**, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

Recommended separation distances between portable and mobile RF communication equipment and the 2863 <b>AIR+</b> pump			
The 2863 <b>AIR+</b> pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 2863 <b>AIR+</b> pump can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 2863 <b>AIR+</b> pump as recommended below, according to the maximum output power of the communications equipment.			
Band (MHz)	Service	Maximum Power (W)	Minimum Separation Distance (m)

Recommended separation distances between portable and mobile RF communication equipment and the 2863 AIR+ pump			
380-390	TETRA 400	1.8	0.3
430-470	GMRS 460; FRS 460	2.0	0.3
704-787	LTE Band 13, 17	0.2	0.3
800-960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	2.0	0.3
1,700-1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	2.0	0.3
2,400-2,570	Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE Band 7	2.0	0.3
5,100-5,800	WLAN 802.11 a/n	0.2	0.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note** - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The 2863 AIR+ pump is suitable for use in the electromagnetic environment specified below. The customer or the user of the 2863 AIR+ pump should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	$\pm 8$ kV contact $\pm 15$ kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV for input/output lines	$\pm 1$ kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.

**Guidance and manufacturer's declaration - electromagnetic immunity**

<p>Voltage dips, voltage variations and short interruptions on power supply input lines</p> <p>IEC 61000-4-11</p>	<p>0%U<sub>T</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0%U<sub>T</sub> for 1 cycle</p> <p>70%U<sub>T</sub> (30% dip in U<sub>T</sub>) for 25/30 cycles</p> <p>0% U<sub>T</sub> for 250/300 cycles</p>	<p>0%U<sub>T</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0%U<sub>T</sub> for 1 cycle</p> <p>70%U<sub>T</sub> (30% dip in U<sub>T</sub>) for 25/30 cycles</p> <p>0% U<sub>T</sub> for 250/300 cycles</p>	<p>Main power quality should be that of a typical commercial or hospital environment. If the user of the 2863 <b>AIR+</b> pump requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field</p> <p>IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

**Note** - U<sub>T</sub> is the a.c. mains voltage before applications of the test level.

<p>Conducted RF IEC 61000- 4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.7 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the 2863 <b>AIR+</b> pump, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>D=(1.2) (\sqrt{P})</math> 80 MHz to 800 MHz</p> <p><math>D=(2.3) (\sqrt{P})</math> 800 MHz to 2.7 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
---	--	----------------------------	--

**Note**

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.
- The ISM (Industrial, Scientific, and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 2863 **AIR+** pump is used exceeds the applicable RF compliance level above, the 2863 **AIR+** pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 2863 **AIR+** pump.



<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 Vrms.

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