AIR+ Pump

Operations Manual

REF 2863

Global symbol glossary

See the Global Symbol Glossary at ifu.stryker.com for symbol definitions.

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR)
	Catalogue number
	Serial number
	Manufacturer
	Date of manufacture
	Unique device identifier
	Mass of product
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provisions for protective earthing or reliance upon installation conditions.
	Type B applied part – SPR Plus ® 2790-100-000 and IsoFlex LAL ® Support Surface 2860
	Fuse

 E347196 MEDICAL ELECTRICAL EQUIPMENT	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1:2012 and CAN/CSA-C22.2 No. 60601-1:14.
IP21	Solids: Protection from touch by fingers and objects greater than 12 mm Liquids: Protection from the fall of vertical water drops
US Patents	For US Patents see www.stryker.com/patents
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	3
Product description	3
Intended use	3
Contraindications	3
Expected service life	3
Specifications	4
Environmental conditions	4
Product illustration	5
Button and icon descriptions	5
Features	6
Contact information	6
Serial number location	7
Date of manufacture	7
Setup	8
Connecting the power cord	8
Connecting the hose	8
Operation	10
Powering on the pump	10
Increasing or decreasing pressure for SPR Plus	10
Activating CPR release	10
Turning the lock on or off	10
Powering off the pump	11
Accessories and parts	12
Troubleshooting	13
Preventive maintenance	14
Product alarms	15
Cleaning and disinfecting with SideKick	16
Cleaning	17
Disinfecting	18
EMC information	19

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always consult the operations manual for the support surface before you attach the pump.
 - Always connect this product to a hospital grade receptacle. Failure to provide proper grounding could result in electrical shock.
 - Never plug the power cord for this product into an extension cord.
 - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
 - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.
 - Always unplug the product before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports or the hose ports.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
 - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
 - Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **AIR+**, including cables specified by the manufacturer.
 - Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
 - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
 - Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 2863 **AIR+** Pump is an electric pump for use with Model 2860 **IsoFlex LAL** support surface or Model 2790-100-000 **SPR Plus** overlay. This pump provides Low Air Loss (LAL) for **IsoFlex LAL** support surface and continuous air to the **SPR Plus** overlay. LAL provides air flow to help manage microclimate of the skin. The hooks on the back of this product allows you to attach the pump to the footboard of the compatible bed frames.

Intended use

The **AIR+** pump when used with the **IsoFlex LAL** support surface or the **SPR Plus** overlay assists in the prevention and treatment of pressure injuries or pressure ulcers (all stages, Unstageable injury, and Deep tissue injury). We recommend you implement this product in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional.

This pump is for use in an acute care setting. This may include critical care, step down, progressive care, med/surg, sub-acute care, and post anesthesia care unit (PACU), or other locations as prescribed by a physician. Operators of this pump include healthcare professionals (such as caregivers, nurses, physicians, patient transporters, and other healthcare providers).

This product is not intended to be sterile, does not include a measuring function, and is not for use in a home healthcare environment.

Contraindications

AIR+ pump when used with **SPR Plus** overlay, air support therapy is not recommended when spinal stability is a concern.







Expected service life

The **AIR+** pump has a five year expected service life under normal use conditions and with appropriate periodic maintenance.

Specifications


Height		17.25 in.	43.8 cm
Width		9.25 in.	23.5 cm
Depth		7.25 in.	18.4 cm
Weight		13.2 lb	6.0 kg
Interface control panel		LED	
Input voltage		120 V	
Electric current		0.6 A	
Input frequency		60 Hz	
Peak noise level		<50 dba	
Mode of operation		Continuous	
Power cord	Short	3.2 ft	1 m
	Long	16.4 ft	5 m
Compliance		IEC 60601-1 3.1 Edition; RoHS Directive 2002/95/EC Reach; IEC 60601-1-8:2012, IEC 60601-1-2:2014; CAN/ CSA C22.2 No. 601.1	
Compatible with bed frame		3002 S3 , 3005 S3 , InTouch , Spirit Select	

Environmental conditions

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

Product illustration

A	User interface
B	Front pump housing
C	Bumper
D	Transport handle
E	Hook
F	Power inlet

G	Anti-vibration foot
H	Air filter
I	Flex grip
J	Fuse
K	Back pump housing
L	Hose port

Button and icon descriptions

Icon	Name	Description
	Power button on, green	Press once to power on
	Power button off	Press and hold for two seconds to power off
	Lock off button	Lock off unlocks the screen and allows you to change settings
	Lock on button	Lock on locks the screen and does not allow you to change settings

Icon	Name	Description
	Decrease or increase firmness buttons	Press down to decrease or press up to increase firmness for SPR Plus
18, 21, 24, 27, 30	Firmness level icons	Five levels to set firmness in mmHg. Icons are green when active
	Pressure alarm on icon	Pressure alarm on when pressure is too low or too high
	Pressure alarm off icon	Pressure alarm off
	SPR Plus button	Activates SPR Plus and allows you to select the firmness levels when selected
	IsoFlex LAL button	Activates Low Air Loss for IsoFlex LAL when selected

Features

Low Air Loss (LAL)

Provides air flow to help manage microclimate of the skin.

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location


Figure 1 – Serial number location

Date of manufacture

The year of manufacture is the first four digits of the serial number.

Setup

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

WARNING

- Always consult the operations manual for the support surface before you attach the pump.
 - Always connect this product to a hospital grade receptacle. Failure to provide proper grounding could result in electrical shock.
 - Never plug the power cord for this product into an extension cord.
 - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
-

CAUTION

- Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
 - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
-

Connecting the power cord

To connect the power cord:

1. Bring the product to room temperature.
2. Attach the power cord to the pump.
3. Insert the power cord under the power cord retention flex grip (A) (Figure 2).


Figure 2 – Cord retention flex grip

4. Plug the pump power cord into an accessible wall outlet.
5. Place the power cord out of the path of foot traffic to avoid trip hazards.
6. Place the pump in a secure location. Allow space around the pump to connect and disconnect the hoses and power cord.
 - a. Option 1: Hang the pump on a compatible bed. Make sure that the pump is secure.
 - b. Option 2: Place the pump on the floor.

Connecting the hose

To connect the hose:

1. Attach the hose to the pump hose port.
2. Attach the hose to the support surface hose port.

Operation

Powering on the pump

To power on the pump:

1. See *Setup* (page 8) to prepare the pump for operation.
2. Press the **Power** button once to turn on the product (Figure 3).

Figure 3 – Status view, operator view

Note - The pump will power on with either **SPR Plus** or **IsoFlex LAL** selected based on the last use.

Increasing or decreasing pressure for SPR Plus

You can increase or decrease the pressure by 3 mmHg increments when you select the **SPR Plus** mode (Figure 4).

To increase pressure, press the **Increase firmness** button.

Figure 4 – Increase or decrease

To decrease pressure, press the **Decrease firmness** button.

Activating CPR release

Unscrew the hose from the **SPR Plus** overlay to deflate the support surface for CPR.

Note - You do not need to remove the hose from **IsoFlex LAL** support surface to perform CPR.

Turning the lock on or off

To lock or unlock product functions:

1. Press and hold the **Lock on** button for one second to turn on the lock (Figure 5).

Figure 5 – Lock on

Note - The **Lock on** button locks the screen and does not allow you to change settings.

2. Press and hold the **Lock off** button for two seconds to turn off the lock (Figure 6).

Figure 6 – Lock off

Powering off the pump

Press and hold the **Power** button for two seconds to power off the pump.

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Part	Number
Power cord, type B, 1 meter	2874-007-001
Power cord, type B, 5 meter	2874-007-002
Internal hose lengths	2863-007-007
Wire harness	2874-007-035
Fuse	2874-007-021
Air filter	2874-007-026

Troubleshooting

Problem	Screen	Cause	Recommended action
Power loss, product does not turn on		Power cord not seated, power cord unplugged from the wall outlet, or possible internal damage	1. Make sure that the power cord is plugged into the product and a wall outlet. 2. Make sure that the power cord is routed through the flex grip.
Button not responsive		Lock is on	Press and hold the Lock button for two seconds to turn lock off.
Alarm icon is on		Pressure alert	1. Make sure there are no kinks in the air hose. 2. Make sure that the hose is connected.
		Mode selected does not match the connected product. For example, if you selected SPR Plus mode but you are connected to an IsoFlex LAL support surface.	Make sure that you are connected to the correct product for the mode selected.
		Mode or pressure setting changed	1. Wait for the pressure to stabilize. The error will clear when the pump reaches the new target pressure. 2. Power cycle the unit. 3. If not resolved, call service.

Preventive maintenance

WARNING - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.

Remove the product from service before you perform preventive maintenance. At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Note - Consult your local regulations to dispose of electronic equipment.

Inspect the following items:

- _____ All fasteners are secure
- _____ Pump housing or components (hose, power cords, or case) are free of cracks, holes, or damage
- _____ Bed frame pump hooks are free of damage
- _____ Hose connections are secure
- _____ User interface are free of cracks or damage
- _____ Fuse
- _____ Air filter (replace each year)
- _____ Test functionality (all tests must pass)

Product serial number:
Completed by:
Date:

Product alarms

The product alarm for this product is display only. There are no audible alarms.

Icon	Priority alarm	Name	Delay	Therapy interrupted
	Low	Pressure alarm	15 minutes	Yes

Note

- Alarm displays on the interface control panel of the pump.
- Alarm resets if there is power loss to the pump.

Cleaning and disinfecting with SideKick

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred 6" x 10" wipes (2060-000-001) include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface with SideKick wipe:

To clean:

1. Wipe down the external product surface with a fresh, clean wipe to remove all visible soils.
2. Repeat as necessary until the external product surface is visibly clean.
3. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Note - Use as many wipes as necessary.

To disinfect:

1. Clean first.
2. Wipe down the external product surface with a fresh, clean wipe until wet.
3. Allow the external product surface to remain wet for two minutes at room temperature.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Cleaning

WARNING

- Always unplug the product before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports or the hose ports.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
-

CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
-

To clean the external side of the pump:

1. Unplug the power cord from the wall outlet.
2. Wipe the outside of the product with a clean, soft cloth moistened with a mild soap and water solution to remove foreign material.
3. Follow specified contact time in accordance with the manufacturer's instructions for use.
4. Using a clean, soft cloth, wipe the outside of the product with water.
5. Dry the external surface with a clean, dry cloth.

Disinfecting

WARNING

- Always unplug the product before you clean or disinfect.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not spray or allow liquids to enter the power cord ports or the hose ports.
 - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
-

CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
 - Failure to follow manufacturing instructions may also affect useful life of the product.
-

Suggested disinfectants:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Accelerated hydrogen peroxide
- Potassium peroxymonosulfate
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients to avoid the risk of cross-contamination and infection.

To disinfect the pump:

1. Unplug the power cord from the wall outlet.
2. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes. Do not soak the product.

Note - Follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.

3. Rinse the product with a clean, damp cloth to remove excess disinfectant.
4. Dry the product with a clean, dry cloth to remove excess liquid or disinfectant.
5. Allow the product to dry before you return to service.
6. Insert the power cord under the power cord retention flex grip.
7. Plug the pump power cord into a wall outlet.

EMC information

WARNING

- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **AIR+**, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

Note

- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- This equipment is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Guidance and manufacturer's declaration - electromagnetic emissions		
The 2863 AIR+ pump is intended for use in the electromagnetic environment specified below. The customer or the user of the 2863 AIR+ pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The 2863 AIR+ pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The 2863 AIR+ pump is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	N/A
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Recommended separation distances between portable and mobile RF communication equipment and the 2863 AIR+ pump			
The 2863 AIR+ pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 2863 AIR+ pump can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 2863 AIR+ pump as recommended below, according to the maximum output power of the communications equipment.			
Band (MHz)	Service	Maximum Power (W)	Minimum Separation Distance (m)
380-390	TETRA 400	1.8	0.3
430-470	GMRS 460; FRS 460	2.0	0.3
704-787	LTE Band 13, 17	0.2	0.3

Recommended separation distances between portable and mobile RF communication equipment and the 2863 AIR+ pump			
800-960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	2.0	0.3
1,700-1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	2.0	0.3
2,400-2,570	Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE Band 7	2.0	0.3
5,100-5,800	WLAN 802.11 a/n	0.2	0.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The 2863 AIR+ pump is suitable for use in the electromagnetic environment specified below. The customer or the user of the 2863 AIR+ pump should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/ output lines	± 2 kV for power supply lines ± 1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV for input/output lines	± 1 kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T for 1 cycle 70% U_T (30% dip in U_T) for 25/30 cycles 0% U_T for 250/300 cycles	0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T for 1 cycle 70% U_T (30% dip in U_T) for 25/30 cycles 0% U_T for 250/300 cycles	Main power quality should be that of a typical commercial or hospital environment. If the user of the 2863 AIR+ pump requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note - U_T is the a.c. mains voltage before applications of the test level.			

Conducted RF IEC 61000- 4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 2863 AIR+ pump, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance $D=(1.2) (\sqrt{P})$ 80 MHz to 800 MHz $D=(2.3) (\sqrt{P})$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note • At 80 MHz and 800 MHz, the higher frequency range applies. • These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people. • The ISM (Industrial, Scientific, and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.			
^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 2863 AIR+ pump is used exceeds the applicable RF compliance level above, the 2863 AIR+ pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 2863 AIR+ pump. ^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 Vrms.			

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **AIR+**, **InTouch**, **IsoFlex**, **LAL**, **S3**, **Sidekick**, **Spirit Select**, **SPR Plus**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.