
IsoTour Gel Support Surface

Operations Manual

Standard and premium

REF 2872


2872-009-001

EN
FR

Global symbol glossary

See the Global Symbol Glossary at ifu.stryker.com for symbol definitions.

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Catalogue number
	Serial number
	Batch code
	Unique device identifier
	Quantity
	Authorized representative in the European Community
	European medical device
	Manufacturer
	Date of manufacture
	Importer
	Type B applied part

	Mass of equipment
	Safe working load
	Wash by hand
	Do not tumble dry
	Do not dry-clean
	Do not iron
	Allow to air dry
	Chlorinated bleach
	CE Mark
	For US Patents see www.stryker.com/patents
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package
	Pressure redistribution
	Shear
	TruTurn


	Low air loss
	Sacral off-loading

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Intended use	4
Indications for use	4
Intended users	4
Clinical benefits	4
Clinical assessment	5
Contraindications	5
Expected service life	5
Disposal/recycle	5
Specifications	6
Features	7
Contact information	7
Serial number location	7
Setup	8
Set up the six inch foot extender pad option	8
Set up the twelve inch foot extender pad option	9
Operation	10
Applying linens	10
Transferring a patient from one patient support platform to another	10
Positioning a patient for TruTurn	11
Accessories and parts	12
Date of manufacture for medical device accessories	12
Preventive maintenance	13
Cleaning and disinfecting with SideKick	14
Care and maintenance	15
Care of bottom cover	16
Top cover replacement	18
Bottom cover replacement	19
Fire barrier replacement	20
Spacer fabric replacement, premium option only	21

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always check the patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
- Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
- Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
- Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
- Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
- Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness of the support surface. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
- Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
- Never transfer or transport a patient while **TruTurn** is in use.
- Do not use the support surface as a transfer device.
- Do not use the support surface handles to lift or move the support surface with a patient on board.
- Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
- Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
- Always center the patient on the support surface. Align the patient's head toward the headboard before you start TruTurn. Check the patient frequently to make sure that you maintain the proper position.
- Always make sure that tubes or wires connected to the patient are long enough, stable, and secure when you use TruTurn.
- Do not exceed the safe working load of this product. Excess weight could cause unpredictable safety and performance of this system.

- Do not wash the internal components of this support surface. Discard the support surface if contamination is found inside.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not allow liquid to pool on the product.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you use cleaning agents. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these care and maintenance instructions may void your warranty.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
 - Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
 - Do not allow sharp edges from the X-ray plate to come into contact with the support surface cover. Stryker recommends that you cover the X-ray plate with a pillow case or other protective material before you place the X-ray plate under the patient. If the support surface cover is damaged, remove it from service immediately to prevent cross-contamination.
 - Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 2872 IsoTour is a gel-filled support surface that is for use non-powered or powered. The support surface is available in a standard version and a premium version.

The Stryker Model 2874 IsoTour Pump is an electrical pump for use with the Stryker Model 2872 IsoTour Support Surface. The pump connects to the support surface with an air hose. Air is delivered through the integrated perforated tubing system within the support surface.

The TruTurn function is available when the pump and premium support surface are connected. The pump inflates the integrated turn bladders in the bottom cover. When the pump is turned on, the Low Air Loss (LAL) function is present for the standard version or premium version of the support surface.

Intended use

The IsoTour gel support surface is intended to assist in the prevention and treatment of pressure injuries [stages (1, 2, 3, 4), unstageable, and deep tissue injury or all pressure injuries] and is recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a health care professional.

Indications for use

The IsoTour gel support surface is indicated for use with human patients to assist in the prevention and treatment of pressure injuries (stages 1, 2, 3, 4, unstageable, and deep tissue injury or all pressure injuries) in an acute care setting within a healthcare facility and is recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by a health care professional.

Intended users

Operators of this product include healthcare professionals (such as nurses, nurse aids, or doctors).

This product is not intended to be sterile, include a measuring function, or be used in a home healthcare environment.

Clinical benefits

Assists in the prevention and treatment of all pressure ulcers or pressure injuries

Clinical assessment

The TruTurn feature is for use only after you complete a clinical assessment per your hospital protocol. This is to make sure that no adverse reaction would result for patients with the following example conditions:

- Spinal cord injury, except for unstable spine, see Contraindications
- Skeletal traction
- Significant hemoptysis
- Bleeding disorders
- Fractures
- Hemodynamically unstable
- Increased intracranial pressure

Contraindications

The TruTurn feature is contraindicated for use with patients that have unstable spine.

Note - TruTurn requires the IsoTour pump to be connected for use. The TruTurn feature is visually identified by the icons on the IsoTour pump (Figure 1).


Figure 1 – TruTurn icon

Expected service life

The product expected life is for normal use, conditions, and with appropriate periodic maintenance.

- Standard top cover (END406) - 5 years
- Bottom cover, standard or premium - 5 years
- Mattress core assembly, spacer fabric, and fire barrier - 10 years






Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

 Safe working load Note - Safe working load indicates the sum of the patient, mattress, and accessory weight.	500 lb	227 kg
Minimum patient weight	50 lb	22.7 kg

Standard	Standard top cover (END406)	2872-000-012		2872-000-016	
		2872-000-020			
Premium		2872-000-011		2872-000-015	
		2872-000-019			
Length		84.25 in.	214 cm	78.75 in.	200 cm
Width		35.5 in.	90.2 cm	35.5 in.	90.2 cm
Thickness		7 - 9.5 in.	17.8 - 24.1 cm	7 - 9.5 in.	17.8 - 24.1 cm
Product weight		85 lb	38.6 kg	79.9 lb	36.2 kg
Compatible bed frames for standard and premium support surfaces		InTouch, ProCuity, 3002 S3, 3005 S3, Spirit Select		N/A	
IsoTour pump		2874-000-001			
Flammability compliance	USA	CALIFORNIA TB 129, BOSTON-BFD IX-11, 16 CFR 1632			
	International	EN 597-1, EN 597-2, CANADA-Method 27.7 of CAN 2-4.2 M77, ITALY UNI 9175			

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

Features

TruTurn	Available for the premium support surface. Air bladders inflate to rotate up to 30 degrees. Helps to off-load the sacrum. Turn angle can be customized per patient tolerance and comfort.
Moisture management	The top surface cover is designed to draw excess moisture away from the patient's skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to help manage microclimate of the skin.
Standard support surface	Orange bottom cover, gel and foam core, LAL capable when connected to pump, handles, and D-rings
Premium support surface	Blue bottom cover, gel and foam core, spacer fabric for enhanced air flow, handles, D-rings, anchor straps, TruTurn and LAL capable when connected to pump

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

Setup

WARNING

- Always check the patient's skin regularly. Consult a physician if erythema or skin breakdown occurs. Serious injury could result if the patient's skin condition is left untreated.
 - Always use extra caution and supervision to help reduce the risk of a patient fall. Patient stability and siderail coverage may be compromised with the use of an overlay.
 - Always inspect for foreign objects between the support surface and the bed frame. Foreign objects may cause the support surface to slide on the support platform.
 - Always consider the use of siderails. The safe use of the support surface is maximized when used in conjunction with siderails. There may be an increased risk of falls when siderails are not present. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of siderails or other restraints. Consider local policies regarding the use of siderails. The physician, operator, or responsible parties should determine whether and how to use siderails based on each patient's individual needs.
 - Always use extra caution with a patient at risk of a fall (such as agitated or confused) to help reduce the likelihood of a fall.
 - Do not use the support surface on a larger or smaller bed frame that does not fit the width, length, or thickness of the support surface. This is to avoid the risk of the support surface sliding, patient injury, or interference with moving parts of the bed.
 - Do not use the support surface when gaps are present. The risk of entrapment can develop when the support surface is placed on incompatible bed frames.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
-

CAUTION

- Always be aware of devices or equipment that are placed on the top of the support surface. Damage to the surface may occur due to the weight of the equipment, heat generated by the equipment, or sharp edges on the equipment.
 - Do not put overlays or accessories inside the cover to avoid the risk of reducing pressure redistribution performance.
-

1. Remove all shipping and packaging materials from the product before use.
2. Place the support surface with the artwork toward the foot end of the bed frame. The artwork is located on the side of the top cover.
3. Place the support surface between the mattress retainers on the bed frame.
 - a. If you have a standard support surface without anchor straps, go to step 4.
 - b. For a premium support surface, set up the anchor strap.
 - c. Remove the headboard and raise the backrest to 30 degrees or a workable height. You need access to the underside of the deck.
 - d. Start with one side of the anchor strap and wrap it around the top edge of the litter deck. Push the buckle up through the slot on the same side.
 - e. Fasten the buckle and pull the strap to tighten the strap to the bed frame. Keep the straps loose enough to allow for the support surface to turn.
 - f. Repeat steps 3d to 3e to attach the anchor strap on the other side.
4. Place linens on the support surface per hospital protocols.

Set up the six inch foot extender pad option

See the bed operations manual for instructions about how to set up the foot extender.

1. Remove all shipping and packaging materials from the product before use.
2. Place the foot extender pad option onto the bed extender. Make sure that you have the appropriate size foot extender pad.

Set up the twelve inch foot extender pad option

See the bed operations manual for instructions about how to set up the foot extender.

1. Remove all shipping and packaging materials from the product before use.
2. Place the foot extender pad option onto the bed extender. Make sure that you have the appropriate size foot extender pad.
3. Align the hook and loop fastener on the foot extender pad with the hook and loop fastener on the bed frame extenders.

Note - The bed frame extender may not have hook and loop fastener applied. Apply or replace hook and loop fastener if necessary.

Operation

Applying linens

To apply linens:

1. Part the D-rings (C) (Figure 2).
2. Thread the four linen corners through the D-rings (C) attached to the bottom cover to secure the linens (A) to the support surface (B).


Figure 2 – Apply linens

Note - Do not pull the linens tight. Keep the linens loose and as smooth as possible on top of the support surface for effective use of the **TruTurn**.

Transferring a patient from one patient support platform to another

WARNING

- Never transfer or transport a patient while **TruTurn** is in use.
 - Do not use the support surface as a transfer device.
 - Do not use the support surface handles to lift or move the support surface with a patient on board.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface.
 - Do not exceed the safe working load of the hospital bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
 - Always make sure that the opposite siderail is raised when placing a patient on the support surface to reduce the risk of patient fall.
-

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes to **ON** for both patient support platforms.
3. Adjust the patient support platform heights so that they are level with one another.
4. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.

Positioning a patient for TruTurn

Before you start TruTurn, you need to center the patient on the support surface.

WARNING

- Always center the patient on the support surface. Align the patient's head toward the headboard before you start TruTurn. Check the patient frequently to make sure that you maintain the proper position.
 - Always make sure that tubes or wires connected to the patient are long enough, stable, and secure when you use TruTurn.
 - Do not exceed the safe working load of this product. Excess weight could cause unpredictable safety and performance of this system.
-

CAUTION

- Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
 - Do not allow sharp edges from the X-ray plate to come into contact with the support surface cover. Stryker recommends that you cover the X-ray plate with a pillow case or other protective material before you place the X-ray plate under the patient. If the support surface cover is damaged, remove it from service immediately to prevent cross-contamination.
-

Note - When you take radiological images of the patient on this product, the internal components of this product may cause artifacts or distort images.

To position the patient:

1. Position the patient in the center of the support surface.
2. Align the patient's head toward the headboard.
3. Check the patient as you use TruTurn for the correct position (Figure 3).

Figure 3 – Correct patient position

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Name	Number	Date of manufacture
Bed extender pad, gel, 6 in. (15.2 cm)	2872-000-800	SN
Bed extender pad, gel, 12 in. (30.4 cm)	2872-000-810	SN
Kit, premium upgrade, 84.25 in. (214 cm)	2872-700-005	
Kit, premium upgrade, 78.75 in. (200 cm)	2872-700-006	
IsoTour pump	2874-000-001	SN
Fire barrier	2872-030-701	N/A
Bottom cover assembly, premium 78.75 in. (200 cm)	2872-001-002	N/A
Bottom cover assembly, premium 84.25 in. (214 cm)	2872-001-001	N/A
Bottom cover assembly, standard 78.75 in. (200 cm)	2872-001-004	N/A
Bottom cover assembly, standard 84.25 in. (214 cm)	2872-001-003	N/A
Spacer fabric 78.75 in. (200 cm), premium option only	2872-035-702	N/A
Spacer fabric 84.25 in. (214 cm), premium option only	2872-030-702	N/A
Standard top cover assembly, 78.75 in. (200 cm)	2872-045-500	N/A
Standard top cover assembly, 84.25 in. (214 cm)	2872-040-500	N/A

Date of manufacture for medical device accessories

The year of manufacture is the first four digits of the serial number.

Preventive maintenance

Remove product from service before you perform preventive maintenance. Check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Note - Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

- _____ Zipper and covers (top and bottom) are free of tears, cuts, holes, or other openings
- _____ Unzip the covers to inspect the internal components for signs of stains from fluid ingress or contamination
- _____ Labels for legibility, proper adherence, and integrity
- _____ Handles are free of rips, cracks, and are intact
- _____ Foam has not degraded or come apart
- _____ Fire barrier for rips or other visible signs of damage
- _____ Spacer fabric, premium option only, for rips or other visible signs of damage
- _____ Run a diagnostic test if paired with a pump - see the **IsoTour** Pump Operations Manual

Note - Only qualified personnel should complete the diagnostic test.

Product serial number:
Completed by:
Date:

Cleaning and disinfecting with SideKick

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred 6" x 10" wipes (2060-000-001) include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface with SideKick wipe:

To clean:

1. Wipe down the external product surface with a fresh, clean wipe to remove all visible soils.
2. Repeat as necessary until the external product surface is visibly clean.
3. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Note - Use as many wipes as necessary.

To disinfect:

1. Clean first.
2. Wipe down the external product surface with a fresh, clean wipe until wet.
3. Allow the external product surface to remain wet for two minutes at room temperature.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

Care and maintenance

WARNING

- Do not wash the internal components of this support surface. Discard the support surface if contamination is found inside.
 - Do not immerse the product in cleaning or disinfectant solutions.
 - Do not allow liquid to pool on the product.
 - Always inspect support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers before each use. If compromised immediately remove the support surface from service.
 - Always make sure that you wipe each product with clean water and dry each product after you use cleaning agents. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow these care and maintenance instructions may void your warranty.
-

CAUTION

- Do not allow liquid to seep into the zipper area or watershed cover barrier when you clean the underside of the support surface. Fluids allowed to come in contact with the zipper may leak into the support surface core.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface as this may damage the product.
 - Always dry the support surface covers before you store, add linens, or place a patient on the surface. A dry product helps to prevent impaired performance.
 - Do not over expose the covers to higher concentration chemical solutions as these may degrade the covers.
 - Failure to follow manufacturing instructions may also affect useful life of the support surface cover.
-

The support surface cover is resistant to the following chemical solutions:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Accelerated Hydrogen Peroxide (AHP) 10,000 ppm hydrogen peroxide plus additional additives $\leq 0.13\%$ peroxyacetic acid
- Potassium peroxymonosulfate
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients, to avoid the risk of cross-contamination and infection.

Care of bottom cover

Always follow your hospital protocol for cleaning and disinfecting. This is a guide to help show you how to access the bottom surface cover for cleaning.

Note - This product is compatible with different bed frames. The bed features listed here may not be available for all compatible bed frames.

1. Remove the headboard and footboard from the bed and set aside.
2. Lower the siderails (Figure 4).

Figure 4 – Siderails lowered

3. Raise the head of bed angle and lower the foot end (Figure 5).

Figure 5 – Head of bed angle

4. With a helper, use the handles on the support surface to help fold the support surface to the foot end of the litter (Figure 6).

Figure 6 – Fold head end

5. Care for the underside of the support surface and the litter deck.
6. Make sure the surface and litter deck are dry.
7. Unfold the support surface and return the bed to a flat position (Figure 4).
8. Reverse steps 4 to 7 for the care of the foot end of the support surface and litter deck (Figure 7).

Figure 7 – Folded foot end

9. Replace the headboard and footboard before you place the product into service.

Top cover replacement

Tools required:

- None

Procedure:

1. Raise the bed height to the highest height position.
2. Lower the Fowler and Gatch sections to the lowest height positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.

Note - Use caution not to damage the fire barrier.

4. Remove and discard the top cover per hospital protocol.
5. Reverse to install.
6. Verify proper operation before you return the product to service.

Bottom cover replacement

A minimum of two operators is recommended for this task.

Tools required:

- None

Procedure:

1. Raise the bed height to the highest height position.
2. Lower the Fowler and Gatch sections to the lowest height positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.
4. Remove and save the top cover.
5. Disconnect the LAL hose.
6. Remove the foam crib assembly from the bed. Save the foam crib assembly.
7. Remove and discard the bottom cover per hospital protocol.
8. Place the supplied bottom cover on the bed with the bottom cover on the litter.
9. Place the foam crib assembly on top of the bottom cover. Align the foam crib to the cover.
10. Reconnect the LAL hose (A) and **TruTurn** hoses (B, C) (Figure 8). Connect the right hose with the red strip to the hose connector with the raised rib (D).


Figure 8 – Reconnect hoses

A	Low air loss (LAL)
B	TruTurn patient left
C	TruTurn patient right
D	Raised rib

11. Place the top cover over the top of the foam crib assembly. Align the top cover with the foam crib assembly.
12. Zip the cover to close. The zipper starts in the middle of the patient right side of the support surface.
13. Run a diagnostic test if paired with a pump. See the **IsoTour** Pump Operations Manual.

Note - Only qualified personnel should complete the diagnostic test.

14. Verify proper operation before you return the product to service.

Fire barrier replacement

Tools required:

- None

Procedure:

1. Raise the bed height to the highest height position.
2. Lower the Fowler and Gatch sections to the lowest height positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.
4. Remove and save the top cover.
5. Starting at the foot end, roll up the fire barrier on the foam crib assembly.

Note - Work from side to side, a little at a time, to move the fire barrier to the top of the support surface.

6. Discard the fire barrier.
7. Starting at the head end, roll the new fire barrier down and slide the fire barrier over the foam crib assembly.

Note - Align the fire barrier on the foam crib assembly before you slide the fire barrier over the foam crib assembly.

8. Slide the fire barrier down the foam crib assembly, work from side to side. Make sure that the fire barrier is tight on the foam crib assembly.
9. Align the foam crib assembly on top of the bottom part of the cover.

Note - Spread the excess fire barrier material equally below the foam crib assembly at the foot end.

10. Fold and align the top cover over the top of the foam crib assembly.
11. Zip the cover to close.
12. Verify proper operation before you return the product to service.

Spacer fabric replacement, premium option only

Tools required:

- None

Procedure:


1. Raise the bed height to the highest height position.
2. Lower the Fowler and Gatch sections to the lowest height positions.
3. Unzip the cover. The zipper starts in the middle of the patient right side of the support surface.
4. Remove and save the top cover.
5. Starting at the foot end, roll up the fire barrier on the foam crib assembly.
6. Remove and save the fire barrier.
7. Remove and discard the spacer fabric per hospital protocol.
8. Place the supplied spacer fabric on the foam core. Pull the fitted ends over the corners of the foam core corners.
9. Replace the fire barrier. Slide the fire barrier down the foam crib assembly, work from side to side. Make sure the fire barrier is tight on the foam crib assembly.
10. Replace the top cover. Align the top cover with the foam crib assembly.
11. Zip the cover to close. The zipper starts in the middle of the patient right side of the support surface.
12. Verify proper operation before you return the product to service.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **InTouch, IsoTour, ProCuity, S3, SideKick, Spirit Select, Stryker, SV2, TruTurn**. All other trademarks are trademarks of their respective owners or holders.