## IsoTour® Pump

Operations Manual

Standard and premium


**REF** 2874

# Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Cardiopulmonary resuscitation (CPR), disconnection point
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
<b>REF</b>	Catalogue number
<b>SN</b>	Serial number
<b>LOT</b>	Batch code
<b>CH REP</b>	Authorized representative in Switzerland
<b>EC REP</b>	Authorized representative in the European Community
<b>MD</b>	European medical device
<b>US Patents</b>	For US Patents see <a href="http://www.stryker.com/patents">www.stryker.com/patents</a>
	Manufacturer
	Date of manufacture
<b>CE</b> 2797	CE mark
<b>UK CA</b>	UK Conformity Assessment mark
	Importer

	Unique device identifier
	Mass of product
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provisions for protective earthing or reliance upon installation conditions.
	Type B applied part - 2872 <b>IsoTour</b> support surface
	Fuse
	Medical Equipment Classified by Underwriters Laboratories Inc. With Respect to Electric Shock, Fire, and Mechanical Hazards Only in Accordance with ANSI/AAMI ES60601-1:2012, IEC 60601-1-8:2012, and CAN/CSA-C22.2 No. 60601-1:16.
IP21	Solids: Protection from touch by fingers and objects greater than 12 mm Liquids: Protection from the fall of vertical water drops
	Keep dry
	Stacking limit by number
	This side up
	Fragile
	Do not use sharp objects to open the package

# Table of Contents

Warning/Caution/Note Definition .....	2
Summary of safety precautions .....	3
Introduction .....	4
Product description .....	4
Intended use .....	4
Contraindications .....	4
Expected service life .....	5
Disposal/recycle .....	5
Specifications .....	5
Environmental conditions .....	6
Product illustration .....	6
Button and icon descriptions .....	7
Features .....	8
Contact information .....	8
Serial number .....	8
Setup .....	9
Connecting the power cord .....	9
Connecting the hose .....	10
Operation .....	11
Powering on the pump .....	11
Activating CPR release .....	12
Powering on or off Low Air Loss (LAL) .....	12
Selecting <b>TruTurn</b> .....	12
Increasing or decreasing <b>TruTurn</b> .....	13
Pausing or canceling <b>TruTurn</b> .....	14
Lock or unlock product functions .....	14
Powering off the pump .....	15
Accessories and parts .....	16
Troubleshooting .....	17
Product alarms .....	19
Preventive maintenance .....	20
Cleaning and disinfecting with SideKick® .....	21
Cleaning .....	22
Disinfecting .....	23
EMC information .....	24

# Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

## **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

## **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

**Note** - Provides special information to make maintenance easier or important instructions clearer.

# Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

## WARNING

- Consult the operations manual for the support surface before you attach the pump.
  - Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
  - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
  - Always setup and use in accordance with the instructions to avoid harmful interference to other devices in the vicinity. This product radiates radio frequency energy.
  - Always make sure the hose is routed down and away from the patient. Avoid accidental patient contact at the foot end of the support surface. Hose temperature may reach 48° C (118.4° F).
  - Never transfer or transport a patient while **TruTurn** is in use.
  - Always raise the bed siderails before you begin **TruTurn**, to avoid the risk of patient fall.
  - Do not leave the patient unattended until **TruTurn** of the patient position is complete.
  - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.
  - Always unplug the product before you begin cleaning or disinfecting.
  - Do not immerse the product in cleaning or disinfectant solutions.
  - Do not spray or allow liquids to enter the power cord ports or the hose ports.
  - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
  - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
  - Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **IsoTour**, including cables specified by the manufacturer.
  - Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.
- 

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
  - Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
  - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
  - Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
-

# Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

---

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- 

## Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

## Product description

The Stryker Model 2874 **IsoTour**® Pump is an electric pump for use with the Stryker Model 2872 **IsoTour** Support Surface. This pump provides Low Air Loss for the **IsoTour** standard or the premium support surface.

The pump when connected to the premium support surface delivers the air to inflate and deflate the air bladders for the **TruTurn**® function.

## Intended use

The **IsoTour** Pump is for use with the **IsoTour** support surface and assists in the prevention and treatment of pressure injuries or pressure ulcers (all stages, unstageable injury, and deep tissue injury). We recommend you implement this product in combination with clinical evaluation of risk factors and skin assessments made by a healthcare professional.

This pump is for use with the **IsoTour** support surface in an acute care setting. This may include critical care, step down, progressive care, med/surg, sub-acute care, and post anesthesia care unit (PACU), or other locations as prescribed by a physician. Operators of this pump include healthcare professionals (such as nurses, nurse aids, or doctors).

This product is not intended to be sterile, does not include a measuring function, or be used in a home healthcare environment.

## Contraindications

The **TruTurn** feature is contraindicated for use with patients that have unstable spine.

**Note** - **TruTurn** requires the **IsoTour** pump to be connected for use. The **TruTurn** feature is visually identified by the icons on the **IsoTour** pump (Figure 1).


Figure 1 – TruTurn icon

## Expected service life

The **IsoTour Pump** has a 5 year expected service life under normal use, conditions, and with appropriate periodic maintenance.

The air hose has a 2 year expected service life under normal use, conditions, and with appropriate periodic maintenance.

## Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

## Specifications


Height		17.25 in.	43.8 cm
Width		9.25 in.	23.5 cm
Depth		7.25 in.	18.4 cm
Weight		16.3 lb	7.4 kg
Interface control panel		LCD	
Input voltage		100-240V	
Current		0.7 Amps	
Input frequency		50/60 Hz	
Noise level		<50 dba	
Mode of operation		Continuous	
Power cord	Short	3.2 ft	1 m
	Long	16.4 ft	5 m
Compliance		IEC 60601-1 3.1 Edition; RoHS Directive 2002/95/EC Reach; IEC 60601-1-8:2012, IEC 60601-1-2:2014	
Compatible with bed frame	IsoTour standard and premium support surface with an IsoTour pump	ProCuity™, 3002 S3®, 3005 S3®, InTouch®, Spirit Select®	

## Environmental conditions

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

Stryker reserves the right to change specifications without notice.

## Product illustration

A	Graphical user interface	H	Power inlet
B	Status indicators	I	Anti-vibration foot
C	Bumpers	J	Fuse
D	Transport handle	K	Flex grip
E	Hooks	L	Hose port
F	HEPA filter	M	Hose
G	Back pump housing	N	Front pump housing

## Button and icon descriptions

The buttons are located on the outside of the graphical user interface.

Icon	Name	Description
	Power button	See <i>Powering on the pump</i> (page 11) or <i>Powering off the pump</i> (page 15)
	Cancel icon	Cancel current function, returns to the main menu
	Pause icon	Pause current function
	Action button, right	User selects a function
	Action button, left	User selects a function
	Lock or unlock icons	See <i>Lock or unlock product functions</i> (page 14)
	Back icon	Returns to the previous screen
	Low Air Loss (LAL) icon	Low Air Loss function active, on by default, see <i>Powering on or off Low Air Loss (LAL)</i> (page 12)
	Fill default segment 2 icon	Select patient turn to the right or to the left
	Maximum fill	Maximum fill is 3 bars, see <i>Increasing or decreasing TruTurn</i> (page 13)
	Decrease or increase	Decrease or increase patient turn
	Check siderails animated icon	Reminder to check siderails before <b>TruTurn</b> is activated
	Check air connector alert icon	Product detects air connector error
	Kinked air hose alert icon	Product detects a bend or twist in the air hose
	Alarm off	Closes and resets pump for kinked hose or check connector

## Features

TruTurn	Available for the premium support surface. Air bladders inflate to rotate up to 30 degrees. Helps to off-load the sacrum. Turn angle can be customized per patient tolerance and comfort.
Moisture management	The top surface cover is designed to draw excess moisture away from the patient's skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to help manage microclimate of the skin.
Standard support surface	Orange bottom cover, gel and foam core, LAL capable when connected to pump, handles, and D-rings
Premium support surface	Blue bottom cover, gel and foam core, spacer fabric for enhanced air flow, handles, D-rings, anchor straps, <b>TruTurn</b> and LAL capable when connected to pump

## Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

**Note** - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

## Serial number

# Setup

Unpack the cartons and check all items. Make sure that the product operates before you place into service.

---

## WARNING


- Consult the operations manual for the support surface before you attach the pump.
  - Do not use in the presence of flammable anesthetics, nitrous oxide, or oxygen-rich environments. Risk of explosion can result.
  - Do not expose this product to liquids while powered as this could result in a severe electrical hazard.
  - Always setup and use in accordance with the instructions to avoid harmful interference to other devices in the vicinity. This product radiates radio frequency energy.
- 

## CAUTION

- Always allow the product to come to room temperature after exposure to extreme high or low temperatures.
  - Always check the power cord for cuts, exposed wires, worn insulation, or other damage that could result in electrical hazards.
- 

## Connecting the power cord

1. Bring the product to room temperature.
2. Attach the power cord to the pump.
3. Insert the power cord under the power cord retention flex grip (A) (Figure 2).

**Figure 2 – Cord retention flex grip**

4. Plug the pump power cord into an accessible outlet.
5. Place the pump power cord out of the path of foot traffic to avoid trip hazards.
6. Place the pump in a secure location. Allow space around the pump to connect and disconnect the hoses and power cord.
  - a. Option 1: Hang the pump on a compatible bed. Make sure that the pump is secure.
  - b. Option 2: Place the pump on the floor.


## Connecting the hose

---

**WARNING** - Always make sure the hose is routed down and away from the patient. Avoid accidental patient contact at the foot end of the support surface. Hose temperature may reach 48° C (118.4° F).

---

1. Attach the hose to the pump hose ports.
2. Open the hose port (A) on the support surface (Figure 3).

**Figure 3 – Hose connection**

3. Attach the hose to the support surface hose port. Make sure that rib (B) lines up with the red line on the hose (C).

# Operation

## Powering on the pump

1. See *Setup* (page 9) to prepare the pump for operation.
2. Press the power button once. The startup screen appears (Figure 4).


Figure 4 – Startup screen - operator view

**Note** - The pump detects either the standard or premium support surface.

3. The pump, on startup, pumps air into the left bladder to determine which compatible mattress is connected.
4. If the pump detects there are no turn bladders, the pump shows the standard software support surface display (Figure 5).

Figure 5 – Standard software support surface display

**Note** - TruTurn is available for the premium support surface only.

5. If the pump detects there are turn bladders, the pump shows the premium software support surface display (Figure 6).

Figure 6 – Premium software support surface display

**Note** - TruTurn is available for the premium support surface only.

## Activating CPR release

Remove the hose from the pump to deflate the support surface for CPR. Follow your hospital protocol for CPR.

## Powering on or off Low Air Loss (LAL)

LAL is on by default. This function can be turned off or on.

1. Press the action button next to the active (green) LAL icon to turn off this function (Figure 7).

Figure 7 – LAL on

2. Press the action button next to the inactive (white) LAL icon to turn on this function (Figure 8).

Figure 8 – LAL off

## Selecting TruTurn

---

### WARNING

- Never transfer or transport a patient while **TruTurn** is in use.
  - Always raise the bed siderails before you begin **TruTurn**, to avoid the risk of patient fall.
  - Do not leave the patient unattended until **TruTurn** of the patient position is complete.
- 

1. Read the **IsoTour** support surface Operations Manual before you start **TruTurn**.
2. Make sure that the head of bed is at zero degrees before you start **TruTurn**.
3. Press the action button next to the patient left or patient right **TruTurn** (Figure 9).

Figure 9 – TruTurn

4. Raise the siderails (Figure 10).


Figure 10 – Check siderail

5. Press the action button next to the time to set the turn duration.

- 02:00:00 = 2 hours
- 00:30:00 = 30 minutes
- 00:00:15 = 15 seconds

**Note**

- Prior to transfer, transport, or other patient handling activity, cancel **TruTurn** function. You can either cancel the **TruTurn** function, turn off the pump, or pull the CPR.
- The support surface will turn up to 30 degrees.
- The count down for the time selected will start after the pump has completed the fill.

## Increasing or decreasing TruTurn

After **TruTurn** has reached the default angle fill, you can increase or decrease the turn.

To increase the turn, press and hold the action button next to the up arrow.

Figure 11 – Increase or decrease turn

To decrease the turn, press and hold the action button next to the down arrow.

**Note** - The increase and decrease arrows turn green when active.

## Pausing or canceling TruTurn

You can only pause **TruTurn** when the pump fills the support surface with air. You can cancel **TruTurn** at any time.

To pause **TruTurn**, press the action button next to the pause icon .


Figure 12 – TruTurn fill

To cancel **TruTurn**, press the action button next to the cancel icon . The turn bladder deflates and returns the surface to a flat position.


Figure 13 – Active TruTurn

## Lock or unlock product functions

The lock icon is orange and closed when active.


To turn the lock on or off:

1. Press and hold the unlock icon for 1 second to turn on the lock.

**Figure 14 – Unlocked**

2. Press and hold the lock icon for 2 seconds to turn off the lock.

**Figure 15 – Locked**

**Note** - When the lock is active, the other functions are not available.

## Powering off the pump

Press the power button and hold for 2 seconds to power off the pump.

### **Note**

- Allow the support surface to deflate before you unplug the power cord.
- For optional deflation, remove the CPR plug from the support surface.

## Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region.

Name	Part number
IsoTour air hose	2874-001-200
Power cord, type B, 1 meter	2874-007-001
Power cord, type B, 5 meter	2874-007-002
Power cord, type E/F, 1 meter	2874-007-003
Power cord, type E/F, 5 meter	2874-007-004
Power cord, type B, 220V, 1 meter	2874-007-005
Power cord, type B, 220V, 5 meter	2874-007-006
Power cord, type J, 1 meter	2874-007-007
Power cord, type J, 5 meter	2874-007-008
Power cord, type I, 1 meter	2874-007-009
Power cord, type I, 5 meter	2874-007-010
Power cord, type G, 1 meter	2874-007-011
Power cord, type G, 5 meter	2874-007-012
Power cord, type L, 1 meter	2874-007-013
Power cord, type L, 5 meter	2874-007-014
Power cord, type N, 1 meter	2874-007-015
Power cord, type N, 5 meter	2874-007-016
Power cord, type D, 1 meter	2874-007-037
Power cord, type D, 5 meter	2874-007-038
Power cord, type M, 1 meter	2874-007-041
Power cord, type M, 5 meter	2874-007-042

# Troubleshooting

Problem	Screen	Cause	Recommended action
Kinked air hoses		Air hose is bent or an obstruction in the hose may cause air flow to be compromised  <b>TruTurn</b> was attempted but you have a standard model (bottom cover is orange)	<ol style="list-style-type: none"> <li>1. Make sure that the hoses are straight and air flow is not obstructed.</li> <li>2. Press the action button next to the Alarm off icon to reset the pump.</li> </ol>  <p>Contact sales support for <b>TruTurn</b> option.</p>
Missed connection with air hoses		Air hoses are not connected to the pump or the support surface	<ol style="list-style-type: none"> <li>1. Make sure that the hoses are seated all the way onto the hose ports on the pump or the support surface.</li> <li>2. Press the action button next to the Alarm off icon to reset the pump.</li> </ol> 
Power loss, product does not turn on		Power cord not seated, power cord unplugged from outlet, or possible internal damage, button is stuck	<ol style="list-style-type: none"> <li>1. Make sure that the power cord is plugged into the product and the outlet.</li> <li>2. Secure the power cord to the flex grip. See <i>Connecting the power cord</i>.</li> <li>3. Replace the LCD screen.</li> </ol>
Power loss, while <b>TruTurn</b> is in use		Unintentional or intentional power loss	<ol style="list-style-type: none"> <li>1. Pull the hose from the pump or the support surface to activate CPR. The bladder deflates and returns the support surface to a flat position.</li> </ol>
Button not responsive		Lock function active, button is stuck	<ol style="list-style-type: none"> <li>1. Press and hold the action button next to the locked icon for 2 seconds to turn off the lock.</li> <li>2. Replace the LCD screen.</li> <li>3. Power cycle the pump.</li> </ol>

Problem	Screen	Cause	Recommended action
<p><b>TruTurn</b> not shown on screen</p>		<p>Kinked hose during setup or <b>IsoTour</b> standard model is connected</p>	<ol style="list-style-type: none"> <li>1. Check for kinked hoses.</li> <li>2. Verify that you have premium bottom cover (black).</li> </ol> <p>Contact sales support for <b>TruTurn</b> option.</p>
<p><b>TruTurn</b> shown on screen</p>		<p>Graphical user interface shows <b>TruTurn</b> but you have a standard model (bottom cover is orange)</p>	<ol style="list-style-type: none"> <li>1. Turn off the pump to reset the pump features.</li> <li>2. Make sure that the air hose is connected before you turn the pump on.</li> </ol>

# Product alarms

The product alarms for this product are display only, no audible alarms.

Icon	Priority alarm	Name	Delay	Therapy interrupted
	Low	Check air connector	Less than 45 seconds	Yes
	Low	Kinked air hose	Less than 5 seconds	Yes

## Note

- Check the air connector alarm is active when **TruTurn** is active.
- Alarms display on the interface control panel of the pump.
- Alarms reset if there is power loss to the pump.
- The operator must press the action button next to the Alarm off icon  to reset the alarm.

# Preventive maintenance

**WARNING** - Do not modify or change this device. Service should only be completed by qualified personnel. Failure could result in injury and void your warranty.

---

Remove product from service before you perform preventive maintenance. At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel

**Note** - Consult your local regulations to dispose of electronic equipment.

Inspect the following items:

- All fasteners are secure
- Pump housing or components (hoses, power cords, or case) for cracks, holes, or damage
- Bed frame pump hooks are not damaged
- No air leaks from the pump or the attached connectors or hoses
- Graphical user interface is not cracked or damaged
- HEPA filter (replace each year)
- Fuse
- All functions on graphical user interface operate
- Run diagnostic test (All okay)

Product serial number:
Completed by:
Date:

# Cleaning and disinfecting with SideKick®

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred spray (2070-000-001) and 6" x 10" wipes (2060-000-001) include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

**Note** - For safety information, read the product label.

## To clean or disinfect the external product surface with SideKick wipe:

### To clean:

1. Wipe down the external product surface with a fresh, clean wipe to remove all visible soils.
2. Repeat as necessary until the external product surface is visibly clean.
3. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

**Note** - Use as many wipes as necessary.

### To disinfect:

1. Clean first.
2. Wipe down the external product surface with a fresh, clean wipe until wet.
3. Allow the external product surface to remain wet for two minutes at room temperature.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

## To clean or disinfect the external product surface with SideKick spray:

### To clean:

1. Spray **SideKick** on a mop, sponge, paper cloth, or cloth towel.
2. Wipe down the external product surface with a fresh, clean mop, sponge, paper cloth, or cloth towel to remove all visible soils.
3. Repeat as necessary until the external product surface is visibly clean.
4. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

### To disinfect:

1. Clean first.
2. Spray **SideKick** on a mop, sponge, paper cloth, or cloth towel.
3. Wipe down the external product surface with a fresh, clean mop, sponge, paper cloth, or cloth towel until wet.
4. Allow the external product surface to remain wet for two minutes at room temperature.
5. Wipe dry with a cloth or allow the external product surface to air dry before you return the product to service.

**Note** - Follow your hospital protocols to launder cloth towels or dispose of wipes or paper cloths.

# Cleaning

---

## WARNING

- Always unplug the product before you begin cleaning or disinfecting.
  - Do not immerse the product in cleaning or disinfectant solutions.
  - Do not spray or allow liquids to enter the power cord ports or the hose ports.
  - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- 

## CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
- 

Recommended cleaning agents include:

- A clean damp cloth and mild soapy water

To clean the external side of the pump:

1. Unplug the controller.
2. Wipe the outside of the product with a clean, soft cloth moistened with a mild soap and water solution to remove foreign material.
3. Follow specified contact time in accordance with the manufacturer's instructions for use.
4. Using a clean, soft cloth, wipe the outside of the product with water.
5. Dry the external surface with a clean, dry cloth.

# Disinfecting

---

## WARNING

- Always unplug the product before you begin cleaning or disinfecting.
  - Do not immerse the product in cleaning or disinfectant solutions.
  - Do not spray or allow liquids to enter the power cord ports or the hose ports.
  - Always make sure that you wipe each product with clean water and dry each product after you clean. Some cleaning agents are corrosive in nature and may cause damage to the product. Failure to follow the cleaning and disinfecting instructions may void your warranty.
- 

## CAUTION

- Do not over expose the product to higher concentration chemical solutions as these may degrade the product.
  - Failure to follow manufacturing instructions may also affect useful life of the product.
- 

Recommended chemical solutions:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution (6500 ppm)
- Accelerated Hydrogen Peroxide (AHP) 10,000 ppm hydrogen peroxide plus additional additives  $\leq$  0.13% peroxyacetic acid
- Potassium peroxymonosulfate
- 70% isopropyl alcohol

Follow hospital protocol for support surface care between patients to avoid the risk of cross-contamination and infection.

**Note** - Make sure you unplug the product before you begin cleaning and disinfecting.

# EMC information

**WARNING** - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

**Note**

- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- This equipment is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Guidance and manufacturer's declaration - electromagnetic emissions		
The 2874 <b>IsoTour</b> pump is intended for use in the electromagnetic environment specified below. The customer or the user of the 2874 <b>IsoTour</b> pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	The 2874 <b>IsoTour</b> pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The 2874 <b>IsoTour</b> pump is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

**WARNING**

- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **IsoTour**, including cables specified by the manufacturer.
- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the products. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they are operating properly.

Recommended separations distances between portable and mobile RF communication equipment and the 2874 IsoTour pump			
The 2874 <b>IsoTour</b> pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 2874 <b>IsoTour</b> pump can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 2874 <b>IsoTour</b> pump as recommended below, according to the maximum output power of the communications equipment.			
Band (MHz)	Service	Maximum Power (W)	Minimum Separation Distance (m)

Recommended separations distances between portable and mobile RF communication equipment and the 2874 IsoTour pump			
380-390	TETRA 400	1.8	0.3
430-470	GMRS 460; FRS 460	2.0	0.3
704-787	LTE Band 13, 17	0.2	0.3
800-960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	2.0	0.3
1,700-1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	2.0	0.3
2,400-2,570	Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE Band 7	2.0	0.3
5,100-5,800	WLAN 802.11 a/n	0.2	0.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note** - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The 2874 IsoTour pump is suitable for use in the electromagnetic environment specified below. The customer or the user of the 2874 IsoTour pump should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	$\pm 8$ kV contact $\pm 15$ kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV for input/output lines	$\pm 1$ kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.

**Guidance and manufacturer's declaration - electromagnetic immunity**

<p>Voltage dips, voltage variations and short interruptions on power supply input lines</p> <p>IEC 61000-4-11</p>	<p>0%U<sub>T</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0%U<sub>T</sub> for 1 cycle</p> <p>70%U<sub>T</sub> (30% dip in U<sub>T</sub>) for 25/30 cycles</p> <p>0% U<sub>T</sub> for 250/300 cycles</p>	<p>0%U<sub>T</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0%U<sub>T</sub> for 1 cycle</p> <p>70%U<sub>T</sub> (30% dip in U<sub>T</sub>) for 25/30 cycles</p> <p>0% U<sub>T</sub> for 250/300 cycles</p>	<p>Main power quality should be that of a typical commercial or hospital environment. If the user of the 2874 <b>IsoTour</b> pump requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field</p> <p>IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

**Note** - U<sub>T</sub> is the a.c. mains voltage before applications of the test level.

<p>Conducted RF IEC 61000- 4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.7 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the 2874 <b>IsoTour</b> pump, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>D=(1.2) (\sqrt{P})</math> 80 MHz to 800 MHz</p> <p><math>D=(2.3) (\sqrt{P})</math> 800 MHz to 2.7 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
---	--	----------------------------	---

**Note**

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- The ISM (Industrial, Scientific, and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 2874 **IsoTour** pump is used exceeds the applicable RF compliance level above, the 2874 **IsoTour** pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the 2874 **IsoTour** pump.


<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 Vrms.

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