Isolibrium™ Support Surface

REF 2971

Integrated with InTouch® FL27 (2131
/ 2141 / 2151 / 2152) version 4.0

stryker®

Operations/Maintenance Manual

Symbols

EN

	See instruction manual/booklet
	Operating instructions
	General warning
	Caution
	Warning; electricity
	Alternating current
	Model
	Serial number
	Manufacturer
	Date of manufacture
	Safe working load
	Mass of equipment
	Protective earth terminal
IPX4	Protection from liquid splash
	Defibrillation proof type B applied part
	Medical Equipment recognized by UL LLC with respect to electric shock, fire, and mechanical hazards in accordance with ANSI/AAMI ES60601-1:2005/ (R)2012, A1:2012, C1:2009/(R)2012, A2:2010/(R)2012, and CAN/CSA-C22.2 No. 60601-1 (2008)
	In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.

Symbols

EN

	Damp wipe only
	Do not tumble dry
	Do not dry-clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
	Pressure redistribution
	Microclimate, Low Air Loss
	Lateral rotation
	Shear
	Turn assist

Isolibrium™ Support Surface

REF 2971

*Integrated with InTouch® FL27 (2131
/ 2141 / 2151 / 2152) version 4.0*

stryker®

Operations/Maintenance Manual

Symbols

EN

	See instruction manual/booklet
	Operating instructions
	General warning
	Caution
	Warning; electricity
	Alternating current
	Model
	Serial number
	Manufacturer
	Date of manufacture
	Safe working load
	Mass of equipment
	Protective earth terminal
IPX4	Protection from liquid splash
	Defibrillation proof type B applied part
	Medical Equipment recognized by UL LLC with respect to electric shock, fire, and mechanical hazards in accordance with ANSI/AAMI ES60601-1:2005/ (R)2012, A1:2012, C1:2009/(R)2012, A2:2010/(R)2012, and CAN/CSA-C22.2 No. 60601-1 (2008)
	In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Refer to your local distributor for return and/or collection systems available in your country.

Symbols

EN

	Damp wipe only
	Do not tumble dry
	Do not dry-clean
	Do not iron
	Allow to completely air dry
	Chlorinated bleach
	Pressure redistribution
	Microclimate, Low Air Loss
	Lateral rotation
	Shear
	Turn assist

Table of Contents

EN

Warning/Caution/Note Definition	1-3
Introduction	1-4
Product description	1-4
Intended use	1-4
Expected service life	1-4
Contraindications.....	1-5
Specifications	1-6
Contact information.....	1-7
Serial number location.....	1-7
Product illustration	1-8
Isolibrium features	1-9
Summary of safety precautions	1-10
Installation	1-12
Operation	1-14
Applying the linens.....	1-14
Positioning a patient on the support surface	1-15
Managing incontinence and drainage.....	1-16
Activating and resetting CPR	1-16
Relocating the support surface	1-17
Transferring of patient from one patient support platform to another.....	1-17
Storing the product	1-18
Cleaning	1-19
Disinfecting	1-20
Preventive maintenance.....	1-21
Quick reference replacement parts	1-23
Block diagram.....	1-24
Isolibrium wiring.....	1-24
Pneumatic	1-25
Circuit boards	1-26
Power supply assembly	1-26
Main power board assembly	1-28
Cables	1-31
Fans and cable assembly.....	1-31
SPI cables (main board to foot box) assembly	1-31
Integration cable assembly 2971-022-133	1-32
Service	1-33
Protecting against Electrostatic Discharge (ESD)	1-33
Cover replacement, top.....	1-34
Cover replacement, bottom.....	1-35
Pod assembly replacement.....	1-37
Foot box cover access.....	1-40
Solenoid valve replacement	1-43
Power supply replacement	1-44

Table of Contents

EN

Main board replacement	1-46
Power cord replacement	1-48
Bed integration cable connector replacement	1-50
Pump assembly replacement	1-51
Low Air Loss (LAL) fan replacement	1-53
Foot box cooling fan replacement	1-56
Sensor board replacement	1-57
Resonator replacement.....	1-59
Turn bladder assembly replacement	1-61
Product labels	1-63
Troubleshooting	1-66
Isolibrium assembly.....	1-71
Footbox assembly.....	1-76
Pump assembly kit 2971-700-001	1-85
Valve manifold assembly kit 2971-700-007	1-86
Recycling passports	1-87
Recycling passports.....	1-88
Recycling passports.....	1-89
EMC Information	1-90
Warranty	1-94
Warranty exclusion and damage limitations.....	1-94
To obtain parts and service.....	1-94
Return authorization.....	1-94
Damaged product.....	1-95
International warranty clause	1-95
Patent information	1-95

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note: Provides special information to make maintenance easier or important instructions clearer.

Introduction

EN

This manual assists you with the operation or maintenance of the Stryker Model 2971 **Isolibrium™** support surface. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

WARNING

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not change the product or any components of the product. Modifying the product can cause unpredictable operation that results in injury to the patient or operator. Modifying the product also voids its warranty.
-

Notes

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. Therefore, while this manual contains the most current product information available at the time of printing, there may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

Isolibrium, Stryker Model 2971, is a powered support surface with features that provide Pressure Redistribution, Low Air Loss, and Lateral Rotation. **Isolibrium** is for use with the **InTouch®** (revisions 4.0 or higher software) bed frame. **Isolibrium** and **InTouch** are an integrated bed and surface system that caregivers use to assist in patient care.

Intended use

The **Isolibrium** support surface assists in the prevention of pressure ulcers. All categories or stages of pressure ulcers including I, II, III, IV, unstageable, and deep tissue injury. Before implementing **Isolibrium**, a healthcare professional should complete a clinical evaluation of risk factors and skin assessments. The **Isolibrium** patient-specific immersion to assist in managing microclimate of the patient skin, patient repositioning, early mobility, and the prevention of pulmonary complications.

The **InTouch** and **Isolibrium** system is for use by patients in an acute care setting. The acute care settings may include critical care, step down, progressive care, medical / surgical, sub-acute care, and post anesthesia care unit (PACU), or other locations as prescribed. The patient, accessories, and any extra loading must not exceed the safe working load of **InTouch** or **Isolibrium**.

The operators for the system are healthcare professionals (nurses, nurse aids, and doctors) who can use all system operations (for example: bed motion functions, nurse call, siderail communications, bed exit, or therapy options). Patients and bystanders can use the bed motion functions, nurse call, and siderail communications but are restricted from the use of the powered support surface functions. Trained professionals install, service, and calibrate the **InTouch** and **Isolibrium** system.

Expected service life

Isolibrium support surface has a five year expected service life under normal use, conditions, and with appropriate periodic maintenance.

Isolibrium covers (top and bottom), air pods, and turn bladders have a two year expected service life under normal use, and conditions.

Contraindications

Stryker promotes the clinical assessment of each patient and appropriate usage by the operator.

The **InTouch** and **Isolibrium** system is not intended:

- to support a patient in a prone position
- for use with patients 35 inches or less in height
- for use with patients that weigh less than 50 pounds
- for use with a full oxygen tent
- to support more than one individual at a time
- for use with a patient less than two years old
- in a home healthcare environment
- as a sterile product
- for use in the presence of flammable anesthetic mixture with air rich in oxygen or nitrous oxide

Note: When **Isolibrium** is on the **InTouch** bed, a half-oxygen tent within the sidrails is allowable. The half-oxygen tent must be above the area of the support surface that contains the foot box.

Isolibrium is not intended for use on patients with the following:

- unstable fractures
- unstable spinal cord injuries

Isolibrium Lateral Rotation and Turn Assist features are for use only after you complete a clinical assessment. This is to make sure that no adverse reaction would result in patients with the following:

- spinal cord injury
- skeletal traction
- significant hemoptysis
- bleeding disorders
- rib fractures
- for whom the techniques cause increased dyspnea or wheezing
- hemodynamic instable
- increased intracranial pressure
- skeletal traction







Introduction

EN

Specifications

Model	2971 Isolibrium	
Length	84 in.	213,4 cm
Width	35 in.	88,9 cm
Thickness	8 - 10.5 in.	20.3 - 26.7 cm
Weight	90 lb	40.8 kg
Therapeutic weight	50 lb to 350 lb	22.7 kg to 158.7 kg
 Safe working load Isolibrium (Sum of the patient and accessory weight)	460 lb	208.6 kg
Electrical requirements	ANSI/AAMI ES60601-1:2005/ (R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012 CAN/CSA-C22.2 No. 60601-1 (2008)	120VAC 50/60Hz, 2.5A
Power cord	3 foot, 16 AWG cord with hospital grade plug for use with the InTouch foot end auxiliary power outlet	
Low Air Loss flow rate	100 LPM	
Leakage current	300 uA Maximum	
Noise level	50 dBA maximum	

Stryker reserves the right to change specifications without notice.

Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

See the **InTouch** manual for bed specifications.

Contact information


Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location


The serial number (A) for your support surface is located on the patient right side of the foot box on the specification label. Unzip the support surface top cover (see [Cover replacement, top on page 1-34](#)).

Introduction

EN

Product illustration

A	Head end
B	Foot end
C	D-rings
D	CPR straps
E	Power cord

F	Integration cable
G	Handles
H	Top support surface cover
I	Bottom support surface cover

Note: To view the power cord label, the operator must position the cord by rotating it until the label is legible from a distance of up to 1m.

Isolibrium features

Physical Concepts Pressure Redistribution	The ability of a support surface to distribute load over the contact areas of the human body.
Pressure Redistribution setpoint	Allows the operator to input patient weight range resulting in automatic change of zoned sensor technology allowing immersion and envelopment.
Lateral Rotation Therapy	A feature of a support surface that provides rotation about a longitudinal axis as characterized by degree of patient turn, duration, and frequency.
Turn Assist	A feature of a support surface that provides a single position turn about a longitudinal axis.
Max Inflate	Allows maximum firmness inflation for use when boosting a patient and on lateral transfer, Promotes patient ease with dangling, egress, and ingress.
Microclimate	The top surface cover draws excess moisture away from the patients skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to assist in the management of (microclimate) the skin.

Summary of safety precautions

EN

Carefully read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not change the product or any components of the product. Modifying the product can cause unpredictable operation that results in injury to the patient or operator. Modifying the product also voids its warranty.
 - Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
 - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
 - Do not move the support surface without assistance.
 - Do not operate the support surface without the **Isolibrium** top and bottom covers and linens.
 - Always position the support surface on the **InTouch** bed frame to that the support surface does not contact the headboard.
 - Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the **Isolibrium** support surface with **InTouch** bed frames.
 - Do not attach the power cord to any moving parts of the **InTouch** bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
 - Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
 - Do not use extension cords with **Isolibrium**. It is only intended to be powered by **InTouch** with the power cord supplied.
 - Always use the power source provided by the **InTouch** foot end auxiliary power outlet with **Isolibrium**.
 - Always make sure that the operator has access to the CPR straps.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not use fitted sheets with this support surface.
 - Always center the patient on the support surface. Align the patient's head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
 - Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
 - Always raise all of the **InTouch** bed siderails before starting Turn Assist or Lateral Rotation functions.
 - Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
 - Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.
 - Always monitor the patient condition at regular intervals for patient safety.
 - Do not use the support surface as a transfer device to avoid the risk of patient injury.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
-

Summary of safety precautions

EN

WARNING (CONTINUED)

- Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
 - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.
 - Do not service or perform maintenance while the product is in use.
-

CAUTION

- Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
 - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
 - Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
 - Always use ESD protective equipment before opening antistatic bags and servicing electronic parts.
 - Do not place unprotected circuit boards on the floor.
-

Installation

EN

The product must reach the ambient operating temperature before installation or operation to avoid the risk of permanent damage to this product. Upon initial installation, reset the CPR.

WARNING

- Risk of Shock. This product must only be connected to a supply mains with protective earth. Grounding reliability can be achieved only when a hospital grade receptacle is used.
 - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
 - Do not move the support surface without assistance.
 - Do not operate the support surface without the **Isolibrium** top and bottom covers and linens.
 - Always position the support surface on the **InTouch** bed frame so that the support surface does not contact the headboard.
 - Do not use when gaps between the support surface and the bed frame are present to prevent the risk of entrapment. Only use the **Isolibrium** support surface with **InTouch** bed frames.
 - Do not use fitted sheets with this support surface.
 - Do not attach the power cord to any moving parts of the **InTouch** bed frame. Improper routing of the power cord could cause the risk of electric shock, pinching the power cord, or interfere with bed functions.
 - Always confirm proper functionality of the scale system after you install the support surface to prevent product malfunction.
 - Do not use extension cords with **Isolibrium**. It is only intended to be powered by **InTouch** with the power cord supplied.
 - Always use the power source provided by the **InTouch** foot end auxiliary power outlet with **Isolibrium**.
-

Note: Do not plug into the power source until the assembly is complete.

Prerequisite: Minimum of two operators are required to install the support surface.

1. Put the **InTouch** bed into the flat position before installing the support surface.
2. If powered, unplug the **InTouch** from the wall outlet and turn the battery switch to off.
3. Remove the **InTouch** headboard and footboard from the bed.
4. Position an operator at each end of the support surface and lift onto the **InTouch** bed frame.
5. Reinstall the **InTouch** headboard and footboard.
6. Position the support surface against the footboard. Allow a small air gap at the headboard.
7. Connect the integration cable from the support surface foot end to the **InTouch** cable adapter (B) ([Figure 1-1 on page 1-13](#)).
8. Turn the locking collars clockwise to secure the connection.
9. Place the linens on the support surface as loose as possible to avoid creating a **hammocking** effect and increasing interface pressures.
10. Plug the power cord for the support surface into the **InTouch** foot end auxiliary power outlet (A) ([Figure 1-1 on page 1-13](#)).

Installation

EN


Figure 1-1: InTouch foot end auxiliary power outlet and mattress connector

11. Make sure that the power cord is routed to the side of the bed and not out the foot end of the **InTouch**.
12. Plug in the **InTouch** power cord into a properly grounded hospital grade three-prong power outlet.
13. Calibrate the bed following the procedures in the **InTouch** maintenance manual.

Notes

- **Isolibrium** is only operational when the power cord for the **InTouch** bed is in a 120v wall outlet.
- When the support surface is plugged into the **InTouch** bed, the **InTouch** footboard will respond. See the **InTouch** manual for more information.

Applying the linens

⚠ WARNING

- Always make sure that the operator has access to the CPR straps.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Do not use fitted sheets with this support surface.

To apply the linens:

1. Apply the linens using the “D” rings for the flat sheet.
2. Thread the four linen corners through the “D” rings (C) attached to the support that are surface to secure the linens (A) to the support surface (B).


Figure 1-2: Apply linens

Note: Do not pull the linens tight. Keep the linens loose and as smooth as possible on top of the support surface to effectively use the Turn Assist and Lateral Rotation functions.

Positioning a patient on the support surface

WARNING

- Always center the patient on the support surface. Align the patient's head toward the headboard before starting functions. Check the patient frequently to make sure that you maintain the proper positioning.
- Always make sure that the tubing and wiring that is connected to the patient is long enough, stable, and secure during Lateral Rotation or Turn Assist.
- Always raise all of the **InTouch** bed siderails before starting Turn Assist or Lateral Rotation functions.
- Do not exceed the safe working load of the **Isolibrium** support surface. Excess weight could cause unpredictable safety and performance of this system.
- Always use extra caution when reading radiology images taken of a patient on a support surface because internal components can cause artifacts and distort readings.

CAUTION

- Do not allow sharp objects to come into contact with the support surface that could puncture, tear, or cut the cover.
- Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillow case or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.

To position the patient:

1. Max Inflate the support surface.
2. Position the patient in the center of the support surface, align the patient's head toward the head board ([Figure 1-3 on page 1-15](#)).
3. Check the patient frequently during Lateral Rotation for proper positioning and support surface inflation ([Figure 1-4 on page 1-15](#) and [Figure 1-5 on page 1-15](#)).


Figure 1-3: Center the patient

Figure 1-4: Correct patient positioning

Figure 1-5: Incorrect patient positioning

Managing incontinence and drainage

WARNING

Always monitor the patient condition at regular intervals for patient safety.

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

Activating and resetting CPR

When you activate for CPR, the system deflates and the **InTouch** control display will change.

WARNING

Always make sure that the operator has access to the CPR straps.

To activate for **CPR**, pull the red CPR strap straight out from the support surface ([Figure 1-6 on page 1-16](#)). The CPR straps are located on the left and right side of the foot end of the support surface. You can activate either one or both CPR straps.

Figure 1-6: CPR activation

Figure 1-7: CPR reset

To reset for **CPR**, pull down on red CPR strap ([Figure 1-7 on page 1-16](#)). The support surface will begin to reinflate.

Relocating the support surface

Relocating the support surface is for the purpose of cleaning or moving to another InTouch bed frame.

WARNING

- Do not move the support surface without assistance.
 - Do not use the support surface as a transfer device to avoid the risk of patient injury.
-

Prerequisite: Minimum of two operators are required to relocate the support surface and setup this product.

1. Make sure that there is not a patient on the support surface.
2. Remove the linens.
3. Unplug the power cord for the support surface from the foot end auxiliary power outlet.
4. Unplug the support surface connector.
5. Remove the headboard and footboard
6. Position an operator at each end of the support surface and remove the support surface and relocate.

Note: Use all four handles when relocating the support surface.

Note: See ([Installation on page 1-12](#)) of the support surface.

7. Reinstall the **InTouch** headboard and footboard.

Transferring of patient from one patient support platform to another

WARNING

- Do not use the support surface as a transfer device to avoid the risk of patient injury.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not exceed the safe working load of the **InTouch** bed frame when supporting both the patient and the support surface. Excess weight could cause unpredictable safety and performance of this product.
-

To transfer the patient from one patient support surface to another:

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

1. Position one patient support platform alongside the other patient support platform while minimizing the gap between the two platforms.
2. Set the brakes to on for both patient support platforms.
3. Activate maximum inflation of the support surface until the inflation level is reached. See the **InTouch** Starting and Stopping Max Inflate section of the manual for use instructions.
4. Adjust the patient support platform heights so that they are level with one another.
5. Transfer the patient following all applicable safety rules and institution protocols for patient and operator safety.
6. Turn off Max Inflate.

Storing the product

Before storing the product, deflate the support surface by activating CPR.

Cleaning

WARNING

- Always unplug the support surface power cord before cleaning or disinfecting to avoid the risk of shock.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - Always perform preventive maintenance based on the usage level of the product. An increase in usage which may include more frequent cleaning and disinfection, may adversely affect the life of the support surface.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
-

CAUTION

- Do not iron, dry-clean, or tumble dry the support surface covers.
 - Do not power wash the support surface.
 - Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
-

Prerequisite: Minimum of two operators are required to clean the bottom support surface cover.

Always follow hospital protocol for cleaning and disinfecting.

To clean the support surface covers between patient uses, follow these steps in order:

1. If powered, unplug the **InTouch** from the wall outlet and turn the battery switch to off.
2. Unplug the support surface power cord.
3. Disconnect the support surface integration cable.
4. Using a clean, soft, damp cloth, wipe the support surface covers with a mild soap and water solution to remove foreign material.
5. Wipe the support surface covers with a clean, dry cloth to remove any excess liquid or cleaning agent.
Note: To clean the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.
6. Thoroughly rinse and dry the support surface covers.

Disinfecting

EN

WARNING

- Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
 - Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning or disinfecting. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, a corrosive residue may be left on the surface of the product that could cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always unplug the product power cord before cleaning or disinfecting to avoid the risk of shock.
-

CAUTION

- Always completely dry the support surface covers before storing, adding linens, or placing a patient on the surface.
 - Do not over expose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers as they may damage the cover or the CPR strap print.
-

Prerequisite: Minimum of two operators are required to disinfect the support surface.

Suggested disinfectants:

- Quaternaries
- Phenolic Disinfectants
- Chlorinated Bleach Solution (5.25% bleach diluted 1 part bleach to 100 parts water)
- 70% Isopropyl Alcohol

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Unplug the InTouch power cord.
2. Unplug the support surface.
3. Disconnect the integration cable.
4. Thoroughly clean and dry the support surface covers (see [Cleaning on page 1-19](#)) before disinfectants are applied.
5. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes (do not soak the support surface).

Notes

- Make sure that you follow the disinfectant manufacturer's instructions for appropriate contact time and rinsing requirements.
 - To disinfect the bottom cover of the support surface, raise the head of bed and foot end to help fold the head end and foot end of the support surface.
6. Rinse the support surface covers with a clean, damp cloth to remove excess disinfectant.
 7. Wipe the support surface covers with a clean, dry cloth to remove any excess liquid or disinfectant.
 8. Allow the support surface covers to dry completely before returning to service.

Preventive maintenance

WARNING

Do not service or perform maintenance while the product is in use.

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage. Service only by qualified personnel.

Remove product from service before you perform preventive maintenance inspection.

Note: Clean and disinfect the exterior of the support surface before inspection, if applicable.

Inspect the following items:

_____ Zipper, snaps, and covers (top, bottom, and fire barrier) are free of tears, cuts, holes, or any other damage
Note: If excessive wear is observed on the cover or the fire barrier, it is strongly recommended to replace the cover.

_____ Support surface cover labels are legible, have proper adherence, and free of damage.

_____ Handles are free of rips or cracks

_____ Power cord, integration cable, and plug are free of damage

_____ **Perform a functional test of the Isolibrium functions

_____ **Perform system diagnostics and confirm there are no errors

_____ *Internal components for signs of staining from fluid ingress or contamination by fully unzipping the cover

_____ Pods are not leaking or cracking

_____ Hose connections to the manifold and sensor tubes (both ends) are seated fully

_____ Foam is free from large tears or large gouges

_____ Turn bladders are free of excessive wear and function properly

_____ Low Air Loss fans rotate freely and are free of debris and dust

_____ Left and Right CPR releases function properly

_____ Both foot box cooling fans rotate freely and are free of debris and dust

_____ Pump spring mount has visible clearance from the foot box bottom

_____ Listen for abnormal pump vibration and noise level

_____ Clean the foot box intake metal screen

_____ Clean the mesh on the bottom cover

_____ Clean the fan filter screen on the foot box

_____ Standoffs (feet of the foot box) are present

_____ Cable tie is present on the power cord

_____ Current leakage not more than 200 μ A

_____ Ground impedance not more than 100 m Ω ([Figure 1-8 on page 1-21](#))

Note: The ground lug can be accessed through the slit in the cover for the CPR, patient left.


Figure 1-8: Ground lug (patient left side)

Preventive maintenance

EN

Notes

- *A zipper pull tool or equivalent is required to access the internal components of the support surface.
- **You cannot perform these preventive maintenance checks without the **InTouch** footboard. See the **InTouch** integrated Models 2131 and 2141, software version 4.0 or higher manual.
- The zipper overlaps at the end and may appear to be misaligned.

Replace the following items every two years:

- Cover assembly, top and bottom

Product serial number:		
Completed by:		Date:

Quick reference replacement parts


These parts are currently available for purchase. Call Stryker Customer Service: 1-800-327-0770 for availability and pricing.

EN

Part name	Part number
Board, main	2971-022-140
Board, sensor	2971-021-033
Cover, bottom	2971-700-008
Cover, top	2971-021-004
Fan, foot box	2971-022-188
Fan, LAL	2971-021-189
Integration cable	2971-022-133
Kit, power cord	2971-700-003
Kit, pump assembly	2971-700-001
Kit, resonator assembly	2971-700-002
Kit, valve manifold	2971-700-007
Pods assembly	2971-021-003
Power supply	2971-022-171
Kit, solenoid valve, SLI, BLI, BLD, BRI, BRD, SRD, SRI, SLD	2971-700-004
Kit, solenoid valve, TRI, TRD, EX, TRI, TU	2971-700-005
Kit, solenoid valve, IN	2971-700-006

Isolbrium wiring

Block diagram


Block diagram

Pneumatic

EN

ISOLIBRIUM PNEUMATIC DIAGRAM


Circuit boards

EN

Power supply assembly


2971-022-171 Rev J (Reference only)

Circuit boards

EN

Power supply assembly (Continued)


Cable location	Voltage	Positive lead	Negative lead	Description
SK1	120AC	Pin 1 brown	Pin 2 blue	120VAC power in to the power supply
SK2 (A)	+5VDC	Pin 1 red	Pin 3 green	+5VDC from the power supply for the main board
SK2	+12VDC	Pin 2 black	Pin 3 green	+12VDC from the power supply for the main board
SK2	+24VDC	Pin 7 blue	Pin 3 green	+24VDC from the power supply
SK2	-12VDC	Pin 8 blue	Pin 6 red	-12VDC from the power supply for the foot box fans

Circuit boards

EN

Main power board assembly

2971-022-140 Rev G (Reference only)

Circuit boards

EN

Main power board assembly (Continued)

Cable and test point (LED) location	Voltage	Positive lead	Negative lead	Description
D43	+5VDC	N/A	N/A	+5VDC in from power supply to main board
D60	+12VDC	N/A	N/A	+12VDC in from power supply to main board
D56	+24VDC	N/A	N/A	+24VDC in from power supply to main board
CPR (left)	+5VDC	Pin 2 red	Pin 3 black	+5VDC power out to CPR switch
CPR (right)	+5VDC	Pin 2 red	Pin 3 black	+5VDC power out to CPR switch
TLI (D30)	+24VDC	TP11	TP24 GND	+24VDC from the main board to solenoid valve (Turn Left Inflate)
TLD (D31)	+24VDC	TP18	TP24 GND	+24VDC from the main board to solenoid valve (Turn Left Deflate)
TRI (D26)	+24VDC	TP9	TP24 GND	+24VDC from the main board to solenoid valve (Turn Right Inflate)
TRD (D27)	+24VDC	TP16	TP24 GND	+24VDC from the main board to solenoid valve (Turn Right Deflate)
BLI (D34)	+24VDC	TP13	TP24 GND	+24VDC from the main board to solenoid valve (Back Left Inflate)
BLD (D36)	+24VDC	TP20	TP24 GND	+24VDC from the main board to solenoid valve (Back Left Deflate)
BRI (D22)	+24VDC	TP7	TP24 GND	+24VDC from the main board to solenoid valve (Back Right Inflate)
BRD (D23)	+24VDC	TP14	TP24 GND	+24VDC from the main board to solenoid valve (Back Right Deflate)

Circuit boards

EN

Main power board assembly (Continued)


Cable and test point (LED) location	Voltage	Positive lead	Negative lead	Description
SLI (D33)	+24VDC	TP12	TP24 GND	+24VDC from the main board to solenoid valve (Seat Left Inflate)
SLD (D32)	+24VDC	TP19	TP24 GND	+24VDC from the main board to solenoid valve (Seat Left Deflate)
SRI (D25)	+24VDC	TP8	TP24 GND	+24VDC from the main board to solenoid valve (Seat Right Inflate)
SRD (D24)	+24VDC	TP15	TP24 GND	+24VDC from the main board to solenoid valve (Seat Right Deflate)
IN (D28)	+24VDC	TP10	TP24 GND	+24VDC from the main board to solenoid valve (Air Intake)
EX (D29)	+24VDC	TP17	TP24 GND	+24VDC from the main board to solenoid valve (Air Exhaust)
J14	+12VDC	Pin 1 red	Pin 6 black	+12VDC out to sensor board from main board
J6	+24VDC	Pin 1 yellow	Pin 2 black	+24VDC out to compressor from the main board (variable voltage) - D65 = 24.4% - D66 = 48.8% - D67 = 73.2% - D68 = 97.6% - D71 = 110%

Cables

EN

Fans and cable assembly

2971-021-189 Rev C (Reference only)

Item	Color
1	Red
2	Blue

SPI cables (main board to foot box) assembly

2971-022-132 Rev D (Reference only)


Item	Color
1	Red
6	Black
2	Black
3	Green
4	Black
5	White
7	Blue
8	Black
9	Brown
10	Black

Cables

EN

Integration cable assembly 2971-022-133

Rev C

Item	Color
1	Black
2	Red
4	Green
6	open

Protecting against Electrostatic Discharge (ESD)

 **CAUTION**

- Always use ESD protective equipment before opening antistatic bags and servicing electronic parts.
 - Do not place unprotected circuit boards on the floor.
-

Note: Always ship back circuit boards to Stryker in the same antistatic bags that the new boards were originally shipped in.

The electronic circuits in the product are completely protected from static electricity damage when factory assembled. Always use adequate static protection when servicing the electronic systems of the product. All service personnel must use static protection whenever they are touching wires.

Sample antistatic protection equipment includes:

- 1 antistatic wrist strap
- 1 grounding plug
- 1 test lead with a banana plug on one end and an alligator clip on the other end

Make sure that you follow the ESD manufacturer's instructions for appropriate protection against static discharge.

Cover replacement, top

Tools required:

- Zipper pull tool or equivalent

Procedure:

1. Apply the **InTouch** brakes.
2. Unplug the support surface power cord.
3. Remove the headboard and footboards.
4. Lower all the siderails.
5. Unsnap the two secure snaps (C) at the foot end corners (Figure 1-9 on page 1-34).
6. Unsnap the four corner retainers (B) (Figure 1-9 on page 1-34).
7. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
8. Unzip the support surface top cover from the bottom using the zipper pull tool or equivalent (A) (Figure 1-9 on page 1-34).

Note: The zipper starts on the patient right side near the foot end.


Figure 1-9: Top cover

9. Remove and discard the top cover.
10. Reverse to install the new cover.

Notes

- After installation, remove the zipper pull tool or equivalent from the zipper.
 - Connect the corner retainers by placing the strap below each D-ring and snaps to the cover. Repeat for the other three corners.
 - Make sure that the watershed covers the zipper.
11. Verify proper operation before returning the product to service.

Cover replacement, bottom

Tools required:

- Zipper pull tool or equivalent

Prerequisite: Minimum of two operators are required to replace the bottom support surface cover.

Procedure:

1. Apply the **InTouch** brakes.
2. Unplug the support surface power cord (D) (Figure 1-10 on page 1-35).
3. Disconnect the support surface integration cable (E) (Figure 1-10 on page 1-35) from the bed.
4. Make sure that the **InTouch** platform is flat.
5. If powered, unplug the **InTouch** from the wall outlet and turn the battery switch to off.
6. Remove the headboard and footboards.
7. Lower all the siderails.
8. Before turning, place a protective sheet on the surface that the support surface top cover will rest.
9. Using two operators, turn the support surface over (Figure 1-10 on page 1-35).


Figure 1-10: Bottom support surface cover

10. Unsnap the four corner cover snaps.
11. Unsnap the two secure snaps at the right and left foot end corners.
12. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
13. Grasp the zipper pull tool or equivalent and pull to unzip the bottom cover from the top cover (B) (Figure 1-10 on page 1-35).
14. Feed the CPR release straps (C), power cord (D), and integration cable (E) through the bottom cover (Figure 1-10 on page 1-35).
15. Lift the foot box assembly slightly toward the head end and separate the two foot box Velcro® retainer straps from the top of the foot box cover.

Note: During reinstall, make sure to align the foot box with the lines on the bottom cover.
16. Starting at the head end of the support surface, separate the six Velcro® patches from the bottom cover (F) (Figure 1-11 on page 1-36).
17. In a gentle motion, unsnap the twelve snaps and separate the Velcro® from the left and right sides of the turn bladder (G) (Figure 1-11 on page 1-36).

Cover replacement, bottom (Continued)


Figure 1-11: Pod Velcro® and snap retainers

18. Detach the foot end foam crib from the bottom cover.
19. Disconnect the exhaust tube quick disconnect (H) from the bladder assembly ([Figure 1-11 on page 1-36](#)).
20. Remove and discard the bottom cover.
21. Reverse to install the new bottom cover.

Notes

- During reinstall, make sure to align the Velcro® patches, snaps and the foot box cover.
 - Connect the corner retainers by placing the strap below each D-ring and snaps to the cover. Repeat for the other three corners.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
22. Verify proper operation before returning the product to service.

Pod assembly replacement

Tools required:

- Zipper pull tool or equivalent

Procedure:

1. Remove the top cover. See [Cover replacement, top on page 1-34](#) for removal.
2. Unsnap the pillow (H) from the left and right side of the pod assembly ([Figure 1-12 on page 1-37](#)). Save the pillow.
3. Separate the six Velcro® connectors holding the foot section gel crib to the pod assembly.
4. Remove the foot section gel crib (F) by lifting upward ([Figure 1-12 on page 1-37](#)). Save the gel crib.


Figure 1-12: Gel crib, pillow, and quick disconnect fitting

5. Disconnect the quick disconnect fitting (Z) ([Figure 1-12 on page 1-37](#)).
6. Remove the two retaining pins (A) from the hose quick connection from the foot box assembly ([Figure 1-13 on page 1-38](#)). Save the pins.

Pod assembly replacement (Continued)


Figure 1-13: Retaining pins

7. Grasp both sides of the hose quick connection from the pod assembly and carefully pull outward toward the head end of the support surface to disconnect from the foot box manifold.
8. Disconnect the four pod sensor hoses (Y) (1 blue, 2 yellow, 4 green, and 5 orange) from the foot box. Push inward on the coupling and pull outward on each of the hoses ([Figure 1-14 on page 1-38](#)).

Notes

- When handling the pod sensor hoses, do not bend or kink the hoses.
- Pay attention to the sensor hose position and insertion color and number labels.
- During re-install, make sure to insert the hose and once you hit a stop, continue pushing until you reach a hard stop.


Figure 1-14: Pod sensor hoses and purse lock wire tie

9. Remove sensor hoses from the purse lock wire ties (M) ([Figure 1-14 on page 1-38](#)).
10. Reach under the pod assembly (C) from the patient's right side near the middle to access the turn bladder and disconnect the turn bladder quick disconnect. Repeat on the patient left side ([Figure 1-15 on page 1-39](#)).

Pod assembly replacement (Continued)


Figure 1-15: Pod assembly

11. Separate the Velcro® from the head end of the pod assembly (C) and feed the turn bladder hose quick connection through the foam crib.
12. Remove and discard the pod assembly (C)([Figure 1-15 on page 1-39](#)).
13. Reverse to install the new pod assembly.

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
14. Run the leak diagnostic test (see *Accessing the Isolibrium Diagnostic Menu* in the **InTouch** maintenance manual).
 15. Verify proper operation before returning the product to service.

Foot box cover access

Tools required:

- Zipper pull tool or equivalent
- T20 Torx driver
- T15 Torx driver

Procedure:

1. Apply the **InTouch** brakes.
2. Unplug the support surface power cord.
3. Remove the headboard and footboards.
4. Lower all the siderails.
5. Unsnap the two secure snaps (C) at the foot end corners (Figure 1-16 on page 1-40).
6. Unsnap the two foot end corner retainers (B) (Figure 1-16 on page 1-40).
7. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
8. Unzip the support surface top cover from the bottom cover using the zipper pull tool or equivalent (A) (Figure 1-16 on page 1-40) and stop unzipping at the middle of the opposite side (D).

Note: The zipper starts on the patient right side near the foot end.


Figure 1-16: Top cover

9. Fold the cover up toward the head end.
10. Separate the six Velcro® connectors holding the foot section gel crib to the pod assembly.
11. Remove the foot section gel crib (F) by lifting upward (Figure 1-17 on page 1-41). Save gel crib.
12. Remove the foot box cover Velcro® retainer straps (E) from the foot box cover (Figure 1-18 on page 1-42).

Foot box cover access (Continued)


Figure 1-17: Gel crib

13. Using a T15 Torx driver, remove the nine ground screws (B) and nine star washers (D) (Figure 1-16 on page 1-40). Save the ground screws and star washers.

Notes

- During reinstall, tighten the screws equally and do not over tighten. Torque (B) to 22 ± 4 in-lb.
- During reinstall, check the LAL fan cable (G) and grommet (H) are in place before replacing the foot box cover (1-18 Foot box cover on page 1-42).

14. Remove the foot box cover (AF) (Figure 1-18 on page 1-42). Save the foot box cover.

Foot box cover access (Continued)


Figure 1-18: Foot box cover

15. Reverse to reinstall.

Solenoid valve replacement

Tools required:

- Zipper pull tool or equivalent
- Torque driver
- ESD system

Procedure:

1. See ([Foot box cover access on page 1-40](#)).
2. Properly ground yourself (see [Protecting against Electrostatic Discharge \(ESD\) on page 1-33](#)).
3. Using the supplied torque driver, remove the two screws that secure the solenoid valve (H, J, or K) to the manifold assembly ([Figure 1-19 on page 1-43](#)). Save the screws.

Note:

During installation tighten the screws equally. Use the supplied torque driver and torque the screws to 4.4 in-lb.


Figure 1-19: Solenoid valve

4. Remove the solenoid valve.
5. Disconnect the connector for the solenoid valve.

Note: Inspect the replacement solenoid valve to make sure that the O-ring is in place. If the O-ring is not in place, transfer the O-ring from the removed solenoid valve to the supplied solenoid valve.
6. Discard the solenoid valve.
7. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
8. Verify proper operation before returning the product to service.

Power supply replacement

Tools required:

- Zipper pull tool or equivalent
- T10 Torx driver
- ESD system

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Properly ground yourself (see [Protecting against Electrostatic Discharge \(ESD\) on page 1-33](#)).
3. Remove the wires for the fans and DC power from the wire clip (A).
4. Disconnect the DC power connector from the main board (AR) (connector R) ([Figure 1-20 on page 1-44](#)).


Figure 1-20: Power supply

5. Disconnect the AC power connector (B) from the power supply ([Figure 1-20 on page 1-44](#)).
6. Disconnect the two foot box fans (D).
7. Disconnect the ground wire connector from the power supply.
8. Using a T10 Torx driver, remove the four screws (C) that secure the power supply (AY) to the foot box ([Figure 1-20 on page 1-44](#)). Save the screws.

Note: During reinstall, torque (C) to 11 ± 2 in-lb.

Power supply replacement (Continued)

9. Remove and discard the power supply (AY) ([Figure 1-20 on page 1-44](#)).
10. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
 - Do not dispose of as unsorted municipal waste. See your local distributor for return or collection systems available in your country.
11. Verify proper operation before returning the product to service.

Main board replacement

Tools required:

- Zipper pull tool or equivalent
- Needle nose pliers
- ESD system

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Properly ground yourself (see [Protecting against Electrostatic Discharge \(ESD\) on page 1-33](#)).
3. Remove all wiring connections from the main board (AR) ([Figure 1-21 on page 1-46](#)).

Note: Pay attention to all connection points before disconnecting.


Figure 1-21: Main board

4. Using needle nose pliers, apply a slight upward pressure on the main board to unclip the seven plastic standoffs.
Note: During reinstall, align the standoffs to the main board and press near the standoffs to seat the main board. Do not bend the board.
5. Remove and discard the main board.
6. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
- After installation, remove the zipper pull tool or equivalent from the zipper.
- Make sure that the watershed covers the zipper.

Main board replacement (Continued)

- Do not dispose of as unsorted municipal waste. See your local distributor for return or collection systems available in your country.
7. Verify proper operation before returning the product to service.

Power cord replacement

Tools required:

- Zipper pull tool or equivalent
- 5/16" Nut driver
- 1" Combination wrench
- Diagonal pliers

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Using diagonal pliers, cut the cable tie (BG) ([Figure 1-22 on page 1-48](#)) and discard.

Notes

- Pay attention to the orientation of the cable tie for installation of the new power cord.
- From the inside of the foot box, install the supplied cable tie through the top hole with the cable tie lock facing upward. Then feed the cable tie around the new power cord and through the bottom hole in the foot box to secure the power cord to the foot box.
- During installation, before cutting the end off of the cable tie, use the diagonal pliers to pry the slack out of the cable tie (one to two clicks).


Figure 1-22: Power cord

3. Using a 5/16" Nut driver, remove the two nuts (H) and four star washers (D) that secures the power cable ground connector ([Figure 1-22 on page 1-48](#)). Save the nuts and star washers.
Note: During reinstall, torque (H) to 22 ± 4 in-lb
4. Using a 1" combination wrench, remove the nut (E) securing the power cord strain relief to the foot box.
5. Feed the ground through the plastic retaining nut (E) ([Figure 1-22 on page 1-48](#)).

Power cord replacement (Continued)

6. Disconnect the power cord quick connect (F) from the power supply ([Figure 1-22 on page 1-48](#)).
7. Feed the two pin connector through the nut (E) ([Figure 1-22 on page 1-48](#)).
8. Remove and discard the power cord (AV) ([Figure 1-22 on page 1-48](#)).
9. Reverse steps to reinstall.

Notes

- During reinstall, first assemble the power cord ground star washers and nut, and then the power supply ground star washers and nut.
 - During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
10. Verify proper operation before returning the product to service.

Bed integration cable connector replacement

Tools required:

- Zipper pull tool or equivalent
- 1" Combination wrench or socket

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Using your fingers, disconnect the bed integration cable connector (AP) from the foot box by rotating the turn-lock counterclockwise ([Figure 1-23 on page 1-50](#)).
3. Using a 1" combination wrench, remove the plastic retaining nut (F) on the bed integration cable connector securing it to the foot box. Save the nut.
4. Disconnect the can inner cable (AT) from the main board (connector J) ([Figure 1-23 on page 1-50](#)).


Figure 1-23: Bed integration cable connector

5. Remove and discard the bed integration cable connector.
6. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
7. Verify proper operation before returning the product to service.

Pump assembly replacement

Tools required:

- Zipper pull tool or equivalent
- Diagonal pliers
- T15 Torx driver

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Disconnect the pump power connector from the main board (connector F).
3. Disconnect the cable from the three wire retainers.
4. Using diagonal pliers, carefully cut the cable ties where the inlet hose connects to the manifold (BG) and the exhaust hose connects to the resonator ([Figure 1-24 on page 1-51](#)). Do not cut the manifold barb.

Note: During installation, before cutting the end off of the cable tie, use the diagonal pliers to pry the slack out of the cable tie (one to two clicks).


Figure 1-24: Pump assembly

5. Using a T15 Torx driver, remove the four screws (A) that secure the pump assembly (AJ) to the foot box ([Figure 1-24 on page 1-51](#)). Save the screws.

Pump assembly replacement (Continued)

Note: During reinstall, torque (A) to 22 ± 4 in-lb.

6. Move the pump assembly to the side so the area in front of the manifold is clear to remove the hose from the manifold.
7. With a firm grasp on the manifold hose, pull straight outward while twisting.
Note: Be careful when removing and installing to not damage to the manifold barb.
8. Grasp the resonator hose and hold the resonator while pulling outward on the hose.
9. Remove and discard the pump assembly.
10. Reverse steps to reinstall.

Notes

- During reinstall, make sure that you hold the pump springs while tightening so the springs do not turn.
 - During reinstall, make sure that you check the pump hose is not kinked or bent.
 - During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
 - Do not dispose of as unsorted municipal waste. See your local distributor for return or collection systems available in your country.
11. Verify proper operation before returning the product to service.

Low Air Loss (LAL) fan replacement

Tools required:

- Zipper pull tool or equivalent

Procedure:

1. Apply the **InTouch** brakes.
2. Unplug the support surface power cord.
3. Remove the headboard and footboards.
4. Lower all the siderails.
5. Unsnap the two secure snaps (C) at the foot end corners ([Figure 1-25 on page 1-53](#)).
6. Unsnap the two foot end corner retainers (B) ([Figure 1-25 on page 1-53](#)).
7. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
8. Unzip the support surface top cover from the bottom cover using the zipper pull tool or equivalent (A) ([Figure 1-25 on page 1-53](#)) and stop unzipping at the middle of the opposite side (D).

Note: The zipper starts on the patient right side near the foot end.


Figure 1-25: Top cover

9. Fold the cover up toward the head end.
10. Separate the six Velcro® connectors holding the foot section gel crib to the pod assembly.
11. Remove the foot section gel crib (F) by lifting upward ([Figure 1-26 on page 1-54](#)). Save gel crib.

Low Air Loss (LAL) fan replacement (Continued)


Figure 1-26: Gel crib

12. Remove the two purse lock wire tie (M) (Figure 1-27 on page 1-54). Save the purse lock wire ties.


Figure 1-27: Purse lock wire tie

13. Disconnect the LAL fan from the quick connection.
14. Remove the LAL fan from the LAL fan nest and discard the LAL fan (J) (Figure 1-28 on page 1-55).
Note: During reinstall, route the cable back through the foam and make sure that the new fan is seated into the LAL fan nest.

Low Air Loss (LAL) fan replacement (Continued)


Figure 1-28: LAL fan

15. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
- After installation, remove the zipper pull tool or equivalent from the zipper.
- Make sure that the watershed covers the zipper.

16. Verify proper operation before returning the product to service.

Sensor board replacement

Tools required:

- Zipper pull tool or equivalent
- Needle nose pliers
- ESD system

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Properly ground yourself (see [Protecting against Electrostatic Discharge \(ESD\) on page 1-33](#)).
3. Disconnect all six of the pod sensor hoses from the foot box by gently pushing inward on the coupling and pulling outward on each of the hoses.

Notes

- When handling the pod sensor hoses, do not bend or kink the hoses.
 - Pay attention to the sensor hose position and insertion color and number labels.
 - During re-install, make sure to insert the hose and once you hit a stop, continue pushing until you reach a hard stop.
4. Disconnect the cable (AN) from the sensor board (AB) ([Figure 1-30 on page 1-57](#)).
 5. Using needle nose pliers, lifting up slightly, disconnect the sensor board (AB) from the three plastic standoffs ([Figure 1-30 on page 1-57](#)).

Note: During reinstall, align the standoffs to the main board and press near the standoffs to seat the main board. Do not bend the board.


Figure 1-30: Sensor board

6. Remove and discard the sensor board.
7. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
- After installation, remove the zipper pull tool or equivalent from the zipper.
- Make sure that the watershed covers the zipper.

Sensor board replacement (Continued)

- Do not dispose of as unsorted municipal waste. See your local distributor for return or collection systems available in your country.
8. Run the leak diagnostic test (see *Accessing the Isolibrium Diagnostic Menu* in the **InTouch** maintenance manual).
 9. Verify proper operation before returning the product to service.

Resonator replacement

Tools required:

- Zipper pull tool or equivalent
- Diagonal pliers

Procedure:

1. See [Foot box cover access on page 1-40](#).
2. Using diagonal pliers, cut the cable ties (BG) from the inlet and outlet hoses of the resonator and discard the cable ties ([Figure 1-31 on page 1-59](#)).

Note: During installation, before cutting the end off of the cable tie, use the diagonal pliers to pry the slack out of the cable tie (one to two clicks).


Figure 1-31: Resonator

3. Pull to disconnect each of the hoses from the resonator (AG) and the manifold ([Figure 1-31 on page 1-59](#)).
4. Remove and discard the resonator.

Note: During reinstall, make sure that the compressor supply hose is not kinked or bent.

5. Reverse steps to reinstall.

Resonator replacement (Continued)

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
6. Verify proper operation before returning the product to service.

Turn bladder assembly replacement

Tools required:

- Zipper pull tool or equivalent
- Diagonal pliers

Procedure:

1. Remove the top cover. See [Cover replacement, top on page 1-34](#) for removal.
2. Separate the six Velcro® connectors holding the foot section gel crib to the pod assembly.
3. Remove the foot section gel crib (F) ([Figure 1-32 on page 1-61](#)). Save the gel crib.


Figure 1-32: Gel crib

4. Lift up on the bottom corner of the pods assembly and foam crib and disconnect the turn bladder quick disconnect. Repeat for the opposite side.
5. Remove the turn bladder red and black sensor hoses from the foot box. Push in on the sensor hose quick connector and then pull outward on the hose.
Note: Do not kink or bend the hoses.
6. Remove the hoses from the purse clips.
7. Lift up on the foam crib, reach between the foam crib and the bottom cover to unhook the six Velcro® retainers (A) and unsnap the 12 snap retainers (B) ([Figure 1-33 on page 1-62](#)). Repeat for the opposite side.

Turn bladder assembly replacement (Continued)


Figure 1-33: Pod Velcro® and snap retainers

8. Remove and discard the turn bladder assembly.
9. Reverse steps to reinstall.

Notes

- During reinstall, make sure to align the Velcro® patches.
 - After installation, remove the zipper pull tool or equivalent from the zipper.
 - Make sure that the watershed covers the zipper.
10. Run the leak diagnostic test (see *Accessing the Isolibrium Diagnostic Menu* in the **InTouch** maintenance manual).
 11. Verify proper operation before returning the product to service.

Product labels

All the product labels listed (Figure 1-34 on page 1-63, Figure 1-35 on page 1-64, Figure 1-36 on page 1-65, and Figure 1-37 on page 1-65) are attached to the **Isolibrium** support surface.

Note: These labels (Figure 1-34 on page 1-63) are required for standards in the United States and Canada and are in English and French Canadian. They are not required to be translated into any other language.

①
P/N 2971-021-902 REV D.0

**UNDER PENALTY OF LAW
THIS TAG NOT TO BE REMOVED
EXCEPT BY THE CONSUMER**

ALL NEW MATERIAL
Consisting of

FOAM.....56%

DRY POLYMER GEL.....26%

POLYURETHANE FILM.....18%

REG. NO. UT-3215 (MI)

Certification is made by the manufacturer that the materials in this article are described in accordance with law.

STRYKER MEDICAL
3800 E. Centre Avenue
Portage, MI 49002-5826 USA

Finished Size: 84 in x 35 in x 8 in
(213.4 cm x 88.9 cm x 20.3 cm)
Net Weight of Filling Mat'l: 50 lb, 0 oz

MADE IN U.S.A.

Model: 2971-ISOLIBRIUM
Prototype ID: 2390428
This mattress meets the requirements of 16 CFR Part 1633 (federal Flammability (open flame) standard for mattress sets) when used without a foundation.

**THIS MATTRESS IS
INTENDED TO BE USED
WITHOUT A FOUNDATION**

Not to be removed until delivered to the consumer.

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This article contains
NEW MATERIAL ONLY

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REF 2971 - ISOLIBRIUM™

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460 lb/
208.6 kg

VALID WITH FLAME BARRIER ONLY
CONFORMS WITH FLAMMABILITY STANDARD USA-16
CFR 1633, 1632 BOSTON-BFD IX-11, CANADA-
METHOD 27.7-1979 of CAN 2-4.2 M77, CALIFORNIA-TB
129 WHEN USED WITHOUT A FOUNDATION.

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Information presented in this section is not to be translated because they are designed for standards in the United States

This product is made with Intelli-Gel®* hollow column configuration and material.*

*Intelli-Gel® is a registered trademark of EdiZONE, LLC of Alpine, UT
DuraGel™ is a trademark of EdiZONE, LLC of Alpine, UT

Covered by one or more of the following:
U.S. Patents Nos.
5,749,111; 6,026,527; 7,076,822; 7,964,664

Patent Pending

WARNING
AVERTISSEMENT

PATIENT STABILITY AND SIDERAILS COVERAGE MAY BE COMPROMISED WITH THE USE OF AN OVERLAY

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WARNING

DO NOT TRANSFER PATIENT FROM ONE BED TO ANOTHER USING THE SUPPORT SURFACE WITH A PATIENT ON IT

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SERIAL NUMBER IS LOCATED INSIDE THE MATTRESS AT THE FOOT END

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Troubleshooting

EN

Problem	Possible Cause	Solution
No power	The support surface power cord is not plugged into the InTouch bed auxiliary outlet	Plug the support surface power cord into the InTouch bed auxiliary outlet
	The support surface power cord is plugged into the InTouch bed auxiliary outlet	Check that the InTouch bed auxiliary outlet breaker is not tripped
	The InTouch bed power cord is not plugged into a hospital grade outlet	Plug the InTouch power cord into a hospital grade outlet
		Using a voltmeter, check the auxiliary outlet on the InTouch bed for 120VAC
	The support surface power supply is not supplying DC power	Voltage (120VAC) is not present on the AC connector on the support surface power supply (see Power supply assembly on page 1-26)
<p>Check the VDC voltage on the support surface terminal block on the power supply:</p> <ol style="list-style-type: none"> 1. Pin 1 (red) and pin 3 (green) = 5VDC 2. Pin 2 (black) and pin 3 (green) = 12VDC 3. Pin 7 (blue) and pin 3 (green) = 24VDC <p>Note: If no low voltage is present, replace the power supply.</p>		

Troubleshooting

EN

Support surface will not inflate	The support surface has no power	See no power
	A CPR pull is active	Pull down the CPR to reset
	The support surface menu is available on the InTouch bed screen	Check the connection of the support surface integration cable (see Installation on page 1-12)
	Pod manifold connector is not set	Set the pod manifold connector and secure
	The pump is not running when you start Max Inflate	<p>Check for 24VDC (variable) is going to the pump on connector F on the main board pin 1 (yellow) and pin 2 (black) (see Power supply assembly on page 1-26)</p> <ol style="list-style-type: none"> If voltage is present, replace pump assembly If no voltage present, replace the control board
	Possible disconnect of hose connections or solenoids or pod assembly	Check the hose connections, solenoids, and pod assembly
No left turn inflate	Loose or disconnection of the left turn bladder hose quick connection coupler	Make sure that the left turn bladder hose quick connection coupler is together and locked
	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and set
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
		If you find a leak, replace the turn bladder assembly
	No air entering the turn bladder	Make sure that the LED (D30) is lit to show the TLI solenoid is energizing (see Main power board assembly on page 1-28)
If energized, replace the TLI inflate solenoid		
No left turn deflate	Loose or disconnection of the left turn bladder hose quick connection coupler	Make sure that the left turn bladder hose quick connection coupler is together and locked
	No air is exiting the turn bladder	Make sure that the LED (D31) is lit to show the TLD solenoid is energizing (see Main power board assembly on page 1-28)
		If energized, replace the TLD inflate solenoid

Troubleshooting

EN

No right turn inflate	Loose or disconnection of the right turn bladder hose quick connection coupler	Make sure that the right turn bladder hose quick connection coupler is together and locked
	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and set
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
		If you find a leak, replace the turn bladder assembly
No air entering the turn bladder	Make sure that the LED (D26) is lit to show the TLI solenoid is energizing (see Main power board assembly on page 1-28)	
	If energized, replace the TLI inflate solenoid	
No right turn deflate	Loose or disconnection of the right turn bladder hose quick connection coupler	Make sure that the right turn bladder hose quick connection coupler is together and locked
	No air is exiting the turn bladder	Make sure that the LED (D27) is lit to show the TLD solenoid is energizing (see Main power board assembly on page 1-28)
		If energized, replace the TLD inflate solenoid
No back right inflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
		If you find a leak, replace the turn bladder assembly
	No air entering the turn bladder	Make sure that the LED (D22) is lit to show the BRI solenoid is energizing (see Main power board assembly on page 1-28)
If energized, replace the BRI inflate solenoid		

Troubleshooting

EN

No back right deflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	No air is exiting the turn bladder	Make sure that the LED (D23) is lit to show the BRD solenoid is energizing (see Main power board assembly on page 1-28)
		If energized, replace the BRD inflate solenoid
No back left inflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
		If you find a leak, replace the turn bladder assembly
	No air entering the turn bladder	Make sure that the LED (D34) is lit to show the BLI solenoid is energizing (see Main power board assembly on page 1-28)
If energized, replace the BLI inflate solenoid		
No back left deflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	No air is exiting the turn bladder	Make sure that the LED (D36) is lit to show the BLD solenoid is energizing (see Main power board assembly on page 1-28)
		If energized, replace the BLD inflate solenoid
No seat right inflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
		If you find a leak, replace the turn bladder assembly
	No air entering the turn bladder	Make sure that the LED (D25) is lit to show the SRI solenoid is energizing (see Main power board assembly on page 1-28)
If energized, replace the SRI inflate solenoid		

Troubleshooting


EN

No seat right deflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	No air is exiting the turn bladder	Make sure that the LED (D24) is lit to show the SRD solenoid is energizing (see Main power board assembly on page 1-28) If energized, replace the SRD inflate solenoid
No seat left inflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	Possible leak in the turn bladder	Listen for air flowing into the turn bladder
		Look for a leak in the turn bladder
	No air entering the turn bladder	If you find a leak, replace the turn bladder assembly
Make sure that the LED (D33) is lit to show the SLI solenoid is energizing (see Main power board assembly on page 1-28) If energized, replace the SLI inflate solenoid		
No seat left deflate	Loose or disconnection of the pod manifold connector	Make sure that the pod manifold connector is set into the manifold and secured
	No air is exiting the turn bladder	Make sure that the LED (D32) is lit to show the SLD solenoid is energizing (see Main power board assembly on page 1-28) If energized, replace the SLD inflate solenoid
LAL fan will not run	LAL icon on the InTouch bed footboard is not orange	Tap the LAL icon to on
LAL fan is noisy	LAL fan is not secure in the frame or obstructed	If you find no obstruction and the fan is set in the frame, replace the LAL fan assembly
Foot box fan will not run	Loose or disconnection of the power supply quick connector	Check the power supply quick connector is together and locked
	Neither foot box fan is running	Use a voltmeter to check for -12VDC on pin 8 (+/ blue) and pin 6 (-/red) If no 12VDC, replace the power supply
Foot box fan is noisy	Fan cage is not secure or obstructed	If you find no obstruction and the fan is set in the cage, replace the foot box fan assembly

Isolibrium assembly


2971-021-001 Rev C (Reference only)

EN



Isolibrium assembly

EN


Isolibrium assembly

EN

Isolibrium assembly

EN

Isolibrium assembly


EN

Item	Part No.	Part Name	Qty
A	0026-672-000	Quick release pin	2
B	2971-021-002	Foam crib assembly	1
C	2971-021-003	Pods layer assembly	1
D	2971-021-004	Top cover assembly	1
E	2971-021-007	Turn bladder assembly	1
F	2971-021-012	Foot foam crib assembly	1
G	2971-021-017	Bottom cover assembly	1
H	2971-021-045	Pillow assembly	1
J	2971-021-189	Fan cable assembly	2
K	2971-021-901	Serial number label	1
L	2971-022-008	Foot box assembly	1
M	0058-383-000	Purse Lock Wire Tie	2

Footbox assembly


EN

2971-022-008 Rev D (Reference only)


Footbox assembly

EN


Footbox assembly

EN


Footbox assembly

EN



Footbox assembly

EN


Footbox assembly

EN


Footbox assembly

EN

Torque item B
to 22 ± 4 in-lb

Footbox assembly

EN

Item	Part No.	Part Name	Qty
A	0004-880-000	BHCS #8-32 x 3/8, Torx with star washer	4
B	0007-094-000	THMS, #8-32 x 3/4", 6 lobe T15	14
C	0004-883-000	BHCS #6-32 x 1/4 , 6 lobe	4
D	0013-018-000	Tooth lock washer	15
E	0015-091-000	1/2 NPT Heyco nylon locknut	1
F	0015-092-000	CAN outer nut	1
G	0015-093-000	Nut, 5/32	6
H	0015-094-000	Hex nut, #8-32, x-small	4
J	0016-033-000	Kep nut	1
K	0027-041-000	Spring detent cotter pin	2
L	0029-028-000	Push pin	4
M	0029-029-000	Wire clip	9
P	0036-046-000	Ground label	1
R	0038-330-000	Extension spring	2
T	0052-916-000	Hex standoff	4
U	0058-380-000	Push in coupler, 5/32	6
V	0058-381-000	Fan filter screen, 3.62 in.	1
W	0025-650-000	Dome head blind rivet	11
Y	0946-001-155	Bumper	5
AA	2011-001-215	Lug grounding	1
AB	2971-021-033	Sensor assembly	1
AC	2971-022-006	Valve manifold assembly	1
AD	2971-022-007	CPR puller assembly	2
AE	2971-022-009	Foot box bottom weldment	1
AF	2971-022-015	Foot box top cover assembly	1
AG	2971-022-012	Resonator assembly	1
AH	2971-022-013	Fan bracket assembly	2
AJ	2971-022-014	Mounted pump assembly	1
AK	2971-022-113	CPR activation cable	1
AL	2971-022-117	CPR switch cable assembly	2
AM	2971-022-129	Foot box power supply cable assembly	1
AN	2971-022-132	Main board to foot box SPI cable	1
AP	2971-022-133	CAN outside cable assembly	1
AR	2971-022-140	Main board PCB assembly	1
AT	2971-022-149	CAN inner cable assembly	1
AU	2971-022-150	Fan box cable assembly	1
AV	2971-022-151	Power cable assembly	1
AW	2971-022-165	Base CPR	2
AY	2971-022-171	Power supply cable assembly	1
BA	2971-022-188	Fan foot box cable assembly	2
BB	2971-022-192	Resonator foam	1
BC	2971-022-903	Color foot box label	1
BD	2971-022-904	Serial number foot box label	1
BE	2971-022-905	Pass foot box label	1
BF	3000-300-115	Standoff	10

Footbox assembly


EN

Item	Part No.	Part Name	Qty
BG	8815-029-200	Cable tie 7.4 in.	4
BJ	0029-003-000	Brass eyelet	2

Pump assembly kit 2971-700-001

EN

Rev B


Item	Part No.	Part Name	Qty
A	2971-022-014	Pump mounted assembly	1
B	8815-029-200	Cable tie 7.4 in.	2

Valve manifold assembly kit 2971-700-007

EN

Rev A


Item	Part No.	Part Name	Qty
A	2971-022-006	Valve manifold assembly	1
B	8815-029-200	Cable tie 7.4 in.	2

Recycling passports

Assembly part number: 2971-022-008 Rev D (Reference only)

EN


Item	Recycling/Material Code	Important Information	Qty
AB	2971-021 -033 Sensor Assembly	Printed circuit board	1

Recycling passports

Assembly part number: 2971-022-008 Rev D (Reference only)

EN


Item	Recycling/Material Code	Important Information	Qty
AJ	2971-022-014 Mounted Pump Assembly	Printed circuit board	1

Recycling passports

EN

Assembly part number: 2971-022-008 Rev D (Reference only)

Item	Recycling/Material Code	Important Information	Qty
AR	2971-022 -140 Main Board PCB Assembly	Printed circuit board	1
AY	2971-022-171 Power Supply Cable Assembly	Printed circuit board	1

EMC Information

EN

Guidance and Manufacturer's declaration - Electromagnetic Immunity			
<p>Isolibrium is suitable for use in the electromagnetic environment specified below. The customer or the user of Isolibrium should make sure that it is used in such an environment.</p>			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+6 kV contact +8 kV air	+6 kV contact +8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/burst IEC 61000-4-4 *	+2 kV for power supply lines +1 kV for input/ output lines	+2 kV for power supply lines +1 kV for input/ output lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5 *	+8 kV differential mode +2 kV common mode	+8 kV differential mode +2 kV common mode	Main power quality is that of typical commercial and/or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11 *	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Main power quality should be that of a typical commercial and/or hospital environment. If the user of Isolibrium requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.
<p>Note: U_T is the a.c. mains voltage before applications of the test level.</p>			

EMC Information

EN

Recommended separation distances between portable and mobile RF communications equipment and Isolibrium.

Isolibrium is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of **Isolibrium** can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and **Isolibrium** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=[3.5/\sqrt{f}] \sqrt{P}$	80 MHz to 800 MHz $d=[3.5/E1] \sqrt{P}$	800 MHz to 2.5 GHz $d=[7/E1] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

EMC Information

EN

Isolibrium is suited for use in the electromagnetic environment specified below. The customer or the user of **Isolibrium** should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6 *</p> <p>Radiated RF IEC 61000-4-3</p>	<p style="text-align: center;">3 Vrms 150 kHz to 80 MHz</p> <p style="text-align: center;">3 V/m 80 MHz to 2.5 GHz</p>	<p style="text-align: center;">3 Vrms</p> <p style="text-align: center;">3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of Isolibrium, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p style="text-align: center;">Recommended Separation Distance</p> <p style="text-align: center;">$d=1.2\sqrt{P}$</p> <p style="text-align: center;">$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p style="text-align: center;">$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **Isolibrium** is used exceeds the applicable RF compliance level above, **Isolibrium** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **Isolibrium**.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

EMC Information

EN

(Continued)

Guidance and Manufacturer's declaration - Electromagnetic Emissions		
Isolibrium is intended for use in an electromagnetic environment specified below. The customer or the user of Isolibrium should make sure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	Isolibrium uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	Isolibrium is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

 **WARNING**

This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating **Isolibrium** or shielding the location.

Warranty

EN

Stryker Medical, a division of Stryker Corporation (“Stryker”), warrants that its Model 2971 **Isolibrium™** product will be free from defects in material and workmanship. This Stryker warranty covers only the following items of the Stryker **Isolibrium** product (each known individually as a “Part” and collectively as the “Product” or “**Isolibrium** product”) during normal use* as follows:

- Support Surface (Air Pods, Turning Bladders and Foam Crib) Warranty Period: 2 years
- Cover and Fire Barrier Assembly Warranty Period: 2 years
- Pump Warranty Period: 1 year

The above noted Warranty periods apply only to the original purchaser of the **Isolibrium** product and begin on the date of delivery to such original purchaser.

If Stryker determines, in its sole discretion, that one or more parts is defective within the above noted warranty periods, then Stryker may, at its option, either repair or replace the **Isolibrium** product or part.

The warranty set forth above does not include or cover the following:

- Abnormal wear and tear on the product, or wear which indicates that the product was not properly maintained in accordance with this Operations/Maintenance Manual, or which product has been subject to unusual stress; or
- Product that has been misused, modified, refurbished or repaired without the prior written consent of Stryker; or
- Damage or product failure due to causes beyond Stryker’s control, including but not limited to, abuse, theft, fire, flood, wind, lightning, freezing, clogging of support surface pores due to tobacco smoke, unusual atmosphere conditions, or material degradation due to exposure to moisture; or
- Damage which is determined to have resulted through the use of the product for patient transfer or transport; or
- Product which serial numbers or other identification marks have been removed or destroyed.

*“Normal use” is defined as use of the product in typical or normal use settings in a hospital or medical facility under normal conditions. Damage to the product which arises from abnormal use, such as damage to the product that may be caused by needle punctures, burns, chemicals, negligent use, improper care, or improper cleaning (as such is detailed in this Operations/Maintenance Manual) or staining resulting from such abnormal uses are exempt from the above note warranty coverage.

Warranty exclusion and damage limitations

The express warranty set forth herein is the only warranty applicable to the product. **Any and all other warranties, whether express or implied, including any implied warranty of merchantability or fitness for a particular purpose are expressly excluded by Stryker.** In no event shall Stryker be liable for incidental or consequential damages.

To obtain parts and service

Stryker products are supported by a nationwide network of dedicated Stryker Field Service Representatives. These representatives are factory trained, available locally, and carry a substantial spare parts inventory to minimize repair time. Simply call your local representative or call Stryker Customer Service at 1-800-327-0770.

Return authorization

Product cannot be returned without prior approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Warranty

EN

Damaged product

ICC Regulations require that claims for damaged product must be made with within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the product, or the damage was not noted on the delivery receipt at the time of receipt, the customer will be responsible for payment of the original invoice in full within thirty (30) days of receipt. Claims for any incomplete shipments must be made within thirty (30) days of invoice.

International warranty clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Contact your local Stryker Medical representative for additional information.

Patent information

This product is made with Intelli-Gel®* hollow column configuration and DurageI™ elastomeric material.

*Intelli-Gel® is a registered trademark of EdiZONE, LLC of Alpine, UT

DurageI™ is a trademark of EdiZONE, LLC of Alpine, UT

Covered by one or more of the following patent numbers:

United States	5,749,111	6,026,527	7,076,822	7,964,664
---------------	-----------	-----------	-----------	-----------

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*Intelli-Gel® est une marque déposée d'EdiZONE, LLC d'Alpine, Utah, États-Unis
DuraGel™ est une marque commerciale d'EdiZONE, LLC d'Alpine, Utah, États-Unis