Isolibrium® PE support surface

Operations Manual

REF 29730000000

2797

EN
CS
DA
DE

EL
ES
ET
FI
FR
HR
HU
IT

JA
KO
LT
LV
NL
NO
PL
PT

RO
RU
SK
SL
SV
TH
TR
ZH

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Patient stability and siderails coverage may be compromised with the use of an overlay
	Warning; electricity
	Indicates that this product does not contain toxic and hazardous substances or elements above the maximum concentration of all 6 values defined by the China RoHS legislation. This product is an environmentally friendly product which can be recycled and reused.
	Catalogue number
	Serial number
	European medical device
	CE mark
	Authorized representative in the European Community
	For US Patents see www.stryker.com/patents
	Manufacturer
	Safe working load
	Mass of equipment
	Protection from liquid splash
	Class II electrical equipment: equipment in which protection against electric shock does not rely on basic insulation only, but in which additional safety precautions such as double insulation or reinforced insulation are provided, there being no provision for protective earthing or reliance upon installation conditions.
	Defibrillation proof type B applied part


	<p>Medical Equipment recognized by UL LLC with respect to electric shock, fire, and mechanical hazards in accordance with ANSI/AAMI ES60601-1:2005/(R)2012, A1:2012, C1:2009/(R)2012, A2:2010/(R)2012, and CAN/CSA-C22.2 No. 60601-1 (2008)</p>
	<p>In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.</p>
	<p>Wash by hand</p>
	<p>Do not tumble dry</p>
	<p>Do not dry-clean</p>
	<p>Do not iron</p>
	<p>Chlorinated bleach</p>
	<p>Keep dry</p>
	<p>Stacking limit by number</p>
	<p>This side up</p>
	<p>Fragile</p>
	<p>Do not use sharp objects to open the package</p>
	<p>Center of gravity/weight distribution</p>
	<p>Two person lift</p>

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	4
Product description	4
Indications for use	4
Clinical benefits	4
Contraindications	4
Expected service life	5
Disposal/recycle	5
Specifications	5
European REACH	6
Product illustration	6
Contact information	7
Serial number	7
Isolibrium PE features	8
Setup	9
Operation	12
Applying the linens	12
Positioning a patient on the support surface	12
Managing incontinence and drainage	13
Activating CPR	13
Resetting CPR	15
Relocating the support surface	16
Transferring of patient from one patient support platform to another	17
Storing the product	17
Accessories and parts	18
Caring for the bottom cover	19
Cleaning and disinfecting with wipes	21
Cleaning	22
Disinfecting	23
Preventive maintenance	24
Cover replacement, top	25
EMC information	27

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Always confirm the scale system functions on the **ProCuity** bed after you setup the support surface to prevent product malfunction. Always make sure that the bed is zeroed/tared with the support surface on the frame for the appropriate patient weight for therapy.
 - Do not use fitted sheets with this support surface.
 - Always follow your hospital CPR protocols. Make sure that the operator has access to the CPR straps.
 - Do not exceed the safe working load of the product.
 - Always monitor the patient condition at regular intervals for patient safety.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Only use **Isolibrium** PE with the compatible **ProCuity** bed frame. Use of this device on an incompatible frame will cause improper or no function of the support surface. Failure to use a compatible frame may void warranty if damages incur.
- Risk of Shock. This product must only be connected to the appropriate terminal on the **ProCuity** bed frame.
- Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
- Do not move the support surface without assistance.
- Always use the bottom cover retainer straps to secure the support surface to the **ProCuity** bed frame.
- Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
- Always center the patient on the support surface. Align the patient's head toward the headboard before you start functions. Check the patient to maintain the proper position.
- Always make sure that the tubes and wires connected to the patient are long enough, stable, and secure while lateral rotation or turn assist is active.
- Always raise the bed siderails before you start turn assist or lateral rotation functions.
- Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillowcase or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
- Always confirm that all persons and equipment are away from the area below and around the Fowler before you activate the CPR release handle. The CPR release handle is for emergency use only.

- Do not use the support surface as a transfer device to avoid the risk of patient injury.
 - Always unplug the support surface power/data cable before you clean or disinfect.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - Always perform preventive maintenance based on the use of the product. An increase in the use of the product may include more frequent cleaning and disinfection and may adversely affect the life of the support surface.
 - Always wipe each product with clean water and dry after you clean or disinfect. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not iron, dry-clean, tumble dry, or power wash the support surface or support surface covers.
 - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
 - Do not overexpose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers greater than 3% glycol ether as they may damage the cover or the CPR strap print.
 - Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the product. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they operate properly.
 - The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The Stryker Model 297300000000 **Isolibrium** PE is a powered support surface with features that provide pressure redistribution, low air loss, turn assist, max inflate, and lateral rotation. The **Isolibrium** PE powered support surface is for use with the Stryker Model 300900000000 **ProCuity**™ bed frame.

Indications for use

The Stryker Model 297300000000 **Isolibrium** PE support surface is intended to assist in the prevention and treatment of all pressure injury or pressure ulcers (includes all stages, unstageable and deep tissue pressure injury). As a recommendation, healthcare professional should complete a clinical evaluation of risk factors and skin assessments. The intended users is both healthcare HCPs (nurses, nurse aids, and medical doctors), transporters, and human patients. This product is for use by human patients in a health delivery organization (HDO).

The **Isolibrium** PE patient-specific immersion is intended to assist in the management of microclimate of the patient's skin, patient repositioning, and the support in prevention of early mobility and pulmonary complications.

The **Isolibrium** PE powered support surface is not intended for use in or with:

- Oxygen rich environments
- Sterile environments
- Home care or long-term care facility settings
- Behavioral health populations

Clinical benefits

Assists in the prevention and treatment of all pressure ulcers or pressure injuries

Contraindications

Isolibrium PE is contraindicated for use on patients with:

- Unstable fractures
- Unstable spinal cord injuries

Expected service life

Isolibrium PE support surface has a five year expected service life under normal use, conditions, and with appropriate periodic maintenance.

Isolibrium PE covers (top and bottom) have a two year expected service life under normal use, and conditions.







Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

 <p>Note - Do not exceed the safe working load of the product.</p>	460 lb	208.6 kg
Minimum patient weight	60 lb	27.2 kg
Length	84 in.	213.4 cm
Width	35 in.	88.9 cm
Thickness	8 - 10.5 in.	20.3 - 26.7 cm
Weight	< 90 lb	< 40.8 kg
Electrical safety	IEC 60601-1-2 Ed. 4.0b:2014, BS EN 60601-1-2:2015, IEC 60601-1-2 Ed. 3.0 b:2007. CAN/CSA C22.2 NO. 60601-1-2:16, AIM 7351731	
DC Voltage	36V +/- 10%	
Current	3.35A	
Maximum noise level	50 dBA	
Flammability	16CFR 1632, 16CFR 1633, California Technical Bulletin 129, BFD IX-11, BS 6807:2006 Clause 9, Method 27.7-1979 of CAN 2-4.2 M77, BS EN 597-1:2015, BS EN 597-2:2015, Italy UNI 9175:2008 Clause 9	
Compatible bed frame	Model 300900000000 ProCuity bed frame	

Stryker reserves the right to change specifications without notice.


Environmental conditions	Operation	Storage and transportation
Ambient temperature		
Relative humidity (non-condensing)		
Atmospheric pressure		

European REACH

In accordance with the European REACH regulation and other environmental regulatory requirements, the components that contain declarable substances are listed.

Description	Number	Substance of very high concern (SVHC) chemical name
Sensor assembly	2971-021-033	1,3,5-Tris(oxiran-2-ylmethyl)-1,3,5-triazinane-2,4,6-trione (TGIC)
		2-Ethoxyethyl acetate
		Cadmium
		Diboron-trioxide
		Lead-monoxide
Main board PCB assembly	2971-022-140	Cadmium
		Diboron-trioxide
		Lead-monoxide
Footbox bottom weldment	2971-022-009	Tris(2,3-epoxypropyl) isocyanurate
Footbox top	2971-022-200	Tris(2,3-epoxypropyl) isocyanurate
Pump assembly	297300220010	Lead-steel
		Lead-brass
		Lead-aluminum
Power supply	297100220910	Lead (Pb) Cas#7439-92-1

Product illustration

A	Head end	F	Top support surface cover
B	Foot end	G	Bottom support cover
C	D-ring	H	Retainer strap (D-ring not shown)
D	CPR strap	I	Four handles on the support surface (not shown)
E	Power/data cable		

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA


Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number

The serial number (A) for your support surface is located on the patient right side of the foot box on the specification label. Unzip the support surface top cover to access the specification label. The serial number (B) is for the support surface cover.

Isolibrium PE features

Pressure redistribution	The ability of a support surface to distribute load over the contact areas of the human body.
Pressure redistribution setpoint	Allows the operator to input patient weight range that results in automatic change of zoned sensor technology for immersion and envelopment to occur.
Lateral Rotation Therapy	A feature of a support surface that provides rotation about a longitudinal axis as characterized by degree of patient turn, duration, and frequency.
BackSmart® Turn Assist	A feature of a support surface that provides a single position turn about a longitudinal axis.
Max Inflate	Allows maximum firmness inflation for use when you boost a patient and on lateral transfer. Promotes patient ease with egress and ingress.
Moisture management	The top surface cover draws excess moisture away from the patients skin and support surface interface.
Low Air Loss (LAL)	Provides air flow to assist in the management of (microclimate) the skin.
HeelGel®¹ Technology	A feature of the support surface to address pressure redistribution in the heel area.
Patient Specific Immersion™	Assists in the management of microclimate of the patient skin, patient repositioning, early mobility, and prevention of pulmonary complications.

¹**Intelli-Gel®** hollow column configuration licensed by EdiZONE LLC of Alpine, UT

Note - See the **ProCuity** bed Operations Manual for operating instructions related to **Isolibrium** PE features.

Setup

The product must reach the ambient operating temperature before setup or operation to avoid the risk of permanent damage to this product.

WARNING

- Always confirm the scale system functions on the **ProCuity** bed after you setup the support surface to prevent product malfunction. Always make sure that the bed is zeroed/tared with the support surface on the frame for the appropriate patient weight for therapy.
 - Do not use fitted sheets with this support surface.
-

CAUTION

- Only use **Isolibrium** PE with the compatible **ProCuity** bed frame. Use of this device on an incompatible frame will cause improper or no function of the support surface. Failure to use a compatible frame may void warranty if damages incur.
 - Risk of Shock. This product must only be connected to the appropriate terminal on the **ProCuity** bed frame.
 - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal, mask type, or half bed-length tent type.
 - Do not move the support surface without assistance.
 - Always use the bottom cover retainer straps to secure the support surface to the **ProCuity** bed frame.
-

Note - Do not operate the support surface without the top cover, bottom cover, and linens.

Prerequisite: Minimum of two operators are required to setup the support surface.

To setup the support surface:

1. Remove all shipping and packaging materials from the product before use.
2. Allow the product to reach room temperature.
3. Apply the brakes on the **ProCuity** bed frame.
4. Lower the siderails to the lowest position of the **ProCuity** bed frame.
5. Put the bed into the flat position.
6. Pull the bed extender out on the **ProCuity** bed frame.
7. If powered, unplug the bed from the wall outlet. Turn the battery switch to off.
8. Remove the headboard and footboard from the bed and set aside.
9. See *Resetting CPR* (page 15).
10. Position an operator at each end of the support surface and lift onto the bed frame.
11. Insert the power/data cable (A) at the foot end of the support surface through the strain relief (B) (Figure 1).


Figure 1 – Secure connection

12. Connect the power/data cable (A) to the **ProCuity** cable adapter (C) (Figure 1).
13. Align the locking collar of the cable adapter (C) (Figure 1). Push to seat the connector. The locking collar will rotate into place when the connector is fully seated.
14. Reinstall the headboard and footboard removed in step 7.
15. Turn the battery switch to on. Plug the bed power cord into a hospital grade protective earthed outlet.
16. Raise the Fowler to access the Fowler retainer holes (D) (Figure 2).


Figure 2 – Retainer holes

17. Thread the retainer strap (E) down through the Fowler retainer hole (D) both sides (Figure 3).


Figure 3 – Insert retainer straps

18. Put the retainer strap (E) through the D-rings (F), (Figure 4).


Figure 4 – Retainer strap

19. Loop the retainer strap between the D-rings (G) (Figure 4). Pull the retainer strap until secure (Figure 5).


Figure 5 – Pull retainer strap

20. Insert the loose end of the retainer strap behind the Fowler.

21. Repeat steps 15 to 18 for the other retainer strap.

22. Place the linens on the support surface as loose as possible. This is to avoid a hammock effect that may increase interface pressures.

Operation

Applying the linens

WARNING

- Always follow your hospital CPR protocols. Make sure that the operator has access to the CPR straps.
 - Do not use fitted sheets with this support surface.
-

CAUTION - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.

To apply the linens:

1. Use the D-rings to apply the flat sheet linens.
2. Thread the four linen corners through the D-rings (C) located on the bottom cover. Secure the linens (A) to the support surface (B).


Figure 6 – Apply linens

Note

- Do not pull the linens tight. Keep the linens loose and as smooth as possible on top of the support surface. This is to effectively use the Turn Assist and Lateral Rotation functions.
- Always make sure that the CPR straps are accessible.

Positioning a patient on the support surface

WARNING

- Do not exceed the safe working load of the product.
-

CAUTION

- Always center the patient on the support surface. Align the patient's head toward the headboard before you start functions. Check the patient to maintain the proper position.
 - Always make sure that the tubes and wires connected to the patient are long enough, stable, and secure while lateral rotation or turn assist is active.
 - Always raise the bed siderails before you start turn assist or lateral rotation functions.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
 - Do not allow sharp edges from the X-ray plate to come in contact with the support surface cover. The recommendation is you cover the X-ray plate with a pillowcase or other device before placement under the patient. If damaged, remove the support surface cover from service immediately to prevent cross contamination.
-

To position the patient:

1. Max Inflate the support surface.
2. Position the patient in the center of the support surface.
3. Align the patient's head toward the headboard (Figure 7).


Figure 7 – Center the patient

4. Check the patient during lateral rotation for proper positioning and support surface inflation (Figure 8).

Figure 8 – Correct patient positioning

Managing incontinence and drainage

WARNING - Always monitor the patient condition at regular intervals for patient safety.

You can use disposable diapers or incontinence pads to manage incontinence. Always provide appropriate skin care after each incontinence episode.

Activating CPR

WARNING - Always follow your hospital CPR protocols. Make sure that the operator has access to the CPR straps.

CAUTION - Always confirm that all persons and equipment are away from the area below and around the Fowler before you activate the CPR release handle. The CPR release handle is for emergency use only.

To prepare the bed frame and support surface for CPR, these steps must be followed:

1. Pull and hold the **ProCuity** CPR release handle (A) to lower the head of bed to the flat position. The two CPR release handles are on the left and right side of the litter Gatch section of the bed (Figure 9).

Figure 9 – Activating the CPR release handle

Note - If **ProCuity** is connected to AC power, the support surface will automatically deflate. If the bed is not connected to AC power, the support surface will not automatically deflate. The support surface will remain in the prior therapy setting until the CPR strap is pulled, go to step 2.

2. A message will appear on the ProCuity bed graphical user interface (Figure 10). Pull the strap on the support surface to deflate faster.

Figure 10 – CPR Activated - Surface Deflating

3. Raise or lower the **ProCuity** foot siderail to the full up or full down position.
4. Pull one of the red CPR support surface straps straight out to manually deflate the surface. The CPR straps are located on the left and right side of the foot end of the support surface. The CPR strap pulls outward about one inch and locks in place to release any air that remains in the mattress (Figure 11).

Figure 11 – CPR strap pulled

Resetting CPR

To reset the **ProCuity** bed CPR status and reinflate the support surface:

1. Pull straight down on the red CPR support surface strap to reset the surface. In a slow motion, release the strap to its original position. The support surface will reinflate when the CPR strap is reset. The status on the **ProCuity** CPR activated display will disappear (Figure 12).

Figure 12 – Display when CPR strap pulled

2. If the support surface CPR strap was not pulled for CPR, the surface needs to be reset. The **ProCuity** CPR Activated - Surface Deflating display (Figure 13).

Figure 13 – Display when CPR strap not pulled

Note - ProCuity display function and therapy control is prevented until the support surface has been reinflated and the display notification resolved. Reset minimizes the occurrence of a patient exposure to a deflated surface.

Relocating the support surface

Relocate the support surface to clean or move to another **ProCuity** bed frame.

CAUTION

- Do not move the support surface without assistance.
- Do not use the support surface as a transfer device to avoid the risk of patient injury.

Prerequisite: Minimum of two operators are required to relocate and setup this product.

To relocate the support surface:

1. Apply the brakes on the **ProCuity** bed frame.
2. Lower the siderails to the lowest position of the **ProCuity** bed frame.
3. Put the bed into the flat position.
4. Remove the patient from the support surface. Follow your hospital protocols.
5. Pull the bed extender out on the **ProCuity** bed frame.
6. Remove the linens.
7. Raise the Fowler to access the Fowler retainer holes.
8. Release the support surface retainer straps from the bed.
9. Lower the Fowler to the flat position.
10. If powered, unplug the bed from the wall outlet. Turn the battery switch to off.
11. Unplug the support surface power/data cable from the bed.
12. Remove the power/data cable from the strain relief bracket.
13. Remove the headboard and footboard from the bed and set aside.
14. Position an operator at each end of the support surface and remove the support surface and relocate.

Note - Use all four handles when you lift to relocate the support surface.

15. Reinstall the headboard and footboard removed in step 9.
16. Turn the battery switch to on. Plug the bed power cord into a hospital grade protective earthed outlet.
17. See *Setup* (page 9) to complete the transfer of the support surface.

Transferring of patient from one patient support platform to another

Do not exceed the safe working load of the product.

CAUTION

- Do not use the support surface as a transfer device to avoid the risk of patient injury.
 - Do not stick needles into a support surface through the support surface cover. Holes may allow body fluids to enter the inside (inner core) of the support surface and could cause cross-contamination, product damage, or product malfunction.
-

Prerequisite: Follow hospital protocols required to transfer a patient from one surface to another.

To transfer the patient from one patient support surface to another:

1. Position one patient support platform next to the other patient support platform. Minimize the space between the two platforms.
2. Apply the brakes for both patient support platforms.
3. Inflate the support surface with max inflate, recommended. See the **ProCuity** Operations Manual to start and stop max inflate.
4. Adjust the patient support platform heights so that they are level with one another.
5. Transfer the patient. Follow hospital protocols for applicable safety rules for patient and operator.
6. Stop max inflate.

Storing the product

Before you store the product, deflate the support surface by activating CPR.

Accessories and parts

These parts are currently available for purchase. Call Stryker Customer Service: 1-800-327-0770 for availability and pricing.

Part name	Part number
Board, main	2971-022-140
Board, sensor	2971-021-033
Cover, bottom - Endurance	297300210017
Cover, top - Endurance	2972-021-004
Fan, foot box	2971-022-188
Fan, LAL	2971-021-189
Power/data cable	297100560801
Power supply PCBA	297100220910
Kit, pump assembly	297307000001
Kit, resonator assembly	2971-700-002
Kit, valve manifold	2971-700-007
Pods assembly	297300210003
Kit, solenoid valve, SLI, BLI, BLD, BRI, BRD, SRD, SRI, SLD	2971-700-004
Kit, solenoid valve, TRI, TRD, EX, TRI, TU	2971-700-005
Kit, solenoid valve, IN	2971-700-006

Caring for the bottom cover

To care for the bottom cover:

1. Apply the brakes on the **ProCuity** bed frame.
2. Unplug the support surface power/data cable from the bed.
3. Remove the power/data cable from the strain relief bracket.
4. Pull one of the CPR straps to deflate the support surface.
5. Lower the siderails on the bed.
6. Extend the **ProCuity** bed extender (Figure 14).

Figure 14 – Bed extended

7. Remove the headboard and footboard from the bed and set aside.
8. Raise the head of bed angle.
9. Place the bed in reverse Trendelenburg (Figure 15).

Figure 15 – Reverse Trendelenburg

10. Remove the retaining straps from the bed.
11. With a helper, use the handles on the support surface to help fold the support surface to the foot end of the litter.
12. Care for the underside of the support surface and litter deck (Figure 16).

Figure 16 – Care of head end

13. Dry the support surface and litter deck.
14. Unfold the support surface and place the bed in Trendelenburg.
15. Reverse steps 9-11 for the care of the foot end of the support surface and litter deck (Figure 17).

Figure 17 – Care of foot end

16. Unfold the support surface.
17. Reattach the retaining straps. See *Setup* (page 9).
18. Replace the headboard and footboard removed in step 6.
19. Plug the support surface power/data cable in through the strain relief.
20. See *Resetting CPR* (page 15).

Cleaning and disinfecting with wipes

For United States only. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Stryker's preferred wipes (2060-000-001 6" x 10" or 2060-000-002 9" x 12") include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

To clean or disinfect the external product surface:

1. To clean, wipe external surfaces with a fresh, clean wipe to remove all visible soils. Repeat as necessary until the product is clean.

Note

- Use as many wipes as necessary.
 - Complete step 1 before you disinfect.
2. To disinfect, wipe external surfaces with a fresh, clean wipe until wet. Allow the external surface to remain wet for two minutes at room temperature.
 3. Allow the product to dry before you return it to service.

Cleaning

CAUTION

- Always unplug the support surface power/data cable before you clean or disinfect.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always inspect the support surface covers (top and bottom) for tears, punctures, excessive wear, and misaligned zippers as you clean the covers. If the cover is compromised, remove the support surface from service immediately and replace the cover to prevent cross-contamination.
 - Always perform preventive maintenance based on the use of the product. An increase in the use of the product may include more frequent cleaning and disinfection and may adversely affect the life of the support surface.
 - Always wipe each product with clean water and dry after you clean or disinfect. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not iron, dry-clean, tumble dry, or power wash the support surface or support surface covers.
-

Prerequisite: A minimum of two operators are required to clean the bottom support surface cover.

Always follow hospital protocol for cleaning and disinfecting.

To clean the support surface covers between patient use, follow these steps in order:

1. If powered, unplug the bed from the wall outlet. Turn the battery switch to off.
2. Unplug the support surface power/data cable from the bed.
3. Remove the power/data cable from the strain relief bracket.
4. Using a clean, soft, damp cloth, wipe the support surface covers with a mild soap and water solution to remove foreign material.
5. Wipe the support surface covers with a clean, dry cloth to remove any excess liquid or cleaning agent.
6. Rinse and dry the support surface covers.
7. Turn the bed battery switch to on.
8. Plug the bed power cord into a hospital grade protective earthed outlet.
9. Plug the support surface power/data cable in through the strain relief.

Note - Dry the support surface covers before you store, add linens, or place a patient on the surface.

Disinfecting

CAUTION

- Always unplug the support surface power/data cable before you clean or disinfect.
 - Do not immerse the support surface or foot box in cleaning or disinfectant solutions.
 - Always disinfect the support surface between patients, to avoid the risk of cross-contamination and infection.
 - Do not allow fluid to pool on the support surface. Fluids may cause degradation of components and may cause unpredictable safety and performance of this product.
 - Always wipe each product with clean water and dry after you clean or disinfect. Some cleaning and disinfecting agents are corrosive in nature and may cause damage to the product. Failure to properly rinse and dry the product leaves a corrosive residue on the surface of the product and may cause premature degradation of critical components. Failure to follow these cleaning or disinfecting instructions may void your warranty.
 - Do not overexpose the covers to higher concentration disinfectant solutions as these may degrade the covers.
 - Do not use accelerated hydrogen peroxides or quaternaries that contain glycol ethers greater than 3% glycol ether as they may damage the cover or the CPR strap print.
-

Prerequisite: Minimum of two operators are required to disinfect the support surface.

Suggested disinfectants:

- Quaternaries (active ingredient - ammonium chloride) that contain less than 3% glycol ether
- Phenolic solution (Matar)
- Chlorinated bleach solution for the Endurance cover is 10,000ppm
- 70% Isopropyl Alcohol

To disinfect the support surface covers after each patient use, follow these steps in order:

1. Unplug the bed from the wall outlet. Turn the battery switch to off.
2. Unplug the support surface power/data cable from the bed.
3. Remove the power/data cable from the strain relief bracket.
4. Clean and dry the support surface covers (see *Cleaning* (page 22)) before disinfectants are applied.
5. Apply recommended disinfectant solution with a damp cloth or pre-soaked wipes. Do not soak the support surface.

Note

- Follow the disinfectant manufacturer's instructions for appropriate contact time and rinse requirements.
 - Raise the head of bed and foot end to help fold the head end and foot end of the support surface to disinfect the bottom cover of the support surface.
6. Rinse the support surface covers with a clean, damp cloth to remove excess disinfectant.
 7. Wipe the support surface covers with a clean, dry cloth to remove excess liquid or disinfectant.
 8. Allow the support surface covers to dry before you return the product to service.
 9. Turn the bed battery switch to on.
 10. Plug the bed power cord into a hospital grade protective earthed outlet.
 11. Plug the support surface power/data cable in through the strain relief.

Note - Dry the support surface covers before you store, add linens, or place a patient on the surface.

Preventive maintenance

Remove product from service before you perform preventive maintenance inspection. Check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage. Service only by qualified personnel.

Note

- Clean and disinfect the exterior of the support surface before inspection, if applicable.
- If excessive wear is observed on the cover or the fire barrier, the recommendation is to replace the cover.

Inspect the following items:

- _____ All fasteners are secure
- _____ Zipper, snaps, and covers (top, bottom, and fire barrier) are free of tears, cuts, holes, or any other damage
- _____ Support surface cover labels are legible, adhere, and free of damage
- _____ Handles are free of rips or cracks
- _____ Power/data cable is free of damage
- _____ **Perform a functional test of the **Isolibrium** PE functions
- _____ **Perform system diagnostics and confirm there are no errors
- _____ *Internal components for signs of stains from fluid ingress or contamination by fully unzipping the cover
- _____ Pods do not leak or are not cracked
- _____ Hose connections to the manifold and sensor tubes (both ends) are seated
- _____ Foam is free from large tears or large gouges
- _____ Turn bladders are free of excessive wear and function
- _____ Low Air Loss fans rotate freely and are free of debris and dust
- _____ Left and right CPR releases function
- _____ Both foot box cooling fans rotate and are free of debris and dust
- _____ Pump spring mount has visible clearance from the foot box bottom
- _____ Listen for abnormal pump vibration and noise level
- _____ Clean the foot box intake metal screen
- _____ Clean the mesh on the bottom cover
- _____ Clean the fan filter screen on the foot box
- _____ Standoffs (feet of the foot box) are present
- _____ Cable tie is present on the power/data cable

Note

- *A zipper pull tool or equivalent is required to access the internal components of the support surface.
- **You cannot perform these preventive maintenance checks without the footboard, see the **ProCuity** Maintenance Manual.
- The zipper overlaps at the end and may appear to be misaligned per design.

Replace the following items every two years:

- Top cover assembly
- Bottom cover assembly

Product serial number:
Completed by:
Date:

Cover replacement, top

Tools required:

- Zipper pull tool or equivalent

Procedure:

1. Apply the brakes on the **ProCuity** bed frame.
2. Unplug the support surface power/data cable from the bed.
3. Remove the power/data cable from the strain relief bracket.
4. Remove the headboard and footboards.
5. Lower all the siderails.
6. Unsnap the two secure snaps (C) at the foot end corners (Figure 18).
7. Unsnap the four corner retainers (B) (Figure 18).

Note

- On install, reconnect the corner retainers. Place the strap below each D-ring at the foot end corners and snap the cover.
- The zipper starts on the patient right side near the foot end.


Figure 18 – Top cover

8. Insert a zipper pull tool or equivalent through the hole in the zipper slider to unlock.
9. Use the zipper pull tool or equivalent, to unzip the support surface top cover (A) from the bottom (Figure 18).
10. Remove and discard the top cover.
11. Reverse to install the new cover.

Note

- After installation, remove the zipper pull tool or equivalent from the zipper.
- Cover the zipper with the support surface cover watershed.

12. Verify proper operation before you return the product to service.

EMC information

CAUTION

- Avoid stacking or placing equipment adjacent with other equipment to prevent improper operation of the product. If such use is necessary, carefully observe stacked or adjacent equipment to make sure that they operate properly.
- The use of accessories, transducers, and cables, other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

Note - Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be no closer than 12 inches (30 cm) to any part of **Isolibrium** PE support surface, including cables specified by the manufacturer.

The Model 297300000000 **Isolibrium** PE support surface was evaluated using the following cables:

Cable	Length (m)
Isolibrium to bed	1.0

Guidance and manufacturer's declaration - electromagnetic emissions

The Model 297300000000 **Isolibrium** PE support surface is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 297300000000 **Isolibrium** PE support surface should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment
RF Emissions CISPR 11	Group 1	Note - The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
RF Emissions CISPR 11	Class A	
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations Flicker Emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity

The Model 297300000000 **Isolibrium** PE support surface is suitable for use in a professional healthcare facility environment and not in environments exceeding immunity test conditions that the product was evaluated to, such as near high frequency (HF) surgical equipment and inside of the radio frequency (RF) shielded room of magnetic resonance imaging (MRI) equipment. The customer or the user of the Model 297300000000 **Isolibrium** PE support surface should assure that it is used in such an environment and that the electromagnetic environment guidance listed below is followed.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Main power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV lines to lines ± 0.5 kV, ± 1 kV, ± 2 kV lines to earth	± 0.5 kV, ± 1 kV lines to lines ± 0.5 kV, ± 1 kV, ± 2 kV lines to earth	Main power quality should be that of a typical commercial or hospital environment.
Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11	<p>0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0% U_T for 1 cycle</p> <p>70% U_T (30% dip in U_T) for 25/30 cycles</p> <p>0% U_T for 250/300 cycles</p>	<p>0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°</p> <p>0% U_T for 1 cycle</p> <p>70% U_T (30% dip in U_T) for 25/30 cycles</p> <p>0% U_T for 250/300 cycles</p>	Main power quality should be that of a typical commercial or hospital environment. If the user of the Model 297300000000 Isolibrium PE support surface requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note - U_T is the a.c. mains voltage before applications of the test level.			

<p>Conducted RF IEC 61000- 4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz</p>	<p>3 Vrms 3 V/m</p>	<p>Portable and mobile RF communications equipment should follow the guidance in the table titled “Recommended separation distances between portable and mobile RF communication equipment and the Model 297300000000 Isolibrium PE support surface.” If the mobile service is not listed in the table, the recommended separation distance should be calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended separation distance $D=(2)(\sqrt{P})$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
---	---	-------------------------	--

Note

- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- The ISM (Industrial, Scientific, and Medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 297300000000 **Isolibrium** PE support surface is used exceeds the applicable RF compliance level above, the Model 297300000000 **Isolibrium** PE support surface should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Model 297300000000 **Isolibrium** PE support surface.

^bOver the frequency range 150 kHz to 80 MHz, field strengths are less than 3 Vrms.

Recommended separation distances between portable and mobile RF communication equipment and the Model 297300000000 Isolibrium PE support surface

The Model 297300000000 **Isolibrium** PE support surface is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model 297300000000 **Isolibrium** PE support surface can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 297300000000 **Isolibrium** PE support surface, including cables, as recommended below, according to the maximum output power of the communications equipment.

Band (MHz)	Service	Maximum power (W)	Minimum separation distance (m)
380-390	TETRA 400	1.8	0.3
430-470	GMRS 460; FRS 460	2.0	0.3
704-787	LTE Band 13, 17	0.2	0.3
800-960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	2.0	0.3
1,700-1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	2.0	0.3
2,400-2,570	Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE Band 7	2.0	0.3
5,100-5,800	WLAN 802.11 a/n	0.2	0.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


Note - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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