stryker

RapR-Round®

REF 8001-061-530

REF 8001-061-535

REF 8001-061-540

REF 8001-061-550

REF 8001-061-560

Operations Manual

Manual de uso

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Manuale d'uso

Bedieningshandleiding

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Benjningsmanual

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Brükerhandbok

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8001-009-111 REV E

www.stryker.com

This manual assists you with the operation of the Stryker Model 8001-061-530, 8001-061-535, and 8001-061-540 RapR-Round® vest leg wrap. Read this manual before operating this product and keep a copy on file. Set methods and procedures to educate and train your staff on the safe operation of this product.

Intended use and Contraindications

Always consult the controller operations manual before using this product.

Product description
The RapR-Round is a vest or leg wrap that wraps around the chest or leg. The Click-Tite® connectors provide the means to connect to the temperature controller.

Expected life
These RapR-Round products have a 30 day single patient use expected life upon first use and under wrap into use conditions.

Specifications

Material	polypropylene, ethylene vinyl acetate, silicone, polyvinyl chloride (PVC)				
Model	8001-061-530	8001-061-535	8001-061-540		
Size	small/medium chest	large chest	one size for the left or right leg - thigh circumference		
Length	32 to 46 in.	81 to 117 in.	46 to 54 in.	117 to 137 in.	52 to 72 cm
Width	20.5 to 28.5 in.				
Compatible controllers	Medi-Therm series temperature controller Altra® Precision Temperature Management System				
Kit	8001-061-550 8001-061-560 Note: Consult your local representative for controller availability.				
Environmental conditions	Operation			Storage and transportation	
Temperature					
Relative humidity					

Stryker reserves the right to change specifications without notice.

Symbols

	General warning		Caution
	Refer to instruction manual/booklet		Consult instructions for use
	Temperature controller		
	Do not puncture		Single patient use
	Batch code		Catalogue number
	Quantity		CE mark
	Date of manufacture		Manufacturer
	Authorized representative in the European Community		

US Federal Law restricts this device to sale by or on the order of a licensed physician.

Contact Information

Contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002 USA

To view your product manual online, see <https://techweb.stryker.com>.

Have the Lot batch code (A) (Figure 1) available when you call Customer Service. Include the Lot batch code in all written communication.

Safety Precautions

- WARNING**
 - Always check the product for leaks, tears, cuts, holes, stains, or any other damage prior to and during use.
 - Always check the patient's skin condition in areas of contact with the blanket according to hospital protocol or as directed by a physician.
 - Always turn or reposition the patient over the duration of therapy according to hospital protocol.
 - Always keep the area between the patient and the thermal transfer device dry.
 - Application of heat or cold are most effective when you apply the body wrap directly to the skin.
 - Do not allow hoses to fold, kink, or to wrap around a patient.

CAUTION

- Do not use if the package is damaged.
- Always use RapR-Round products with an approved Stryker temperature controller. Using this product with other controllers has not been tested.
- Do not use this product on another patient to avoid the risk of cross-contamination and infection.
- Do not use this product to position or transfer a patient.
- Always pre-fit the vest or leg wrap to lower the risk of pressure ulcers.
- Always follow hospital protocol when regulating temperature with external devices. Frequency of assessment and documentation will vary upon the individual response of the patient.
- Always consult the controller operations manual before connecting this product.
- Always complete a finger check, after applying a vest or wrap, to make sure that it is not on the patient's skin. This check is to avoid the risk of pressure ulcers.
- Do not clean this product. If soiled, dispose of product according to hospital protocol. This is a single patient use product.
- Do not use sharp objects or pins with this product.
- Do not block hard pathways of the thermal transfer device with restraints.

Return authorization

Product cannot be returned without prior approval from the Stryker Customer Service Department. An authorization number will be provided which must be printed on the returned product. Stryker reserves the right to charge shipping and restocking fees on returned product. Special, modified, or discontinued products are not subject to return.

Damaged merchandise

ICC Regulations require that claims for damaged product must be made with within fifteen (15) days of receipt of the product. Do not accept damaged shipments unless such damage is noted on the delivery receipt at the time of receipt. Upon prompt notification, Stryker will file a freight claim with the appropriate carrier for damages incurred. Claims will be limited in amount to the actual replacement cost. In the event that this information is not received by Stryker within the fifteen (15) day period following the delivery of the product, the damage was not noted on the delivery receipt at the time of receipt, the invoice will be responsible for payment for the original purchase price. Stryker does not accept claims for any incomplete shipments must be made within thirty (30) days of receipt.

International warranty clause

This warranty reflects U.S. domestic policy. Warranty outside the U.S. may vary by country. Please contact your local Stryker Medical representative for extra information.

Operations
Assessing the patient's skin condition
Note: Follow your facilities skin care protocol. The following instructions are recommendations.

- Document initial skin assessment.
- Check for any foreign items such as medical patches, IV's, sutures.
- Check for fluids such as moisture or gel.
- Make sure that the vest or leg wrap size is correct and pre-filled before applying to the patient.

Connecting and pre-filling the vest or leg wrap (Click-Tite connectors)

- Close the pinch clamps on the connector hose and vest or leg wrap (Figure 2).
- Attach the connector hoses to the Click-Tite connectors.
- Attach the connector hoses to the vest or leg wrap. Insert the male coupling of vest or leg wrap into the female coupling of hose with a twisting motion (Figure 3 and Figure 4).
- Attach the connector hoses to the vest or leg wrap (Figure 7).
- Insert the power cord into a wall outlet.
- Press the power button to turn on the controller. See the controller instructions for use.

9. Check that the vest or leg wrap fills with water.
Note: Isotonic based disinfectants will stain the product. Staining does not affect the patient nor the use of the product.

Applying the vest or leg wrap to the patient

- Apply the pre-filled vest or leg wrap directly to the patient's skin with the logo on the outside.
- Close the vest or leg wrap by fastening the Velcro® closures.
- Slide a finger along both sides of the vest or leg wrap (top and bottom). Check for constriction of any soft tissue.

Notes:

- Application of heat and cold are most effective when you apply body wraps directly to the skin.
- Make sure there are no folds or creases that may prevent the flow of water.
- Locate the control part of the leg wrap over the knee.
- Avoid positioning the vest wrap edges above the bony prominence or areas of sensitive skin.
- Place the hoses on the patient's left side. Do not wrap the hoses around or under the patient.
- Following the application of the vest or leg wraps, you may see the pattern of these products on the patient's skin. This pattern on the skin may redden (hyperemia). This occurrence is normal. The reddening will disappear within the same length of time that you applied the thermal transfer device. This phenomenon is often mistaken for injury, but is in fact a reactive hyperemia.^{1,2}

Rechecking the patient's skin condition
Check the patient at regular intervals as directed by hospital protocol. Open the Velcro® closures to check the entire area that is in contact with the vest or leg wrap. Note any change in the skin integrity that relates to:

- Excessive moisture - dry the skin surface by wiping away the moisture
- Color of the epidermis
- Skin texture
- Patent's skin condition is acceptable to continue therapy
- Make sure there are no folds or creases that may prevent the flow of water.
- Locate the control part of the leg wrap over the knee.
- Avoid positioning the vest wrap edges above the bony prominence or areas of sensitive skin.
- Place the hoses on the patient's left side. Do not wrap the hoses around or under the patient.
- Following the application of the vest or leg wraps, you may see the pattern of these products on the patient's skin. This pattern on the skin may redden (hyperemia). This occurrence is normal. The reddening will disappear within the same length of time that you applied the thermal transfer device. This phenomenon is often mistaken for injury, but is in fact a reactive hyperemia.^{1,2}

Discontinuing the vest or leg wrap from the hoses

- Power off the controller. Allow most of the water to drain back into the machine. (Approximately 10 minutes).
- Close vest or leg wrap clamps (Figure 2).
- Close the connector hose clamps.
- Unsnap the fitting. Pull the male coupling out to disconnect.
- Discard vest or leg wrap after use per facility protocol.

Note: Do not discard the controller connector hose after use, save for next use.

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