## Rapr-Round

### Operations Manual

<b>REF</b>	8001-061-630
<b>REF</b>	8001-061-635
<b>REF</b>	8001-061-640
<b>REF</b>	8001-061-650
<b>REF</b>	8001-061-660

EN	NL
CS	NO
DA	PL
DE	PT
ES	RU
FR	SV
IT	TR
KO	ZH

## R<sub>x</sub> ONLY

8001-009-151 Rev AB.1

2024-10

EN

### Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Single patient multiple reuse
	Do not puncture
<b>REF</b>	Catalogue number
<b>LOT</b>	Batch code
<b>CE</b>	CE mark
<b>UK CA</b>	UK Conformity Assessment mark
	Importer
<b>UDI</b>	Unique device identifier
<b>QTY</b>	Quantity
	Manufacturer
	Date of manufacture
<b>CH REP</b>	Authorized representative in Switzerland
<b>EC REP</b>	Authorized representative in the European Community
<b>MD</b>	European medical device
<b>R<sub>x</sub> ONLY</b>	<b>CAUTION</b> - US Federal Law restricts this device to sale by or on the order of a licensed physician.
	Temperature controller

### Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

#### WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

#### CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

**Note** - Provides special information to make maintenance easier or important instructions clearer.

### Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

#### WARNING

- Always check the product for leaks, tears, cuts, holes, stains, or any other damage prior to and during use.
- Always check the patient's skin condition in areas of contact with the product according to hospital protocol or as directed by physician.
- Always turn or reposition the patient over the duration of therapy according to hospital protocol.
- Always keep the area between the patient and the product dry.
- Always apply the product directly to the skin so application of heat or cold are most effective.
- Do not allow hoses to fold, kink, or to wrap around a patient.

#### CAUTION

- US Federal Law restricts this device to sale by or on the order of a licensed physician.
- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Do not use if the package is damaged.
- Do not connect two or more thermal transfer devices in a series on a single port.
- Always use **Rapr-Round** products with an approved temperature controller. The use of this product with other controllers has not been tested.
- Do not reuse this product on another patient to avoid the risk of cross-contamination and infection.
- Do not use this product to position or transfer a patient.
- Always prefill the product before you apply it to the patient to avoid the risk of pressure injuries or pressure ulcers.
- Always follow hospital protocol when regulating temperature with external devices. Frequency of assessment and documentation will vary upon the individual responses of the patient.
- Always consult the **Altrix** Operations Manual before you connect this product.
- Always check the placement and tightness of the vest or leg wrap according to hospital protocol.
- Do not clean this product. If soiled, dispose of product according to hospital protocol.
- Do not use sharp objects or pins with this product.
- Do not block fluid pathways of the thermal transfer device with restraints.

### Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

#### CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

#### Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

### Product description

The **Rapr-Round** hyper/hypothermia vest or leg wrap is applied around a patient chest or leg. The **Colder**® style connectors provide the means to connect to the **Altrix** Precision Temperature Management System (referred to as the controller).

### Intended use, Contraindications, Clinical benefits, and Environmental conditions

Always consult the **Altrix** Operations Manual before you use this product.

### Expected life

The **Rapr-Round** products have a 30 day single patient use expected life upon first use and under normal use conditions.

### Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

### Specifications

Material	Polypropylene, ethylene vinyl acetate, silicone, polyvinyl chloride (PVC)		
Model	8001-061-630	8001-061-635	8001-061-640

Size	Small/medium chest		Large chest		One size for the left or right leg - thigh circumference	
Chest size	32 in. to 46 in.	81 cm to 117 cm	46 in. to 54 in.	117 cm to 137 cm	20.5 in. to 28.5 in.	52 cm to 72 cm
Compatible controllers	Model 8001 <b>Altrix</b> Precision Temperature Management System					

Kit	Content	Quantity
8001-061-650	8001-061-630	5
	8001-061-640	10
8001-061-660	8001-061-635	5
	8001-061-640	10

Stryker reserves the right to change specifications without notice.

### Contact information


Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical  
3800 E. Centre Avenue  
Portage, MI 49002  
USA

**Note** - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established. To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the batch lot code number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

### Batch lot code

### Operation

#### Assessing the patient's skin condition

**Note** - Follow your hospital protocol for skin care. The following instructions are recommendations.

- Document initial skin assessment.
- Check for any foreign items such as medical patches, IVs, or sutures.
- Check for fluids such as moisture or gels.
- Verify the wrap size is correct and pre-filled before you apply to the patient.

#### Connecting and prefilling a vest or leg wrap

- Attach the connector hoses to the controller.
- Attach the controller hoses to the wrap. Insert the male coupling from the wrap into the female coupling of the hose (Figure 1).


Figure 1 - Connect couplings

- With a light pull, make sure that you lock all the connectors.
- Plug the power cord into a wall outlet.
- Press the power button to turn on the controller. See the **Altrix** Operations Manual.
- Check that the wrap fills with water.

**Note** - Iodine based disinfectants will stain the product. Staining does not affect the patient or the use of the product.

#### Applying the vest or leg wrap to the patient

- Apply the pre-filled vest or leg wrap directly to the patient's skin with the logo on the outside.
- Close the vest or leg wrap by fastening the hook and loop closures.
- Slide a finger along both sides of the vest or leg wrap (top and bottom). Check for constriction of any soft tissue.

#### Note

- Application of heat and cold are most effective when you apply body wraps directly to the skin.
- Make sure there are no folds or creases as those may prevent the flow of water.
- Locate the cutout part of the leg wrap over the knee.
- Avoid positioning the vest wrap edges above the bony prominence or areas of sensitive skin.
- Place the hoses on the patient's left side. Do not wrap the hoses around or under the patient.
- Following the application of the vest or leg wraps, you may see the pattern of these products on the patient's skin. This pattern on the skin may redden (hyperemic). This occurrence is normal. The reddening will disappear within the same length of time that you applied the thermal transfer device. This phenomenon is often mistaken for injury, but is in fact reactive hyperemia.<sup>1,2</sup>

#### References:

<sup>1</sup>Stewart, T.P., Magnano, S.J. Burns or Pressure Ulcers in the Surgical Patient? Decubitus 1988 | 3640

<sup>2</sup>Stewart, T.P., When a burn is not a burn. International Journal of Dermatology 1998 37: 729732

### Rechecking the patient's skin condition

Check the patient at regular intervals as directed by hospital protocol. Open the hook and loop closures to check the entire area of skin that is in contact with the vest or leg wrap. Note any change in the skin integrity that relates to:

- Excessive moisture - wipe away the moisture to dry the skin surface
- Color of the epidermis
- Skin texture
- Patient's skin condition is acceptable to continue therapy

**Note** - Check the vest or leg wrap position. Readjust the hook and loop closures, as necessary.

### Rechecking the vest or leg wrap

Recheck the vest or leg wrap for the following:

- Position
- Water status
- Dry surface
- Leaks
- Cracks

### Defibrillation considerations

- Open the vest to expose the patient's chest.
- Remove excess moisture.

### Disconnecting the vest or leg wrap from the hoses

- Power off the controller. Allow most of the water to drain from the wrap back into the controller (approximately 10 minutes).
- Press down on the thumb tab of the female coupling. Pull the male coupling out to disconnect (Figure 2).


Figure 2 - Disconnect couplings

- Discard the vest or leg wrap after use per hospital protocol.

**Note** - Do not discard the controller connector hose after use, save for next use.

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Altrix**, **Rapr-Round**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.