Altrix® Temperature Management Wraps

Operations Manual

REF	8003-002-001
REF	8003-002-002
REF	8003-002-003
REF	8003-002-004
REF	8003-002-005
REF	8003-002-006
REF	8003-003-001
REF	8003-003-002
REF	8003-003-003
REF	8003-003-004

R_x ONLY


EN
CS
DA
DE
ES
FR
IT
KO

NL
NO
PL
PT
RU
SV
TR
ZH

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	WARNING magnetic field ≤ 3T
	Caution
	Do not puncture
	Single patient multiple reuse
	Catalogue number
	Batch code
	MR conditional under the following conditions: 3T static magnetic field
	CE mark
	UK Conformity Assessment mark
	Importer

	Unique device identifier
	Manufacturer
	Date of manufacture
	Authorized representative in Switzerland
	Authorized representative in the European Community
	European medical device
	For US Patents see www.stryker.com/patents
	CAUTION - US Federal Law restricts this device to sale by or on the order of a licensed physician.
	Use-by date
	Quantity
	Altrix controller

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

WARNING

- Always check the product for leaks, tears, cuts, holes, stains, or any other damage prior to and during use.
- Always check the patient's skin condition in areas of contact with the torso or thigh wrap according to hospital protocol or as directed by a physician.
- Always turn or reposition the patient over the duration of therapy according to hospital protocol.
- Always apply the product directly to the skin so application of heat or cold are most effective.
- Do not allow hoses to fold, kink, or to wrap around a patient.

- Always keep the area between the patient and the torso or thigh wrap dry.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always follow hospital protocol when regulating temperature with external devices.
 - Do not use if the package is damaged.
 - Always use these products with the **Altrix** controller. Consult the **Altrix** controller Operations Manual before you connect this product.
 - Do not use this product if the patient has a transdermal medication (patch).
 - Avoid the use of materials of good thermal conductivity, such as water, gel, or similar substances with the **Altrix** controller not powered on.
 - Do not apply thermal transfer devices to patients with ischemic limbs. This may result in harm to the patient.
 - Do not reuse this product on another patient.
 - Do not use this product to position or transfer a patient.
 - Always wipe soil or contaminants from the wrap for patient safety.
 - Always check the placement and tightness of the torso or leg wrap according to hospital protocol.
 - Do not use sharp objects or pins with this product.
 - Do not block fluid pathways of the thermal transfer device with restraints.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The **Altrix** Temperature Management Wraps is a torso or thigh wrap that you can use individually or in combination. The **Colder®** style connectors provide the means to connect to the **Altrix** Precision Temperature Management System (referred to as the **Altrix** controller).

Indications for use, Contraindications, Clinical benefits, and Environmental conditions

Always consult the **Altrix** Operations Manual before you apply the wrap.

Expected service life

Each **Altrix** Temperature Management Wrap has a 14 day single patient multiple use expected life under normal use conditions.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

8003-002-001 small torso		
Size	26 in. - 47 in.	66 cm – 119 cm
Weight range	95–240 lb	43–109 kg
8003-002-002 large torso		
Size	35 in. - 60 in.	89 cm – 152 cm
Weight range	150–350 lb	68–159 kg
8003-002-005 large left thigh		
8003-002-006 large right thigh		
Size	16 in. - 28 in.	41 cm – 71 cm
Weight range	115–250 lb	52–113 kg
8003-002-003 small left thigh		
8003-002-004 small right thigh		
Size	19.5 in. - 32 in.	50 cm – 81 cm
Weight range	157–313 lb	71–142 kg
Compatible controller	Model 8001 Altrix Precision Temperature Management System	

European REACH

In accordance with the European REACH regulation and other environmental regulatory requirements the components that contain declarable substances are listed.

Description	Number	Substance of very high concern (SVHC) chemical name
Altrix wrap specification label	8003-001-908	Nonylphenol Ethoxylate

Material

Polyurethane (ether-based), Polyester Polyurethane Foam, Non-Woven Spun-Bonded Polyester Fabric (patient side)

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.


Stryker Medical
3800 E. Centre Avenue
Portage, MI 49002
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the lot batch code (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Lot batch code location

Operation

Assessing the skin condition of the patient

Follow your hospital protocol.

Applying the torso wrap to the patient


Figure 1 – Torso wrap

1. Position the torso wrap under the supine patient's rib cage (Figure 1).
Note - Avoid positioning the torso wrap edges above the bony prominences or areas of sensitive skin.
2. Close the wrap by fastening the hook and loop closures.
3. Route the hose to the lateral side of the patient.

Applying the thigh wrap to the patient


Figure 2 – Thigh wrap right and left

1. Position the thigh wrap just above the knee (Figure 2).
2. Close the wrap by fastening the hook and loop closures.

3. Route the hose to lateral side of the patient's knee.

²Stewart, T.P., When a burn is not a burn. International Journal of Dermatology 1998 37: 730-732

Checking for constriction of soft tissue

1. Slide a finger along both sides of the wrap (top and bottom).
2. Adjust the hook and loop closures, as necessary.

Connecting the hose to the Altrix Controller

1. Attach the connector hoses to the torso or thigh wrap. Insert the male coupling of torso or thigh wrap into the female coupling of hose. Press until you hear a click.
2. With a light pull, make sure that you lock all the connectors.
3. Check that the torso or thigh wrap fills with water.

Note

- When reading radiology images taken of a patient with the torso or thigh wrap, internal components may cause artifacts or distort readings.
- After the application of the wrap, the pattern of the wrap may appear on the skin. The reddening will disappear within the same length of time that you applied the wrap. This phenomenon is often mistaken for injury, but is in fact reactive hyperemia.^{1 2}

¹Stewart, T.P., Magnano, S.J. Burns or Pressure Ulcers in the Surgical Patient? Decubitus 1988; 1: 36-40

Reassessing the patient's skin

Follow your hospital protocol.

Defibrillation considerations

1. Open the torso wrap to expose the patient's chest.
2. Remove excess moisture.

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region.

Name	Part number
Altrix wraps, complete kit, small	8003-003-001 (8003-002-003, 8003-002-004, 8003-002-001)
Altrix wraps, thigh kit, small	8003-003-003 (8003-002-003, 8003-002-004)
Altrix wraps, complete kit, large	8003-003-002 (8003-002-005, 8003-002-006, 8003-002-002)
Altrix wraps, thigh kit, large	8003-003-004 (8003-002-005, 8003-002-006)

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Altrix**, **Stryker**. All other trademarks are trademarks of their respective owners or holders.