Altrix® Temperature Management Wraps™

Operations Manual
 Manual de uso
 Bedienungsanleitung
 Manuel d'utilisation
 Manuale d'uso
 Bedieningshandleiding
 Användarhandbok
 Betjeningsmanual
 Brukerhåndbok
 Manual de utilização
 Podręcznik użytkowania
 Руководство по эксплуатации
 Üzemeltetési utasítás
 Kullanma Kılavuzu
 Příručka pro obsluhu
 操作手册
 작동 설명서

REF 8003-002-001	REF 8003-003-001
REF 8003-002-002	REF 8003-003-002
REF 8003-002-003	REF 8003-003-003
REF 8003-002-004	REF 8003-003-004
REF 8003-002-005	
REF 8003-002-006	


Figure 1 - Torso wrap
 Figura 1 - Ortesis torácica
 Abbildung 1 - Rumpfwickel
 Figure 1 - Enveloppe-torse
 Figura 1 - Fascia toracica
 Figuur 1 - Torsowikkel
 Figur 1 - Omslag torso
 Figur 1 - Torsoomslag
 Figur 1 - Overkroppsomslag
 Figura 1 - Envolvimento do tronco
 Rys. 1 - Opaska na klatkę piersiową
 Рисунок 1 - Бандаж
 1. ábra - Törzsszorító
 Şekil 1 - Gövde sargısı
 Obrázek 1 - Tělový zábal
 图 1 - 胸部敷带
 그림 1 - 몸체 감싸개


Figure 2 - Thigh wrap (right and left)
 Figura 2 - Ortesis de muslo (derecha e izquierda)
 Abbildung 2 - Schenkelwickel (rechts und links)
 Figure 2 - Enveloppe-cuisse (droite et gauche)
 Figura 2 - Fascia per coscia (destra e sinistra)
 Figuur 2 - Dijwikkkel (rechts en links)
 Figur 2 - Omslag lår (höger och vänster)
 Figur 2 - Låromslag (højre og venstre)
 Figur 2 - Låromslag (høyre og venstre)
 Figura 2 - Envolvimento da coxa (direita e esquerda)
 Rys. 2 - Opaska na udo (prawa i lewa)
 Рисунок 2 - Бедренная манжета (правая и левая)
 2. ábra - Combszorító (jobb és bal)
 Şekil 2 - Uyluk sargısı (sağ ve sol)
 Obrázek 2 - Stehenní zábal (pravý a levý)
 图 2 - 大腿敷带 (右和左)
 그림 2 - 허벅지 감싸개(오른쪽 및 왼쪽)

English (EN)

Symbols

 General warning	 Caution	 Catalogue number	 Quantity
 Refer to instruction manual/ booklet	 Consult instructions for use	 CE mark	 Authorized representative in the European Community
 Altrix controller	 Single patient use	 Manufacturer	 Batch code
 Do not puncture	 MR conditional under the following conditions: 3T static magnetic field	 European medical device	 WARNING magnetic field
 CAUTION US Federal Law restricts this device to sale by or on the order of a licensed physician.	 Use-by date	 For US Patents see www.stryker.com/patents	

Introduction

This manual assists you with the operation of your Stryker product. Read this manual before operating this product. Set methods and procedures to educate and train your staff on the safe operation of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note: Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual.

Product description

The **Altrix® Temperature Management Wraps™** is a torso or thigh wrap that you can use individually or in combination. The Colder style connectors provide the means to connect to the **Altrix Precision Temperature Management System** (referred to as the **Altrix controller**).

Indications for use, Contraindications, Clinical benefits, and Environmental conditions

Always consult the **Altrix Precision Temperature Management System Operations Manual** before you apply the wrap.

Expected life

Each **Altrix Temperature Management Wraps** has a 14 day single patient use expected life under normal use conditions.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Material

Polyurethane (ether-based), Polyester Polyurethane Foam, Non-Woven Spun-Bonded Polyester Fabric (patient side)


Contact information

If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your product manual online, see <https://techweb.stryker.com>.

Have the Lot batch code (A) available when you call Customer Service. Include the Lot batch code in all written communication.

Small torso	26 in. - 47 in. (66-119 cm)	95-240 lb (43-109 kg)	8003-002-001
Large torso	35 in. - 60 in. (89-152 cm)	150-350 lb (68-159 kg)	8003-002-002
Small left thigh	16 in. - 28 in. (41-71 cm)	115-250 lb (52-113 kg)	8003-002-003
Small right thigh			8003-002-004
Large left thigh	19.5 in. - 32 in. (50-81 cm)	157-313 lb (71-142 kg)	8003-002-005
Large right thigh			8003-002-006

Small kits	8003-002-003	8003-002-004	8003-002-001
8003-003-001	x	x	x
8003-003-003	x	x	n/a
Large kits	8003-002-005	8003-002-006	8003-002-002
8003-003-002	x	x	x
8003-003-004	x	x	n/a

English (EN)

Safety precautions

WARNING

- Always check the product for leaks, tears, cuts, holes, stains, or any other damage prior to and during use.
- Always check the patient's skin condition in areas of contact with the torso or thigh wrap according to hospital protocol or as directed by a physician.
- Always turn or reposition the patient over the duration of therapy according to hospital protocol.
- Application of heat or cold are most effective when you apply the body wrap directly to the skin.
- Do not allow hoses to fold, kink, or to wrap around a patient.
- Always keep the area between the patient and the torso or thigh wrap dry.

CAUTION

- Always follow hospital protocol when regulating temperature with external devices.
- Do not use if the package is damaged.
- Always use these products with the **Altrix** controller. Consult the **Altrix** controller operations manual before you connect this product.
- Do not use this product if the patient has a transdermal medication (patch).
- Avoid the use of materials of good thermal conductivity, such as water, gel, or similar substances with the **Altrix** controller not powered on.
- Do not apply thermal transfer devices to patients with ischemic limbs. This may result in harm to the patient.
- Do not reuse this product on another patient.
- Do not use this product to position or transfer a patient.
- Always wipe soil or contaminants from the wrap for patient safety.
- Check the placement and tightness of the torso or leg wrap according to hospital protocol.
- Do not use sharp objects or pins with this product.
- Do not block fluid pathways of the thermal transfer device with restraints.

Assessing the patient's skin condition

Follow your hospital protocol.

Applying the torso wrap to the patient (see Figure 1)

1. Position the torso wrap under the supine patient's rib cage.
Note: Avoid positioning the torso wrap edges above the bony prominences or areas of sensitive skin.
2. Close the wrap by fastening the **Velcro®** closures.
3. Route the hose to the lateral side of the patient.

Applying the thigh wrap to the patient (see Figure 2)

1. Position the thigh wrap just above the knee.
2. Close the wrap by fastening the **Velcro®** closures.
3. Route the hose to the lateral side of the patient's knee.

Checking for constriction of soft tissue

1. Slide a finger along both sides of the wrap (top and bottom).
2. Adjust the **Velcro®** closures, if necessary.

Connecting the hose to the Altrix Controller

1. Attach the connector hoses to the torso or thigh wrap. Insert the male coupling of torso or thigh wrap into the female coupling of hose. Press until you hear a click.
2. With a light pull, make sure that you lock all the connectors.
3. Check that the torso or thigh wrap fills with water.

Note:

- When reading radiology images taken of a patient with the torso or thigh wrap, internal components may cause artifacts or distort readings.
- After the application of the wrap, the pattern of the wrap may appear on the skin. The reddening will disappear within the same length of time that you applied the wrap. This phenomenon is often mistaken for injury, but is in fact reactive hyperemia.^{1,2}

¹Stewart, T.P., Magnano, S.J. Burns or Pressure Ulcers in the Surgical Patient? *Decubitus* 1988; 1: 36-40

²Stewart, T.P., When a burn is not a burn. *International Journal of Dermatology* 1998 37: 730-732

Rechecking the patient's skin condition

Follow your hospital protocol.


Defibrillation considerations

1. Open the torso wrap to expose the patient's chest.
2. Remove excess moisture.

In accordance with the European REACH regulation and other environmental regulatory requirements the components that contain declarable substances are listed.

Description	Number	Substance of very high concern (SVHC) chemical name
Altrix wrap specification label	8003-001-908	Nonylphenol Ethoxylate

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