Unity™ Sleep-Eez™ Recliner

Operations/Maintenance Manual

REF 3500

Symbols

	Consult instructions for use
	Caution
	Catalogue number
	Serial number
	Distributed by
	Manufacturer
	For US Patents see www.stryker.com/patents
	Safe working load

Table of Contents

Warning/Caution/Note Definition	2
Summary of safety precautions	2
Introduction	3
Product description	3
Intended use	3
Expected service life	3
Specifications	3
Product illustration	4
Contact information	4
Serial number location	5
Operation	6
Locking and unlocking the casters	6
Reclining	6
Lying flat for a sleep surface	7
Accessories and parts	9
Preventive maintenance	10
Cleaning with SideKick® wipes	11
Cleaning the product	12
Service	13
Caster replacement	13
Ottoman cushion replacement	13
Seat cushion replacement	14
Back cushion replacement	15
Footboard replacement	16
Arm cap replacement	16
Recliner assembly	17

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note - Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - Always confirm stable foot placement and back support or get assistance from another person to avoid the risk of injury or product damage.
-

Introduction

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
-

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at 1-800-327-0770.

Product description

The **Sleep-Eez™** Recliner is a furniture product that functions as a padded seat, capable of reclining. The product is also capable of lying flat to function as a visitor sleep surface.

Intended use

The **Sleep-Eez** Recliner is intended to be used as furniture within the healthcare setting.

Expected service life

The **Sleep-Eez** Recliner has a 10 year expected service life under normal use conditions and with appropriate periodic maintenance.

Specifications


3500-000-691

	Safe working load	300 lb	136 kg
Overall length		37.0 in.	940 mm
Overall width		35.0 in.	889 mm
Overall height		35.0 in.	889 mm
Seat depth (minimum)		21.0 in.	533 mm
Seat width (minimum)		27.5 in.	699 mm
Seat height (minimum)		17.5 in.	445 mm
Back width		27.5 in.	699 mm
Back height (minimum)		18.5 in.	470 mm

Armrest clearance	27.5 in.	699 mm
Arm height (from seat)	9.0 in.	229 mm
Arm height (from floor)	26.2 in.	667 mm
Weight	146 lb	66.4 kg

Stryker reserves the right to change specifications without notice.

Product illustration

A	Strap
B	Footboard
C	Ottoman cushion
D	Arm
E	Seat cushion

F	Back cushion
G	Arm cap
H	Caster
I	Extension frame

Contact information

Contact Stryker Customer Service or Technical Support at: 1-800-327-0770.

Stryker Medical
3800 E. Centre Avenue

Portage, MI 49002
USA

To view your operations or maintenance manual online, see <https://techweb.stryker.com/>.

Have the serial number (A) of your Stryker product available when calling Stryker Customer Service or Technical Support. Include the serial number in all written communication.

Serial number location

The serial number for your product is located behind the seat cushion (A).


Figure 1 – Serial number location

Operation

Locking and unlocking the casters

The locking casters are located at the front of the bench.

To lock the casters, use your foot to push the caster lock (A) down (Figure 2).

To unlock the casters, use your foot to pull the caster lock up.


Figure 2 – Caster lock

Reclining

1. Lock the casters (*Locking and unlocking the casters* (page 6)).
2. Pull the strap (A) on the front of the recliner to slide out the extension frame (B) (Figure 3).
3. Pull the strap up and out to raise the ottoman cushion (C) to its upright position.
4. Reverse steps to store away.

Figure 3 – Recline

Lying flat for a sleep surface

1. After you recline (*Reclining* (page 6)), pull the strap (A) on the front of the recliner to slide out the extension frame (B) until it stops (Figure 4).

Note - The seat cushion should slide out when extended.

2. Grasp the top of the back cushion (C) and fold forward.
3. Reverse steps to store away.


Figure 4 – Lying flat

Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call Stryker Customer Service: 1-800-327-0770.

Name	Part number
Arm cap, black	3500-022-690
Arm cap, kalua	3500-023-690
Arm cap , latte	3500-024-690
Arm cap, mocha	3500-025-690
Arm cap, shadow	3500-026-690
Arm cap, toast	3500-027-690
Caster, lock	3500-100-054
Caster, swivel	3500-100-055
Cushion, back, assembly	3500-020-690
Cushion, ottoman	3500-014-690
Cushion, seat	3500-015-690
Foot kick	3500-017-690
Footrest frame kit, stage 3	3500-013-690
Roller frame kit, stage 1	3500-011-690
Screws, back cushion hinge to chair frame	3500-100-151
Screws, caster	3500-100-056
Screws, cushions	3500-100-077
Screws, foot kick	3500-100-056
Seat frame kit, stage 2	3500-012-690
Spacers, foot kick	3500-100-176

Preventive maintenance

Remove product from service before you perform the preventive maintenance inspection. Check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified personnel.

Inspect the following items:

- All fasteners are secure
- All casters are secure
- All casters roll and swivel freely
- Front casters lock
- No debris in casters
- Back, arms, seat, and foot section are not damaged
- Foot and seat sections slide out and stop freely

Product serial number:
Completed by:
Date:

Cleaning with SideKick® wipes

For the United States only. Call Stryker Customer Service: 1-800-327-0770 to confirm availability for your configuration or region.

Stryker's preferred wipes (2060-000-001 6" x 10" or 2060-000-002 9" x 12") include the following active ingredients:

- n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride - 0.154%
- n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride - 0.154%
- Isopropanol - 21.000%

Non-active ingredient: Ethylene Glycol Monobutyl Ether – < 3%

Note - For safety information, read the product label.

Follow these steps to clean the external product surface using **SideKick** wipes. **SideKick** wipes are not approved to disinfect your product.

1. Rub with a fresh, clean **SideKick** wipe to remove all visible soils. Repeat as necessary until the product is clean.

Note - Use as many wipes as necessary.

2. Allow the external product surface to air dry before you return the product to service.

Cleaning the product

Recommended cleaners:

- Quaternary solutions (active ingredient – ammonium chloride) that contain less than 3% glycol ether
- Chlorinated bleach solution (1 part bleach solution to 100 parts water, 520 ppm available chlorine solution)
- Isopropyl alcohol (active ingredient \leq 21%)

Recommended cleaning method:

1. Follow the cleaning solution manufacturer's dilution recommendations.
2. Hand wash all surfaces of the product with a cleaning solution.
3. Avoid oversaturation. Do not allow the product to stay wet longer than the chemical manufacturer's guidelines for cleaning.
4. Dry the product.
5. Verify proper operation before you return the product to service.

Service

Caster replacement

Tools required:

- #3 Phillips screwdriver
- Blanket or equivalent

Procedure:

1. Pull the recliner away from the wall.
2. Lay a blanket, or equivalent, on the floor behind the recliner.
3. Using two hands, lift up on the bottom of the recliner to flip the recliner onto its back.

CAUTION - Always confirm stable foot placement and back support or get assistance from another person to avoid the risk of injury or product damage.

4. Using a #3 Phillips screwdriver, remove the four screws that secure the caster to the recliner frame and remove the caster. Save the screws.
5. Reverse steps to reinstall.
6. Verify proper operation before you return the product to service.

Ottoman cushion replacement

Tools required:

- #3 Phillips screwdriver

Procedure:

1. Pull the strap (A) on the front of the recliner to slide out the extension frame (B) (Figure 5).
2. Pull the strap up and out to raise the ottoman cushion (C) to its upright position.
3. Using a #3 Phillips screwdriver, reach through the extension frame to remove the four screws that secure the ottoman cushion to the extension frame. Save the screws.
4. Remove the ottoman cushion.
5. Reverse steps to reinstall.
6. Verify proper operation before you return the product to service.

Figure 5 – Ottoman cushion replacement

Seat cushion replacement

Tools required:

- #3 Phillips screwdriver

Procedure:

1. Pull the strap (A) on the front of the recliner to slide out the extension frame (B) (Figure 6).
2. Pull the strap up and out to raise the ottoman cushion (C) to its upright position.
3. Pull the strap to slide out the extension frame until it stops.

Note - The seat cushion should slide out when extended.

4. Using a #3 Phillips screwdriver, reach through the seat frame (F) to remove the four screws that secure the seat cushion (D) to the seat frame. Save the screws.
5. Remove the cushion.
6. Reverse steps to reinstall.
7. Verify proper operation before you return the product to service.


Figure 6 – Seat cushion replacement

Back cushion replacement

Tools required:

- #3 Phillips screwdriver

Procedure:

1. Pull the strap (A) on the front of the recliner to slide out the extension frame (B) (Figure 6).
2. Pull the strap up and out to raise the ottoman cushion (C) to its upright position.
3. Pull the strap to slide out the extension frame until it stops.

Note - The seat cushion should slide out when extended.

4. Using a #3 Phillips screwdriver, reach through the seat frame (F) to remove the two bottom screws (one on each side) that secure the back cushion hinge mechanism to the recliner frame. Save the screws.
5. Grasp the top of the back cushion (E) and fold forward.
6. Using a #3 Phillips screwdriver, remove the two top screws (one on each side) that secure the back cushion hinge mechanism to the recliner frame. Save the screws.
7. Remove the back cushion assembly.
8. Reverse steps to reinstall.

Note

- The left and right stampings on the hinges are a directional reference as you look at the front of the chair. The hinges have two screw holes on each side. Reinstall the screws through the inner holes.
 - Use caution when you slide the hinges down to install the new back cushion assembly to avoid the risk of damage to the fabric.
9. Verify proper operation before you return the product to service.

Footboard replacement

Tools required:

- #3 Phillips screwdriver

Procedure:

1. Pull the strap (A) on the front of the recliner to slide out the extension frame (B) (Figure 5).
2. Pull the strap up and out to raise the ottoman cushion (C) to its upright position.
3. Using a #3 Phillips screwdriver, reach through the extension frame to remove the four screws that secure the footboard to the extension frame. Save the screws.
4. Remove the footboard.
5. Reverse steps to reinstall.
6. Verify proper operation before you return the product to service.

Arm cap replacement

Tools required:

- Rubber hammer
- Large flat screwdriver

Procedure:

1. Using a rubber hammer, hit the back of the arm cap (A) toward the front of the recliner until it stops moving forward (Figure 7).

Note - Use caution when you hit the arm cap to avoid the risk of damage to the fabric.

2. Position a large flat screwdriver between the arm cap and arm to pull up on the arm cap and loosen the arm cap from the recliner. Remove and discard the arm cap.
3. Align the new arm cap bolts with the holes at the front of the arm slots.
4. Using a rubber hammer, hit the arm cap down until it is tight to the arm.
5. Using a rubber hammer, hit the arm cap toward the back of the recliner until it stops moving back.

Note - Use caution when you hit the arm cap to avoid the risk of damage to the fabric.


6. Verify proper operation before you return the product to service.

Figure 7 – Arm cap

Recliner assembly

3500-000-691 (Reference only)

Item	Name	Quantity
A	Caster, swivel	2
B	Arm cap with wood insert	2
C	Extension frame	1
D	Back cushion, with hinge, upholstered	1
E	Ottoman cushion with strap	1
F	Caster with brake	2

Item	Name	Quantity
G	Metal retainer plate	8
H	Arm, upholstered, left	1
I	Arm, upholstered, right	1
J	Back frame	1
K	Footboard	1
L	Seat cushion, upholstered	1

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