



1. PURPOSE

This document provides the list of service parts for the LIFEPAK 15 monitor/defibrillator (V4 hardware configuration only) and provides cross references between reference designators (ex. A01) used in the LIFEPAK 15 V4 Service Manual (3316925-008), and part and catalog numbers.

2. SCOPE

Scope is limited to individual parts and service repair kits which are available for service purposes for the LP15 V4 Service Manual (3316925-008). Kit contents are not listed in this document and exist in electronic bills of materials (eBOMs) in the Product Lifecycle Management (PLM) system.

3. SERVICE BOM PART

Item Number	Description
3343060-000	Bill of Material, Service, LIFEPAK 15 V4

See 3343060-000 EBOM for parts list or reference the tables in Sections 4 and 5.

4. SERVICE INDIVIDUAL PARTS LIST

Note: Reference Designators are unique identifiers for the part used to cross-reference between the Service Manual or Service Process Work Instruction and the Service BOM. The reference designators have been assigned based on the following prefix system:

- A prefix- Printed Circuit Board Assemblies (PCBAs), keypads, and modules
- W prefix – Cables, subassemblies
- F prefix – Fasteners (e.g. screws, nuts, and standoffs)
- M prefix – Plastic parts and mechanical components
- E prefix – Electrical components
- L prefix – Labels and adhesives
- K prefix – Kits, or parts assembled into a kit for a specific repair.

Note: For some reference designators, there are multiple language options (e.g. labels and keypads). Refer to the tables below to select the version that applies.

Click on a link to navigate to a parts list:

- [External Parts List](#)
- [Carrying Case Parts List](#)
- [Front Case Parts List](#)
- [System/Therapy Assembly Parts List](#)
- [Parameter Bezel Parts List](#)
- [Rear Case Parts List](#)
- [SpO2 Internal Parts List](#)
- [NIBP Parts List](#)
- [CO2 Parts List](#)
- [Printer Control Keypad Parts List](#)
- [Main Keypad Parts List](#)
- [LIFEPAK 15 Label Set Languages Parts List](#)
- [3rd Edition Label Set Parts List](#)
- [Service Repair Kits](#)



External Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07 or L08 (see notes)	Multiple	Multiple	LABEL SET, language and feature-specific	Refer to table for 3rd Edition Label Sets Parts List (L08) for multiple language options.
A10	Multiple	Multiple	A10 - KEYPAD ASSY - MAIN CNTRL	Refer to Main Keypad Parts List for multiple language options.
A09	Multiple	Multiple	A09 - KEYPAD ASSY - PRINTER CONTROL	Refer to Printer Control Keypads Parts List for multiple language options.
M01	3012120-004	21300-007445	COVER - C02 CONNECTOR	
A12	3208290-001	11241-000009	PRINTER-100MM, GRAY, WITH FERRITES	Part replaced with 11241-000017
A12	3344311-000	11241-000017	PRINTER-100MM, GRAY, WITH FERRITES	
M02	3206959-001	21330-001169	BRACKET - COVER, OEM MODULES	Door - Blank, USB
M03	3207307-001	21300-007356	GUARD - LOWER RIGHT	Bumper (guard), lower right corner
M04	3206968-000	21300-007010	GUARD - CORNER, LOWER LEFT	Bumper (guard), lower left corner
M05	3207247-001	21300-007319	Foot - Mounting	
F01	3207361-375	21300-007253	SCREW-CAP, SCH, REC T15 TORX 6-32, 0.375L, SS, NYLOK	
L09	3207131-001	21501-001997	LABEL - SERIAL NUMBER, V4	Orderable only through Service Depot
L01	3207359-001	21501-001993	LABEL-SPO2 MASIMO PATENT	
F02	201407-069	21300-000777	SCREW, SELF-SEAL, SELF-LOCK, 6-32X0.375	



L02	3320539-000	21501-002882	Label - UDI	Orderable only through Service Depot
-----	-------------	--------------	-------------	--------------------------------------

Carrying Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F03	3207492-001	21300-007429	SCREW-PHILLIPS,TRUSSHEAD,8-32,.312 LONG,SS	Mounts carry case at device bottom
F04	200476-708	21300-007966	SCREW-PH,P,SS,PV,1/4-20 X 1.00 L	Mounts carry case at top holes
F05	202253-576	21300-001033	SCREW,MACHINE, PANHEAD,NYLOK,6-32 X .750L	Mounts on rear pouch
K28	3207824-001	11577-000002	KIT - CARRY BAG, MAIN BAG	Standard Carry Case with Shoulder Strap
K29	3207827-003	11260-000039	KIT - CARRY BAG, REAR POUCH, 3RD EDITION	Back Pouch
K30	3207823-000	11577-000001	KIT - CARRY BAG, SHOULDER STRAP	Shoulder Strap
K31	3010267-00	11220-000028	POUCH, TOP, CARRYING BAG, LP12	Top Pouch

Front Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F06	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
F07	202253-764	21300-004599	SCREW,MACHINE,PANHEAD,NY LOK,4-40 X .500	
M06	3006809-00	21300-001299	GASKET, PRINTER, SUB-D	
M07	3006810-01	21300-006100	BRACKET, CONNECTOR, SUBD	



Front Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
W16	3009724-00	21330-000165	CABLE ASSY, RIBBON,INTERFACE PCB/CHART RCDR	Part # 3009724-00 to be used up before using part # 3009724-002
	3009724-002			
W19	3009726-01	21300-001370	WIRE HARNESS -GROUND/PRINTER BRKT	
W17	3009726-03	21300-001371	WIRE HARNESS-SPEAKER INTERFACE PCB	Part # 3009726-03 to be used up before using part # 3009726-000
	3009726-000			
W15	3011128-002	21300-003713	SWITCH-ROTARY	Part # 3011128-002 to be used up before using part # 3349211-000
	3349211-000	21300-008390		
M08	3012693-00	21300-001575	SPRING CLAMP, STAINLESS STEEL	
A11	3203221-002	21300-007358	DISPLAY-ACTIVE MATRIX, 8.4 ,COLOR	Part replaced with 21300-008140
A11	3319977-000	21300-008140	DISPLAY, ACTIVE MATRIX, 8.4 IN, 262K COLOR	See Display Repair Kit (K15)
M09	3206724-000	21330-001173	LENS-DISPLAY	See Front Case Repair Kit (K08)
M10	3206725-001	21330-001174	SHIELD-DISPLAY	See Display Shield Repair Kit (K04)
M11	3206726-001	21330-001168	FRAME - DISPLAY SHIELD	Part replaced with 21330-001697
	3206726-002	21330-001641	FRAME - DISPLAY SHIELD	Part replaced with 21330-001697
	3206726-003	21330-001697	FRAME - DISPLAY SHIELD, V5	See Display Shield Repair Kit (K04)
M12	3305427-000	21330-001423	ASSEMBLY - DISPLAY SHIELD	



Front Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A05	3347646-000	21330-001652	PCB ASSY – INTERFACE, LP15, V4	See Interface PCBA Repair Kit (K23)
M13	3206936-002	21330-001189	BRACKET - LCD DISPLAY MOUNTING	
W12	3206989-000	21330-001214	CABLE ASSY - SM KEYPAD, INTERFACE PCB	
W18	3206990-001	21330-001417	CABLE ASSY - FLEX, LCD/INTERFACE PCB	
W04	3206991-003	21330-001226	CABLE ASSY - FLEX, SYSTEM PCB/INTERFACE PCB	Part # 3206991-003 to be used up before using part # 3206991-004
	3206991-004			
W06	3318862-000	21330-001525	CABLE ASSY - BACKINV, INTERFACE PCB, ROHS	
M14	3009277-01	21300-001340	KNOB, ENERGY SELECT,PLAIN	Replaced with 21300-007051
M14	3207030-000	21300-007051	KNOB-ROTARY SWITCH	
W11	3312554-004	21330-001510	WIRE HARNESS-THERAPY RECEPTACLE	See Therapy Connector Repair Kit (K13)
	3312554-005			Part # 3312554-004 to be used up before using part # 3312554-005
M15	3207252-001	21300-007079	SHIELD,INVERTER	
F08	3207337-312	21300-007297	SCREW-MACH,PNH,PHH,NYLOCK,4-40,0.312L,WSHR,CS,ZN	
F01	3207361-375	21300-007253	SCREW-CAP,SCH,REC T15 TORX 6-32,0.375L,SS,NYLOK	
F09	3207367-312	21300-007196	SCREW-MACH,TRH,T10 TORX, 4-40,0.312L,SS	
W13	3207388-001	21330-001323	CABLE-ASSY,MAIN KEYPAD,INTERFACE PCB	



Front Case Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
M16	3207701-000	21300-007464	SEAL-THERAPY RECEPTACLE, LP15	See Therapy Connector Repair Kit (K13)
M17	3207725-002	21330-001273	ASSEMBLY-ENCLOSURE,FRONT	See Front Case Repair Kit (K08)
A08	3319108-000	21300-008138	DRIVER,LED BACKLIGHT,12VDC,85C	See Display Driver Repair Kit (K14)
L03	3208003-000	21501-002038	LABEL-PHYSIO-CONTROL ICON	Replaced with 21501-003065
	3345260-000	21501-003065	LABEL, V4+, LP15	
L04	800943-09	21501-000248	LABEL-SYMBOL INTL	
M18	804234-03	21300-002394	SEAL-PERIMETER,CASE	

System/Therapy Assembly Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F06	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
A04	3311744-011	21330-001538	PCB ASSY - THERAPY, LP15 V4	See Therapy PCBA Repair Kit (K12)
W03	3009878-002	21300-001379	CONN-HDR, SQUAREPIN DUAL ROW AU PLATED, 30 PIN	
A01	3342432-001	21330-001564	PCB ASSY - SYSTEM, LP15 V4	See System PCBA Repair Kit V4 (K26)
N/A	3207575-001	21300-007553	SHIELD - SYSTEM PCBA, SBC, TOP	
M19	3312315-001	21300-008205	BRACKET ASSY-MOUNTING, SYSTEM-THERAPY PCBA	See System PCBA Repair Kit V4 (K26)
L05	804447-27	21300-003172	SPACER-FOAM, SYS PCBCAPS	



System/Therapy Assembly Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
E01	202305-000	21300-001052	BATTERY-LITH,3V,COIN	

Parameter Bezel Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L06	3207407-000	21300-007362	SPACER - CONNECTOR, CO2	This spacer is used in devices without CO2 feature.
M20	3206723-005	21330-001211	BEZEL - PARAMETER	
W07	3318313-001	21330-001527	W07 - RECEPTACLE-ECG,12 CONTACT, WIRED, ROHS	See ECG Connector Repair Kit (K18)
M21	805915-01	21300-002683	SEAL-ECG RECEPTACLE	
M22	3205680-001	21300-006304	BRACKET - CABLE SUPPORT	
F10	202253-763	21300-006965	SCREW,MACHINE,PANHEAD,NYLOK,4-40 X .437	
W22	3206995-004	21330-001290	CABLE ASSY - FLEX, SPO2, MASIMO	See OEM PCBA Repair Kit (K22)
	3206995-006			Part # 3206995-004 to be used up before using part # 3206995-006
M23	3207009-000	21300-007112	HOUSING-FERRITE,SPO2	
F11	3205311-001	21300-007557	SCREW-MACHINE, PAN, TORX , NYLOK, 4-40, .6875L	QTY 4 for SpO2
W28	3317964-004	21330-001529	FRS ASSEMBLY-CO2, MINI MODULE, ROHS	See CO2 Connector Repair Kit (K19)
M24	3007997-01	21300-001314	SEAL-CO2 CONNECTOR	MiniMedi



Parameter Bezel Parts List				
M25	3012119-02	21300-001548	ADAPTER-CO2 CONNECTOR	MiniMedi
M26	3012121-002	21300-007998	RETAINER-CO2 CONN	MiniMedi
M70	3340287-000	21300-008212	RETAINER CLIP, NIBP COUPLER, LP15	
F12	202253-732	21300-007353	SCREW-MACH,NYLOK,PNH,PHH ,2-56X0.500,CS,ZN	QTY 4 for CO2
M27	3207033-001	21300-007425	CONNECTOR - PNEUMATIC COUPLER, NIBP	See NIBP Connector Repair Kit (K07)
M28	3012180-02	21300-001560	TUBING-NIBP	
M29	3012333-001	21300-001566	CONN-PNEU,STR,0.125IDTBG ,NYL,WHT	For NIBP
W33	3318652-000	21330-001526	WIRE HARNESS-INVASIVE PRESSURE 1 AND 2, ROHS	See Invasive Pressure Connector Repair Kit (K16)
M30	3007998-00	21300-001315	GASKET-CONNECTOR,INVASIVE PRESSURE	QTY 2
F08	3207337-312	21300-007297	SCREW-MACH,PNH,PHH,NYLOC K,4-40,0.312L,WSHR,CS,ZN	QTY 8 for IP
W35	3319355-000	21330-001522	CABLE ASSEMBLY - TEMPERATURE, ROHS	Part # 3319355-000 to be used up before using part # 3319355-001
	3319355-001			
F18	3340265-000	21300-008213	WASHER, NIBP COUPLER, LP15	

Rear Case Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A03	3336515-000	21330-001555	PCB ASSY - POWER, LP15	Order kit K32
A05	3317231-001	21330-001537	PCB ASSY - INTERFACE, LP15 V4	Order kit K23
A06	3317232-000	21330-001532	PCB ASSY - OEM, LP15	Order kit K22
A07	3207037-002	21330-001288	CONTACT ASSY-BATTERY	



Rear Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A15	3200846-002	21300-004232	CAPACITOR-ENERGY STORAGE , 2.25 IN DIA	Part # 3200846-002 to be used up before using part # 3200846-003
	3200846-003			
M31	3207726-000	21330-001272	ASSEMBLY - ENCLOSURE, REAR	Order kit K09
M32	3207881-000	21330-001276	ASSEMBLY-RETAINER, BATTERY	
M33	3206965-001	21330-001195	BRACKET - COVER, OEM MODULES	
M34	3206964-001	21330-001194	BRACKET - MOUNTING, OEM MODULES	
M35	3206961-003	21330-001257	BRACKET ASSY - POWER PCB	
M36	3207708-000	21300-007472	BRACKET-SHOULDER STRAP	
M37	3318613-000	21300-008135	MICRO WIRE SADDLE, TEAR DROP,NYLON 6/6,ROHS	
M38	802278-02	21300-002023	CONN-PLUG,BANANA PIN,LP10	Part replaced with 21300-008193
M38	3335400-000	21300-008193	PLUG, MINI-BANANA, NYLONPATCH	
M39	3006766-003	21330-001241	COVER - LATCH, PADDLE, GRAY	Part of Rear Case Assembly. See Paddle Retainer Repair Kit (K03)
M40	3207746-000	21300-007470	COVER PLATE - HANDLE, RIGHT	
M41	3207707-000	21300-007469	COVER PLATE-HANDLE, LEFT	
M42	3206960-002	21330-001190	GROUND PLANE - FLEXIBLE REAR CASE	
M43	3207706-000	21300-007471	HANDLE	Order kit K06
F13	806091-00	21300-002691	NUT-AUX,SYSTEM,CONN,STAINLESS STEEL	
F14	201508-000	21300-000804	NUT-LOCK,CARBON STEEL,KEP,#4-40 THREAD	Replaced with 21300-000584



Rear Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F14	200805-000	21300-000584	NUT-HEX,SS,LOCK 4-40X.250W	
M44	805487-00	21300-002650	PLATE-SEAL,CONNECTOR,REAR	
M45	200536-011	21300-000504	RETNR-CABLE TIE,NYLON, . 35W X 21.0 L	
F06	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
F15	202253-760	21300-006251	SCREW-M,PH,NYLOK,CS,4-40 ,.250L	
F08	3207337-312	21300-007297	SCREW-MACH,PNH,PHH,NYLOCK,4-40,0.312L,WSHR,CS,ZN	
F02	201407-069	21300-000777	SCREW,SELF-SEAL,SELF-LOCK,6-32X0.375	
F10	202253-763	21300-006965	SCREW,MACHINE,PANHEAD,NYLOK,4-40 X .437	
M46	200060-011	21300-000203	SEAL-O-RING,RUBBER, .551 ID	
M47	3006291-00	21300-001284	SEAL, DRAIN	
M59	3009787-003	21300-004110	SHIELD - DIELECTRIC, CAPACITOR	
M48	3318900-001	21300-008127	SHIELD-DIELECTRIC, OEM, POWER	
M08	3012693-00	21300-001575	SPRING CLAMP, STAINLESS STEEL	
F16	3010569-02	21300-001402	STANDOFF, HEX, NYLON	
M49	3012345-00	21300-001567	STRAP-GROUND,BTE PCB	
W01	3207692-000	21300-007456	CABLE ASSY - POWER,SYSTEM, PCBA	
W02	3009726-05	21300-001373	WIRE HARNESS-POWER-THERAPY-PCB	
W05	3319354-000	21330-001316	CABLE ASSEMBLY- CONTACT PCB/POWER PCB,LP15 ROHS	
W08	3009652-005	21330-001554	HARNESS, CONNECTOR, SYSTEM	Replaced with 21330-000156
W08	3009652-01	21330-000156	CABLE ASSY- SYSTEM CONNECTOR	
W09	3303848-002	21330-001485	CABLE ASSEMBLY - AUXILIARY POWER	



Rear Case Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
W10	3303863-000	21330-001394	CABLE ASSEMBLY - BATTERY POWER	
W14	3318946-000	21330-001534	FLEX ASSY - USB, ROHS	Part # 3318946-000 to be used up before using part # 3318946-001
	3318946-001			
M50	3312617-002	21300-008120	RETAINING END CAP-THERAPY CAPACITOR	
A24	3314519-001	21330-001536	PCB ASSY - CAP DISCHARGE , LP15 V4	
M51	3312615-001	21300-008122	RETAINING BASE-THERAPY CAPACITOR	
M52	3312614-001	21300-008121	RETAINING BRACKET-THERAPY CAPACITOR	
W24	3312616-003	21330-001500	WIRE HARNESS, THERAPY CAP	
M53	3317079-000	21300-008128	FOAM WITH ADHESIVE	
M54	3318064-001	21300-008126	DIELECTRIC SHIELD, SYSTEM PCBA	Part number 3318064-001 to be used up before using part number 3349630-000.
	3349630-000	21330-001661	ASSEMBLY, DIELECTRIC SHIELD,SYSTEM PCBA	
M55	200536-001	21300-000499	RETAINR-CABLE TIE,NYLON, .10 W, X 4.0 L	
M60	3336625-000	21300-008250	MOUNTING BRACKET- OEM MODULES	
M61	3338476-000	21300-008249	ASSEMBLY- COVER, OEM MODULES	



SpO2 Internal Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F17	202253-760	21300-006251	SCREW-M,PH,NYLOK,CS,4-40 ,.250L	
F07	202253-764	21300-004599	SCREW,MACHINE,PANHEAD,NYLOK,4-40 X .500	
M56	3206962-001	21300-007550	BRACKET - MOUNTING, SPO2 MODULE	
W21	3318916-000	21330-001524	CABLE, ASSY-FLEX, SPO2, OEM PCB, ROHS	Part number 3318916-000 to be used up before using part number 3318916-001.
	3318916-001			
A16	3317849-001	21330-001530	MODULE-OEM,PULSE OXIMETER, MX-51, ROHS	Order kit K20

NIBP Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F15	202253-760	21300-006251	SCREW-M,PH,NYLOK,CS,4-40 ,.250L	
M57	3010805-000	21300-001416	SOCKET RETAINER CLIP -10 PIN	
W27	3012181-00	21330-000235	CABLE ASSY-RIBBON,NIBP/OEM PCB	
A21	3317965-002	21330-001528	MODULE-NIBP, MAXIQ, ROHSROHS	Replaced with 21330-001660
A21	3317965-003	21330-001556	MODULE-NIBP, ORION, ROHS	Replaced with 21330-001660
A21	3346604-000	21330-001660	MODULE-NIBP, ORION, ROHS	See NIBP Module Repair Kit (K37). Part number 3346604-000 to be used up before using part number 3346604-001.
	3346604-001			



MiniMedi CO2 Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
F15	202253-760	21300-006251	SCREW-M,PH,NYLOK,CS,4-40 ,.250L	
M58	3012140-002	21300-001557	EXHAUST TUBING-CO2	
A23	3012140-006	21300-006948	A23 - MODULE, CO2, MINI	Part replaced with 21330-001535
A23	3317964-003	21330-001535	MODULE - MINIMEDI CO2, ROHS	Order kit K34 & K35
W30	3012397-01	21330-000237	CABLE ASSY-CO2, ADAPTER	
W26	3206998-002	21330-001217	CABLE ASSY-CO2 MODULE, OEM PCB	

NanoMedi CO2 Parts List				
Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
M62	3336628-000	21300-008224	RETAINER, CONNECTOR, CO2 NANOMEDI	
M63	3345277-000	21300-008245	ADAPTER ASSEMBLY - CO2 CONNECTOR	
M64	3345235-000	21300-008247	SEAL - CO2 CONNECTOR, NANOMEDI	
W29	3346405-000	21330-001650	CONNECTOR ASSY, FLR, NANOMEDICO2, MEDTRONIC, LP15	Order kit K34
A25	3336624-000	21300-008225	MODULE, CO2, NANOMEDICO2, MEDTRONIC, LP15	Order kit K34
M60	3336625-000	21300-008250	MOUNTING BRACKET - OEM MODULES	
M61	3338476-000	21300-008249	ASSEMBLY - COVER, OEM MODULES	
M65	3345604-000	21300-008251	INTAKE TUBE COVER, CO2, NANOMEDICO2 LP15	
M66	3342395-000	21300-008248	FERRITE, CLAMP ON,OD.465IN, ID.169IN, LG.913IN	



NanoMedi CO2 Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
W31	3336626-000	21330-001649	ASSY, CABLE, CO2 MODULE TO PCBA, LP15	
M67	3345815-000	21300-008252	EXHAUST TUBING CO2, NANOMEDI, INNER	
M68	3345815-001	21300-008253	EXHAUST TUBING CO2, NANOMEDI, OUTER	
M69	3345815-002	21300-008254	EXHAUST TUBING CO2, NANOMEDI, FITTING	

Printer Control Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A09	3207080-000	21330-001238	KEYPAD ASSY - PRINTER, 12 LEAD CONTROL, ENGLISH	
A09	3207080-001	21330-001264	KEYPAD ASSY-PRINTER CONTROL, ENGLISH	
A09	3207080-040	21330-001326	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, GERMAN	
A09	3207080-041	21330-001327	KEYPAD ASSY-PRINTER CONTROL, GERMAN	
A09	3207080-080	21330-001330	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, FRENCH	
A09	3207080-081	21330-001331	KEYPAD ASSY-PRINTER CONTROL, FRENCH	
A09	3207080-120	21330-001334	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, SPANISH	
A09	3207080-121	21330-001335	KEYPAD ASSY-PRINTER CONTROL, SPANISH	
A09	3207080-060	21330-001328	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, ITALIAN	
A09	3207080-061	21330-001329	KEYPAD ASSY-PRINTER CONTROL, ITALIAN	



Printer Control Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A09	3207080-140	21330-001383	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, PORTUGUESE	
A09	3207080-141	21330-001395	KEYPAD ASSY-PRINTER CONTROL, PORTUGUESE	
A09	3207080-180	21330-001336	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, SWEDISH	
A09	3207080-181	21330-001337	KEYPAD ASSY-PRINTER CONTROL, SWEDISH	
A09	3207080-200	21330-001338	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, DANISH	
A09	3207080-201	21330-001339	KEYPAD ASSY-PRINTER CONTROL, DANISH	
A09	3207080-220	21330-001384	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, FINNISH	
A09	3207080-221	21330-001402	KEYPAD ASSY-PRINTER CONTROL, FINNISH	
A09	3207080-240	21330-001340	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, NORWEGIAN	
A09	3207080-241	21330-001341	KEYPAD ASSY-PRINTER CONTROL, NORWEGIAN	
A09	3207080-100	21330-001332	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, DUTCH	
A09	3207080-101	21330-001333	KEYPAD ASSY-PRINTER CONTROL, DUTCH	
A09	3207080-260	21330-001342	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, POLISH	
A09	3207080-261	21330-001343	KEYPAD ASSY-PRINTER CONTROL, POLISH	
A09	3207080-280	21330-001385	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, HUNGARIAN	
A09	3207080-281	21330-001398	KEYPAD ASSY-PRINTER CONTROL, HUNGARIAN	
A09	3207080-300	21330-001386	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, CZECH	
A09	3207080-301	21330-001401	KEYPAD ASSY-PRINTER CONTROL, CZECH	



Printer Control Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A09	3207080-160	21330-001441	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, BRAZILIAN	
A09	3207080-161	21330-001442	KEYPAD ASSY-PRINTER CONTROL, BRAZILIAN	
A09	3207080-320	21330-001443	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, RUSSIAN	
A09	3207080-321	21330-001444	KEYPAD ASSY-PRINTER CONTROL, RUSSIAN	
A09	3207080-380	21330-001445	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, KOREAN	
A09	3207080-381	21330-001446	KEYPAD ASSY-PRINTER CONTROL, KOREAN	
A09	3207080-400	21330-001504	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, JAPANESE	
A09	3207080-401	21330-001505	KEYPAD ASSY-PRINTER CONTROL, JAPANESE	

Main Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A10	3302470-004	21330-001496	KEYPAD ASSY-MAIN CONTROL, NIBP, ENGLISH	Replaced with 21330-001603
A10	3302470-006	21330-001603	KEYPAD ASSY-MAIN CONTROL, NIBP, ENGLISH	
A10	3302470-005	21330-001497	KEYPAD ASSY-MAIN CONTROL, ENGLISH	Replaced with 21330-001604
A10	3302470-007	21330-001604	KEYPAD ASSY-MAIN CONTROL, ENGLISH	
A10	3302470-044	21330-001450	KEYPAD ASSY-MAIN CONTROL, NIBP, GERMAN	Replaced with 21330-001605
A10	3302470-046	21330-001605	KEYPAD ASSY-MAIN CONTROL, NIBP, GERMAN	
A10	3302470-045	21330-001451	KEYPAD ASSY-MAIN CONTROL, GERMAN	Replaced with 21330-001606



Main Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A10	3302470-047	21330-001606	KEYPAD ASSY-MAIN CONTROL, GERMAN	
A10	3302470-064	21330-001452	KEYPAD ASSY-MAIN CONTROL, NIBP, ITALIAN	Replaced with 21330-001607
A10	3302470-066	21330-001607	KEYPAD ASSY-MAIN CONTROL, NIBP, ITALIAN	
A10	3302470-065	21330-001453	KEYPAD ASSY-MAIN CONTROL, ITALIAN	Replaced with 21330-001608
A10	3302470-067	21330-001608	KEYPAD ASSY-MAIN CONTROL, ITALIAN	
A10	3302470-084	21330-001454	KEYPAD ASSY-MAIN CONTROL, NIBP, FRENCH	Replaced with 21330-001609
A10	3302470-086	21330-001609	KEYPAD ASSY-MAIN CONTROL, NIBP, FRENCH	
A10	3302470-085	21330-001455	KEYPAD ASSY-MAIN CONTROL, FRENCH	Replaced with 21330-001610
A10	3302470-087	21330-001610	KEYPAD ASSY-MAIN CONTROL, FRENCH	
A10	3302470-104	21330-001456	KEYPAD ASSY-MAIN CONTROL, NIBP, DUTCH	Replaced with 21330-001611
A10	3302470-106	21330-001611	KEYPAD ASSY-MAIN CONTROL, NIBP, DUTCH	
A10	3302470-105	21330-001457	KEYPAD ASSY-MAIN CONTROL, DUTCH	Replaced with 21330-001612
A10	3302470-107	21330-001612	KEYPAD ASSY-MAIN CONTROL, DUTCH	
A10	3302470-124	21330-001458	KEYPAD ASSY-MAIN CONTROL, NIBP, SPANISH	Replaced with 21330-001613
A10	3302470-126	21330-001613	KEYPAD ASSY-MAIN CONTROL, NIBP, SPANISH	
A10	3302470-125	21330-001459	KEYPAD ASSY-MAIN CONTROL, SPANISH	Replaced with 21330-001614
A10	3302470-127	21330-001614	KEYPAD ASSY-MAIN CONTROL, SPANISH	
A10	3302470-144	21330-001460	KEYPAD ASSY-MAIN CONTROL, NIBP, PORTUGUESE	Replaced with 21330-001615



Main Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A10	3302470-146	21330-001615	KEYPAD ASSY-MAIN CONTROL, NIBP, PORTUGUESE	
A10	3302470-145	21330-001461	KEYPAD ASSY-MAIN CONTROL, PORTUGUESE	Replaced with 21330-001616
A10	3302470-147	21330-001616	KEYPAD ASSY-MAIN CONTROL, PORTUGUESE	
A10	3302470-164	21330-001462	KEYPAD ASSY-MAIN CONTROL, NIBP, BRAZILIAN	Replaced with 21330-001617
A10	3302470-166	21330-001617	KEYPAD ASSY-MAIN CONTROL, NIBP, BRAZILIAN	
A10	3302470-165	21330-001463	KEYPAD ASSY-MAIN CONTROL, BRAZILIAN	Replaced with 21330-001618
A10	3302470-167	21330-001618	KEYPAD ASSY-MAIN CONTROL, BRAZILIAN	
A10	3302470-184	21330-001464	KEYPAD ASSY-MAIN CONTROL, NIBP, SWEDISH	Replaced with 21330-001619
A10	3302470-186	21330-001619	KEYPAD ASSY-MAIN CONTROL, NIBP, SWEDISH	
A10	3302470-185	21330-001465	KEYPAD ASSY-MAIN CONTROL, SWEDISH	Replaced with 21330-001620
A10	3302470-187	21330-001620	KEYPAD ASSY-MAIN CONTROL, SWEDISH	
A10	3302470-204	21330-001466	KEYPAD ASSY-MAIN CONTROL, NIBP, DANISH	Replaced with 21330-001621
A10	3302470-206	21330-001621	KEYPAD ASSY-MAIN CONTROL, NIBP, DANISH	
A10	3302470-205	21330-001469	KEYPAD ASSY-MAIN CONTROL, DANISH	Replaced with 21330-001622
A10	3302470-207	21330-001622	KEYPAD ASSY-MAIN CONTROL, DANISH	
A10	3302470-224	21330-001467	KEYPAD ASSY-MAIN CONTROL, NIBP, FINNISH	Replaced with 21330-001623
A10	3302470-226	21330-001623	KEYPAD ASSY-MAIN CONTROL, NIBP, FINNISH	
A10	3302470-225	21330-001468	KEYPAD ASSY-MAIN CONTROL, FINNISH	Replaced with 21330-001624



Main Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A10	3302470-227	21330-001624	KEYPAD ASSY-MAIN CONTROL, FINNISH	
A10	3302470-244	21330-001474	KEYPAD ASSY-MAIN CONTROL, NIBP, NORWEGIAN	Replaced with 21330-001625
A10	3302470-246	21330-001625	KEYPAD ASSY-MAIN CONTROL, NIBP, NORWEGIAN	
A10	3302470-245	21330-001470	KEYPAD ASSY-MAIN CONTROL, NORWEGIAN	Replaced with 21330-001626
A10	3302470-247	21330-001626	KEYPAD ASSY-MAIN CONTROL, NORWEGIAN	
A10	3302470-264	21330-001471	KEYPAD ASSY-MAIN CONTROL, NIBP, POLISH	Replaced with 21330-001627
A10	3302470-266	21330-001627	KEYPAD ASSY-MAIN CONTROL, NIBP, POLISH	
A10	3302470-265	21330-001472	KEYPAD ASSY-MAIN CONTROL, POLISH	Replaced with 21330-001628
A10	3302470-267	21330-001628	KEYPAD ASSY-MAIN CONTROL, POLISH	
A10	3302470-284	21330-001473	KEYPAD ASSY-MAIN CONTROL, NIBP, HUNGARIAN	Replaced with 21330-001629
A10	3302470-286	21330-001629	KEYPAD ASSY-MAIN CONTROL, NIBP, HUNGARIAN	
A10	3302470-285	21330-001476	KEYPAD ASSY-MAIN CONTROL, HUNGARIAN	Replaced with 21330-001630
A10	3302470-287	21330-001630	KEYPAD ASSY-MAIN CONTROL, HUNGARIAN	
A10	3302470-304	21330-001475	KEYPAD ASSY-MAIN CONTROL, NIBP, CZECH	Replaced with 21330-001631
A10	3302470-306	21330-001631	KEYPAD ASSY-MAIN CONTROL, NIBP, CZECH	
A10	3302470-305	21330-001477	KEYPAD ASSY-MAIN CONTROL, CZECH	Replaced with 21330-001632
A10	3302470-307	21330-001632	KEYPAD ASSY-MAIN CONTROL, CZECH	
A10	3302470-324	21330-001478	KEYPAD ASSY-MAIN CONTROL, NIBP, RUSSIAN	Replaced with 21330-001633



Main Keypad Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
A10	3302470-326	21330-001633	KEYPAD ASSY-MAIN CONTROL, NIBP, RUSSIAN	
A10	3302470-325	21330-001479	KEYPAD ASSY-MAIN CONTROL, RUSSIAN	Replaced with 21330-001634
A10	3302470-327	21330-001634	KEYPAD ASSY-MAIN CONTROL, RUSSIAN	
A10	3302470-384	21330-001480	KEYPAD ASSY-MAIN CONTROL, NIBP, KOREAN	Replaced with 21330-001635
A10	3302470-386	21330-001635	KEYPAD ASSY-MAIN CONTROL, NIBP, KOREAN	
A10	3302470-385	21330-001481	KEYPAD ASSY-MAIN CONTROL, KOREAN	Replaced with 21330-001636
A10	3302470-387	21330-001636	KEYPAD ASSY-MAIN CONTROL, KOREAN	
A10	3302470-402	21330-001506	KEYPAD ASSY-MAIN CONTROL, NIBP, JAPANESE	Replaced with 21330-001637
A10	3302470-404	21330-001637	KEYPAD ASSY-MAIN CONTROL, NIBP, JAPANESE	
A10	3302470-403	21330-001507	KEYPAD ASSY-MAIN CONTROL, JAPANESE	Replaced with 21330-001638
A10	3302470-405	21330-001638	KEYPAD ASSY-MAIN CONTROL, JAPANESE	

LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-000	21501-002408	LABEL SET-ECG,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-001	21501-002409	LABEL SET-ECG,SPO2,ENGLISH	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-002	21501-002410	LABEL SET-ECG,NIBP,SPO2,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-003	21501-002411	LABEL SET-ECG,SPO2,CO2,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-004	21501-002412	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-005	21501-002413	LABEL SET-ECG,SPO2,NIBP,IP,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-006	21501-002414	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-007	21501-002415	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ENGLISH	Replaced with L08 3 rd edition labels
L07	3305642-040	21501-002475	LABEL SET-ECG,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-041	21501-002476	LABEL SET-ECG,SPO2,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-042	21501-002477	LABEL SET-ECG,NIBP,SPO2,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-043	21501-002478	LABEL SET-ECG,SPO2,CO2,GERMAN	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-044	21501-002479	LABEL SET-ECG,SPO2,NIBP,CO2,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-045	21501-002545	LABEL SET-ECG,SPO2,NIBP,IP,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-046	21501-002546	LABEL SET-ECG,SPO2,NIBP,CO2,IP,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-047	21501-002547	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,GERMAN	Replaced with L08 3 rd edition labels
L07	3305642-060	21501-002548	LABEL SET-ECG,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-061	21501-002549	LABEL SET-ECG,SPO2,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-062	21501-002550	LABEL SET-ECG,NIBP,SPO2,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-063	21501-002551	LABEL SET-ECG,SPO2,CO2,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-064	21501-002552	LABEL SET-ECG,SPO2,NIBP,CO2,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-065	21501-002553	LABEL SET-ECG,SPO2,NIBP,IP,ITALIAN	Replaced with L08 3 rd edition labels

**LIFEPAK 15 Label Set Languages Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-066	21501-002554	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-067	21501-002555	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ITALIAN	Replaced with L08 3 rd edition labels
L07	3305642-080	21501-002480	LABEL SET-ECG,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-081	21501-002481	LABEL SET-ECG,SPO2,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-082	21501-002482	LABEL SET-ECG,NIBP,SPO2,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-083	21501-002483	LABEL SET-ECG,SPO2,CO2,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-084	21501-002484	LABEL SET-ECG,SPO2,NIBP,CO2,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-085	21501-002485	LABEL SET-ECG,SPO2,NIBP,IP,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-086	21501-002486	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FRENCH	Replaced with L08 3 rd edition labels
L07	3305642-087	21501-002487	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FRENCH	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-100	21501-002488	LABEL SET-ECG,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-101	21501-002489	LABEL SET-ECG,SPO2,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-102	21501-002527	LABEL SET-ECG,NIBP,SPO2,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-103	21501-002528	LABEL SET-ECG,SPO2,CO2,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-104	21501-002529	LABEL SET-ECG,SPO2,NIBP,CO2,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-105	21501-002530	LABEL SET-ECG,SPO2,NIBP,IP,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-106	21501-002536	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-107	21501-002500	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DUTCH	Replaced with L08 3 rd edition labels
L07	3305642-120	21501-002531	LABEL SET-ECG,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-121	21501-002532	LABEL SET-ECG,SPO2,SPANISH	Replaced with L08 3 rd edition labels

**LIFEPAK 15 Label Set Languages Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-122	21501-002533	LABEL SET-ECG,NIBP,SPO2,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-123	21501-002534	LABEL SET-ECG,SPO2,CO2,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-124	21501-002535	LABEL SET-ECG,SPO2,NIBP,CO2,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-125	21501-002491	LABEL SET-ECG,SPO2,NIBP,IP,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-126	21501-002493	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-127	21501-002494	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SPANISH	Replaced with L08 3 rd edition labels
L07	3305642-140	21501-002501	LABEL SET-ECG,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-141	21501-002502	LABEL SET-ECG,SPO2,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-142	21501-002503	LABEL SET-ECG,NIBP,SPO2,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-143	21501-002504	LABEL SET-ECG,SPO2,CO2,PORTUGESE	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-144	21501-002505	LABEL SET-ECG,SPO2,NIBP,CO2,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-145	21501-002506	LABEL SET-ECG,SPO2,NIBP,IP,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-146	21501-002507	LABEL SET-ECG,SPO2,NIBP,CO2,IP,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-147	21501-002508	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,PORTUGESE	Replaced with L08 3 rd edition labels
L07	3305642-160	21501-002509	LABEL SET-ECG,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-161	21501-002557	LABEL SET-ECG,SPO2,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-162	21501-002558	LABEL SET-ECG,NIBP,SPO2,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-163	21501-002559	LABEL SET-ECG,SPO2,CO2,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-164	21501-002560	LABEL SET-ECG,SPO2,NIBP,CO2,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-165	21501-002561	LABEL SET-ECG,SPO2,NIBP,IP,BRAZILIAN	Replaced with L08 3 rd edition labels

**LIFEPAK 15 Label Set Languages Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-166	21501-002562	LABEL SET-ECG,SPO2,NIBP,CO2,IP,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-167	21501-002556	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,BRAZILIAN	Replaced with L08 3 rd edition labels
L07	3305642-180	21501-002510	LABEL SET-ECG,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-181	21501-002511	LABEL SET-ECG,SPO2,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-182	21501-002512	LABEL SET-ECG,NIBP,SPO2,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-183	21501-002513	LABEL SET-ECG,SPO2,CO2,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-184	21501-002490	LABEL SET-ECG,SPO2,NIBP,CO2,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-185	21501-002492	LABEL SET-ECG,SPO2,NIBP,IP,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-186	21501-002495	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SWEDISH	Replaced with L08 3 rd edition labels
L07	3305642-187	21501-002496	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SWEDISH	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-200	21501-002497	LABEL SET-ECG,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-201	21501-002498	LABEL SET-ECG,SPO2,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-202	21501-002499	LABEL SET-ECG,NIBP,SPO2,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-203	21501-002514	LABEL SET-ECG,SPO2,CO2,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-204	21501-002515	LABEL SET-ECG,SPO2,NIBP,CO2,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-205	21501-002516	LABEL SET-ECG,SPO2,NIBP,IP,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-206	21501-002517	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-207	21501-002518	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DANISH	Replaced with L08 3 rd edition labels
L07	3305642-220	21501-002519	LABEL SET-ECG,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-221	21501-002520	LABEL SET-ECG,SPO2,FINNISH	Replaced with L08 3 rd edition labels

**LIFEPAK 15 Label Set Languages Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-222	21501-002521	LABEL SET-ECG,NIBP,SPO2,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-223	21501-002522	LABEL SET-ECG,SPO2,CO2,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-224	21501-002523	LABEL SET-ECG,SPO2,NIBP,CO2,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-225	21501-002524	LABEL SET-ECG,SPO2,NIBP,IP,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-226	21501-002525	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-227	21501-002526	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FINNISH	Replaced with L08 3 rd edition labels
L07	3305642-240	21501-002537	LABEL SET-ECG,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-241	21501-002538	LABEL SET-ECG,SPO2,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-242	21501-002540	LABEL SET-ECG,NIBP,SPO2,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-243	21501-002539	LABEL SET-ECG,SPO2,CO2,NORWEGIAN	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-244	21501-002541	LABEL SET-ECG,SPO2,NIBP,CO2,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-245	21501-002542	LABEL SET-ECG,SPO2,NIBP,IP,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-246	21501-002543	LABEL SET-ECG,SPO2,NIBP,CO2,IP,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-247	21501-002544	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,NORWEGIAN	Replaced with L08 3 rd edition labels
L07	3305642-260	21501-002563	LABEL SET-ECG,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-261	21501-002564	LABEL SET-ECG,SPO2,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-262	21501-002565	LABEL SET-ECG,NIBP,SPO2,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-263	21501-002566	LABEL SET-ECG,SPO2,CO2,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-264	21501-002567	LABEL SET-ECG,SPO2,NIBP,CO2,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-265	21501-002568	LABEL SET-ECG,SPO2,NIBP,IP,POLISH	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-266	21501-002569	LABEL SET-ECG,SPO2,NIBP,CO2,IP,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-267	21501-002570	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,POLISH	Replaced with L08 3 rd edition labels
L07	3305642-280	21501-002571	LABEL SET-ECG,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-281	21501-002572	LABEL SET-ECG,SPO2,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-282	21501-002573	LABEL SET-ECG,NIBP,SPO2,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-283	21501-002574	LABEL SET-ECG,SPO2,CO2,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-284	21501-002575	LABEL SET-ECG,SPO2,NIBP,CO2,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-285	21501-002576	LABEL SET-ECG,SPO2,NIBP,IP,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-286	21501-002577	LABEL SET-ECG,SPO2,NIBP,CO2,IP,HUNGARIAN	Replaced with L08 3 rd edition labels
L07	3305642-287	21501-002578	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,HUNGARIAN	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-300	21501-002579	LABEL SET-ECG,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-301	21501-002580	LABEL SET-ECG,SPO2,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-302	21501-002581	LABEL SET-ECG,NIBP,SPO2,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-303	21501-002582	LABEL SET-ECG,SPO2,CO2,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-304	21501-002583	LABEL SET-ECG,SPO2,NIBP,CO2,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-305	21501-002584	LABEL SET-ECG,SPO2,NIBP,IP,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-306	21501-002585	LABEL SET-ECG,SPO2,NIBP,CO2,IP,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-307	21501-002586	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,CZECH	Replaced with L08 3 rd edition labels
L07	3305642-320	21501-002587	LABEL SET-ECG,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-321	21501-002588	LABEL SET-ECG,SPO2,RUSSIAN	Replaced with L08 3 rd edition labels

**LIFEPAK 15 Label Set Languages Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-322	21501-002589	LABEL SET-ECG,NIBP,SPO2,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-323	21501-002590	LABEL SET-ECG,SPO2,CO2,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-324	21501-002591	LABEL SET-ECG,SPO2,NIBP,CO2,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-325	21501-002592	LABEL SET-ECG,SPO2,NIBP,IP,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-326	21501-002593	LABEL SET-ECG,SPO2,NIBP,CO2,IP,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-327	21501-002594	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,RUSSIAN	Replaced with L08 3 rd edition labels
L07	3305642-380	21501-002595	LABEL SET-ECG,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-381	21501-002596	LABEL SET-ECG,SPO2,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-382	21501-002597	LABEL SET-ECG,NIBP,SPO2,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-383	21501-002598	LABEL SET-ECG,SPO2,CO2,KOREAN	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-384	21501-002599	LABEL SET-ECG,SPO2,NIBP,CO2,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-385	21501-002600	LABEL SET-ECG,SPO2,NIBP,IP,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-386	21501-002601	LABEL SET-ECG,SPO2,NIBP,CO2,IP,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-387	21501-002602	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,KOREAN	Replaced with L08 3 rd edition labels
L07	3305642-400	21501-002649	LABEL SET-ECG,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-401	21501-002650	LABEL SET-ECG,SPO2,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-402	21501-002651	LABEL SET-ECG,NIBP, SPO2 ,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-403	21501-002652	LABEL SET-ECG,SPO2,CO2,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-404	21501-002653	LABEL SET-ECG,SPO2,NIBP,CO2,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-405	21501-002654	LABEL SET-ECG,SPO2,NIBP,IP,JAPANESE	Replaced with L08 3 rd edition labels



LIFEPAK 15 Label Set Languages Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L07	3305642-406	21501-002655	LABEL SET-ECG,SPO2,NIBP,CO2,IP,JAPANESE	Replaced with L08 3 rd edition labels
L07	3305642-407	21501-002656	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,JAPANESE	Replaced with L08 3 rd edition labels

3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-000	21501-002700	LABEL SET-ECG,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-001	21501-002701	LABEL SET-ECG,SPO2,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-002	21501-002702	LABEL SET-ECG,NIBP,SPO2,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-003	21501-002703	LABEL SET-ECG,SPO2,CO2,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-004	21501-002704	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-005	21501-002705	LABEL SET-ECG,SPO2,NIBP,IP,ENGLISH, 3RD EDITION	Bluetooth
L08	3349906-006	21501-002706	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ENGLISH, 3RD ED	Bluetooth
L08	3349906-007	21501-002707	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ENGLISH, 3RD ED	Bluetooth
L08	3349906-040	21501-002708	LABEL SET-ECG,GERMAN, 3RD EDITION	Bluetooth
L08	3349906-041	21501-002709	LABEL SET-ECG,SPO2,GERMAN, 3RD EDITION	Bluetooth
L08	3349906-042	21501-002710	LABEL SET-ECG,NIBP,SPO2,GERMAN, 3RD EDITION	Bluetooth
L08	3349906-043	21501-002711	LABEL SET-ECG,SPO2,CO2,GERMAN, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-044	21501-002712	LABEL SET-ECG,SPO2,NIBP,CO2,GERMAN, 3RD EDITION	Bluetooth
L08	3349906-045	21501-002713	LABEL SET-ECG,SPO2,NIBP,IP,GERMAN, 3RD EDITION	Bluetooth
L08	3349906-046	21501-002714	LABEL SET-ECG,SPO2,NIBP,CO2,IP,GERMAN, 3RD ED	Bluetooth
L08	3349906-047	21501-002715	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,GERMAN, 3RD ED	Bluetooth
L08	3349906-060	21501-002716	LABEL SET-ECG,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-061	21501-002717	LABEL SET-ECG,SPO2,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-062	21501-002718	LABEL SET-ECG,NIBP,SPO2,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-063	21501-002719	LABEL SET-ECG,SPO2,CO2,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-064	21501-002720	LABEL SET-ECG,SPO2,NIBP,CO2,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-065	21501-002721	LABEL SET-ECG,SPO2,NIBP,IP,ITALIAN, 3RD EDITION	Bluetooth
L08	3349906-066	21501-002722	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ITALIAN, 3RD ED	Bluetooth
L08	3349906-067	21501-002723	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ITALIAN, 3RD ED	Bluetooth
L08	3349906-080	21501-002724	LABEL SET-ECG,FRENCH, 3RD EDITION	Bluetooth
L08	3349906-081	21501-002725	LABEL SET-ECG,SPO2,FRENCH, 3RD EDITION	Bluetooth
L08	3349906-082	21501-002726	LABEL SET-ECG,NIBP,SPO2,FRENCH, 3RD EDITION	Bluetooth
L08	3349906-083	21501-002727	LABEL SET-ECG,SPO2,CO2,FRENCH, 3RD EDITION	Bluetooth
L08	3349906-084	21501-002728	LABEL SET-ECG,SPO2,NIBP,CO2,FRENCH, 3RD EDITION	Bluetooth
L08	3349906-085	21501-002729	LABEL SET-ECG,SPO2,NIBP,IP,FRENCH, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-086	21501-002730	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FRENCH, 3RD ED	Bluetooth
L08	3349906-087	21501-002731	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FRENCH, 3RD ED	Bluetooth
L08	3349906-100	21501-002732	LABEL SET-ECG,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-101	21501-002733	LABEL SET-ECG,SPO2,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-102	21501-002734	LABEL SET-ECG,NIBP,SPO2,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-103	21501-002735	LABEL SET-ECG,SPO2,CO2,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-104	21501-002736	LABEL SET-ECG,SPO2,NIBP,CO2,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-105	21501-002737	LABEL SET-ECG,SPO2,NIBP,IP,DUTCH, 3RD EDITION	Bluetooth
L08	3349906-106	21501-002738	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DUTCH, 3RD ED	Bluetooth
L08	3349906-107	21501-002739	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DUTCH, 3RD ED	Bluetooth
L08	3349906-120	21501-002740	LABEL SET-ECG,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-121	21501-002741	LABEL SET-ECG,SPO2,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-122	21501-002742	LABEL SET-ECG,NIBP,SPO2,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-123	21501-002743	LABEL SET-ECG,SPO2,CO2,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-124	21501-002744	LABEL SET-ECG,SPO2,NIBP,CO2,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-125	21501-002745	LABEL SET-ECG,SPO2,NIBP,IP,SPANISH, 3RD EDITION	Bluetooth
L08	3349906-126	21501-002746	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SPANISH, 3RD ED	Bluetooth
L08	3349906-127	21501-002747	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SPANISH, 3RD ED	Bluetooth
L08	3349906-140	21501-002748	LABEL SET-ECG,PORTUGUESE, 3RD EDITION	Bluetooth
L08	3349906-141	21501-002749	LABEL SET-ECG,SPO2,PORTUGUESE, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-142	21501-002750	LABEL SET-ECG,NIBP,SPO2,PORTUGUESE, 3RD EDITION	Bluetooth
L08	3349906-143	21501-002751	LABEL SET-ECG,SPO2,CO2,PORTUGUESE, 3RD EDITION	Bluetooth
L08	3349906-144	21501-002752	LABEL SET-ECG,SPO2,NIBP,CO2,PORTUGUESE, 3RD ED	Bluetooth
L08	3349906-145	21501-002753	LABEL SET-ECG,SPO2,NIBP,IP,PORTUGUESE, 3RD ED	Bluetooth
L08	3349906-146	21501-002754	LABEL SET-ECG,SPO2,NIBP,CO2,IP,PORTUGUESE, 3RDED	Bluetooth
L08	3349906-147	21501-002755	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,PORTUGUESE, 3ED	Bluetooth
L08	3349906-160	21501-002756	LABEL SET-ECG,BRAZILIAN, 3RD EDITION	Bluetooth
L08	3349906-161	21501-002757	LABEL SET-ECG,SPO2,BRAZILIAN, 3RD EDITION	Bluetooth
L08	3349906-162	21501-002758	LABEL SET-ECG,NIBP,SPO2,BRAZILIAN, 3RD EDITION	Bluetooth
L08	3349906-163	21501-002759	LABEL SET-ECG,SPO2,CO2,BRAZILIAN, 3RD EDITION	Bluetooth
L08	3349906-164	21501-002760	LABEL SET-ECG,SPO2,NIBP,CO2,BRAZILIAN, 3RD ED	Bluetooth
L08	3349906-165	21501-002761	LABEL SET-ECG,SPO2,NIBP,IP,BRAZILIAN, 3RD ED	Bluetooth
L08	3349906-166	21501-002762	LABEL SET-ECG,SPO2,NIBP,CO2,IP,BRAZILIAN, 3RD ED	Bluetooth
L08	3349906-167	21501-002763	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,BRAZILIAN, 3ED	Bluetooth
L08	3349906-180	21501-002764	LABEL SET-ECG,SWEDISH, 3RD EDITION	Bluetooth
L08	3349906-181	21501-002765	LABEL SET-ECG,SPO2,SWEDISH, 3RD EDITION	Bluetooth
L08	3349906-182	21501-002766	LABEL SET-ECG,NIBP,SPO2,SWEDISH, 3RD EDITION	Bluetooth
L08	3349906-183	21501-002767	LABEL SET-ECG,SPO2,CO2,SWEDISH, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-184	21501-002768	LABEL SET-ECG,SPO2,NIBP,CO2,SWEDISH, 3RD EDITION	Bluetooth
L08	3349906-185	21501-002769	LABEL SET-ECG,SPO2,NIBP,IP,SWEDISH, 3RD EDITION	Bluetooth
L08	3349906-186	21501-002770	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SWEDISH, 3RD ED	Bluetooth
L08	3349906-187	21501-002771	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SWEDISH, 3RD ED	Bluetooth
L08	3349906-200	21501-002772	LABEL SET-ECG,DANISH, 3RD EDITION	Bluetooth
L08	3349906-201	21501-002773	LABEL SET-ECG,SPO2,DANISH, 3RD EDITION	Bluetooth
L08	3349906-202	21501-002774	LABEL SET-ECG,NIBP,SPO2,DANISH, 3RD EDITION	Bluetooth
L08	3349906-203	21501-002775	LABEL SET-ECG,SPO2,CO2,DANISH, 3RD EDITION	Bluetooth
L08	3349906-204	21501-002776	LABEL SET-ECG,SPO2,NIBP,CO2,DANISH, 3RD EDITION	Bluetooth
L08	3349906-205	21501-002777	LABEL SET-ECG,SPO2,NIBP,IP,DANISH, 3RD EDITION	Bluetooth
L08	3349906-206	21501-002778	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DANISH, 3RD ED	Bluetooth
L08	3349906-207	21501-002779	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DANISH, 3RD ED	Bluetooth
L08	3349906-220	21501-002780	LABEL SET-ECG,FINNISH, 3RD EDITION	Bluetooth
L08	3349906-221	21501-002781	LABEL SET-ECG,SPO2,FINNISH, 3RD EDITION	Bluetooth
L08	3349906-222	21501-002782	LABEL SET-ECG,NIBP,SPO2,FINNISH, 3RD EDITION	Bluetooth
L08	3349906-223	21501-002783	LABEL SET-ECG,SPO2,CO2,FINNISH, 3RD EDITION	Bluetooth
L08	3349906-224	21501-002784	LABEL SET-ECG,SPO2,NIBP,CO2,FINNISH, 3RD EDITION	Bluetooth
L08	3349906-225	21501-002785	LABEL SET-ECG,SPO2,NIBP,IP,FINNISH, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-226	21501-002786	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FINNISH, 3RD ED	Bluetooth
L08	3349906-227	21501-002787	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FINNISH, 3RD ED	Bluetooth
L08	3349906-240	21501-002788	LABEL SET-ECG,NORWEGIAN, 3RD EDITION	Bluetooth
L08	3349906-241	21501-002789	LABEL SET-ECG,SPO2,NORWEGIAN, 3RD EDITION	Bluetooth
L08	3349906-242	21501-002790	LABEL SET-ECG,NIBP,SPO2,NORWEGIAN, 3RD EDITION	Bluetooth
L08	3349906-243	21501-002791	LABEL SET-ECG,SPO2,CO2,NORWEGIAN, 3RD EDITION	Bluetooth
L08	3349906-244	21501-002792	LABEL SET-ECG,SPO2,NIBP,CO2,NORWEGIAN, 3RD ED	Bluetooth
L08	3349906-245	21501-002793	LABEL SET-ECG,SPO2,NIBP,IP,NORWEGIAN, 3RD ED	Bluetooth
L08	3349906-246	21501-002794	LABEL SET-ECG,SPO2,NIBP,CO2,IP,NORWEGIAN, 3RD ED	Bluetooth
L08	3349906-247	21501-002795	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,NORWEGIAN, 3ED	Bluetooth
L08	3349906-260	21501-002796	LABEL SET-ECG,POLISH, 3RD EDITION	Bluetooth
L08	3349906-261	21501-002797	LABEL SET-ECG,SPO2,POLISH, 3RD EDITION	Bluetooth
L08	3349906-262	21501-002798	LABEL SET-ECG,NIBP,SPO2,POLISH, 3RD EDITION	Bluetooth
L08	3349906-263	21501-002799	LABEL SET-ECG,SPO2,CO2,POLISH, 3RD EDITION	Bluetooth
L08	3349906-264	21501-002800	LABEL SET-ECG,SPO2,NIBP,CO2,POLISH, 3RD EDITION	Bluetooth
L08	3349906-265	21501-002801	LABEL SET-ECG,SPO2,NIBP,IP,POLISH, 3RD EDITION	Bluetooth
L08	3349906-266	21501-002802	LABEL SET-ECG,SPO2,NIBP,CO2,IP,POLISH, 3RD ED	Bluetooth
L08	3349906-267	21501-002803	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,POLISH, 3RD ED	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-280	21501-002804	LABEL SET-ECG,HUNGARIAN, 3RD EDITION	Bluetooth
L08	3349906-281	21501-002805	LABEL SET-ECG,SPO2,HUNGARIAN, 3RD EDITION	Bluetooth
L08	3349906-282	21501-002806	LABEL SET-ECG,NIBP,SPO2,HUNGARIAN, 3RD EDITION	Bluetooth
L08	3349906-283	21501-002807	LABEL SET-ECG,SPO2,CO2,HUNGARIAN, 3RD EDITION	Bluetooth
L08	3349906-284	21501-002808	LABEL SET-ECG,SPO2,NIBP,CO2,HUNGARIAN, 3RD ED	Bluetooth
L08	3349906-285	21501-002809	LABEL SET-ECG,SPO2,NIBP,IP,HUNGARIAN, 3RD ED	Bluetooth
L08	3349906-286	21501-002810	LABEL SET-ECG,SPO2,NIBP,CO2,IP,HUNGARIAN, 3RD ED	Bluetooth
L08	3349906-287	21501-002811	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,HUNGARIAN,3RDED	Bluetooth
L08	3349906-300	21501-002812	LABEL SET-ECG,CZECH, 3RD EDITION	Bluetooth
L08	3349906-301	21501-002813	LABEL SET-ECG,SPO2,CZECH , 3RD EDITION	Bluetooth
L08	3349906-302	21501-002814	LABEL SET-ECG,NIBP,SPO2,CZECH, 3RD EDITION	Bluetooth
L08	3349906-303	21501-002815	LABEL SET-ECG,SPO2,CO2,CZECH, 3RD EDITION	Bluetooth
L08	3349906-304	21501-002816	LABEL SET-ECG,SPO2,NIBP,CO2,CZECH, 3RD EDITION	Bluetooth
L08	3349906-305	21501-002817	LABEL SET-ECG,SPO2,NIBP,IP,CZECH, 3RD EDITION	Bluetooth
L08	3349906-306	21501-002818	LABEL SET-ECG,SPO2,NIBP,CO2,IP,CZECH, 3RD ED	Bluetooth
L08	3349906-307	21501-002819	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,CZECH, 3RD ED	Bluetooth
L08	3349906-320	21501-002820	LABEL SET-ECG,RUSSIAN, 3RD EDITION	Bluetooth
L08	3349906-321	21501-002821	LABEL SET-ECG,SPO2,RUSSIAN, 3RD EDITION	Bluetooth
L08	3349906-322	21501-002822	LABEL SET-ECG,NIBP,SPO2,RUSSIAN, 3RD EDITION	Bluetooth
L08	3349906-323	21501-002823	LABEL SET-ECG,SPO2,CO2,RUSSIAN, 3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-324	21501-002824	LABEL SET-ECG,SPO2,NIBP,CO2,RUSSIAN, 3RD EDITION	Bluetooth
L08	3349906-325	21501-002825	LABEL SET-ECG,SPO2,NIBP, IP,RUSSIAN, 3RD EDITION	Bluetooth
L08	3349906-326	21501-002826	LABEL SET-ECG,SPO2,NIBP,CO2,IP,RUSSIAN, 3RD ED	Bluetooth
L08	3349906-327	21501-002827	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,RUSSIAN, 3RD ED	Bluetooth
L08	3349906-380	21501-002828	LABEL SET-ECG,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-381	21501-002829	LABEL SET-ECG,SPO2,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-382	21501-002830	LABEL SET-ECG,NIBP,SPO2,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-383	21501-002831	LABEL SET-ECG,SPO2,CO2,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-384	21501-002832	LABEL SET-ECG,SPO2,NIBP,CO2,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-385	21501-002833	LABEL SET-ECG,SPO2,NIBP,IP,KOREAN, 3RD EDITION	Bluetooth
L08	3349906-386	21501-002834	LABEL SET-ECG,SPO2,NIBP,CO2,IP,KOREAN, 3RD ED	Bluetooth
L08	3349906-387	21501-002835	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,KOREAN, 3RD ED	Bluetooth
L08	3349906-400	21501-002836	LABEL SET-ECG, JAPANESE, 3RD EDITION	Bluetooth
L08	3349906-401	21501-002837	LABEL SET-ECG,SPO2,JAPANESE, 3RD EDITION	Bluetooth
L08	3349906-402	21501-002838	LABEL SET-ECG,NIBP, SPO2, JAPANESE, 3RD EDITION	Bluetooth
L08	3349906-403	21501-002839	LABEL SET-ECG,SPO2,CO2,JAPANESE, 3RD EDITION	Bluetooth
L08	3349906-404	21501-002840	LABEL SET-ECG,SPO2,NIBP,CO2,JAPANESE,3RD EDITION	Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3349906-405	21501-002841	LABEL SET-ECG,SPO2,NIBP,IP,JAPANESE, 3RD EDITION	Bluetooth
L08	3349906-406	21501-002842	LABEL SET-ECG,SPO2,NIBP,CO2,IP,JAPANESE, 3RD ED	Bluetooth
L08	3349906-407	21501-002843	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,JAPANESE,3RD ED	Bluetooth
L08	3348613-000	21501-003074	LABEL SET-ECG,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-001	21501-003075	LABEL SET-ECG,SPO2,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-002	21501-003076	LABEL SET-ECG,NIBP,SPO2,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-003	21501-003077	LABEL SET-ECG,SPO2,CO2,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-004	21501-003078	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-005	21501-003079	LABEL SET-ECG,SPO2,NIBP,IP,ENGLISH,3RD EDITION	Non-Bluetooth
L08	3348613-006	21501-003080	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ENGLISH,3RD ED	Non-Bluetooth
L08	3348613-007	21501-003081	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ENGLISH,3RD ED	Non-Bluetooth
L08	3348613-040	21501-003082	LABEL SET-ECG,GERMAN,3RD EDITION	Non-Bluetooth
L08	3348613-041	21501-003083	LABEL SET-ECG,SPO2,GERMAN,3RD EDITION	Non-Bluetooth
L08	3348613-042	21501-003084	LABEL SET-ECG,NIBP,SPO2,GERMAN,3RD EDITION	Non-Bluetooth
L08	3348613-043	21501-003085	LABEL SET-ECG,SPO2,CO2,GERMAN, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-044	21501-003086	LABEL SET-ECG,SPO2,NIBP,CO2,GERMAN, 3RD EDITION	Non-Bluetooth
L08	3348613-045	21501-003087	LABEL SET-ECG,SPO2,NIBP,IP,GERMAN, 3RD EDITION	Non-Bluetooth
L08	3348613-046	21501-003088	LABEL SET-ECG,SPO2,NIBP,CO2,IP,GERMAN, 3RD ED	Non-Bluetooth
L08	3348613-047	21501-003089	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,GERMAN, 3RD ED	Non-Bluetooth
L08	3348613-060	21501-003090	LABEL SET-ECG,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-061	21501-003091	LABEL SET-ECG,SPO2,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-062	21501-003092	LABEL SET-ECG,NIBP,SPO2,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-063	21501-003093	LABEL SET-ECG,SPO2,CO2,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-064	21501-003094	LABEL SET-ECG,SPO2,NIBP,CO2,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-065	21501-003095	LABEL SET-ECG,SPO2,NIBP,IP,ITALIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-066	21501-003096	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ITALIAN, 3RD ED	Non-Bluetooth
L08	3348613-067	21501-003097	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ITALIAN, 3RD ED	Non-Bluetooth
L08	3348613-080	21501-003098	LABEL SET-ECG,FRENCH, 3RD EDITION	Non-Bluetooth
L08	3348613-081	21501-003099	LABEL SET-ECG,SPO2,FRENCH, 3RD EDITION	Non-Bluetooth
L08	3348613-082	21501-003100	LABEL SET-ECG,NIBP,SPO2,FRENCH, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-083	21501-003101	LABEL SET-ECG,SPO2,CO2,FRENCH, 3RD EDITION	Non-Bluetooth
L08	3348613-084	21501-003102	LABEL SET-ECG,SPO2,NIBP,CO2,FRENCH, 3RD EDITION	Non-Bluetooth
L08	3348613-085	21501-003103	LABEL SET-ECG,SPO2,NIBP,IP,FRENCH, 3RD EDITION	Non-Bluetooth
L08	3348613-086	21501-003104	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FRENCH, 3RD ED	Non-Bluetooth
L08	3348613-087	21501-003105	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FRENCH, 3RD ED	Non-Bluetooth
L08	3348613-100	21501-003106	LABEL SET-ECG,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-101	21501-003107	LABEL SET-ECG,SPO2,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-102	21501-003108	LABEL SET-ECG,NIBP,SPO2,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-103	21501-003109	LABEL SET-ECG,SPO2,CO2,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-104	21501-003110	LABEL SET-ECG,SPO2,NIBP,CO2,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-105	21501-003111	LABEL SET-ECG,SPO2,NIBP,IP,DUTCH, 3RD EDITION	Non-Bluetooth
L08	3348613-106	21501-003112	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DUTCH, 3RD ED	Non-Bluetooth
L08	3348613-107	21501-003113	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DUTCH, 3RD ED	Non-Bluetooth
L08	3348613-120	21501-003114	LABEL SET-ECG,SPANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-121	21501-003115	LABEL SET-ECG,SPO2,SPANISH, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-122	21501-003116	LABEL SET-ECG,NIBP,SPO2,SPANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-123	21501-003117	LABEL SET-ECG,SPO2,CO2,SPANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-124	21501-003118	LABEL SET-ECG,SPO2,NIBP,CO2,SPANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-125	21501-003119	LABEL SET-ECG,SPO2,NIBP,IP,SPANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-126	21501-003120	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SPANISH, 3RD ED	Non-Bluetooth
L08	3348613-127	21501-003121	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SPANISH, 3RD ED	Non-Bluetooth
L08	3348613-140	21501-003122	LABEL SET-ECG,PORTUGUESE, 3RD EDITION	Non-Bluetooth
L08	3348613-141	21501-003123	LABEL SET-ECG,SPO2,PORTUGUESE, 3RD EDITION	Non-Bluetooth
L08	3348613-142	21501-003124	LABEL SET-ECG,NIBP,SPO2,PORTUGUESE, 3RD EDITION	Non-Bluetooth
L08	3348613-143	21501-003125	LABEL SET-ECG,SPO2,CO2,PORTUGUESE, 3RD EDITION	Non-Bluetooth
L08	3348613-144	21501-003126	LABEL SET-ECG,SPO2,NIBP,CO2,PORTUGUESE, 3RD ED	Non-Bluetooth
L08	3348613-145	21501-003127	LABEL SET-ECG,SPO2,NIBP,IP,PORTUGUESE, 3RD ED	Non-Bluetooth
L08	3348613-146	21501-003128	LABEL SET-ECG,SPO2,NIBP,CO2,IP,PORTUGUESE, 3RDED	Non-Bluetooth
L08	3348613-147	21501-003129	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,PORTUGUESE, 3ED	Non-Bluetooth
L08	3348613-160	21501-003130	LABEL SET-ECG,BRAZILIAN, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-161	21501-003131	LABEL SET-ECG,SPO2,BRAZILIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-162	21501-003132	LABEL SET-ECG,NIBP,SPO2,BRAZILIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-163	21501-003133	LABEL SET-ECG,SPO2,CO2,BRAZILIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-164	21501-003134	LABEL SET-ECG,SPO2,NIBP,CO2,BRAZILIAN, 3RD ED	Non-Bluetooth
L08	3348613-165	21501-003135	LABEL SET-ECG,SPO2,NIBP,IP,BRAZILIAN, 3RD ED	Non-Bluetooth
L08	3348613-166	21501-003136	LABEL SET-ECG,SPO2,NIBP,CO2,IP,BRAZILIAN, 3RD ED	Non-Bluetooth
L08	3348613-167	21501-003137	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,BRAZILIAN, 3ED	Non-Bluetooth
L08	3348613-180	21501-003138	LABEL SET-ECG,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-181	21501-003139	LABEL SET-ECG,SPO2,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-182	21501-003140	LABEL SET-ECG,NIBP,SPO2,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-183	21501-003141	LABEL SET-ECG,SPO2,CO2,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-184	21501-003142	LABEL SET-ECG,SPO2,NIBP,CO2,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-185	21501-003143	LABEL SET-ECG,SPO2,NIBP,IP,SWEDISH, 3RD EDITION	Non-Bluetooth
L08	3348613-186	21501-003144	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SWEDISH, 3RD ED	Non-Bluetooth
L08	3348613-187	21501-003145	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SWEDISH, 3RD ED	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-200	21501-003146	LABEL SET-ECG,DANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-201	21501-003147	LABEL SET-ECG,SPO2,DANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-202	21501-003148	LABEL SET-ECG,NIBP,SPO2,DANISH,3RD EDITION	Non-Bluetooth
L08	3348613-203	21501-003149	LABEL SET-ECG,SPO2,CO2,DANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-204	21501-003150	LABEL SET-ECG,SPO2,NIBP,CO2,DANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-205	21501-003151	LABEL SET-ECG,SPO2,NIBP,IP,DANISH, 3RD EDITION	Non-Bluetooth
L08	3348613-206	21501-003152	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DANISH, 3RD ED	Non-Bluetooth
L08	3348613-207	21501-003153	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DANISH, 3RD ED	Non-Bluetooth
L08	3348613-220	21501-003154	LABEL SET-ECG,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-221	21501-003155	LABEL SET-ECG,SPO2,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-222	21501-003156	LABEL SET-ECG,NIBP,SPO2,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-223	21501-003157	LABEL SET-ECG,SPO2,CO2,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-224	21501-003158	LABEL SET-ECG,SPO2,NIBP,CO2,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-225	21501-003159	LABEL SET-ECG,SPO2,NIBP,IP,FINNISH, 3RD EDITION	Non-Bluetooth
L08	3348613-226	21501-003160	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FINNISH, 3RD ED	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-227	21501-003161	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FINNISH, 3RD ED	Non-Bluetooth
L08	3348613-240	21501-003162	LABEL SET-ECG,NORWEGIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-241	21501-003163	LABEL SET-ECG,SPO2,NORWEGIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-242	21501-003164	LABEL SET-ECG,NIBP,SPO2,NORWEGIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-243	21501-003165	LABEL SET-ECG,SPO2,CO2,NORWEGIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-244	21501-003166	LABEL SET-ECG,SPO2,NIBP,CO2,NORWEGIAN, 3RD ED	Non-Bluetooth
L08	3348613-245	21501-003167	LABEL SET-ECG,SPO2,NIBP,IP,NORWEGIAN, 3RD ED	Non-Bluetooth
L08	3348613-246	21501-003168	LABEL SET-ECG,SPO2,NIBP,CO2,IP,NORWEGIAN, 3RD ED	Non-Bluetooth
L08	3348613-247	21501-003169	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,NORWEGIAN, 3ED	Non-Bluetooth
L08	3348613-260	21501-003170	LABEL SET-ECG,POLISH, 3RD EDITION	Non-Bluetooth
L08	3348613-261	21501-003171	LABEL SET-ECG,SPO2,POLISH, 3RD EDITION	Non-Bluetooth
L08	3348613-262	21501-003172	LABEL SET-ECG,NIBP,SPO2,POLISH, 3RD EDITION	Non-Bluetooth
L08	3348613-263	21501-003173	LABEL SET-ECG,SPO2,CO2,POLISH, 3RD EDITION	Non-Bluetooth
L08	3348613-264	21501-003174	LABEL SET-ECG,SPO2,NIBP,CO2,POLISH, 3RD EDITION	Non-Bluetooth
L08	3348613-265	21501-003175	LABEL SET-ECG,SPO2,NIBP,IP,POLISH, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-266	21501-003176	LABEL SET-ECG,SPO2,NIBP,CO2,IP,POLISH, 3RD ED	Non-Bluetooth
L08	3348613-267	21501-003177	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,POLISH, 3RD ED	Non-Bluetooth
L08	3348613-280	21501-003178	LABEL SET-ECG,HUNGARIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-281	21501-003179	LABEL SET-ECG,SPO2,HUNGARIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-282	21501-003180	LABEL SET-ECG,NIBP,SPO2,HUNGARIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-283	21501-003181	LABEL SET-ECG,SPO2,CO2,HUNGARIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-284	21501-003182	LABEL SET-ECG,SPO2,NIBP,CO2,HUNGARIAN, 3RD ED	Non-Bluetooth
L08	3348613-285	21501-003183	LABEL SET-ECG,SPO2,NIBP,IP,HUNGARIAN, 3RD ED	Non-Bluetooth
L08	3348613-286	21501-003184	LABEL SET-ECG,SPO2,NIBP,CO2,IP,HUNGARIAN, 3RD ED	Non-Bluetooth
L08	3348613-287	21501-003185	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,HUNGARIAN,3RDED	Non-Bluetooth
L08	3348613-300	21501-003186	LABEL SET-ECG,CZECH, 3RD EDITION	Non-Bluetooth
L08	3348613-301	21501-003187	LABEL SET-ECG,SPO2,CZECH, 3RD EDITION	Non-Bluetooth
L08	3348613-302	21501-003188	LABEL SET-ECG,NIBP,SPO2,CZECH, 3RD EDITION	Non-Bluetooth
L08	3348613-303	21501-003189	LABEL SET-ECG,SPO2,CO2,CZECH, 3RD EDITION	Non-Bluetooth
L08	3348613-304	21501-003190	LABEL SET-ECG,SPO2,NIBP,CO2,CZECH, 3RD EDITION	Non-Bluetooth

**3rd Edition Label Set Parts List**

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-305	21501-003191	LABEL SET-ECG,SPO2,NIBP,IP,CZECH, 3RD EDITION	Non-Bluetooth
L08	3348613-306	21501-003192	LABEL SET-ECG,SPO2,NIBP,CO2,IP,CZECH, 3RD ED	Non-Bluetooth
L08	3348613-307	21501-003193	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,CZECH, 3RD ED	Non-Bluetooth
L08	3348613-320	21501-003194	LABEL SET-ECG,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-321	21501-003195	LABEL SET-ECG,SPO2,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-322	21501-003196	LABEL SET-ECG,NIBP,SPO2,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-323	21501-003197	LABEL SET-ECG,SPO2,CO2,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-324	21501-003198	LABEL SET-ECG,SPO2,NIBP,CO2,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-325	21501-003199	LABEL SET-ECG,SPO2,NIBP,IP,RUSSIAN, 3RD EDITION	Non-Bluetooth
L08	3348613-326	21501-003200	LABEL SET-ECG,SPO2,NIBP,CO2,IP,RUSSIAN, 3RD ED	Non-Bluetooth
L08	3348613-327	21501-003201	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,RUSSIAN, 3RD ED	Non-Bluetooth
L08	3348613-380	21501-003202	LABEL SET-ECG,KOREAN, 3RD EDITION	Non-Bluetooth
L08	3348613-381	21501-003203	LABEL SET-ECG,SPO2,KOREAN, 3RD EDITION	Non-Bluetooth
L08	3348613-382	21501-003204	LABEL SET-ECG,NIBP,SPO2,KOREAN, 3RD EDITION	Non-Bluetooth
L08	3348613-383	21501-003205	LABEL SET-ECG,SPO2,CO2,KOREAN, 3RD EDITION	Non-Bluetooth



3rd Edition Label Set Parts List

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
L08	3348613-384	21501-003206	LABEL SET-ECG,SPO2,NIBP,CO2,KOREAN, 3RD EDITION	Non-Bluetooth
L08	3348613-385	21501-003207	LABEL SET-ECG,SPO2,NIBP,IP,KOREAN, 3RD EDITION	Non-Bluetooth
L08	3348613-386	21501-003208	LABEL SET-ECG,SPO2,NIBP,CO2,IP,KOREAN, 3RD ED	Non-Bluetooth
L08	3348613-387	21501-003209	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,KOREAN, 3RD ED	Non-Bluetooth
L08	3348613-400	21501-003210	LABEL SET-ECG, JAPANESE, 3RD EDITION	Non-Bluetooth
L08	3348613-401	21501-003211	LABEL SET-ECG,SPO2,JAPANESE, 3RD EDITION	Non-Bluetooth
L08	3348613-402	21501-003212	LABEL SET-ECG,NIBP, SPO2, JAPANESE, 3RD EDITION	Non-Bluetooth
L08	3348613-403	21501-003213	LABEL SET-ECG,SPO2,CO2,JAPANESE, 3RD EDITION	Non-Bluetooth
L08	3348613-404	21501-003214	LABEL SET-ECG,SPO2,NIBP,CO2,JAPANESE,3RD EDITION	Non-Bluetooth
L08	3348613-405	21501-003215	LABEL SET-ECG,SPO2,NIBP,IP,JAPANESE, 3RD EDITION	Non-Bluetooth
L08	3348613-406	21501-003216	LABEL SET-ECG,SPO2,NIBP,CO2,IP,JAPANESE, 3RD ED	Non-Bluetooth
L08	3348613-407	21501-003217	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,JAPANESE,3RD ED	Non-Bluetooth

5. SERVICE REPAIR KITS

LIFEPAK 15 Service Repair Kits include components that support a particular replacement activity. These repair kits are also listed in 3305431, LP15 Service Repair Kit. Reference the electronic PLM system Bill of Materials for part numbers and quantity for each assembly.



Service Repair Kits

Reference Designator	Part Number	Catalog Number	Description	Note (Optional)
K01	3305431-015	40577-000003	KIT - REPAIR, INTERNAL HARDWARE, LP15	
K02	3305431-016	40577-000004	KIT - REPAIR, EXTERNAL HARDWARE, LP15	
K03	3305431-017	40577-000005	KIT - REPAIR, PADDLE RETAINER, LP15	Replaced with 40577-000054
K03	3305431-539	40577-000054	KIT - REPAIR, PADDLE RETAINER, LP15	
K04	3305431-018	40577-000006	KIT - REPAIR, DISPLAY SHIELD, LP15	Replaced with 40577-000063
K04	3305431-548	40577-000063	KIT - REPAIR, DISPLAY SHIELD, BLK, LP15	
K05	3305431-019	40577-000007	KIT - REPAIR, GUARD AND FEET, LP15	
K06	3305431-020	40577-000008	KIT - REPAIR, HANDLE, LP15	
K07	3305431-022	40577-000009	KIT - REPAIR, NIBP CONNECTOR, LP15	Replaced with 40577-000055
K07	3305431-540	40577-000055	KIT - REPAIR, NIBP CONNECTOR, LP15	
K08	3305431-023	40577-000010	KIT - REPAIR, FRONT CASE , LP15	
K09	3305431-024	40577-000011	KIT - REPAIR, REAR CASE,LP15	Replaced with 40577-000051
K09	3305431-536	40577-000051	KIT - REPAIR, REAR CASE, LP15	
K10	3305431-026	40577-000025	KIT-REPAIR, MASIMO SPO2 PANEL MOUNT CABLE, LP15	
K11	3305431-538	40577-000053	KIT - REPAIR, SYSTEM PCBA V4	
K12	3305431-033	40577-000032	KIT - REPAIR, THERAPY PCBA V4	Replaced with 40577-000050
K12	3305431-535	40577-000050	KIT - REPAIR, THERAPY PCBA V4	
K13	3305431-034	40577-000033	KIT - REPAIR, THERAPY CONNECTOR V4	
K14	3305431-035	40577-000001	KIT - REPAIR, LED DISPLAY DRIVER, LP15	
K15	3305431-036	40577-000044	KIT - REPAIR, LED DISPLAY, LP15	



Service Repair Kits

K16	3305431-001	40577-000013	KIT - REPAIR, IP CONNECTORS, LP15	Replaced with 40577-000034
K16	3305431-501	40577-000034	KIT - REPAIR, IP CONNECTORS, LP15 ROHS	
K18	3305431-510	40577-000037	KIT - REPAIR, ECG CONNECTOR, LP15 ROHS	
K19	3305431-514	40577-000038	KIT - REPAIR, CO2 CONNECTOR, LP15 ROHS	MiniMedi Module
K20	3305431-525	40577-000039	KIT - REPAIR, MASIMO SPO2, LP15 ROHS	
K21	3305431-027	40577-000026	KIT - REPAIR, TEMPERATURE CABLE, LP15	Replaced with 40577-000040
K21	3305431-527	40577-000040	KIT - REPAIR, TEMPERATURE CABLE, LP15 ROHS	
K22	3305431-528	40577-000041	KIT - REPAIR, OEM PCBA, LP15 ROHS	
K23	3305431-529	40577-000042	KIT - REPAIR, INTERFACEPCBA, LP15 ROHS	Replaced with 40577-000062
K23	3305431-545	40577-000060	KIT – REPAIR, INTERFACE PCBA, NO BT	Replaced with 40577-000062
K23	3305431-547	40577-000062	KIT- REPAIR, INTERFACE PCBA	Requires software -110 or greater
K24	3305431-534	40577-000049	KIT - REPAIR, POWER PCBA , V2/V4 LP15 ROHS	
K25	3305431-531	40577-000036	KIT - REPAIR, NIBP MAXIQLP15 ROHS	Replaced with 40577-000061
K25	3305431-532	40577-000045	KIT - REPAIR, NIBP ORION LP15 ROHS	Replaced with 40577-000061
K25	3305431-546	40577-000061	KIT - REPAIR, NIBP LP15 ROHS	
K26	3305431-032	40577-000031	KIT - REPAIR, SYSTEM PCBA V4	Replaced with 40577-000048
K26	3305431-533	40577-000048	KIT - REPAIR, SYSTEM PCBA V4 ROHS	
K27	3305431-535	40577-000050	KIT – REPAIR, THERAPY PCBA V4	
K28	3207824-001	11577-000002	KIT - CARRY BAG, MAIN BAG	From Front Case Parts List



Service Repair Kits

Service Repair Kits				
K29	3207827-003	11260-000039	KIT - CARRY BAG, REAR POUCH, 3RD EDITION	From Front Case Parts List
K30	3207823-000	11577-000001	KIT - CARRY BAG, SHOULDER STRAP	From Front Case Parts List
K31	3010267-00	11220-000028	POUCH, TOP, CARRYING BAG, LP12	From Front Case Parts List
K32	3305431-530	40577-000043	KIT - REPAIR, POWER PCBA, V2/V4 LP15 ROHS	Replaced with 40577-000052
K32	3305431-537	40577-000052	KIT - REPAIR, POWER PCBA, V4 ROHS	
K33	3305431-542	40577-000057	KIT, REPAIR, NANOMEDI CO2 CONNECTOR, LP15, ROHS	Replaced with 40577-000058
K34	3305431-002	40577-000012	KIT - REPAIR, ETCO2 MODULE, LP15	Replaced with 40577-000058
K34	3305431-502	40577-000035	KIT - REPAIR, ETCO2 MODULE, LP15 ROHS	Replaced with 40577-000058
K34	3305431-543	40577-000058	KIT, REPAIR, NANOMEDI CO2 MODULE, LP15 ROHS	Requires software -109 or greater
K35	3305431-544	40577-000059	KIT, REPAIR, NANOMEDI CO2 HARDWARE, LP15 ROHS	