

Title	LP15 V2 Service BOM			Page	1 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

This document is electronically signed in the Physio-Control, Inc. Product Lifecycle Management (PLM) System. Approvals can be obtained from the PLM system which displays the selected approvers, their approval roles, and approval dates.

PROPRIETARY AND CONFIDENTIAL DOCUMENT FOR PHYSIO-CONTROL, INC. USE ONLY. This document is property of Physio-Control and may not be used, reproduced, published or disclosed to others without authorization from Physio-Control.



Physio-Control, Inc.
 11811 Willows Road NE, P.O. Box 97006, Redmond, WA 98052
 Tel 425.867.4000 Toll-free 800.442.1142
www.physio-control.com

Owner

Corporate

Attachments (Optional):	Qty	Pages	Page Sequencing
3347648_A_LP15_V1_V2_Service_BOM.pdf	1	42	1-42

Change Description
Initial Release. The LP15 V1/V2 Service BOM defines the individual replaceable parts or service kits for servicing the LP15 V1/V2 device.

Title	LP15 V2 Service BOM			Page	2 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

1. PURPOSE

This document provides the list of service parts for the LIFEPAK 15 V2 monitor/defibrillator and provides cross references between reference designators (ex. A01) used in the LIFEPAK 15 V1/V2 Service Manual (3300117), and part and catalog numbers, per the requirements in Service Support Document Development Work Instruction 7000742.

2. SCOPE

Scope is limited to individual parts and service repair kits which are available for service purposes for the LP15 V1/V2 Service Manual (3300117). Kit contents are not listed in this document and exist in electronic bills of materials (eBOMs) in the Product Lifecycle Management (PLM) system.

3. SERVICE BOM PART

Item Number	Description
3347648-000	Bill of Material, Service, LIFEPAK 15 V2

See 3347648-000 EBOM for parts list or reference the tables in Sections 4 and 5.

4. SERVICE INDIVIDUAL PARTS LIST

Note: Reference Designators are unique identifiers for the part used to cross-reference between the Service Manual or Service Process Work Instruction and the Service BOM. The reference designators have been assigned based on the following prefix system:

- A prefix- Printed Circuit Board Assemblies (PCBAs), keypads, and modules
- W prefix – Cables, subassemblies
- F prefix – Fasteners (e.g. screws, nuts, and standoffs)
- M prefix – Plastic parts and mechanical components
- E prefix – Electrical components
- L prefix – Labels and adhesives
- K prefix – Kits, or parts assembled into a kit for a specific repair.

Note: For some reference designators, there are multiple language options (e.g. labels and keypads). Refer to the tables below to select the version that applies.

[Click on a link to navigate to a parts list:](#)

- [External Parts List](#)
- [Carrying Case Parts List](#)
- [Front Case Parts List](#)
- [System/Therapy Assembly Parts List](#)
- [Parameter Bezel Parts List](#)
- [Rear Case Parts List](#)
- [SpO2 Internal Parts List](#)
- [NIBP Parts List](#)
- [CO2 Parts List](#)
- [Printer Control Keypad Parts List](#)
- [A10 Main Keypad -Languages Parts](#)
- [LIFEPAK 15 Label Set - Languages](#)
- [LIFEPAK 15 Label Set - Languages \(3rd Edition\)](#)
- [Service Repair Kits](#)

Title	LP15 V2 Service BOM		Page	3 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

External Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	Multiple	Multiple	V1	LABEL SET, language and feature-specific (V2- Auxiliary Power feature)	Refer to tables for LIFEPAK 15 Label Set Languages Parts List (L07)
L08	Multiple	Multiple	V2	LABEL SET, language and feature-specific (V2- Auxiliary Power feature) 3rd Edition	Refer to tables for 3rd Edition Label Sets Parts List (L08) for multiple language options.
A10	Multiple	Multiple	V1, V2	A10-KEYPAD ASSY -MAIN CNTRL (V2- Auxiliary Power feature)	Refer to Main Keypad Parts List for multiple language options.
A09	Multiple	Multiple	V1, V2	A09- KEYPAD ASSY - PRINTER CONTROL	Refer to Printer Control Keypads Parts List for multiple language options.
M01	3012120-004	21300-007445	V1, V2	COVER - CO2 CONNECTOR	
A12	3208290-001	11241-000017	V1, V2	PRINTER-100MM, GRAY,WITH FERRITES	Use up 11241-000009 before switching to this part.
M02	3206959-001	21330-001169	V1, V2	DOOR - BLANK, USB	
M03	3207307-001	21300-007356	V1, V2	GUARD - CORNER, LOWER RIGHT	Bumper (guard), lower right corner.
M04	3206968-000	21300-007010	V1, V2	GUARD - CORNER, LOWER LEFT	Bumper (guard), lower left corner
F01	3207361-375	21300-007253	V1, V2	SCREW-CAP,SCH,REC T15 TORX 6-32,0.375L,SS,NYLOK	
L01	3207131-011	21501-003005	V1, V2	LABEL-SERIAL NUMBER, V1, V2, LP15	Not orderable.
L01	3207131-001	21501-001997	V1, V2	LABEL-SERIAL NUMBER	Orderable only through Service Depot.
M60	3207247-001	21300-007319	V1, V2	FOOT - MOUNTING	
M61	3207575-001	21300-007553	V1, V2	SHIELD - SYSTEM PCBA, SBC, TOP	
M62	3318900-001	21300-008127	V1, V2	SHIELD-DIELECTRIC, OEM, POWER	
F19	3335400-000	21300-008193	V1, V2	PLUG, MINI-BANANA, NYLONPATCH	
L02	3207359-001	21501-001993	V1, V2	LABEL-SPO2 MASIMO PATENT	

Title	LP15 V2 Service BOM		Page	4 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

External Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F02	201407-069	21300-000777	V1, V2	SCREW,SELF-SEAL,SELF-LOCK,6-32X0.375	

Carrying Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F03	3207492-001	21300-007429	V1, V2	SCREW-PHILLIPS, TRUSS HEAD,8-32, 0.312IN LONG,SS	Mounts carry case at device bottom
F04	200476-708	21300-007966	V1, V2	SCREW-PH,P,SS,PV,1/4-20 X 1.00 L	Mounts carry case at top holes
F05	202253-576	21300-001033	V1, V2	SCREW,MACHINE, PANHEAD,NYLOK,6-32 X .750L	Mounts on rear pouch
K27	3207824-001	11577-000002	V1, V2	KIT - CARRY BAG, MAIN BAG	Standard Carry Case with Shoulder Strap
K28	3207827-003	11260-000039	V1, V2	KIT - CARRY BAG, REAR POUCH, 3RD EDITION	Back Pouch
K29	3207823-000	11577-000001	V1, V2	KIT - CARRY BAG, SHOULDER STRAP	Shoulder Strap
K30	3010267-00	11220-000028	V1, V2	POUCH,TOP,CARRYINGBAG,LP12	Top Pouch

Front Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F06	202253-761	21300-001038	V1, V2	SCREW-M,CS,Z,PH, NYLOK,4-40 X .312L	
F07	202253-764	21300-004599	V1, V2	SCREW,MACHINE,PANHEAD,NYLOK,4-40 X .500	
M05	3006809-00	21300-001299	V1, V2	GASKET, PRINTER, SUB-D	
M06	3006810-01	21300-006100	V1, V2	BRACKET, CONNECTOR, SUBD	
W16	3009724-00	21330-000165	V1, V2	CABLE ASSY, RIBBON,INTERFACE PCB/CHART RCDR	
W19	3009726-01	21300-001370	V1, V2	WIRE HARNESS - GROUND/PRINTER BRKT	

Title	LP15 V2 Service BOM		Page	5 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Front Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
W17	3009726-03	21300-001371	V1, V2	WIRE HARNESS-SPKRINTERFACE PCB	
W15	3011128-002	21300-003713	V1, V2	SWITCH-ROTARY	
M07	3012693-00	21300-001575	V1, V2	SPRING CLAMP, STAINLESS STEEL	
A11	3319977-000	21300-008140	V1, V2	DISPLAY, ACTIVE MATRIX, 8.4 IN, 262K COLOR, ANTI GLARE	See Display Repair Kit (K14)
M08	3206724-000	21330-001173	V1, V2	LENS-DISPLAY	See Front Case Repair Kit (K10)
M09	3206725-001	21330-001174	V1, V2	SHIELD-DISPLAY	See Display Shield Repair Kit (K07)
M10	3206726-003	21330-001641	V1, V2	FRAME - DISPLAY SHIELD, V5	See Display Shield Repair Kit (K07)
M11	3305427-000	21330-001423	V1, V2	ASSEMBLY, SHIELD, DISPLAY, CLEAR	See Display Shield Repair Kit (K07)
A05	3206815-003	21330-001236	V1	PCB ASSY - INTERFACE	No longer available for order. Order V1 Interface PCB Kit 40577-000023.
A05	3317231-002	21300-001533	V2	PCB ASSY - INTERFACE, LP15 V4	Interface PCBA and Repair Kit no longer available for order.
M12	3206936-002	21330-001189	V1, V2	BRACKET - LCD DISPLAY MOUNTING	
W12	3206989-000	21330-001214	V1, V2	CABLE ASSY - SM KEYPAD, INTERFACE PCB	
W18	3206990-001	21330-001417	V1, V2	CABLE ASSY - FLEX, LCD/INTERFACE PCB	
W04	3206991-003	21330-001226	V1, V2	CABLE ASSY - FLEX, SYSTEM PCB/INTERFACE PCB	
W06	3318862-000	21330-001525	V1, V2	CABLE ASSY - BACKINV, INTERFACE PCB, ROHS	
M13	3207030-000	21300-007051	V1, V2	KNOB-ROTARY SWITCH	
W11	3207044-002	21300-007463	V1, V2	WIRE HARNESS-THERAPY, INTERNAL	See Therapy Connector Repair Kit (K03)

Title	LP15 V2 Service BOM			Page	6 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

Front Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
M14	3207252-001	21300-007079	V1, V2	SHIELD,INVERTER	
F08	3207337-312	21300-007297	V1, V2	SCREW-MACH,PNH,PHH,NYLOCK,4-40,0.312L,WSHR,CS,ZN	
F01	3207361-375	21300-007253	V1, V2	SCREW-CAP,SCH,REC T15 TORX 6-32,0.375L,SS,NYLOK	
F09	3207367-312	21300-007196	V1, V2	SCREW-MACH,TRH,T10 TORX,4-40,0.312L,SS	
W13	3207388-001	21330-001323	V1, V2	CABLE-ASSY,MAIN KEYPAD,INTERFACE PCB	
M15	3207701-000	21300-007464	V1, V2	SEAL-THERAPY RECEPTACLE, LP15	
A25	3319108-000	21300-008138	V1, V2	DRIVER,LED BACKLIGHT,12VDC,85C	Not orderable. See LED Backlight Repair Kit K13.
M16	3207725-002	21330-001273	V1, V2	ASSEMBLY-ENCLOSURE,FRONT	No longer available for order. Order Front Case Repair Kit 40577-000010 (K10)
A08	3207933-000	21300-007555	V1, V2	INVERTER-BACKLIGHT,DUAL TUBE,12VDC,85C	
L03	3208003-000	21501-002038	V1, V2	LABEL-PHYSIO-CONTROL ICON	
L04	800943-09	21501-000248	V1, V2	LABEL-SYMBOL INTL	
M17	804234-03	21501-002394	V1, V2	SEAL-PERIMETER,CASE	

System/Therapy Assembly Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F10	202253-550	21300-001030	V1, V2	SCREW-PH,NYLOCK,PHIL,SS,#4-40X.1.125L	
F06	202253-761	21300-001038	V1, V2	SCREW-M,CS,Z,PH, NYLOCK,4-40 X .312L	
A04	3306311-000	21330-001449	V1, V2	PCB ASSY- THERAPY, LP15	See Therapy PCBA Repair Kit (K02)
M18	3009331-00	21300-001343	V1, V2	SHIELD EMI, THERAPY PCB	

Title	LP15 V2 Service BOM		Page	7 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

System/Therapy Assembly Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
W03	3009878-002	21300-001379	V1, V2	CONN-HDR, SQUAREPIN DUAL ROW AU PLATED, 30 PIN	
F11	3011629-00	21300-001478	V1, V2	INSERT-HEX,LP12	
F12	3011630-00	21300-001479	V1, V2	SPACER-PCB,LP12	
A01	3206834-013	21330-001448	V1	PCB ASSY - SYSTEM, LP15	V1 System PCBA Repair Kit no longer available for order. No substitute part.
M19	3207010-004	21330-001318	V1, V2	ASSY-MOUNTING,BRACKET,SYSTEM, THERAPY PCB	
L05	804447-27	21300-003172	V1, V2	SPACER-FOAM, SYS PCBCAPS	
E01	202305-000	21300-001052	V1, V2	BATTERY-LITH,3V,COIN	

Parameter Bezel Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L06	3207407-000	21300-007362	V1, V2	SPACER - CONNECTOR, CO2	
M20	3206723-005	21330-001211	V1, V2	BEZEL - PARAMETER	
W07	3318313-001	21330-001527	V1, V2	RECEPTACLE-ECG,12 CONTACT,WIRED, ROHS	See ECG Connector Repair Kit (K17)
M21	805915-01	21300-002683	V1, V2	SEAL-ECG RECEPTACLE	
M22	3205680-001	21300-006304	V1, V2	BRACKET - CABLE SUPPORT, ASSY	
F13	202253-763	21300-006965	V1, V2	SCREW,MACHINE,PANHEAD,NY LOK,4-40 X .437	
W22	3206995-004	21330-001290	V1, V2	CABLE ASSY - FLEX, SPO2, MASIMO	See MASIMO SPO2 PANEL MOUNT CABLE Repair Kit (K11)
M23	3207009-000	21300-007112	V1, V2	HOUSING-FERRITE,SPO2	
F14	3205311-001	21300-007557	V1, V2	SCREW-MACHINE, PAN, TORX, NYLOK, 4-40, .6875L	
W28	3317964-004	21330-001529	V1, V2	FRS ASSEMBLY-CO2, MINI MODULE, ROHS	See Connector Repair Kit (K18)

Title	LP15 V2 Service BOM			Page	8 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

Parameter Bezel Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
M24	3007997-01	21300-001314	V1, V2	SEAL-CO2 CONNECTOR	
M25	3012119-02	21300-001548	V1, V2	ADAPTER-CO2 CONNECTOR	
M26	3012121-002	21300-007998	V1, V2	RETAINER-CO2 CONN	
F15	202253-732	21300-007353	V1, V2	SCREW- MACH,NYLOK,PNH,PHH,2- 56X0.500,CS,ZN	
M27	3207033-001	21300-007425	V1, V2	CONNECTOR - PNEUMATIC COUPLER, NIBP	See NIBP Connector Repair Kit (K25)
M28	3012180-02	21300-01560	V1, V2	TUBING-NIBP	
M29	3012333-001	21300-001566	V1, V2	CONN- PNEU,STR,0.125IDTBG,NYL,WH T	
W30	3318652-000	21330-001526	V1, V2	WIRE HARNESS-INVASIVE PRESSURE 1 AND 2, ROHS	See Invasive Pressure Connector Repair Kit (K15)
M30	3007998-00	21300-001315	V1, V2	GASKET-CONNECTOR,INVASIVE PRESURE	
F08	3207337-312	21300-007297	V1, V2	SCREW- MACH,PNH,PHH,NYLOCK,4- 40,0.312L,WSHR,CS,ZN	
W30	3200466-01	21300-001583	V1, V2	WIRE HARNESS- INVASIVEPRESSURE 1 / 2	Part no longer available for order. See Repair Kit K15.
W31	3319355-000	21330-001522	V1, V2	CABLE ASSEMBLY - TEMPERATURE, ROHS	See Temperature Cable Repair Kit (K20)

Rear Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A03	3302519-002	21330-001408	V2	PCB ASSY - POWER,LP15	Part and kit no longer available for order. No substitute part.

Title	LP15 V2 Service BOM		Page	9 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Rear Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A03	3206749-003	21330-001404	V1	PCB ASSY - POWER, LP15	Part and kit no longer available for order. No substitute part.
A06	3317232-000	21330-001532	V2	PCB ASSY - OEM, LP15	See OEM PCB Repair Kit (K21)
A06	3206813-003	21330-001285	V1	PCB ASSY - OEM, LP15	Part and kit no longer available. No substitute.
A07	3207037-003	21330-001288	V1, V2	CONTACT ASSY-BATTERY	
A13	3201583-001	21330-000920	V1, V2	RELAY ASSEMBLY - ISOLATION, DPDT, BTE, LP12	
A14	3010212-02	21300-001392	V1, V2	INDUCTIVE RESISTOR-5 OHM	
A15	3008164-002	21300-001320	V1, V2	CAPACITOR-ENERGY STORAGE,BTE,LP12	
A17	3008897-002	21300-001337	V1, V2	BRACKET-INTERCONNECT	
A22	3010178-010	21330-000176	V1, V2	PCB ASSY - BIPHASIC MODULE, CEDAR, LP12	Order Biphasic PCB Kit (K04).
M31	3207726-000	21330-001272	V1, V2	ASSEMBLY - ENCLOSURE, REAR	Order Rear Case Kit (K23).
M32	3207881-000	21330-001276	V1, V2	ASSEMBLY-RETAINER, BATTERY	
M33	3207031-001	21300-007317	V1, V2	BRACKET-CAPACITOR	
M34	3206965-001	21330-001195	V1, V2	BRACKET - COVER, OEM MODULES	
M35	3206964-001	21330-001194	V1, V2	BRACKET - MOUNTING, OEM MODULES	
M36	3206961-003	21330-001257	V1, V2	BRACKET ASSY - POWER PCB	
M37	3011589-004	21300-007546	V1, V2	BRACKET-MOUNTING, BTE PCB	
M38	3207708-000	21300-007472	V1, V2	BRACKET-SHOULDER STRAP	
F22	802278-02	21300-002023	V1, V2	CONN-PLUG,BANANA PIN,LP10	Part replaced with 21300-008193 (F19).
M40	3006766-003	21330-001241	V1, V2	COVER - LATCH, PADDLE, GRAY	See Retainer Repair Kit (K24).

Title	LP15 V2 Service BOM		Page	10 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Rear Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
M41	3207746-000	21300-007470	V1, V2	COVER PLATE - HANDLE, RIGHT	
M42	3207707-000	21300-007469	V1, V2	COVER PLATE-HANDLE, LEFT	
M43	3206960-002	21330-001190	V1, V2	GROUND PLANE - FLEXIBLE REAR CASE	
M44	3207706-000	21300-007471	V1, V2	HANDLE	See Handle Repair Kit (K09).
L04	800943-09	21501-000248	V1, V2	LABEL-SYMBOL INTL	
L09	3320539-000	21501-002882	V1, V2	LABEL - UDI	
L10	3011980-00	21501-000152	V1, V2	LABEL-SHOCK HAZARD	Physio order only.
F16	806091-00	21300-002691	V1, V2	NUT-AUX,SYSTEM,CONN,STAINLESS STEEL	
F17	201508-000	21300-000804	V1, V2	NUT-LOCK,CARBON STEEL,KEP,#4-40 THREAD	Replaced by F21.
M45	805487-00	21300-002650	V1, V2	PLATE-SEAL,CONNECTOR,REAR	
M46	200536-011	21300-000504	V1, V2	RETNR-CABLE TIE,NYLON, .35W X 21.0 L	
M47	200536-001	21300-000499	V1, V2	RETAINR-CABLE TIE,NYLON,.10 W, X 4.0 L	
F06	202253-761	21300-001038	V1, V2	SCREW-M,CS,Z,PH, NYLOCK,4-40 X .312L	
F18	202253-760	21300-006251	V1, V2	SCREW-M,PH,NYLOK,CS,4-40,.250L	
F08	3207337-312	21300-007297	V1, V2	SCREW-MACH,PNH,PHH,NYLOCK,4-40,0.312L,WSHR,CS,ZN	
F02	201407-069	21300-000777	V1, V2	SCREW,SELF-SEAL,SELF-LOCK,6-32X0.375	
F13	202253-763	21300-006965	V1, V2	SCREW,MACHINE,PANHEAD,N YLOK,4-40 X .437	
M48	200060-011	21300-000203	V1, V2	SEAL-O-RING,RUBBER, .551ID	
M49	3006291-00	21300-001284	V1, V2	SEAL, DRAIN	
M51	3208298-000	21300-007744	V1, V2	SHIELD-DIELECTRIC, OEM, POWER	
M52	805542-00	21300-002660	V1, V2	SHIELD-	

Title	LP15 V2 Service BOM		Page	11 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Rear Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				END,CAPACITOR,SLOTTED	
M53	3010593-00	21300-001407	V1, V2	SHIELD-HIGH VOLTAGE	
M07	3012693-00	21300-001575	V1, V2	SPRING CLAMP, STAINLESS STEEL	
F19	3010569-02	21300-001402	V1, V2	STANDOFF, HEX, NYLON	
F20	3340265-000	21300-008213	V1, V2	WASHER, NIBP COUPLER, LP15	
F21	200805-000	21300-000584	V1, V2	NUT-HEX,SS,LOCK 4-40X.250W	
M54	3012345-00	21300-001567	V1, V2	STRAP-GROUND,BTE PCB	
M58	3009787-003	21300-001377	V1, V2	SHIELD - DIELECTRIC,CAPACITOR	
M59	804234-03	21300-002394	V1, V2	SEAL-PERIMETER,CASE	
M63	3318613-000	21300-008135	V1, V2	MICRO WIRE SADDLE, TEAR DROP,NYLON 6/6,ROHS	
M64	3312316-001	21300-008205	V1, V2	BRACKET ASSY-MOUNTING, SYSTEM-THERAPY PCBA	
M65	3340287-000	21300-008212	V1, V2	RETAINER CLIP, NIBP COUPLER, LP15	
W01	3207692-000	21300-007456	V1, V2	CABLE ASSY - POWER,SYSTEM, PCBA	
W02	3009726-05	21300-001373	V1, V2	WIRE HARNESS-POWER-THERAPY-PCB	
W05	3319354-000	21330-001316	V1, V2	CABLE ASSEMBLY- CONTACT PCB/POWER PCB,LP15 ROHS	
W08	3009652-01	21330-000156	V1, V2	CABLE ASSY-SYSTEM CONNECTOR	Order catalog number 21330-000156 in FDA regions.
W08	3009652-005	21330-001554	V1, V2	HARNESS, CONNECTOR, SYSTEM	Order catalog number 21330-000156 in FDA regions.
W09	3303848-003	21330-001485	V2	CABLE ASSEMBLY - AUXILIARY POWER	
W10	3303863-000	21330-001394	V2	CABLE ASSEMBLY - BATTERY POWER	
W10	3009726-08	21300-001374	V1	WIRE HARNESS, BATTERY POWER	
W14	3318946-000	21330-001534	V1, V2	FLEX ASSY - USB, ROHS	
W20	3011792-005	21300-007670	V1, V2	CABLE - FLEX, 12 CONTACT,	This cable

Title	LP15 V2 Service BOM			Page	12 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

Rear Case Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				MTE	replaced 21300-001528.
W20	3206999-003	21330-001286	V1, V2	CABLE ASSY-FLEX, BTE PCB, SYSTEM,THERAPY PCB	No longer available. Order 21300-007670.
W23	3011979-00	21300-001541	V1, V2	WIRE HARNESS-BTE,LP12, J103-PIN9	
W24	3011979-01	21300-001542	V1, V2	WIRE HARNESS-BTE,LP12, J101- PIN3	
W25	3011979-02	21300-001543	V1, V2	WIRE HARNESS-BTE,LP12, J105 TO PIN6	
W36	3207000-003	21330-001321	V1, V2	W21 - ASSY-CABLE, FLEX, Sp02, OEM, PCB	Replaced by W21.

SpO2 Internal Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F18	202253-760	21300-006251	V1, V2	SCREW-M,PH,NYLOK,CS,4-40,.250L	
F07	202253-764	21300-004599	V1, V2	SCREW,MACHINE,PANHEAD,N YLOK,4-40 X .500	
M55	3206962-001	21300-007550	V1, V2	BRACKET - MOUNTING, SPO2 MODULE	
W21	3318916-000	21330-001524	V1, V2	CABLE, ASSY-FLEX, SPO2, OEM PCB, ROHS	
A16	3317849-001	21330-001530	V1, V2	MODULE-OEM,PULSE OXIMETER, MX-51, ROHS	No longer available for order. See Masimo SPO2 Repair Kit (K19).

NIBP Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F18	202253-760	21300-006251	V1, V2	SCREW-M,PH,NYLOK,CS,4-40,.250L	

Title	LP15 V2 Service BOM			Page	13 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

NIBP Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
M56	3010805-000	21300-001416	V1, V2	SOCKET RETAINER CLIP -10 PIN	
W27	3012181-00	21330-000235	V1, V2	CABLE ASSY-RIBBON,NIBP/OEM PCB	
A21	3317965-003	21330-001556	V1, V2	MODULE-NIBP, ORION, ROHS	See NIBP Module Repair Kit (K22).
A24	3317965-003	21330-001528	V1, V2	MODULE-NIBP, MAXIQ, ROHS	No longer available for order. See NIBP Repair Kit K22.

CO2 Parts List					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
F18	202253-760	21300-006251	V1, V2	SCREW-M,PH,NYLOK,CS,4-40,.250L	
M57	3012140-002	21300-001557	V1, V2	EXHAUST TUBING-CO2	
A23	3317964-003	21330-001535	V1, V2	MODULE - MINIMEDI CO2, ROHS	No longer available for order. See CO2 Module Repair Kit (K16).
W29	3012397-01	21330-000237	V1, V2	CABLE ASSY-CO2, ADAPTER	
W26	3206998-002	21330-001217	V1, V2	CABLE ASSY-CO2 MODULE, OEM PCB	

Printer Control Keypad Parts List (V1 and V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A09	3207080-000	21330-001238	V1, V2	KEYPAD ASSY - PRINTER,12 LEAD CONTROL, MEMBRANE	
A09	3207080-001	21330-001264	V1, V2	KEYPAD ASSY-PRINTER, ENGLISH	
A09	3207080-040	21330-001326	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, GERMAN	

Title	LP15 V2 Service BOM		Page	14 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Printer Control Keypad Parts List (V1 and V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A09	3207080-041	21330-001327	V1, V2	KEYPAD ASSY-PRINTER CONTROL, GERMAN	
A09	3207080-060	21330-001328	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, ITALIAN	
A09	3207080-061	21330-001329	V1, V2	KEYPAD ASSY-PRINTER CONTROL, ITALIAN	
A09	3207080-080	21330-001330	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, FRENCH	
A09	3207080-081	21330-001331	V1, V2	KEYPAD ASSY-PRINTER CONTROL, FRENCH	
A09	3207080-100	21330-001332	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, DUTCH	
A09	3207080-101	21330-001333	V1, V2	KEYPAD ASSY-PRINTER CONTROL, DUTCH	
A09	3207080-120	21330-001334	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, SPANISH	
A09	3207080-121	21330-001335	V1, V2	KEYPAD ASSY-PRINTER CONTROL, SPANISH	
A09	3207080-140	21330-001383	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, PORTUGUESE	
A09	3207080-141	21330-001395	V1, V2	KEYPAD ASSY-PRINTER CONTROL, PORTUGUESE	
A09	3207080-160	21330-001441	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, BRAZILIAN	
A09	3207080-161	21330-001442	V1, V2	KEYPAD ASSY-PRINTER CONTROL, BRAZILIAN	
A09	3207080-180	21330-001336	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, SWEDISH	
A09	3207080-181	21330-001337	V1, V2	KEYPAD ASSY-PRINTER CONTROL, SWEDISH	
A09	3207080-200	21330-001338	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, DANISH	
A09	3207080-201	21330-001339	V1, V2	KEYPAD ASSY-PRINTER CONTROL, DANISH	
A09	3207080-220	21330-001384	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, FINNISH	
A09	3207080-221	21330-001402	V1, V2	KEYPAD ASSY-PRINTER CONTROL, FINNISH	

Title	LP15 V2 Service BOM		Page	15 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Printer Control Keypad Parts List (V1 and V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A09	3207080-240	21330-001340	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, NORWEGIAN	
A09	3207080-241	21330-001341	V1, V2	KEYPAD ASSY-PRINTER CONTROL, NORWEGIAN	
A09	3207080-260	21330-001342	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, POLISH	
A09	3207080-261	21330-001343	V1, V2	KEYPAD ASSY-PRINTER CONTROL, POLISH	
A09	3207080-280	21330-001385	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, HUNGARIAN	
A09	3207080-281	21330-001398	V1, V2	KEYPAD ASSY-PRINTER CONTROL, HUNGARIAN	
A09	3207080-300	21330-001386	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, CZECH	
A09	3207080-301	21330-001401	V1, V2	KEYPAD ASSY-PRINTER CONTROL, CZECH	
A09	3207080-320	21330-001443	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, RUSSIAN	
A09	3207080-321	21330-001444	V1, V2	KEYPAD ASSY-PRINTER CONTROL, RUSSIAN	
A09	3207080-380	21330-001445	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, KOREAN	
A09	3207080-381	21330-001446	V1, V2	KEYPAD ASSY-PRINTER CONTROL, KOREAN	
A09	3207080-400	21330-001504	V1, V2	KEYPAD ASSY-PRINTER CONTROL, 12-LEAD, JAPANESE	
A09	3207080-401	21330-001505	V1, V2	KEYPAD ASSY-PRINTER CONTROL, JAPANESE	

Title	LP15 V2 Service BOM		Page	16 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Main Keypad Parts List (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A10	3207079-007	21330-001575	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, ENGLISH	
A10	3207079-008	21330-001576	V1	KEYPAD ASSY-MAIN CONTROL, CPR, ENG	
A10	3207079-045	21330-001577	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, GERMAN	
A10	3207079-046	21330-001578	V1	KEYPAD ASSY-MAIN CONTROL, CPR, GERMAN	
A10	3207079-065	21330-001579	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, ITALIAN	
A10	3207079-066	21330-001580	V1	KEYPAD ASSY-MAIN CONTROL, CPR, ITALIAN	
A10	3207079-085	21330-001581	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, FRENCH	
A10	3207079-086	21330-001582	V1	KEYPAD ASSY-MAIN CONTROL, CPR, FRENCH	
A10	3207079-105	21330-001583	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, DUTCH	
A10	3207079-106	21330-001584	V1	KEYPAD ASSY-MAIN CONTROL, CPR, DUTCH	
A10	3207079-125	21330-001585	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, SPANISH	
A10	3207079-126	21330-001586	V1	KEYPAD ASSY-MAIN CONTROL, CPR, SPANISH	
A10	3207079-145	21330-001587	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, PORTUGUESE	
A10	3207079-146	21330-001588	V1	KEYPAD ASSY-MAIN CONTROL, CPR, PORTUGUESE	
A10	3207079-185	21330-001589	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, SWEDISH	
A10	3207079-186	21330-001590	V1	KEYPAD ASSY-MAIN CONTROL, CPR, SWEDISH	
A10	3207079-205	21330-001591	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, DANISH	
A10	3207079-206	21330-001592	V1	KEYPAD ASSY-MAIN CONTROL, CPR, DANISH	
A10	3207079-225	21330-001593	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, FINNISH	
A10	3207079-226	21330-001594	V1	KEYPAD ASSY-MAIN CONTROL,	

Title	LP15 V2 Service BOM		Page	17 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Main Keypad Parts List (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				CPR, FINNISH	
A10	3207079-245	21330-001595	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, NORWEGIAN	
A10	3207079-246	21330-001596	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NORWEGIAN	
A10	3207079-265	21330-001597	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, POLISH	
A10	3207079-266	21330-001598	V1	KEYPAD ASSY-MAIN CONTROL, CPR, POLISH	
A10	3207079-285	21330-001599	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, HUNGARIAN	
A10	3207079-286	21330-001600	V1	KEYPAD ASSY-MAIN CONTROL, CPR, HUNGARIAN	
A10	3207079-305	21330-001601	V1	KEYPAD ASSY-MAIN CONTROL, CPR, NIBP, CZECH	
A10	3207079-306	21330-001602	V1	KEYPAD ASSY-MAIN CONTROL, CPR, CZECH	

Main Keypad Parts List (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A10	3302470-006	21330-001603	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, ENGLISH	
A10	3302470-007	21330-001604	V2	KEYPAD ASSY-MAIN CONTROL, ENGLISH	
A10	3302470-046	21330-001605	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, GERMAN	
A10	3302470-047	21330-001606	V2	KEYPAD ASSY-MAIN CONTROL, GERMAN	
A10	3302470-066	21330-001607	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, ITALIAN	
A10	3302470-067	21330-001608	V2	KEYPAD ASSY-MAIN CONTROL, ITALIAN	
A10	3302470-086	21330-001609	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, FRENCH	
A10	3302470-087	21330-001610	V2	KEYPAD ASSY-MAIN CONTROL,	

Title	LP15 V2 Service BOM		Page	18 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Main Keypad Parts List (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				FRENCH	
A10	3302470-106	21330-001611	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, DUTCH	
A10	3302470-107	21330-001612	V2	KEYPAD ASSY-MAIN CONTROL, DUTCH	
A10	3302470-126	21330-001613	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, SPANISH	
A10	3302470-127	21330-001614	V2	KEYPAD ASSY-MAIN CONTROL, SPANISH	
A10	3302470-146	21330-001615	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, PORTUGUESE	
A10	3302470-147	21330-001616	V2	KEYPAD ASSY-MAIN CONTROL, PORTUGUESE	
A10	3302470-166	21330-001617	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, BRAZILIAN	
A10	3302470-167	21330-001618	V2	KEYPAD ASSY-MAIN CONTROL, BRAZILIAN	
A10	3302470-186	21330-001619	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, SWEDISH	
A10	3302470-187	21330-001620	V2	KEYPAD ASSY-MAIN CONTROL, SWEDISH	
A10	3302470-206	21330-001621	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, DANISH	
A10	3302470-207	21330-001622	V2	KEYPAD ASSY-MAIN CONTROL, DANISH	
A10	3302470-226	21330-001623	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, FINNISH	
A10	3302470-227	21330-001624	V2	KEYPAD ASSY-MAIN CONTROL, FINNISH	
A10	3302470-246	21330-001625	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, NORWEGIAN	
A10	3302470-247	21330-001626	V2	KEYPAD ASSY-MAIN CONTROL, NORWEGIAN	
A10	3302470-266	21330-001627	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, POLISH	
A10	3302470-267	21330-001628	V2	KEYPAD ASSY-MAIN CONTROL, POLISH	

Title	LP15 V2 Service BOM		Page	19 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Main Keypad Parts List (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
A10	3302470-286	21330-001629	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, HUNGARIAN	
A10	3302470-287	21330-001630	V2	KEYPAD ASSY-MAIN CONTROL, HUNGARIAN	
A10	3302470-306	21330-001631	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, CZECH	
A10	3302470-307	21330-001632	V2	KEYPAD ASSY-MAIN CONTROL, CZECH	
A10	3302470-326	21330-001633	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, RUSSIAN	
A10	3302470-327	21330-001634	V2	KEYPAD ASSY-MAIN CONTROL, RUSSIAN	
A10	3302470-386	21330-001635	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, KOREAN	
A10	3302470-387	21330-001636	V2	KEYPAD ASSY-MAIN CONTROL, KOREAN	
A10	3302470-404	21330-001637	V2	KEYPAD ASSY-MAIN CONTROL, NIBP, JAPANESE	
A10	3302470-405	21330-001638	V2	KEYPAD ASSY-MAIN CONTROL, JAPANESE	

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3207318-013	21501-002039	V1	LABEL SET-ECG, ENGLISH	
L07	3207318-014	21501-002040	V1	LABEL SET-ECG,SPO2, ENGLISH	
L07	3207318-015	21501-002041	V1	LABEL SET-ECG,SPO2,CO2,ENGLISH	
L07	3207318-016	21501-002042	V1	LABEL SET-ECG,SPO2,NIBP,IP, ENGLISH	
L07	3207318-017	21501-002043	V1	LABEL SET-ECG,SPO2,NIBP,CO2,IP, ENGLISH	
L07	3207318-018	21501-002044	V1	LABEL SET-ECG,NIBP,SPO2,ENGLISH	

Title	LP15 V2 Service BOM		Page	20 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3207318-019	21501-002045	V1	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH	
L07	3207318-040	21501-002050	V1	LABEL SET-ECG, GERMAN	
L07	3207318-041	21501-002051	V1	LABEL SET-ECG, SPO2, GERMAN	
L07	3207318-042	21501-002052	V1	LABEL SET-ECG, SPO2, CO2, GERMAN	
L07	3207318-043	21501-002053	V1	LABEL SET-ECG, SPO2, NIBP, IP, GERMAN	
L07	3207318-044	21501-002054	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, GERMAN	
L07	3207318-045	21501-002055	V1	LABEL SET-ECG, NIBP, SPO2, GERMAN	
L07	3207318-046	21501-002056	V1	LABEL SET-ECG, SPO2, NIBP, CO2, GERMAN	
L07	3207318-060	21501-002057	V1	LABEL SET-ECG, ITALIAN	
L07	3207318-061	21501-002058	V1	LABEL SET-ECG, SPO2, ITALIAN	
L07	3207318-062	21501-002059	V1	LABEL SET-ECG, SPO2, CO2, ITALIAN	
L07	3207318-063	21501-002060	V1	LABEL SET-ECG, SPO2, NIBP, IP, ITALIAN	
L07	3207318-064	21501-002061	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, ITALIAN	
L07	3207318-065	21501-002062	V1	LABEL SET-ECG, NIBP, SPO2, ITALIAN	
L07	3207318-066	21501-002063	V1	LABEL SET-ECG, SPO2, NIBP, CO2, ITALIAN	
L07	3207318-080	21501-002064	V1	LABEL SET-ECG, FRENCH	
L07	3207318-081	21501-002065	V1	LABEL SET-ECG, SPO2, FRENCH	
L07	3207318-082	21501-002066	V1	LABEL SET-ECG, SPO2, CO2, FRENCH	
L07	3207318-083	21501-002067	V1	LABEL SET-ECG, SPO2, NIBP, IP, FRENCH	
L07	3207318-084	21501-002068	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, FRENCH	
L07	3207318-085	21501-002069	V1	LABEL SET-ECG, NIBP, SPO2, FRENCH	
L07	3207318-086	21501-002070	V1	LABEL SET-ECG, SPO2, NIBP, CO2, FRENCH	

Title	LP15 V2 Service BOM		Page	21 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3207318-100	21501-002162	V1	LABEL SET-ECG, DUTCH	
L07	3207318-101	21501-002163	V1	LABEL SET-ECG, SPO2, DUTCH	
L07	3207318-102	21501-002164	V1	LABEL SET-ECG, SPO2, CO2, DUTCH	
L07	3207318-103	21501-002165	V1	LABEL SET-ECG, SPO2, NIBP, IP, DUTCH	
L07	3207318-104	21501-002166	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, DUTCH	
L07	3207318-105	21501-002167	V1	LABEL SET-ECG, NIBP, SPO2, DUTCH	
L07	3207318-106	21501-002168	V1	LABEL SET-ECG, SPO2, NIBP, CO2, DUTCH	
L07	3207318-120	21501-002169	V1	LABEL SET-ECG, SPANISH	
L07	3207318-121	21501-002170	V1	LABEL SET-ECG, SPO2, SPANISH	
L07	3207318-122	21501-002171	V1	LABEL SET-ECG, SPO2, CO2, SPANISH	
L07	3207318-123	21501-002172	V1	LABEL SET-ECG, SPO2, NIBP, IP, SPANISH	
L07	3207318-124	21501-002173	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, SPANISH	
L07	3207318-125	21501-002174	V1	LABEL SET-ECG, NIBP, SPO2, SPANISH	
L07	3207318-126	21501-002175	V1	LABEL SET-ECG, SPO2, NIBP, CO2, SPANISH	
L07	3207318-140	21501-002297	V1	LABEL SET-ECG, PORTUGUESE	
L07	3207318-141	21501-002298	V1	LABEL SET-ECG, SPO2, PORTUGUESE	
L07	3207318-142	21501-002299	V1	LABEL SET-ECG, SPO2, CO2, PORTUGUESE	
L07	3207318-143	21501-002300	V1	LABEL SET-ECG, SPO2, NIBP, IP, PORTUGUESE	
L07	3207318-144	21501-002301	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, PORTUGUESE	
L07	3207318-145	21501-002302	V1	LABEL SET-ECG, NIBP, SPO2, PORTUGUESE	
L07	3207318-146	21501-002303	V1	LABEL SET-ECG, SPO2, NIBP, CO2, PORTUGUESE	

Title	LP15 V2 Service BOM		Page	22 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3207318-180	21501-002176	V1	LABEL SET-ECG, SWEDISH	
L07	3207318-181	21501-002177	V1	LABEL SET-ECG, SPO2, SWEDISH	
L07	3207318-182	21501-002178	V1	LABEL SET-ECG, SPO2, CO2, SWEDISH	
L07	3207318-183	21501-002179	V1	LABEL SET-ECG, SPO2, NIBP, IP, SWEDISH	
L07	3207318-184	21501-002180	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, SWEDISH	
L07	3207318-185	21501-002181	V1	LABEL SET-ECG, NIBP, SPO2, SWEDISH	
L07	3207318-186	21501-002182	V1	LABEL SET-ECG, SPO2, NIBP, CO2, SWEDISH	
L07	3207318-200	21501-002183	V1	LABEL SET-ECG, DANISH	
L07	3207318-201	21501-002184	V1	LABEL SET-ECG, SPO2, DANISH	
L07	3207318-202	21501-002185	V1	LABEL SET-ECG, SPO2, CO2, DANISH	
L07	3207318-203	21501-002186	V1	LABEL SET-ECG, SPO2, NIBP, IP, DANISH	
L07	3207318-204	21501-002187	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, DANISH	
L07	3207318-205	21501-002188	V1	LABEL SET-ECG, NIBP, SPO2, DANISH	
L07	3207318-206	21501-002189	V1	LABEL SET-ECG, SPO2, NIBP, CO2, DANISH	
L07	3207318-220	21501-002304	V1	LABEL SET-ECG, FINNISH	
L07	3207318-221	21501-002306	V1	LABEL SET-ECG, SPO2, FINNISH	
L07	3207318-222	21501-002305	V1	LABEL SET-ECG, SPO2, CO2, FINNISH	
L07	3207318-223	21501-002307	V1	LABEL SET-ECG, SPO2, NIBP, IP, FINNISH	
L07	3207318-224	21501-002308	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, FINNISH	
L07	3207318-225	21501-002309	V1	LABEL SET-ECG, NIBP, SPO2, FINNISH	
L07	3207318-226	21501-002310	V1	LABEL SET-ECG, SPO2, NIBP, CO2, FINNISH	
L07	3207318-240	21501-002190	V1	LABEL SET-ECG, NORWEGIAN	
L07	3207318-241	21501-002191	V1	LABEL SET-ECG, SPO2,	

Title	LP15 V2 Service BOM		Page	23 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				NORWEGIAN	
L07	3207318-242	21501-002192	V1	LABEL SET-ECG, SPO2, CO2, NORWEGIAN	
L07	3207318-243	21501-002193	V1	LABEL SET-ECG, SPO2, NIBP, IP, NORWEGIAN	
L07	3207318-244	21501-002194	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, NORWEGIAN	
L07	3207318-245	21501-002195	V1	LABEL SET-ECG, NIBP, SPO2, NORWEGIAN	
L07	3207318-246	21501-002196	V1	LABEL SET-ECG, SPO2, NIBP, CO2, NORWEGIAN	
L07	3207318-260	21501-002197	V1	LABEL SET-ECG, POLISH	
L07	3207318-261	21501-002198	V1	LABEL SET-ECG, SPO2, POLISH	
L07	3207318-262	21501-002199	V1	LABEL SET-ECG, SPO2, CO2, POLISH	
L07	3207318-263	21501-002200	V1	LABEL SET-ECG, SPO2, NIBP, IP, POLISH	
L07	3207318-264	21501-002201	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, POLISH	
L07	3207318-265	21501-002202	V1	LABEL SET-ECG, NIBP, SPO2, POLISH	
L07	3207318-266	21501-002203	V1	LABEL SET-ECG, SPO2, NIBP, CO2, POLISH	
L07	3207318-280	21501-002283	V1	LABEL SET-ECG, HUNGARIAN	
L07	3207318-281	21501-002284	V1	LABEL SET-ECG, SPO2, HUNGARIAN	
L07	3207318-282	21501-002285	V1	LABEL SET-ECG, SPO2, CO2, HUNGARIAN	
L07	3207318-283	21501-002286	V1	LABEL SET-ECG, SPO2, NIBP, IP, HUNGARIAN	
L07	3207318-284	21501-002287	V1	LABEL SET-ECG, SPO2, NIBP, CO2, IP, HUNGARIAN	
L07	3207318-285	21501-002288	V1	LABEL SET-ECG, NIBP, SPO2, HUNGARIAN	
L07	3207318-286	21501-002289	V1	LABEL SET-ECG, SPO2, NIBP, CO2, HUNGARIAN	
L07	3207318-300	21501-002290	V1	LABEL SET-ECG, CZECH	
L07	3207318-301	21501-002291	V1	LABEL SET-ECG, SPO2, CZECH	
L07	3207318-302	21501-002292	V1	ECG, SpO2, CO2, Czech	

Title	LP15 V2 Service BOM			Page	24 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V1)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3207318-303	21501-002293	V1	ECG, SpO2, NIBP, IP, Czech	
L07	3207318-305	21501-002294	V1	ECG, SpO2, NIBP, CO2, IP, Czech	
L07	3207318-304	21501-002295	V1	ECG, SpO2, NIBP, Czech	
L07	3207318-306	21501-002296	V1	ECG, SpO2, NIBP, CO2, Czech	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-000	21501-002408	V2	LABEL SET-ECG,ENGLISH	
L07	3305642-001	21501-002409	V2	LABEL SET-ECG,SPO2,ENGLISH	
L07	3305642-002	21501-002410	V2	LABEL SET-ECG,NIBP,SPO2,ENGLISH	
L07	3305642-003	21501-002411	V2	LABEL SET-ECG,SPO2,CO2,ENGLISH	
L07	3305642-004	21501-002412	V2	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH	
L07	3305642-005	21501-002413	V2	LABEL SET-ECG,SPO2,NIBP,IP,ENGLISH	
L07	3305642-006	21501-002414	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ENGLISH	
L07	3305642-007	21501-002415	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ENGLISH	
L07	3305642-040	21501-002475	V2	LABEL SET-ECG,GERMAN	
L07	3305642-041	21501-002476	V2	LABEL SET-ECG,SPO2,GERMAN	
L07	3305642-042	21501-002477	V2	LABEL SET-ECG,NIBP,SPO2,GERMAN	
L07	3305642-043	21501-002478	V2	LABEL SET-ECG,SPO2,CO2,GERMAN	
L07	3305642-044	21501-002479	V2	LABEL SET-ECG,SPO2,NIBP,CO2,GERMAN	
L07	3305642-045	21501-002545	V2	LABEL SET-ECG,SPO2,NIBP,IP,GERMAN	
L07	3305642-046	21501-002546	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,GERMAN	
L07	3305642-047	21501-002547	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,GERMAN	
L07	3305642-060	21501-002548	V2	LABEL SET-ECG,ITALIAN	

Title	LP15 V2 Service BOM		Page	25 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-061	21501-002549	V2	LABEL SET-ECG,SPO2,ITALIAN	
L07	3305642-062	21501-002550	V2	LABEL SET-ECG,NIBP,SPO2,ITALIAN	
L07	3305642-063	21501-002551	V2	LABEL SET-ECG,SPO2,CO2,ITALIAN	
L07	3305642-064	21501-002552	V2	LABEL SET-ECG,SPO2,NIBP,CO2,ITALIAN	
L07	3305642-065	21501-002553	V2	LABEL SET-ECG,SPO2,NIBP,IP,ITALIAN	
L07	3305642-066	21501-002554	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ITALIAN	
L07	3305642-067	21501-002555	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ITALIAN	
L07	3305642-080	21501-002480	V2	LABEL SET-ECG,FRENCH	
L07	3305642-081	21501-002481	V2	LABEL SET-ECG,SPO2,FRENCH	
L07	3305642-082	21501-002482	V2	LABEL SET-ECG,NIBP,SPO2,FRENCH	
L07	3305642-083	21501-002483	V2	LABEL SET-ECG,SPO2,CO2,FRENCH	
L07	3305642-084	21501-002484	V2	LABEL SET-ECG,SPO2,NIBP,CO2,FRENCH	
L07	3305642-085	21501-002485	V2	LABEL SET-ECG,SPO2,NIBP,IP,FRENCH	
L07	3305642-086	21501-002486	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FRENCH	
L07	3305642-087	21501-002487	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FRENCH	
L07	3305642-100	21501-002488	V2	LABEL SET-ECG,DUTCH	
L07	3305642-101	21501-002489	V2	LABEL SET-ECG,SPO2,DUTCH	
L07	3305642-102	21501-002527	V2	LABEL SET-ECG,NIBP,SPO2,DUTCH	
L07	3305642-103	21501-002528	V2	LABEL SET-ECG,SPO2,CO2,DUTCH	
L07	3305642-104	21501-002529	V2	LABEL SET-ECG,SPO2,NIBP,CO2,DUTCH	
L07	3305642-105	21501-002530	V2	LABEL SET-ECG,SPO2,NIBP,IP,DUTCH	
L07	3305642-106	21501-002536	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DUTCH	
L07	3305642-107	21501-002500	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DUTCH	
L07	3305642-120	21501-002531	V2	LABEL SET-ECG,SPANISH	
L07	3305642-121	21501-002532	V2	LABEL SET-ECG,SPO2,SPANISH	
L07	3305642-122	21501-002533	V2	LABEL SET-ECG,NIBP,SPO2,SPANISH	
L07	3305642-123	21501-002534	V2	LABEL SET-ECG,SPO2,CO2,SPANISH	

Title	LP15 V2 Service BOM		Page	26 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-124	21501-002535	V2	LABEL SET-ECG,SPO2,NIBP,CO2,SPANISH	
L07	3305642-125	21501-002491	V2	LABEL SET-ECG,SPO2,NIBP,IP,SPANISH	
L07	3305642-126	21501-002493	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SPANISH	
L07	3305642-127	21501-002494	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SPANISH	
L07	3305642-140	21501-002501	V2	LABEL SET-ECG,PORTUGESE	
L07	3305642-141	21501-002502	V2	LABEL SET-ECG,SPO2,PORTUGESE	
L07	3305642-142	21501-002503	V2	LABEL SET-ECG,NIBP,SPO2,PORTUGESE	
L07	3305642-143	21501-002504	V2	LABEL SET-ECG,SPO2,CO2,PORTUGESE	
L07	3305642-144	21501-002505	V2	LABEL SET-ECG,SPO2,NIBP,CO2,PORTUGESE	
L07	3305642-145	21501-002506	V2	LABEL SET-ECG,SPO2,NIBP,IP,PORTUGESE	
L07	3305642-146	21501-002507	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,PORTUGESE	
L07	3305642-147	21501-002508	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,PORTUGESE	
L07	3305642-160	21501-002509	V2	LABEL SET-ECG,BRAZILIAN	
L07	3305642-161	21501-002557	V2	LABEL SET-ECG,SPO2,BRAZILIAN	
L07	3305642-162	21501-002558	V2	LABEL SET-ECG,NIBP,SPO2,BRAZILIAN	
L07	3305642-163	21501-002559	V2	LABEL SET-ECG,SPO2,CO2,BRAZILIAN	
L07	3305642-164	21501-002560	V2	LABEL SET-ECG,SPO2,NIBP,CO2,BRAZILIAN	
L07	3305642-165	21501-002561	V2	LABEL SET-ECG,SPO2,NIBP,IP,BRAZILIAN	
L07	3305642-166	21501-002562	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,BRAZILIAN	
L07	3305642-167	21501-002556	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,BRAZILIAN	
L07	3305642-180	21501-002510	V2	LABEL SET-ECG,SWEDISH	
L07	3305642-181	21501-002511	V2	LABEL SET-ECG,SPO2,SWEDISH	

Title	LP15 V2 Service BOM		Page	27 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-182	21501-002512	V2	LABEL SET-ECG,NIBP,SPO2,SWEDISH	
L07	3305642-183	21501-002513	V2	LABEL SET-ECG,SPO2,CO2,SWEDISH	
L07	3305642-184	21501-002490	V2	LABEL SET-ECG,SPO2,NIBP,CO2,SWEDISH	
L07	3305642-185	21501-002492	V2	LABEL SET-ECG,SPO2,NIBP,IP,SWEDISH	
L07	3305642-186	21501-002495	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SWEDISH	
L07	3305642-187	21501-002496	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SWEDISH	
L07	3305642-200	21501-002497	V2	LABEL SET-ECG,DANISH	
L07	3305642-201	21501-002498	V2	LABEL SET-ECG,SPO2,DANISH	
L07	3305642-202	21501-002499	V2	LABEL SET-ECG,NIBP,SPO2,DANISH	
L07	3305642-203	21501-002514	V2	LABEL SET-ECG,SPO2,CO2,DANISH	
L07	3305642-204	21501-002515	V2	LABEL SET-ECG,SPO2,NIBP,CO2,DANISH	
L07	3305642-205	21501-002516	V2	LABEL SET-ECG,SPO2,NIBP,IP,DANISH	
L07	3305642-206	21501-002517	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DANISH	
L07	3305642-207	21501-002518	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DANISH	
L07	3305642-220	21501-002519	V2	LABEL SET-ECG,FINNISH	
L07	3305642-221	21501-002520	V2	LABEL SET-ECG,SPO2,FINNISH	
L07	3305642-222	21501-002521	V2	LABEL SET-ECG,NIBP,SPO2,FINNISH	
L07	3305642-223	21501-002522	V2	LABEL SET-ECG,SPO2,CO2,FINNISH	
L07	3305642-224	21501-002523	V2	LABEL SET-ECG,SPO2,NIBP,CO2,FINNISH	
L07	3305642-225	21501-002524	V2	LABEL SET-ECG,SPO2,NIBP,IP,FINNISH	
L07	3305642-226	21501-002525	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FINNISH	
L07	3305642-227	21501-002526	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FINNISH	
L07	3305642-240	21501-002537	V2	LABEL SET-ECG,NORWEGIAN	
L07	3305642-241	21501-002538	V2	LABEL SET-ECG,SPO2,NORWEGIAN	
L07	3305642-242	21501-002540	V2	LABEL SET-ECG,NIBP,SPO2,NORWEGIAN	

Title	LP15 V2 Service BOM		Page	28 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-243	21501-002539	V2	LABEL SET-ECG,SPO2,CO2,NORWEGIAN	
L07	3305642-244	21501-002541	V2	LABEL SET-ECG,SPO2,NIBP,CO2,NORWEGIAN	
L07	3305642-245	21501-002542	V2	LABEL SET-ECG,SPO2,NIBP,IP,NORWEGIAN	
L07	3305642-246	21501-002543	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,NORWEGIAN	
L07	3305642-247	21501-002544	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,NORWEGIAN	
L07	3305642-260	21501-002563	V2	LABEL SET-ECG,POLISH	
L07	3305642-261	21501-002564	V2	LABEL SET-ECG,SPO2,POLISH	
L07	3305642-262	21501-002565	V2	LABEL SET-ECG,NIBP,SPO2,POLISH	
L07	3305642-263	21501-002566	V2	LABEL SET-ECG,SPO2,CO2,POLISH	
L07	3305642-264	21501-002567	V2	LABEL SET-ECG,SPO2,NIBP,CO2,POLISH	
L07	3305642-265	21501-002568	V2	LABEL SET-ECG,SPO2,NIBP,IP,POLISH	
L07	3305642-266	21501-002569	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,POLISH	
L07	3305642-267	21501-002570	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,POLISH	
L07	3305642-280	21501-002571	V2	LABEL SET-ECG,HUNGARIAN	
L07	3305642-281	21501-002572	V2	LABEL SET-ECG,SPO2,HUNGARIAN	
L07	3305642-282	21501-002573	V2	LABEL SET-ECG,NIBP,SPO2,HUNGARIAN	
L07	3305642-283	21501-002574	V2	LABEL SET-ECG,SPO2,CO2,HUNGARIAN	
L07	3305642-284	21501-002575	V2	LABEL SET-ECG,SPO2,NIBP,CO2,HUNGARIAN	
L07	3305642-285	21501-002576	V2	LABEL SET-ECG,SPO2,NIBP,IP,HUNGARIAN	
L07	3305642-286	21501-002577	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,HUNGARIAN	
L07	3305642-287	21501-002578	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,HUNGARIAN	
L07	3305642-300	21501-002579	V2	LABEL SET-ECG,CZECH	

Title	LP15 V2 Service BOM		Page	29 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-301	21501-002580	V2	LABEL SET-ECG,SPO2,CZECH	
L07	3305642-302	21501-002581	V2	LABEL SET-ECG,NIBP,SPO2,CZECH	
L07	3305642-303	21501-002582	V2	LABEL SET-ECG,SPO2,CO2,CZECH	
L07	3305642-304	21501-002583	V2	LABEL SET-ECG,SPO2,NIBP,CO2,CZECH	
L07	3305642-305	21501-002584	V2	LABEL SET-ECG,SPO2,NIBP,IP,CZECH	
L07	3305642-306	21501-002585	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,CZECH	
L07	3305642-307	21501-002586	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,CZECH	
L07	3305642-320	21501-002587	V2	LABEL SET-ECG,RUSSIAN	
L07	3305642-321	21501-002588	V2	LABEL SET-ECG,SPO2,RUSSIAN	
L07	3305642-322	21501-002589	V2	LABEL SET-ECG,NIBP,SPO2,RUSSIAN	
L07	3305642-323	21501-002590	V2	LABEL SET-ECG,SPO2,CO2,RUSSIAN	
L07	3305642-324	21501-002591	V2	LABEL SET-ECG,SPO2,NIBP,CO2,RUSSIAN	
L07	3305642-325	21501-002592	V2	LABEL SET-ECG,SPO2,NIBP,IP,RUSSIAN	
L07	3305642-326	21501-002593	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,RUSSIAN	
L07	3305642-327	21501-002594	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,RUSSIAN	
L07	3305642-380	21501-002595	V2	LABEL SET-ECG,KOREAN	
L07	3305642-381	21501-002596	V2	LABEL SET-ECG,SPO2,KOREAN	
L07	3305642-382	21501-002597	V2	LABEL SET-ECG,NIBP,SPO2,KOREAN	
L07	3305642-383	21501-002598	V2	LABEL SET-ECG,SPO2,CO2,KOREAN	
L07	3305642-384	21501-002599	V2	LABEL SET-ECG,SPO2,NIBP,CO2,KOREAN	
L07	3305642-385	21501-002600	V2	LABEL SET-ECG,SPO2,NIBP,IP,KOREAN	
L07	3305642-386	21501-002601	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,KOREAN	
L07	3305642-387	21501-002602	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,KOREAN	
L07	3305642-400	21501-002649	V2	LABEL SET-ECG, JAPANESE	
L07	3305642-401	21501-002650	V2	LABEL SET-ECG,SPO2,JAPANESE	
L07	3305642-402	21501-002651	V2	LABEL SET-ECG,NIBP, SPO2, JAPANESE	

Title	LP15 V2 Service BOM			Page	30 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L07	3305642-403	21501-002652	V2	LABEL SET-ECG,SPO2,CO2,JAPANESE	
L07	3305642-404	21501-002653	V2	LABEL SET-ECG,SPO2,NIBP,CO2,JAPANESE	
L07	3305642-405	21501-002654	V2	LABEL SET-ECG,SPO2,NIBP,IP,JAPANESE	
L07	3305642-406	21501-002655	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,JAPANESE	
L07	3305642-407	21501-002656	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,JAPANESE	

LIFEPAK 15 Label Set Languages (V2, 3 rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-008	21501-002700	V2	LABEL SET-ECG,ENGLISH, 3RD EDITION	
L08	3305642-009	21501-002701	V2	LABEL SET-ECG,SPO2,ENGLISH, 3RD EDITION	
L08	3305642-010	21501-002702	V2	LABEL SET-ECG,NIBP,SPO2,ENGLISH, 3RD EDITION	
L08	3305642-011	21501-002703	V2	LABEL SET-ECG,SPO2,CO2,ENGLISH, 3RD EDITION	
L08	3305642-012	21501-002704	V2	LABEL SET-ECG,SPO2,NIBP,CO2,ENGLISH, 3RD EDITION	
L08	3305642-013	21501-002705	V2	LABEL SET-ECG,SPO2,NIBP,IP,ENGLISH, 3RD EDITION	
L08	3305642-014	21501-002706	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ENGLISH, 3RD ED	

Title	LP15 V2 Service BOM		Page	31 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-015	21501-002707	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ENGLISH, 3RD ED	
L08	3305642-048	21501-002708	V2	LABEL SET-ECG,GERMAN, 3RD EDITION	
L08	3305642-049	21501-002709	V2	LABEL SET-ECG,SPO2,GERMAN, 3RD EDITION	
L08	3305642-050	21501-002710	V2	LABEL SET-ECG,NIBP,SPO2,GERMAN, 3RD EDITION	
L08	3305642-051	21501-002711	V2	LABEL SET-ECG,SPO2,CO2,GERMAN, 3RD EDITION	
L08	3305642-052	21501-002712	V2	LABEL SET-ECG,SPO2,NIBP,CO2,GERMAN, 3RD EDITION	
L08	3305642-053	21501-002713	V2	LABEL SET-ECG,SPO2,NIBP,IP,GERMAN, 3RD EDITION	
L08	3305642-054	21501-002714	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,GERMAN, 3RD ED	
L08	3305642-055	21501-002715	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,GERMAN, 3RD ED	
L08	3305642-068	21501-002716	V2	LABEL SET-ECG,ITALIAN, 3RD EDITION	
L08	3305642-069	21501-002717	V2	LABEL SET-ECG,SPO2,ITALIAN, 3RD EDITION	
L08	3305642-070	21501-002718	V2	LABEL SET-ECG,NIBP,SPO2,ITALIAN, 3RD EDITION	
L08	3305642-071	21501-002719	V2	LABEL SET-ECG,SPO2,CO2,ITALIAN, 3RD EDITION	
L08	3305642-072	21501-002720	V2	LABEL SET-ECG,SPO2,NIBP,CO2,ITALIAN, 3RD EDITION	

Title	LP15 V2 Service BOM		Page	32 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-073	21501-002721	V2	LABEL SET-ECG,SPO2,NIBP,IP,ITALIAN, 3RD EDITION	
L08	3305642-074	21501-002722	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,ITALIAN , 3RD ED	
L08	3305642-075	21501-002723	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,ITALIAN, 3RD ED	
L08	3305642-088	21501-002724	V2	LABEL SET-ECG,FRENCH, 3RD EDITION	
L08	3305642-089	21501-002725	V2	LABEL SET-ECG,SPO2,FRENCH, 3RD EDITION	
L08	3305642-090	21501-002726	V2	LABEL SET-ECG,NIBP,SPO2,FRENCH, 3RD EDITION	
L08	3305642-091	21501-002727	V2	LABEL SET-ECG,SPO2,CO2,FRENCH, 3RD EDITION	
L08	3305642-092	21501-002728	V2	LABEL SET-ECG,SPO2,NIBP,CO2,FRENCH, 3RD EDITION	
L08	3305642-093	21501-002729	V2	LABEL SET-ECG,SPO2,NIBP,IP,FRENCH, 3RD EDITION	
L08	3305642-094	21501-002730	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FRENCH , 3RD ED	
L08	3305642-095	21501-002731	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FRENCH, 3RD ED	
L08	3305642-108	21501-002732	V2	LABEL SET-ECG,DUTCH, 3RD EDITION	
L08	3305642-109	21501-002733	V2	LABEL SET-ECG,SPO2,DUTCH, 3RD EDITION	
L08	3305642-110	21501-002734	V2	LABEL SET-ECG,NIBP,SPO2,DUTCH, 3RD EDITION	
L08	3305642-111	21501-002735	V2	LABEL SET-ECG,SPO2,CO2,DUTCH, 3RD	

Title	LP15 V2 Service BOM		Page	33 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				EDITION	
L08	3305642-112	21501-002736	V2	LABEL SET- ECG,SPO2,NIBP,CO2,DUTCH, 3RD EDITION	
L08	3305642-113	21501-002737	V2	LABEL SET- ECG,SPO2,NIBP,IP,DUTCH, 3RD EDITION	
L08	3305642-114	21501-002738	V2	LABEL SET- ECG,SPO2,NIBP,CO2,IP,DUTCH, 3RD ED	
L08	3305642-115	21501-002739	V2	LABEL SET- ECG,SPO2,NIBP,CO2,TEMP,DUT CH, 3RD ED	
L08	3305642-128	21501-002740	V2	LABEL SET-ECG,SPANISH, 3RD EDITION	
L08	3305642-129	21501-002741	V2	LABEL SET-ECG,SPO2,SPANISH, 3RD EDITION	
L08	3305642-130	21501-002742	V2	LABEL SET- ECG,NIBP,SPO2,SPANISH, 3RD EDITION	
L08	3305642-131	21501-002743	V2	LABEL SET- ECG,SPO2,CO2,SPANISH, 3RD EDITION	
L08	3305642-132	21501-002744	V2	LABEL SET- ECG,SPO2,NIBP,CO2,SPANISH, 3RD EDITION	
L08	3305642-133	21501-002745	V2	LABEL SET- ECG,SPO2,NIBP,IP,SPANISH, 3RD EDITION	
L08	3305642-134	21501-002746	V2	LABEL SET- ECG,SPO2,NIBP,CO2,IP,SPANIS H, 3RD ED	
L08	3305642-135	21501-002747	V2	LABEL SET- ECG,SPO2,NIBP,CO2,TEMP,SPA NISH, 3RD ED	
L08	3305642-148	21501-002748	V2	LABEL SET-ECG,PORTUGUESE, 3RD EDITION	

Title	LP15 V2 Service BOM			Page	34 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-149	21501-002749	V2	LABEL SET- ECG,SPO2,PORTUGUESE, 3RD EDITION	
L08	3305642-150	21501-002750	V2	LABEL SET- ECG,NIBP,SPO2,PORTUGUESE, 3RD EDITION	
L08	3305642-151	21501-002751	V2	LABEL SET- ECG,SPO2,CO2,PORTUGUESE, 3RD EDITION	
L08	3305642-152	21501-002752	V2	LABEL SET- ECG,SPO2,NIBP,CO2,PORTUGU ESE, 3RD ED	
L08	3305642-153	21501-002753	V2	LABEL SET- ECG,SPO2,NIBP,IP,PORTUGUES E, 3RD ED	
L08	3305642-154	21501-002754	V2	LABEL SET- ECG,SPO2,NIBP,CO2,IP,PORTU GUESE, 3RDED	
L08	3305642-155	21501-002755	V2	LABEL SET- ECG,SPO2,NIBP,CO2,TEMP,POR TUGUESE, 3ED	
L08	3305642-168	21501-002756	V2	LABEL SET-ECG,BRAZILIAN, 3RD EDITION	
L08	3305642-169	21501-002757	V2	LABEL SET- ECG,SPO2,BRAZILIAN, 3RD EDITION	
L08	3305642-170	21501-002758	V2	LABEL SET- ECG,NIBP,SPO2,BRAZILIAN, 3RD EDITION	
L08	3305642-171	21501-002759	V2	LABEL SET- ECG,SPO2,CO2,BRAZILIAN, 3RD EDITION	
L08	3305642-172	21501-002760	V2	LABEL SET- ECG,SPO2,NIBP,CO2,BRAZILIAN , 3RD ED	
L08	3305642-173	21501-002761	V2	LABEL SET- ECG,SPO2,NIBP,IP,BRAZILIAN, 3RD ED	

Title	LP15 V2 Service BOM			Page	35 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-174	21501-002762	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,BRAZILIAN, 3RD ED	
L08	3305642-175	21501-002763	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,BRAZILIAN, 3ED	
L08	3305642-188	21501-002764	V2	LABEL SET-ECG,SWEDISH, 3RD EDITION	
L08	3305642-189	21501-002765	V2	LABEL SET-ECG,SPO2,SWEDISH, 3RD EDITION	
L08	3305642-190	21501-002766	V2	LABEL SET-ECG,NIBP,SPO2,SWEDISH, 3RD EDITION	
L08	3305642-191	21501-002767	V2	LABEL SET-ECG,SPO2,CO2,SWEDISH, 3RD EDITION	
L08	3305642-192	21501-002768	V2	LABEL SET-ECG,SPO2,NIBP,CO2,SWEDISH, 3RD EDITION	
L08	3305642-193	21501-002769	V2	LABEL SET-ECG,SPO2,NIBP,IP,SWEDISH, 3RD EDITION	
L08	3305642-194	21501-002770	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,SWEDISH, 3RD ED	
L08	3305642-195	21501-002771	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,SWEDISH, 3RD ED	
L08	3305642-208	21501-002772	V2	LABEL SET-ECG,DANISH, 3RD EDITION	
L08	3305642-209	21501-002773	V2	LABEL SET-ECG,SPO2,DANISH, 3RD EDITION	
L08	3305642-210	21501-002774	V2	LABEL SET-ECG,NIBP,SPO2,DANISH, 3RD EDITION	
L08	3305642-211	21501-002775	V2	LABEL SET-ECG,SPO2,CO2,DANISH, 3RD EDITION	

Title	LP15 V2 Service BOM		Page	36 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-212	21501-002776	V2	LABEL SET-ECG,SPO2,NIBP,CO2,DANISH, 3RD EDITION	
L08	3305642-213	21501-002777	V2	LABEL SET-ECG,SPO2,NIBP,IP,DANISH, 3RD EDITION	
L08	3305642-214	21501-002778	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,DANISH, 3RD ED	
L08	3305642-215	21501-002779	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,DANISH, 3RD ED	
L08	3305642-228	21501-002780	V2	LABEL SET-ECG,FINNISH, 3RD EDITION	
L08	3305642-229	21501-002781	V2	LABEL SET-ECG,SPO2,FINNISH, 3RD EDITION	
L08	3305642-230	21501-002782	V2	LABEL SET-ECG,NIBP,SPO2,FINNISH, 3RD EDITION	
L08	3305642-231	21501-002783	V2	LABEL SET-ECG,SPO2,CO2,FINNISH, 3RD EDITION	
L08	3305642-232	21501-002784	V2	LABEL SET-ECG,SPO2,NIBP,CO2,FINNISH, 3RD EDITION	
L08	3305642-233	21501-002785	V2	LABEL SET-ECG,SPO2,NIBP,IP,FINNISH, 3RD EDITION	
L08	3305642-234	21501-002786	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,FINNISH , 3RD ED	
L08	3305642-235	21501-002787	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,FINNISH, 3RD ED	
L08	3305642-248	21501-002788	V2	LABEL SET-ECG,NORWEGIAN, 3RD EDITION	
L08	3305642-249	21501-002789	V2	LABEL SET-ECG,SPO2,NORWEGIAN, 3RD EDITION	

Title	LP15 V2 Service BOM		Page	37 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-250	21501-002790	V2	LABEL SET-ECG,NIBP,SPO2,NORWEGIAN, 3RD EDITION	
L08	3305642-251	21501-002791	V2	LABEL SET-ECG,SPO2,CO2,NORWEGIAN, 3RD EDITION	
L08	3305642-252	21501-002792	V2	LABEL SET-ECG,SPO2,NIBP,CO2,NORWEGIAN, 3RD ED	
L08	3305642-253	21501-002793	V2	LABEL SET-ECG,SPO2,NIBP,IP,NORWEGIAN , 3RD ED	
L08	3305642-254	21501-002794	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,NORWEGIAN, 3RD ED	
L08	3305642-255	21501-002795	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,NORWEGIAN3ED	
L08	3305642-268	21501-002796	V2	LABEL SET-ECG,POLISH, 3RD EDITION	
L08	3305642-269	21501-002797	V2	LABEL SET-ECG,SPO2,POLISH, 3RD EDITION	
L08	3305642-270	21501-002798	V2	LABEL SET-ECG,NIBP,SPO2,POLISH, 3RD EDITION	
L08	3305642-271	21501-002799	V2	LABEL SET-ECG,SPO2,CO2,POLISH, 3RD EDITION	
L08	3305642-272	21501-002800	V2	LABEL SET-ECG,SPO2,NIBP,CO2,POLISH, 3RD EDITION	
L08	3305642-273	21501-002801	V2	LABEL SET-ECG,SPO2,NIBP,IP,POLISH, 3RD EDITION	
L08	3305642-274	21501-002802	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,POLISH, 3RD ED	
L08	3305642-275	21501-002803	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,POLI	

Title	LP15 V2 Service BOM			Page	38 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
				SH, 3RD ED	
L08	3305642-288	21501-002804	V2	LABEL SET-ECG,HUNGARIAN, 3RD EDITION	
L08	3305642-289	21501-002805	V2	LABEL SET-ECG,SPO2,HUNGARIAN, 3RD EDITION	
L08	3305642-290	21501-002806	V2	LABEL SET-ECG,NIBP,SPO2,HUNGARIAN, 3RD EDITION	
L08	3305642-291	21501-002807	V2	LABEL SET-ECG,SPO2,CO2,HUNGARIAN, 3RD EDITION	
L08	3305642-292	21501-002808	V2	LABEL SET-ECG,SPO2,NIBP,CO2,HUNGARIAN, 3RD ED	
L08	3305642-293	21501-002809	V2	LABEL SET-ECG,SPO2,NIBP,IP,HUNGARIAN, 3RD ED	
L08	3305642-294	21501-002810	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,HUNGARIAN, 3RD ED	
L08	3305642-295	21501-002811	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,HUNGARIAN,3RDED	
L08	3305642-308	21501-002812	V2	LABEL SET-ECG,CZECH, 3RD EDITION	
L08	3305642-309	21501-002813	V2	LABEL SET-ECG,SPO2,CZECH, 3RD EDITION	
L08	3305642-310	21501-002814	V2	LABEL SET-ECG,NIBP,SPO2,CZECH, 3RD EDITION	
L08	3305642-311	21501-002815	V2	LABEL SET-ECG,SPO2,CO2,CZECH, 3RD EDITION	
L08	3305642-312	21501-002816	V2	LABEL SET-ECG,SPO2,NIBP,CO2,CZECH, 3RD EDITION	

Title	LP15 V2 Service BOM			Page	39 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-313	21501-002817	V2	LABEL SET-ECG,SPO2,NIBP,IP,CZECH, 3RD EDITION	
L08	3305642-314	21501-002818	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,CZECH, 3RD ED	
L08	3305642-315	21501-002819	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,CZECH, 3RD ED	
L08	3305642-328	21501-002820	V2	LABEL SET-ECG,RUSSIAN, 3RD EDITION	
L08	3305642-329	21501-002821	V2	LABEL SET-ECG,SPO2,RUSSIAN, 3RD EDITION	
L08	3305642-330	21501-002822	V2	LABEL SET-ECG,NIBP,SPO2,RUSSIAN, 3RD EDITION	
L08	3305642-331	21501-002823	V2	LABEL SET-ECG,SPO2,CO2,RUSSIAN, 3RD EDITION	
L08	3305642-332	21501-002824	V2	LABEL SET-ECG,SPO2,NIBP,CO2,RUSSIAN, 3RD EDITION	
L08	3305642-333	21501-002825	V2	LABEL SET-ECG,SPO2,NIBP,IP,RUSSIAN, 3RD EDITION	
L08	3305642-334	21501-002826	V2	LABEL SET-ECG,SPO2,NIBP,CO2,IP,RUSSIAN, 3RD ED	
L08	3305642-335	21501-002827	V2	LABEL SET-ECG,SPO2,NIBP,CO2,TEMP,RUSSIAN, 3RD ED	
L08	3305642-388	21501-002828	V2	LABEL SET-ECG,KOREAN, 3RD EDITION	
L08	3305642-389	21501-002829	V2	LABEL SET-ECG,SPO2,KOREAN, 3RD EDITION	
L08	3305642-390	21501-002830	V2	LABEL SET-ECG,NIBP,SPO2,KOREAN, 3RD EDITION	

Title	LP15 V2 Service BOM			Page	40 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev	A
Doc State	Released	State Date	3/1/2024		

LIFEPAK 15 Label Set Languages (V2, 3rd Edition)					
Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
L08	3305642-391	21501-002831	V2	LABEL SET- ECG,SPO2,CO2,KOREAN, 3RD EDITION	
L08	3305642-392	21501-002832	V2	LABEL SET- ECG,SPO2,NIBP,CO2,KOREAN, 3RD EDITION	
L08	3305642-393	21501-002833	V2	LABEL SET- ECG,SPO2,NIBP,IP,KOREAN, 3RD EDITION	
L08	3305642-394	21501-002834	V2	LABEL SET- ECG,SPO2,NIBP,CO2,IP,KOREAN , 3RD ED	
L08	3305642-395	21501-002835	V2	LABEL SET- ECG,SPO2,NIBP,CO2,TEMP,KOR EAN, 3RD ED	
L08	3305642-408	21501-002836	V2	LABEL SET-ECG, JAPANESE, 3RD EDITION	
L08	3305642-409	21501-002837	V2	LABEL SET- ECG,SPO2,JAPANESE, 3RD EDITION	
L08	3305642-410	21501-002838	V2	LABEL SET-ECG,NIBP, SPO2, JAPANESE, 3RD EDITION	
L08	3305642-411	21501-002839	V2	LABEL SET- ECG,SPO2,CO2,JAPANESE, 3RD EDITION	
L08	3305642-412	21501-002840	V2	LABEL SET- ECG,SPO2,NIBP,CO2,JAPANESE, 3RD EDITION	
L08	3305642-413	21501-002841	V2	LABEL SET- ECG,SPO2,NIBP,IP,JAPANESE, 3RD EDITION	
L08	3305642-414	21501-002842	V2	LABEL SET- ECG,SPO2,NIBP,CO2,IP,JAPANE SE, 3RD ED	
L08	3305642-415	21501-002843	V2	LABEL SET- ECG,SPO2,NIBP,CO2,TEMP,JAP ANESE,3RD ED	

Title	LP15 V2 Service BOM		Page	41 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

5. SERVICE REPAIR KITS

LIFEPAK 15 V1 and V2 Service Repair Kits include components that support a particular replacement activity. Due to the need for special tooling and processes, the front and rear case repair kits come partially assembled from the factory. These repair kits are also listed in 3305431, LP15 Service Repair Kit. Reference the electronic PLM system Bill of Materials for part numbers and quantity for each assembly.

Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
K01	3305431-005	40577-000023	V1	KIT - REPAIR, INTERFACE PCBA, FOR V1 ONLY	
K02	3305431-008	40577-000016	V1, V2	KIT - REPAIR, THERAPY PCBA	
K03	3305431-009	40577-000019	V1, V2	KIT - REPAIR, THERAPY CONNECTOR, LP15	
K04	3305431-011	40577-000022	V1, V2	KIT - REPAIR, BIPHASIC MODULE, LP15	
K05	3305431-015	40577-000003	V1, V2	KIT - REPAIR, INTERNAL HARDWARE, LP15	
K06	3305431-016	40577-000004	V1, V2	KIT - REPAIR, EXTERNAL HARDWARE, LP15	
K07	3305431-018	40577-000006	V1, V2	KIT - REPAIR, DISPLAY SHIELD, LP15	
K08	3305431-019	40577-000007	V1, V2	KIT - REPAIR, GUARD AND FEET, LP15	
K09	3305431-020	40577-000008	V1, V2	KIT - REPAIR, HANDLE, LP15	
K10	3305431-023	40577-000010	V1, V2	KIT - REPAIR, FRONT CASE , LP15	
K11	3305431-026	40577-000025	V1, V2	KIT-REPAIR, MASIMO SPO2 PANEL MOUNT CABLE, LP15	
K12	3305431-027	40577-000026	V1, V2	KIT - REPAIR, TEMPERATURE CABLE, LP15	
K13	3305431-035	40577-000001	V1, V2	KIT - REPAIR, LED DISPLAY DRIVER, LP15	
K14	3305431-036	40577-000044	V1, V2	KIT - REPAIR, LED DISPLAY, LP15	If replacing older display (Board PN 3203221-002), also order new backlight inverter (40577-000001).
K15	3305431-501	40577-000034	V1, V2	KIT - REPAIR, IP CONNECTORS, LP15 ROHS	

Title	LP15 V2 Service BOM		Page	42 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
K16	3305431-502	40577-000035	V1, V2	KIT - REPAIR, ETCO2 MODULE, LP15 ROHS	
K17	3305431-510	40577-000037	V1, V2	KIT - REPAIR, ECG CONNECTOR, LP15 ROHS	
K18	3305431-514	40577-000038	V1, V2	KIT - REPAIR, CO2 CONNECTOR, LP15 ROHS	
K19	3305431-525	40577-000039	V1, V2	KIT - REPAIR, MASIMO SPO2, LP15 ROHS	
K20	3305431-527	40577-000040	V1, V2	KIT - REPAIR, TEMPERATURE CABLE, LP15 ROHS	
K21	3305431-528	40577-000041	V2	KIT - REPAIR, OEM PCBA, LP15 ROHS	
K22	3305431-532	40577-000045	V1, V2	KIT - REPAIR, NIBP ORION LP15 ROHS	
K23	3305431-536	40577-000051	V1, V2	KIT - REPAIR, REAR CASE, LP15	
K24	3305431-539	40577-000054	V1, V2	KIT - REPAIR, PADDLE RETAINER, LP15	
K25	3305431-540	40577-000055	V1, V2	KIT - REPAIR, NIBP CONNECTOR, LP15	
K26	3207824-001	11577-000002	V1, V2	KIT - CARRY BAG, MAIN BAG	From Front Case Parts List.
K27	3207827-003	11260-000039	V1, V2	KIT - CARRY BAG, REAR POUCH, 3RD EDITION	From Front Case Parts List.
K28	3207823-000	11577-000001	V1, V2	KIT - CARRY BAG, SHOULDER STRAP	From Front Case Parts List.
K29	3010267-00	11220-000028	V1, V2	POUCH, TOP, CARRYING BAG, L P12	From Front Case Parts List.
K31	3305431-006	40577-000015	V1	KIT - REPAIR, POWER PCBA, FOR V1 ONLY	Not orderable.
K32	3305431-004	40577-000014	V1	KIT - REPAIR, OEM PCBA, FOR V1 ONLY	Not orderable.
K33	3305431-007	40577-000020	V1	KIT - REPAIR, SYSTEM PCBA, FOR V1 ONLY	Not orderable.
K34	3305431-029	40577-000028	V1	KIT - REPAIR, INTERFACE PCBA, V2	Not orderable.
K35	3305431-529	40577-000042	V2	KIT - REPAIR, INTERFACE PCBA, LP15 ROHS	Not orderable.
K37	3305431-031	40577-000030	V2	KIT - REPAIR, SYSTEM PCBA	Not orderable.

Title	LP15 V2 Service BOM		Page	43 of 43
Doc Type	Service Product Drawing	Doc ID	3347648	Doc Rev A
Doc State	Released	State Date	3/1/2024	

Reference Designator	Part Number	Catalog Number	Version	Description	Note (Optional)
K38	3305431-030	40577-000029	V1, V2	KIT - REPAIR, POWER PCBA, V2	Not orderable.
K39	3305431-001	40577-000013	V1, V2	KIT - REPAIR, IP CONNECTORS, LP15	Not orderable.
K40	3305431-025	40577-000024	V1, V2	KIT - REPAIR, MASIMO SPO2, LP15	Not orderable.
K41	3305431-037	40577-000047	V1, V2	KIT - REPAIR, SYSTEM BOARD, V2	Not orderable.